

ORGANIC STANDARDS

Pesticide and Metabolite Standards

ISO methods ■

EN standards ■

ASTM methods ■

EPA methods ■

Contaminant standards ■

***Single Components Solutions available
in multiple solvents***

General Information

CPAchem is a French-Bulgarian company, located in Bulgaria.

Since 1998, we have been manufacturing a higher class of standard solutions and reagents in inorganic and organic.

CPAchem. is well recognized as a leader in production of **custom-made standard solutions according to customers specific requirements**. More than 4000 laboratories worldwide trust our brand based on the experience gained through the years.

CAPABILITIES

Complete range of analytical standards

- **CPA-spectr** : Inorganic calibration standards for Atomic Spectrometry (AAS, Graphite Furnace and ICP), Mass Spectrometry (ICP-MS) and Ion Chromatography.
- **CPA-pharm** : European Pharmacopoeia Reagents and Standards (Primary Opalescence Suspension & Colour solutions, Conductivity standards & Reference Buffer solutions, Reagents & Volumetric solutions, Ready to use Volumetric solutions)
- **CPA-org** : Organic analytical standards for residue analysis and environmental analysis (pesticides, pesticides metabolites, Deuterium and C13 labelled compounds, other...).

Custom manufacturing

- **Custom standard solutions in Organic**
- **Custom mono & multi-element standard in Inorganic**
- **Custom standards solution & reagents in Pharmacopoeia**

OUR TECHNOLOGY

The control of all the processes is fully automated by our dedicated Computer Aided Manufacturing (CAM) software, in addition to the modern network SQL –based data collecting system. The software system performs and controls the following :

- Automated calculations
- Barcode-driven movement
- Computer control of balances and other hardware.
- Incoming control of the raw materials
- Manufacturing and control of the intermediate solutions (bulks)
- Preparation of a custom-made solution (the program determines the needed weights and controls the gravimetric process on the analytical balances)
- Control of the final product (AAS, ICP, IC, GC, GC-MS, HPLC, LC-MS)
- Evaluating the final date, received by both gravimetric and instrumental analysis and calculating the certified values and uncertainties
- Automatic printing of labels, certificates and MSDS
- Automatic printing of Delivery Note for the products intended to be exported in order to obtain a full traceability every product and raw material is marked with a barcode.

OUR QUALITY

Each standard is rigorously tested and analyzed as a finished product and accompanied with a certificate of analysis and material safety data documentation.

Our products are manufactured under quality management system which has been introduced to the requirements of **ISO 9001:2000** and certified by Lloyd's Register Quality Assurance. We are currently in-process of implementing **ISO 17025** quality standards.

Custom Standards Information

GRAVIMETRIC DATA

A gravimetric certificate is provided with all standard mixtures, with the following informations :

- purity of the compound,
- source of raw material with lot number,
- concentration.

EXPIRY DATES

An expiry date is printed on both solution label and on the certificate.

INDEPENDENT LOTS

Differents lots are available . please specify " 2 differents lots".

SOLVENTS

Acetone	Dichloromethane	Iso-Octane	Petroleum Ether
Acetonitrile	Dimethylsulfoxide	Methanol	2- Propanol
Benzene	Ethanol	Methyl-tert-butylether	Tetrahydrofuran (THF)
Carbon disulfide	Ethyl acetate	n-Nonane	Toluene
Cyclohexane	n-Heptane	n-Octane	Water
Chloroform	n-Hexane	n-Pentane	Xylene

According to your specific needs other solvents are available.

Packaging information



Ampoule

1 ml – 2 ml
5 ml – 10 ml



Certan© bottle

1, 5 ml – 4,5 ml
10 ml



Amber bottle

20 ml – 40 ml
100 ml



Headspace Ampoule

10 ml – 20 ml

You can offer other bottles according your needs.

About the catalogue

ABBREVIATION TERMS

AC / ACs	Aromatic Compounds	PCB / PCBs	Polychlorinated Biphenyls
NCC / NCCs	Nitrogen-Containing Compounds	PAH / PAHs	Polycyclic Aromatic Hydrocarbons
OC / OCs	Organochlorine pesticides	VOC / VOCs	Volatile Organic Compounds
OPP / OPPs	Organophosphorus Pesticides	P & T	Purge & Trap

CONCENTRATION

The stated concentration in the mixtures refers to each of the components separately and it is not a total concentration.

Take advantage of our CUSTOM MULTI-ELEMENT SOLUTIONS

- *Quotations supplied quickly*
- *Mixtures made for your special needs*
- *Analysis Certificates according with the requirements of ISO 17025*
- *Available in ampoules, Certan[®], headspace bottle or other*

For our Custom Made certified Mixtures Request Form, see inside the page 149 or visit our web site at www.cpachem.com



Certificate of Analysis

REFERENCE MATERIAL CALIBRATION
Reference Materials Calibration Laboratory
CPA Ltd

Products for Analysis

(*)

Nº S24760

Date: 06.04.2010

Description of the RM
subject to calibration:

n-Alkane Standard Solution 6 components acc. ISO 9377-2-Mod
in n-Hexane

Catalog number:

9BAD.K1.ZT.1

Calibration method:

CRM's calibration procedure (MK-003)

Result from calibration
(certified value)
/Uncertainty:

Component	Chem. Formula	CAS No.	Certified Value and uncertainty (mg/l)*
n-Decane	C10H22	124-18-5	50.27 +/- 0.70
n-Hexadecane	C16H34	544-76-3	50.71 +/- 1.31
n-Eicosane	C20H42	112-95-8	49.97 +/- 1.68
n-Tetracosane	C24H50	646-31-1	50.06 +/- 1.18
n-Triacontane	C30H62	638-68-6	50.45 +/- 0.94
n-Tetracontane	C40H82	4181-95-7	50.09 +/- 2.01

* The expanded traceability of measurement is expressed as a mean quadratic uncertainty of measurement multiplied by the coverage factor of confidence interval $k = 2$ that corresponds to a confidence level of 95 % at a normal distribution. The mean quadratic uncertainty is determined in compliance with EA 4/02.

Metrological traceability:

n-Decane	ChemService, O-729
n-Hexadecane	ChemService, O-734
n-Eicosane	ChemService, F2183
n-Tetracosane	Dr Ehrenstorfer, C 17395800
n-Triacontane	ChemService, F2297
n-Tetracontane	Dr Ehrenstorfer, C 17395500

The metrological traceability is assured through calibration on GC/MS. The calibration curve is drawn using a series of standard solutions prepared from a certified reference material of Dr Ehrenstorfer, traceable to SI. All contributions in relation to the preparation of standard solutions are considered when evaluating the uncertainty

Contents

ISO Methods

Petroleum products

ISO 3924:2010	Determination of boiling range distribution - Gas chromatography method	5
---------------	---	---

Water quality

ISO 6468:1996	Determination of certain organochlorine insecticides, polychlorinated biphenyls and chlorobenzenes - Gas chromatographic method after liquid-liquid extraction	5
ISO 7981-1:2005	Determination of polycyclic aromatic hydrocarbons (PAH) - Part 1: Determination of six PAH by high-performance thin-layer chromatography with fluorescence detection after liquid-liquid extraction	6
ISO 7981-2:2005	Determination of polycyclic aromatic hydrocarbons (PAH) - Part 2: Determination of six PAH by high-performance liquid chromatography with fluorescence detection after liquid-liquid extraction	6
ISO 8165-1:1992	Determination of selected monovalent phenols - Part 1: Gas-chromatographic method after enrichment by extraction	6
ISO 8165-2:1999	Determination of selected monovalent phenols - Part 2: Method by derivatization and gas chromatography	7
ISO 9377-2:2000	Determination of hydrocarbon oil index - Part 2: Method using solvent extraction and gas chromatography	7
ISO 9377-2-Mod:	Determination of hydrocarbon oil index down to C7	7-9
ISO 9562:2004	Determination of absorbable organically bound halogens (AOX)	9
ISO 10301:1997	Determination of highly volatile halogenated hydrocarbons - Gas-chromatographic methods	9
ISO 10695:2000	Determination of selected organic nitrogen and phosphorus compounds - Gas chromatographic methods	10
ISO 11369:1997	Determination of selected plant treatment agents - Method using high performance liquid chromatography with UV detection after solid-liquid extraction	10
ISO 11423-2:1997	Determination of benzene and some derivatives - Part 2: Method using extraction and gas chromatography	10
ISO 15680:2003	Gas-chromatographic determination of a number of monocyclic aromatic hydrocarbons, naphthalene and several chlorinated compounds using purge-and-trap and thermal desorption	11,12
ISO 15913:2000	Determination of selected phenoxyalkanoic herbicides, including bentazones and hydroxybenzonitriles by gas chromatography and mass spectrometry after solid phase extraction and derivatization	12,13
ISO 17495:2001	Determination of selected nitrophenols - Method by solid-phase extraction and gas chromatography with mass spectrometric detection	13
ISO 17858:2007	Determination of dioxin-like polychlorinated biphenyls - Method using gas chromatography/mass spectrometry	13
ISO 17993:2002	Determination of 15 polycyclic aromatic hydrocarbons (PAH) in water by HPLC with fluorescence detection after liquid-liquid extraction	13
ISO 18856:2004	Determination of selected phthalates using gas chromatography/mass spectrometry	14
ISO 18857-1:2005	Determination of selected alkylphenols - Part 1: Method for non-filtered samples using liquid-liquid extraction and gas chromatography with mass selective detection	14
ISO 23631:2006	Determination of dalapon, trichloroacetic acid and selected haloacetic acids - Method using gas chromatography (GC-ECD and/or GC-MS detection) after liquid-liquid extraction and derivatization	14

Soil quality

ISO 10382:2002	Determination of organochlorine pesticides and polychlorinated biphenyls - Gas-chromatographic method with electron capture detection	9
ISO 13877:1998	Determination of polynuclear aromatic hydrocarbons - Method using high performance liquid chromatography	10
ISO 14154:2005	Determination of some selected chlorophenols - Gas-chromatographic method with electron-capture detection	10,11
ISO 15009:2002	Gas chromatographic determination of the content of volatile aromatic hydrocarbons, naphthalene and volatile halogenated hydrocarbons - Purge-and-trap method with thermal desorption	11
ISO 18287:2006	Determination of polycyclic aromatic hydrocarbons (PAH) - Gas chromatographic method with mass spectrometric detection (GC-MS)	14

Animal and vegetable fats and oils

ISO 15753:2006	Determination of polycyclic aromatic hydrocarbons	12
----------------	---	----

EN Standards

EN 12673	Water quality - Gas chromatographic determination of some selected chlorophenols in water	15
EN 12766	Petroleum products and used oils - Determination of PCBs and related products	15
EN 12916	Petroleum products - Determination of aromatic hydrocarbon types in middle distillates - High performance liquid chromatography method with refractive index detection	16
EN 12918	Water quality - Determination of parathion, parathion-methyl and some other organophosphorus compounds in water by dichloromethane extraction and gas chromatographic analysis	16
EN 13628	Packaging - Flexible packaging material - Determination of residual solvents by static headspace gas chromatography - Part 2: Industrial methods	16
EN 14039	Characterization of waste - Determination of hydrocarbon content in the range of C10 to C40 by gas chromatography	16,17
EN 15308:2008	Characterization of waste - Determination of selected polychlorinated biphenyls (PCB) in solid waste by using capillary gas chromatography with electron capture or mass spectrometric detection	17
EN 61619	Insulating liquids - Contamination by polychlorinated biphenyls (PCBs) - Method of determination by capillary column gas chromatography (IEC 61619:1997)	17

ASTM Methods

ASTM D664 -07	Standard Test Method for Acid Number of Petroleum Products by Potentiometric Titration	18
ASTM D2887-06a	Standard Test Method for Boiling Range Distribution of Petroleum Fractions by Gas Chromatography	18
ASTM D3606-07	Standard Test Method for Determination of Benzene and Toluene in Finished Motor and Aviation Gasoline by Gas Chromatography	19
ASTM D3710-95	Standard Test Method for Boiling Range Distribution of Gasoline and Gasoline Fractions by Gas Chromatography	19
ASTM D3798 -03	Standard Test Method for Analysis of p-Xylene by Gas Chromatography	19
ASTM D4815 -04	Standard Test Method for Determination of ethers and alcohols in gasoline by Gas Chromatography	19,20
ASTM D5186 -03	Standard Test Method which measures the total aromatic and PNA content in Diesel and Jet Fuel	21
ASTM D5307 -97	Standard Test Method for Determination of Boiling Range Distribution of Crude Petroleum by Gas Chromatography	21
ASTM D5442 -93	Standard Test Method for Analysis of Petroleum Waxes by Gas Chromatography.....	21
ASTM D5580 -02	Standard Test Method for Determination of Benzene, Toluene, Ethylbenzene, p/m-Xylene, o-Xylene, C9 and Heavier Aromatics, and Total Aromatics in Finished Gasoline by Gas Chromatography	22,23
ASTM E1387 -01	Standard Test Method for Ignitable Liquid Residues in Extracts from Fire Debris Samples by Gas Chromatography	23

EPA Methods

EPA Methods 500	Series for Drinking Water	24-39
EPA Methods 600	Series for monitoring organic pollutants in industrial and municipal waste discharges	40-49
EPA Methods 1300	SeriesMethods	49-50
EPA Methods 8000	Series for monitoring organic pollutants in ground water	50-67

Contaminant Standards

Aromatic Hydrocarbons	68,69
Carbonyl and Carbonyl -DNPH standards	69-71
Ethers and Alcohols	71
Haloorganic Acids and Esters	71,72
Hydrocarbon Standards	72-76
Nitrogen containing compounds (Amines, Aniline derivates, nitroaromatics)	77
Pesticides	77-90
Phenols	90-93
Phthalates	93,94
Polychlorinated Biphenyls (PCBs) Congeners Mixtures	94-97
Polycyclic aromatic hydrocarbons (PAH)	97-100
Residual Solvents in Flexible packaging material	100
Residual Solvents in Medicinal Products	100,101
Revised USP 467 Residual Solvent Method	101
Sulfur in Petroleum Products	102,103
Volatile and Semi- Volatile Organic Compounds	103-110

Single component Solutions

Compounds	111-148
-----------------	---------

Custom-made Solution Standards

Custom Standard Request Form	149
Certificate of Analysis	150,151

Index

Cas Number	152,153
Reference Number	154-161

ISO Methods

ISO 3924:2010

Petroleum products - Determination of boiling range distribution - Gas chromatography method

ISO 3924:2010 specifies a method for the determination of the boiling-range distribution of petroleum products.
ISO 3924:2010 has successfully been applied to samples containing biodiesel up to 10%.

Calibration Mix Solution 1 - 17 components

n-Hexane (C6)	CAS:110-54-3	600 ug/ml	n-Octadecane (C18)	CAS:593-45-3	500 ug/ml
n-Heptane (C7)	CAS:142-82-5	600 ug/ml	n-Eicosane (C20)	CAS:112-95-8	200 ug/ml
n-Octane (C8)	CAS:111-65-9	800 ug/ml	n-Tetracosane (C24)	CAS:646-31-1	200 ug/ml
n-Nonane (C9)	CAS:111-84-2	800 ug/ml	n-Octacosane (C28)	CAS:630-02-4	50 ug/ml
n-Decane (C10)	CAS:124-18-5	1 200 ug/ml	n-Dotriacontane (C32)	CAS:544-85-4	100 ug/ml
n-Undecane (C11)	CAS:1120-21-4	1 200 ug/ml	n-Hexatriacontane (C36)	CAS:630-06-8	100 ug/ml
n-Dodecane (C12)	CAS:112-40-3	1 200 ug/ml	n-Tetracontane (C40)	CAS:4181-95-7	100 ug/ml
n-Tetradecane (C14)	CAS:629-59-4	1 200 ug/ml	n-Tetratetracontane (C44)	CAS:7098-22-8	100 ug/ml
n-Hexadecane (C16)	CAS:544-76-3	1 000 ug/ml			
Solvent : Chloroform		ampoule 1 ml		Ref : F110061	Price : GF

Calibration Mix Solution 2 – 20 components

n-Pentane (C5)	CAS:109-66-0	n-Dodecane (C12)	CAS:112-40-3	n-Tetracosane (C24)	CAS:646-31-1
n-Hexane (C6)	CAS:110-54-3	n-Tetradecane (C14)	CAS:629-59-4	n-Octacosane (C28)	CAS:630-02-4
n-Heptane (C7)	CAS:142-82-5	n-Pentadecane (C15)	CAS:629-62-9	n-Dotriacontane (C32)	CAS:544-85-4
n-Octane (C8)	CAS:111-65-9	n-Hexadecane (C16)	CAS:544-76-3	n-Hexatriacontane (C36)	CAS:630-06-8
n-Nonane (C9)	CAS:111-84-2	n-Heptadecane (C17)	CAS:629-78-7	n-Tetracontane (C40)	CAS:4181-95-7
n-Decane (C10)	CAS:124-18-5	n-Octadecane (C18)	CAS:593-45-3	n-Tetratetracontane (C44)	CAS:7098-22-8
n-Undecane (C11)	CAS:1120-21-4	n-Eicosane (C20)	CAS:112-95-8		
Solvent : n-Octane		ampoule 1 ml	0.1% (w/w)	Ref : F110071	Price : HF

Column Test Mix -2 components

n-Hexadecane (C16)	CAS:544-76-3	n-Octadecane (C18)	CAS:593-45-3		
Solvent : n-Octane		ampoule 1 ml	1% (w/v)	Ref : F110051	Price : FA

Fuel Oil Degradation/Retention Time Mix C17-C20 – 4 components

n-Heptadecane (C17)	CAS:629-78-7	Phytane	CAS:638-36-8	Pristane	CAS:1921-70-6
n-Octadecane (C18)	CAS:593-45-3				
Solvent : Dichloromethane/Carbon disulfide (1/1)		ampoule 1 ml	2 000 µg/ml	Ref : F110081	Price : DG

Hydrocarbon Window Defining Standard C8-C40 – 35 components

n-Octane (C8)	CAS:111-65-9	Pristane (C19)	CAS:1921-70-6	n-Triacontane (C30)	CAS:638-68-6
n-Nonane (C9)	CAS:111-84-2	n-Eicosane (C20)	CAS:112-95-8	n-Hentriacontane (C31)	CAS:630-04-6
n-Decane (C10)	CAS:124-18-5	Phytane (C20)	CAS:638-36-8	n-Dotriacontane (C32)	CAS:544-85-4
n-Undecane (C11)	CAS:1120-21-4	n-Heneicosane (C21)	CAS:629-94-7	n-Tritriacontane (C33)	CAS:630-05-7
n-Dodecane (C12)	CAS:112-40-3	n-Docosane (C22)	CAS:629-97-0	n-Tettracontane (C34)	CAS:14167-59-0
n-Tridecane (C13)	CAS:629-50-5	n-Tricosane (C23)	CAS:638-67-5	n-Pentatriacontane (C35)	CAS:630-07-9
n-Tetradecane (C14)	CAS:629-59-4	n-Tetracosane (C24)	CAS:646-31-1	n-Hexatriacontane (C36)	CAS:630-06-8
n-Pentadecane (C15)	CAS:629-62-9	n-Pentacosane (C25)	CAS:629-99-2	n-Heptatriacontane (C37)	CAS:7194-84-5
n-Hexadecane (C16)	CAS:544-76-3	n-Hexacosane (C26)	CAS:630-01-3	n-Octatriacontane (C38)	CAS:7194-85-6
n-Heptadecane (C17)	CAS:629-78-7	n-Heptacosane (C27)	CAS:593-49-7	n-Nonatriacontane (C39)	CAS:7194-86-7
n-Octadecane (C18)	CAS:593-45-3	n-Octacosane (C28)	CAS:630-02-4	n-Tetracontane (C40)	CAS:4181-95-7
n-Nonadecane (C19)	CAS:629-92-5	n-Nonacosane (C29)	CAS:630-03-5		
Solvent : Carbon disulfide/Dichloromethane (3/1)		ampoule 1 ml	500 µg/ml	Ref : F061935	Price : KF
		ampoule 5 ml	500 µg/ml	Ref : F110241	Price : DJA

ISO 6488 :1996

Water quality - Determination of certain organochlorine insecticides, polychlorinated biphenyls and chlorobenzenes - Gas chromatographic method after liquid-liquid extraction.

Describes a method for determining certain organochlorine insecticides, polychlorinated biphenyls (PCBs) and chlorobenzenes (except the mono- and dichlorobenzenes) in drinking water, ground water, surface waters and waste waters. The method is applicable to samples containing up to 0,05 g/l of suspended solids.

OCs Pesticide Standard Solution - 20 components

Aldrin	CAS:309-00-2	Dieldrin	CAS:60-57-1	Delta-HCH	CAS:319-86-8
2,4'-DDE	CAS:3424-82-6	Endosulfan-alpha	CAS:959-98-8	epsilon-HCH	CAS:6108-10-7
4,4'-DDE	CAS:72-55-9	Endosulfan-beta	CAS:33213-65-9	Heptachlor	CAS:76-44-8
2,4'-DDD (o,p'-TDE)	CAS:53-19-0	Endrin	CAS:72-20-8	Heptachlor-endo-epoxide	CAS:28044-83-9
4,4'-DDD (TDE)	CAS:72-54-8	Alpha-HCH	CAS:319-84-6	Heptachlor-exo-epoxide	CAS:1024-57-3
2,4'-DDT	CAS:789-02-6	Beta-HCH	CAS:319-85-7	Methoxychlor (DMTD)	CAS:72-43-5
4,4'-DDT	CAS:50-29-3	Gamma-HCH (Lindane)	CAS:58-89-9		
Solvent : Iso-octane		ampoule 1 ml	10 µg/ml	Ref : F147410	Price : EG

OCs/PCBs/ACs Standard Solution - 36 components

Aldrin	CAS:309-00-2	delta -HCH	CAS:319-86-8	PCB 153	CAS:35065-27-1
2,4'-DDE	CAS:3424-82-6	gamma -HCH (Lindane)	CAS:58-89-9	PCB 180	CAS:35065-29-3
4,4'-DDE	CAS:72-55-9	epsilon-HCH	CAS:6108-10-7	PCB 194	CAS:35694-08-7
4,4'-DDD (TDE)	CAS:72-54-8	Heptachlor	CAS:6108-10-7	Pentachlorobenzene	CAS:608-93-5
2,4'-DDT	CAS:789-02-6	Heptachlor-endo-epoxide	CAS:76-44-8	Pentachloronitrobenzene	CAS:82-68-8
4,4'-DDT	CAS:50-29-3	Heptachlor-exo-epoxide	CAS:28044-83-9	o,p'-TDE	CAS:53-19-0
Dieldrin	CAS:60-57-1	Hexachlorobenzene	CAS:118-74-1	1,2,3,4-Tetrachlorobenzene	CAS:634-66-2
Endosulfan-alpha	CAS:959-98-8	Methoxychlor (DMTD)	CAS:72-43-5	1,2,3,5-Tetrachlorobenzene	CAS:634-90-2
Endosulfan-beta	CAS:33213-65-9	PCB 28	CAS:7012-37-5	1,2,4,5-Tetrachlorobenzene	CAS:95-94-3
Endrin	CAS:72-20-8	PCB 52	CAS:35693-99-3	1,2,3-Trichlorobenzene	CAS:87-61-6
alpha -HCH	CAS:319-84-6	PCB 101	CAS:37680-73-2	1,2,4-Trichlorobenzene	CAS:120-82-1
beta -HCH	CAS:319-85-7	PCB 138	CAS:35065-28-2	1,3,5-Trichlorobenzene	CAS:108-70-3
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F109230	Price : BEF

ISO 7981-1:2005**Water quality - Determination of polycyclic aromatic hydrocarbons (PAH) - Part 1: Determination of six PAH by high-performance thin-layer chromatography with fluorescence detection after liquid-liquid extraction**

Specifies the determination of six selected polycyclic aromatic hydrocarbons (PAH) in drinking water. The 6 PAH are: fluoranthene, benzo[b]fluoranthene, benzo[a]pyrene, benzo[k]fluoranthene, indeno[1,2,3-cd]pyrene, and benzo[ghi]perylene.

PAH Standard Solution - 6 components

Fluoranthene	CAS:206-44-0	10 µg/ml	Benzo(g,h,i)perylene	CAS:191-24-2	2 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	2 µg/ml	Benzo(a)pyrene	CAS:50-32-8	2 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	2 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	2 µg/ml
Solvent : Acetonitrile		ampoule 5 ml		Ref : F109250	Price : EF

ISO 7981-2:2005**Water quality - Determination of polycyclic aromatic hydrocarbons (PAH) - Part 2: Determination of six PAH by high-performance liquid chromatography with fluorescence detection after liquid-liquid extraction.**

ISO 7981-2:2005 specifies the determination of six selected PAH in drinking, mineral and table waters and ground and surface waters in mass concentrations above 0,005 micrograms per litre by high-performance liquid chromatography with fluorescence detection after liquid-liquid extraction. The six PAH are: fluoranthene, benzo[b]fluoranthene, benzo[a]pyrene, benzo[k]fluoranthene, indeno[1,2,3-cd]pyrene, and benzo[ghi]perylene. With some modification, this method is also applicable for the analysis of moderately polluted waste waters.

PAH Standard Solution - 6 components

Fluoranthene	CAS:206-44-0	Benzo(k)fluoranthene	CAS:207-08-9	Benzo(a)pyrene	CAS:50-32-8
Benzo(b)fluoranthene	CAS:205-99-2	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Solvent : Acetonitrile		ampoule 5 ml	10 µg/ml	Ref : F109260	Price : EF

ISO 8165-1:1992**Water quality - Determination of selected monovalent phenols -- Part 1: Gas-chromatographic method after enrichment by extraction.**

This international method specifies a method for determining phenols by a gas chromatographic method. The determination of selected phenols and nitrophenols in drinking water, ground water and surface water. The method provides a procedure for solid-phase extraction of nitrophenols, followed by solvent elution, derivatization with diazomethane and determination by gas chromatography and mass spectrometry.

Phenol Internal Standard Stock Solution - 1 component

2,4-Dibromophenol	CAS:615-58-7				
Solvent : Acetone		ampoule 5 ml	1000 µg/ml	Ref : F109280	Price : CF

Phenol Standard Stock Solution - 33 components

2-Benzylphenol	CAS:28994-41-4	2,4-Dichloro-3,5-dimethylphenol	CAS:133-53-9	2-Naphthol	CAS:135-19-3
2-Benzyl-4-chlorophenol	CAS:120-32-1	2,3-Dichlorophenol	CAS:576-24-9	Pentachlorophenol	CAS:87-86-5
4-tert-Butyl-2-chlorophenol	CAS:98-28-2	2,4-Dichlorophenol	CAS:120-83-2	Phenol	CAS:108-95-2
2-Chlorophenol	CAS:95-57-8	2,5-Dichlorophenol	CAS:583-78-8	2-Phenylphenol	CAS:90-43-7
3-Chlorophenol	CAS:108-43-0	2,6-Dichlorophenol	CAS:87-65-0	2,3,4,5-Tetrachlorophenol	CAS:4901-51-3
4-Chlorophenol	CAS:106-48-9	2,4-Dimethylphenol	CAS:105-67-9	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
4-Chloro-2-isopropyl-5-methylphenol	CAS:89-68-9	4-Ethylphenol	CAS:123-07-9	2,3,5,6-Tetrachlorophenol	CAS:935-95-5
4-Chloro-2-methylphenol	CAS:1570-64-5	2-Methylphenol	CAS:95-48-7	2,3,5-Trichlorophenol	CAS:933-78-8
4-Chloro-3-methylphenol	CAS:59-50-7	3-Methylphenol	CAS:108-39-4	2,3,6-Trichlorophenol	CAS:933-75-5
2-Chloro-5-methylphenol	CAS:615-74-7	4-Methylphenol	CAS:106-44-5	2,4,5-Trichlorophenol	CAS:95-95-4
2,6-Di-tert-butyl-4-methylphenol	CAS:128-37-0	1-Naphthol	CAS:90-15-3	2,4,6-Trichlorophenol	CAS:88-06-2
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F109270	Price : KA

ISO 8165-2 :1999**Water quality - Determination of selected monovalent phenols - Part 2: Method by derivatization and gas chromatography**

This international method specifies a method for the determination of phenols by gas chromatography, following pentafluorobenzoyl chloride derivatization. It may in particular be applied to the examination of drinking water, ground water and moderately contaminated surface water. With this method, lower limits of detection may be obtained compared with extraction procedures.

Internal Control Stock Solution - 1 component

2,4-Dibromophenol CAS:615-58-7

Solvent : Methanol purge & trap ampoule 5 ml 10 µg/ml Ref : F109310 Price : BJ

Phenol Standard Solution - 26 components

Phenol	CAS:108-95-2	2,4-Dichlorophenol	CAS:120-83-2	2-Phenylphenol	CAS:90-43-7
2-Benzylphenol	CAS:28994-41-4	2,5-Dichlorophenol	CAS:583-78-8	2,3,4,5-Tetrachlorophenol	CAS:4901-51-3
2-Chlorophenol	CAS:95-57-8	2,6-Dichlorophenol	CAS:87-65-0	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
3-Chlorophenol	CAS:108-43-0	2,4-Dimethylphenol	CAS:105-67-9	2,3,5,6-Tetrachlorophenol	CAS:935-95-5
4-Chlorophenol	CAS:106-48-9	4-Ethylphenol	CAS:123-07-9	2,3,5-Trichlorophenol	CAS:933-78-8
4-Chloro-2-methylphenol	CAS:1570-64-5	2-Methylphenol	CAS:95-48-7	2,3,6-Trichlorophenol	CAS:933-75-5
4-Chloro-3-methylphenol	CAS:59-50-7	3-Methylphenol	CAS:108-39-4	2,4,5-Trichlorophenol	CAS:95-95-4
2,4-Dichloro-3,5-dimethylphenol	CAS:133-53-9	4-Methylphenol	CAS:106-44-5	2,4,6-Trichlorophenol	CAS:88-06-2
2,3-Dichlorophenol	CAS:576-24-9	Pentachlorophenol	CAS:87-86-5		

Solvent : Methanol purge & trap ampoule 1 ml 30 µg/ml Ref : F129181 Price : HF

Phenol Standard Stock Solution - 31 components

Phenol	CAS:108-95-2	2,6-Di-tert-butyl-4-methylphenol	CAS:128-37-0	Pentachlorophenol	CAS:87-86-5
2-Benzyl-4-chlorophenol	CAS:120-32-1	2,4-Dichloro-3,5-dimethylphenol	CAS:133-53-9	2-Phenylphenol	CAS:90-43-7
2-Benzylphenol	CAS:28994-41-4	2,4-Dimethylphenol	CAS:105-67-9	2,3,4,5-Tetrachlorophenol	CAS:4901-51-3
4-tert-Butyl-2-chlorophenol	CAS:98-28-2	2,3-Dichlorophenol	CAS:576-24-9	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
4-Chloro-2-isopropyl-5-methylphenol	CAS:89-68-9	2,4-Dichlorophenol	CAS:120-83-2	2,3,5,6-Tetrachlorophenol	CAS:935-95-5
4-Chloro-2-methylphenol	CAS:1570-64-5	2,5-Dichlorophenol	CAS:583-78-8	2,3,5-Trichlorophenol	CAS:933-78-8
4-Chloro-3-methylphenol	CAS:59-50-7	2,6-Dichlorophenol	CAS:87-65-0	2,3,6-Trichlorophenol	CAS:933-75-5
2-Chloro-5-methylphenol	CAS:615-74-7	4-Ethylphenol	CAS:123-07-9	2,4,5-Trichlorophenol	CAS:95-95-4
2-Chlorophenol	CAS:95-57-8	2-Methylphenol	CAS:95-48-7	2,4,6-Trichlorophenol	CAS:88-06-2
3-Chlorophenol	CAS:108-43-0	3-Methylphenol	CAS:108-39-4		
4-Chlorophenol	CAS:106-48-9	4-Methylphenol	CAS:106-44-5		

Solvent : Methanol purge & trap ampoule 1 ml 300 µg/ml Ref : F109300 Price : BBA

ISO 9377-2 :2000**Water quality - Determination of hydrocarbon oil index - Part 2: Method using solvent extraction and gas chromatography**

The method of choice for the determination of oil and grease in water is based on extraction with hydrocarbon solvent like pentane or hexane.

This test determines the hydrocarbon oil index in water by means of gas chromatography. The method is suitable for surface water, waste water, and water from sewage treatment plants, and allows the determination of the hydrocarbon oil index in concentrations above 0.1 mg/l.

The index is the sum of compounds with retention times between n-decane and n-tetracontane. Substances complying with this definition are long-chain or branched aliphatic, alicyclic, aromatic or alkylsubstituted aromatic hydrocarbons.

ISO 9377-2 :Mod**Water quality -Determination of hydrocarbon oil index down to C7**

A modified method for oil in water analysis of produced water from offshore petroleum production installations has recently been taken into force.

ISO 9377-2 is not applicable for volatile hydrocarbons, and a modification of the method is proposed by the OSPAR commission in order to include the determination of certain hydrocarbons with boiling points between 98 and 174 °C from produced water.

The OSPAR Reference Method of Analysis for the Determination of the Dispersed Oil Content in Produced Water.

Mixture of Mineral Oils (A+B) (Diesel Fuel / Lubricating Oil) : 1/1(neat) - 2 components

Mineral Oil 50% CAS:8042-47-5 Diesel Oil 50% N/A (id : 381)

ampoule 1 ml Ref : F109631 Price : DF

ampoule 5 ml Ref : F061945 Price : KJ

Calibration Mixture of Mineral Oils (A+B) Set - 2 components

Ampoule 1	Blank	n-Hexane			Ampoule 4	Mineral Oil	CAS : 8042-47-5 + Diesel Oil	0.6 mg/ml each
Ampoule 2	Mineral Oil	CAS : 8042-47-5	+ Diesel Oil	0.2 mg/ml each	Ampoule 5	Mineral Oil	CAS : 8042-47-5 + Diesel Oil	0.8 mg/ml each
Ampoule 3	Mineral Oil	CAS : 8042-47-5	+ Diesel Oil	0.4 mg/ml each	Ampoule 6	Mineral Oil	CAS : 8042-47-5 + Diesel Oil	1.0 mg/ml each
Solvent : n-Hexane					6 ampoules X 1 ml			
					6 ampoules X 5 ml			
					Ref : F061955			
					Price : EF			
					Ref : F109681			
					Price : BDA			

Quality Control Standard of Mineral Oils - 2 components

Mineral Oil	CAS:8042-47-5	Diesel Oil	N/A (id : 381)		
Solvent : Acetone		ampoule 1 ml	500 µg/ml each	Ref : P858050	Price : DF
		ampoule 5 ml	500 µg/ml each	Ref : F109641	Price : KJ
		ampoule 1 ml	1 000 µg/ml each	Ref : F061965	Price : DF
		ampoule 5 ml	1 000 µg/ml each	Ref : F145091	Price : KJ
Solvent : n-Hexane		ampoule 1 ml	5 000 µg/ml each	Ref : P858010	Price : DF
		ampoule 5 ml	5 000 µg/ml each	Ref : F109661	Price : KJ
		ampoule 1 ml	10 000 µg/ml each	Ref : F061975	Price : DF
		ampoule 5 ml	10 000 µg/ml each	Ref : F109671	Price : KJ

Florisol Cartridge QC Standard Mixture - 2 components

Mineral Oil	CAS:8042-47-5	Diesel Oil	N/A (id : 381)		
Solvent : n-Hexane		ampoule 1 ml	1 000 µg/ml each	Ref : P858090	Price : DA
		ampoule 5 ml	1 000 µg/ml each	Ref : F109651	Price : FF

Extraction Solvent Stock Solution - 2 components

n-Decane (C10)	CAS:124-18-5	20 µl/l	n-Tetracontane (C40)	CAS:4181-95-7	20 mg/l
Solvent : n-Hexane		ampoule 1 ml		Ref : P858570	Price : DA
		ampoule 10 ml		Ref : F109741	Price : FE
Solvent : n-Hexane/Petroleum ether (1/1)		ampoule 1 ml		Ref : P858530	Price : DA
		ampoule 10 ml		Ref : F109751	Price : FE

Stearyl Stearate Test Solution – 1 component

Stearyl stearate	CAS:2778-96-3				
Solvent : Cyclohexane		ampoule 1 ml	2 000 µg/ml	Ref : P858490	Price : CF
		ampoule 10 ml	2 000 µg/ml	Ref : P858500	Price : JD
Solvent : n-Hexane		ampoule 1 ml	2 000 µg/ml	Ref : P858450	Price : CF
		ampoule 10 ml	2 000 µg/ml	Ref : F109761	Price : JD

BTEX Standard Solution - 6 components

Benzene	CAS:71-43-2	Toluene	CAS:108-88-3	m-Xylene	CAS:108-38-3
Ethylbenzene	CAS:100-41-4	o-Xylene	CAS:95-47-6	p-Xylene	CAS:106-42-3
Solvent : n-Hexane		ampoule 1 ml	200 µg/ml	Ref : F061985	Price : DF
		ampoule 5 ml	200 µg/ml	Ref : F109801	Price : KJ
		ampoule 1 ml	2 000 µg/ml	Ref : F061995	Price : DF
		ampoule 5 ml	2 000 µg/ml	Ref : F109811	Price : KJ

n-Alkane Standard Solution - 6 components

n-Decane (C10)	CAS:124-18-5	n-Eicosane (C20)	CAS:112-95-8	n-Triacontane (C30)	CAS:638-68-6
n-Hexadecane (C16)	CAS:544-76-3	n-Tetracosane (C24)	CAS:646-31-1	n-Tetracontane (C40)	CAS:4181-95-7
Solvent : n-Hexane		ampoule 1 ml	100 µg/ml	Ref : F062005	Price : DA
		ampoule 5 ml	100 µg/ml	Ref : F109791	Price : KJ

n-Alkane Standard Solution - 16 components C10÷C40 (all even)

n-Decane (C10)	CAS:124-18-5	n-Docosane (C22)	CAS:629-97-0	n-Tetratriacontane (C34)	CAS:14167-59-0
n-Dodecane (C12)	CAS:112-40-3	n-Tetracosane (C24)	CAS:646-31-1	n-Hexatriacontane (C36)	CAS:630-06-8
n-Tetradecane (C14)	CAS:629-59-4	n-Hexacosane (C26)	CAS:630-01-3	n-Octatriacontane (C38)	CAS:7194-85-6
n-Hexadecane (C16)	CAS:544-76-3	n-Octacosane (C28)	CAS:630-02-4	n-Tetracontane (C40)	CAS:4181-95-7
n-Octadecane (C18)	CAS:593-45-3	n-Triacontane (C30)	CAS:638-68-6		
n-Eicosane (C20)	CAS:112-95-8	n-Dotriacontane (C32)	CAS:544-85-4		
Solvent : n-Hexane/Petroleum ether (1/1)		ampoule 1 ml	100 µg/ml	Ref : P858130	Price : GA
		ampoule 5 ml	100 µg/ml	Ref : F109771	Price : BGA
Solvent : n-Hexane		ampoule 1 ml	50 µg/ml	Ref : P858170	Price : GA
		ampoule 5 ml	50 µg/ml	Ref : F109781	Price : BGA

n-Alkane Standard Solution - 19 components C7+C8+C9+C10÷C40(all even)

n-Heptane (C7)	CAS:142-82-5	n-Octadecane (C18)	CAS:593-45-3	n-Dotriacontane (C32)	CAS:544-85-4
n-Octane (C8)	CAS:111-65-9	n-Eicosane (C20)	CAS:112-95-8	n-Tetratriacontane (C34)	CAS:14167-59-0
n-Nonane (C9)	CAS:111-84-2	n-Docosane (C22)	CAS:629-97-0	n-Hexatriacontane (C36)	CAS:630-06-8
n-Decane (C10)	CAS:124-18-5	n-Tetracosane (C24)	CAS:646-31-1	n-Octatriacontane (C38)	CAS:7194-85-6
n-Dodecane (C12)	CAS:112-40-3	n-Hexacosane (C26)	CAS:630-01-3	n-Tetracontane (C40)	CAS:4181-95-7
n-Tetradecane (C14)	CAS:629-59-4	n-Octacosane (C28)	CAS:630-02-4		
n-Hexadecane (C16)	CAS:544-76-3	n-Triacontane (C30)	CAS:638-68-6		
Solvent : Iso-Octane		ampoule 1 ml	50 µg/ml	Ref : F062015	Price : GF
		ampoule 5 ml	50 µg/ml	Ref : F109821	Price : BJA

n-Alkane Standard Solution - 25 components C7+C8+C9+C10+C40+ BTEX

n-Heptane (C7)	CAS:142-82-5	n-Docosane (C22)	CAS:629-97-0	n-Tetracontane (C40)	CAS:4181-95-7
n-Octane (C8)	CAS:111-65-9	n-Tetracosane (C24)	CAS:646-31-1	Benzene	CAS:71-43-2
n-Nonane (C9)	CAS:111-84-2	n-Hexacosane (C26)	CAS:630-01-3	Ethylbenzene	CAS:100-41-4
n-Decane (C10)	CAS:124-18-5	n-Octacosane (C28)	CAS:630-02-4	Toluene	CAS:108-88-3
n-Dodecane (C12)	CAS:112-40-3	n-Triacontane (C30)	CAS:638-68-6	m-Xylene	CAS:108-38-3
n-Tetradecane (C14)	CAS:629-59-4	n-Dotriacontane (C32)	CAS:544-85-4	o-Xylene	CAS:95-47-6
n-Hexadecane (C16)	CAS:544-76-3	n-Tetatriacontane (C34)	CAS:14167-59-0	p-Xylene	CAS:106-42-3
n-Octadecane (C18)	CAS:593-45-3	n-Hexatriacontane (C36)	CAS:630-06-8		
n-Eicosane (C20)	CAS:112-95-8	n-Octatriacontane (C38)	CAS:7194-85-6		
Solvent : n-Hexane		ampoule 1 ml	50 µg/ml	Ref : F062025	Price : GA
		ampoule 5 ml	50 µg/ml	Ref : F109851	Price : BHA

ISO 9562 :2004**Water quality - Determination of absorbable organically bound halogens (AOX)**

This European standard specifies a method for the direct determination of an amount of more than 10 microgram/l in water of organically bound chlorine, bromine and iodine (expressed as chloride) absorbable on activated carbon. The concentration of inorganic chloride ions in the test sample (see clause 8) needs to be less than 1 g/l. Samples with higher concentrations need to be diluted prior to analysis. For samples containing suspended solids, halogens adsorbed onto the solid matter are also included.

Stock Solution - 1 component

4-Chlorophenol	CAS:106-48-9				
Solvent : Water		ampoule 1 ml	725 µg/ml	Ref : F128981	Price : CF

Stock Solution - 1 component

2-Chlorobenzoic acid	CAS:118-91-2				
Solvent : Water		ampoule 1 ml	1 104 µg/ml	Ref : F128991	Price : BH

ISO 10301 :1997**Water quality - Determination of highly volatile halogenated hydrocarbons - Gas-chromatographic methods**

This International Standard specifies two methods for the determination of highly volatile halogenated hydrocarbons in drinking water, ground water, swimming pool water, most rivers and lakes and many sewage and industrial effluents using gas chromatography.

VOC Standard Solution - 22 components

Bromochloromethane	CAS:74-97-5	1,1-Dichloroethene	CAS:75-35-4	Tetrachloroethene	CAS:127-18-4
Bromodichloromethane	CAS:75-27-4	cis-1,2-Dichloroethene	CAS:156-59-2	Tetrachloromethane	CAS:56-23-5
Chloroform	CAS:67-66-3	trans-1,2-Dichloroethene	CAS:156-60-5	Tribromomethane	CAS:75-25-2
Dibromochloromethane	CAS:124-48-1	Dichloromethane	CAS:75-09-2	1,1,1-Trichloroethane	CAS:71-55-6
1,2-Dibromoethane	CAS:106-93-4	1,2-Dichloropropane	CAS:78-87-5	1,1,2-Trichloroethane	CAS:79-00-5
Dibromomethane	CAS:74-95-3	1,3-Dichloropropane	CAS:142-28-9	Trichloroethene	CAS:79-01-6
1,1-Dichloroethane	CAS:75-34-3	cis-1,3-Dichloropropene	CAS:10061-01-5		
1,2-Dichloroethane	CAS:107-06-2	trans-1,3-Dichloropropene	CAS:10061-02-6		
Solvent : Iso-Octane		ampoule 10 ml	10 µg/ml	Ref : F128961	Price : GF

ISO 10382:2002**Soil quality - Determination of organochlorine pesticides and polychlorinated biphenyls - Gas-chromatographic method with electron capture detection**

Specifies a method for quantitative determination of seven polychlorinated biphenyls and seventeen organochlorine pesticides in soil. ISO 10382:2002 is applicable to all types of soil.

Multicomponent Stock Solution of PCBs - 7 components

PCB 28	CAS:7012-37-5	PCB 118	CAS:31508-00-6	PCB 180	CAS:35065-29-3
PCB 52	CAS:35693-99-3	PCB 138	CAS:35065-28-2		
PCB 101	CAS:37680-73-2	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 5 ml	100 µg/ml	Ref : F109831	Price : BEA

Multicomponent Stock Solution of PCBs & Organochlorine Pesticides - 24 components

PCB 28	CAS:7012-37-5	Alpha-HCH	CAS:319-84-6	Heptachlor-exo-epoxide	CAS:1024-57-3
PCB 52	CAS:35693-99-3	Beta-HCH	CAS:319-85-7	Endosulfan-alpha	CAS:959-98-8
PCB 101	CAS:37680-73-2	Gamma-HCH (Lindane)	CAS:58-89-9	4,4'-DDE	CAS:72-55-9
PCB 118	CAS:31508-00-6	Aldrin	CAS:309-00-2	2,4'-DDD	CAS:53-19-0
PCB 138	CAS:35065-28-2	Dieldrin	CAS:60-57-1	2,4'-DDT	CAS:789-02-6
PCB 153	CAS:35065-27-1	Endrin	CAS:72-20-8	4,4'-DDD (TDE)	CAS:72-54-8
PCB 180	CAS:35065-29-3	Heptachlor	CAS:76-44-8	2,4'-DDE	CAS:3424-82-6
Hexachlorobenzene	CAS:118-74-1	Heptachlor-endo-epoxide	CAS:28044-83-9	4,4'-DDT	CAS:50-29-3
Solvent : Iso-Octane		ampoule 1 ml	100 µg/ml	Ref : F109841	Price : KF

ISO 10695 :2000

Water quality - Determination of selected organic nitrogen and phosphorus compounds - Gas chromatographic methods

This International Standard specifies two methods for the determination of certain organic nitrogen and phosphorus compounds in waters by gas chromatography.

NCC / OPP Standard Solution - 12 components

Atrazine	CAS:1912-24-9	Parathion-methyl	CAS:298-00-0	Simazine (CAT)	CAS:122-34-9
Cyanazine	CAS:21725-46-2	Pendimethalin	CAS:40487-42-1	Terbutylazine	CAS:5915-41-3
Metazachlor	CAS:67129-08-2	Propazine	CAS:139-40-2	Trifluralin	CAS:1582-09-8
Parathion-ethyl	CAS:56-38-2	Sebutylazine	CAS:7286-69-3	Vinclozolin	CAS:50471-44-8
Solvent : Acetone		ampoule 10 ml	10 µg/ml	Ref : F128971	Price : EF

ISO 11369:1997

Water quality - Determination of selected plant treatment agents - Method using high performance liquid chromatography with UV detection after solid-liquid extraction

This International Standard describes a method for the determination of organic plant treatment agents in drinking and ground water using high performance liquid chromatography (HPLC) with UV detection after solid-liquid extraction.

NCC Standard Solution - 17 components

Atrazine	CAS:1912-24-9	Isoproturon	CAS:34123-59-6	Metoxuron	CAS:19937-59-8
Atrazine-desethyl	CAS:6190-65-4	Linuron	CAS:330-55-2	Monolinuron	CAS:1746-81-2
Chlortoluron	CAS:15545-48-9	Metazachlor	CAS:67129-08-2	Sebutylazine	CAS:7286-69-3
Cyanazine	CAS:21725-46-2	Methabenzthiazuron	CAS:18691-97-9	Simazine (CAT)	CAS:122-34-9
Diuron	CAS:330-54-1	Metobromuron	CAS:3060-89-7	Terbutylazine	CAS:5915-41-3
Hexazinone	CAS:51235-04-2	Metolachlor	CAS:51218-45-2		
Solvent : Acetonitrile		ampoule 10 ml	10 µg/ml	Ref : F129021	Price : EF

ISO 11423-2:1997

Water quality - Determination of benzene and some derivatives - Part 2: Method using extraction and gas chromatography

The method described is applicable to the determination of benzene, methylbenzene (toluene), dimethylbenzenes (xylenes) and ethylbenzene (abbreviated hereafter to BTX) in water and waste water in concentrations above 5 µg/l. High concentrations may be determined by diluting the extract. A number of further derivatives and nonpolar compounds with similar boiling points may also be determined by this method.

ACs Internal Standard Solution - 1 component

Toluene D8	CAS:2037-26-5				
Solvent : Acetone		ampoule 1 ml	50 µg/ml	Ref : F129011	Price : BK

ACs Standard Solution - 6 components

Benzene	CAS:71-43-2	Toluene	CAS:108-88-3	m-Xylene	CAS:108-38-3
Ethylbenzene	CAS:100-41-4	o-Xylene	CAS:95-47-6	p-Xylene	CAS:106-42-3
Solvent : Acetone		ampoule 5 ml	50 ml/l	Ref : F129001	Price : EF
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F129231	Price : CD

ISO 13877:1998

Soil quality - Determination of polynuclear aromatic hydrocarbons - Method using high performance liquid chromatography

PAH Standard Stock Solution - 16 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Acenaphthylene	CAS:208-96-8	Benzo(k)fluoranthene	CAS:207-08-9	Naphthalene	CAS:91-20-3
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Phenanthrene	CAS:85-01-8
Benzo(a)anthracene	CAS:56-55-3	Dibenzo(a,h)anthracene	CAS:53-70-3	Pyrene	CAS:129-00-0
Benzo(a)pyrene	CAS:50-32-8	Fluoranthene	CAS:206-44-0		
Benzo(b)fluoranthene	CAS:205-99-2	Fluorene	CAS:86-73-7		
Solvent : Acetonitrile		ampoule 5 ml	10 µg/ml	Ref : F109871	Price : GF
		ampoule 1 ml	100 µg/ml	Ref : F109861	Price : GF

ISO 14154:2005

Soil quality - Determination of some selected chlorophenols - Gas chromatographic method with electron-capture detection

ISO 14154:2005 describes the gas chromatographic determination of 15 chlorophenols in soil samples. This method can also be applied to other solid samples such as sediments and solid wastes.

ISO 14154:2005 describes an acid-base liquid extraction, followed by acetylation and then liquid/liquid extraction. Determination of mass concentration is then carried out by gas chromatography and electron-capture detection.

Internal Standard Stock Solution - 1 component

2,4,6-Tribromophenol	CAS:118-79-6	ampoule 5 ml	1.148 mg/ml	Ref : F109891	Price : CC
Solvent : Ethanol					

Chlorophenol Standard, Stock Solution - 15 components

2,3-Dichlorophenol	CAS:576-24-9	400 µg/ml	2,3,4,6-Tetrachlorophenol	CAS:58-90-2	600 µg/ml
2,4-Dichlorophenol	CAS:120-83-2	400 µg/ml	2,3,4-Trichlorophenol	CAS:15950-66-0	400 µg/ml
2,5-Dichlorophenol	CAS:583-78-8	400 µg/ml	2,3,5-Trichlorophenol	CAS:933-78-8	400 µg/ml
2,6-Dichlorophenol	CAS:87-65-0	400 µg/ml	2,3,6-Trichlorophenol	CAS:933-75-5	400 µg/ml
3,4-Dichlorophenol	CAS:95-77-2	400 µg/ml	2,4,5-Trichlorophenol	CAS:95-95-4	400 µg/ml
3,5-Dichlorophenol	CAS:591-35-5	400 µg/ml	2,4,6-Trichlorophenol	CAS:88-06-2	600 µg/ml
Pentachlorophenol	CAS:87-86-5	1 000 µg/ml	3,4,5-Trichlorophenol	CAS:609-19-8	200 µg/ml
2,3,4,5-Tetrachlorophenol	CAS:4901-51-3	200 µg/ml			
Solvent : Ethanol		ampoule 1 ml		Ref : F109881	Price : HF

Chlorophenyl Acetate Standard, Stock Solution - 15 components

2,3-Dichlorophenol acetate	CAS:61925-85-7	503 µg/ml	2,3,5-Trichlorophenol acetate	CAS:61925-88-0	485 µg/ml
2,4-Dichlorophenol acetate	CAS:6341-97-5	503 µg/ml	2,3,6-Trichlorophenol acetate	CAS:61925-87-9	485 µg/ml
2,5-Dichlorophenol acetate	CAS:30124-46-0	503 µg/ml	2,4,5-Trichlorophenol acetate	CAS:5393-75-9	485 µg/ml
2,6-Dichlorophenol acetate	CAS:28165-71-1	503 µg/ml	2,4,6-Trichlorophenol acetate	CAS:23399-90-8	728 µg/ml
3,4-Dichlorophenol acetate	CAS:17847-51-7	503 µg/ml	3,4,5-Trichlorophenol acetate	id : 15591	243 µg/ml
3,5-Dichlorophenol acetate	CAS:61925-86-8	503 µg/ml	2,3,4,5-Tetrachlorophenol acetate	id : 14700	250 µg/ml
Pentachlorophenol acetate	CAS:1441-02-7	1 158 µg/ml	2,3,4,6-Tetrachlorophenol acetate	id : 14720	750 µg/ml
2,3,4-Trichlorophenol acetate	id : 15561	485 µg/ml			
Solvent : Ethanol		ampoule 1 ml		Ref : F109901	Price : CBA

ISO 15009:2002

Soil quality - Gas chromatographic determination of the content of volatile aromatic hydrocarbons, naphthalene and volatile halogenated hydrocarbons - Purge-and-trap method with thermal desorption

Internal Stock Standard Solution for aromatic hydrocarbons - 2 components

Ethylbenzene D10	CAS:25837-05-2	Toluene D8	CAS:2037-26-5		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062035	Price : CF
		ampoule 5 ml	2 000 µg/ml	Ref : F109911	Price : HF

Internal Stock Standard Solution for volatile halogenated hydrocarbons - 3 components

2-Bromofluorobenzene	CAS:1072-85-1	1,4-Dichlorobutane	CAS:110-56-5	α, α, α -Trifluorotoluene	CAS:98-08-8
Solvent : Methanol purge-&-trap		ampoule 1 ml	2 000 µg/ml	Ref : F062045	Price : CA
		ampoule 5 ml	2 000 µg/ml	Ref : F109921	Price : GF

Volatile Aromatic Hydrocarbons - 8 components

Benzene	CAS:71-43-2	Styrene	CAS:100-42-5	m-Xylene	CAS:108-38-3
Ethylbenzene	CAS:100-41-4	Toluene	CAS:108-88-3	p-Xylene	CAS:106-42-3
Naphthalene	CAS:91-20-3	o-Xylene	CAS:95-47-6		
Solvent : Methanol purge & trap		ampoule 1 ml	4 000 µg/ml	Ref : F062055	Price : CF
		ampoule 5 ml	4 000 µg/ml	Ref : F109931	Price : JA

Volatile Halogenated Hydrocarbons - 18 components

3-Chloropropene	CAS:107-05-1	Dichloromethane	CAS:75-09-2	Tetrachloromethane	CAS:56-23-5
1,2-Dichlorobenzene	CAS:95-50-1	1,2-Dichloropropane	CAS:78-87-5	1,1,1-Trichloroethane	CAS:71-55-6
1,1-Dichloroethane	CAS:75-34-3	cis-1,3-Dichloropropene	CAS:10061-01-5	1,1,2-Trichloroethane	CAS:79-00-5
1,2-Dichloroethane	CAS:107-06-2	trans-1,3-Dichloropropene	CAS:10061-02-6	Trichloroethene	CAS:79-01-6
cis-1,2-Dichloroethene	CAS:156-59-2	Monochlorobenzene	CAS:108-90-7	Trichloromethane	CAS:67-66-3
trans-1,2-Dichloroethene	CAS:156-60-5	Tetrachloroethene	CAS:127-18-4	1,2,3-Trichloropropane	CAS:96-18-4
Solvent : Methanol purge & trap		ampoule 1 ml	4 000 µg/ml	Ref : F062065	Price : BCA
		ampoule 5 ml	4 000 µg/ml	Ref : F109941	Price : EJA

ISO 15680:2003

Water quality - Gas chromatographic determination of a number of monocyclic aromatic hydrocarbons, naphthalene and several chlorinated compounds using purge-and-trap and thermal desorption

ISO 15680:2003 specifies a general method for the determination of volatile organic compounds (VOCs) in water by purge-and-trap isolation and gas chromatography (GC). ISO 15680:2003 is applicable to drinking water, ground water, surface water, seawater and to (diluted) waste water.

ACs Internal Stock Standard Solution - 4 components

1,4-Difluorobenzene	CAS:540-36-3	Fluorobenzene	CAS:462-06-6	Monochlorobenzene D5	CAS:3114-55-4
1,4-Dichlorobenzene D4	CAS:3855-82-1				
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062075	Price : DF
		ampoule 5 ml	2 000 µg/ml	Ref : F109951	Price : BCA

VOC Stock Calibration Standard Solution 1- 14 components

Bromochloromethane	CAS:74-97-5	1,1-Dichloroethane	CAS:75-34-3	2,2-Dichloropropane	CAS:594-20-7
Bromomethane	CAS:74-83-9	1,1-Dichloroethene	CAS:75-35-4	1,1,1-Trichloroethane	CAS:71-55-6
Chloroethane	CAS:75-00-3	cis-1,2-Dichloroethene	CAS:156-59-2	Trichlorofluoromethane	CAS:75-69-4
Chloromethane	CAS:74-87-3	trans-1,2-Dichloroethene	CAS:156-60-5	Trichloromethane	CAS:67-66-3
Dichlorodifluoromethane	CAS:75-71-8	Dichloromethane	CAS:75-09-2		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062085	Price : JF
		ampoule 5 ml	2 000 µg/ml	Ref : F109961	Price : DCA

VOC Stock Calibration Standard Solution 2 - 13 components

Benzene	CAS:71-43-2	1,2-Dichloropropane	CAS:78-87-5	Toluene	CAS:108-88-3
Bromodichloromethane	CAS:75-27-4	1,1-Dichloropropene	CAS:563-58-6	1,1,2-Trichloroethane	CAS:79-00-5
1,2-Dibromoethane	CAS:106-93-4	cis-1,3-Dichloropropene	CAS:10061-01-5	Trichloroethene	CAS:79-01-6
Dibromomethane	CAS:74-95-3	trans-1,3-Dichloropropene	CAS:10061-02-6		
1,2-Dichloroethane	CAS:107-06-2	Tetrachloromethane	CAS:56-23-5		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062095	Price : JA
		ampoule 5 ml	2 000 µg/ml	Ref : F109971	Price : CJA

VOC /ACs Stock Calibration Standard Solution 3 - 11 components

Dibromochloromethane	CAS:124-48-1	Styrene	CAS:100-42-5	m-Xylene	CAS:108-38-3
1,3-Dichloropropane	CAS:142-28-9	1,1,1,2-Tetrachloroethane	CAS:630-20-6	o-Xylene	CAS:95-47-6
Ethylbenzene	CAS:100-41-4	Tetrachloroethene	CAS:127-18-4	p-Xylene	CAS:106-42-3
Monochlorobenzene	CAS:108-90-7	Tribromomethane	CAS:75-25-2		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062105	Price : CA
		ampoule 5 ml	2 000 µg/ml	Ref : F109981	Price : JA

Stock Calibration Standard Solution 4 - 22 components

Bromobenzene	CAS:108-86-1	1,3-Dichlorobenzene	CAS:541-73-1	1,2,3-Trichlorobenzene	CAS:87-61-6
n-Butylbenzene	CAS:104-51-8	1,4-Dichlorobenzene	CAS:106-46-7	1,2,4-Trichlorobenzene	CAS:120-82-1
sec-Butylbenzene	CAS:135-98-8	Hexachlorobutadiene	CAS:87-68-3	1,3,5-Trichlorobenzene	CAS:108-70-3
tert-Butylbenzene	CAS:98-06-6	Isopropylbenzene	CAS:98-82-8	1,2,3-Trichloropropane	CAS:96-18-4
2-Chlorotoluene	CAS:95-49-8	4-Isopropyltoluene	CAS:99-87-6	1,2,4-Trimethylbenzene	CAS:95-63-6
4-Chlorotoluene	CAS:106-43-4	Naphthalene	CAS:91-20-3	1,3,5-Trimethylbenzene	CAS:108-67-8
1,2-Dibromo-3-chloropropane	CAS:96-12-8	n-Propylbenzene	CAS:103-65-1		
1,2-Dichlorobenzene	CAS:95-50-1	1,1,2,2-Tetrabromoethane	CAS:79-27-6		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062155	Price : DF
		ampoule 5 ml	2 000 µg/ml	Ref : F109991	Price : BDA

ISO 15753:2006**Animal and vegetable fats and oils - Determination of polycyclic aromatic hydrocarbons**

ISO 15753:2006 describes two methods for the determination of 15 polycyclic aromatic hydrocarbons (PAHs) in animal and vegetable fats and oils: a general method, and a method specific for coconut oil and short-chain vegetable oils.

These methods are not quantitative for the very volatile compounds such as naphthalene, acenaphthene and fluorene. Due to interferences provided by the matrix itself, palm oil and olive pomace oil cannot be analysed using this method.

Multiple Compound Standard Solution - 15 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Fluorene	CAS:86-73-7
Anthracene	CAS:120-12-7	Benzo(k)fluoranthene	CAS:207-08-9	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Benzo(a)anthracene	CAS:56-55-3	Chrysene	CAS:218-01-9	Naphthalene	CAS:91-20-3
Benzo(a)pyrene	CAS:50-32-8	Dibenzo(a,h)anthracene	CAS:53-70-3	Phenanthrene	CAS:85-01-8
Benzo(b)fluoranthene	CAS:205-99-2	Fluoranthene	CAS:206-44-0	Pyrene	CAS:129-00-0
Solvent : Acetonitrile		ampoule 1 ml	200 µg/ml	Ref : F110011	Price : HF
Solvent : Acetonitrile/THF/Methanol (1/1/1)		ampoule 10 ml	200 µg/ml	Ref : F146591	Price : CGA
Solvent : Toluene		ampoule 1 ml	100 µg/ml	Ref : F110001	Price : GF

ISO 15913:2000**Water quality - Determination of selected phenoxyalkanoic herbicides, including bentazones and hydroxybenzonitriles by gas chromatography and mass spectrometry after solid phase extraction and derivatization.**

This standard specifies a method for the determination of phenoxyalkanoic acids in ground and drinking water in mass concentrations > 0,05 Bµg/l. The single substances are listed in the above mentioned standard. The method includes a solid-liquid extraction, derivatization and gas chromatography with MS detection.

Phenoxyacetic Herbicides Standard Solution - 10 components

Bentazone	CAS:25057-89-0	Dichlorprop	CAS:120-36-5	Mecoprop (MCP)	CAS:7085-19-0
Bromoxynil	CAS:1689-84-5	Fenoprop	CAS:93-72-1	2,4,5-Trichlorophenoxy acetic acid	CAS:93-76-5
2,4-D	CAS:94-75-7	MCPA	CAS:94-74-6		
2,4-DB	CAS:94-82-6	MCPB	CAS:94-81-5		
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F129031	Price : EF

Phenoxyalkanoic Carbonic Esters Standard Solution - 10 components

Bentazone methyl ester	id : 17751	Dichlorprop methyl ester	CAS:57153-17-0	Mecoprop methyl ester	CAS:23844-56-6
Bromoxynil methyl ester	id : 17761	Fenoprop methyl ester	CAS:4841-20-7	2,4,5-Trichlorophenoxy	CAS:1928-37-6
2,4-D methyl ester	CAS:1928-38-7	MCPA methyl ester	CAS:2436-73-9	acetic acid	
2,4-DB methyl ester	CAS:18625-12-2	MCPB methyl ester	CAS:57153-18-1		
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F129041	Price : EE

ISO 17495:2001**Water quality - Determination of selected nitrophenols - Method by solid-phase extraction and gas chromatography with mass spectrometric detection**

This standard specifies a method for the determination of selected nitrophenols in drinking, ground and surface water in mass concentrations > 0,5 Bµg/l. The method includes extraction, derivatization and gas chromatography. Exact figures on the application ranges are given in the table listing the result from the interlaboratory trial. The method includes the determination of mononitrophenols, methylnitrophenols, and selected chloronitrophenols and methylated chloronitrophenols.

Phenols Internal Standard Solution - 2 components

2,4-Dibromophenol	CAS:615-58-7	2,3,6-Trichlorophenol	CAS:933-75-5		
Solvent : Ethyl Acetate		ampoule 1 ml	10 µg/ml	Ref : F129061	Price : CF

Phenols Standard Solution - 14 components

2,4-Dichloro-6-nitrophenol	CAS:609-89-2	2,5-Dinitrophenol	CAS:329-71-5	5-Methyl-2-nitrophenol	CAS:700-38-9
2,6-Dichloro-4-nitrophenol	CAS:618-80-4	2,6-Dinitrophenol	CAS:573-56-8	2-Nitrophenol	CAS:88-75-5
2,6-Dimethyl-4-nitrophenol	CAS:2423-71-4	3-Methyl-2-nitrophenol	CAS:4920-77-8	3-Nitrophenol	CAS:554-84-7
2,4-Dinitro-6-methylphenol	CAS:534-52-1	3-Methyl-4-nitrophenol	CAS:2581-34-2	4-Nitrophenol	CAS:100-02-7
2,4-Dinitrophenol	CAS:51-28-5	4-Methyl-2-nitrophenol	CAS:119-33-5		
Solvent : Ethyl Acetate		ampoule 1 ml	500 µg/ml	Ref : F129051	Price : DJ

ISO 17858:2007**Water quality - Determination of dioxin-like polychlorinated biphenyls - Method using gas chromatography/mass spectrometry**

ISO 17858:2007 specifies a method for the determination of dioxin-like tetra- to hepta-chlorinated biphenyls (PCBs) in waters and wastewaters (containing less than 1 % suspended solids) using high-resolution gas chromatography/high-resolution mass spectrometry (HRGC/HRMS). The minimum levels at which the dioxin-like PCBs can typically be determined with no interferences present are specified. This method is "performance based". The analyst is permitted to modify the method to overcome interferences or lower the cost of measurements, provided that all performance criteria in this method are met. The requirements for establishing method equivalency are given.

Precision and Recovery PCBs Standard Solution -13 components

PCB 77	CAS:32598-13-3	PCB 126	CAS:57465-28-8	PCB 170	CAS:35065-30-6
PCB 81	CAS:70362-50-4	PCB 156	CAS:38380-08-4	PCB 180	CAS:35065-29-3
PCB 105	CAS:32598-14-4	PCB 157	CAS:69782-90-7	PCB 189	CAS:39635-31-9
PCB 114	CAS:74472-37-0	PCB 167	CAS:52663-72-6		
PCB 123	CAS:65510-44-3	PCB 169	CAS:32774-16-6		
Solvent : n-Nonane		ampoule 1 ml	20 µg/l	Ref : F062125	Price : DF
		ampoule 5 ml	20 µg/l	Ref : F110031	Price : BCA

ISO 17993:2002**Water quality - Determination of 15 polycyclic aromatic hydrocarbons (PAH) in water by HPLC with fluorescence detection after liquid-liquid extraction**

ISO 17993:2002 specifies a method using high performance liquid chromatography (HPLC) with fluorescence detection for the determination of 15 selected polycyclic aromatic hydrocarbons (PAH) in drinking and ground water in mass concentrations greater than 0,005 microgram/litre (for each single compound) and surface waters in mass concentrations above 0,01 microgram/litre.

This method is, with some modification, also suitable for the analysis of waste water.

Multiple Compound Stock Solution - 15 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Fluorene	CAS:86-73-7
Anthracene	CAS:120-12-7	Benzo(k)fluoranthene	CAS:207-08-9	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Benzo(a)anthracene	CAS:56-55-3	Chrysene	CAS:218-01-9	Naphthalene	CAS:91-20-3
Benzo(a)pyrene	CAS:50-32-8	Dibenzo(a,h)anthracene	CAS:53-70-3	Phenanthrene	CAS:85-01-8
Benzo(b)fluoranthene	CAS:205-99-2	Fluoranthene	CAS:206-44-0	Pyrene	CAS:129-00-0
Solvent : Acetonitrile		ampoule 5 ml	10 µg/ml	Ref : F110041	Price : GF

ISO 18287:2006

Soil quality - Determination of polycyclic aromatic hydrocarbons (PAH) - Gas chromatographic method with mass spectrometric detection (GC-MS)

ISO 18287:2006 specifies the quantitative determination of 16 polycyclic aromatic hydrocarbons (PAH) according to the priority list of the Environmental Protection Agency, USA (EPA, 1982). ISO 18287:2006 is applicable to all types of soil (field-moist or chemically dried samples), covering a wide range of PAH contamination levels. Under the conditions specified in ISO 18287:2006, a lower limit of application of 0,01 mg/kg (expressed as dry matter) can be ensured for each individual PAH.

PAH Standard Solution -16 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Acenaphthylene	CAS:208-96-8	Benzo(k)fluoranthene	CAS:207-08-9	Naphthalene	CAS:91-20-3
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Phenanthrene	CAS:85-01-8
Benzo(a)anthracene	CAS:56-55-3	Dibenzo(a,h)anthracene	CAS:53-70-3	Pyrene	CAS:129-00-0
Benzo(a)pyrene	CAS:50-32-8	Fluoranthene	CAS:206-44-0		
Benzo(b)fluoranthene	CAS:205-99-2	Fluorene	CAS:86-73-7		
Solvent : Cyclohexane		ampoule 1 ml	100 µg/ml	Ref : F062135	Price : FF
		ampoule 1 ml	200 µg/ml	Ref : F131381	Price : JA

ISO 18856:2004

Water quality - Determination of selected phthalates using gas chromatography/mass spectrometry

This future standard will specify a method for the determination of phthalates in water after solid phase extraction and gas chromatography - mass spectrometry. The method is applicable to the determination of phthalates (a list of common phthalates is included in ISO/CD 18856) in ground water, surface water, waste water and drinking water in mass concentrations between 0,02 Bµg/l and 0,150 Bµg/l, depending on the individual substance. The applicability of the method to other phthalates not given in the table is not excluded, but shall be verified in each case.

Phthalates Standard Solution -11 components

Butyl benzyl phthalate	CAS:85-68-7	Di-2-ethylhexylphthalate	CAS:117-81-7	Di-n-octylphthalate	CAS:117-84-0
Dibutylphthalate	CAS:84-74-2	Diethyl phthalate	CAS:84-66-2	Diundecyl phthalate	CAS:3648-20-2
Dicyclohexyl phthalate	CAS:84-61-7	Di-iso-butylphthalate	CAS:84-69-5	Phthalic acid, bis-propyl ester	CAS:131-16-8
Didecyl phthalate	CAS:84-77-5	Dimethyl phthalate	CAS:131-11-3		
Solvent : Ethyl Acetate		ampoule 1 ml	1 000 µg/ml	Ref : F129191	Price : DC

ISO 18857-1:2005

Water quality - Determination of selected alkylphenols - Part 1: Method for non-filtered samples using liquid-liquid extraction and gas chromatography with mass selective detection

This part of ISO 18857 specifies a method for the determination of 4-nonylphenol (mixture of isomers) and 4-(1,1,3,3-tetramethylbutyl)phenol in non-filtered samples of drinking water, ground water and surface water. The method is applicable in a concentration range from 0,005 Bµg/l to 0,2 Bµg/l for 4-(1,1,3,3-tetramethylbutyl)phenol and from 0,02 Bµg/l to 0,2 Bµg/l for 4-nonylphenol (mixture of isomers). Depending on the matrix, the method is also applicable to waste water in a concentration range from 0,1 Bµg/l to 50 Bµg/l. Higher concentrations may be measured after appropriate dilution of the sample.

Phenol Calibration Standard - 1 component

Solvent : Toluene

4-Nonylphenol (technical)	CAS:84852-15-3	ampoule 1 ml	100 µg/ml	Ref : F129201	Price : BG
4-tert- Octylphenol	CAS:140-66-9	ampoule 1 ml	100 µg/ml	Ref : F129211	Price : BG

ISO 23631 :2006

Water quality - Determination of dalapon, trichloroacetic acid and selected haloacetic acids - Method using gas chromatography (GC-ECD and/or GC-MS detection) after liquid-liquid extraction and derivatization

ISO 23631:2006 specifies a method for the determination of dalapon, trichloroacetic acid (TCA) and selected haloacetic acids in ground water and drinking water by gas chromatography (GC-ECD and/or GC-MS detection) after liquid-liquid-extraction and derivatization using diazomethane. Depending on the matrix, the method is applicable to a concentration range from 0,5 to 10 micrograms per litre.

Standard Solution - 7 components

Bromochloroacetic Acid	CAS:5589-96-8	Dichloroacetic Acid	CAS:79-43-6	Trichloroacetic acid	CAS:76-03-9
Dalapon	CAS:75-99-0	Monobromoacetic Acid	CAS:79-08-3		
Dibromoacetic acid	CAS:631-64-1	Monochloroacetic Acid	CAS:79-11-8		
Solvent : Methyl-tert.butylether		ampoule 10 ml	10 µg/ml	Ref : F129500	Price : EF

Standard Solution - 7 components

Bromochloroacetic acid-methyl ester	CAS:20428-74-4	Dichloroacetic acid-methyl ester	CAS:116-54-1	Trichloroacetic acid-methyl ester	CAS:598-99-2
Dalapon methyl ester	CAS:17640-02-7	Monobromoacetic Acid	CAS:96-32-2		
Dibromoacetic acid-methyl ester	CAS:6482-26-4	Monochloroacetic Acid	CAS:96-34-4		
Solvent : Methyl-tert.butylether		ampoule 10 ml	10 µg/ml	Ref : F129510	Price : EF

EN 12673

Water quality - Gas chromatographic determination of some selected chlorophenols in water

This European Standard describes the gas chromatographic determination of 19 chlorophenols (2-, 3-, and 4-chlorophenol, 2,3-, 2,4-, 2,5-, 2,6-, 3,4- and 3,5-dichlorophenol, 2,3,4-, 2,3,5-, 2,3,6-, 2,4,5-, 2,4,6- and 3,4,5-trichlorophenol, 2,3,4,5-, 2,3,4,6-, and 2,3,5,6-tetrachlorophenol and pentachlorophenol) in drinking water, groundwater, rainwater, waste water, sea water and surface water.

Internal Standard Solution - 1 component			Solvent : Ethanol		
2,4-Dibromophenol	CAS:615-58-7	ampoule 1 ml	1 000 µg/ml	Ref : F131331	Price : BJ
2,6-Dibromophenol	CAS:608-33-3	ampoule 1 ml	1 000 µg/ml	Ref : F131341	Price : BH
2,3,6-Trichlorophenol	CAS:933-75-5	ampoule 1 ml	1 000 µg/ml	Ref : F131351	Price : BH
2,4,6-Tribromophenol	CAS:118-79-6	ampoule 1 ml	1 000 µg/ml	Ref : F131361	Price : BK

Phenols Standard Solution - 19 components

2-Chlorophenol	CAS:95-57-8	3000 µg/l	2,3,4,5-Tetrachlorophenol acetate	id : 14700	200 µg/l
3-Chlorophenol	CAS:108-43-0	3000 µg/l	2,3,4,6-Tetrachlorophenol acetate	id : 14 720	200 µg/l
4-Chlorophenol	CAS:106-48-9	3000 µg/l	2,3,5,6-Tetrachlorophenol acetate	id : 14730	200 µg/l
2,3-Dichlorophenol	CAS:576-24-9	400 µg/l	2,3,4-Trichlorophenol acetate	id : 15561	300 µg/l
2,4-Dichlorophenol	CAS:120-83-2	400 µg/l	2,3,5-Trichlorophenol acetate	CAS:61925-88-0	300 µg/l
2,5-Dichlorophenol	CAS:583-78-8	400 µg/l	2,3,6-Trichlorophenol acetate	CAS:61925-87-9	300 µg/l
2,6-Dichlorophenol	CAS:87-65-0	400 µg/l	2,4,5-Trichlorophenol acetate	CAS:5393-75-9	300 µg/l
3,4-Dichlorophenol	CAS:95-77-2	400 µg/l	2,4,6-Trichlorophenol acetate	CAS:23399-90-8	300 µg/l
3,5-Dichlorophenol	CAS:591-35-5	400 µg/l	3,4,5-Trichlorophenol acetate	id : 15591	300 µg/l
Pentachlorophenol	CAS:87-86-5	100 µg/l			
Solvent : Ethanol		ampoule 1 ml		Ref : F131321	Price : BCA

EN 12766

Petroleum products and used oils - Determination of PCBs and related products

The document specifies a method for the determination of the content of PCB congeners in petroleum products by gas chromatography using an electron capture detector.

PCBs Standard Solution - 1 component			Solvent : Iso-octane		
PCB 30	CAS:35693-92-6	ampoule 1 ml	10 µg/ml	Ref : F131480	Price : BG
PCB 209	CAS:2051-24-3	ampoule 1 ml	10 µg/ml	Ref : F131500	Price : BG

PCBs Standard Solution - 8 components

PCB 1	CAS:2051-60-7	PCB 47	CAS:2437-79-8	PCB 171	CAS:52663-71-5
PCB 5	CAS:16605-91-7	PCB 89	CAS:73575-57-2	PCB 200	CAS:40186-71-8
PCB 29	CAS:15862-07-4	PCB 154	CAS:60145-22-4		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F127921	Price : JJ

Quality Control Mix – 5 components

Arochlor 1242	CAS : 53469-21-9	0.50 µg/l	PCB 30	CAS : 35693-92-6	0.02 µg/l
Arochlor 1254	CAS : 11097-69-1	0.25 µg/l	PCB 209	CAS : 2051-24-3	0.02 µg/l
Arochlor 1260	CAS : 11096-82-5	0.25 µg/l			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F131520	Price : EA

Calibration Mixture without PCB 170 – 13 components

PCB 18	CAS:37680-65-2	PCB 101	CAS:37680-73-2	PCB180	CAS:35065-29-3
PCB 28	CAS:7012-37-5	PCB 118	CAS:31508-00-6	PCB194	CAS:35694-08-7
PCB 31	CAS:16606-02-3	PCB 138	CAS:35065-28-2	PCB 209	CAS:2051-24-3
PCB 44	CAS:41464-39-5	PCB 149	CAS:38380-04-0		
PCB 52	CAS:35693-99-3	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F131470	Price : JA

PCBs Calibration Mixture - 14 components

PCB 18	CAS:37680-65-2	PCB 101	CAS:37680-73-2	PCB 170	CAS:35065-30-6
PCB 28	CAS:7012-37-5	PCB 118	CAS:31508-00-6	PCB 180	CAS:35065-29-3
PCB 31	CAS:16606-02-3	PCB 138	CAS:35065-28-2	PCB 194	CAS:35694-08-7
PCB 44	CAS:41464-39-5	PCB 149	CAS:38380-04-0	PCB 209	CAS:2051-24-3
PCB 52	CAS:35693-99-3	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F127901	Price : JC

Petroleum products - Determination of aromatic hydrocarbon types in middle distillates - High performance liquid chromatography method with refractive index detection

This European Standard specifies a test method for the determination of the content of mono-aromatic, di-aromatic and tri-aromatic hydrocarbons in diesel fuels that may contain fatty acid methyl esters (FAME) up to 5 % (V/V) and petroleum distillates in the boiling range from 150 B°C to 400 B°C.

System Calibration Standard 1 - 7 components

Cyclohexane	CAS:110-82-7	1.00 % w/v	9-Methylanthracene	CAS:779-02-2	0.05 % w/v
Dibenzothiophene	CAS:132-65-0	0.05 % w/v	Naphthalene	CAS:91-20-3	0.01 % w/v
1,2-Dimethylbenzene	CAS:95-47-6	0.50 % w/v	1-Phenyldodecane	CAS:123-01-3	0.10 % w/v
Hexamethylbenzene	CAS:87-85-4	0.01 % w/v			

Solvent : n-Heptane

ampoule 5 ml

Ref : F110101

Price : EF

Water quality - Determination of parathion, parathion-methyl and some other organophosphorus compounds in water by dichloromethane extraction and gas chromatographic analysis

This European Standard specifies the extraction processes and gas chromatographic (GC) methods for determining parathion, parathion-methyl and some other organophosphorus compounds in drinking waters, surface waters and waste waters. This standard may also be suitable for the determination of other organic compounds. The range is dependent on the compound and the source of water and is typically up to 1 µg/l with a reporting limit of 0,001 µg/l for drinking waters involving a 1 000 to 1 extraction ratio.

OPP Standard Solution - 19 components

Azinphos-ethyl	CAS:2642-71-9	Dichlorvos	CAS:62-73-7	Parathion-methyl	CAS:298-00-0
Azinphos-methyl	CAS:86-50-0	Dimethoate	CAS:60-51-5	Phosalone	CAS:2310-17-0
Bromophos	CAS:2104-96-3	Fenitrothion	CAS:122-14-5	Propetamphos	CAS:31218-83-4
Chlorfenvinphos	CAS:470-90-6	Fenthion	CAS:55-38-9	Triadimefon	CAS:43121-43-3
Chlorpyrifos-ethyl	CAS:2921-88-2	Malathion	CAS:121-75-5	Triazophos	CAS:24017-47-8
Chlorpyrifos methyl	CAS:5598-13-0	Mevinphos	CAS:7786-34-7		
Diazinon	CAS:333-41-5	Parathion-ethyl	CAS:56-38-2		

Solvent : Acetone

ampoule 1 ml

100 µg/ml

Ref : F128951

Price : GF

Solvent : Methanol purge & trap

ampoule 10 ml

10 µg/ml

Ref : F131371

Price : GF

Packaging - Flexible packaging material - Determination of residual solvents by static headspace gas chromatography - Part 2: Industrial methods

This part 2 should be read in conjunction with prEN (261190) part 1 and describes rapid methods as commonly used in quality control for example for monitoring the level of residual solvents used in the production of flexible packaging by static headspace chromatography. The procedure described in this part involves one single injection of the headspace which implies an incomplete extraction of the solvent. The values obtained may be lower than the absolute content which should be determined according to Part 1. During the analysis there may be interferences from possible products of thermal decomposition.

VOC Standard Solution - 10 components

Acetic acid-isobutyl ester	CAS:110-19-0	Ethyl Acetate	CAS:141-78-6	Trichloroethene	CAS:79-01-6
Acetone	CAS:67-64-1	Methanol	CAS:67-56-1	Xylenes	CAS:1330-20-7
1-Butanol	CAS:71-36-3	Methylisobutylketone	CAS:108-10-1		
Ethanol	CAS:64-17-5	Toluene	CAS:108-88-3		

Solvent : Iso-Octane

ampoule 1 ml

200 µg/ml

Ref : F129221

Price : CJ

Characterization of waste - Determination of hydrocarbon content in the range of C10 to C40 by gas chromatography

This document specifies a method for the quantitative determination of the hydrocarbon content (C10 to C40) in solid waste by gas chromatography.

It is applicable to hydrocarbon content between 100 mg/kg and 10 000 mg/kg expressed as dry matter basis.

Using this standard all hydrocarbons with a boiling range of approximately 175 B°C to 525 B°C, e.g. n-alkanes from C10H22 to C40H82, isoalkanes, cycloalkanes, alkyl benzenes, alkyl naphthalenes and polycyclic aromatic compounds are determined as hydrocarbons, provided they do not adsorb on the Florisil column during clean-up.

n-Alkane Standard Solution - 16 components C10÷C40 (all even)

n-Decane (C10)	CAS:124-18-5	n-Docosane (C22)	CAS:629-97-0	n-Tetatriacontane (C34)	CAS:14167-59-0
n-Dodecane (C12)	CAS:112-40-3	n-Tetracosane (C24)	CAS:646-31-1	n-Hexatriacontane (C36)	CAS:630-06-8
n-Tetradecane (C14)	CAS:629-59-4	n-Hexacosane (C26)	CAS:630-01-3	n-Octatriacontane (C38)	CAS:7194-85-6
n-Hexadecane (C16)	CAS:544-76-3	n-Octacosane (C28)	CAS:630-02-4	n-Tetracontane (C40)	CAS:4181-95-7
n-Octadecane (C18)	CAS:593-45-3	n-Triacontane (C30)	CAS:638-68-6		
n-Eicosane (C20)	CAS:112-95-8	n-Dotriacontane (C32)	CAS:544-85-4		

Solvent : n-Hexane

ampoule 1 ml

50 µg/ml

Ref : F062145

Price : GA

ampoule 5 ml

50 µg/ml

Ref : F109781

Price : BGA

Standard Solution - 22 components

Bromobenzene	CAS:108-86-1	1,3-Dichlorobenzene	CAS:541-73-1	1,2,3-Trichlorobenzene	CAS:87-61-6
n-Butylbenzene	CAS:104-51-8	1,4-Dichlorobenzene	CAS:106-46-7	1,2,4-Trichlorobenzene	CAS:120-82-1
sec-Butylbenzene	CAS:135-98-8	Hexachlorobutadiene	CAS:87-68-3	1,3,5-Trichlorobenzene	CAS:108-70-3
tert-Butylbenzene	CAS:98-06-6	Isopropylbenzene	CAS:98-82-8	1,2,3-Trichloropropane	CAS:96-18-4
2-Chlorotoluene	CAS:95-49-8	4-Isopropyltoluene	CAS:99-87-6	1,2,4-Trimethylbenzene	CAS:95-63-6
4-Chlorotoluene	CAS:106-43-4	Naphthalene	CAS:91-20-3	1,3,5-Trimethylbenzene	CAS:108-67-8
1,2-Dibromo-3-chloropropane	CAS:96-12-8	n-Propylbenzene	CAS:103-65-1		
1,2-Dichlorobenzene	CAS:95-50-1	1,1,2,2-Tetrabromoethane	CAS:79-27-6		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062155	Price : DF
		ampoule 5 ml	2 000 µg/ml	Ref : F109991	Price : BDA

EN 15308 :2008**Characterization of waste - Determination of selected polychlorinated biphenyls (PCB) in solid waste by using capillary gas chromatography with electron capture or mass spectrometric detection**

This document specifies a method for quantitative determination of seven polychlorinated biphenyl congeners (PCB-28, PCB-52, PCB-101, PCB-118, PCB-138, PCB-153 and PCB-180) in solid waste using high-resolution gas chromatography with electron capture or mass spectrometric detection.

Internal and injection Standard - 1 component

Solvent : n-Hexane

PCB 29	CAS:15862-07-4	ampoule 1 ml	400 µg/ml	Ref : F131401	Price : CE
PCB 30	CAS:35693-92-6	ampoule 1 ml	400 µg/ml	Ref : F131411	Price : BG
PCB 143	CAS:68194-15-0	ampoule 1 ml	400 µg/ml	Ref : F131421	Price : FA
PCB 155	CAS:33979-03-2	ampoule 1 ml	400 µg/ml	Ref : F131431	Price : CE
PCB 198	CAS:68194-17-2	ampoule 1 ml	400 µg/ml	Ref : F131441	Price : EA
PCB 207	CAS:52663-79-3	ampoule 1 ml	400 µg/ml	Ref : F131451	Price : FA
PCB 209	CAS:2051-24-3	ampoule 1 ml	400 µg/ml	Ref : F131461	Price : BG

PCBs Standard Solution - 7 components

PCB 28	CAS:7012-37-5	PCB 118	CAS:31508-00-6	PCB 180	CAS:35065-29-3
PCB 52	CAS:35693-99-3	PCB 138	CAS:35065-28-2		
PCB 101	CAS:37680-73-2	PCB 153	CAS:35065-27-1		
Solvent : n-Pentane		ampoule 1 ml	100 µg/ml	Ref : F131391	Price : GA

EN 61619**Insulating liquids - Contamination by polychlorinated biphenyls (PCBs) - Method of determination by capillary column gas chromatography (IEC 61619:1997)**

Specifies a method for the determination of polychlorinated biphenyl (PCB) concentration in non-halogenated insulating liquids by high-resolution capillary column gas chromatography using an electron capture detector (ECD). Gives total PCB content; especially useful when a detailed analysis of PCB congeners is necessary.

Calibration PCB Solution - 14 components

PCB 18	CAS:37680-65-2	PCB 101	CAS:37680-73-2	PCB 170	CAS:35065-30-6
PCB 28	CAS:7012-37-5	PCB 118	CAS:31508-00-6	PCB 180	CAS:35065-29-3
PCB 31	CAS:16606-02-3	PCB 138	CAS:35065-28-2	PCB 194	CAS:35694-08-7
PCB 44	CAS:41464-39-5	PCB 149	CAS:38380-04-0	PCB 209	CAS:2051-24-3
PCB 52	CAS:35693-99-3	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F127901	Price : JC

ASTM D664 - 07

Standard Test Method for Acid Number of Petroleum Products by Potentiometric Titration

This test method covers procedures for the determination of acidic constituents in petroleum products and lubricants soluble or nearly soluble in mixtures of toluene and propan-2-ol. It is applicable for the determination of acids whose dissociation constants in water are larger than 10⁻⁹; extremely weak acids whose dissociation constants are smaller than 10⁻⁹ do not interfere. Salts react if their hydrolysis constants are larger than 10⁻⁹. The range of acid numbers included in the precision statement is 0.1 mg/g KOH to 150 mg/g KOH.

Solvent - 3 components		1 000 ml		Ref : F110211	Price : CDE
Isopropanol (2-propanol) 49.5%	CAS:67-63-0	Toluene	50%	Water, ion-free	0.5%
			CAS:108-88-3		CAS:7732-18-5
Titrant - Solvent : Isopropanol		1 000 ml		Ref : F110221	Price : DB
Potassium Hydroxide		0.1 M			

ASTM D2887 – 06a

Standard Test Method for Boiling Range Distribution of Petroleum Fractions by Gas Chromatography

This test method covers the determination of the boiling range distribution of petroleum products. The test method is applicable to petroleum products and fractions having a final boiling point of 538°C (1000°F) or lower at atmospheric pressure as measured by this test method. This test method is limited to samples having a boiling range greater than 55°C (100°F), and having a vapor pressure sufficiently low to permit sampling at ambient temperature.

Fuel Oil Degradation/Retention Time Mix C17-C20 – 4 components

n-Heptadecane (C17)	CAS:629-78-7	Phytane (C20)	CAS:638-36-8	Pristane (C19)	CAS:1921-70-6
n-Octadecane (C18)	CAS:593-45-3				
Solvent : Dichloromethane/Carbon disulfide (1/1)		ampoule 1 ml	2 000 µg/ml	Ref : F110261	Price : DF

Column Test – 2 components

n-Hexadecane (C16)	CAS:544-76-3	1 % w/v	n-Octadecane (C18)	CAS:593-45-3	1 % w/v
Solvent : n-Octane		ampoule 1 ml		Ref : F110051	Price : FA

Column Test Mixture - 17 components

n-Hexane (C6)	CAS:110-54-3	6 % wt	n-Octadecane (C18)	CAS:593-45-3	5 % wt
n-Heptane (C7)	CAS:142-82-5	6 % wt	n-Eicosane (C20)	CAS:112-95-8	2 % wt
n-Octane (C8)	CAS:111-65-9	8 % wt	n-Tetracosane (C24)	CAS:646-31-1	2 % wt
n-Nonane (C9)	CAS:111-84-2	8 % wt	n-Octacosane (C28)	CAS:630-02-4	1 % wt
n-Decane (C10)	CAS:124-18-5	12 % wt	n-Dotriacontane (C32)	CAS:544-85-4	1 % wt
n-Undecane (C11)	CAS:1120-21-4	12 % wt	n-Hexatriacontane (C36)	CAS:630-06-8	1 % wt
n-Dodecane (C12)	CAS:112-40-3	12 % wt	n-Tetracontane (C40)	CAS:4181-95-7	1 % wt
n-Tetradecane (C14)	CAS:629-59-4	12 % wt	n-Tetratetracontane (C44)	CAS:7098-22-8	1 % wt
n-Hexadecane (C16)	CAS:544-76-3	10 % wt			
		ampoule 1 ml		Ref : F112171	Price : EEA

Calibration Mixture C5-C44 – 18 components

n-Pentane (C5)	CAS:109-66-0	n-Undecane (C11)	CAS:1120-21-4	n-Tetracosane (C24)	CAS:646-31-1
n-Hexane (C6)	CAS:110-54-3	n-Dodecane (C12)	CAS:112-40-3	n-Octacosane (C28)	CAS:630-02-4
n-Heptane (C7)	CAS:142-82-5	n-Tetradecane (C14)	CAS:629-59-4	n-Dotriacontane (C32)	CAS:544-85-4
n-Octane (C8)	CAS:111-65-9	n-Hexadecane (C16)	CAS:544-76-3	n-Hexatriacontane (C36)	CAS:630-06-8
n-Nonane (C9)	CAS:111-84-2	n-Octadecane (C18)	CAS:593-45-3	n-Tetracontane (C40)	CAS:4181-95-7
n-Decane (C10)	CAS:124-18-5	n-Eicosane (C20)	CAS:112-95-8	n-Tetratetracontane (C44)	CAS:7098-22-8
Solvent : Carbon disulfide		ampoule 1 ml	100 µg/ml	Ref : F110231	Price : GF

Hydrocarbon Window Defining Standard C8-C40 – 35 components

n-Octane (C8)	CAS:111-65-9	Pristane (C19)	CAS:1921-70-6	n-Triacontane (C30)	CAS:638-68-6
n-Nonane (C9)	CAS:111-84-2	n-Eicosane (C20)	CAS:112-95-8	n-Hentriacontane (C31)	CAS:630-04-6
n-Decane (C10)	CAS:124-18-5	Phytane (C20)	CAS:638-36-8	n-Dotriacontane (C32)	CAS:544-85-4
n-Undecane (C11)	CAS:1120-21-4	n-Heneicosane (C21)	CAS:629-94-7	n-Tritriacontane (C33)	CAS:630-05-7
n-Dodecane (C12)	CAS:112-40-3	n-Docosane (C22)	CAS:629-97-0	n-Tettriacontane (C34)	CAS:14167-59-0
n-Tridecane (C13)	CAS:629-50-5	n-Tricosane (C23)	CAS:638-67-5	n-Pentatriacontane (C35)	CAS:630-07-9
n-Tetradecane (C14)	CAS:629-59-4	n-Tetracosane (C24)	CAS:646-31-1	n-Hexatriacontane (C36)	CAS:630-06-8
n-Pentadecane (C15)	CAS:629-62-9	n-Pentacosane (C25)	CAS:629-99-2	n-Heptatriacontane (C37)	CAS:7194-84-5
n-Hexadecane (C16)	CAS:544-76-3	n-Hexacosane (C26)	CAS:630-01-3	n-Octatriacontane (C38)	CAS:7194-85-6
n-Heptadecane (C17)	CAS:629-78-7	n-Heptacosane (C27)	CAS:593-49-7	n-Nonatriacontane (C39)	CAS:7194-86-7
n-Octadecane (C18)	CAS:593-45-3	n-Octacosane (C28)	CAS:630-02-4	n-Tetracontane (C40)	CAS:4181-95-7
n-Nonadecane (C19)	CAS:629-92-5	n-Nonacosane (C29)	CAS:630-03-5		
Solvent : Carbon disulfide/Dichloromethane (3/1)		ampoule 1 ml	500 µg/ml	Ref : F061935	Price : KF
		ampoule 5 ml	500 µg/ml	Ref : F110241	Price : DJA

ASTM D3606 - 07**Standard Test Method for Determination of Benzene and Toluene in Finished Motor and Aviation Gasoline by Gas Chromatography**

This test method covers the determination of benzene and toluene in finished motor and aviation gasolines by gas chromatography.

Aromatics Quantitative Calibration Standards (without Internal Standard)**Kit of 2 Analytes - 7 ampoules X 10 ml**

Solvent : Iso-octane [CAS 540-84-1]

Calibration range : Benzene [CAS 71-43-2]: 0.06-5.0 Toluene [CAS 108-88-3] : 0.5-20

Ampoule A	Benzene 5.00 %Vol + Toluene 20.00 %Vol + IO 75.00 %Vol	Ampoule E	Benzene 0.33 %Vol + Toluene 2.50 %Vol + IO 97.17 %Vol
Ampoule B	Benzene 2.50 %Vol + Toluene 15.00 %Vol + IO 82.50 %Vol	Ampoule F	Benzene 0.12 %Vol + Toluene 1.00 %Vol + IO 98.88 %Vol
Ampoule C	Benzene 1.25 %Vol + Toluene 10.00 %Vol + IO 88.75 %Vol	Ampoule G	Benzene 0.06 %Vol + Toluene 0.50 %Vol + IO 99.44 %Vol
Ampoule D	Benzene 0.67 %Vol + Toluene 5.00 %Vol + IO 94.33 %Vol		

Ref : F111805 Price : BCL

Aromatics Quantitative Calibration Standards (with Internal Standard)**Kit of 2 Analytes - 7 ampoules X 1 ml**

Solvent : Iso-octane [CAS 540-84-1]

Calibration range : Benzene [CAS 71-43-2]: 0.06-5.0 Toluene [CAS 108-88-3] : 0.5-20

Internal Standard : Methyl ethyl ketone : 4% Vol.

Ampoule A	Benzene 4.80 %Vol + Toluene 19.20 %Vol + IO 72.00 %Vol	Ampoule E	Benzene 0.32 %Vol + Toluene 2.40 %Vol + IO 93.28 %Vol
Ampoule B	Benzene 2.40 %Vol + Toluene 14.40 %Vol + IO 79.20 %Vol	Ampoule F	Benzene 0.12 %Vol + Toluene 0.95 %Vol + IO 94.93 %Vol
Ampoule C	Benzene 1.20 %Vol + Toluene 9.60 %Vol + IO 85.20 %Vol	Ampoule G	Benzene 0.06 %Vol + Toluene 0.48 %Vol + IO 95.46 %Vol
Ampoule D	Benzene 0.64 %Vol + Toluene 4.80 %Vol + IO 90.56 %Vol		

Ref : F111806 Price :BAJ

ASTM D3710 - 95(2004)**Standard Test Method for Boiling Range Distribution of Gasoline and Gasoline Fractions by Gas Chromatography**

This test method covers the determination of the boiling range distribution of gasoline and gasoline components. This test method is applicable to petroleum products and fractions with a final boiling point of 500°F (260°C) or lower as measured by this test method.

Hydrocarbons Calibration Standard - 17 components

n-Pentane (C5)	CAS:109-66-0	8.00 % wt	n-Pentadecane (C15)	CAS:629-62-9	2.50 % wt
n-Hexane (C6)	CAS:110-54-3	6.00 % wt	n-Butylbenzene	CAS:104-51-8	4.00 % wt
n-Heptane (C7)	CAS:142-82-5	10.00 % wt	2,4-Dimethylpentane	CAS:108-08-7	6.00 % wt
n-Octane (C8)	CAS:111-65-9	6.00 % wt	2-Methylbutane	CAS:78-78-4	10.00 % wt
n-Decane (C10)	CAS:124-18-5	4.00 % wt	2-Methylpentane	CAS:107-83-5	6.00 % wt
n-Dodecane (C12)	CAS:112-40-3	3.50 % wt	n-Propylbenzene	CAS:103-65-1	5.00 % wt
n-Tridecane (C13)	CAS:629-50-5	2.50 % wt	Toluene	CAS:108-88-3	11.00 % wt
n-Tetradecane (C14)	CAS:629-59-4	2.50 % wt	p-Xylene	CAS:106-42-3	13.00 % wt

ampoule 1 ml

Ref : F110421 Price : BJA

ASTM D3798 - 03**Standard Test Method for Analysis of p-Xylene by Gas Chromatography**

This test method covers the determination of known hydrocarbon impurities in, and the purity of p-xylene by gas chromatography (GC). It is generally meant for the analysis of p-xylene of 99 % or greater purity. Impurity concentrations that can be measured range from 0.001 to 1.000 weight %.

p-Xylene Impurity Standard - 10 components

Benzene	CAS:71-43-2	0.15 % wt	n-Propylbenzene	CAS:103-65-1	0.15 % wt
Cumene	CAS:98-82-8	0.15 % wt	Toluene	CAS:108-88-3	0.15 % wt
Ethylbenzene	CAS :100-41-4	0.15 % wt	m-Xylene	CAS:108-38-3	0.15 % wt
n-Octane (C8)	CAS:111-65-9	0.15 % wt	o-Xylene	CAS:95-47-6	0.15 % wt
n-Pentane (C5)	CAS:109-66-0	0.15 % wt	p-Xylene	CAS:106-42-3	98.65 % wt

ampoule 1 ml

Ref : F110431 Price : HF

ASTM D4815 - 04**Standard Test Method for determination of ethers and alcohols in gasolines by gas chromatography**

This test method covers the determination of ethers and alcohols in gasolines by gas chromatography. Specific compounds determined are methyl tert-butylether (MTBE), ethyl tert-butylether (ETBE), tert-amylmethylether (TAME), diisopropylether (DIPE), methanol, ethanol, isopropanol, n-propanol, isobutanol, tert-butanol, sec-butanol, n-butanol, and tert-pentanol (tert-amylalcohol)

Valve Timing Mixture - 4 components

Diisopropylether	CAS:108-20-3	Methyl-tert.butylether	CAS:1634-04-4	Methylcyclopentane	CAS:96-37-7
Ethyl-tert-butylether (ETBE)	CAS:637-92-3				

Solvent : n-Hexane (C6) 60 %wt

ampoule 1 ml

10 % wt

Ref : F111911

Price : BFA

Standard Solution - 5 components

Ethanol	CAS:64-17-5	2.85 % wt	Methyl-tert.butylether	CAS:1634-04-4	19.00 % wt
tert.-Butanol	CAS:75-65-0	0.10 % wt	tert-Pentanol	CAS:75-85-4	1.19 % wt
1,2-Dimethoxyethane	CAS:110-71-4	5.00 % wt			
Solvent : Iso-octane /Xylene (46.72/25.15)		ampoule 1 ml	Ref : F062175	Price : FA	
		ampoule 10 ml	Ref : F111851	Price : BJA	

Standard Solution - 5 components

tert.-Butanol	CAS:75-65-0	2.85 % wt	Methyl-tert.butylether	CAS:1634-04-4	14.25 % wt
Ethanol	CAS:64-17-5	0.10 % wt	tert-Pentanol	CAS:75-85-4	4.75 % wt
1,2-Dimethoxyethane	CAS:110-71-4	5.00 % wt			
Solvent : Iso-octane /Xylene (47.49/25. 57)		ampoule 1 ml	Ref : F062185	Price : GF	
		ampoule 10 ml	Ref : F111861	Price : CBA	

Standard Solution - 5 components

tert.-Butanol	CAS:75-65-0	5.70 % wt	Methyl-tert.butylether	CAS:1634-04-4	9.50 % wt
Ethanol	CAS:64-17-5	5.70 % wt	tert-Pentanol	CAS:75-85-4	2.38 % wt
1,2-Dimethoxyethane	CAS:110-71-4	5.00 % wt			
Solvent : Iso-octane /Xylene (46.62/25.11)		ampoule 1 ml	Ref : LM12085	Price : GF	
		ampoule 10 ml	Ref : F111871	Price : CBA	

Standard Solution - 5 components

tert.-Butanol	CAS:75-65-0	7.60 % wt	Methyl-tert.butylether	CAS:1634-04-4	4.75 % wt
Ethanol	CAS:64-17-5	8.55 % wt	tert-Pentanol	CAS:75-85-4	3.56 % wt
1,2-Dimethoxyethane	CAS:110-71-4	5.00 % wt			
Solvent : Iso-octane /Xylene (45.85/24.69)		ampoule 1 ml	Ref : F062205	Price : GF	
		ampoule 10 ml	Ref : F111881	Price : CCA	

Standard Solution - 5 components

tert.-Butanol	CAS:75-65-0	11.40 % wt	Methyl-tert.butylether	CAS:1634-04-4	0.10 % wt
Ethanol	CAS:64-17-5	11.40 % wt	tert-Pentanol	CAS:75-85-4	0.10 % wt
1,2-Dimethoxyethane	CAS:110-71-4	5.00 % wt			
Solvent : Iso-octane /Xylene (46.8/25.2)		ampoule 1 ml	Ref : LM12095	Price : GF	
		ampoule 10 ml	Ref : F111891	Price : CCA	

Qualitative Peak ID Mixture - 16 components

tert.-Amyl-methylether	CAS:994-05-8	7.30 % wt	Ethyl-tert-butylether	CAS:637-92-3	4.00 % wt
Benzene	CAS:71-43-2	5.00 % wt	Isobutanol	CAS:78-83-1	7.30 % wt
1-Butanol	CAS:71-36-3	7.30 % wt	Isopropanol	CAS:67-63-0	7.30 % wt
tert.-Butanol	CAS:75-65-0	7.30 % wt	Methanol	CAS:67-56-1	7.30 % wt
sec. Butylalcohol	CAS:78-92-2	4.00 % wt	Methyl-tert.butylether	CAS:1634-04-4	4.00 % wt
Diisopropylether	CAS:108-20-3	4.00 % wt	Methylcyclopentane	CAS:96-37-7	4.00 % wt
1,2-Dimethoxyethane	CAS:110-71-4	6.00 % wt	1-Propanol	CAS:71-23-8	7.30 % wt
Ethanol	CAS:64-17-5	7.30 % wt	tert-Pentanol	CAS:75-85-4	7.30 % wt
		ampoule 1 ml	Ref : F111901	Price : BEA	

Oxygenated Quantitative Calibration Mixtures**Kit of 4 Analytes – 5 ampoules**

Calibration range : Ethanol [CAS 64-17-5] : 0.10-11.40 tert-Butanol [CAS 75-65-0] : 0.10 -11.40
Methyl tert butylether [CAS 1634 04 4] : 0.10-19.00 tert -Pentanol [CAS 75-85-4] : 0.10 - 4.79

Ampoule A	Ethanol	3.00 % wt	+ tert-Butanol	0.10 % wt	+ Methyl-tert.butyl ether	20.00 %wt	+ tert-Pentanol	1.25 % wt
Ampoule B	Ethanol	0.10 % wt	+ tert-Butanol	3.00 % wt	+ Methyl-tert.butyl ether	15.00 %wt	+ tert-Pentanol	5.00 % wt
Ampoule C	Ethanol	6.00 % wt	+ tert-Butanol	6.00 % wt	+ Methyl-tert.butyl ether	10.00 %wt	+ tert-Pentanol	2.50 % wt
Ampoule D	Ethanol	9.00 % wt	+ tert-Butanol	8.00 % wt	+ Methyl-tert.butyl ether	5.00 %wt	+ tert-Pentanol	3.75 % wt
Ampoule E	Ethanol	12.00 % wt	+ tert-Butanol	12.00 % wt	+ Methyl-tert.butyl ether	0.10 %wt	+ tert-Pentanol	0.10 % wt
Solvent : Iso-Octane /Xylene (65/35)		5 ampoules X 1 ml	Ref : F062215	Price : DF				
		5 ampoules X 10 ml	Ref : F111801	Price : BCA				

Oxygenated Quantitative Calibration Mixture (with Internal Standard)**Kit of 4 Analytes – 5 ampoules**

Calibration range : Ethanol : 0.10-11.40 Butanol : 0.10-11.40 Methyl t-butylether : 0.10–19.00 Pentanol : 0.10-4.79
Internal Standard : 1,2-Dimethoxyethane: 5 % wt

Ampoule	Ethanol	Butanol	Methyl t-butylether	Pentanol	Total Oxygenates & int. standard
Ampoule A	2.85 % wt	0.10 %wt	19.00 % wt	1.19 % wt	28.14 % wt
Ampoule B	0.10 % wt	2.85 %wt	14.25 % wt	4.75 % wt	26.95 % wt
Ampoule C	5.70 % wt	5.70 %wt	9.50 % wt	2.38 % wt	28.28 % wt
Ampoule D	8.55 % wt	7.60 %wt	4.75 % wt	3.56 % wt	29.46 % wt
Ampoule E	11.40 % wt	11.40 %wt	0.10 % wt	0.10 % wt	28.00 % wt
Solvent : Iso-Octane /Xylene (65/35)		5 ampoules X 1 ml	Ref : F062225	Price : BFE	
		5 ampoules X 10 ml	Ref : F111807	Price : FDK	

ASTM D5186 – 03**Standard Test Method which measures the total aromatic and PNA content in Diesel and jet Fuel**

Standard Test Method for Determination of Aromatic Content and Polynuclear Aromatic Content of Diesel Fuels and Aviation Turbine Fuels by Supercritical Fluid Chromatography. This test method covers the determination of the total amounts of monoaromatic and polynuclear aromatic hydrocarbon compounds in motor diesel fuels, aviation turbine fuels, and blend stocks by supercritical fluid chromatography (SFC).

Performance Control Mixture – 4 components

n-Hexadecane (C16)	CAS:544-76-3	75%	Tetralin	CAS:119-64-2	3%
Naphthalene	CAS:91-20-3	2%	Toluene	CAS:108-88-3	20%
		ampoule 1 ml	Ref: F111921		Price : BFF

ASTM D5307 – 97(2007)**Standard Test Method for Determination of Boiling Range Distribution of Crude Petroleum by Gas Chromatography**

This test method covers the determination of the boiling range distribution of water-free crude petroleum through 538°C (1000°F). Material boiling above 538°C is reported as residue. This test method is applicable to whole crude samples, that can be solubilized in a solvent to permit sampling by means of a microsyringe

Column Resolution Test Mix C16-C18 – 2 components

n-Hexadecane (C16)	CAS:544-76-3	n-Octadecane (C18)	CAS:593-45-3		
Solvent : n-Octane		ampoule 1 ml	1 % wt	Ref: F111951	Price : EF

Petrochemical Internal Standard C14-C17 – 4 components

1-Tetradecene (C14)	CAS:1120-36-1	n-Hexadecane (C16)	CAS:544-76-3	n-Heptadecane (C17)	CAS:629-78-7
n-Pentadecane (C15)	CAS:629-62-9				
		bottle 40 ml	25 % wt	Ref: F111941	Price : EEA

Petrochemical Quantitative Standard C10-C44 – 16 components

n-Decane (C10)	CAS:124-18-5	n-Hexadecane (C16)	CAS:544-76-3	n-Dotriacontane (C32)	CAS:544-85-4
n-Undecane (C11)	CAS:1120-21-4	n-Heptadecane (C17)	CAS:629-78-7	n-Hexatriacontane (C36)	CAS:630-06-8
n-Dodecane (C12)	CAS:112-40-3	n-Octadecane (C18)	CAS:593-45-3	n-Tetracontane (C40)	CAS:4181-95-7
n-Tridecane (C13)	CAS:629-50-5	n-Eicosane (C20)	CAS:112-95-8	n-Tetratetracontane (C44)	CAS:7098-22-8
n-Tetradecane (C14)	CAS:629-59-4	n-Tetracosane (C24)	CAS:646-31-1		
n-Pentadecane (C15)	CAS:629-62-9	n-Octacosane (C28)	CAS:630-02-4		
Solvent : Carbon disulfide		ampoule 1 ml	5 000 µg/ml	Ref: F111931	Price : BGA

ASTM D5442 – 93(2003)e1**Standard Test Method for Analysis of Petroleum Waxes by Gas Chromatography**

This test method covers the quantitative determination of the carbon number distribution of petroleum waxes in the range from n-C17 through n-C44 by gas chromatography using internal standardization. In addition, the content of normal and non-normal hydrocarbons for each carbon number is also determined.

Column Resolution Test Mix C20-C24 - 2 components

n-Eicosane (C20)	CAS:112-95-8	n-Tetracosane (C24)	CAS:646-31-1		
Solvent : Cyclohexane		ampoule 1 ml	0.05 % wt	Ref: F112011	Price : EF

Quantitative Linearity Wax Mixture C16-C44 – 12 components

n-Hexadecane (C16)	CAS:544-76-3	n-Tetracosane (C24)	CAS:646-31-1	n-Dotriacontane (C32)	CAS:544-85-4
n-Octadecane (C18)	CAS:593-45-3	n-Hexacosane (C26)	CAS:630-01-3	n-Hexatriacontane (C36)	CAS:630-06-8
n-Eicosane (C20)	CAS:112-95-8	n-Octacosane (C28)	CAS:630-02-4	n-Tetracontane (C40)	CAS:4181-95-7
n-Docosane (C22)	CAS:629-97-0	n-Triacontane (C30)	CAS:638-68-6	n-Tetratetracontane (C44)	CAS:7098-22-8
Solvent : Cyclohexane		ampoule 1 ml	0.01 % wt	Ref: F111991	Price : FA

Retention Time Wax Mixture C16-C44 – 12 components

n-Hexadecane (C16)	CAS:544-76-3	n-Tetracosane (C24)	CAS:646-31-1	n-Dotriacontane (C32)	CAS:544-85-4
n-Octadecane (C18)	CAS:593-45-3	n-Hexacosane (C26)	CAS:630-01-3	n-Hexatriacontane (C36)	CAS:630-06-8
n-Eicosane (C20)	CAS:112-95-8	n-Octacosane (C28)	CAS:630-02-4	n-Tetracontane (C40)	CAS:4181-95-7
n-Docosane (C22)	CAS:629-97-0	n-Triacontane (C30)	CAS:638-68-6	n-Tetratetracontane (C44)	CAS:7098-22-8
Solvent : Cyclohexane		ampoule 1 ml	8.3 % wt	Ref: F112001	Price : DEA

Standard Test Method for Determination of Benzene, Toluene, Ethylbenzene, p/m-Xylene, o-Xylene, C9 and Heavier Aromatics, and Total Aromatics in Finished Gasoline by Gas Chromatography

This test method covers the determination of benzene, toluene, ethylbenzene, the xylenes, C9 and heavier aromatics, and total aromatics in finished motor gasoline by gas chromatography

Selectivity Check Standard – 1 component

n-Dodecane (C12) 1.5% wt CAS:112-40-3

Solvent : Iso-octane

ampoule 1 ml

Ref : F112141

Price : DF

Standard Solution - 2 components

2-Bromo-1-chloropropane CAS:3017-95-6

Fluorobenzene

CAS:462-06-6

Solvent : Methanol purge & trap

ampoule 1 ml

2 000 µg/ml

Ref : F112761

Price : DF

Standard Solution - 2 components

2-Bromo-1-chloropropane CAS:3017-95-6

1-Chloro-2-fluorobenzene

CAS:348-51-6

Solvent : Methanol purge & trap

ampoule 1 ml

2 000 µg/ml

Ref : F112771

Price : DF

Standard Solution - 6 components

Bromomethane

CAS:74-83-9

Chloromethane

CAS:74-87-3

Fluorotrichloromethane

CAS:75-69-4

Chloroethane

CAS:75-00-3

Dichlorodifluoromethane

CAS:75-71-8

Vinylchloride

CAS:75-01-4

Solvent : Methanol purge & trap

ampoule 1 ml

200 µg/ml

Ref : F112901

Price : FA

ampoule 1 ml

2000 µg/ml

Ref : F112911

Price : HF

Valve Timing Calibration Mix (without Internal Standard) – 4 components

Benzene

CAS:71-43-2

5.00 % wt

Toluene

CAS:108-88-3

5.00 % wt

Ethylbenzene

CAS:100-41-4

10.00 % wt

o-Xylene

CAS:95-47-6

10.00 % wt

Solvent : Iso-octane 70% wt

ampoule 1 ml

Ref : F062235

Price : EF

ampoule 10 ml

Ref : F112121

Price : BEA

Valve Timing Calibration Mix (with Internal Standard) – 5 components

Benzene

CAS:71-43-2

4.50 % wt

Toluene

CAS:108-88-3

4.50 % wt

Ethylbenzene

CAS:100-41-4

9.00 % wt

o-Xylene

CAS:95-47-6

9.00 % wt

2-Hexanone

CAS:591-78-6

10.00 %wt

Solvent : Iso-octane 63 %wt

ampoule 1 ml

Ref : F062245

Price : EF

ampoule 10 ml

Ref : F112131

Price : BEA

Aromatics Quantitative Calibration Mixes (without Internal Standard)**5 Analytes – 5 ampoules**

Calibration range : Benzene [CAS 71-43-2] : 0.10-5.0 Toluene [CAS 108-88-3] : 1.00-15.00 o-Xylene [CAS 95-47-6] : 0.50-10.00
Ethylbenzene [CAS 100-41-4] : 0.50-10.00 1,2,4- Trimethylbenzene [CAS 95-63-6] : 0.50-10.00

Ampoule A	Benzene	0.10 % wt	+ Toluene	15.00 % wt	+ Ethyl Benzene	0.50 % wt	+ o-Xylene	1.00 % wt	+ 1,2,4-Trimethylbenzene	1.00 % wt
Ampoule B	Benzene	0.50 % wt	+ Toluene	10.00 % wt	+ Ethyl Benzene	1.00 % wt	+ o-Xylene	2.50 % wt	+ 1,2,4-Trimethylbenzene	10.00 % wt
Ampoule C	Benzene	1.00 % wt	+ Toluene	5.00 % wt	+ Ethyl Benzene	2.50 % wt	+ o-Xylene	10.00 % wt	+ 1,2,4-Trimethylbenzene	0.50 % wt
Ampoule D	Benzene	2.00 % wt	+ Toluene	2.50 % wt	+ Ethyl Benzene	5.00 % wt	+ o-Xylene	5.00 % wt	+ 1,2,4-Trimethylbenzene	5.00 % wt
Ampoule E	Benzene	5.00 % wt	+ Toluene	1.00 % wt	+ Ethyl Benzene	10.00 % wt	+ o-Xylene	0.50 % wt	+ 1,2,4-Trimethylbenzene	2.50 % wt

Solvent : Iso-Octane [CAS 540-84-1]

5 ampoules X 1 ml

Ref : F062255

Price : JA

5 ampoules X 10 ml

Ref : F112021

Price : DCA

Aromatics Quantitative Calibration Mixes (with Internal Standard)**5 Analytes + 2-hexanone as internal standard – 5 ampoules**

Internal standard : 2-hexanone [CAS 591-78-6] 10.00 % wt

Calibration range : Benzene [CAS 71-43-2] : 0.10-5.0 Toluene [CAS 108-88-3] : 1.00-15.00 o-Xylene [CAS 95-47-6] : 0.50-10.00
Ethylbenzene [CAS 100-41-4] : 0.50-10.00 1,2,4- Trimethylbenzene [CAS 95-63-6] : 0.50-10.00

Ampoule A	Benzene	0.09 %wt	+ Toluene	13.50 %wt	+ Ethyl Benzene	0.45 %wt	+ o-Xylene	0.90 %wt	+ 1,2,4-Trimethylbenzene	9.00 % wt
Ampoule B	Benzene	0.45 %wt	+ Toluene	9.00 %wt	+ Ethyl Benzene	0.90 %wt	+ o-Xylene	2.25 %wt	+ 1,2,4-Trimethylbenzene	0.45 % wt
Ampoule C	Benzene	0.90 %wt	+ Toluene	4.50 %wt	+ Ethyl Benzene	2.25 %wt	+ o-Xylene	9.00 %wt	+ 1,2,4-Trimethylbenzene	4.50 % wt
Ampoule D	Benzene	1.80 %wt	+ Toluene	2.25 %wt	+ Ethyl Benzene	4.50 %wt	+ o-Xylene	4.50 %wt	+ 1,2,4-Trimethylbenzene	2.25 % wt
Ampoule E	Benzene	4.50 %wt	+ Toluene	0.90 %wt	+ Ethyl Benzene	9.00 %wt	+ o-Xylene	0.45 %wt	+ 1,2,4-Trimethylbenzene	0.90 % wt

Solvent : Iso-Octane [CAS 540-84-1]

5 ampoules X 1 ml

Ref : F112071

Price : BCA

Daily Quality Control Standard (without Internal Standard) - 14 components

n-Hexane (C6)	CAS:110-54-3	12 %wt	Iso-octane	CAS:540-84-1	20 %wt
n-Heptane (C7)	CAS:142-82-5	20 %wt	Naphthalene	CAS:91-20-3	1 %wt
n-Octane (C8)	CAS:111-65-9	14 %wt	Toluene	CAS:108-88-3	9 %wt
n-Decane (C10)	CAS:124-18-5	10 %wt	1,2,4,5-Tetramethylbenzene	CAS:95-93-2	1 %wt
n-Dodecane (C12)	CAS:112-40-3	1 %wt	1,2,4-Trimethylbenzene	CAS:95-63-6	3 %wt
Benzene	CAS:71-43-2	1 %wt	o-Xylene	CAS:95-47-6	2 %wt
Ethylbenzene	CAS:100-41-4	2 %wt	p-Xylene	CAS:106-42-3	3 %wt
ampoule 1 ml			Ref : F062265	Price : JF	
ampoule 10 ml			Ref : F112151	Price : DEA	

Daily Quality Control Standard (with Internal Standard) – 15 components

n-Hexane (C6)	CAS:110-54-3	11 %wt	Naphthalene	CAS:91-20-3	1 %wt
n-Heptane (C7)	CAS:142-82-5	19 %wt	Toluene	CAS:108-88-3	8 %wt
n-Octane (C8)	CAS:111-65-9	14 %wt	1,2,4,5-Tetramethylbenzene	CAS:95-93-2	1 %wt
n-Decane (C10)	CAS:124-18-5	9 %wt	1,2,4-Trimethylbenzene	CAS:95-63-6	2 %wt
n-Tridecane (C13)	CAS:629-50-5	1 %wt	o-Xylene	CAS:95-47-6	1 %wt
Benzene	CAS:71-43-2	1 %wt	p-Xylene	CAS:106-42-3	2 %wt
Ethylbenzene	CAS:100-41-4	1 %wt	2-Hexanone	CAS:591-78-6	10 %wt
Iso-octane	CAS:540-84-1	19 %wt			
ampoule 1 ml			Ref : F062275	Price : JF	
ampoule 10 ml			Ref : F112161	Price : DEA	

Standard Solution - 54 components

Benzene	CAS:71-43-2	1,3-Dichlorobenzene	CAS:541-73-1	Naphthalene	CAS:91-20-3
Bromobenzene	CAS:108-86-1	1,4-Dichlorobenzene	CAS:106-46-7	n-Propylbenzene	CAS:103-65-1
Bromochloromethane	CAS:74-97-5	1,1-Dichloroethane	CAS:75-34-3	Styrene	CAS:100-42-5
Bromodichloromethane	CAS:75-27-4	1,2-Dichloroethane	CAS:107-06-2	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromoform	CAS:75-25-2	1,1-Dichloroethene	CAS:75-35-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
n-Butylbenzene	CAS:104-51-8	cis-1,2-Dichloroethene	CAS:156-59-2	Tetrachloroethene	CAS:127-18-4
sec-Butylbenzene	CAS:135-98-8	trans-1,2-Dichloroethene	CAS:156-60-5	Toluene	CAS:108-88-3
tert-Butylbenzene	CAS:98-06-6	1,2-Dichloropropane	CAS:78-87-5	1,2,3-Trichlorobenzene	CAS:87-61-6
Carbon tetrachloride	CAS:56-23-5	1,3-Dichloropropane	CAS:142-28-9	1,2,4-Trichlorobenzene	CAS:120-82-1
Chlorobenzene	CAS:108-90-7	2,2-Dichloropropane	CAS:594-20-7	1,1,1-Trichloroethane	CAS:71-55-6
Chloroform	CAS:67-66-3	1,1-Dichloropropene	CAS:563-58-6	1,1,2-Trichloroethane	CAS:79-00-5
2-Chlorotoluene	CAS:95-49-8	cis-1,3-Dichloropropene	CAS:10061-01-5	Trichloroethene	CAS:79-01-6
4-Chlorotoluene	CAS:106-43-4	trans-1,3-Dichloropropene	CAS:10061-02-6	1,2,3-Trichloropropane	CAS:96-18-4
Dibromochloromethane	CAS:124-48-1	Ethylbenzene	CAS:100-41-4	1,2,4-Trimethylbenzene	CAS:95-63-6
1,2-Dibromo-3-chloropropane	CAS:96-12-8	Hexachlorobutadiene	CAS:87-68-3	1,3,5-Trimethylbenzene	CAS:108-67-8
1,2-Dibromoethane	CAS:106-93-4	Isopropylbenzene	CAS:98-82-8	o-Xylene	CAS:95-47-6
Dibromomethane	CAS:74-95-3	4-Isopropyltoluene	CAS:99-87-6	m-Xylene	CAS:108-38-3
1,2-Dichlorobenzene	CAS:95-50-1	Methylene chloride	CAS:75-09-2	p-Xylene	CAS:106-42-3
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112781	Price : BDA

ASTM E1387 - 01**Standard Test Method for Ignitable Liquid Residues in Extracts from Fire Debris Samples by Gas Chromatography**

This test method covers the identification of residues of ignitable liquids in extracts from fire debris samples. These materials also can be used for underground storage tank monitoring.

Column Resolution Check Mix - 13 components

n-Hexane (C6)	CAS:110-54-3	n-Hexadecane (C16)	CAS:544-76-3	Toluene	CAS:108-88-3
n-Octane (C8)	CAS:111-65-9	n-Octadecane (C18)	CAS:593-45-3	1,2,4-Trimethylbenzene	CAS:95-63-6
n-Decane (C10)	CAS:124-18-5	n-Eicosane (C20)	CAS:112-95-8	p-Xylene	CAS:106-42-3
n-Dodecane (C12)	CAS:112-40-3	2-Ethyltoluene	CAS:611-14-3		
n-Tetradecane (C14)	CAS:629-59-4	3-Ethyltoluene	CAS:620-14-4		
Solvent : Dichloromethane		ampoule 1 ml	2 000 µg/ml	Ref : F128941	Price : FA

EPA methods 500 for Drinking Water

The 500 Series product line contains standards used in proposed and promulgated methods for the identification and quantification of organic compounds in drinking water. The organic compounds listed in the various methods include volatile organic compounds (VOCs), pesticides, synthetic organic compounds (SOCs), and trihalomethane disinfection by-products. Analytical techniques used in the identification and quantification of the above compounds include gas chromatography with selective detectors (PID, ELCD, ECD, FID, NPD, FPD), gas chromatography/mass spectrometry and high performance liquid chromatography.

EPA method 501

Trihalomethane Analysis

Method 501 is a purge and trap method for measurement of total trihalomethanes in raw source water, or drinking water in any stage of treatment, and in finished drinking water by gas chromatography.

Trihalomethanes Standard Solution - 4 components

Bromodichloromethane	CAS:75-27-4	Chloroform	CAS:67-66-3	Dibromochloromethane	CAS:124-48-1
Bromoform	CAS:75-25-2				
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F112711	Price : CF
		ampoule 1 ml	200 µg/ml	Ref : F112721	Price : CF
		ampoule 1 ml	2 000 µg/ml	Ref : F112741	Price : CF
		ampoule 1 ml	5 000 µg/ml	Ref : F112731	Price : DA

EPA method 502

Volatile Organic Compounds in Water

EPA Method 502.1 identifies volatile halogenated organic compounds, in water, by purge and trap gas chromatography.

EPA Method 502.2 is an enhanced and expanded version 502.1. It is a general purpose method for the identification and measurement of purgeable volatile organic compounds in finished drinking water, raw source water, or drinking water in any treatment stage.

Aromatic Halocarbons Mixture - 9 components

Bromobenzene	CAS:108-86-1	4-Chlorotoluene	CAS:106-43-4	1,4-Dichlorobenzene	CAS:106-46-7
Chlorobenzene	CAS:108-90-7	1,2-Dichlorobenzene	CAS:95-50-1	1,2,3-Trichlorobenzene	CAS:87-61-6
2-Chlorotoluene	CAS:95-49-8	1,3-Dichlorobenzene	CAS:541-73-1	1,2,4-Trichlorobenzene	CAS:120-82-1
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113011	Price : DF

Aromatic Hydrocarbons Mixture - 16 components

Benzene	CAS:71-43-2	4-Isopropyltoluene	CAS:99-87-6	1,3,5-Trimethylbenzene	CAS:108-67-8
n-Butylbenzene	CAS:104-51-8	Naphthalene	CAS:91-20-3	m-Xylene	CAS:108-38-3
sec-Butylbenzene	CAS:135-98-8	n-Propylbenzene	CAS:103-65-1	o-Xylene	CAS:95-47-6
tert-Butylbenzene	CAS:98-06-6	Styrene	CAS:100-42-5	p-Xylene	CAS:106-42-3
Ethylbenzene	CAS:100-41-4	Toluene	CAS:108-88-3		
Isopropylbenzene	CAS:98-82-8	1,2,4-Trimethylbenzene	CAS:95-63-6		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113001	Price : FA

Aromatics Volatile Organics Mixture - 25 components

Benzene	CAS:71-43-2	1,3-Dichlorobenzene	CAS:541-73-1	1,2,3-Trichlorobenzene	CAS:87-61-6
Bromobenzene	CAS:108-86-1	1,4-Dichlorobenzene	CAS:106-46-7	1,2,4-Trichlorobenzene	CAS:120-82-1
n-Butylbenzene	CAS:104-51-8	Ethylbenzene	CAS:100-41-4	1,2,4-Trimethylbenzene	CAS:95-63-6
sec-Butylbenzene	CAS:135-98-8	Isopropylbenzene	CAS:98-82-8	1,3,5-Trimethylbenzene	CAS:108-67-8
tert-Butylbenzene	CAS:98-06-6	4-Isopropyltoluene	CAS:99-87-6	o-Xylene	CAS:95-47-6
Chlorobenzene	CAS:108-90-7	Naphthalene	CAS:91-20-3	m-Xylene	CAS:108-38-3
2-Chlorotoluene	CAS:95-49-8	n-Propylbenzene	CAS:103-65-1	p-Xylene	CAS:106-42-3
4-Chlorotoluene	CAS:106-43-4	Styrene	CAS:100-42-5		
1,2-Dichlorobenzene	CAS:95-50-1	Toluene	CAS:108-88-3		
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F113021	Price : EF
		ampoule 1 ml	200 µg/ml	Ref : F113031	Price : EF
		ampoule 1 ml	2 000 µg/ml	Ref : F113051	Price : HA

Haloalkanes Mixture - 35 components

Bromochloromethane	CAS:74-97-5	Dichlorodifluoromethane	CAS:75-71-8	Hexachlorobutadiene	CAS:87-68-3
Bromodichloromethane	CAS:75-27-4	1,1-Dichloroethane	CAS:75-34-3	Methylene chloride	CAS:75-09-2
Bromoform	CAS:75-25-2	1,2-Dichloroethane	CAS:107-06-2	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromomethane	CAS:74-83-9	1,1-Dichloroethene	CAS:75-35-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Carbon tetrachloride	CAS:56-23-5	cis-1,2-Dichloroethene	CAS:156-59-2	Tetrachloroethene	CAS:127-18-4
Chloroethane	CAS:75-00-3	trans-1,2-Dichloroethene	CAS:156-60-5	1,1,1-Trichloroethane	CAS:71-55-6
Chloroform	CAS:67-66-3	1,2-Dichloropropane	CAS:78-87-5	1,1,2-Trichloroethane	CAS:79-00-5
Chloromethane	CAS:74-87-3	1,3-Dichloropropane	CAS:142-28-9	Trichloroethene	CAS:79-01-6
Dibromochloromethane	CAS:124-48-1	2,2-Dichloropropane	CAS:594-20-7	Trichlorofluoromethane	CAS:75-69-4
1,2-Dibromo-3-chloropropane	CAS:96-12-8	1,1-Dichloropropene	CAS:563-58-6	1,2,3-Trichloropropane	CAS:96-18-4
1,2-Dibromoethane	CAS:106-93-4	cis-1,3-Dichloropropene	CAS:10061-01-5	Vinyl chloride	CAS:75-01-4
Dibromomethane	CAS:74-95-3	trans-1,3-Dichloropropene	CAS:10061-02-6		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112991	Price : KF

Haloethanes Mixture - 14 components

Chloroethane	CAS:75-00-3	cis-1,2-Dichloroethene	CAS:156-59-2	1,1,1-Trichloroethane	CAS:71-55-6
1,2-Dibromoethane	CAS:106-93-4	trans-1,2-Dichloroethene	CAS:156-60-5	1,1,2-Trichloroethane	CAS:79-00-5
1,1-Dichloroethane	CAS:75-34-3	1,1,1,2-Tetrachloroethane	CAS:630-20-6	Trichloroethene	CAS:79-01-6
1,2-Dichloroethane	CAS:107-06-2	1,1,2,2-Tetrachloroethane	CAS:79-34-5	Vinylchloride	CAS:75-01-4
1,1-Dichloroethene	CAS:75-35-4	Tetrachloroethene	CAS:127-18-4		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112971	Price : FA

Halomethanes Mixture - 12 components

Bromochloromethane	CAS:74-97-5	Carbon tetrachloride	CAS:56-23-5	Dibromomethane	CAS:74-95-3
Bromodichloromethane	CAS:75-27-4	Chloroform	CAS:67-66-3	Dichlorodifluoromethane	CAS:75-71-8
Bromoform	CAS:75-25-2	Chloromethane	CAS:74-87-3	Methylene chloride	CAS:75-09-2
Bromomethane	CAS:74-83-9	Dibromochloromethane	CAS:124-48-1	Trichlorofluoromethane	CAS:75-69-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112961	Price : FA

Halopropanes Mixture - 9 components

1,2-Dibromo-3-chloropropane	CAS:96-12-8	2,2-Dichloropropane	CAS:594-20-7	trans-1,3-Dichloropropene	CAS:10061-02-6
1,2-Dichloropropane	CAS:78-87-5	1,1-Dichloropropene	CAS:563-58-6	Hexachlorobutadiene	CAS:87-68-3
1,3-Dichloropropane	CAS:142-28-9	cis-1,3-Dichloropropene	CAS:10061-01-5	1,2,3-Trichloropropane	CAS:96-18-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112981	Price : FA

VOC Standard Solution - 1 component

2-Bromo-1-chloropropane	CAS:3017-95-6	ampoule 1 ml	2 000 µg/ml	Ref : F112931	Price : DA
1-Chloro-2-fluorobenzene	CAS:348-51-6	ampoule 1 ml	2 000 µg/ml	Ref : F112951	Price : DA
1,4-Dichlorobutane	CAS:110-56-5	ampoule 1 ml	2 000 µg/ml	Ref : F112941	Price : DA
Fluorobenzene	CAS:462-06-6	ampoule 1 ml	2 000 µg/ml	Ref : F112921	Price : DA

VOC Gas Mixture - 6 components

Bromomethane	CAS:74-83-9	Chloromethane	CAS:74-87-3	Fluorotrichloromethane	CAS:75-69-4
Chloroethane	CAS:75-00-3	Dichlorodifluoromethane	CAS:75-71-8	Vinylchloride	CAS:75-01-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112901	Price : FA
		ampoule 1 ml	2 000 µg/ml	Ref : F112911	Price : HF

VOC Internal Standard Mixture - 2 components

2-Bromo-1-chloropropane	CAS:3017-95-6	Fluorobenzene	CAS:462-06-6		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112761	Price : DF

VOC Internal Standard Mixture - 2 components

2-Bromo-1-chloropropane	CAS:3017-95-6	1-Chloro-2-fluorobenzene	CAS:348-51-6		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112771	Price : DF

VOC Standard Solution - 54 components

Benzene	CAS:71-43-2	1,3-Dichlorobenzene	CAS:541-73-1	Naphthalene	CAS:91-20-3
Bromobenzene	CAS:108-86-1	1,4-Dichlorobenzene	CAS:106-46-7	n-Propylbenzene	CAS:103-65-1
Bromochloromethane	CAS:74-97-5	1,1-Dichloroethane	CAS:75-34-3	Styrene	CAS:100-42-5
Bromodichloromethane	CAS:75-27-4	1,2-Dichloroethane	CAS:107-06-2	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromoform	CAS:75-25-2	1,1-Dichloroethene	CAS:75-35-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
n-Butylbenzene	CAS:104-51-8	cis-1,2-Dichloroethene	CAS:156-59-2	Tetrachloroethene	CAS:127-18-4
sec-Butylbenzene	CAS:135-98-8	trans-1,2-Dichloroethene	CAS:156-60-5	Toluene	CAS:108-88-3
tert-Butylbenzene	CAS:98-06-6	1,2-Dichloropropane	CAS:78-87-5	1,2,3-Trichlorobenzene	CAS:87-61-6
Carbon tetrachloride	CAS:56-23-5	1,3-Dichloropropane	CAS:142-28-9	1,2,4-Trichlorobenzene	CAS:120-82-1
Chlorobenzene	CAS:108-90-7	2,2-Dichloropropane	CAS:594-20-7	1,1,1-Trichloroethane	CAS:71-55-6
Chloroform	CAS:67-66-3	1,1-Dichloropropene	CAS:563-58-6	1,1,2-Trichloroethane	CAS:79-00-5
2-Chlorotoluene	CAS:95-49-8	cis-1,3-Dichloropropene	CAS:10061-01-5	Trichloroethene	CAS:79-01-6
4-Chlorotoluene	CAS:106-43-4	trans-1,3-Dichloropropene	CAS:10061-02-6	1,2,3-Trichloropropane	CAS:96-18-4
Dibromochloromethane	CAS:124-48-1	Ethylbenzene	CAS:100-41-4	1,2,4-Trimethylbenzene	CAS:95-63-6
1,2-Dibromo-3-chloropropane	CAS:96-12-8	Hexachlorobutadiene	CAS:87-68-3	1,3,5-Trimethylbenzene	CAS:108-67-8
1,2-Dibromoethane	CAS:106-93-4	Isopropylbenzene	CAS:98-82-8	o-Xylene	CAS:95-47-6
Dibromomethane	CAS:74-95-3	4-Isopropyltoluene	CAS:99-87-6	m-Xylene	CAS:108-38-3
1,2-Dichlorobenzene	CAS:95-50-1	Methylene chloride	CAS:75-09-2	p-Xylene	CAS:106-42-3
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112781	Price : BDA

VOC Standard Solution - 60 components

Benzene	CAS:71-43-2	1,2-Dichlorobenzene	CAS:95-50-1	Naphthalene	CAS:91-20-3
Bromobenzene	CAS:108-86-1	1,3-Dichlorobenzene	CAS:541-73-1	n-Propylbenzene	CAS:103-65-1
Bromochloromethane	CAS:74-97-5	1,4-Dichlorobenzene	CAS:106-46-7	Styrene	CAS:100-42-5
Bromodichloromethane	CAS:75-27-4	Dichlorodifluoromethane	CAS:75-71-8	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromoform	CAS:75-25-2	1,1-Dichloroethane	CAS:75-34-3	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Bromomethane	CAS:74-83-9	1,2-Dichloroethane	CAS:107-06-2	Tetrachloroethene	CAS:127-18-4
n-Butylbenzene	CAS:104-51-8	1,1-Dichloroethene	CAS:75-35-4	Toluene	CAS:108-88-3
sec-Butylbenzene	CAS:135-98-8	cis-1,2-Dichloroethene	CAS:156-59-2	1,2,3-Trichlorobenzene	CAS:87-61-6
tert-Butylbenzene	CAS:98-06-6	trans-1,2-Dichloroethene	CAS:156-60-5	1,2,4-Trichlorobenzene	CAS:120-82-1
Carbon tetrachloride	CAS:56-23-5	1,2-Dichloropropane	CAS:78-87-5	1,1,1-Trichloroethane	CAS:71-55-6
Chlorobenzene	CAS:108-90-7	1,3-Dichloropropane	CAS:142-28-9	1,1,2-Trichloroethane	CAS:79-00-5
Chloroethane	CAS:75-00-3	2,2-Dichloropropane	CAS:594-20-7	Trichloroethene	CAS:79-01-6
Chloroform	CAS:67-66-3	1,1-Dichloropropene	CAS:563-58-6	Trichlorofluoromethane	CAS:75-69-4
Chloromethane	CAS:74-87-3	cis-1,3-Dichloropropene	CAS:10061-01-5	1,2,3-Trichloropropane	CAS:96-18-4
2-Chlorotoluene	CAS:95-49-8	trans-1,3-Dichloropropene	CAS:10061-02-6	1,2,4-Trimethylbenzene	CAS:95-63-6
4-Chlorotoluene	CAS:106-43-4	Ethylbenzene	CAS:100-41-4	1,3,5-Trimethylbenzene	CAS:108-67-8
Dibromochloromethane	CAS:124-48-1	Hexachlorobutadiene	CAS:87-68-3	Vinylchloride	CAS:75-01-4
1,2-Dibromo-3-chloropropane	CAS:96-12-8	Isopropylbenzene	CAS:98-82-8	o-Xylene	CAS:95-47-6
1,2-Dibromoethane	CAS:106-93-4	4-Isopropyltoluene	CAS:99-87-6	m-Xylene	CAS:108-38-3
Dibromomethane	CAS:74-95-3	Methylene chloride	CAS:75-09-2	p-Xylene	CAS:106-42-3
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112751	Price : BGA

EPA method 503.1

Method 503.1 is applicable to the determination of various volatile aromatic and unsaturated organic compounds by purge and trap gas chromatography in finished drinking water, any raw source water, or drinking water in any treatment stage.

ACs Standard Solution - 1 component

α, α, α-Trifluorotoluene	CAS:98-08-8				
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113071	Price : CF
		ampoule 1 ml	2 000 µg/ml	Ref : F113061	Price : CF

Aromatics & Alkenes Mixture - 28 components

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	1,2,3-Trichlorobenzene	CAS:87-61-6
Bromobenzene	CAS:108-86-1	Ethylbenzene	CAS:100-41-4	1,2,4-Trichlorobenzene	CAS:120-82-1
n-Butylbenzene	CAS:104-51-8	Hexachlorobutadiene	CAS:87-68-3	Trichloroethene	CAS:79-01-6
sec-Butylbenzene	CAS:135-98-8	Isopropylbenzene	CAS:98-82-8	1,2,4-Trimethylbenzene	CAS:95-63-6
tert-Butylbenzene	CAS:98-06-6	4-Isopropyltoluene	CAS:99-87-6	1,3,5-Trimethylbenzene	CAS:108-67-8
Chlorobenzene	CAS:108-90-7	Naphthalene	CAS:91-20-3	o-Xylene	CAS:95-47-6
2-Chlorotoluene	CAS:95-49-8	n-Propylbenzene	CAS:103-65-1	m-Xylene	CAS:108-38-3
4-Chlorotoluene	CAS:106-43-4	Styrene	CAS:100-42-5	p-Xylene	CAS:106-42-3
1,2-Dichlorobenzene	CAS:95-50-1	Tetrachloroethene	CAS:127-18-4		
1,3-Dichlorobenzene	CAS:541-73-1	Toluene	CAS:108-88-3		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113041	Price : DF
		ampoule 1 ml	2 000 µg/ml	Ref : F113081	Price : HA

Aromatic Hydrocarbons Mixture - 16 components

Benzene	CAS:71-43-2	4-Isopropyltoluene	CAS:99-87-6	1,3,5-Trimethylbenzene	CAS:108-67-8
n-Butylbenzene	CAS:104-51-8	Naphthalene	CAS:91-20-3	m-Xylene	CAS:108-38-3
sec-Butylbenzene	CAS:135-98-8	n-Propylbenzene	CAS:103-65-1	o-Xylene	CAS:95-47-6
tert-Butylbenzene	CAS:98-06-6	Styrene	CAS:100-42-5	p-Xylene	CAS:106-42-3
Ethylbenzene	CAS:100-41-4	Toluene	CAS:108-88-3		
Isopropylbenzene	CAS:98-82-8	1,2,4-Trimethylbenzene	CAS:95-63-6		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113001	Price : FA

Halocarbons Mixture - 12 components

Bromobenzene	CAS:108-86-1	1,2-Dichlorobenzene	CAS:95-50-1	Tetrachloroethene	CAS:127-18-4
Chlorobenzene	CAS:108-90-7	1,3-Dichlorobenzene	CAS:541-73-1	1,2,3-Trichlorobenzene	CAS:87-61-6
2-Chlorotoluene	CAS:95-49-8	1,4-Dichlorobenzene	CAS:106-46-7	1,2,4-Trichlorobenzene	CAS:120-82-1
4-Chlorotoluene	CAS:106-43-4	Hexachlorobutadiene	CAS:87-68-3	Trichloroethene	CAS:79-01-6
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113101	Price : EF

EPA method 504. 504.1

Method 504.1 is applicable to the determination of 1,2-Dibromoethane (EDB) and 1,2-Dibromo-3-chloropropane (DBCP) in finished drinking water, raw source water, or drinking water in any treatment stage by microextraction and gas chromatography.

Method 504.1 adds 1,2,3-trichloropropane to the analyte list.

EDB/DBCP Mixture - 2 components

1,2-Dibromo-3-chloropropane	CAS:96-12-8	1,2-Dibromoethane	CAS:106-93-4		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113111	Price : DA
		ampoule 1 ml	2 000 µg/ml	Ref : F113121	Price : DA

EPA method 504.1 Mixtures**Standard Solution - 3 components**

1,2-Dibromo-3-chloropropane	CAS:96-12-8	1,2-Dibromoethane	CAS:106-93-4	1,2,3-Trichloropropane	CAS:96-18-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113131	Price : DF
		ampoule 1 ml	2 000 µg/ml	Ref : F113141	Price : DF

EPA method 505

Method 505 is applicable to the determination of organohalide pesticides and commercial polychlorinated biphenyl (PCB) products in finished drinking water, during intermediate stages of treatment, and raw source water, by microextraction and gas chromatography.

NCC Standard Solution - 8 components

Atrazine	CAS:1912-24-9	Heptachlor	CAS:76-44-8	trans-Nonachlor	CAS:39765-80-5
cis-Chlordane	CAS:5103-71-9	Hexachlorocyclopentadiene	CAS:77-47-4	Simazine (CAT)	CAS:122-34-9
trans-Chlordane	CAS:5103-74-2	cis-Nonachlor	CAS:5103-73-1		
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F113171	Price : FF

NCC Standard Solution - 16 components

Alachlor	CAS:15972-60-8	10 µg/ml	Heptachlor	CAS:76-44-8	1 µg/ml
Aldrin	CAS:309-00-2	1 µg/ml	Hexachlorobenzene	CAS:118-74-1	1 µg/ml
Atrazine	CAS:1912-24-9	250 µg/ml	Hexachlorocyclopentadiene	CAS:77-47-4	1 µg/ml
alpha-Chlordane	CAS:5103-71-9	1 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	5 µg/ml
Gamma-Chlordane	CAS:5103-74-2	1 µg/ml	cis-Nonachlor	CAS:5103-73-1	1 µg/ml
Dieldrin	CAS:60-57-1	1 µg/ml	trans-Nonachlor	CAS:39765-80-5	1 µg/ml
Endrin	CAS:72-20-8	1 µg/ml	Simazine (CAT)	CAS:122-34-9	250 µg/ml
Gamma-HCH (Lindane)	CAS:58-89-9	1 µg/ml	cis-Heptachlorepoxyde (cis-, exo-, isomer B)	CAS:1024-57-3	1 µg/ml
Solvent : Acetone		ampoule 1 ml		Ref : F113181	Price : GF

OCs Standard Solution - 1 component

Toxaphene	CAS:8001-35-2				
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F120031	Price : BG

OCs Instrument Performance Check Solution - 1 component

			Solvent : Methanol		
Endrin	CAS:72-20-8	ampoule 1 ml	100 µg/ml	Ref : F113291	Price : BG
Endrin ketone	CAS:53494-70-5	ampoule 1 ml	100 µg/ml	Ref : F113301	Price : CE

OCs / NCC Standard Solution - 12 components

Alachlor	CAS:15972-60-8	50 µg/ml	Gamma-HCH (Lindane)	CAS:58-89-9	20 µg/ml
Aldrin	CAS:309-00-2	20 µg/ml	Heptachlor	CAS:76-44-8	20 µg/ml
Atrazine	CAS:1912-24-9	500 µg/ml	Hexachlorobenzene	CAS:118-74-1	10 µg/ml
Dieldrin	CAS:60-57-1	20 µg/ml	Hexachlorocyclopentadiene	CAS:77-47-4	20 µg/ml
Endrin	CAS:72-20-8	20 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	200 µg/ml
cis-Heptachlorepoxyde (cis-, exo-, isomer B)	CAS:1024-57-3	20 µg/ml	Simazine (CAT)	CAS:122-34-9	100 µg/ml
Solvent : Acetone		ampoule 1 ml	200 µg/ml	Ref : F113161	Price : EF

Organochlorine Pesticides Mixture - 9 components

Alachlor	CAS:15972-60-8	Endrin	CAS:72-20-8	Hexachlorobenzene	CAS:118-74-1
Aldrin	CAS:309-00-2	Heptachlor	CAS:76-44-8	Gamma-HCH (Lindane)	CAS:58-89-9
Dieldrin	CAS:60-57-1	Heptachlor-exo-epoxyde	CAS:1024-57-3	Methoxychlor (DMTD)	CAS:72-43-5
Solvent : Toluene		ampoule 1 ml	100 µg/ml	Ref : F113151	Price : EF

PCBs Standard Solution - 1 component

Solvent : Methanol

Aroclor 1016	CAS:12674-11-2	ampoule 1 ml	100 µg/ml	Ref : F113211	Price : BH
Aroclor 1221	CAS:11104-28-2	ampoule 1 ml	100 µg/ml	Ref : F113221	Price : BG
Aroclor 1232	CAS:11141-16-5	ampoule 1 ml	100 µg/ml	Ref : F113231	Price : CC
Aroclor 1242	CAS:53469-21-9	ampoule 1 ml	100 µg/ml	Ref : F113241	Price : BH
Aroclor 1248	CAS:12672-29-6	ampoule 1 ml	100 µg/ml	Ref : F113251	Price : BJ
Aroclor 1254	CAS:11097-69-1	ampoule 1 ml	100 µg/ml	Ref : F113261	Price : BJ
Aroclor 1260	CAS:11096-82-5	ampoule 1 ml	100 µg/ml	Ref : F113271	Price : BJ

Standard Solution - 1 component

Chlordane (technical)	CAS:57-74-9				
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F113191	Price : BG

Standard Solution - 16 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Acenaphthylene	CAS:208-96-8	Benzo(k)fluoranthene	CAS:207-08-9	Naphthalene	CAS:91-20-3
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Phenanthrene	CAS:85-01-8
Benzo(a)anthracene	CAS:56-55-3	Dibenzo(a,c)anthracene	CAS:215-58-7	Pyrene	CAS:129-00-0
Benzo(a)pyrene	CAS:50-32-8	Fluoranthene	CAS:206-44-0		
Benzo(b)fluoranthene	CAS:205-99-2	Fluorene	CAS:86-73-7		
Solvent : Toluene		ampoule 1 ml	100 µg/ml	Ref : F062285	Price : DF
		ampoule 5 ml	100 µg/ml	Ref : F112441	Price : BCA

EPA method 506**Phthalate and Adipate**

Method 506 describes a procedure for the determination of certain phthalate and adipate esters in drinking waters by liquid/liquid or liquid/solid extraction.

Phthalates Standard Solution - 7 components

Adipic acid-bis-2-ethylhexyl ester	CAS:103-23-1	Phthalic acid, bis-ethyl ester	CAS:84-66-2	Phthalic acid, bis-1-octyl ester	CAS:117-84-0
Phthalic acid, benzylbutyl ester	CAS:85-68-7	Phthalic acid,bis-2-ethylhexylester	CAS:117-81-7		
Phthalic acid, bis-butyl ester	CAS:84-74-2	Phthalic acid, bis-methyl ester	CAS:131-11-3		
Solvent : Iso-octane		ampoule 1 ml	1 000 µg/ml	Ref : F113281	Price : EA

Phthalates Standard Solution - 7 components

Adipic acid-bis-2-ethylhexyl ester	CAS:103-23-1	1200 µg/ml	Phthalic acid,bis-2-ethylhexylester	CAS:117-81-7	250 µg/ml
Phthalic acid, benzylbutyl ester	CAS:85-68-7	250 µg/ml	Phthalic acid, bis-methyl ester	CAS:131-11-3	100 µg/ml
Phthalic acid, bis-butyl ester	CAS:84-74-2	100 µg/ml	Phthalic acid, bis-1-octyl ester	CAS:117-84-0	650 µg/ml
Phthalic acid, bis-ethyl ester	CAS:84-66-2	100 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F113311	Price : EF

EPA method 507**Nitrogen and Phosphorus containing pesticides**

Method 507 is applicable to the determination of Nitrogen and Phosphorous containing pesticides in ground water and finishing drinking water by gas chromatography.

NCC Standard Solution - 1 component

Simazine (CAT)	CAS:122-34-9				
Solvent : Acetone		ampoule 1 ml	1 000 µg/ml	Ref : F113411	Price : BK

NCC Surrogate Standard Solution - 1 component

1,3-Dimethyl-2-nitrobenzene CAS:81-20-9

Solvent : Methyl-tert butylether ampoule 1 ml 250 µg/ml Ref : F113321 Price : BH

NCC / OPP Laboratory performance check mixture - 6 components

Atrazine	CAS:1912-24-9	15 µg/ml	Prometon	CAS:1610-18-0	30 µg/ml
Bromacil	CAS:314-40-9	500 µg/ml	Triphenyl phosphate	CAS:115-86-6	250 µg/ml
1,3-Dimethyl-2-nitrobenzene	CAS:81-20-9	250 µg/ml	Vernolate	CAS:1929-77-7	5 µg/ml

Solvent : Methyl-tert butylether ampoule 1 ml Ref : F113431 Price : DF

NCC / OPP Standard Solution - 2 components

Simazine (CAT)	CAS:122-34-9	Paraoxon-methyl	CAS:950-35-6		
----------------	--------------	-----------------	--------------	--	--

Solvent : Acetone ampoule 1 ml 1 000 µg/ml Ref : F113401 Price : CF

NCC / OPP Standard Solution - 6 components

Ametryn	CAS:834-12-8	Disulfoton	CAS:298-04-4	Merphos	CAS:150-50-5
Cycloate	CAS:1134-23-2	Fenamiphos	CAS:22224-92-6	Prometon	CAS:1610-18-0

Solvent : Methyl-tert butylether ampoule 1 ml 1 000 µg/ml Ref : F113341 Price : GF

NCC / OPP Standard Solution - 8 components

Dichlorvos	CAS:62-73-7	Napropamide	CAS:15299-99-7	Tebuthiuron	CAS:34014-18-1
Fenarimol	CAS:60168-88-9	Pebulate	CAS:1114-71-2	Terbacil	CAS:5902-51-2
Fluridone	CAS:59756-60-4	Simetryn	CAS:1014-70-6		

Solvent : Methyl-tert butylether ampoule 1 ml 1 000 µg/ml Ref : F113381 Price : EF

NCC / OPP Standard Solution - 8 components

Benfluralin	CAS:1861-40-1	Oxyfluorfen	CAS:42874-03-3	Propachlor	CAS:1918-16-7
Isopropalin	CAS:33820-53-0	Pendimethalin	CAS:40487-42-1	Trifluralin	CAS:1582-09-8
Oxadiazon	CAS:19666-30-9	Profluralin	CAS:26399-36-0		

Solvent : Methyl-tert butylether ampoule 1 ml 1 000 µg/ml Ref : F113391 Price : EF

NCC / OPP Standard Solution - 9 components

Atrazine	CAS:1912-24-9	Ethoprophos	CAS:13194-48-4	Propazine	CAS:139-40-2
Diphenamid	CAS:957-51-7	Mevinphos	CAS:7786-34-7	Terbutryn	CAS:886-50-0
EPTC	CAS:759-94-4	Prometryn	CAS:7287-19-6	Triadimefon	CAS:43121-43-3

Solvent : Methyl-tert butylether ampoule 1 ml 1 000 µg/ml Ref : F113351 Price : EF

NCC / OPP Standard Solution 9 components

Butachlor	CAS:23184-66-9	Metolachlor	CAS:51218-45-2	Norflurazon	CAS:27314-13-2
Carboxin	CAS:5234-68-4	Metribuzin	CAS:21087-64-9	Terbufos	CAS:13071-79-9
Diazinon	CAS:333-41-5	MGK 264	CAS:113-48-4	Vernolate	CAS:1929-77-7

Solvent : Methyl-tert butylether ampoule 1 ml 1 000 µg/ml Ref : F113361 Price : FA

NCC / OPP Standard Solution 10 components

Alachlor	CAS:15972-60-8	Chlorpropham	CAS:101-21-3	trans-Tetrachlorvinphos	CAS:22248-79-9
Atraton	CAS:1610-17-9	Hexazinone	CAS:51235-04-2	Tricyclazole	CAS:41814-78-2
Bromacil	CAS:314-40-9	Molinat	CAS:2212-67-1		
Butylate	CAS:2008-41-5	Pronamide	CAS:23950-58-5		

Solvent : Methyl-tert butylether ampoule 1 ml 1 000 µg/ml Ref : F113371 Price : FA

OPP Standard Solution - 1 component

Tribufos	CAS:78-48-8				
----------	-------------	--	--	--	--

Solvent : Methyl-tert butylether ampoule 1 ml 1 000 µg/ml Ref : F113421 Price : BK

OPP Standard Solution - 1 component

Triphenyl phosphate	CAS:115-86-6				
---------------------	--------------	--	--	--	--

Solvent : Methyl-tert butylether ampoule 1 ml 500 µg/ml Ref : F113331 Price : BG

EPA method 508, 508.1, 508A

Chlorinated pesticides and PCBs

Method 508 is applicable to the determination of chlorinated pesticides in ground water and finished water by gas chromatography with an electron capture detector.

Method 508.1 is applicable to the determination of chlorinated pesticides, herbicides, and organohalides in ground water and drinking water in any treatment stage using extraction and electron-capture detector gas chromatography.

Method 508A is applicable for screening of Polychlorinated Biphenyls (PCBs) by gas chromatography in raw source water or drinking water in any treatment stage.

EPA method 508

OCs Internal Standard Solution - 1 component

Quintozene	CAS:82-68-8				
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F120131	Price : CA

OCs Laboratory Performance Check Solution - 4 components

Chlorothalonil	CAS:1897-45-6	50 µg/l	Chlorthal-dimethyl	CAS:1861-32-1	50 µg/l
Chlorpyrifos	CAS:2921-88-2	2 µg/l	Delta-HCH	CAS:319-86-8	40 µg/l
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113531	Price : DF

OCs Standard Solution - 1 component

Solvent : Methyl-tert butylether

Chlordane (technical)	CAS:57-74-9	ampoule 1 ml	2 500 µg/ml	Ref : F120051	Price : CA
Toxaphene	CAS:8001-35-2	ampoule 1 ml	2 500 µg/ml	Ref : F113491	Price : CA

OCs Standard Solution - 7 components

Alachlor	CAS:15972-60-8	Hexachlorocyclopentadiene	CAS:77-47-4	Simazine (CAT)	CAS:122-34-9
Atrazine	CAS:1912-24-9	Metolachlor	CAS:51218-45-2		
Cyanazine	CAS:21725-46-2	Metribuzin	CAS:21087-64-9		
Solvent : Ethyl Acetate		ampoule 1 ml	500 µg/ml	Ref : F113461	Price : EA

OCs Standard Solution - 9 components

Chlorobenzilate	CAS:510-15-6	Chlorthal-dimethyl	CAS:1861-32-1	Permethrin	CAS:52645-53-1
Chloroneb	CAS:2675-77-6	Etridiazole	CAS:2593-15-9	Propachlor	CAS:1918-16-7
Chlorothalonil	CAS:1897-45-6	Hexachlorobenzene	CAS:118-74-1	Trifluralin	CAS:1582-09-8
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113481	Price : FA

OCs Standard Solution - 12 components

alpha-Chlordane	CAS:5103-71-9	Chlorothalonil	CAS:1897-45-6	Hexachlorobenzene	CAS:118-74-1
Gamma-Chlordane	CAS:5103-74-2	Chlorpyrifos	CAS:2921-88-2	Permethrin	CAS:52645-53-1
Chlorobenzilate	CAS:510-15-6	DCPA	CAS:1861-32-1	Propachlor	CAS:1918-16-7
Chloroneb	CAS:2675-77-6	Etridiazole	CAS:2593-15-9	Trifluralin	CAS:1582-09-8
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113451	Price : GH
		(Etridiazole : 2 000 µg/ml)			

OCs Standard Solution - 16 components

4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	Delta-HCH	CAS:319-86-8
4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8
4,4'-DDT	CAS:50-29-3	Endrin aldehyde	CAS:7421-93-4	Heptachlor-exo-epoxide	CAS:1024-57-3
Dieldrin	CAS:60-57-1	Alpha-HCH	CAS:319-84-6	Methoxychlor (DMTD)	CAS:72-43-5
Endosulfan-alpha	CAS:959-98-8	Beta-HCH	CAS:319-85-7		
Endosulfan-beta	CAS:33213-65-9	Gamma-HCH (Lindane)	CAS:58-89-9		
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113441	Price : JD

PCBs Standard Solution - 1 component

Solvent : Methyl-tert butylether

Aroclor 1016	CAS:12674-11-2	ampoule 1 ml	1 000 µg/ml	Ref : F120061	Price : CA
Aroclor 1221	CAS:11104-28-2	ampoule 1 ml	1 000 µg/ml	Ref : F120071	Price : BK
Aroclor 1232	CAS:11141-16-5	ampoule 1 ml	1 000 µg/ml	Ref : F120081	Price : EF
Aroclor 1242	CAS:53469-21-9	ampoule 1 ml	1 000 µg/ml	Ref : F120091	Price : CA
Aroclor 1248	CAS:12672-29-6	ampoule 1 ml	1 000 µg/ml	Ref : F120101	Price : CC
Aroclor 1254	CAS:11097-69-1	ampoule 1 ml	1 000 µg/ml	Ref : F120111	Price : CC
Aroclor 1260	CAS:11096-82-5	ampoule 1 ml	1 000 µg/ml	Ref : F120121	Price : CC

PCBs Surrogate Standard Solution - 1 component

Decachlorobiphenyl

CAS:2051-24-3

Solvent : Methyl-tert butylether

ampoule 1 ml

5 000 µg/ml

Ref : F120141

Price : FA

EPA method 508.1**Organochlorine pesticides mixtures****Degradation Calibration Mixture - 2 components**

4,4'-DDT

CAS:50-29-3

200 µg/ml

Endrin

CAS:72-20-8

100 µg/ml

Solvent : Methyl-tert butylether

ampoule 1 ml

Ref : F113561

Price : CF

Degradation Products Mixture - 4 components

4,4'-DDD (TDE)

CAS:72-54-8

200 µg/ml

Endrin aldehyde

CAS:7421-93-4

100 µg/ml

4,4'-DDE

CAS:72-55-9

200 µg/ml

Endrin ketone

CAS:53494-70-5

100 µg/ml

Solvent : Methyl-tert butylether

ampoule 1 ml

Ref : F113571

Price : EA

OCs Internal Standard Solution - 1 component

Pentachloronitrobenzene

CAS:82-68-8

Solvent : Ethyl Acetate

ampoule 1 ml

1 000 µg/ml

Ref : F113611

Price : BH

Solvent : Iso-Octane

ampoule 1 ml

1 000 µg/ml

Ref : F113601

Price : BH

Solvent : Methanol

ampoule 1 ml

100 µg/ml

Ref : F113591

Price : BH

Surrogate Standard Solution - 1 component

4,4'-Dibromobiphenyl

CAS:92-86-4

Solvent : Ethyl Acetate

ampoule 1 ml

1 000 µg/ml

Ref : F113651

Price : CF

Solvent : n-Hexane

ampoule 1 ml

100 µg/ml

Ref : F113691

Price : CF

ampoule 1 ml

1 000 µg/ml

Ref : F147420

Price : CF

OCs Standard Solution - 17 components

Aldrin

CAS:309-00-2

Endosulfan-beta

CAS:33213-65-9

Gamma-HCH (Lindane)

CAS:58-89-9

4,4'-DDD (TDE)

CAS:72-54-8

Endosulfan-total (sulfate)

CAS:1031-07-8

Delta-HCH

CAS:319-86-8

4,4'-DDE

CAS:72-55-9

Endrin

CAS:72-20-8

Heptachlor

CAS:76-44-8

4,4'-DDT

CAS:50-29-3

Endrin aldehyde

CAS:7421-93-4

Heptachlor-exo-epoxide

CAS:1024-57-3

Dieldrin

CAS:60-57-1

Alpha-HCH

CAS:319-84-6

Methoxychlor (DMTD)

CAS:72-43-5

Endosulfan-alpha

CAS:959-98-8

Beta-HCH

CAS:319-85-7

Solvent : Ethyl Acetate

ampoule 1 ml

100 µg/ml

Ref : F113541

Price : GF

OCs /NCS Standard Solution - 19 components

Alachlor

CAS:15972-60-8

Chlorothalonil

CAS:1897-45-6

Metribuzin

CAS:21087-64-9

Atrazine

CAS:1912-24-9

Cyanazine

CAS:21725-46-2

Permethrin

CAS:52645-53-1

Butachlor

CAS:23184-66-9

DCPA

CAS:1861-32-1

Propachlor

CAS:1918-16-7

alpha-Chlordane

CAS:5103-71-9

Etridiazole

CAS:2593-15-9

Simazine (CAT)

CAS:122-34-9

Gamma-Chlordane

CAS:5103-74-2

Hexachlorobenzene

CAS:118-74-1

Trifluralin

CAS:1582-09-8

Chlorobenzilate

CAS:510-15-6

Hexachlorocyclopentadiene

CAS:77-47-4

Chloroneb

CAS:2675-77-6

Metolachlor

CAS:51218-45-2

Solvent : Ethyl Acetate

ampoule 1 ml

100 µg/ml

Ref : F113551

Price : GF

(Permethrin: 200 µg/ml)

Pesticide Degradation Check Solution - 2 components

4,4'-DDT

CAS:50-29-3

Endrin

CAS:72-20-8

Solvent : Ethyl Acetate

ampoule 1 ml

1 µg/ml

Ref : F113581

Price : CF

EPA method 508.A**Aroclor 1260 Standard Solution - 1 component**

Aroclor 1260

CAS:11096-82-5

Solvent : Methanol purge & trap

ampoule 1 ml

100 µg/ml

Ref : F113271

Price : BJ

ampoule 1 ml

1 000 µg/ml

Ref : F113631

Price : CC

ampoule 1 ml

5 000 µg/ml

Ref : F113641

Price : CH

Decachlorobiphenyl Standard Solution - 1 component

PCB 209

CAS:2051-24-3

Solvent : Acetone	ampoule 1 ml	500 µg/ml	Ref : F113681	Price : BG
Solvent : n-Hexane	ampoule 1 ml	100 µg/ml	Ref : F113671	Price : BG
Solvent Toluene	ampoule 1 ml	1 000 µg/ml	Ref : F113661	Price : CE

EPA method 515.1 & 515.2**Chlorinated Acids**

Method 515.1 is applicable to the determination of certain chlorinated acids in ground water and finished drinking water.

Method 515.2 is applicable to the determination of chlorinated acids in ground water, and finished drinking water, using extraction, derivation, and gas chromatography with electron capture detector.

Method 515.3 is applicable to the determination of chlorinated acids in drinking water, ground water, raw source water, and water at any intermediate stage by extraction, derivation, and gas chromatography with electron capture detector.

Chlorinated Herbicides Mixture - 6 components

3,5-Dichlorobenzoic acid	CAS:51-36-5	500 µg/ml	Pentachlorophenol	CAS:87-86-5	100 µg/ml
Dichlorprop	CAS:120-36-5	100 µg/ml	2,4,5-T	CAS:93-76-5	100 µg/ml
Dinoseb	CAS:88-85-7	200 µg/ml	Tetrachloroterephthalic acid	CAS:2136-79-0	100 µg/ml
Solvent : Methanol purge & trap	ampoule 1 ml		Ref : F113831	Price : EA	

Chlorinated Herbicides Mixture - 6 components

3,5-Dichlorobenzoic acid methyl ester	CAS:2905-67-1	500 µg/ml	Dinoseb methyl ether	CAS:6099-79-2	200 µg/ml
Dichlorprop methyl ester	CAS:57153-17-0	100 µg/ml	Pentachloroanisole	CAS:1825-21-4	100 µg/ml
DCPA methyl ester	CAS:1861-32-1	100 µg/ml	2,4,5-T methyl ester	CAS:1928-37-6	100 µg/ml
Solvent : Methanol purge & trap	ampoule 1 ml		Ref : F113841	Price : FA	

Chlorinated Herbicides Acids Mixture - 7 components

Acifluorfen	CAS:50594-66-6	100 µg/ml	Dicamba	CAS:1918-00-9	300 µg/ml
Bentazone	CAS:25057-89-0	1 000 µg/ml	Picloram	CAS:1918-02-1	300 µg/ml
2,4-D	CAS:94-75-7	100 µg/ml	2,4,5-TP	CAS:93-72-1	100 µg/ml
2,4-DB	CAS:94-82-6	1 000 µg/ml			
Solvent : Methanol purge & trap	ampoule 1 ml		Ref : F113851	Price : EA	

Chlorinated Herbicides Mixture - 7 components

Acifluorfen methyl ester	id : 15801	200 µg/ml	Dicamba-methyl ester	CAS:6597-78-0	300 µg/ml
Bentazon methyl derivative	CAS:2225-40-3	1 000 µg/ml	Picloram methyl ester	id : 15841	300 µg/ml
2,4-D methyl ester	CAS:1928-38-7	100 µg/ml	Silvex methyl ester	CAS:4841-20-7	100 µg/ml
2,4-DB methyl ester	CAS:18625-12-2	1 000 µg/ml			
Solvent : Methanol purge & trap	ampoule 1 ml		Ref : F113861	Price : GA	

Chlorinated Herbicides Standard Solution - 16 components

Acifluorfen	CAS:50594-66-6	100 µg/ml	Dichlorprop	CAS:120-36-5	300 µg/ml
Bentazone	CAS:25057-89-0	200 µg/ml	Dinoseb	CAS:88-85-7	200 µg/ml
Chloramben	CAS:133-90-4	100 µg/ml	4-Nitrophenol	CAS:100-02-7	100 µg/ml
2,4-D	CAS:94-75-7	200 µg/ml	Pentachlorophenol	CAS:87-86-5	100 µg/ml
Dalapon	CAS:75-99-0	1 300 µg/ml	Picloram	CAS:1918-02-1	100 µg/ml
2,4-DB	CAS:94-82-6	800 µg/ml	2,4,5-T	CAS:93-76-5	100 µg/ml
Dicamba	CAS:1918-00-9	100 µg/ml	Tetrachloroterephthalic acid	CAS:2136-79-0	100 µg/ml
3,5-Dichlorobenzoic acid	CAS:51-36-5	100 µg/ml	2,4,5-TP	CAS:93-72-1	100 µg/ml
Solvent : Methyl-tert butylether	ampoule 1 ml		Ref : F113711	Price : GA	

Internal Standard Solution - 1 component

4,4'-Dibromooctafluorobiphenyl	CAS:10386-84-2				
Solvent : Methyl-tert butylether	ampoule 1 ml	100 µg/ml	Ref : F113781	Price : BG	
Solvent : Methanol purge & trap	ampoule 1 ml	5 000 µg/ml	Ref : F113911	Price : CE	

Laboratory Performance Check Solution - 5 components

4,4'-Dibromooctafluorobiphenyl	CAS:10386-84-2	2 500 µg/l	Dinoseb	CAS:88-85-7	40 µg/l
3,5-Dichlorobenzoic acid	CAS:51-36-5	6 000 µg/l	4-Nitrophenol	CAS:100-02-7	16 000 µg/l
2,4-Dichlorophenylacetic acid	CAS:19719-28-9	5 000 µg/l			
Solvent : Methyl-tert butylether	ampoule 1 ml		Ref : F113791	Price : FA	

Methyl Derivatives Solution - 5 components

4,4'-Dibromooctafluorobiphenyl	CAS:10386-84-2	250 µg/l	Methyl 3-amino-2,5-dichlorobenzoate	CAS:7286-84-2	600 µg/l
2,4-Dichlorophenylacetic acid methyl ester	id : 15851	500 µg/l	4-Nitroanisole	CAS:100-17-4	1 600 µg/l
Dinoseb methyl ether	CAS:6099-79-2	4 µg/l			
Solvent : Methyl-tert butylether	ampoule 1 ml		Ref : F113811	Price : EA	

Methyl Derivatives Solution - 5 components

4,4'-Dibromooctafluorobiphenyl	CAS:10386-84-2	250 µg/ml	Dinoseb	CAS:88-85-7	4 µg/ml
3,5-Dichlorobenzoic acid	CAS:51-36-5	600 µg/ml	4-Nitrophenol	CAS:100-02-7	16 000 µg/ml
2,4-Dichlorophenylacetic acid	CAS:19719-28-9	500 µg/ml			
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113821	Price : DF

Methyl Derivatives Standard Solution - 16 components

Acifluorfen methyl ester	id : 15801	DCPA	CAS:1861-32-1	Pentachloroanisole	CAS:1825-21-4
Bentazon methyl derivative	CAS:2225-40-3	Dicamba-methyl ester	CAS:6597-78-0	Picloram methyl ester	id : 15841
Chloramben methyl ester	CAS:7286-84-2	3,5-Dichlorobenzoic acid methyl ester	CAS:2905-67-1	Silvex methyl ester	CAS:4841-20-7
2,4-D methyl ester	CAS:1928-38-7	Dichlorprop methyl ester	CAS:57153-17-0	2,4,5-T methyl ester	CAS:1928-37-6
Dalapon methyl ester	CAS:17640-02-7	Dinoseb methyl ether	CAS:6099-79-2		
2,4-DB methyl ester	CAS:18625-12-2	4-Nitroanisole	CAS:100-17-4		
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113751	Price : BCA

Methyl Derivatives Standard Solution - 16 components

Acifluorfen methyl ester	id : 15801	100 µg/ml	3,5-Dichlorobenzoic acid methyl ester	CAS:2905-67-1	100 µg/ml
Bentazon methyl derivative	CAS:2225-40-3	200 µg/ml	Dichlorprop methyl ester	CAS:57153-17-0	300 µg/ml
Chloramben methyl ester	CAS:7286-84-2	100 µg/ml	Dinoseb methyl ether	CAS:6099-79-2	200 µg/ml
2,4-D methyl ester	CAS:1928-38-7	200 µg/ml	4-Nitroanisole	CAS:100-17-4	100 µg/ml
Dalapon methyl ester	CAS:17640-02-7	1 300 µg/ml	Pentachloroanisole	CAS:1825-21-4	100 µg/ml
2,4-DB methyl ester	CAS:18625-12-2	800 µg/ml	Picloram methyl ester	id : 15841	100 µg/ml
DCPA	CAS:1861-32-1	100 µg/ml	Silvex methyl ester	CAS:4841-20-7	100 µg/ml
Dicamba-methyl ester	CAS:6597-78-0	100 µg/ml	2,4,5-T methyl ester	CAS:1928-37-6	100 µg/ml
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113721	Price : HA

SDWA Herbicides Mixture - 6 components

2,4-D	CAS:94-75-7	1 300 µg/ml	Pentachlorophenol	CAS:87-86-5	100 µg/ml
Dalapon	CAS:75-99-0	1 300 µg/ml	Picloram	CAS:1918-02-1	100 µg/ml
Dinoseb	CAS:88-85-7	200 µg/ml	Silvex	CAS:93-72-1	100 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F113871	Price : EA

SDWA Methylated Herbicides Mixture - 6 components

2,4-D methyl ester	CAS:1928-38-7	200 µg/ml	Pentachloroanisole	CAS:1825-21-4	100 µg/ml
Dalapon methyl ester	CAS:17640-02-7	1 300 µg/ml	Picloram methyl ester	id : 15841	100 µg/ml
Dinoseb methyl ether	CAS:6099-79-2	200 µg/ml	Silvex methyl ester	CAS:4841-20-7	100 µg/ml
Solvent : Methanol purge & trap		ampoule 1 ml		Ref : F113881	Price : FA

Surrogate Standard Solution - 1 component

2,4-Dichlorophenylacetic acid	CAS:19719-28-9				
Solvent : Methanol purge & trap		ampoule 1 ml	5 000 µg/ml	Ref : F113891	Price : CE
Solvent : Methyl-tert butylether		ampoule 1 ml	100 µg/ml	Ref : F113761	Price : BG

Surrogate Standard Solution - 1 component

2,4-Dichlorophenylacetic acid methyl ester	id : 15851				
Solvent : Methanol purge & trap		ampoule 1 ml	5 000 µg/ml	Ref : F113901	Price : CJ
Solvent : Methyl-tert butylether		ampoule 1 ml	100 µg/ml	Ref : F113771	Price : BG

Underivatized Analyte - 5 components

4,4'-Dibromooctafluorobiphenyl	CAS:10386-84-2	250 µg/ml	Dinoseb	CAS:88-85-7	4 µg/ml
3,5-Dichlorobenzoic acid	CAS:51-36-5	600 µg/ml	4-Nitrophenol	CAS:100-02-7	1600 µg/ml
2,4-Dichlorophenylacetic acid	CAS:19719-28-9	500 µg/ml			
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113801	Price : DF

EPA method 524.2**Purgeable organic compounds**

Method 524.2 is a purge and trap GC/MS method allowing determination of all VOCs using capillary column.

Haloalkanes Mixture - 35 components

Bromochloromethane	CAS:74-97-5	Dichlorodifluoromethane	CAS:75-71-8	Hexachlorobutadiene	CAS:87-68-3
Bromodichloromethane	CAS:75-27-4	1,1-Dichloroethane	CAS:75-34-3	Methylene chloride	CAS:75-09-2
Bromoform	CAS:75-25-2	1,2-Dichloroethane	CAS:107-06-2	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromomethane	CAS:74-83-9	1,1-Dichloroethene	CAS:75-35-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Carbon tetrachloride	CAS:56-23-5	cis-1,2-Dichloroethene	CAS:156-59-2	Tetrachloroethene	CAS:127-18-4
Chloroethane	CAS:75-00-3	trans-1,2-Dichloroethene	CAS:156-60-5	1,1,1-Trichloroethane	CAS:71-55-6
Chloroform	CAS:67-66-3	1,2-Dichloropropane	CAS:78-87-5	1,1,2-Trichloroethane	CAS:79-00-5
Chloromethane	CAS:74-87-3	1,3-Dichloropropane	CAS:142-28-9	Trichloroethene	CAS:79-01-6
Dibromochloromethane	CAS:124-48-1	2,2-Dichloropropane	CAS:594-20-7	Trichlorofluoromethane	CAS:75-69-4
1,2-Dibromo-3-chloropropane	CAS:96-12-8	1,1-Dichloropropene	CAS:563-58-6	1,2,3-Trichloropropane	CAS:96-18-4
1,2-Dibromoethane	CAS:106-93-4	cis-1,3-Dichloropropene	CAS:10061-01-5	Vinyl chloride	CAS:75-01-4
Dibromomethane	CAS:74-95-3	trans-1,3-Dichloropropene	CAS:10061-02-6		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112991	Price : KF

Haloethanes Mixture - 14 components

Chloroethane	CAS:75-00-3	cis-1,2-Dichloroethene	CAS:156-59-2	1,1,1-Trichloroethane	CAS:71-55-6
1,2-Dibromoethane	CAS:106-93-4	trans-1,2-Dichloroethene	CAS:156-60-5	1,1,2-Trichloroethane	CAS:79-00-5
1,1-Dichloroethane	CAS:75-34-3	1,1,1,2-Tetrachloroethane	CAS:630-20-6	Trichloroethene	CAS:79-01-6
1,2-Dichloroethane	CAS:107-06-2	1,1,2,2-Tetrachloroethane	CAS:79-34-5	Vinylchloride	CAS:75-01-4
1,1-Dichloroethene	CAS:75-35-4	Tetrachloroethene	CAS:127-18-4		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112971	Price : FA

Halopropanes Mixture - 9 components

1,2-Dibromo-3-chloropropane	CAS:96-12-8	2,2-Dichloropropane	CAS:594-20-7	trans-1,3-Dichloropropene	CAS:10061-02-6
1,2-Dichloropropane	CAS:78-87-5	1,1-Dichloropropene	CAS:563-58-6	Hexachlorobutadiene	CAS:87-68-3
1,3-Dichloropropane	CAS:142-28-9	cis-1,3-Dichloropropene	CAS:10061-01-5	1,2,3-Trichloropropane	CAS:96-18-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112981	Price : FA

Internal Standard Solution - 1 component

1,2-Dichlorobenzene D4	CAS:2199-69-1				
Solvent : Methanol		ampoule 1 ml	2 000 µg/ml	Ref : F113961	Price : CJ

Internal Standard Solution - 2 components

1,2-Dichlorobenzene D4	CAS:2199-69-1	Fluorobenzene	CAS:462-06-6		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F113931	Price : CF

Internal Standard Solution - 3 components

4-Bromofluorobenzene	CAS:460-00-4	1,2-Dichlorobenzene D4	CAS:2199-69-1	Fluorobenzene	CAS:462-06-6
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F113941	Price : DF

Surrogate Standard - 1 component

4-Bromofluorobenzene	CAS:460-00-4				
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F113921	Price : CA

VOC Internal Standard Solution - 1 component

Fluorobenzene	CAS:462-06-6				
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112921	Price : DA

VOC Standard Solution - 6 components

Bromomethane	CAS:74-83-9	Chloromethane	CAS:74-87-3	Fluorotrichloromethane	CAS:75-69-4
Chloroethane	CAS:75-00-3	Dichlorodifluoromethane	CAS:75-71-8	Vinylchloride	CAS:75-01-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112901	Price : FA
		ampoule 1 ml	2 000 µg/ml	Ref : F112911	Price : HF

VOC Standard Solution - 12 components

Bromochloromethane	CAS:74-97-5	Carbon tetrachloride	CAS:56-23-5	Dibromomethane	CAS:74-95-3
Bromodichloromethane	CAS:75-27-4	Chloroform	CAS:67-66-3	Dichlorodifluoromethane	CAS:75-71-8
Bromoform	CAS:75-25-2	Chloromethane	CAS:74-87-3	Methylene chloride	CAS:75-09-2
Bromomethane	CAS:74-83-9	Dibromochloromethane	CAS:124-48-1	Trichlorofluoromethane	CAS:75-69-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112961	Price : FA

VOC Standard Solution - 24 components

Acetone	CAS:67-64-1	Diethylether	CAS:60-29-7	Methyl methacrylate	CAS:80-62-6
Acrylonitrile	CAS:107-13-1	Ethyl methacrylate	CAS:97-63-2	MIBK	CAS:108-10-1
Allyl chloride	CAS:107-05-1	Hexachloroethane	CAS:67-72-1	MTBE	CAS:1634-04-4
Carbon disulfide	CAS:75-15-0	2-Hexanone	CAS:591-78-6	Nitrobenzene	CAS:98-95-3
Chloroacetone	CAS:107-14-2	MEK	CAS:78-93-3	2-Nitropropane	CAS:79-46-9
1-Chlorobutane	CAS:109-69-3	Methacrylonitrile	CAS:126-98-7	Pentachloroethane	CAS:76-01-7
trans-1,4-Dichloro-2-butene	CAS:110-57-6	Methyl acrylate	CAS:96-33-3	Propionitrile	CAS:107-12-0
1,1-Dichloro-2-propanone	CAS:513-88-2	Methyl iodide	CAS:74-88-4	THF	CAS:109-99-9
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F113951	Price : GF

VOC Standard Solution - 54 components

Benzene	CAS:71-43-2	1,3-Dichlorobenzene	CAS:541-73-1	Naphthalene	CAS:91-20-3
Bromobenzene	CAS:108-86-1	1,4-Dichlorobenzene	CAS:106-46-7	n-Propylbenzene	CAS:103-65-1
Bromochloromethane	CAS:74-97-5	1,1-Dichloroethane	CAS:75-34-3	Styrene	CAS:100-42-5
Bromodichloromethane	CAS:75-27-4	1,2-Dichloroethane	CAS:107-06-2	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromoform	CAS:75-25-2	1,1-Dichloroethene	CAS:75-35-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
n-Butylbenzene	CAS:104-51-8	cis-1,2-Dichloroethene	CAS:156-59-2	Tetrachloroethene	CAS:127-18-4
sec-Butylbenzene	CAS:135-98-8	trans-1,2-Dichloroethene	CAS:156-60-5	Toluene	CAS:108-88-3
tert-Butylbenzene	CAS:98-06-6	1,2-Dichloropropane	CAS:78-87-5	1,2,3-Trichlorobenzene	CAS:87-61-6
Carbon tetrachloride	CAS:56-23-5	1,3-Dichloropropane	CAS:142-28-9	1,2,4-Trichlorobenzene	CAS:120-82-1
Chlorobenzene	CAS:108-90-7	2,2-Dichloropropane	CAS:594-20-7	1,1,1-Trichloroethane	CAS:71-55-6
Chloroform	CAS:67-66-3	1,1-Dichloropropene	CAS:563-58-6	1,1,2-Trichloroethane	CAS:79-00-5
2-Chlorotoluene	CAS:95-49-8	cis-1,3-Dichloropropene	CAS:10061-01-5	Trichloroethene	CAS:79-01-6
4-Chlorotoluene	CAS:106-43-4	trans-1,3-Dichloropropene	CAS:10061-02-6	1,2,3-Trichloropropane	CAS:96-18-4

Dibromochloromethane	CAS:124-48-1	Ethylbenzene	CAS:100-41-4	1,2,4-Trimethylbenzene	CAS:95-63-6
1,2-Dibromo-3-chloropropane	CAS:96-12-8	Hexachlorobutadiene	CAS:87-68-3	1,3,5-Trimethylbenzene	CAS:108-67-8
1,2-Dibromoethane	CAS:106-93-4	Isopropylbenzene	CAS:98-82-8	o-Xylene	CAS:95-47-6
Dibromomethane	CAS:74-95-3	4-Isopropyltoluene	CAS:99-87-6	m-Xylene	CAS:108-38-3
1,2-Dichlorobenzene	CAS:95-50-1	Methylene chloride	CAS:75-09-2	p-Xylene	CAS:106-42-3
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112781	Price : BDA

VOC Standard Solution - 60 components

Benzene	CAS:71-43-2	1,2-Dichlorobenzene	CAS:95-50-1	Naphthalene	CAS:91-20-3
Bromobenzene	CAS:108-86-1	1,3-Dichlorobenzene	CAS:541-73-1	n-Propylbenzene	CAS:103-65-1
Bromochloromethane	CAS:74-97-5	1,4-Dichlorobenzene	CAS:106-46-7	Styrene	CAS:100-42-5
Bromodichloromethane	CAS:75-27-4	Dichlorodifluoromethane	CAS:75-71-8	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromoform	CAS:75-25-2	1,1-Dichloroethane	CAS:75-34-3	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Bromomethane	CAS:74-83-9	1,2-Dichloroethane	CAS:107-06-2	Tetrachloroethene	CAS:127-18-4
n-Butylbenzene	CAS:104-51-8	1,1-Dichloroethene	CAS:75-35-4	Toluene	CAS:108-88-3
sec-Butylbenzene	CAS:135-98-8	cis-1,2-Dichloroethene	CAS:156-59-2	1,2,3-Trichlorobenzene	CAS:87-61-6
tert-Butylbenzene	CAS:98-06-6	trans-1,2-Dichloroethene	CAS:156-60-5	1,2,4-Trichlorobenzene	CAS:120-82-1
Carbon tetrachloride	CAS:56-23-5	1,2-Dichloropropane	CAS:78-87-5	1,1,1-Trichloroethane	CAS:71-55-6
Chlorobenzene	CAS:108-90-7	1,3-Dichloropropane	CAS:142-28-9	1,1,2-Trichloroethane	CAS:79-00-5
Chloroethane	CAS:75-00-3	2,2-Dichloropropane	CAS:594-20-7	Trichloroethene	CAS:79-01-6
Chloroform	CAS:67-66-3	1,1-Dichloropropene	CAS:563-58-6	Trichlorofluoromethane	CAS:75-69-4
Chloromethane	CAS:74-87-3	cis-1,3-Dichloropropene	CAS:10061-01-5	1,2,3-Trichloropropane	CAS:96-18-4
2-Chlorotoluene	CAS:95-49-8	trans-1,3-Dichloropropene	CAS:10061-02-6	1,2,4-Trimethylbenzene	CAS:95-63-6
4-Chlorotoluene	CAS:106-43-4	Ethylbenzene	CAS:100-41-4	1,3,5-Trimethylbenzene	CAS:108-67-8
Dibromochloromethane	CAS:124-48-1	Hexachlorobutadiene	CAS:87-68-3	Vinylchloride	CAS:75-01-4
1,2-Dibromoethane	CAS:106-93-4	Isopropylbenzene	CAS:98-82-8	m-Xylene	CAS:108-38-3
1,2-Dibromo-3-chloropropane	CAS:96-12-8	4-Isopropyltoluene	CAS:99-87-6	o-Xylene	CAS:95-47-6
Dibromomethane	CAS:74-95-3	Methylene chloride	CAS:75-09-2	p-Xylene	CAS:106-42-3
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112751	Price : BGA

EPA method 525.1

Determination of Organic Compounds in Drinking Water by Liquid-Solid Extraction and Capillary Column Gas Chromatography/Mass Spectrometry

Method 525.1 is applicable to the determination of purgeable volatile organic compounds in surface water, ground water, and drinking water in any treatment stage.

Extractables Mixture - 9 components

Bis(2-ethylhexyl)adipate	CAS:103-23-1	Di-n-butylphthalate	CAS:84-74-2	HCB	CAS:118-74-1
Bis(2-ethylhexyl)phthalate	CAS:117-81-7	Diethyl phthalate	CAS:84-66-2	Hexachlorocyclopentadiene	CAS:77-47-4
Butyl benzyl phthalate	CAS:85-68-7	Dimethyl phthalate	CAS:131-11-3	Pentachlorophenol	CAS:87-86-5
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F114011	Price : EF
			(Pentachlorophenol: 400 µg/ml)		
		ampoule 1 ml	500 µg/ml	Ref : F114021	Price : EF
			(Pentachlorophenol : 2 000 µg/ml)		

Fortification Solution for Optional Recovery Standard -1 component

p-Terphenyl D14	CAS:1718-51-0				
Solvent : Dichloromethane		ampoule 1 ml	500 µg/ml	Ref : F114101	Price : BK

Fortification Solution of Internal Standards and Surrogates - 4 components

Acenaphthene D10	CAS:15067-26-2	Perylene D12	CAS:1520-96-3	Phenanthrene D10	CAS:1517-22-2
Chrysene D12	CAS:1719-03-5				
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F114081	Price : DF

MS Performance Check Solution - 1 component

Decafluorotriphenylphosphine	CAS:5074-71-5				
Solvent : Dichloromethane		ampoule 1 ml	500 µg/ml	Ref : F114071	Price : BG

OCs Standard Solution - 12 components

Alachlor	CAS:15972-60-8	Gamma-Chlordane	CAS:5103-74-2	Heptachlor-exo-epoxide	CAS:1024-57-3
Aldrin	CAS:309-00-2	Gamma-HCH (Lindane)	CAS:58-89-9	Methoxychlor (DMTD)	CAS:72-43-5
Atrazine	CAS:1912-24-9	Endrin	CAS:72-20-8	trans-Nonachlor	CAS:39765-80-5
alpha-Chlordane	CAS:5103-71-9	Heptachlor	CAS:76-44-8	Simazine (CAT)	CAS:122-34-9
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F113991	Price : GA
		ampoule 1 ml	500 µg/ml	Ref : F114001	Price : BCA

PAH Surrogate Standard Solution - 1 component

Pyrene D10	CAS:1718-52-1				
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F114091	Price : CE

PAH Standard Solution - 13 components

Acenaphthylene	CAS:208-96-8	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Anthracene	CAS:120-12-7	Benzo(a)pyrene	CAS:50-32-8	Phenanthrene	CAS:85-01-8
Benzo(a)anthracene	CAS:56-55-3	Chrysene	CAS:218-01-9	Pyrene	CAS:129-00-0
Benzo(j)fluoranthene	CAS:205-82-3	Dibenzo(a,h)anthracene	CAS:53-70-3		
Benzo(k)fluoranthene	CAS:207-08-9	Fluorene	CAS:86-73-7		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F113971	Price : EF
Solvent : Toluene		ampoule 1 ml	500 µg/ml	Ref : F113981	Price : KF

PCBs Standard Solution - 8 components

PCB 1	CAS:2051-60-7	PCB 47	CAS:2437-79-8	PCB 171	CAS:52663-71-5
PCB 5	CAS:16605-91-7	PCB 98	CAS:60233-25-2	PCB 200	CAS:40186-71-8
PCB 29	CAS:15862-07-4	PCB 154	CAS:60145-22-4		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F114031	Price : GH
		ampoule 1 ml	500 µg/ml	Ref : F114041	Price : BKA

SDWA SOC's Mixture - 6 components

Benzo(a)pyrene	CAS:50-32-8	Bis(2-ethylhexyl)phthalate	CAS:117-81-7	Hexachlorocyclopentadiene	CAS:77-47-4
Bis(2-ethylhexyl)adipate	CAS:103-23-1	Hexachlorobenzene	CAS:118-74-1	Pentachlorophenol	CAS:87-86-5
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F114051	Price : EF
		(Pentachlorophenol: 2 000 µg/ml)			

Toxaphene Calibration Standard Solution - 1 component

Toxaphene	CAS:8001-35-2				
Solvent : Acetone		ampoule 1 ml	2 500 µg/ml	Ref : F113495	Price : BH

EPA method 525.2**Determination of organic compounds in drinking water by liquid-solid extraction and capillary column Gas Chromatography/ Mass Spectrometry**

This is a general purpose method that provides procedures for determination of organic compounds in finished drinking water, source water, or drinking water in any treatment stage.

GC/MS Performance Check Solution - 3 components

4,4'-DDT	CAS:50-29-3	Decafluorotriphenylphosphine	CAS:5074-71-5	Endrin	CAS:72-20-8
Solvent : Dichloromethane		ampoule 1 ml	5 µg/ml	Ref : F114141	Price : CF

PAH Internal & Surrogate Standard Fortification Solution - 3 components

Acenaphthene D10	CAS:15067-26-2	Chrysene D12	CAS:1719-03-5	Phenanthrene D10	CAS:1517-22-2
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F114111	Price : DF

PAH Standard Solution - 6 components

Acenaphthene D10	CAS:15067-26-2	1,3-Dimethyl-2-nitrobenzene	CAS:81-20-9	Phenanthrene D10	CAS:1517-22-2
Chrysene D12	CAS:1719-03-5	Perylene D12	CAS:1520-96-3	Triphenylphosphate	CAS:115-86-6
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F114121	Price : EA

PAH Standard Solution - 7 components

Acenaphthene D10	CAS:15067-26-2	Perylene D12	CAS:1520-96-3	Pyrene D10	CAS:1718-52-1
Chrysene D12	CAS:1719-03-5	Phenanthrene D10	CAS:1517-22-2	Triphenylphosphate	CAS:115-86-6
1,3-Dimethyl-2-nitrobenzene	CAS:81-20-9				
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F114131	Price : FF

Standard Solution - 3 components

4,4'-DDT	CAS:50-29-3	Decafluorotriphenylphosphine	CAS:5074-71-5	Endrin	CAS:72-20-8
Solvent : Dichloromethane		ampoule 1 ml	500 µg/ml	Ref : F114151	Price : DA

Toxaphene Calibration Standard Solution - 1 component

Toxaphene	CAS:8001-35-2				
Solvent : Acetone		ampoule 1 ml	2 500 µg/ml	Ref : F113495	Price : CA

EPA method 531.1**Measurement of N-methylcarbamoyloximes and N-methylcarbamates in water by direct aqueous injection HPLC with pst column derivatization**

This is a high performance liquid chromatographic (HPLC) method applicable to the determinations of certain N-methylcarbamoyloximes and N-methylcarbamates in ground water and finished drinking water.

Carbamates Internal Standard Solution - 1 component

BDMC	CAS:672-99-1				
Solvent : Methanol		ampoule 10 ml	100 µg/ml	Ref : F114171	Price : BG

Carbamates Laboratory Performances Check Solution - 4 components

Aldicarb-sulfoxide	CAS:1646-87-3	100 µg/ml	Methiocarb	CAS:2032-65-7	20 µg/ml
Carbofuran-3-hydroxy	CAS:16655-82-6	2 µg/ml	BDMC	CAS:672-99-1	10 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F114181	Price : DF

Carbamates Standard Solution - 10 components

Aldicarb	CAS:116-06-3	Carbaryl	CAS:63-25-2	Methomyl	CAS:16752-77-5
Aldicarb-sulfone	CAS:1646-88-4	Carbofuran	CAS:1563-66-2	Oxamyl	CAS:23135-22-0
Aldicarb-sulfoxide	CAS:1646-87-3	Carbofuran-3-hydroxy	CAS:16655-82-6		
Baygon	CAS:114-26-1	Methiocarb	CAS:2032-65-7		
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F114161	Price : FA

EPA method 547**Determination of glyphosate in drinking water by direct aqueous injection HPLC, post-column derivatization, and fluorescence detection**

This method describes a procedure for the identification and measurement of Glyphosate (N-phosphonomethyl glycine) in drinking water matrices.

Standard Solution - 1 component

Glyphosate	CAS:1071-83-6				
Solvent : Water		ampoule 1 ml	100 µg/ml	Ref : F114201	Price : BG

EPA method 549.2**Determination of Diquat and Paraquat in drinking water by liquid-solid extraction and HPLC**

... This is a high performance liquid chromatography (HPLC) method for the determination of diquat (1,1'-ethylene-2,2'-bipyridilium dibromide salt) and paraquat (1,1'-dimethyl-4,4'-bipyridilium dichloride salt) in drinking water sources and finished drinking water.

Standard Solution - 2 components

Diquat (corrected from dibromide)	CAS:83-32-9	Paraquat (corrected from dichloride)	CAS:208-96-8		
Solvent : Water		ampoule 1 ml	1 000 µg/ml	Ref : F114211	Price : DA

EPA method 550 & 550.1**Determination of polycyclic aromatic hydrocarbons in drinking water by extraction and HPLC with coupled ultraviolet and fluorescence detection**

These methods describe a procedure for determination of certain polycyclic aromatic hydrocarbons (PAH) in drinking water sources and finished drinking water. Method 550 uses liquid-liquid, and 550.1 - liquid-solid extraction.

PAH Standard Solution - 16 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Acenaphthylene	CAS:208-96-8	Benzo(k)fluoranthene	CAS:207-08-9	Naphthalene	CAS:91-20-3
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Phenanthrene	CAS:85-01-8
Benzo(a)anthracene	CAS:56-55-3	Dibenzo(a,h)anthracene	CAS:53-70-3	Pyrene	CAS:129-00-0
Benzo(a)pyrene	CAS:50-32-8	Fluoranthene	CAS:206-44-0		
Benzo(b)fluoranthene	CAS:205-99-2	Fluorene	CAS:86-73-7		
Solvent : Acetonitrile		ampoule 1 ml	1 000 µg/ml	Ref : F114231	Price : BDF

PAH Standard Solution - 16 components

Acenaphthene	CAS:83-32-9	1 000 µg/ml	Chrysene	CAS:218-01-9	625 µg/ml
Acenaphthylene	CAS:208-96-8	1 000 µg/ml	Dibenzo(a,h)anthracene	CAS:53-70-3	1.25 µg/ml
Anthracene	CAS:120-12-7	625 µg/ml	Fluoranthene	CAS:206-44-0	2.5 µg/ml
Benzo(a)anthracene	CAS:56-55-3	1 µg/ml	Fluorene	CAS:86-73-7	100 µg/ml
Benzo(a)pyrene	CAS:50-32-8	5 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	12.5 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	1 µg/ml	Naphthalene	CAS:91-20-3	1 000 µg/ml
Benzo(g,h,i)perylene	CAS:191-24-2	5 µg/ml	Phenanthrene	CAS:85-01-8	50 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	1.25 µg/ml	Pyrene	CAS:129-00-0	62.5 µg/ml
Solvent : Acetonitrile		ampoule 1 ml		Ref : F114221	Price : GF

EPA method 552, 552.1, 552.2, 552.3

Determination of haloacetic acids and dalapon in drinking water by extraction, derivatization, and Gas Chromatography with electron capture detection

This is a gas chromatographic (GC) method applicable to the determination of halogenated acetic acids and Dalapon in drinking water, ground water, raw water and any intermediate treatment stage.

Haloacetic Acids Mixture - 8 components

Bromochloroacetic acid	CAS:5589-96-8	2,4-Dichlorophenol	CAS:120-83-2	Trichloroacetic acid	CAS:76-03-9
Dibromoacetic acid	CAS:631-64-1	Monobromoacetic acid	CAS:79-08-3	2,4,6-Trichlorophenol	CAS:88-06-2
Dichloroacetic acid	CAS:79-43-6	Monochloroacetic acid	CAS:79-11-8		
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F114241	Price : FA

Haloacetic Acids Mixture without Surrogate - 10 components

Bromochloroacetic acid	CAS:5589-96-8	40 µg/ml	Dichloroacetic acid	CAS:79-43-6	60 µg/ml
Bromodichloroacetic acid	CAS:71133-14-7	40 µg/ml	Monobromoacetic acid	CAS:79-08-3	40 µg/ml
Chlorodibromoacetic acid	CAS:5278-95-5	100 µg/ml	Monochloroacetic acid	CAS:79-11-8	60 µg/ml
Dalapon	CAS:75-99-0	40 µg/ml	Tribromoacetic acid	CAS:75-96-7	200 µg/ml
Dibromoacetic acid	CAS:631-64-1	20 µg/ml	Trichloroacetic acid	CAS:76-03-9	20 µg/ml
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F114571	Price : FA

Haloacetic Acids Mixture with Surrogate - 11 components

Bromochloroacetic acid	CAS:5589-96-8	40 µg/ml	2,3-Dibromopropionic acid	CAS:600-05-5	100 µg/ml
Bromodichloroacetic acid	CAS:71133-14-7	40 µg/ml	Monobromoacetic acid	CAS:79-08-3	40 µg/ml
Chlorodibromoacetic acid	CAS:5278-95-5	100 µg/ml	Monochloroacetic acid	CAS:79-11-8	60 µg/ml
Dalapon	CAS:75-99-0	40 µg/ml	Tribromoacetic acid	CAS:75-96-7	200 µg/ml
Dibromoacetic acid	CAS:631-64-1	20 µg/ml	Trichloroacetic acid	CAS:76-03-9	20 µg/ml
Dichloroacetic Acid	CAS:79-43-6	60 µg/ml			
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F114431	Price : FA

Haloorganic Acids Surrogate Standard - 1 component

Solvent : Methyl-tert butylether

2,3-Dibromopropionic acid	CAS:600-05-5	ampoule 1 ml	1 000 µg/ml	Ref : F114611	Price : BH
2,3-Dibromopropionic acid methyl ester	CAS:1729-67-5	ampoule 1 ml	1 000 µg/ml	Ref : F114601	Price : BH
3,5-Dichlorobenzoic acid	CAS:51-36-5	ampoule 1 ml	1 000 µg/ml	Ref : F114621	Price : BK
3,5-Dichlorobenzoic acid methyl ester	CAS:2905-67-1	ampoule 1 ml	1 000 µg/ml	Ref : F114631	Price : EC

Methylated Haloacetic Acids Mixture - 8 components

Bromoacetic acid-methyl ester	CAS:96-32-2	Dibromoacetic acid-methyl ester	CAS:6482-26-4	Trichloroacetic acid-methyl ester	CAS:598-99-2
Bromochloroacetic acid-methyl ester	CAS:20428-74-4	Dichloroacetic acid-methyl ester	CAS:116-54-1	2,4,6-Trichloroanisole	CAS:87-40-1
Chloroacetic acid methyl ester	CAS:96-34-4	2,4-Dichloroanisole	CAS:553-82-2		
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F114251	Price : HA

Methylated Haloacetic Acids Mixture without Surrogate - 10 components

Bromodichloroacetic acid-methyl ester	CAS:20428-76-6	40 µg/ml	Methyl bromochloroacetate	CAS:20428-74-4	40 µg/ml
Dalapon methyl ester	CAS:17640-02-7	40 µg/ml	Methyl chloroacetate	CAS:96-34-4	60 µg/ml
Bromochloroacetic acid-methyl ester	CAS:20428-75-5	100 µg/ml	Methyl dibromoacetate	CAS:6482-26-4	20 µg/ml
Dichloroacetic acid-methyl ester	CAS:116-54-1	60 µg/ml	Methyl tribromoacetate	CAS:3222-05-7	200 µg/ml
Methyl bromoacetate	CAS:96-32-2	40 µg/ml	Methyl trichloroacetate	CAS:598-99-2	20 µg/ml
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F114591	Price : FF

Methylated Haloacetic Acids Mixture with Surrogate - 11 components

Bromodichloroacetic acid-methyl ester	CAS:20428-76-6	40 µg/ml	Methyl chloroacetate	CAS:96-34-4	60 µg/ml
Dalapon methyl ester	CAS:17640-02-7	40 µg/ml	Methyl dibromoacetate	CAS:6482-26-4	20 µg/ml
Dibromochloroacetic acid-methyl ester	CAS:20428-75-5	100 µg/ml	Methyl tribromoacetate	CAS:3222-05-7	200 µg/ml
Dichloroacetic acid-methyl ester	CAS:116-54-1	60 µg/ml	Methyl trichloroacetate	CAS:598-99-2	20 µg/ml
Methyl bromoacetate	CAS:96-32-2	40 µg/ml	Methyl 2,3-dibromopropionate	CAS:1729-67-5	100 µg/ml
Methyl bromochloroacetate	CAS:20428-74-4	40 µg/ml			
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F114581	Price : FF

EPA method 554

Determination of carbonyl compounds in drinking water by dinitrophenylhydrazine derivatization and HPLC

This is a high performance liquid chromatographic (HPLC) method optimized for the determination of carbonyl compounds in finished drinking water and raw source water.

Carbonyl Compounds Standard Solution - 12 components

Acetaldehyde	CAS:75-07-0	Decanal	CAS:112-31-2	Nonanal	CAS:124-19-6
Butanal	CAS:123-72-8	Formaldehyde	CAS:50-00-0	Octanal	CAS:124-13-0
Crotonaldehyde	CAS:123-73-9	Heptanal	CAS:111-71-7	Pentanal	CAS:110-62-3
Cyclohexanone	CAS:108-94-1	Hexanal	CAS:66-25-1	Propanal	CAS:123-38-6
Solvent : Acetonitrile		ampoule 1 ml	1 000 µg/ml	Ref : F114641	Price : EF

Derivatised Carbonyl Compounds Mixture - 12 components

Acetaldehyde-DNPH	CAS:1019-57-4	Decanal-DNPH	id : 4110	Nonanal-DNPH	id : 4100
Butanal-DNPH	CAS:1527-98-6	Formaldehyde-DNPH	CAS:1081-15-8	Octanal-DNPH	CAS:1726-77-8
Crotonaldehyde-DNPH	CAS:1527-96-4	Heptanal-DNPH	CAS:2074-05-7	Pentanal-DNPH	CAS:2057-84-3
Cyclohexanone-DNPH	CAS:1589-62-4	Hexaldehyde-DNPH	CAS:1527-97-5	Propanal-DNPH	CAS:725-00-8
Solvent : Methanol/Acetonitrile (4/1)		ampoule 1 ml	1 000 µg/ml	Ref : F114651	Price : BCA

EPA method 555

Determination of chlorinated acids in water by HPLC with a photodiode array ultraviolet detector

This is a high performance liquid chromatographic (HPLC) method for the determination of certain chlorinated acids in ground water and finished drinking water.

Chlorinated Acids Mixture A - 8 components

Acifluorfen	CAS:50594-66-6	2,4-D	CAS:94-75-7	Picloram	CAS:1918-02-1
Bentazone	CAS:25057-89-0	Dicamba	CAS:1918-00-9	Silvex	CAS:93-72-1
Chloramben	CAS:133-90-4	Dichlorprop	CAS:120-36-5		
Solvent : Acetonitrile		ampoule 1 ml	1 000 µg/ml	Ref : F114661	Price : EC

Chlorinated Acids Mixture B - 8 components

2,4-DB	CAS:94-82-6	MCPA	CAS:94-74-6	Pentachlorophenol	CAS:87-86-5
3,5-Dichlorobenzoic acid	CAS:51-36-5	Mecoprop (MCP)	CAS:7085-19-0	2,4,5-T	CAS:93-76-5
Dinoseb	CAS:88-85-7	4-Nitrophenol	CAS:100-02-7		
Solvent : Acetonitrile		ampoule 1 ml	1 000 µg/ml	Ref : F114671	Price : EC

EPA method 556 & 556.1

Determination of carbonyl compounds in drinking water by pentafluorobenzylhydroxylamine derivatization and capillary Gas chromatography with electron capture detection

This is a gas chromatographic method optimized for the determination of carbonyl compounds in finished drinking water and raw source water. The primary difference between method 556 and 556.1 is that crotonaldehyde has been removed from the analyte list in the 556.1 method.

Aldehydes Mixtures (556 & 556.1) - 14 components

Acetaldehyde	CAS:75-07-0	Ethanedial	CAS:107-22-2	Nonanal	CAS:124-19-6
Benzaldehyde	CAS:100-52-7	Formaldehyde	CAS:50-00-0	Octanal	CAS:124-13-0
Butanal	CAS:123-72-8	Heptanal	CAS:111-71-7	Pentanal	CAS:110-62-3
Cyclohexanone	CAS:108-94-1	Hexanal	CAS:66-25-1	Propanal	CAS:123-38-6
Decanal	CAS:112-31-2	Methylglyoxal	CAS:78-98-8		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F114771	Price : EF

Aldehydes Mixtures (556) - 15 components

Acetaldehyde	CAS:75-07-0	Decanal	CAS:112-31-2	Methylglyoxal	CAS:78-98-8
Benzaldehyde	CAS:100-52-7	Ethanedial	CAS:107-22-2	Nonanal	CAS:124-19-6
Butanal	CAS:123-72-8	Formaldehyde	CAS:50-00-0	Octanal	CAS:124-13-0
Crotonaldehyde	CAS:123-73-9	Heptanal	CAS:111-71-7	Pentanal	CAS:110-62-3
Cyclohexanone	CAS:108-94-1	Hexanal	CAS:66-25-1	Propanal	CAS:123-38-6
Solvent : Acetonitrile		ampoule 1 ml	1 000 µg/ml	Ref : F114681	Price : EF

Internal Standard Solution - 1 component

1,2-Dibromopropane	CAS:78-75-1				
Solvent : n-Hexane		ampoule 10 ml	10 µg/ml	Ref : F114791	Price : BG
		ampoule 1 ml	1 000 µg/ml	Ref : F114781	Price : BH

PFBHA Reagent - 1 component

O-(2,3,4,5,6-Pentafluorobenzyl) hydroxyl-amine hydrochloride	CAS:57981-02-9				
Solvent : Water		ampoule 1 ml	15 µg/ml	Ref : F114821	Price : CD

Surrogate Standard Solution - 1 component

2',4',5'-Trifluoroacetophenone	CAS:129322-83-4				
Solvent : Acetonitrile		ampoule 1 ml	2 000 µg/ml	Ref : F114811	Price : CB
		ampoule 1 ml	10 000 µg/ml	Ref : F114801	Price : CK

Working level (Internal Standard) - 1 component

1,2-Dibromopropane	CAS:78-75-1				
Solvent : n-Hexane		ampoule 10 ml	400 µg/ml	Ref : F114831	Price : BG

EPA methods 600 Series

The EPA 600 methods are designed for monitoring organic pollutants in industrial and municipal waste discharges. These methods are cited under the Clean Water Act (CWA).

EPA methods 601

Purgeable halocarbons

This is a purge and trap gas chromatographic (GC/MS) method applicable to the determination of purgeable halocarbons.

Purgeable Halocarbon Mix - 28 components

Bromodichloromethane	CAS:75-27-4	1,3-Dichlorobenzene	CAS:541-73-1	trans-1,3-Dichloropropene	CAS:10061-02-6
Bromoform	CAS:75-25-2	1,4-Dichlorobenzene	CAS:106-46-7	Methylene chloride	CAS:75-09-2
Bromomethane	CAS:74-83-9	Dichlorodifluoromethane	CAS:75-71-8	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Carbon tetrachloride	CAS:56-23-5	1,1-Dichloroethane	CAS:75-34-3	Tetrachloroethene	CAS:127-18-4
Chlorobenzene	CAS:108-90-7	1,2-Dichloroethane	CAS:107-06-2	1,1,1-Trichloroethane	CAS:71-55-6
Chloroethane	CAS:75-00-3	1,1-Dichloroethene	CAS:75-35-4	1,1,2-Trichloroethane	CAS:79-00-5
Chloroform	CAS:67-66-3	trans-1,2-Dichloroethene	CAS:156-60-5	Trichloroethene	CAS:79-01-6
Chloromethane	CAS:74-87-3	1,2-Dichloropropane	CAS:78-87-5	Trichlorofluoromethane	CAS:75-69-4
Dibromochloromethane	CAS:124-48-1	cis-1,3-Dichloropropene	CAS:10061-01-5	Vinylchloride	CAS:75-01-4
1,2-Dichlorobenzene	CAS:95-50-1				
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F114841	Price : HF

VOC Purgeable Internal Standards (Surrogate Standard Mixture) - 3 components

Bromochloromethane	CAS:74-97-5	2-Bromo-1-chloropropane	CAS:3017-95-6	1,4-Dichlorobutane	CAS:110-56-5
Solvent : Methanol purge & trap		ampoule 1 ml	20 µg/ml	Ref : F114861	Price : CF

VOC Purgeable Gas Mixture - 5 components

Bromomethane	CAS:74-83-9	Chloromethane	CAS:74-87-3	Vinylchloride	CAS:75-01-4
Chloroethane	CAS:75-00-3	Dichlorodifluoromethane	CAS:75-71-8		
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F114851	Price : EA

VOC Performance Check Solution - 8 components

Benzene	CAS:71-43-2	1,2-Dichloroethane	CAS:107-06-2	Trichloroethene	CAS:79-01-6
Carbontetrachloride	CAS:56-23-5	1,1-Dichloroethene	CAS:75-35-4	Vinylchloride	CAS:75-01-4
1,4-Dichlorobenzene	CAS:106-46-7	1,1,1-Trichloroethane	CAS:71-55-6		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F114871	Price : EA

EPA method 602

Purgeable aromatics

This is a purge and trap gas chromatographic (GC/MS) method applicable to the determination of various purgeable aromatics.

ACs Surrogate Standard Solution - 1 component

α,α,α -Trifluorotoluene	CAS:98-08-8				
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113071	Price : CF
		ampoule 1 ml	2 000 µg/ml	Ref : F113061	Price : CF

Purgeable Aromatics Mixture - 7 components

Benzene	CAS:71-43-2	1,3-Dichlorobenzene	CAS:541-73-1	Toluene	CAS:108-88-3
Chlorobenzene	CAS:108-90-7	1,4-Dichlorobenzene	CAS:106-46-7		
1,2-Dichlorobenzene	CAS:95-50-1	Ethylbenzene	CAS:100-41-4		
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F114881	Price : DA

EPA method 603

Acrolein and Acrylonitrile

This is a purge and trap gas chromatographic (GC/FID) method applicable to the determination of acrolein and acrylonitrile in municipal and industrial discharges as provided under 40 CFR Part 136.1.

Carbonyl Compounds Standard Solution - 2 components

Acrolein (2-Propenal)	CAS:107-02-8	Acrylonitrile	CAS:107-13-1		
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F114931	Price : CF
		ampoule 1 ml	2 000 µg/ml	Ref : F114941	Price : CF
Solvent : Water		ampoule 1 ml	1 000 µg/ml	Ref : F114951	Price : CF
		ampoule 1 ml	10 000 µg/ml	Ref : F114961	Price : CF

EPA method 604**Phenols**

This method is a flame ionization detector gas chromatographic (GC/FID) method applicable to the determination of phenol and certain substituted phenols in municipal and industrial discharges as provided under 40 CFR Part 136.1. This method describes analytical conditions for derivatization, cleanup, and electron capture detector gas chromatography (ECDGC) that can be used to confirm measurements made by FIDGC.

Phenols Standard Solution - 11 components

4-Chloro-3-methylphenol	CAS:59-50-7	2,4-Dinitrophenol	CAS:51-28-5	Pentachlorophenol	CAS:87-86-5
2-Chlorophenol	CAS:95-57-8	2-Methyl-4,6-dinitrophenol	CAS:534-52-1	Phenol	CAS:108-95-2
2,4-Dichlorophenol	CAS:120-83-2	2-Nitrophenol	CAS:88-75-5	2,4,6-Trichlorophenol	CAS:88-06-2
2,4-Dimethylphenol	CAS:105-67-9	4-Nitrophenol	CAS:100-02-7		
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F114911	Price : EA

Phenols Surrogate Standard Solution - 1 component

2,4,6-Tribromophenol	CAS:118-79-6				
Solvent : Methanol		ampoule 1 ml	2 000 µg/ml	Ref : F114921	Price : CA

EPA method 605**Benzidines**

This is a high performance liquid chromatography (HPLC) method applicable to the determination of certain benzidines in municipal and industrial discharges as provided under 40 CFR Part 136.1.

Standard Solution - 2 components

Benzidine	CAS:92-87-5	3,3'-Dichlorobenzidine	CAS:91-94-1		
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F114971	Price : CF

EPA method 606**Phthalate ester**

This is a gas chromatographic (GC/ECD) method applicable to the determination of certain phthalate esters in municipal and industrial discharges as provided under 40 CFR Part 136.1.

Phthalates QC Check Sample Concentrate - 6 components

Bis(2-ethylhexyl)phthalate	CAS:117-81-7	50 µg/ml	Diethyl phthalate	CAS:84-66-2	25 µg/ml
Butyl benzyl phthalate	CAS:85-68-7	10 µg/ml	Dimethyl phthalate	CAS:131-11-3	25 µg/ml
di-n-butyl phthalate	CAS:84-74-2	25 µg/ml	di-n-octyl phthalate	CAS:117-84-0	50 µg/ml
Solvent : Acetone		ampoule 1 ml		Ref : F115031	Price : DF

Phthalates Standard Solution - 6 components

Bis(2-ethylhexyl)phthalate	CAS:117-81-7	di-n-butyl phthalate	CAS:84-74-2	Dimethyl phthalate	CAS:131-11-3
Butyl benzyl phthalate	CAS:85-68-7	Diethyl phthalate	CAS:84-66-2	di-n-octyl phthalate	CAS:117-84-0
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F115021	Price : DF

EPA method 607**Nitrosamines**

This is a gas chromatographic (GC/NPD) method applicable to the determination of certain nitrosamines in municipal and industrial discharges as provided under 40 CFR 136.1.

NCC Standard Solution - 3 components

N-Nitroso-dimethylamine	CAS:62-75-9	N-Nitroso-diphenylamine	CAS:86-30-6	N-Nitroso-di-n-propylamine	CAS:621-64-7
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F115041	Price : DF

EPA methods 608,608.1 & 608.2**Organochlorine pesticides and PCBs**

Method 608 is a gas chromatographic (GC/ECD) method applicable to the determination of certain organochlorine pesticides and PCBs in municipal and industrial discharges as provided under 40 CFR Part 136.1. Methods 608.1 and 608.2 include additional analytes.

OCs Solution - 1 component

Solvent : Iso-Octane

Chlordane (technical)	CAS:57-74-9	ampoule 1 ml	100 µg/ml	Ref : F115071	Price : BG
Toxaphene	CAS:8001-35-2	ampoule 1 ml	100 µg/ml	Ref : F115061	Price : BG

OCs Standard Solution - 5 components

Chlorothalonil	CAS:1897-45-6	100 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	100 µg/ml
DCPA	CAS:1861-32-1	100 µg/ml	Permethrin	CAS:52645-53-1	200 µg/ml
Dicloran (CNA)	CAS:99-30-9	100 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F115191	Price : DF

OCs Standard Solution - 5 components

Chlorothalonil	CAS:1897-45-6	1 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	40 µg/ml
DCPA	CAS:1861-32-1	3 µg/ml	Permethrin	CAS:52645-53-1	400 µg/ml
Dicloran (CNA)	CAS:99-30-9	2 µg/ml			
Solvent : n-Hexane		ampoule 1 ml		Ref : F115211	Price : DF

OCs Standard Solution - 7 components

Chlorobenzilate	CAS:510-15-6	1,2-Dibromo-3-chloropropane	CAS:96-12-8	Propachlor	CAS:1918-16-7
Chloroneb	CAS:2675-77-6	Etridiazole	CAS:2593-15-9		
Chloropropylate	CAS:5836-10-2	Pentachloronitrobenzene	CAS:82-68-8		
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F115181	Price : FA

OCs Standard Solution - 7 components

Chlorobenzilate	CAS:510-15-6	200 µg/ml	Etridiazole	CAS:2593-15-9	40 µg/ml
Chloroneb	CAS:2675-77-6	40 µg/ml	Pentachloronitrobenzene	CAS:82-68-8	60 µg/ml
Chloropropylate	CAS:5836-10-2	200 µg/ml	Propachlor	CAS:1918-16-7	1000 µg/ml
1,2-Dibromo-3-chloropropane	CAS:96-12-8	40 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F115201	Price : EA

OCs Standard Solution - 16 components

Aldrin	CAS:309-00-2	Endosulfan II	CAS:33213-65-9	Gamma-HCH (Lindane)	CAS:58-89-9
4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	Delta-HCH	CAS:319-86-8
4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8
4,4'-DDT	CAS:50-29-3	Endrin aldehyde	CAS:7421-93-4	Heptachlor epoxide	CAS:1024-57-3
Dieldrin	CAS:60-57-1	Alpha-HCH	CAS:319-84-6		
Endosulfan I	CAS:959-98-8	Beta-HCH	CAS:319-85-7		
Solvent : Iso-Octane		ampoule 1 ml	20 µg/ml	Ref : F115051	Price : FC

OCs QC Check Sample Concentrate - 16 components

Aldrin	CAS:309-00-2	20 µg/ml	Endrin	CAS:72-20-8	100 µg/ml
4,4'-DDD (TDE)	CAS:72-54-8	100 µg/ml	Endrin aldehyde	CAS:7421-93-4	20 µg/ml
4,4'-DDE	CAS:72-55-9	20 µg/ml	Alpha-HCH	CAS:319-84-6	20 µg/ml
4,4'-DDT	CAS:50-29-3	100 µg/ml	Beta-HCH	CAS:319-85-7	20 µg/ml
Dieldrin	CAS:60-57-1	20 µg/ml	Gamma-HCH (Lindane)	CAS:58-89-9	20 µg/ml
Endosulfan I	CAS:959-98-8	20 µg/ml	Delta-HCH	CAS:319-86-8	20 µg/ml
Endosulfan II	CAS:33213-65-9	100 µg/ml	Heptachlor	CAS:76-44-8	20 µg/ml
Endosulfan-total (sulfate)	CAS:1031-07-8	100 µg/ml	Heptachlor epoxide	CAS:1024-57-3	20 µg/ml
Solvent : Iso-Octane		ampoule 1 ml		Ref : F115171	Price : FF

PCBs Standard Solution - 1 component

Solvent : Iso-Octane

PCB 1016	CAS:12674-11-2	ampoule 1 ml	100 µg/ml	Ref : F115081	Price : BH
PCB 1221	CAS:11104-28-2	ampoule 1 ml	100 µg/ml	Ref : F115091	Price : BG
PCB 1232	CAS:11141-16-5	ampoule 1 ml	100 µg/ml	Ref : F115121	Price : CC
PCB 1242	CAS:53469-21-9	ampoule 1 ml	100 µg/ml	Ref : F115131	Price : BH
PCB 1248	CAS:12672-29-6	ampoule 1 ml	100 µg/ml	Ref : F115141	Price : BJ
PCB 1254	CAS:11097-69-1	ampoule 1 ml	100 µg/ml	Ref : F115151	Price : BJ
PCB 1260	CAS:11096-82-5	ampoule 1 ml	100 µg/ml	Ref : F115161	Price : BJ

EPA method 609**Nitroaromatics and isophorone**

This is a gas chromatographic (GC/ECD/FID) method applicable to the determination of certain nitroaromatics and isophorone in municipal and industrial discharges as provided under 40 CFR Part 136.1.

ACs Standard Solution - 4 components

2,4-Dinitrotoluene	CAS:121-14-2	Isophorone	CAS:78-59-1	Nitrobenzene	CAS:98-95-3
2,6-Dinitrotoluene	CAS:606-20-2				
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F115221	Price : DF

EPA method 610

Polynuclear aromatic hydrocarbons

This is a chromatographic method (GC/FID or HPLC) applicable to the determination of certain polynuclear aromatic hydrocarbons (PAH) in municipal and industrial discharges as provided under 40 CFR Part 136.1.

PAH Standard Solution - 16 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Acenaphthylene	CAS:208-96-8	Benzo(k)fluoranthene	CAS:207-08-9	Naphthalene	CAS:91-20-3
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Phenanthrene	CAS:85-01-8
Benzo(a)anthracene	CAS:56-55-3	Dibenzo(a,h)anthracene	CAS:53-70-3	Pyrene	CAS:129-00-0
Benzo(a)pyrene	CAS:50-32-8	Fluoranthene	CAS:206-44-0		
Benzo(b)fluoranthene	CAS:205-99-2	Fluorene	CAS:86-73-7		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F109861	Price : GF
Solvent : Methanol / Dichloromethane (1/1)		ampoule 1 ml	20 µg/ml	Ref : F115521	Price : GF
		ampoule 1 ml	100 µg/ml	Ref : F115511	Price : GF

PAH Standard Solution - 16 components

Acenaphthene	CAS:83-32-9	100 µg/ml	Chrysene	CAS:218-01-9	100 µg/ml
Acenaphthylene	CAS:208-96-8	200 µg/ml	Dibenzo(a,h)anthracene	CAS:53-70-3	200 µg/ml
Anthracene	CAS:120-12-7	100 µg/ml	Fluoranthene	CAS:206-44-0	200 µg/ml
Benzo(a)anthracene	CAS:56-55-3	100 µg/ml	Fluorene	CAS:86-73-7	200 µg/ml
Benzo(a)pyrene	CAS:50-32-8	100 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	100 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	200 µg/ml	Naphthalene	CAS:91-20-3	1000 µg/ml
Benzo(g,h,i)perylene	CAS:191-24-2	200 µg/ml	Phenanthrene	CAS:85-01-8	100 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	100 µg/ml	Pyrene	CAS:129-00-0	100 µg/ml
Solvent : Methanol / Dichloromethane (1/1)		ampoule 1 ml		Ref : F115531	Price : HF

QC PAH Mixture - 16 components

Acenaphthene	CAS:83-32-9	100 µg/ml	Chrysene	CAS:218-01-9	10 µg/ml
Acenaphthylene	CAS:208-96-8	100 µg/ml	Dibenzo(a,h)anthracene	CAS:53-70-3	10 µg/ml
Anthracene	CAS:120-12-7	100 µg/ml	Fluoranthene	CAS:206-44-0	10 µg/ml
Benzo(a)anthracene	CAS:56-55-3	10 µg/ml	Fluorene	CAS:86-73-7	100 µg/ml
Benzo(a)pyrene	CAS:50-32-8	10 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	10 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	10 µg/ml	Naphthalene	CAS:91-20-3	100 µg/ml
Benzo(g,h,i)perylene	CAS:191-24-2	10 µg/ml	Phenanthrene	CAS:85-01-8	100 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	5 µg/ml	Pyrene	CAS:129-00-0	10 µg/ml
Solvent : Acetonitrile		ampoule 1 ml		Ref : F117531	Price : GF

EPA method 611

Haloethers

This is a gas chromatographic (GC/ECD or ECLD) method applicable to the determination of certain haloethers in municipal and industrial discharges as provided under 40 CFR Part 136.1.

QC Haloethers Mixture - 5 components

Bis-(2-chloroethyl)-ether	CAS:111-44-4	Bis-(2-chloro-1-methylethyl)ether	CAS:108-60-1	4-Chlorophenyl phenyl ether	CAS:7005-72-3
Bis-(2-chloroethoxy)-methane	CAS:111-91-1	4-Bromophenyl phenyl ether	CAS:101-55-3		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F115561	Price : DA

EPA method 612

Chlorinated hydrocarbons

This is a gas chromatographic (GC/ECD or ECLD) method applicable to the determination of certain chlorinated hydrocarbons in municipal and industrial discharges as provided under 40 CFR Part 136.1.

ACs Standard Solution - 9 components

2-Chloronaphthalene	CAS:91-58-7	1,4-Dichlorobenzene	CAS:106-46-7	Hexachlorocyclopentadiene	CAS:77-47-4
1,2-Dichlorobenzene	CAS:95-50-1	Hexachlorobenzene	CAS:118-74-1	Hexachloroethane	CAS:67-72-1
1,3-Dichlorobenzene	CAS:541-73-1	Hexachlorobutadiene	CAS:87-68-3	1,2,4-Trichlorobenzene	CAS:120-82-1
Solvent : Methanol / Dichloromethane (1/1)		ampoule 1 ml	100 µg/ml	Ref : F115571	Price : DF

ACs Standard Solution - 9 components

2-Chloronaphthalene	CAS:91-58-7	400 µg/ml	Hexachlorobutadiene	CAS:87-68-3	1 µg/ml
1,2-Dichlorobenzene	CAS:95-50-1	200 µg/ml	Hexachlorocyclopentadiene	CAS:77-47-4	1 µg/ml
1,3-Dichlorobenzene	CAS:541-73-1	200 µg/ml	Hexachloroethane	CAS:67-72-1	1 µg/ml
1,4-Dichlorobenzene	CAS:106-46-7	400 µg/ml	1,2,4-Trichlorobenzene	CAS:120-82-1	100 µg/ml
Hexachlorobenzene	CAS:118-74-1	1 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F115581	Price : DF

EPA method 614

The Determination of Organophosphorus Pesticides in Municipal and Industrial Wastewater

This is a gas chromatographic (GC/NPD or GC/MS) method applicable to the determination of certain organophosphorus pesticides in industrial and municipal discharges as provided under 40 CFR 136.1.

OPP Standard Solution - 8 components

Azinphos-methyl	CAS:86-50-0	Disulfoton	CAS:298-04-4	Parathion-ethyl	CAS:56-38-2
Demeton (O+S)	CAS:8065-48-3	Ethion	CAS:563-12-2	Parathion-methyl	CAS:298-00-0
Diazinon	CAS:333-41-5	Malathion	CAS:121-75-5		
Solvent : Acetone		ampoule 1 ml	200 µg/ml	Ref : F115591	Price : EA

EPA method 614.1

The Determination of Organophosphorus Pesticides in Municipal and Industrial Wastewater

This method covers the determination by GC/NPD or GC/MS of certain organophosphorus pesticides in municipal and industrial wastewater.

OPP Standard Solution - 4 components

Dioxathion	CAS:78-34-2	Ethion	CAS:563-12-2	Terbufos	CAS:13071-79-9
EPN	CAS:2104-64-5				
Solvent : n-Hexane		ampoule 1 ml	200 µg/ml	Ref : F115601	Price : DF

OPP Standard Solution - 4 components

Dioxathion	CAS:78-34-2	10 µg/ml	Ethion	CAS:563-12-2	100 µg/ml
EPN	CAS:2104-64-5	200 µg/ml	Terbufos	CAS:13071-79-9	4 µg/ml
Solvent : n-Hexane		ampoule 1 ml		Ref : F115611	Price : DF

EPA method 615

The Determination of Chlorinated Herbicides in Municipal and Industrial Wastewater

This is a gas chromatographic (GC/ECD or GC/MS) method applicable to the determination of certain chlorinated herbicides, as well as of salts and esters of these compounds in industrial and municipal discharges as provided under 40 CFR 136.1.

Chlorinated Herbicides Mixture - 10 components

2,4-D	CAS:94-75-7	100 µg/ml	Dinoseb	CAS:88-85-7	100 µg/ml
Dalapon	CAS:75-99-0	100 µg/ml	MCPA	CAS:94-74-6	10 000 µg/ml
2,4-DB	CAS:94-82-6	100 µg/ml	Mecoprop (MCP)	CAS:7085-19-0	10 000 µg/ml
Dicamba	CAS:1918-00-9	100 µg/ml	2,4,5-TP	CAS:93-72-1	10 000 µg/ml
Dichlorprop	CAS:120-36-5	100 µg/ml	2,4,5-T	CAS:93-76-5	10 000 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F115621	Price : HF

Methylated Chlorinated Herbicides Mixture - 10 components

2,4-D methyl ester	CAS:1928-38-7	100 µg/ml	Dinoseb methyl ether	CAS:6099-79-2	100 µg/ml
Dalapon methyl ester	CAS:17640-02-7	100 µg/ml	MCPA-methyl ester	CAS:2436-73-9	10 000 µg/ml
2,4-DB methyl ester	CAS:18625-12-2	100 µg/ml	Mecoprop-methyl ester	CAS:23844-56-6	10 000 µg/ml
Dicamba-methyl ester	CAS:6597-78-0	100 µg/ml	2,4,5-TP	CAS:4841-20-7	10 000 µg/ml
Dichlorprop methyl ester	CAS:57153-17-0	100 µg/ml	2,4,5-T methyl ester	CAS:1928-37-6	10 000 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F115631	Price : JF

EPA method 617

The Determination of Organohalide Pesticides and PCBs in Municipal and Industrial Wastewater

This is a gas chromatographic (GC) method applicable to the determination of certain organohalide pesticides and PCBs in industrial and municipal discharges as provided under 40 CFR 136.1.

OCs Solution - 1 component

Solvent : Iso-Octane

Chlordane (technical)	CAS:57-74-9	ampoule 1 ml	1 000 µg/ml	Ref : F115661	Price : BK
Toxaphene	CAS:8001-35-2	ampoule 1 ml	1 000 µg/ml	Ref : F115671	Price : BK

OCs Standard Solution - 9 components

Captan	CAS:133-06-2	Dicofol	CAS:115-32-2	PCNB	CAS:82-68-8
Carbophenothion	CAS:786-19-6	Isodrin	CAS:465-73-6	Perthane	CAS:72-56-0
Dicloran (CNA)	CAS:99-30-9	Mirex	CAS:2385-85-5	Trifluralin	CAS:1582-09-8
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F115651	Price : EF

OCs Standard Solution - 18 components

Aldrin	CAS:309-00-2	4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8
alpha-BHC	CAS:319-84-6	4,4'-DDT	CAS:50-29-3	Endrin ketone	CAS:53494-70-5
beta-BHC	CAS:319-85-7	Dieldrin	CAS:60-57-1	Endrin aldehyde	CAS:7421-93-4
gamma-BHC	CAS:58-89-9	Endosulfan I	CAS:959-98-8	Heptachlor	CAS:76-44-8
delta-BHC	CAS:319-86-8	Endosulfan II	CAS:33213-65-9	Methoxychlor (DMTD)	CAS:72-43-5
4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	cis-Heptachlorepoxyde (cis-, exo-, isomer B)	CAS:1024-57-3
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F115641	Price : BCF

PCBs Standard Solution - 1 component

Solvent : Iso-Octane

PCB 1016	CAS:12674-11-2	ampoule 1 ml	1 000 µg/ml	Ref : F115681	Price : CA
PCB 1221	CAS:11104-28-2	ampoule 1 ml	1 000 µg/ml	Ref : F115691	Price : BK
PCB 1232	CAS:11141-16-5	ampoule 1 ml	1 000 µg/ml	Ref : F115701	Price : EF
PCB 1242	CAS:53469-21-9	ampoule 1 ml	1 000 µg/ml	Ref : F115711	Price : CA
PCB 1248	CAS:12672-29-6	ampoule 1 ml	1 000 µg/ml	Ref : F115721	Price : CC
PCB 1254	CAS:11097-69-1	ampoule 1 ml	1 000 µg/ml	Ref : F115731	Price : CC
PCB 1260	CAS:11096-82-5	ampoule 1 ml	1 000 µg/ml	Ref : F115741	Price : CC

EPA method 619**The Determination of Triazine Pesticides in Municipal and Industrial Wastewater**

This is a gas chromatographic (GC/NPD or GC/MS) method applicable to the determination of certain triazine pesticides in industrial and municipal discharges as provided under 40 CFR 136.1.

NCC Standard Solution - 11 components

Ametryn	CAS:834-12-8	Prometryn	CAS:7287-19-6	Simazine (CAT)	CAS:122-34-9
Atraton	CAS:1610-17-9	Propazine	CAS:139-40-2	Terbutylazine	CAS:5915-41-3
Atrazine	CAS:1912-24-9	Sebumeton	CAS:26259-45-0	Terbutryn	CAS:886-50-0
Prometon	CAS:1610-18-0	Simetryn	CAS:1014-70-6		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F115751	Price : EF

EPA method 622**The Determination of Organophosphorus Pesticides in Municipal and Industrial Wastewater**

This is a gas chromatographic (GC/NPD or GC/MS) method applicable to the determination of certain organophosphorus pesticides in industrial and municipal discharges as provided under 40 CFR 136.

OPP Standard Solution - 1 component

Naled	CAS:300-76-5				
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F115801	Price : BG

OPP Standard Solution - 3 components

Dichlorvos	CAS:62-73-7	10 µg/ml	cis-Tetrachlorvinphos	CAS:961-11-5	500 µg/ml
Phosdrin	CAS:7786-34-7	30 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F115781	Price : CH

OPP Standard Solution - 7 components

Chlorpyrifos-ethyl	CAS:2921-88-2	30 µg/ml	Merphos	CAS:150-50-5	25 µg/ml
Chlorpyrifos methyl	CAS:5598-13-0	30 µg/ml	Parathion-methyl	CAS:298-00-0	30 µg/ml
Diazinon	CAS:333-41-5	60 µg/ml	Ronnel	CAS:299-84-3	30 µg/ml
Ethoprop	CAS:13194-48-4	25 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F115791	Price : EF

OPP Standard Solution - 9 components

Azinphos-methyl	CAS:86-50-0	150 µg/ml	Fenthion	CAS:55-38-9	10 µg/ml
Coumaphos	CAS:56-72-4	150 µg/ml	Phorate	CAS:298-02-2	15 µg/ml
Demeton (O+S)	CAS:8065-48-3	25 µg/ml	Tokuthion	CAS:34643-46-4	50 µg/ml
Disulfoton	CAS:298-04-4	20 µg/ml	Trichloronat	CAS:327-98-0	15 µg/ml
Fensulfthion	CAS:115-90-2	150 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F115771	Price : EF

OPP Standard Solution - 21 components

Azinphos-methyl	CAS:86-50-0	Dichlorvos	CAS:62-73-7	Naled	CAS:300-76-5
Bolstar	CAS:35400-43-2	Disulfoton	CAS:298-04-4	Parathion-methyl	CAS:298-00-0
Chlorpyrifos-ethyl	CAS:2921-88-2	Ethoprop	CAS:13194-48-4	Phorate	CAS:298-02-2
Chlorpyrifos methyl	CAS:5598-13-0	Fensulfthion	CAS:115-90-2	Ronnel	CAS:299-84-3
Coumaphos	CAS:56-72-4	Fenthion	CAS:55-38-9	Stirofos	CAS:22248-79-9
Demeton	CAS:8065-48-3	Merphos	CAS:150-50-5	Tokuthion	CAS:34643-46-4
Diazinon	CAS:333-41-5	Mevinphos	CAS:7786-34-7	Trichloronat	CAS:327-98-0
Solvent : Iso-Octane		ampoule 1 ml	100 µg/ml	Ref : F115761	Price : HF

EPA method 622.1

The Determination of Thiophosphate Pesticides in Municipal and Industrial Wastewater

This is a gas chromatographic (GC/NPD or GC/MS) method applicable to the determination of certain thiophosphate pesticides in municipal and industrial discharges as provided under 40 CFR 136.1.

OPP Standard Solution - 7 components

Aspon	CAS:3244-90-4	Fenitrothion	CAS:122-14-5	Thionazin	CAS:297-97-2
Dichlofenthion	CAS:97-17-6	Fonofos	CAS:944-22-9		
Famphur	CAS:52-85-7	Phosmet	CAS:732-11-6		
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F115811	Price : EF

EPA method 624

Purgeables

This is a purge and trap gas chromatographic/mass spectrometer (GC/MS) method applicable to the determination of a number of purgeable organics in municipal and industrial discharges as provided under 40 CFR Part 136.1.

VOC Standard Solution - 31 components

Benzene	CAS:71-43-2	1,2-Dichlorobenzene	CAS:95-50-1	Methylene chloride	CAS:75-09-2
Bromodichloromethane	CAS:75-27-4	1,3-Dichlorobenzene	CAS:541-73-1	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Bromoform	CAS:75-25-2	1,4-Dichlorobenzene	CAS:106-46-7	Tetrachloroethene	CAS:127-18-4
Bromomethane	CAS:74-83-9	1,1-Dichloroethane	CAS:75-34-3	Toluene	CAS:108-88-3
Carbon tetrachloride	CAS:56-23-5	1,2-Dichloroethane	CAS:107-06-2	1,1,1-Trichloroethane	CAS:71-55-6
Chlorobenzene	CAS:108-90-7	1,1-Dichloroethene	CAS:75-35-4	1,1,2-Trichloroethane	CAS:79-00-5
Chloroethane	CAS:75-00-3	trans-1,2-Dichloroethene	CAS:156-60-5	Trichloroethene	CAS:79-01-6
2-Chloroethyl-vinylether	CAS:110-75-8	1,2-Dichloropropane	CAS:78-87-5	Trichlorofluoromethane	CAS:75-69-4
Chloroform	CAS:67-66-3	cis-1,3-Dichloropropene	CAS:10061-01-5	Vinylchloride	CAS:75-01-4
Chloromethane	CAS:74-87-3	trans-1,3-Dichloropropene	CAS:10061-02-6		
Dibromochloromethane	CAS:124-48-1	Ethylbenzene	CAS:100-41-4		
Solvent : Methanol purge & trap		ampoule 1 ml	20 µg/ml	Ref : F115831	Price : FF
		ampoule 1 ml	200 µg/ml	Ref : F115821	Price : GF

VOC Surrogate and Internal Mixture - 3 components

Bromochloromethane	CAS:74-97-5	2-Bromo-1-chloropropane	CAS:3017-95-6	1,4-Dichlorobutane	CAS:110-56-5
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F115841	Price : CF
		ampoule 1 ml	20 000 µg/ml	Ref : F115851	Price : FA

VOC Surrogate and Internal Mixture - 3 components

4-Bromofluorobenzene	CAS:460-00-4	Fluorobenzene	CAS:462-06-6	Pentafluorobenzene	CAS:363-72-4
Solvent : Methanol purge & trap		ampoule 1 ml	20 000 µg/ml	Ref : F115861	Price : DF

EPA method 625

Base/neutrals and acids

This method covers the determination of a number of organic compounds that are partitioned into an organic solvent and are amenable to gas chromatography.

Base/Neutral Extractables Mixture - 44 components

Acenaphthene	CAS:83-32-9	2-Chloronaphthalene	CAS:91-58-7	Hexachlorobenzene	CAS:118-74-1
Acenaphthylene	CAS:208-96-8	4-Chlorophenyl phenyl ether	CAS:7005-72-3	Hexachlorobutadiene	CAS:87-68-3
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Hexachlorocyclopentadiene	CAS:77-47-4
Azobenzene	CAS:103-33-3	Dibenzo(a,h)anthracene	CAS:53-70-3	Hexachloroethane	CAS:67-72-1
Benzo(a)anthracene	CAS:56-55-3	Di-n-butylphthalate	CAS:84-74-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Benzo(b)fluoranthene	CAS:205-99-2	1,2-Dichlorobenzene	CAS:95-50-1	Isophorone	CAS:78-59-1
Benzo(k)fluoranthene	CAS:207-08-9	1,3-Dichlorobenzene	CAS:541-73-1	Naphthalene	CAS:91-20-3
Benzo(g,h,i)perylene	CAS:191-24-2	1,4-Dichlorobenzene	CAS:106-46-7	Nitrobenzene	CAS:98-95-3
Benzo(a)pyrene	CAS:50-32-8	Diethyl phthalate	CAS:84-66-2	N-Nitroso-dimethylamine	CAS:62-75-9
Bis-(2-chloroethyl)-ether	CAS:111-44-4	Dimethyl phthalate	CAS:131-11-3	N-Nitrosodi-n-propylamine	CAS:621-64-7
Bis-(2-chloroethoxy)-methane	CAS:111-91-1	2,4-Dinitrotoluene	CAS:121-14-2	N-Nitroso-diphenylamine	CAS:86-30-6
Bis(2-ethylhexyl)phthalate	CAS:117-81-7	2,6-Dinitrotoluene	CAS:606-20-2	Phenanthrene	CAS:85-01-8
Bis(2-Chloroisopropyl)ether	CAS:108-60-1	di-n-octyl phthalate	CAS:117-84-0	Pyrene	CAS:129-00-0
4-Bromophenyl phenyl ether	CAS:101-55-3	Fluoranthene	CAS:206-44-0	1,2,4-Trichlorobenzene	CAS:120-82-1
Butyl benzyl phthalate	CAS:85-68-7	Fluorene	CAS:86-73-7		
Solvent : Dichloromethane/Benzene/Acetonitrile (2/2/1)		ampoule 1 ml	2 000 µg/ml	Ref : F115871	Price : BJA

Organochlorine Pesticide Mixture - 16 components

Aldrin	CAS:309-00-2	4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8
alpha-BHC	CAS:319-84-6	4,4'-DDT	CAS:50-29-3	Endrin aldehyde	CAS:7421-93-4
beta-BHC	CAS:319-85-7	Dieldrin	CAS:60-57-1	Heptachlor	CAS:76-44-8
delta-BHC	CAS:319-86-8	Endosulfan I	CAS:959-98-8	cis-Heptachlorepoxyde	CAS:1024-57-3
gamma-BHC	CAS:58-89-9	Endosulfan II	CAS:33213-65-9	(cis-,exo-, isomer B)	
4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8		
Solvent : Hexane/Toluene (1/1)		ampoule 1 ml	2 000 µg/ml	Ref : F115891	Price : BCF

Phenols Standard Solution - 11 components

4-Chloro-3-methylphenol	CAS:59-50-7	2,4-Dinitrophenol	CAS:51-28-5	Pentachlorophenol	CAS:87-86-5
2-Chlorophenol	CAS:95-57-8	2-Methyl-4,6-dinitrophenol	CAS:534-52-1	Phenol	CAS:108-95-2
2,4-Dichlorophenol	CAS:120-83-2	2-Nitrophenol	CAS:88-75-5	2,4,6-Trichlorophenol	CAS:88-06-2
2,4-Dimethylphenol	CAS:105-67-9	4-Nitrophenol	CAS:100-02-7		
Solvent : Dichloromethane		ampoule 1 ml	2 000 µg/ml	Ref : F115881	Price : FA

Semi-Volatiles Internal Standard Mixture - 6 components

Acenaphthene D10	CAS:15067-26-2	1,4-Dichlorobenzene D4	CAS:3855-82-1	Perylene D12	CAS:1520-96-3
Chrysene D12	CAS:1719-03-5	Naphthalene D8	CAS:1146-65-2	Phenanthrene D10	CAS:1517-22-2
Solvent : Dichloromethane		ampoule 1 ml	2 000 µg/ml	Ref : F115901	Price : HF

Surrogate Standard Solution - 1 component

Solvent : Acetone

Aniline D5	CAS:4165-61-1	ampoule 1 ml	100 µg/ml	Ref : F146943	Price : BJ
Anthracene D10	CAS:1719-06-8	ampoule 1 ml	100 µg/ml	Ref : F146771	Price : BJ
2-Fluorophenol	CAS:367-12-4	ampoule 1 ml	100 µg/ml	Ref : F146761	Price : BH
Pentafluorophenol	CAS:771-61-9	ampoule 1 ml	100 µg/ml	Ref : F146781	Price : BG

EPA method 627

The Determination of Dinitroaniline Pesticides in Municipal and Industrial Wastewater

This is a gas chromatographic (GC) method applicable to the determination of certain dinitroaniline pesticides in industrial and municipal discharges as provided under 40 CFR 136.1.

Standard Solution - 5 components

Benfluralin	CAS:1861-40-1	Isopropalin	CAS:33820-53-0	Trifluralin	CAS:1582-09-8
Ethalfuralin	CAS:55283-68-6	Profluralin	CAS:26399-36-0		
Solvent : n-Hexane		ampoule 1 ml	1 000 µg/ml	Ref : F115911	Price : EF

EPA method 632

The Determination of Carbamate and Urea Pesticides in Municipal and Industrial Wastewater

This is a high-performance liquid chromatographic (HPLC) method applicable to the determination of certain carbamate and urea pesticides in industrial and municipal discharges as provided under 40 CFR 136.1.

Carbamates Standard Solution - 21 components

Aminocarb	CAS:2032-59-9	Fenuron TCA	CAS:4482-55-7	Monuron TCA	CAS:140-41-0
Barban	CAS:101-27-9	Fluometuron	CAS:2164-17-2	Neburon	CAS:555-37-3
Carbaryl	CAS:63-25-2	Linuron	CAS:330-55-2	Oxamyl	CAS:23135-22-0
Carbofuran	CAS:1563-66-2	Methiocarb	CAS:2032-65-7	Propham	CAS:122-42-9
Chlorpropham	CAS:101-21-3	Methomyl	CAS:16752-77-5	Propoxur	CAS:114-26-1
Diuron	CAS:330-54-1	Mexacarbate	CAS:315-18-4	Siduron	CAS:1982-49-6
Fenuron	CAS:101-42-8	Monuron	CAS:150-68-5	Swep	CAS:1918-18-9
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F115921	Price : EE

EPA method 632.1

The Determination of Carbamate and Amide Pesticides in Municipal and Industrial Wastewater

This is a high-performance liquid chromatographic (HPLC) method applicable to the determination of certain carbamate/amide pesticides in municipal and industrial discharges.

Carbamates Standard Solution - 2 components

Napropamide	CAS:15299-99-7	Propanil	CAS:709-98-8		
Solvent : Acetonitrile/Acetone (1/1)		ampoule 1 ml	100 µg/ml	Ref : F115941	Price : DF

EPA method 633

The Determination of Organonitrogen Pesticides in Municipal and Industrial Wastewater

This is a gas chromatographic (GC/NPD or GC/MS) method applicable to the determination of certain organonitrogen pesticides in industrial and municipal discharges as provided under 40 CFR 136.1.

NCC Standard Solution - 7 components

Bromacil	CAS:314-40-9	Metribuzin	CAS:21087-64-9	Tricyclazole	CAS:41814-78-2
Deet	CAS:134-62-3	Terbacil	CAS:5902-51-2		
Hexazinone	CAS:51235-04-2	Triadimefon	CAS:43121-43-3		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F115931	Price : DF

EPA method 634

Thiocarbamate Pesticides

This is a gas chromatographic (GC/NPD or GC/MS) method applicable to the determination of thiocarbamate pesticides in industrial and municipal wastewaters.

Carbamates Standard Solution - 6 components

Butylate	CAS:2008-41-5	EPTC	CAS:759-94-4	Pebulate	CAS:1114-71-2
Cycloate	CAS:1134-23-2	Molinate	CAS:2212-67-1	Vernolate	CAS:1929-77-7
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F115951	Price : EF

Carbamates Internal Standard Solution - 1 component

Carbazole	CAS:86-74-8				
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F115961	Price : BK

EPA method 635

Rotenone

This is a high-performance liquid chromatographic (HPLC) method applicable to the determination of rotenone in industrial and municipal wastewaters Rotenone [CAS:83-79-4] 1000ug/ml in Methanol

Standard Solution - 1 component

Rotenone	CAS:83-79-4				
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F115971	Price : CF

EPA method 636

Bensulide

This is a high-performance liquid chromatographic (HPLC) method applicable to the determination of Bensulide pesticide in industrial and municipal discharges as provided under 40 CFR 136.1.

OPP Standard Solution - 1 component

Bensulide	CAS:741-58-2				
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F115981	Price : BK

EPA method 638

Oryzalin

This is a high-performance liquid chromatographic (HPLC) method applicable to the determination of Oryzalin pesticide in industrial and municipal discharges as provided under 40 CFR 136.1.

NCC Standard Solution - 1 component

Oryzalin	CAS:19044-88-3				
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F115991	Price : BK

EPA method 639

Bendiocarb

This is a high-performance liquid chromatographic (HPLC) method applicable to the determination of Bendiocarb pesticide in industrial and municipal discharges as provided under 40 CFR 136.1.

Carbamates Standard Solution - 1 component

Bendiocarb	CAS:22781-23-3				
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F116001	Price : BK

EPA method 641

Thiabendazole

This is a high-performance liquid chromatographic (HPLC) method applicable to the determination of Thiabendazole in industrial and municipal discharges as provided under 40 CFR 136.1.

NCC Standard Solution - 1 component

Thiabendazole	CAS:148-79-8				
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F116011	Price : CF

EPA method 642**Biphenyl and 2-Biphenylol**

This is a high-performance liquid chromatographic (HPLC) method applicable to the determination of Biphenyl and o-Biphenylol in industrial and municipal wastewaters.

ACs Standard Solution - 1 component

Biphenyl	CAS:92-52-4				
Solvent : n-Hexane		ampoule 1 ml	1 000 µg/ml	Ref : F116021	Price : BK

ACs Standard Solution - 1 component

o-Phenylphenol	CAS:90-43-7				
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F116031	Price : BK

EPA method 643**Bentazon**

This is a high-performance liquid chromatographic (HPLC) method applicable to the determination of Bentazon in industrial and municipal wastewaters.

NCC Standard Solution - 1 component

Bentazon	CAS:25057-89-0				
Solvent : Acetonitrile		ampoule 1 ml	1 000 µg/ml	Ref : F116041	Price : BK

EPA method 644**Picloram**

This is a high-performance liquid chromatographic (HPLC) method applicable to the determination of Picloram in industrial and municipal wastewaters.

Internal Standard Solution - 1 component

Salicylic acid	CAS:69-72-7				
Solvent : 2-Propanol		ampoule 1 ml	2 000 µg/ml	Ref : F116061	Price : CA

NCC Standard Solution - 1 component

Picloram	CAS:1918-02-1				
Solvent : Acetonitrile		ampoule 1 ml	1 000 µg/ml	Ref : F116051	Price : BK

EPA methods 1300 Series**EPA method 1311****Toxicity characteristic leaching procedure**

The TCLP is designed to determine the mobility of both organic and inorganic analytes present in liquid, solid, and multiphasic wastes.

TCLP Acids Mixture - 6 components

2-Methylphenol	CAS:95-48-7	4-Methylphenol	CAS:106-44-5	2,4,5-Trichlorophenol	CAS:95-95-4
3-Methylphenol	CAS:108-39-4	Pentachlorophenol	CAS:87-86-5	2,4,6-Trichlorophenol	CAS:88-06-2
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F116081	Price : EA

TCLP Base/Neutrals Mixture - 7 components

1,4-Dichlorobenzene	CAS:106-46-7	Hexachlorobutadiene	CAS:87-68-3	Pyridine	CAS:110-86-1
2,4-Dinitrotoluene	CAS:121-14-2	Hexachloroethane	CAS:67-72-1		
Hexachlorobenzene	CAS:118-74-1	Nitrobenzene	CAS:98-95-3		
Solvent : Acetone		ampoule 1 ml	1 000 µg/ml	Ref : F116071	Price : EA

TCLP Herbicides Spiking Mixture - 2 components

2,4-D	CAS:94-75-7	Silvex	CAS:93-72-1		
Solvent : Methanol		ampoule 1 ml	2 000 µg/ml	Ref : F116151	Price : DF

TCLP methylated Herbicides Spiking Mixture - 2 components

2,4-D methyl ester	CAS:1928-38-7	Silvex	CAS:4841-20-7		
Solvent : Methanol		ampoule 1 ml	2 000 µg/ml	Ref : F116161	Price : DF

TCLP Pesticides Mixture - 5 components

Endrin	CAS:72-20-8	100 µg/ml	Heptachlor	CAS:76-44-8	100 µg/ml
Gamma-HCH (Lindane)	CAS:58-89-9	100 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	1 000 µg/ml
cis-Heptachlorepoxyde (cis, exo-, isomer B)	CAS:1024-57-3	100 µg/ml			
Solvent : Methanol		ampoule 1 ml		Ref : F116111	Price : DF

TCLP Pesticides Spiking Mixture - 2 components

Chlordane (technical)	CAS:57-74-9	2 000 µg/ml	Toxaphene	CAS:8001-35-2	4 000 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F116141	Price : FA

TCLP Pesticides Spiking Mixture - 5 components

Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8	Methoxychlor (DMTD)	CAS:72-43-5
Gamma-HCH (Lindane)	CAS:58-89-9	cis-Heptachlorepoxyde (cis-, exo-, isomer B)	CAS:1024-57-3		
Solvent : Methanol		ampoule 1 ml	2 000 µg/ml	Ref : F116131	Price : GA

TCLP Pesticides Spiking Mixture - 7 components

Chlordane (technical)	CAS:57-74-9	2 000 µg/ml	Heptachlor	CAS:76-44-8	2 000 µg/ml
Endrin	CAS:72-20-8	2 000 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	2 000 µg/ml
Gamma-HCH (Lindane)	CAS:58-89-9	2 000 µg/ml	Toxaphene	CAS:8001-35-2	4 000 µg/ml
cis-Heptachlorepoxyde (cis, exo-, isomer B)	CAS:1024-57-3	2 000 µg/ml			
Solvent : Methanol		ampoule 1 ml		Ref : F116121	Price : GF

TCLP Semi-Volatile Spiking Mixture - 13 components

1,4-Dichlorobenzene	CAS:106-46-7	2-Methylphenol	CAS:95-48-7	Pyridine	CAS:110-86-1
2,4-Dinitrotoluene	CAS:121-14-2	3-Methylphenol	CAS:108-39-4	2,4,5-Trichlorophenol	CAS:95-95-4
Hexachlorobenzene	CAS:118-74-1	4-Methylphenol	CAS:106-44-5	2,4,6-Trichlorophenol	CAS:88-06-2
Hexachlorobutadiene	CAS:87-68-3	Nitrobenzene	CAS:98-95-3		
Hexachloroethane	CAS:67-72-1	Pentachlorophenol	CAS:87-86-5		
Solvent : Dichloromethane		ampoule 1 ml	1 000 µg/ml	Ref : F116101	Price : FF

TCLP Volatiles Mixture - 11 components

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	Tetrachloroethene	CAS:127-18-4
Carbon tetrachloride	CAS:56-23-5	1,2-Dichloroethane	CAS:107-06-2	Trichloroethene	CAS:79-01-6
Chlorobenzene	CAS:108-90-7	1,1-Dichloroethene	CAS:75-35-4	Vinylchloride	CAS:75-01-4
Chloroform	CAS:67-66-3	MEK	CAS:78-93-3		
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F116091	Price : EF

EPA method 1312**Synthetic precipitation leaching procedure**

Method 1312 is designed to determine the mobility of both organic and inorganic analytes present in liquids, soils, and wastes.

OCs / ACs / VOC Standard Solution - 14 components

Acenaphthene	CAS:83-32-9	1,4-Dichlorobenzene	CAS:106-46-7	Gamma-HCH (Lindane)	CAS:58-89-9
Bis-(2-chloroethyl)-ether	CAS:111-44-4	2,4-Dimethylphenol	CAS:105-67-9	Hexachlorobenzene	CAS:118-74-1
2-Chlorophenol	CAS:95-57-8	2,4-Dinitrophenol	CAS:51-28-5	Hexachlorobutadiene	CAS:87-68-3
o-Cresol	CAS:95-48-7	2,4-Dinitrotoluene	CAS:121-14-2	Nitrobenzene	CAS:98-95-3
1,2-Dichlorobenzene	CAS:95-50-1	Beta-HCH	CAS:319-85-7		
Solvent : Dichloromethane		ampoule 1 ml	2 000 µg/ml	Ref : F116171	Price : GA

EPA methods 8000 Series

The EPA 8000 methods are designed for monitoring organic pollutants in ground water, as prescribed in the Resource Recovery and Conservation Act (RCRA).

EPA method 8010B**Halogenated volatile organics by GC**

Method 8010 provides gas chromatographic conditions for the detection of halogenated volatile organic compounds. Samples can be introduced into the GC using direct injection or purge-and-trap (Method 5030). Ground water samples must be analyzed using Method 5030. A temperature program is used in the gas chromatograph to separate the organic compounds. Detection is achieved by a electrolytic conductivity detector (HECD).

Ethers Standard Solution - 1 component

2-Chloroethyl-vinylether	CAS:110-75-8				
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F116211	Price : BG
		ampoule 1 ml	5 000 µg/ml	Ref : F116221	Price : CE

VOC Internal & Surrogate Standard Solution - 3 components

4-Bromochlorobenzene	CAS:106-39-8	Bromochloromethane	CAS:74-97-5	4-Bromofluorobenzene	CAS:460-00-4
Solvent : Methanol purge & trap		ampoule 1 ml	1 500 µg/ml	Ref : F116241	Price : CF

VOC Standard Solution - 1 component

Chloroprene	CAS:126-99-8				
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F116231	Price : BG

VOC Standard Solution - 33 components

Benzylchloride	CAS:100-44-7	1,2-Dichlorobenzene	CAS:95-50-1	trans-1,3-Dichloropropene	CAS:10061-02-6
Bromobenzene	CAS:108-86-1	1,3-Dichlorobenzene	CAS:541-73-1	Methylene chloride	CAS:75-09-2
Bromoform	CAS:75-25-2	1,4-Dichlorobenzene	CAS:106-46-7	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromomethane	CAS:74-83-9	Dichlorobromomethane	CAS:75-27-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Carbon tetrachloride	CAS:56-23-5	Dichlorodifluoromethane	CAS:75-71-8	Tetrachloroethene	CAS:127-18-4
Chlorobenzene	CAS:108-90-7	1,1-Dichloroethane	CAS:75-34-3	1,1,1-Trichloroethane	CAS:71-55-6
Chloroethane	CAS:75-00-3	1,2-Dichloroethane	CAS:107-06-2	1,1,2-Trichloroethane	CAS:79-00-5
Chloroform	CAS:67-66-3	1,1-Dichloroethene	CAS:75-35-4	Trichloroethene	CAS:79-01-6
Chloromethane	CAS:74-87-3	trans-1,2-Dichloroethene	CAS:156-60-5	Trichlorofluoromethane	CAS:75-69-4
Dibromochloromethane	CAS:124-48-1	1,2-Dichloropropane	CAS:78-87-5	1,2,3-Trichloropropane	CAS:96-18-4
Dibromomethane	CAS:74-95-3	cis-1,3-Dichloropropene	CAS:10061-01-5	Vinylchloride	CAS:75-01-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F116181	Price : BAF

VOC Standard Solution - 40 components

Allyl chloride	CAS:107-05-1	Dibromomethane	CAS:74-95-3	trans-1,3-Dichloropropene	CAS:10061-02-6
Bromobenzene	CAS:108-86-1	1,2-Dichlorobenzene	CAS:95-50-1	Ethylene dibromide	CAS:106-93-4
Bromoform	CAS:75-25-2	1,3-Dichlorobenzene	CAS:541-73-1	Methylene chloride	CAS:75-09-2
Bromomethane	CAS:74-83-9	1,4-Dichlorobenzene	CAS:106-46-7	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Carbon tetrachloride	CAS:56-23-5	Dichlorobromomethane	CAS:75-27-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Chlorobenzene	CAS:108-90-7	1,4-Dichloro-2-butene	CAS:764-41-0	Tetrachloroethene	CAS:127-18-4
Chloroethane	CAS:75-00-3	Dichlorodifluoromethane	CAS:75-71-8	1,1,1-Trichloroethane	CAS:71-55-6
2-Chloroethanol	CAS:107-07-3	1,1-Dichloroethane	CAS:75-34-3	1,1,2-Trichloroethane	CAS:79-00-5
Chloroform	CAS:67-66-3	1,2-Dichloroethane	CAS:107-06-2	Trichloroethene	CAS:79-01-6
1-Chlorohexane	CAS:544-10-5	1,1-Dichloroethene	CAS:75-35-4	Trichlorofluoromethane	CAS:75-69-4
Chloromethane	CAS:74-87-3	trans-1,2-Dichloroethene	CAS:156-60-5	1,2,3-Trichloropropane	CAS:96-18-4
4-Chlorotoluene	CAS:106-43-4	1,2-Dichloropropane	CAS:78-87-5	Vinylchloride	CAS:75-01-4
Dibromochloromethane	CAS:124-48-1	1,3-Dichloro-2-propanol	CAS:96-23-1		
1,2-Dibromo-3-chloropropane	CAS:96-12-8	cis-1,3-Dichloropropene	CAS:10061-01-5		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F116201	Price : BBF

EPA method 8011**1,2-Dibromoethane and 1,2-dibromo-3chloropropane by microextraction and GC**

This method is applicable to the determination of 1,2-Dibromoethane and 1,2-Dibromo-3-chloropropane in drinking water and ground water

Standard Solution - 2 components

1,2-Dibromo-3-chloropropane	CAS:96-12-8	1,2-Dibromoethane	CAS:106-93-4		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F116251	Price : CF
		ampoule 1 ml	2 000 µg/ml	Ref : F116261	Price : CF

EPA method 8015B**Nonhalogenated organics using GC/FID**

Method 8015 is used to determine the concentration of various nonhalogenated volatile organic compounds and semivolatiles organic compounds by gas chromatography.

VOC Standard Solution - 4 components

Diethylether	CAS:60-29-7	MEK	CAS:78-93-3	MIBK	CAS:108-10-1
Ethanol	CAS:64-17-5				
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F116271	Price : DF

VOC Standard Solution - 12 components

Acetonitrile	CAS:75-05-8	Ethyl alcohol	CAS:64-17-5	Methacrylonitrile	CAS:126-98-7
Acrylamide	CAS:79-06-1	Ethyl methacrylate	CAS:97-63-2	Methyl methacrylate	CAS:80-62-6
Diethylether	CAS:60-29-7	Isobutyl alcohol	CAS:78-83-1	MIBK	CAS:108-10-1
1,4-Dioxane	CAS:123-91-1	MEK	CAS:78-93-3	Propionitrile	CAS:107-12-0
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F116281	Price : EF

EPA method 8020A

Aromatic volatile organics by GC

Method 8020 provides chromatographic conditions for the detection of aromatic volatile compounds. Samples can be introduced into the GC using direct injection or purge-and-trap (Method 5030). Ground water samples must be determined using Method 5030. A temperature program is used in the gas chromatograph to separate the organic compounds. Detection is achieved by a photo-ionization detector (PID).

ACs Internal Standard Solution - 2 components

2-Bromofluorobenzene	CAS:1072-85-1	α , α , α -Trifluorotoluene	CAS:98-08-8		
Solvent : Methanol purge & trap		ampoule 5 ml	2 000 µg/ml	Ref : F116321	Price : EA

ACs Standard Solution - 10 components

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	m-Xylene	CAS:108-38-3
Chlorobenzene	CAS:108-90-7	Ethylbenzene	CAS:100-41-4	p-Xylene	CAS:106-42-3
1,2-Dichlorobenzene	CAS:95-50-1	Toluene	CAS:108-88-3		
1,3-Dichlorobenzene	CAS:541-73-1	o-Xylene	CAS:95-47-6		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F116291	Price : EF

ACs Standard Solution - 11 components

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	m-Xylene	CAS:108-38-3
Chlorobenzene	CAS:108-90-7	Ethylbenzene	CAS:100-41-4	o-Xylene	CAS:95-47-6
1,2-Dichlorobenzene	CAS:95-50-1	Toluene	CAS:108-88-3	p-Xylene	CAS:106-42-3
1,3-Dichlorobenzene	CAS:541-73-1	Styrene	CAS:100-42-5		
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F116301	Price : DF

ACs Surrogate Standard Solution - 1 component

α , α , α -Trifluorotoluene	CAS:98-08-8				
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113071	Price : CF
		ampoule 1 ml	2 000 µg/ml	Ref : F113061	Price : CF

ACs Surrogate Standard Solution - 3 components

4-Bromochlorobenzene	CAS:106-39-8	1,4-Difluorobenzene	CAS:540-36-3	Fluorobenzene	CAS:462-06-6
Solvent : Methanol purge & trap		ampoule 5 ml	2 000 µg/ml	Ref : F116331	Price : FA

ACs Surrogate Standard Solution - 5 components

α , α , α -Trifluorotoluene	CAS:98-08-8	4-Bromochlorobenzene	CAS:106-39-8	Fluorobenzene	CAS:462-06-6
2-Bromofluorobenzene	CAS:1072-85-1	1,4-Difluorobenzene	CAS:540-36-3		
Solvent : Methanol purge & trap		ampoule 5 ml	2 000 µg/ml	Ref : F116341	Price : FF

Ethers Performance Check Solution - 1 component

Methyl-tert.butylether	CAS:1634-04-4				
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F116311	Price : CF

EPA method 8021B

Aromatic and halogenated volatiles by GC using photoionization and/or electrolytic conductivity detectors

Method 8021 is used to determine volatile organic compounds in a variety of solid waste matrices. This method is applicable to nearly all types of samples, regardless of water content, including ground water, aqueous sludges, caustic liquors, acid liquors, waste solvents, oily wastes, mousses, tars, fibrous wastes, polymeric emulsions, filter cakes, spent carbons, spent catalysts, soils, and sediments.

VOC Internal Standard Solution - 2 components

2-Bromo-1-chloropropane	CAS:3017-95-6	Fluorobenzene	CAS:462-06-6		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112761	Price : DF

VOC Gas Standard Solution - 1 component

2-Chloroethyl-vinylether	CAS:110-75-8				
Solvent : Methanol purge & trap		ampoule 1 ml	5 000 µg/ml	Ref : F116221	Price : CE

VOC Gas Standard Solution - 1 component

Chloroprene	CAS:126-99-8				
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F116231	Price : BG

VOC Gas Standard Solution - 6 components

Bromomethane	CAS:74-83-9	Chloromethane	CAS:74-87-3	Trichlorofluoromethane	CAS:75-69-4
Chloroethane	CAS:75-00-3	Dichlorodifluoromethane	CAS:75-71-8	Vinylchloride	CAS:75-01-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112901	Price : FA
		ampoule 1 ml	2 000 µg/ml	Ref : F112911	Price : HF

VOC Standard Solution - 54 components

Benzene	CAS:71-43-2	1,3-Dichlorobenzene	CAS:541-73-1	Naphthalene	CAS:91-20-3
Bromobenzene	CAS:108-86-1	1,4-Dichlorobenzene	CAS:106-46-7	n-Propylbenzene	CAS:103-65-1
Bromochloromethane	CAS:74-97-5	1,1-Dichloroethane	CAS:75-34-3	Styrene	CAS:100-42-5
Bromodichloromethane	CAS:75-27-4	1,2-Dichloroethane	CAS:107-06-2	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromoform	CAS:75-25-2	1,1-Dichloroethene	CAS:75-35-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
n-Butylbenzene	CAS:104-51-8	cis-1,2-Dichloroethene	CAS:156-59-2	Tetrachloroethene	CAS:127-18-4
sec-Butylbenzene	CAS:135-98-8	trans-1,2-Dichloroethene	CAS:156-60-5	Toluene	CAS:108-88-3
tert-Butylbenzene	CAS:98-06-6	1,2-Dichloropropane	CAS:78-87-5	1,2,3-Trichlorobenzene	CAS:87-61-6
Carbon tetrachloride	CAS:56-23-5	1,3-Dichloropropane	CAS:142-28-9	1,2,4-Trichlorobenzene	CAS:120-82-1
Chlorobenzene	CAS:108-90-7	2,2-Dichloropropane	CAS:594-20-7	1,1,1-Trichloroethane	CAS:71-55-6
Chloroform	CAS:67-66-3	1,1-Dichloropropene	CAS:563-58-6	1,1,2-Trichloroethane	CAS:79-00-5
2-Chlorotoluene	CAS:95-49-8	cis-1,3-Dichloropropene	CAS:10061-01-5	Trichloroethene	CAS:79-01-6
4-Chlorotoluene	CAS:106-43-4	trans-1,3-Dichloropropene	CAS:10061-02-6	1,2,3-Trichloropropane	CAS:96-18-4
Dibromochloromethane	CAS:124-48-1	Ethylbenzene	CAS:100-41-4	1,2,4-Trimethylbenzene	CAS:95-63-6
1,2-Dibromo-3-chloropropane	CAS:96-12-8	Hexachlorobutadiene	CAS:87-68-3	1,3,5-Trimethylbenzene	CAS:108-67-8
1,2-Dibromoethane	CAS:106-93-4	Isopropylbenzene	CAS:98-82-8	o-Xylene	CAS:95-47-6
Dibromomethane	CAS:74-95-3	4-Isopropyltoluene	CAS:99-87-6	m-Xylene	CAS:108-38-3
1,2-Dichlorobenzene	CAS:95-50-1	Methylene chloride	CAS:75-09-2	p-Xylene	CAS:106-42-3
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112781	Price : BDA

VOC Standard Solution - 60 components

Benzene	CAS:71-43-2	1,2-Dichlorobenzene	CAS:95-50-1	Naphthalene	CAS:91-20-3
Bromobenzene	CAS:108-86-1	1,3-Dichlorobenzene	CAS:541-73-1	n-Propylbenzene	CAS:103-65-1
Bromochloromethane	CAS:74-97-5	1,4-Dichlorobenzene	CAS:106-46-7	Styrene	CAS:100-42-5
Bromodichloromethane	CAS:75-27-4	Dichlorodifluoromethane	CAS:75-71-8	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromoform	CAS:75-25-2	1,1-Dichloroethane	CAS:75-34-3	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Bromomethane	CAS:74-83-9	1,2-Dichloroethane	CAS:107-06-2	Tetrachloroethene	CAS:127-18-4
n-Butylbenzene	CAS:104-51-8	1,1-Dichloroethene	CAS:75-35-4	Toluene	CAS:108-88-3
sec-Butylbenzene	CAS:135-98-8	cis-1,2-Dichloroethene	CAS:156-59-2	1,2,3-Trichlorobenzene	CAS:87-61-6
tert-Butylbenzene	CAS:98-06-6	trans-1,2-Dichloroethene	CAS:156-60-5	1,2,4-Trichlorobenzene	CAS:120-82-1
Carbon tetrachloride	CAS:56-23-5	1,2-Dichloropropane	CAS:78-87-5	1,1,1-Trichloroethane	CAS:71-55-6
Chlorobenzene	CAS:108-90-7	1,3-Dichloropropane	CAS:142-28-9	1,1,2-Trichloroethane	CAS:79-00-5
Chloroethane	CAS:75-00-3	2,2-Dichloropropane	CAS:594-20-7	Trichloroethene	CAS:79-01-6
Chloroform	CAS:67-66-3	1,1-Dichloropropene	CAS:563-58-6	Trichlorofluoromethane	CAS:75-69-4
Chloromethane	CAS:74-87-3	cis-1,3-Dichloropropene	CAS:10061-01-5	1,2,3-Trichloropropane	CAS:96-18-4
2-Chlorotoluene	CAS:95-49-8	trans-1,3-Dichloropropene	CAS:10061-02-6	1,2,4-Trimethylbenzene	CAS:95-63-6
4-Chlorotoluene	CAS:106-43-4	Ethylbenzene	CAS:100-41-4	1,3,5-Trimethylbenzene	CAS:108-67-8
Dibromochloromethane	CAS:124-48-1	Hexachlorobutadiene	CAS:87-68-3	Vinylchloride	CAS:75-01-4
1,2-Dibromoethane	CAS:106-93-4	Isopropylbenzene	CAS:98-82-8	m-Xylene	CAS:108-38-3
1,2-Dibromo-3-chloropropane	CAS:96-12-8	4-Isopropyltoluene	CAS:99-87-6	o-Xylene	CAS:95-47-6
Dibromomethane	CAS:74-95-3	Methylene chloride	CAS:75-09-2	p-Xylene	CAS:106-42-3
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112751	Price : BGA

VOC Surrogate Standard Solution - 2 components

4-Bromochlorobenzene	CAS:106-39-8	1,4-Dichlorobutane	CAS:110-56-5		
Solvent : Methanol purge & trap		ampoule 1 ml	1 500 µg/ml	Ref : F116361	Price : CF

EPA method 8030A**Acrolein and acrylonitrile by GC**

Method 8030 provides gas chromatographic conditions for the detection of the target analytes. Samples can be analyzed using direct injection or purge and-trap (Method 5030). Tenax should be used as the trap packing material. Ground water samples must be analyzed using Method 5030. A temperature program is used in the gas chromatograph to separate the organic compounds. Detection is achieved by a flame ionization detector (FID).

Carbonyl Compounds Standard Solution - 2 components

Acrolein (2-Propenal)	CAS:107-02-8	Acrylonitrile	CAS:107-13-1		
Solvent : Water		ampoule 1 ml	100 µg/ml	Ref : F116441	Price : CF
		ampoule 1 ml	1 000 µg/ml	Ref : F114951	Price : CF
		ampoule 1 ml	2 000 µg/ml	Ref : F116421	Price : CF
		ampoule 1 ml	10 000 µg/ml	Ref : F114961	Price : CF

EPA method 8031**Acrylonitrile**

Method 8031 is applicable for the determination of Acrylonitrile by gas chromatography.

Carbonyl Compound Standard Solution - 1 component

Acrylonitrile	CAS:107-13-1				
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F116461	Price : BK

EPA method 8032

Acrylamide

Method 8032 is applicable for the determination of trace amounts of Acrylamide monomer in aqueous by gas chromatography with electron capture detection.

Phthalates Internal Standard Solution - 1 component

Dimethyl phthalate	CAS:131-11-3				
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F116481	Price : BH

Standard Solution - 1 component

Acrylamide	CAS:79-06-1				
Solvent : Water		ampoule 1 ml	1 000 µg/ml	Ref : F116471	Price : BK

EPA method 8033

Acetonitrile

Method 8033 is applicable for the determination and concentration of Acetonitrile in aqueous by gas chromatography with nitrogen-phosphorous detection.

VOC Standard Solution - 1 component

Acetonitrile	CAS:75-05-8				
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F116491	Price : BG

EPA method 8041

Phenols by GC

Method 8041 describes open-tubular, capillary column gas chromatography procedures for the analysis of phenols, using both single-column and dual-column/dual-detector approaches.

ACs Internal Standard Solution - 2 components

2,5-Dibromotoluene	CAS:615-59-8	2,2',5,5'-Tetrabromobiphenyl	CAS:59080-37-4		
Solvent : 2-Propanol		ampoule 1 ml	1 000 µg/ml	Ref : F116551	Price : FD

Phenols Standard Solution - 1 component

2,4-Dibromophenol	CAS:615-58-7				
Solvent : 2-Propanol		ampoule 1 ml	1 000 µg/ml	Ref : F116531	Price : BJ

Phenols Standard Solution - 1 component

Solvent : Methanol purge & trap

2,4,6-Tribromophenol	CAS:118-79-6	ampoule 1 ml	2 000 µg/ml	Ref : F114921	Price : CA
2-Fluorophenol	CAS:367-12-4	ampoule 1 ml	2 000 µg/ml	Ref : F116591	Price : CC
Pentafluorophenol	CAS:771-61-9	ampoule 1 ml	2 000 µg/ml	Ref : F116601	Price : CA
Phenol D5	CAS:4165-62-2	ampoule 1 ml	2 000 µg/ml	Ref : F116611	Price : CG

Phenols Standard Solution - 3 components

2-Cyclohexyl-4,6-dinitrophenol	CAS:131-89-5	2,3,4,5-Tetrachlorophenol	CAS:4901-51-3	2,3,5,6-Tetrachlorophenol	CAS:935-95-5
Solvent : 2-Propanol		ampoule 1 ml	2 000 µg/ml	Ref : F116511	Price : HJ

Phenols Standard Solution - 9 components

4-Chloro-3-methylphenol	CAS:59-50-7	2-Methyl-4,6-dinitrophenol	CAS:534-52-1	Pentachlorophenol	CAS:87-86-5
o-Cresol	CAS:95-48-7	2-Nitrophenol	CAS:88-75-5	Phenol	CAS:108-95-2
2,4-Dichlorophenol	CAS:120-83-2	4-Nitrophenol	CAS:100-02-7	2,4,6-Trichlorophenol	CAS:88-06-2
Solvent : 2-Propanol		ampoule 1 ml	2 000 µg/ml	Ref : F116501	Price : EF

Phenols Standard Solution - 9 components

2-Chlorophenol	CAS:95-57-8	2,6-Dichlorophenol	CAS:87-65-0	Dinoseb	CAS:88-85-7
m-Cresol	CAS:108-39-4	2,4-Dimethylphenol	CAS:105-67-9	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
p-Cresol	CAS:106-44-5	2,4-Dinitrophenol	CAS:51-28-5	2,4,5-Trichlorophenol	CAS:95-95-4
Solvent : 2-Propanol		ampoule 1 ml	2 000 µg/ml	Ref : F116521	Price : FA

Phenols Standard Solution - 17 components

4-Chloro-3-methylphenol	CAS:59-50-7	2,6-Dichlorophenol	CAS:87-65-0	Pentachlorophenol	CAS:87-86-5
2-Chlorophenol	CAS:95-57-8	2,4-Dimethylphenol	CAS:105-67-9	Phenol	CAS:108-95-2
o-Cresol	CAS:95-48-7	2,4-Dinitrophenol	CAS:51-28-5	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
m-Cresol	CAS:108-39-4	2-Methyl-4,6-dinitrophenol	CAS:534-52-1	2,4,5-Trichlorophenol	CAS:95-95-4
p-Cresol	CAS:106-44-5	2-Nitrophenol	CAS:88-75-5	2,4,6-Trichlorophenol	CAS:88-06-2
2,4-Dichlorophenol	CAS:120-83-2	4-Nitrophenol	CAS:100-02-7		
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F116571	Price : FJ

Phenols Standard Solution - 21 components

4-Chloro-3-methylphenol	CAS:59-50-7	2,6-Dichlorophenol	CAS:87-65-0	Pentachlorophenol	CAS:87-86-5
2-Chlorophenol	CAS:95-57-8	2,4-Dimethylphenol	CAS:105-67-9	Phenol	CAS:108-95-2
o-Cresol	CAS:95-48-7	2,4-Dinitrophenol	CAS:51-28-5	2,3,4,5-Tetrachlorophenol	CAS:4901-51-3
m-Cresol	CAS:108-39-4	Dinoseb	CAS:88-85-7	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
p-Cresol	CAS:106-44-5	2-Methyl-4,6-dinitrophenol	CAS:534-52-1	2,3,5,6-Tetrachlorophenol	CAS:935-95-5
2-Cyclohexyl-4,6-dinitrophenol	CAS:131-89-5	2-Nitrophenol	CAS:88-75-5	2,4,6-Trichlorophenol	CAS:88-06-2
2,4-Dichlorophenol	CAS:120-83-2	4-Nitrophenol	CAS:100-02-7	2,4,5-Trichlorophenol	CAS:95-95-4
Solvent : 2-Propanol		ampoule 1 ml	100 µg/ml	Ref : F116561	Price : HF

Phenols Surrogate Standard Mixture - 2 components

2-Fluorophenol	CAS:367-12-4	2,4,6-Tribromophenol	CAS:118-79-6		
Solvent : 2-Propanol		ampoule 1 ml	2 000 µg/ml	Ref : F116541	Price : DA

EPA method 8061A**Phthalate esters by Gas chromatography with electron capture detection (GC/ECD)**

Method 8061 is used to determine the identities and concentrations of various phthalate esters in aqueous and solid matrices including groundwater, leachate, soil, sludge and sediment.

ACs Internal Standard Solution - 1 component

Benzyl benzoate	CAS:120-51-4				
Solvent : n-Hexane		ampoule 1 ml	5 000 µg/ml	Ref : F116661	Price : CE

Phthalates Matrix Spike Solution - 2 components

Bis(2-ethylhexyl)phthalate	CAS:117-81-7	Butyl benzyl phthalate	CAS:85-68-7		
Solvent : Acetone		ampoule 1 ml	2 000 µg/ml	Ref : F116651	Price : CF

Phthalates Standard Solution - 6 components

di-n-butylphthalate	CAS:84-74-2	Dimethyl phthalate	CAS:131-11-3	Phthalic acid,bis-2-ethylhexylester	CAS:117-81-7
Diethyl phthalate	CAS:84-66-2	di-n-octylphthalate	CAS:117-84-0	Phthalic acid, benzylbutyl ester	CAS:85-68-7
Solvent : n-Hexane		ampoule 1 ml	100 µg/ml	Ref : F116631	Price : EF
		ampoule 1 ml	1 000 µg/ml	Ref : F116621	Price : EF

Phthalates Surrogate Standard Solution - 3 components

Dibenzyl phthalate	CAS:523-31-9	Diphenyl isophthalate	id : 16211	Diphenyl phthalate	CAS:84-62-8
Solvent : Acetone		ampoule 1 ml	50 µg/ml	Ref : F116681	Price : DA
		ampoule 1 ml	500 µg/ml	Ref : F116671	Price : DA

EPA method 8070**Nitrosoamines**

Method 8070 is applicable for the determination of nitrosoamines by gas chromatography (GC/NPD).

Standard Solution - 3 components

n-Nitroso-dimethylamine	CAS:62-75-9	n-Nitroso-diphenylamine	CAS:86-30-6	n-Nitroso-di-n-propylamine	CAS:621-64-7
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F116691	Price : DH

EPA method 8080A**Organochlorine pesticides and polychlorinated biphenyls by GC**

Method 8080 provides gas chromatographic conditions (GC-ECD or GC-HECD) for the detection of ppb concentrations of certain organochlorine pesticides and PCBs.

OCs Decomposition Solution - 2 components

p,p'-DDT	CAS:50-29-3	2 µg/ml	Endrin	CAS:72-20-8	1 µg/ml
Solvent : Acetone		ampoule 1 ml		Ref : F116821	Price : CF

OCS Standard Solution - 1 component

				Solvent : n-Hexane	
Aroclor 1016 (PCB 1016)	CAS:12674-11-2	ampoule 1 ml	100 µg/ml	Ref : F116741	Price : BH
Aroclor 1221 (PCB 1221)	CAS:11104-28-2	ampoule 1 ml	100 µg/ml	Ref : F116751	Price : BG
Aroclor 1232 (PCB 1232)	CAS:11141-16-5	ampoule 1 ml	100 µg/ml	Ref : F116761	Price : CC
Aroclor 1242	CAS:53469-21-9	ampoule 1 ml	100 µg/ml	Ref : F116771	Price : BH
Aroclor 1248	CAS:12672-29-6	ampoule 1 ml	100 µg/ml	Ref : F116781	Price : BJ
Aroclor 1254	CAS:11097-69-1	ampoule 1 ml	100 µg/ml	Ref : F116791	Price : BJ
Aroclor 1260	CAS:11096-82-5	ampoule 1 ml	100 µg/ml	Ref : F116801	Price : BJ
Camphechlor	CAS:8001-35-2	ampoule 1 ml	100 µg/ml	Ref : F116731	Price : BG
Chlordane (technical)	CAS:57-74-9	ampoule 1 ml	100 µg/ml	Ref : F116721	Price : BG

OCs Standard Solution - 3 components

2,4'-DDD	CAS:53-19-0	2,4'-DDE	CAS:3424-82-6	2,4'-DDT	CAS:789-02-6
Solvent : Iso-Octane		ampoule 1 ml	250 µg/ml	Ref : F116811	Price : DF

OCs Standard Solution - 17 components

Aldrin	CAS:309-00-2	Endosulfan II	CAS:33213-65-9	gamma -HCH (Lindane)	CAS:58-89-9
4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	delta -HCH	CAS:319-86-8
4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8
4,4'-DDT	CAS:50-29-3	Endrin aldehyde	CAS:7421-93-4	Heptachlor epoxide	CAS:1024-57-3
Dieldrin	CAS:60-57-1	alpha -HCH	CAS:319-84-6	Methoxychlor (DMTD)	CAS:72-43-5
Endosulfan I	CAS:959-98-8	beta -HCH	CAS:319-85-7		
Solvent : Acetone		ampoule 1 ml	2 000 µg/ml	Ref : F116701	Price : BEA

OCs QC Standard Solution - 17 components

Aldrin	CAS:309-00-2	0.02 µg/ml	Endrin aldehyde	CAS:7421-93-4	0.02 µg/ml
4,4'-DDD (TDE)	CAS:72-54-8	0.1 µg/ml	alpha -HCH	CAS:319-84-6	0.02 µg/ml
4,4'-DDE	CAS:72-55-9	0.02 µg/ml	beta -HCH	CAS:319-85-7	0.02 µg/ml
4,4'-DDT	CAS:50-29-3	0.1 µg/ml	gamma -HCH (Lindane)	CAS:58-89-9	0.02 µg/ml
Dieldrin	CAS:60-57-1	0.02 µg/ml	delta -HCH	CAS:319-86-8	0.02 µg/ml
Endosulfan I	CAS:959-98-8	0.02 µg/ml	Heptachlor	CAS:76-44-8	0.02 µg/ml
Endosulfan II	CAS:33213-65-9	0.1 µg/ml	Heptachlor epoxide	CAS:1024-57-3	0.02 µg/ml
Endosulfan-total (sulfate)	CAS:1031-07-8	0.1 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	0.02 µg/ml
Endrin	CAS:72-20-8	0.1 µg/ml			
Solvent : Acetone		ampoule 1 ml		Ref : F116711	Price : GF

Surrogate Standard Solution - 2 components

Decachlorobiphenyl (PCB 209)	CAS:2051-24-3	2,4,5,6-Tetrachloro-m-xylene	CAS:877-09-8		
Solvent : Acetone		ampoule 1 ml	200 µg/ml	Ref : F116831	Price : DA

EPA method 8081A**Organochlorine pesticides by GC**

Method 8081 is used to determine the concentrations of various organochlorine pesticides in extracts from solid and liquid matrices, using fused-silica, open-tubular, capillary columns with electron capture detectors (ECD).

ACs Internal Standard Solution - 1 component

2-bromo-nitrobenzene	CAS:577-19-5				
Solvent : Acetone		ampoule 1 ml	5 000 µg/ml	Ref : F117261	Price : CE

OCs Internal Standard Solution - 1 component

Pentachloronitrobenzene	CAS:82-68-8				
Solvent : Acetone		ampoule 1 ml	5 000 µg/ml	Ref : F117251	Price : CF

OCS Standard Solution - 1 component

Solvent : n-Hexane

Chlordane (technical)	CAS:57-74-9	ampoule 1 ml	100 µg/ml	Ref : F116721	Price : BG
Toxaphene	CAS:8001-35-2	ampoule 1 ml	100 µg/ml	Ref : F116731	Price : BG

OCs Standard Solution - 6 components

Chlorobenzilate	CAS:510-15-6	1,2-Dibromo-3-chloropropane	CAS:96-12-8	Hexachlorocyclopentadiene	CAS:77-47-4
Diallate	CAS:2303-16-4	Hexachlorobenzene	CAS:118-74-1	Isodrin	CAS:465-73-6
Solvent : n-Hexane/Toluene (1/1)		ampoule 1 ml	1 000 µg/ml	Ref : F117231	Price : EF

OCs Standard Solution - 20 components

Aldrin	CAS:309-00-2	4,4'-DDD (TDE)	CAS:72-54-8	Endrin	CAS:72-20-8
alpha-BHC	CAS:319-84-6	4,4'-DDE	CAS:72-55-9	Endrin aldehyde	CAS:7421-93-4
beta-BHC	CAS:319-85-7	4,4'-DDT	CAS:50-29-3	Endrin ketone	CAS:53494-70-5
gamma-BHC	CAS:58-89-9	Dieldrin	CAS:60-57-1	Heptachlor	CAS:76-44-8
delta-BHC	CAS:319-86-8	Endosulfan-alpha	CAS:959-98-8	Methoxychlor (DMTD)	CAS:72-43-5
alpha-Chlordane	CAS:5103-71-9	Endosulfan-beta	CAS:33213-65-9	cis-Heptachlorepoxyde	CAS:1024-57-3
gamma-Chlordane	CAS:5103-74-2	Endosulfan-total (sulfate)	CAS:1031-07-8	(cis-, exo-, isomer B)	
Solvent : n-Hexane/Toluene (1/1)		ampoule 1 ml	1 000 µg/ml	Ref : F117221	Price : BFA

Pesticides Surrogate Standard Spiking Solution - 2 components

PCB 209	CAS:2051-24-3	2,4,5,6-Tetrachloro-m-xylene	CAS:877-09-8		
Solvent : Acetone		ampoule 1 ml	200 µg/ml	Ref : F117241	Price : DA

EPA method 8082

Polychlorinated biphenyls (PCBs) by GC

Method 8082 is used to determine the concentrations of polychlorinated biphenyls (PCBs) as Aroclors or as individual PCB congeners in extracts from solid and aqueous matrices. Open tubular, capillary columns are employed with electron capture detectors (ECD) or electrolytic conductivity detectors (ELCD).

ACs Surrogate Standard Solution - 1 component

Tetrachloro-m-xylene	CAS:877-09-8				
Solvent : Acetone		ampoule 1 ml	1 000 µg/ml	Ref : F117321	Price : BK

Internal & Surrogate Solution - 2 components

Decachlorobiphenyl (PCB 209)	CAS:2051-24-3	Tetrachloro-m-xylene	CAS : 877-09-8		
Solvent : Acetone		ampoule 1 ml	1 000 µg/ml	Ref : F117331	Price : DA

PCBs Standard Solution - 1 component

Solvent : Iso-Octane

Aroclor 1016 (PCB 1016)	CAS:12674-11-2	ampoule 1 ml	100 µg/ml	Ref : F115081	Price : BH
Aroclor 1221 (PCB 1221)	CAS:11104-28-2	ampoule 1 ml	100 µg/ml	Ref : F115091	Price : BG
Aroclor 1232 (PCB 1232)	CAS:11141-16-5	ampoule 1 ml	100 µg/ml	Ref : F115121	Price : CC
Aroclor 1242	CAS:53469-21-9	ampoule 1 ml	100 µg/ml	Ref : F115131	Price : BH
Aroclor 1248 (PCB 1248)	CAS:12672-29-6	ampoule 1 ml	100 µg/ml	Ref : F115141	Price : BJ
Aroclor 1254	CAS:11097-69-1	ampoule 1 ml	100 µg/ml	Ref : F115151	Price : BJ
Aroclor 1260	CAS:11096-82-5	ampoule 1 ml	100 µg/ml	Ref : F115161	Price : BJ

PCBs Standard Solution - 2 components

Aroclor 1016 (PCB 1016)	CAS:12674-11-2	Aroclor 1260 (PCB 1260)	CAS:11096-82-5		
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F117301	Price : DF

PCBs Standard Solution - 19 components

PCB 1	CAS:2051-60-7	PCB 87	CAS:38380-02-8	PCB 170	CAS:35065-30-6
PCB 5	CAS:16605-91-7	PCB 101	CAS:37680-73-2	PCB 180	CAS:35065-29-3
PCB 18	CAS:37680-65-2	PCB 110	CAS:38380-03-9	PCB 183	CAS:52663-69-1
PCB 31	CAS:16606-02-3	PCB 138	CAS:35065-28-2	PCB 187	CAS:52663-68-0
PCB 44	CAS:41464-39-5	PCB 141	CAS:52712-04-6	PCB 206	CAS:40186-72-9
PCB 52	CAS:35693-99-3	PCB 151	CAS:52663-63-5		
PCB 66	CAS:32598-10-0	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 1 ml	100 µg/ml	Ref : F117291	Price : BDA

PCBs Surrogate Standard Solution - 1 component

Decachlorobiphenyl (PCB 209)	CAS:2051-24-3				
Solvent : Acetone		ampoule 1 ml	1 000 µg/ml	Ref : F117311	Price : BK

EPA method 8091

Nitroaromatics and cyclic ketones by GC

Method 8091 is a gas chromatographic (GC) method used to determine the concentration of nitroaromatics and cyclic ketones. It describes wide-bore, open-tubular, capillary column gas chromatography procedures using either electron capture (ECD) or nitrogen-phosphorous (NPD) detectors.

ACs Internal Standard Solution - 1 component

Hexachlorobenzene	CAS:118-74-1				
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F117491	Price : BK

ACs Surrogate Standard Solution - 1 component

1-Chloro-3-nitrobenzene	CAS:121-73-3				
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F117501	Price : BK

RCRA Analytes - 6 components

1,4-Dinitrobenzene	CAS:100-25-4	2,6-Dinitrotoluene	CAS:606-20-2	Nitrobenzene	CAS:98-95-3
2,4-Dinitrotoluene	CAS:121-14-2	1,4-Naphthoquinone	CAS:130-15-4	Pentachloronitrobenzene	CAS:82-68-8
Solvent : Iso-Octane/Toluene (1/1)		ampoule 1 ml	1 000 µg/ml	Ref : F117371	Price : DJ

Non-RCA Analytes - 30 components

Benefin	CAS:1861-40-1	2,4-Dichloronitrobenzene	CAS:611-06-3	3-Nitrotoluene	CAS:99-08-1
Butralin	CAS:33629-47-9	3,5-Dichloronitrobenzene	CAS:618-62-2	4-Nitrotoluene	CAS:99-99-0
1-Chloro-2,4-dinitrobenzene	CAS:97-00-7	3,4-Dichloronitrobenzene	CAS:99-54-7	Pendimethalin	CAS:40487-42-1
1-Chloro-3,4-dinitrobenzene	CAS:610-40-2	2,5-Dichloronitrobenzene	CAS:89-61-2	Profluralin	CAS:26399-36-0
1-Chloro-2-nitrobenzene	CAS:88-73-3	Dinitramine	CAS:29091-05-2	2,3,4,5-Tetrachloronitrobenzene	CAS:879-39-0
1-Chloro-4-nitrobenzene	CAS:100-00-5	1,2-Dinitrobenzene	CAS:528-29-0	2,3,5,6-Tetrachloronitrobenzene	CAS:117-18-0
2-Chloro-6-nitrotoluene	CAS:83-42-1	1,3-Dinitrobenzene	CAS:99-65-0	2,3,4-Trichloronitrobenzene	CAS:17700-09-3
4-Chloro-2-nitrotoluene	CAS:89-59-8	Isopropalin	CAS:33820-53-0	2,4,5-Trichloronitrobenzene	CAS:89-69-0
4-Chloro-3-nitrotoluene	CAS:89-60-1	1,2-Naphthoquinone	CAS:524-42-5	1,3,5-Trichloro-2-nitrobenzene	CAS:18708-70-8
2,3-Dichloronitrobenzene	CAS:3209-22-1	2-Nitrotoluene	CAS:88-72-2	Trifluralin	CAS:1582-09-8
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F117381	Price : JJ

EPA method 8100

Polynuclear aromatic hydrocarbons

Method 8100 provides gas chromatographic conditions (GC/FID) for the detection of ppb levels of certain polynuclear aromatic hydrocarbons.

PAH Standard Solution - 8 components

Benzo(j)fluoranthene	CAS:205-82-3	7-H-Dibenzo(c,g)carbazole	CAS:194-59-2	Dibenzo(a,i)pyrene	CAS:189-55-9
Dibenz(a,h)acridine	CAS:226-36-8	Dibenzo(a,e)pyrene	CAS:192-65-4	3-Methylcholanthrene	CAS:56-49-5
Dibenz(a,j)acridine	CAS:224-42-0	Dibenzo(a,h)pyrene	CAS:189-64-0		
Solvent : Dichloromethane		ampoule 1 ml	1 000 µg/ml	Ref : F117521	Price : CKA

PAH Standard Solution - 16 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Acenaphthylene	CAS:208-96-8	Benzo(a)pyrene	CAS:50-32-8	Naphthalene	CAS:91-20-3
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Phenanthrene	CAS:85-01-8
Benzo(a)anthracene	CAS:56-55-3	Dibenzo(a,h)anthracene	CAS:53-70-3	Pyrene	CAS:129-00-0
Benzo(b)fluoranthene	CAS:205-99-2	Fluoranthene	CAS:206-44-0		
Benzo(k)fluoranthene	CAS:207-08-9	Fluorene	CAS:86-73-7		
Solvent : Dichloromethane		ampoule 1 ml	1 000 µg/ml	Ref : F117511	Price : BDF

PAH QC Mixture - 16 components

Acenaphthene	CAS:83-32-9	100 µg/ml	Chrysene	CAS:218-01-9	10 µg/ml
Acenaphthylene	CAS:208-96-8	100 µg/ml	Dibenzo(a,h)anthracene	CAS:53-70-3	10 µg/ml
Anthracene	CAS:120-12-7	100 µg/ml	Fluoranthene	CAS:206-44-0	10 µg/ml
Benzo(a)anthracene	CAS:56-55-3	10 µg/ml	Fluorene	CAS:86-73-7	100 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	10 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	10 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	5 µg/ml	Naphthalene	CAS:91-20-3	100 µg/ml
Benzo(g,h,i)perylene	CAS:191-24-2	10 µg/ml	Phenanthrene	CAS:85-01-8	100 µg/ml
Benzo(a)pyrene	CAS:50-32-8	10 µg/ml	Pyrene	CAS:129-00-0	10 µg/ml
Solvent : Acetonitrile		ampoule 1 ml		Ref : F117531	Price : GF

Surrogate Standard Solution - 2 components

2-Fluorobiphenyl	CAS:321-60-8	1-Fluoronaphthalene	CAS:321-38-0		
Solvent : Dichloromethane		ampoule 1 ml	2 000 µg/ml	Ref : F117541	Price : DA

EPA method 8111

Haloethers by GC

Method 8111 is a gas chromatographic (GC) method used to determine the concentration of haloethers. It describes wide-bore open-tubular, capillary column gas chromatography procedures using a dual-column/dual-detector approach, however, a single column/single detector approach is acceptable.

Ethers Standard Solution - 4 components

Bis-(2-chloroethoxy)-methane	CAS:111-91-1	Bis(2-chloroisopropyl) ether	CAS:108-60-1	4-Chlorophenyl phenyl ether	CAS:7005-72-3
Bis-(2-chloroethyl)-ether	CAS:111-44-4				
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F117551	Price : DF

Internal Standard Solution - 1 component

4-bromophenyl phenyl ether	CAS:92-86-4				
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F117561	Price : CC

EPA method 8121

Chlorinated hydrocarbons by GC: capillary column technique

Method 8121 describes the determination of chlorinated hydrocarbons in extracts prepared from environmental samples and RCRA wastes. It describes wide-bore open-tubular, capillary column gas chromatography procedures using both single column/single detector and dual-column/dual-detector approaches.

ACs Internal Standard Solution - 1 component

1,3,5-Tribromobenzene	CAS:626-39-1				
Solvent : n-Hexane		ampoule 1 ml	1 000 µg/ml	Ref : F117581	Price : BK

ACs Internal Standard Solution - 1 component

1,3,5-Tribromobenzene	CAS:626-39-1				
Solvent : Acetone		ampoule 1 ml	50 µg/ml	Ref : F117591	Price : BK

ACs Surrogate Standard Solution - 3 components

1,4-Dichloronaphthalene	CAS:1825-31-6	10 µg/ml	alpha,2,6-Trichlorotoluene	CAS:2014-83-7	1 µg/ml
2,3,4,5,6-Pentachlorotoluene	CAS:877-11-2	1 µg/ml			
Solvent : Acetone		ampoule 1 ml		Ref : F117611	Price : CF

ACs / VOC Standard Solution - 22 components

2-Chloronaphthalene	CAS:91-58-7	2000 µg/ml	1,2,3,4-Tetrachlorobenzene	CAS:634-66-2	100 µg/ml
1,2-Dichlorobenzene	CAS:95-50-1	1000 µg/ml	1,2,4,5-Tetrachlorobenzene	CAS:95-94-3	100 µg/ml
1,3-Dichlorobenzene	CAS:541-73-1	1000 µg/ml	1,2,3,5-Tetrachlorobenzene	CAS:634-90-2	100 µg/ml
1,4-Dichlorobenzene	CAS:106-46-7	1000 µg/ml	1,2,3-Trichlorobenzene	CAS:87-61-6	100 µg/ml
Benzal chloride	CAS:98-87-3	100 µg/ml	1,2,4-Trichlorobenzene	CAS:120-82-1	100 µg/ml
Benzotrichlorid	CAS:98-07-7	100 µg/ml	1,3,5-Trichlorobenzene	CAS:108-70-3	100 µg/ml
Benzyl chloride	CAS:100-44-7	100 µg/ml	Hexachlorobenzene	CAS:118-74-1	10 µg/ml
alpha-BHC	CAS:319-84-6	100 µg/ml	Hexachlorobutadiene	CAS:87-68-3	10 µg/ml
beta-BHC	CAS:319-85-7	100 µg/ml	Hexachlorocyclopentadiene	CAS:77-47-4	10 µg/ml
gamma-BHC	CAS:58-89-9	100 µg/ml	Hexachloroethane	CAS:67-72-1	10 µg/ml
delta-BHC	CAS:319-86-8	100 µg/ml	Pentachlorobenzene	CAS:608-93-5	10 µg/ml
Solvent : n-Hexane		ampoule 1 ml		Ref : F117621	Price : GH

OCs / ACs Standard Solution - 22 components

Benzal chloride	CAS:98-87-3	1,2-Dichlorobenzene	CAS:95-50-1	1,2,3,4-Tetrachlorobenzene	CAS:634-66-2
Benzotrichlorid	CAS:98-07-7	1,3-Dichlorobenzene	CAS:541-73-1	1,2,3,5-Tetrachlorobenzene	CAS:634-90-2
Benzyl chloride	CAS:100-44-7	1,4-Dichlorobenzene	CAS:106-46-7	1,2,4,5-Tetrachlorobenzene	CAS:95-94-3
alpha-BHC	CAS:319-84-6	Hexachlorobenzene	CAS:118-74-1	1,2,3-Trichlorobenzene	CAS:87-61-6
beta-BHC	CAS:319-85-7	Hexachlorobutadiene	CAS:87-68-3	1,2,4-Trichlorobenzene	CAS:120-82-1
gamma-BHC	CAS:58-89-9	Hexachlorocyclopentadiene	CAS:77-47-4	1,3,5-Trichlorobenzene	CAS:108-70-3
delta-BHC	CAS:319-86-8	Hexachloroethane	CAS:67-72-1		
2-Chloronaphthalene	CAS:91-58-7	Pentachlorobenzene	CAS:608-93-5		
Solvent : n-Hexane		ampoule 1 ml	1 000 µg/ml	Ref : F117571	Price : GC

EPA method 8131**Aniline and selected derivatives by GC**

Method 8131 is used to determine (by gas chromatography with a specific detector) the concentration of aniline and certain derivatives of aniline in extracts prepared from environmental samples and RCRA wastes.

NCC Standard Solution - 19 components

Aniline	CAS:62-53-3	2-Chloro-4,6-dinitroaniline	CAS:3531-19-9	2-Nitroaniline	CAS:88-74-4
4-Bromoaniline	CAS:106-40-1	2-Chloro-4-nitroaniline	CAS:121-87-9	3-Nitroaniline	CAS:99-09-2
2-Bromo-6-chloro-4-nitroaniline	CAS:99-29-6	4-Chloro-2-nitroaniline	CAS:89-63-4	4-Nitroaniline	CAS:100-01-6
2-Bromo-4,6-dinitroaniline	CAS:1817-73-8	2,6-Dibromo-4-nitroaniline	CAS:827-94-1	2,4,6-Trichloroaniline	CAS:634-93-5
2-Chloroaniline	CAS:95-51-2	Dicloran (CNA)	CAS:99-30-9	2,4,5-Trichloroaniline	CAS:636-30-6
3-Chloroaniline	CAS:108-42-9	3,4-Dichloroaniline	CAS:95-76-1		
4-Chloroaniline	CAS:106-47-8	2,4-Dinitroaniline	CAS:97-02-9		
Solvent : Toluene		ampoule 1 ml	1 000 µg/ml	Ref : F117631	Price : FJ

EPA method 8141A**Organophosphorus compounds by GC: capillary column technique**

Method 8141 is a capillary gas chromatographic (GC-FPD or NPD) method used to determine the concentration of organophosphorus (OP) compounds.

Internal Standard for NPD - 1 component

1-Bromo-2-nitrobenzene	CAS:577-19-5				
Solvent : Acetone		ampoule 1 ml	1 000 µg/ml	Ref : F119061	Price : BK

Industrial Chemicals & Triazine Herbicides Mixture - 2 components

Hexamethylphosphoric triamide (HMPA)	CAS:680-31-9	Tri-o-cresylphosphate (TOCP)	CAS:78-30-8		
Solvent : n-Hexane		ampoule 1 ml	200 µg/ml	Ref : F119071	Price : DA

Industrial Chemicals & Triazine Herbicides Mixture - 2 components

Atrazine	CAS:1912-24-9	Simazine (CAT)	CAS:122-34-9		
Solvent : Acetone		ampoule 1 ml	200 µg/ml	Ref : F119081	Price : DA

OPP Standard Solution - 7 components

Dimethoate	CAS:60-51-5	Monocrotophos	CAS:6923-22-4	TEPP	CAS:107-49-3
EPN	CAS:2104-64-5	Parathion-ethyl	CAS:56-38-2		
Malathion	CAS:121-75-5	Sulfotep	CAS:3689-24-5		
Solvent : n-Hexane/Acetone (9/1)		ampoule 1 ml	200 µg/ml	Ref : F119031	Price : EF

OPP Standard Solution - 9 components

Aspon	CAS:3244-90-4	Dichlofenthion	CAS:97-17-6	Fonofos	CAS:944-22-9
Chlorpyrifos methyl	CAS:5598-13-0	Diclotophos	CAS:141-66-2	Thionazin	CAS:297-97-2
Crotoxypfos	CAS:7700-17-6	Fenitrothion	CAS:122-14-5	Trichlorfon	CAS:52-68-6
Solvent : n-Hexane		ampoule 1 ml	200 µg/ml	Ref : F119051	Price : FA

OPP Standard Solution - 10 components

Azinphos-ethyl	CAS:2642-71-9	Famphur	CAS:52-85-7	Terbufos	CAS:13071-79-9
Chlorfenvinphos	CAS:470-90-6	Leptophos	CAS:21609-90-5	Trithion	CAS:786-19-6
Dioxathion	CAS:78-34-2	Phosmet	CAS:732-11-6		
Ethion	CAS:563-12-2	Phosphamidon	CAS:13171-21-6		
Solvent : n-Hexane/Acetone (9/1)		ampoule 1 ml	200 µg/ml	Ref : F119041	Price : FA

OPP Standard Solution - 20 components

Azinphos-methyl	CAS:86-50-0	Disulfoton	CAS:298-04-4	Phosdrin	CAS:7786-34-7
Bolstar	CAS:35400-43-2	Ethoprophos	CAS:13194-48-4	Phorate	CAS:298-02-2
Chlorpyrifos	CAS:2921-88-2	Fensulfotthion	CAS:115-90-2	Ronnel	CAS:299-84-3
Coumaphos	CAS:56-72-4	Fenthion	CAS:55-38-9	Stirofos	CAS:22248-79-9
Demeton (O+S)	CAS:8065-48-3	Merphos	CAS:150-50-5	Tokuthion	CAS:34643-46-4
Diazinon	CAS:333-41-5	Naled	CAS:300-76-5	Trichloronat	CAS:327-98-0
Dichlorvos	CAS:62-73-7	Parathion-methyl	CAS:298-00-0		
Solvent : n-Hexane/Acetone (9/1)		ampoule 1 ml	200 µg/ml	Ref : F119021	Price : KG

Surrogate Standard for NPD only - 1 component

4-Chloro-3-nitrobenzotrifluoride	CAS:121-17-5				
Solvent : Acetone		ampoule 1 ml	1 000 µg/ml	Ref : F119101	Price : BH

Surrogate Standard for NPD & FPD - 2 components

Tributyl phosphate	CAS:126-73-8	Triphenylphosphate	CAS:115-86-6		
Solvent : Acetone		ampoule 1 ml	1 000 µg/ml	Ref : F119091	Price : DA

EPA method 8151A**Chlorinated herbicides by GC using methylation or pentafluorobenzoylation derivatization**

Method 8151 is a capillary gas chromatographic (GC) method for determining certain chlorinated acid herbicides and related compounds in aqueous, soil and waste matrices.

Internal Standard Solution - 1 component

DBOB	CAS:10386-84-2				
Solvent : Acetone		ampoule 1 ml	250 µg/ml	Ref : F119150	Price : BG

Methylated Chlorinated Herbicides Mixture - 8 components

2,4-D methyl ester	CAS:1928-38-7	Dicamba-methyl ester	CAS:6597-78-0	Silvex methyl ester	CAS:4841-20-7
Dalapon methyl ester	CAS:17640-02-7	Dichlorprop methyl ester	CAS:57153-17-0	2,4,5-T methyl ester	CAS:1928-37-6
2,4-DB methyl ester	CAS:18625-12-2	Dinoseb methyl ether	CAS:6099-79-2		
Solvent : n-Hexane		ampoule 1 ml	20 µg/ml	Ref : F119161	Price : EA

OCs Laboratory Performance Check Solution - 5 components

DBOB	CAS:10386-84-2	250 µg/ml	Dinoseb	CAS:88-85-7	4 µg/ml
3,5-Dichlorobenzoic acid	CAS:51-36-5	600 µg/ml	4-Nitrophenol	CAS:100-02-7	1 600 µg/ml
2,4-Dichlorophenylacetic acid	CAS:19719-28-9	500 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F119181	Price : DF

OCs / Phenoxyacetic Herbicides Standard Solution - 6 components

Acifluorfen	CAS:50594-66-6	Chloramben	CAS:133-90-4	3,5-Dichlorobenzoic acid	CAS:51-36-5
Bentazone	CAS:25057-89-0	D CPA diacid	CAS:2136-79-0	Picloram	CAS:1918-02-1
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F119120	Price : EF

OCs / Phenoxyacetic Herbicides Standard Solution - 12 components

2,4-D	CAS:94-75-7	Dichlorprop	CAS:120-36-5	4-Nitrophenol	CAS:100-02-7
Dalapon	CAS:75-99-0	Dinoseb	CAS:88-85-7	Pentachlorophenol	CAS:87-86-5
2,4-DB	CAS:94-82-6	MCPA	CAS:94-74-6	Silvex	CAS:93-72-1
Dicamba	CAS:1918-00-9	Mecoprop (MCP)	CAS:7085-19-0	2,4,5-T	CAS:93-76-5
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F119110	Price : FA

Phenoxyacetic Herbicides Standard Solution - 3 components

2,4-D methyl ester	CAS:1928-38-7	Silvex methyl ester	CAS:4841-20-7	2,4,5-T methyl ester	CAS:1928-37-6
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F119171	Price : DA

Phenoxyacetic Herbicides Surrogate Solution - 1 component

2,4-Dichlorophenylacetic acid

CAS:19719-28-9

Solvent : Acetone	ampoule 1 ml	100 µg/ml	Ref : F119130	Price : BG
	ampoule 1 ml	2 000 µg/ml	Ref : F119140	Price : CA

EPA method 8240B

Volatile organic compounds by Gas chromatography/Mass spectrometry (GC/MS)

Method 8240 is used to determine volatile organic compounds in a variety of solid waste matrices. This method is applicable to nearly all types of samples, regardless of water content, including ground water, aqueous sludges, caustic liquors, acid liquors, waste solvents, oily wastes, mousses, tars, fibrous wastes, polymeric emulsions, filter cakes, spent carbons, spent catalysts, soils, and sediments.

Internal Standard Solution - 1 component

Ethylene Oxide

CAS:75-21-8

Solvent : Iso-Octane	ampoule 1 ml	200 µg/ml	Ref : F147440	Price : BG
----------------------	--------------	-----------	---------------	------------

Internal Standard Mixture - 3 components

Bromochloromethane	CAS:74-97-5	Chlorobenzene D5	CAS:3114-55-4	1,4-Difluorobenzene	CAS:540-36-3
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F119221	Price : DF

VOC Standard Solution - 15 components

Acrolein (2-Propenal)	CAS:107-02-8	Chloroethane	CAS:75-00-3	Dichlorodifluoromethane	CAS:75-71-8
Acrylonitrile	CAS:107-13-1	2-Chloroethanol	CAS:107-07-3	Trichlorofluoromethane	CAS:75-69-4
Bromomethane	CAS:74-83-9	2-Chloroethyl-vinylether	CAS:110-75-8	1,2,3-Trichloropropane	CAS:96-18-4
Bromoacetone	CAS:598-31-2	Chloromethane	CAS:74-87-3	Vinylacetate	CAS:108-05-4
Chloralhydrate	CAS:302-17-0	Chloroprene	CAS:126-99-8	Vinylchloride	CAS:75-01-4
Solvent : Methanol purge & trap/Xylene (1/traces)		ampoule 1 ml	200 µg/ml	Ref : F119211	Price : FF

VOC Standard Solution - 28 components

Acetonitrile	CAS:75-05-8	1,2,3,4-Diepoxybutane	CAS:1464-53-5	2-Picoline	CAS:109-06-8
Allyl alcohol	CAS:107-18-6	1,4-Dioxane	CAS:123-91-1	Propargyl alcohol	CAS:107-19-7
Allyl chloride	CAS:107-05-1	Epichlorhydrin	CAS:106-89-8	b-Propiolactone	CAS:57-57-8
Benzyl chloride	CAS:100-44-7	Ethyl methacrylate	CAS:97-63-2	Propionitrile	CAS:107-12-0
2-Chloroethanol	CAS:107-07-3	2-Hydroxypropionitrile	CAS:78-97-7	n- Propyl amine	CAS:107-10-8
3-Chloropropionitrile	CAS:542-76-7	Isobutyl alcohol	CAS:78-83-1	Pyridine	CAS:110-86-1
1,2-Dibromo-3-chloropropane	CAS:96-12-8	Malononitrile	CAS:109-77-3	1,1,1,2-Tetrachloroethane	CAS:630-20-6
1,2-Dibromoethane	CAS:106-93-4	Methacrylonitrile	CAS:126-98-7	1,2,3-Trichloropropane	CAS:96-18-4
Dibromomethane	CAS:74-95-3	Methyl methacrylate	CAS:80-62-6		
1,3-Dichloropropane-2-ol	CAS:96-23-1	Pentachloroethane	CAS:76-01-7		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F119201	Price : GA

VOC Standard Solution - 34 components

Acetone	CAS:67-64-1	1,2-Dichloroethane	CAS:107-06-2	Styrene	CAS:100-42-5
Benzene	CAS:71-43-2	1,1-Dichloroethene	CAS:75-35-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Bromodichloromethane	CAS:75-27-4	trans-1,2-Dichloroethene	CAS:156-60-5	Tetrachloroethene	CAS:127-18-4
Bromoform	CAS:75-25-2	1,2-Dichloropropane	CAS:78-87-5	Toluene	CAS:108-88-3
2-Butanone	CAS:78-93-3	cis-1,3-Dichloropropene	CAS:10061-01-5	1,1,1-Trichloroethane	CAS:71-55-6
Carbon disulfide	CAS:75-15-0	trans-1,3-Dichloropropene	CAS:10061-02-6	1,1,2-Trichloroethane	CAS:79-00-5
Carbon tetrachloride	CAS:56-23-5	Ethanol	CAS:64-17-5	Trichloroethene	CAS:79-01-6
Chlorobenzene	CAS:108-90-7	Ethylbenzene	CAS:100-41-4	o-Xylene	CAS:95-47-6
Chlorodibromomethane	CAS:124-48-1	2-Hexanone	CAS:591-78-6	m-Xylene	CAS:108-38-3
Chloroform	CAS:67-66-3	Methyl iodide	CAS:74-88-4	p-Xylene	CAS:106-42-3
1,4-Dichloro-2-butene	CAS:764-41-0	Methylene chloride	CAS:75-09-2		
1,1-Dichloroethane	CAS:75-34-3	4-Methyl-2-pentanone	CAS:108-10-1		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F119191	Price : HF

Surrogate Standard Mixture - 3 components

4-Bromofluorobenzene	CAS:460-00-4	1,2-Dichloroethane D4	CAS:17060-07-0	Toluene D8	CAS:2037-26-5
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F119231	Price : DF

EPA method 8260B

Volatile organic compounds by Gas chromatography/Mass spectrometry (GC/MS)

Method 8260 is used to determine volatile organic compounds in a variety of solid waste matrices. This method is applicable to nearly all types of samples, regardless of water content, including various air sampling trapping media, ground and surface water, aqueous sludges, caustic liquors, acid liquors, waste solvents, oily wastes, mousses, tars, fibrous wastes, polymeric emulsions, filter cakes, spent carbons, spent catalysts, soils, and sediments.

Internal Standard Solution - 1 component

Ethylene oxide

CAS:75-21-8

Solvent : Water	ampoule 1 ml	5 000 µg/ml	Ref : F147430	Price : BJ
-----------------	--------------	-------------	---------------	------------

Internal Standard Mixture - 4 components

Chlorobenzene D5	CAS:3114-55-4	1,4-Difluorobenzene	CAS:540-36-3	Fluorobenzene	CAS:462-06-6
1,4-Dichlorobenzene D4	CAS:3855-82-1				
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F119291	Price : DH

VOC Standard Solution - 8 components

Acrolein (2-Propenal)	CAS:107-02-8	Chloralhydrate	CAS:302-17-0	Paraldehyde	CAS:123-63-7
Acrylonitrile	CAS:107-13-1	2-Chloroethyl-vinylether	CAS:110-75-8	Vinylacetate	CAS:108-05-4
Bromoacetone	CAS:598-31-2	Chloroprene	CAS:126-99-8		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F119281	Price : EF

VOC Standard Solution - 14 components

Allyl chloride	CAS:107-05-1	Iodomethane	CAS:74-88-4	2-Picoline	CAS:109-06-8
Benzyl chloride	CAS:100-44-7	Nitrobenzene	CAS:98-95-3	n- Propyl amine	CAS:107-10-8
cis-1,4-Dichloro-2-butene	CAS:1476-11-5	2-Nitropropane	CAS:79-46-9	Pyridine	CAS:110-86-1
trans-1,4-Dichloro-2-butene	CAS:110-57-6	N-Nitrosodi-n-butylamine	CAS:924-16-3	o-Toluidine	CAS:95-53-4
Hexachloroethane	CAS:67-72-1	Pentachloroethane	CAS:76-01-7		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F119271	Price : FA

VOC Standard Solution - 15 components

Acetone	CAS:67-64-1	Dioxan	CAS:123-91-1	Methanol	CAS:67-56-1
1-Butanol	CAS:71-36-3	Ethanol	CAS:64-17-5	MIBK	CAS:108-10-1
2-Butanone	CAS:78-93-3	Ethyl Acetate	CAS:141-78-6	2-Pentanone	CAS:107-87-9
tert-butyl alcohol	CAS:75-65-0	2-Hexanone	CAS:591-78-6	1-Propanol	CAS:71-23-8
Diethylether	CAS:60-29-7	Isobutyl alcohol	CAS:78-83-1	2- Propanol	CAS:67-63-0
Solvent : Water		ampoule 1 ml	200 µg/ml	Ref : F119251	Price : FA

VOC Standard Solution - 17 components

Acetonitrile	CAS:75-05-8	1,3-Dichloro-2-propanol	CAS:96-23-1	Methacrylonitrile	CAS:126-98-7
Allyl alcohol	CAS:107-18-6	1,2,3,4-Diepoxybutane	CAS:1464-53-5	Methyl methacrylate	CAS:80-62-6
Carbon disulfide	CAS:75-15-0	Epichlorhydrin	CAS:106-89-8	Propargyl alcohol	CAS:107-19-7
2-Chloroethanol	CAS:107-07-3	Ethyl methacrylate	CAS:97-63-2	b-Propiolactone	CAS:57-57-8
3-Chloropropionitrile	CAS:542-76-7	2-Hydroxypropionitrile	CAS:78-97-7	Ethylcyanid	CAS:107-12-0
Crotonaldehyde	CAS:123-73-9	Malononitrile	CAS:109-77-3		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F119261	Price : GA

VOC Standard Solution - 48 components

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	Naphthalene	CAS:91-20-3
Bromochloromethane	CAS:74-97-5	Dichlorodifluoromethane	CAS:75-71-8	Styrene	CAS:100-42-5
Bromodichloromethane	CAS:75-27-4	1,1-Dichloroethane	CAS:75-34-3	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromoform	CAS:75-25-2	1,2-Dichloroethane	CAS:107-06-2	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Bromomethane	CAS:74-83-9	1,1-Dichloroethene	CAS:75-35-4	Tetrachloroethene	CAS:127-18-4
Carbon tetrachloride	CAS:56-23-5	cis-1,2-Dichloroethene	CAS:156-59-2	Toluene	CAS:108-88-3
Chlorobenzene	CAS:108-90-7	trans-1,2-Dichloroethene	CAS:156-60-5	1,2,4-Trichlorobenzene	CAS:120-82-1
Chlorodibromomethane	CAS:124-48-1	1,2-Dichloropropane	CAS:78-87-5	1,1,1-Trichloroethane	CAS:71-55-6
Chloroethane	CAS:75-00-3	1,3-Dichloropropane	CAS:142-28-9	1,1,2-Trichloroethane	CAS:79-00-5
Chloroform	CAS:67-66-3	2,2-Dichloropropane	CAS:594-20-7	Trichloroethene	CAS:79-01-6
Chloromethane	CAS:74-87-3	cis-1,3-Dichloropropene	CAS:10061-01-5	Trichlorofluoromethane	CAS:75-69-4
1,2-Dibromo-3-chloropropane	CAS:96-12-8	trans-1,3-Dichloropropene	CAS:10061-02-6	1,2,3-Trichloropropane	CAS:96-18-4
1,2-Dibromoethane	CAS:106-93-4	Ethylbenzene	CAS:100-41-4	Vinylchloride	CAS:75-01-4
Dibromomethane	CAS:74-95-3	Hexachlorobutadiene	CAS:87-68-3	o-Xylene	CAS:95-47-6
1,2-Dichlorobenzene	CAS:95-50-1	Isopropylbenzene	CAS:98-82-8	m-Xylene	CAS:108-38-3
1,3-Dichlorobenzene	CAS:541-73-1	Methylene chloride	CAS:75-09-2	p-Xylene	CAS:106-42-3
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F119241	Price : BCL

Surrogate Standard Mixture - 3 components

4-Bromofluorobenzene	CAS:460-00-4	1,2-Dichloroethane D4	CAS:17060-07-0	Toluene D8	CAS:2037-26-5
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F119231	Price : KK

EPA method 8270C**Semivolatile organic compounds by Gas chromatography/Mass spectrometry (GC/MS)**

Method 8270 is used to determine the concentration of semivolatile organic compounds in extracts prepared from many types of solid waste matrices, soils, air sampling media and water samples. Direct injection of a sample may be used in limited applications.

ACs Internal Standard mixture - 6 components

Acenaphthene D10	CAS:15067-26-2	1,4-Dichlorobenzene D4	CAS:3855-82-1	Perylene D12	CAS:1520-96-3
Chrysene D12	CAS:1719-03-5	Naphthalene D8	CAS:1146-65-2	Phenanthrene D10	CAS:1517-22-2
Solvent : Dichloromethane		ampoule 1 ml	2 000 µg/ml	Ref : F119821	Price : GD

ACs Surrogate Standard Mixture - 6 components

2-Fluorobiphenyl	CAS:321-60-8	Nitrobenzene D5	CAS:4165-60-0	p-Terphenyl D14	CAS:1718-51-0
2-Fluorophenol	CAS:367-12-4	Phenol D5	CAS:4165-62-2	2,4,6-Tribromophenol	CAS:118-79-6
Solvent : Dichloromethane		ampoule 1 ml	4 000 µg/ml	Ref : F119831	Price : KK

NCC Standard Solution - 9 components

N-Nitrosodibutylamine	CAS:924-16-3	N-Nitroso-diphenylamine	CAS:86-30-6	N-Nitrosopiperidine	CAS:100-75-4
N-Nitrosodiethylamine	CAS:55-18-5	N-Nitrosomethylethylamine	CAS:10595-95-6	N-Nitroso-di-n-propylamine	CAS:621-64-7
N-Nitroso-dimethylamine	CAS:62-75-9	N-Nitrosomorpholine	CAS:59-89-2	N-Nitrosopyrrolidine	CAS:930-55-2
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119341	Price : EF

NCC Standard Solution - 17 components

Anilazine	CAS:101-05-3	3,3'-Dimethoxybenzidine	CAS:119-90-4	6-Propyl-2-thiouracil	CAS:51-52-5
o-Anisidine	CAS:90-04-0	3,3'-Dimethylbenzidine	CAS:119-93-7	Saffrole	CAS:94-59-7
Barban	CAS:101-27-9	Fluchloralin	CAS:33245-39-5	Strychnine	CAS:57-24-9
Benzidine	CAS:92-87-5	Isosafrole	CAS:120-58-1	Toluene-2,4-diisocyanate	CAS:584-84-9
Bromoxynil	CAS:1689-84-5	Mestranol	CAS:72-33-3	Trifluralin	CAS:1582-09-8
3,3'-Dichlorobenzidine	CAS:91-94-1	Mexacarbate	CAS:315-18-4		
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119721	Price : FF

NCC Standard Solution - 32 components

4-Aminoazobenzene	CAS:60-09-3	p-Dimethylaminoazobenzene	CAS:60-11-7	5-Nitro-o-anisidine	CAS:99-59-2
4-Aminobiphenyl	CAS:92-67-1	Diphenylamine	CAS:122-39-4	5-Nitro-o-toluidine	CAS:99-55-8
3-Amino-9-ethylcarbazole	CAS:132-32-1	1,2-Diphenylhydrazine	CAS:122-66-7	4,4'-Oxydianiline	CAS:101-80-4
Aniline	CAS:62-53-3	4,4'-Methylene-bis(2-chloroanil	CAS:101-14-4	Phenacetin	CAS:62-44-2
4-Chloroaniline	CAS:106-47-8	4,4'-Methylene bis(N,Ndimethylaniline)	CAS:101-61-1	1,4-Phenylenediamine	CAS:106-50-3
5-Chloro-2-methylaniline	CAS:95-79-4	1-Naphthylamine	CAS:134-32-7	2-Picoline	CAS:109-06-8
3-(Chloromethyl) pyridine hydrochloride	CAS:6959-48-4	2-Naphthylamine	CAS:91-59-8	Propyzamide	CAS:23950-58-5
4-Chloro-1,2-phenylenediamine	CAS:95-83-0	Nicotine	CAS:54-11-5	Pyridine	CAS:110-86-1
4-Chloro-1,3-phenylenediamine	CAS:5131-60-2	2-Nitroaniline	CAS:88-74-4	o-Toluidine	CAS:95-53-4
p-Cresidine	CAS:120-71-8	3-Nitroaniline	CAS:99-09-2	2,4,5-Trimethylaniline	CAS:137-17-7
2,4-Diaminotoluene	CAS:95-80-7	4-Nitroaniline	CAS:100-01-6		
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119331	Price : GC

OCs Standard Solution - 1 component

Chlordane (technical)	CAS:57-74-9	ampoule 1 ml	100 µg/ml	Ref : F119801	Price : BG
Toxaphene	CAS:8001-35-2	ampoule 1 ml	100 µg/ml	Ref : F119811	Price : BG

OCs Standard Solution - 20 components

Aldrin	CAS:309-00-2	4,4'-DDT	CAS:50-29-3	Endrin ketone	CAS:53494-70-5
alpha-BHC	CAS:319-84-6	Dieldrin	CAS:60-57-1	Heptachlor	CAS:76-44-8
beta-BHC	CAS:319-85-7	Endosulfan I	CAS:959-98-8	Heptachlor epoxide	CAS:1024-57-3
delta-BHC	CAS:319-86-8	Endosulfan II	CAS:33213-65-9	Isodrin	CAS:465-73-6
gamma-BHC	CAS:58-89-9	Endosulfan-total (sulfate)	CAS:1031-07-8	Methoxychlor (DMTD)	CAS:72-43-5
4,4'-DDD (TDE)	CAS:72-54-8	Endrin	CAS:72-20-8	Mirex	CAS:2385-85-5
4,4'-DDE	CAS:72-55-9	Endrin aldehyde	CAS:7421-93-4		
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119351	Price : FF

OPP Standard Solution - 34 components

Azinphos-methyl	CAS:86-50-0	Ethion	CAS:563-12-2	Sulfoxide	CAS:120-62-7
Benzenethiol	CAS:108-98-5	EPN	CAS:2104-64-5	Phosphamidon	CAS:13171-21-6
Captafol	CAS:2425-06-1	Fensulfothion	CAS:115-90-2	Sulfallate	CAS:95-06-7
Captan	CAS:133-06-2	Fenthion	CAS:55-38-9	Terbufos	CAS:13071-79-9
Chlorfenvinphos	CAS:470-90-6	HMPA	CAS:680-31-9	cis-Tetrachlorvinphos	CAS:961-11-5
Coumaphos	CAS:56-72-4	Leptophos	CAS:21609-90-5	Tetraethyl dithiopyrophosphate	CAS:3689-24-5
Crotoxyphos	CAS:7700-17-6	Malathion	CAS:121-75-5	Tetraethylpyrophosphate	CAS:107-49-3
Demeton-O	CAS:298-03-3	Mevinphos	CAS:7786-34-7	Trimethyl phosphate	CAS:512-56-1
Demeton-S	CAS:126-75-0	Monocrotophos	CAS:6923-22-4	Tris (2,3-dibromopropyl) phosphate	CAS:126-72-7
Dichlorvos	CAS:62-73-7	Naled	CAS:300-76-5	Tri-4-cresyl phosphate	CAS:78-32-0
Dicrotophos	CAS:141-66-2	Phosalone	CAS:2310-17-0		
Dioxathion	CAS:78-34-2	Phosmet	CAS:732-11-6		
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119371	Price : GE

PCBs Standard Solution - 1 component

Arochlor 1016 (pcb 1016)	CAS:12674-11-2	ampoule 1 ml	100 µg/ml	Ref : F119731	Price : BH
Arochlor 1221 (PCB 1221)	CAS:11104-28-2	ampoule 1 ml	100 µg/ml	Ref : F119741	Price : BG
Arochlor 1232 (PCB 1232)	CAS:11141-16-5	ampoule 1 ml	100 µg/ml	Ref : F119751	Price : CC
Arochlor 1242	CAS:53469-21-9	ampoule 1 ml	100 µg/ml	Ref : F119761	Price : BH
Arochlor 1248 (PCB 1248)	CAS:12672-29-6	ampoule 1 ml	100 µg/ml	Ref : F119771	Price : BJ
Arochlor 1254	CAS:11097-69-1	ampoule 1 ml	100 µg/ml	Ref : F119781	Price : BJ
Arochlor 1260	CAS:11096-82-5	ampoule 1 ml	100 µg/ml	Ref : F119791	Price : BJ

Standard Solution - 14 components

Aramite	CAS:140-57-8	Dinoseb	CAS:88-85-7	Parathion-methyl	CAS:298-00-0
Chlorobenzilate	CAS:510-15-6	Disulfoton	CAS:298-04-4	Phorate	CAS:298-02-2
Diallate	CAS:2303-16-4	Famphur	CAS:52-85-7	Thionazin	CAS:297-97-2
Dimethoate	CAS:60-51-5	Kepone	CAS:143-50-0	O.O.O-Triethylphosphorothioate	CAS:126-68-1
Dinocap	CAS:39300-45-3	Parathion	CAS:56-38-2		
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119361	Price : FA

Standard Solution - 18 components

Acetophenone	CAS:98-86-2	Carbofuran	CAS:1563-66-2	Hydroquinone	CAS:123-31-9
2-Acetamidofluorene	CAS:53-96-3	Carbophenothion	CAS:786-19-6	Maleic anhydride	CAS:108-31-6
1-Acetyl-2-thiourea	CAS:591-08-2	Dichloro	CAS:117-80-6	1,4-Naphthoquinone	CAS:130-15-4
2-Aminoanthraquinone	CAS:117-79-3	Diethylstilbestrol	CAS:56-53-1	Nitrofen	CAS:1836-75-5
p-Benzoquinone	CAS:106-51-4	5,5-Diphenylhydantoin	CAS:57-41-0	4-Nitroquinoline-N-oxide	CAS:56-57-5
Carbaryl	CAS:63-25-2	Ethyl carbamate	CAS:51-79-6	Phthalic anhydride	CAS:85-44-9
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119381	Price : FF

VOC / Phenols Standard Solution - 25 components

Benzoic acid	CAS:65-85-0	4,6-Dinitro-2-methylphenol	CAS:534-52-1	4-Nitrophenol	CAS:100-02-7
Benzyl alcohol	CAS:100-51-6	2,4-Dinitrophenol	CAS:51-28-5	Pentachlorophenol	CAS:87-86-5
4-Chloro-3-methylphenol	CAS:59-50-7	Ethyl methanesulfonate	CAS:62-50-0	Phenol	CAS:108-95-2
2-Chlorophenol	CAS:95-57-8	Hexachlorophen	CAS:70-30-4	Resorcinol	CAS:108-46-3
2-Cyclohexyl-4,6-dinitrophenol	CAS:131-89-5	Methyl methanesulfonate	CAS:66-27-3	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
2,4-Dichlorophenol	CAS:120-83-2	2-Methylphenol	CAS:95-48-7	2,4,5-Trichlorophenol	CAS:95-95-4
2,6-Dichlorophenol	CAS:87-65-0	3-Methylphenol	CAS:108-39-4	2,4,6-Trichlorophenol	CAS:88-06-2
Diethyl sulfate	CAS:64-67-5	4-Methylphenol	CAS:106-44-5		
2,4-Dimethylphenol	CAS:105-67-9	2-Nitrophenol	CAS:88-75-5		
Solvent : Dichloromethane/Methanol p & t (1/1)		ampoule 1 ml	100 µg/ml	Ref : F119321	Price : GA

VOC / PAH / Phthalates Standard Solution - 58 components

Acenaphthene	CAS:83-32-9	Dibenzo(a,e)pyrene	CAS:192-65-4	Hexachloroethane	CAS:67-72-1
Acenaphthylene	CAS:208-96-8	1,2-Dibromo-3-chloropropane	CAS:96-12-8	Hexachloropropene	CAS:1888-71-7
Anthracene	CAS:120-12-7	Di-n-butylphthalate	CAS:84-74-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Benzo(a)anthracene	CAS:56-55-3	1,2-Dichlorobenzene	CAS:95-50-1	Isophorone	CAS:78-59-1
Benzo(b)fluoranthene	CAS:205-99-2	1,3-Dichlorobenzene	CAS:541-73-1	3-Methylcholanthrene	CAS:56-49-5
Benzo(k)fluoranthene	CAS:207-08-9	1,4-Dichlorobenzene	CAS:106-46-7	2-Methylnaphthalene	CAS:91-57-6
Benzo(g,h,i)perylene	CAS:191-24-2	Diethyl phthalate	CAS:84-66-2	Naphthalene	CAS:91-20-3
Benzo(a)pyrene	CAS:50-32-8	7,12-Dimethylbenzo(a)anthracene	CAS:57-97-6	5-Nitroacenaphthene	CAS:602-87-9
Bis-(2-chloroethoxy)-methane	CAS:111-91-1	Dimethyl phthalate	CAS:131-11-3	Nitrobenzene	CAS:98-95-3
Bis-(2-chloroethyl)-ether	CAS:111-44-4	1,2-Dinitrobenzene	CAS:528-29-0	4-Nitrobiphenyl	CAS:92-93-3
Bis-(2-chloro-1-methylethyl)ether	CAS:108-60-1	1,3-Dinitrobenzene	CAS:99-65-0	Pentachlorobenzene	CAS:608-93-5
4-Bromophenyl-phenyl ether	CAS:101-55-3	1,4-Dinitrobenzene	CAS:100-25-4	Pentachloronitrobenzene	CAS:82-68-8
Butyl benzyl phthalate	CAS:85-68-7	2,4-Dinitrotoluene	CAS:121-14-2	Phenanthrene	CAS:85-01-8
1-Chloronaphthalene	CAS:90-13-1	2,6-Dinitrotoluene	CAS:606-20-2	Pyrene	CAS:129-00-0
2-Chloronaphthalene	CAS:91-58-7	Di-n-octylphthalate	CAS:117-84-0	1,2,4,5-Tetrachlorobenzene	CAS:95-94-3
4-Chlorophenyl phenyl ether	CAS:7005-72-3	Fluoranthene	CAS:206-44-0	1,2,4-Trichlorobenzene	CAS:120-82-1
Chrysene	CAS:218-01-9	Fluorene	CAS:86-73-7	1,3,5-Trinitrobenzene	CAS:99-35-4
Dibenz(a,j)acridine	CAS:224-42-0	Hexachlorobenzene	CAS:118-74-1	Phthalic acid,bis-2-	CAS:117-81-7
Dibenzo(a,h)anthracene	CAS:53-70-3	Hexachlorobutadiene	CAS:87-68-3	ethylhexylester	
Dibenzofuran	CAS:132-64-9	Hexachlorocyclopentadiene	CAS:77-47-4		
Solvent : Dichloromethane/Benzene (3/1)		ampoule 1 ml	100 µg/ml	Ref : F119311	Price : BCA

EPA method 8275A**Semi-volatile organic compounds (PAHs & PCBs) in soils/ sludges and solid wastes using thermal extraction/ Gas chromatography/Mass spectrometry (TE/GC/MS)**

Method 8275 is a thermal extraction capillary GC/MS procedure for the rapid quantitative determination of targeted PCBs and PAHs in soils, sludges and solid wastes.

ACs internal Standard Mixture - 2 components

2-Fluorobiphenyl	CAS:321-60-8	Phenanthrene D10	CAS:1517-22-2		
Solvent : Dichloromethane		ampoule 1 ml	1 000 µg/ml	Ref : F119861	Price : DF

PAH / ACs Standard Solution - 22 components

Acenaphthene	CAS:83-32-9	4-Bromobiphenyl ether	CAS:101-55-3	Hexachlorobenzene	CAS:118-74-1
Acenaphthylene	CAS:208-96-8	1-Chloronaphthalene	CAS:90-13-1	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Naphthalene	CAS:91-20-3
Benzo(a)anthracene	CAS:56-55-3	Dibenzo(a,h)anthracene	CAS:53-70-3	Phenanthrene	CAS:85-01-8
Benzo(b)fluoranthene	CAS:205-99-2	Dibenzofuran	CAS:132-64-9	Pyrene	CAS:129-00-0
Benzo(k)fluoranthene	CAS:207-08-9	Dibenzothiophene	CAS:132-65-0	1,2,4-Trichlorobenzene	CAS:120-82-1
Benzo(g,h,i)perylene	CAS:191-24-2	Fluoranthene	CAS:206-44-0		
Benzo(a)pyrene	CAS:50-32-8	Fluorene	CAS:86-73-7		
Solvent : Dichloromethane/Toluene (3/1)		ampoule 1 ml	100 µg/ml	Ref : F119841	Price : GA

PCBs / ACs Standard Solution - 19 components

PCB 1	CAS:2051-60-7	PCB 52	CAS:35693-99-3	PCB 180	CAS:35065-29-3
PCB 11	CAS:2050-67-1	PCB 66	CAS:32598-10-0	PCB187	CAS:52663-68-0
PCB 18	CAS:37680-65-2	PCB 101	CAS:37680-73-2	PCB 194	CAS:35694-08-7
PCB 26	CAS:38444-81-4	PCB 118	CAS:31508-00-6	PCB 206	CAS:40186-72-9
PCB 31	CAS:16606-02-3	PCB128	CAS:38380-07-3	PCB 209	CAS:2051-24-3
PCB 44	CAS:41464-39-5	PCB 138	CAS:35065-28-2		
PCB 49	CAS:41464-40-8	PCB170	CAS:35065-30-6		
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119851	Price : BAF

PAH Internal Standard Mixture - 6 components

Acenaphthene D10	CAS:15067-26-2	1,4-Dichlorobenzene D4	CAS:3855-82-1	Perylene D12	CAS:1520-96-3
Chrysene D12	CAS:1719-03-5	Naphthalene D8	CAS:1146-65-2	Phenanthrene D10	CAS:1517-22-2
Solvent : Dichloromethane		ampoule 1 ml	4 000 µg/ml	Ref : F119871	Price : BBC

EPA method 8310**Polynuclear aromatic hydrocarbons**

Method 8310 is used to determine the concentration of certain polynuclear aromatic hydrocarbons (PAH) in ground water and wastes.

PAH Standard Solution - 16 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Acenaphthylene	CAS:208-96-8	Benzo(a)pyrene	CAS:50-32-8	Naphthalene	CAS:91-20-3
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Phenanthrene	CAS:85-01-8
Benzo(a)anthracene	CAS:56-55-3	Dibenzo(a,h)anthracene	CAS:53-70-3	Pyrene	CAS:129-00-0
Benzo(b)fluoranthene	CAS:205-99-2	Fluoranthene	CAS:206-44-0		
Benzo(k)fluoranthene	CAS:207-08-9	Fluorene	CAS:86-73-7		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F109861	Price : GF

PAH Quality Calibration Mixture - 16 components

Acenaphthene	CAS:83-32-9	100 µg/ml	Chrysene	CAS:218-01-9	10 µg/ml
Acenaphthylene	CAS:208-96-8	100 µg/ml	Dibenzo(a,h)anthracene	CAS:53-70-3	10 µg/ml
Anthracene	CAS:120-12-7	100 µg/ml	Fluoranthene	CAS:206-44-0	10 µg/ml
Benzo(a)anthracene	CAS:56-55-3	10 µg/ml	Fluorene	CAS:86-73-7	100 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	10 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	10 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	5 µg/ml	Naphthalene	CAS:91-20-3	100 µg/ml
Benzo(g,h,i)perylene	CAS:191-24-2	10 µg/ml	Phenanthrene	CAS:85-01-8	100 µg/ml
Benzo(a)pyrene	CAS:50-32-8	10 µg/ml	Pyrene	CAS:129-00-0	10 µg/ml
Solvent : Acetonitrile		ampoule 1 ml		Ref : F117531	Price : GF

Surrogate Standard Solution - 1 component

Decafluorobiphenyl	CAS:434-90-2				
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F119901	Price : BG

EPA method 8315A**Determination of carbonyl compounds by HPLC**

This method provides procedures for the determination of free carbonyl compounds in various matrices by derivatization with 2,4-dinitrophenylhydrazine (DNPH). The method utilizes high performance liquid chromatography (HPLC) with ultraviolet/visible (UV/vis) detection to identify and quantitate the target analytes. This method includes two procedures encompassing all aspects of the analysis (extraction to determination of concentration). Procedure 1 is appropriate for the analysis of aqueous, soil and waste samples and stack samples collected by Method 0011. Procedure 2 is appropriate for the analysis of indoor air samples collected by Method 0100.

Carbonyl Compounds Standard Solution - 2 components

Acetaldehyde	CAS:75-07-0	Formaldehyde	CAS: 50-00-0		
Solvent : Water		ampoule 1 ml	1 000 µg/ml	Ref : F120008	Price : BG

Carbonyl Compounds Standard Solution - 6 components

Decanal-DNPH	id : 4110	Hexanal-DNPH	CAS:1527-97-5	Octanal-DNPH	CAS:1726-77-8
Heptanal-DNPH	CAS:2074-05-7	Nonanal-DNPH	id : 4100	Pentanal-DNPH	CAS:2057-84-3
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F120011	Price : DF

Carbonyl Compounds Standard Solution - 16 components

Acetaldehyde-DNPH	CAS:1019-57-4	Cyclohexanone-DNPH	CAS:1589-62-4	Propanal-DNPH	CAS:725-00-8
Acetone-DNPH	CAS:1567-89-1	2,5-Dimethylbenzaldehyde-DNPH	CAS:152477-96-8	m-Tolualdehyde-DNPH	CAS:2880-05-9
Acrolein-DNPH	CAS:888-54-0	Formaldehyde-DNPH	CAS:1081-15-8	o-Tolualdehyde-DNPH	CAS:1773-44-0
Benzaldehyde-DNPH	CAS:1157-84-2	Hexaldehyde-DNPH	CAS:1527-97-5	p-Tolualdehyde-DNPH	CAS:2571-00-8
Butanal-DNPH	CAS:1527-98-6	Isovaleraldehyde-DNPH	CAS:2256-01-1		
Crotonaldehyde-DNPH	CAS:1527-96-4	Pentanal-DNPH	CAS:2057-84-3		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F120001	Price : FH

EPA method 8316

Acrylamide, acrylonitrile and acrolein by HPLC

Water samples are analyzed by high performance liquid chromatography (HPLC). A 200 µL aliquot is injected into a C-18 reverse-phase column, and compounds in the effluent are detected with an ultraviolet (UV) detector.

Carbonyl Compounds / NCC Standard Solution - 3 components

Acrolein (2-Propenal)	CAS:107-02-8	Acrylamide	CAS:79-06-1	Acrylonitrile	CAS:107-13-1
Solvent : Water		ampoule 1 ml	1 000 µg/ml	Ref : F119911	Price : DA

EPA method 8318

N-methylcarbamates by HPLC

Method 8318 is used to determine the concentration of N-methylcarbamates in soil, water and waste matrices.

Carbamates Solution - 10 components

Aldicarb	CAS:116-06-3	Dioxacarb	CAS:6988-21-2	Promecarb	CAS:2631-37-0
Aldicarb-sulfone	CAS:1646-88-4	3-Hydroxycarbofuran	CAS:16655-82-6	Propoxur	CAS:114-26-1
Carbaryl	CAS:63-25-2	Methiocarb	CAS:2032-65-7		
Carbofuran	CAS:1563-66-2	Methomyl	CAS:16752-77-5		
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F119921	Price : GA

EPA method 8321A

Solvent extractable non volatile compounds by HPLC/ Thermospray/ Mass Spectrometry (HPLC/TS/MS) or ultraviolet (UV) detection

This method covers the use of high performance liquid chromatography (HPLC), coupled with either thermospray-mass spectrometry (TS-MS), and/or ultraviolet (UV), for the determination of disperse azo dyes, organophosphorus compounds, and tris(2,3-dibromopropyl)phosphate, chlorinated phenoxyacid compounds and their esters, and carbamates in wastewater, ground water, and soil/sediment matrices. Data are also provided for chlorophenoxy acid herbicides in fly ash, however, recoveries for most compounds are very poor indicating poor extraction efficiency for these analytes using the extraction procedure included in this method. Additionally, it may apply to other non-volatile compounds that are solvent extractable, are amenable to HPLC, and are ionizable under thermospray introduction for mass spectrometric detection or may be determined by a UV detector.

Chlorinated Phenoxyacid Compounds Mixture - 14 components

2,4-D	CAS:94-75-7	Dicamba	CAS:1918-00-9	Silvex	CAS:93-72-1
2,4-D butoxyethyl ester	CAS:1929-73-3	Dichlorprop	CAS:120-36-5	2,4,5-T	CAS:93-76-5
2,4-D 2-ethylhexyl ester	CAS:1928-43-4	Dinoseb	CAS:88-85-7	2,4,5-T n-butyl ester	CAS:93-79-8
Dalapon	CAS:75-99-0	MCPA	CAS:94-74-6	2,4,5-T-butoxyethyl ester	CAS:2545-59-7
2,4-DB	CAS:94-82-6	Mecoprop (MCP)	CAS:7085-19-0		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F119941	Price : FF

Organophosphorus Compounds Mixture - 15 components

Asulam	CAS:3337-71-1	Fensulfothion	CAS:115-90-2	Parathion-methyl	CAS:298-00-0
Dichlorvos	CAS:62-73-7	Merphos	CAS:150-50-5	Phorate	CAS:298-02-2
Dimethoate	CAS:60-51-5	Methomyl	CAS:16752-77-5	Thiofanox	CAS:39196-18-4
Disulfoton	CAS:298-04-4	Monocrotophos	CAS:6923-22-4	Trichlorfon	CAS:52-68-6
Famphur	CAS:52-85-7	Naled	CAS:300-76-5	Tris(2,3-dibromopropyl)phosphate	CAS:126-72-7
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F119931	Price : GF

Carbamates Mixture - 29 components

Aldicarb	CAS:116-06-3	Carbofuran	CAS:1563-66-2	Mexacarbate	CAS:315-18-4
Aldicarb-sulfone	CAS:1646-88-4	Chloroxuron	CAS:1982-47-4	Monuron (CMU)	CAS:150-68-5
Aldicarb-sulfoxide	CAS:1646-87-3	Chlorpropham	CAS:101-21-3	Neburon	CAS:555-37-3
Aminocarb	CAS:2032-59-9	Diuron	CAS:330-54-1	Oxamyl	CAS:23135-22-0
Barban	CAS:101-27-9	Fenuron	CAS:101-42-8	Propachlor	CAS:1918-16-7
Benomyl	CAS:17804-35-2	Fluometuron	CAS:2164-17-2	Propham	CAS:122-42-9
Bendiocarb	CAS:22781-23-3	3-Hydroxycarbofuran	CAS:16655-82-6	Propoxur	CAS:114-26-1
Bromacil	CAS:314-40-9	Linuron	CAS:330-55-2	Siduron	CAS:1982-49-6
Carbaryl	CAS:63-25-2	Methiocarb	CAS:2032-65-7	Tebuthiuron	CAS:34014-18-1
Carbendazim	CAS:10605-21-7	Methomyl	CAS:16752-77-5		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F119951	Price : GF

EPA method 8325

Solvent extractable non volatile compounds by High performance liquid chromatography/ particle beam/ Mass spectrometry (HPLC/PB/MS)

This method describes the use of high performance liquid chromatography (HPLC), coupled with particle beam (PB) mass spectrometry (MS), for the determination of benzidines and nitrogen-containing pesticides in water and wastewater.

NCC Standard Solution - 12 components

Benzidine	CAS:92-87-5	250 µg/ml	3,3'-Dimethylbenzidine	CAS:119-93-7	350 µg/ml
Benzoylprop ethyl	CAS:22212-55-1	350 µg/ml	Diuron	CAS:330-54-1	450 µg/ml
Carbaryl	CAS:63-25-2	1000 µg/ml	Linuron	CAS:330-55-2	1300 µg/ml
o-Chlorophenyl thiourea	CAS:5344-82-1	750 µg/ml	Monuron (CMU)	CAS:150-68-5	400 µg/ml
3,3'-Dichlorobenzidine	CAS:91-94-1	250 µg/ml	Rotenone	CAS:83-79-4	3200 µg/ml
3,3'-Dimethoxybenzidine	CAS:119-90-4	750 µg/ml	Siduron	CAS:1982-49-6	450 µg/ml
Solvent : Acetonitrile/Methanol (1/1)		ampoule 1 ml		Ref : F119961	Price : FH

Performance Check Solution - 1 component

Decafluorotriphenylphosphine oxide	id: 17201				
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F119971	Price : CA

EPA method 8410

Gas chromatography/ Fourier transform infrared (GC/FT-IR) spectrometry for semivolatile organics: capillary column

This method covers the automated identification, or compound class assignment of unidentifiable compounds, of solvent extractable semivolatile organic compounds which are amenable to gas chromatography, by GC/FT-IR.

GC/FT-IR can be a useful complement to GC/MS analysis (Method 8270).

Internal Standard Mixture - 2 components

1-Fluoronaphthalene	CAS:321-38-0	p-Terphenyl D14	CAS:1718-51-0		
Solvent : Dichloromethane		ampoule 1 ml	2 000 µg/ml	Ref : F119981	Price : FA

EPA method 8430

Analysis of bis(2-chloroethyl)ether and hydrolysis products by direct aqueous injection GC/FT-IR

This method provides procedures for the identification and quantitation of bis(2-chloroethyl) ether and its hydrolysis compounds in aqueous matrices by direct aqueous injection (DAI) and gas chromatography with detection by a fourier transform infrared spectrometer (GC/FTIR).

Ethers and Alcohols Standard Solution - 5 components

Bis-(2-chloroethyl)-ether	CAS:111-44-4	2-(2-Chloroethoxy)ethanol	CAS:628-89-7	Ethyleneglycol	CAS:107211
2-Chloroethanol	CAS:107-07-3	Diethylene glycol	CAS:111-46-6		
Solvent : Water		ampoule 1 ml	1 000 µg/ml	Ref : F119991	Price : DF

Contaminants Standards

Aromatic Hydrocarbons

ACs Internal Stock Standard Solution - 4 components

(ISO15680:2003)

1,4-Difluorobenzene	CAS:540-36-3	Fluorobenzene	CAS:462-06-6	Monochlorobenzene D5	CAS:3114-55-4
1,4-Dichlorobenzene D4	CAS:3855-82-1				
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062075	Price : DF
		ampoule 5 ml	2 000 µg/ml	Ref : F109951	Price : BCA

ACs Standard Solution - 9 components

(EPA 612)

2-Chloronaphthalene	CAS:91-58-7	1,4-Dichlorobenzene	CAS:106-46-7	Hexachlorocyclopentadiene	CAS:77-47-4
1,2-Dichlorobenzene	CAS:95-50-1	Hexachlorobenzene	CAS:118-74-1	Hexachloroethane	CAS:67-72-1
1,3-Dichlorobenzene	CAS:541-73-1	Hexachlorobutadiene	CAS:87-68-3	1,2,4-Trichlorobenzene	CAS:120-82-1
Solvent : Methanol / Dichloromethane (1/1)		ampoule 1 ml	100 µg/ml	Ref : F115571	Price : DF

ACs Standard Solution - 9 components

(EPA 612)

2-Chloronaphthalene	CAS:91-58-7	400 µg/ml	Hexachlorobutadiene	CAS:87-68-3	1 µg/ml
1,2-Dichlorobenzene	CAS:95-50-1	200 µg/ml	Hexachlorocyclopentadiene	CAS:77-47-4	1 µg/ml
1,3-Dichlorobenzene	CAS:541-73-1	200 µg/ml	Hexachloroethane	CAS:67-72-1	1 µg/ml
1,4-Dichlorobenzene	CAS:106-46-7	400 µg/ml	1,2,4-Trichlorobenzene	CAS:120-82-1	100 µg/ml
Hexachlorobenzene	CAS:118-74-1	1 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F115581	Price : DF

ACs Standard Solution - 10 components

(EPA 8020A)

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	m-Xylene	CAS:108-38-3
Chlorobenzene	CAS:108-90-7	Ethylbenzene	CAS:100-41-4	p-Xylene	CAS:106-42-3
1,2-Dichlorobenzene	CAS:95-50-1	Toluene	CAS:108-88-3		
1,3-Dichlorobenzene	CAS:541-73-1	o-Xylene	CAS:95-47-6		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F116291	Price : EF

ACs Standard Solution - 11 components

(EPA 8020A)

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	m-Xylene	CAS:108-38-3
Chlorobenzene	CAS:108-90-7	Ethylbenzene	CAS:100-41-4	p-Xylene	CAS:106-42-3
1,2-Dichlorobenzene	CAS:95-50-1	Toluene	CAS:108-88-3	Styrene	CAS:100-42-5
1,3-Dichlorobenzene	CAS:541-73-1	o-Xylene	CAS:95-47-6		
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F116301	Price : DF

Standard Solution - 22 components

(EN 14039)

Bromobenzene	CAS:108-86-1	1,3-Dichlorobenzene	CAS:541-73-1	1,2,3-Trichlorobenzene	CAS:87-61-6
n-Butylbenzene	CAS:104-51-8	1,4-Dichlorobenzene	CAS:106-46-7	1,2,4-Trichlorobenzene	CAS:120-82-1
sec-Butylbenzene	CAS:135-98-8	Hexachlorobutadiene	CAS:87-68-3	1,3,5-Trichlorobenzene	CAS:108-70-3
tert-Butylbenzene	CAS:98-06-6	Isopropylbenzene	CAS:98-82-8	1,2,3-Trichloropropane	CAS:96-18-4
2-Chlorotoluene	CAS:95-49-8	4-Isopropyltoluene	CAS:99-87-6	1,2,4-Trimethylbenzene	CAS:95-63-6
4-Chlorotoluene	CAS:106-43-4	Naphthalene	CAS:91-20-3	1,3,5-Trimethylbenzene	CAS:108-67-8
1,2-Dibromo-3-chloropropane	CAS:96-12-8	n-Propylbenzene	CAS:103-65-1		
1,2-Dichlorobenzene	CAS:95-50-1	1,1,2,2-Tetrabromoethane	CAS:79-27-6		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062155	Price : DF
		ampoule 5 ml	2 000 µg/ml	Ref : F109991	Price : BDA

Aromatics & Alkenes Mixture - 28 components

(EPA 503.1)

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	1,2,3-Trichlorobenzene	CAS:87-61-6
Bromobenzene	CAS:108-86-1	Ethylbenzene	CAS:100-41-4	1,2,4-Trichlorobenzene	CAS:120-82-1
n-Butylbenzene	CAS:104-51-8	Hexachlorobutadiene	CAS:87-68-3	Trichloroethene	CAS:79-01-6
sec-Butylbenzene	CAS:135-98-8	Isopropylbenzene	CAS:98-82-8	1,2,4-Trimethylbenzene	CAS:95-63-6
tert-Butylbenzene	CAS:98-06-6	4-Isopropyltoluene	CAS:99-87-6	1,3,5-Trimethylbenzene	CAS:108-67-8
Chlorobenzene	CAS:108-90-7	Naphthalene	CAS:91-20-3	o-Xylene	CAS:95-47-6
2-Chlorotoluene	CAS:95-49-8	n-Propylbenzene	CAS:103-65-1	m-Xylene	CAS:108-38-3
4-Chlorotoluene	CAS:106-43-4	Styrene	CAS:100-42-5	p-Xylene	CAS:106-42-3
1,2-Dichlorobenzene	CAS:95-50-1	Tetrachloroethene	CAS:127-18-4		
1,3-Dichlorobenzene	CAS:541-73-1	Toluene	CAS:108-88-3		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113041	Price : DF
		ampoule 1 ml	2 000 µg/ml	Ref : F113081	Price : HA

Aromatic Halocarbons Mixture - 9 components

(EPA 502)

Bromobenzene	CAS:108-86-1	4-Chlorotoluene	CAS:106-43-4	1,4-Dichlorobenzene	CAS:106-46-7
Chlorobenzene	CAS:108-90-7	1,2-Dichlorobenzene	CAS:95-50-1	1,2,3-Trichlorobenzene	CAS:87-61-6
2-Chlorotoluene	CAS:95-49-8	1,3-Dichlorobenzene	CAS:541-73-1	1,2,4-Trichlorobenzene	CAS:120-82-1
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113011	Price : DF

Aromatic Halocarbons Mixture - 12 components**(EPA 503.1)**

Bromobenzene	CAS:108-86-1	1,2-Dichlorobenzene	CAS:95-50-1	Tetrachloroethene	CAS:127-18-4
Chlorobenzene	CAS:108-90-7	1,3-Dichlorobenzene	CAS:541-73-1	1,2,3-Trichlorobenzene	CAS:87-61-6
2-Chlorotoluene	CAS:95-49-8	1,4-Dichlorobenzene	CAS:106-46-7	1,2,4-Trichlorobenzene	CAS:120-82-1
4-Chlorotoluene	CAS:106-43-4	Hexachlorobutadiene	CAS:87-68-3	Trichloroethene	CAS:79-01-6
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113101	Price : EF

Aromatic Hydrocarbons Mixture - 16 components**(EPA 503.1; EPA 502)**

Benzene	CAS:71-43-2	4-Isopropyltoluene	CAS:99-87-6	1,3,5-Trimethylbenzene	CAS:108-67-8
n-Butylbenzene	CAS:104-51-8	Naphthalene	CAS:91-20-3	m-Xylene	CAS:108-38-3
sec-Butylbenzene	CAS:135-98-8	n-Propylbenzene	CAS:103-65-1	o-Xylene	CAS:95-47-6
tert-Butylbenzene	CAS:98-06-6	Styrene	CAS:100-42-5	p-Xylene	CAS:106-42-3
Ethylbenzene	CAS:100-41-4	Toluene	CAS:108-88-3		
Isopropylbenzene	CAS:98-82-8	1,2,4-Trimethylbenzene	CAS:95-63-6		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F113001	Price : FA

Base/Neutral Extractables Mixture - 44 components**(EPA 625)**

Acenaphthene	CAS:83-32-9	2-Chloronaphthalene	CAS:91-58-7	Hexachlorobenzene	CAS:118-74-1
Acenaphthylene	CAS:208-96-8	4-Chlorophenyl phenyl ether	CAS:7005-72-3	Hexachlorobutadiene	CAS:87-68-3
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Hexachlorocyclopentadiene	CAS:77-47-4
Azobenzene	CAS:103-33-3	Dibenzo(a,h)anthracene	CAS:53-70-3	Hexachloroethane	CAS:67-72-1
Benzo(a)anthracene	CAS:56-55-3	Di-n-butylphthalate	CAS:84-74-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Benzo(b)fluoranthene	CAS:205-99-2	1,2-Dichlorobenzene	CAS:95-50-1	Isophorone	CAS:78-59-1
Benzo(k)fluoranthene	CAS:207-08-9	1,3-Dichlorobenzene	CAS:541-73-1	Naphthalene	CAS:91-20-3
Benzo(g,h,i)perylene	CAS:191-24-2	1,4-Dichlorobenzene	CAS:106-46-7	Nitrobenzene	CAS:98-95-3
Benzo(a)pyrene	CAS:50-32-8	Diethyl phthalate	CAS:84-66-2	N-Nitroso-dimethylamine	CAS:62-75-9
Bis-(2-chloroethyl)-ether	CAS:111-44-4	Dimethyl phthalate	CAS:131-11-3	N-nitrosodi-n-propylamine	CAS:621-64-7
Bis-(2-chloroethoxy)-methane	CAS:111-91-1	2,4-Dinitrotoluene	CAS:121-14-2	N-Nitroso-diphenylamine	CAS:86-30-6
Bis(2-ethylhexyl)phthalate	CAS:117-81-7	2,6-Dinitrotoluene	CAS:606-20-2	Phenanthrene	CAS:85-01-8
bis(2-Chloroisopropyl)ether	CAS:108-60-1	di-n-octyl phthalate	CAS:117-84-0	Pyrene	CAS:129-00-0
4-bromophenyl phenyl ether	CAS:101-55-3	Fluoranthene	CAS:206-44-0	1,2,4-Trichlorobenzene	CAS:120-82-1
Butyl benzyl phthalate	CAS:85-68-7	Fluorene	CAS:86-73-7		
Solvent : Dichloromethane/Benzene/Acetonitrile (2/2/1) ampoule 1 ml		2 000 µg/ml	Ref : F115871	Price : BJA	

Internal Stock Standard Solution for aromatic hydrocarbons - 2 components**(ISO15009:2002)**

Ethylbenzene D10	CAS:25837-05-2	Toluene D8	CAS:2037-26-5		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062035	Price : CF
		ampoule 5 ml	2 000 µg/ml	Ref : F109911	Price : HF

Purgeable Aromatics Mixture - 7 components**(EPA 602)**

Benzene	CAS:71-43-2	1,3-Dichlorobenzene	CAS:541-73-1	Toluene	CAS:108-88-3
Chlorobenzene	CAS:108-90-7	1,4-Dichlorobenzene	CAS:106-46-7		
1,2-Dichlorobenzene	CAS:95-50-1	Ethylbenzene	CAS:100-41-4		
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F114881	Price : DA

Volatile Aromatic Hydrocarbons - 8 components**(ISO15009:2002)**

Benzene	CAS:71-43-2	Styrene	CAS:100-42-5	m-Xylene	CAS:108-38-3
Ethylbenzene	CAS:100-41-4	Toluene	CAS:108-88-3	p-Xylene	CAS:106-42-3
Naphthalene	CAS:91-20-3	o-Xylene	CAS:95-47-6		
Solvent : Methanol purge & trap		ampoule 1 ml	4 000 µg/ml	Ref : F062055	Price : CF
		ampoule 5 ml	4 000 µg/ml	Ref : F109931	Price : JA

Volatile Aromatic Organic Mixture - 25 components**(EPA 502)**

Benzene	CAS:71-43-2	1,3-Dichlorobenzene	CAS:541-73-1	1,2,3-Trichlorobenzene	CAS:87-61-6
Bromobenzene	CAS:108-86-1	1,4-Dichlorobenzene	CAS:106-46-7	1,2,4-Trichlorobenzene	CAS:120-82-1
n-Butylbenzene	CAS:104-51-8	Ethylbenzene	CAS:100-41-4	1,2,4-Trimethylbenzene	CAS:95-63-6
sec-Butylbenzene	CAS:135-98-8	Isopropylbenzene	CAS:98-82-8	1,3,5-Trimethylbenzene	CAS:108-67-8
tert-Butylbenzene	CAS:98-06-6	4-Isopropyltoluene	CAS:99-87-6	o-Xylene	CAS:95-47-6
Chlorobenzene	CAS:108-90-7	Naphthalene	CAS:91-20-3	m-Xylene	CAS:108-38-3
2-Chlorotoluene	CAS:95-49-8	n-Propylbenzene	CAS:103-65-1	p-Xylene	CAS:106-42-3
4-Chlorotoluene	CAS:106-43-4	Styrene	CAS:100-42-5		
1,2-Dichlorobenzene	CAS:95-50-1	Toluene	CAS:108-88-3		
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F113021	Price : EF
		ampoule 1 ml	200 µg/ml	Ref : F113031	Price : EF
		ampoule 1 ml	2 000 µg/ml	Ref : F113051	Price : HA

Carbonyl & Carbonyl-DNPH standards**Aldehyde/Ketone-DNPH Standard Solution - 8 components**

Acetaldehyde-DNPH	CAS:1019-57-4	2-Butanone-DNPH	CAS:958-60-1	Methyl Isobutyl Ketone-DNPH	CAS:1655-42-1
Acetone-DNPH	CAS:1567-89-1	Formaldehyde-DNPH	CAS:1081-15-8	Propionaldehyde-DNPH	CAS:725-00-8
Acrolein-DNPH	CAS:888-54-0	Methacrolein-DNPH	CAS:5077-73-6		
Solvent : Acetonitrile		ampoule 1 ml	15 µg/ml	Ref : F128671	Price : EF

Aldehyde/Ketone-DNPH Working Standard Solution - 13 components

Acetaldehyde-DNPH	CAS:1019-57-4	n-Butyraldehyde-DNPH	CAS:1527-98-6	Propionaldehyde-DNPH	CAS:725-00-8
Acetone-DNPH	CAS:1567-89-1	Crotonaldehyde-DNPH	CAS:1527-96-4	m-Tolualdehyde-DNPH	CAS:2880-05-9
Acrolein-DNPH	CAS:888-54-0	Formaldehyde-DNPH	CAS:1081-15-8	Valeraldehyde-DNPH	CAS:2057-84-3
Benzaldehyde-DNPH	CAS:1157-84-2	Hexaldehyde-DNPH	CAS:1527-97-5		
2-Butanone-DNPH	CAS:958-60-1	Methacrolein-DNPH	CAS:5077-73-6		
Solvent : Acetonitrile		ampoule 1 ml	3 µg/ml	Ref : F128651	Price : FC

Aldehyde/Ketone-DNPH Check Standard Solution - 13 components

Acetaldehyde-DNPH	CAS:1019-57-4	n-Butyraldehyde-DNPH	CAS:1527-98-6	Propionaldehyde-DNPH	CAS:725-00-8
Acetone-DNPH	CAS:1567-89-1	Crotonaldehyde-DNPH	CAS:1527-96-4	m-Tolualdehyde-DNPH	CAS:2880-05-9
Acrolein-DNPH	CAS:888-54-0	Formaldehyde-DNPH	CAS:1081-15-8	Valeraldehyde-DNPH	CAS:2057-84-3
Benzaldehyde-DNPH	CAS:1157-84-2	Hexaldehyde-DNPH	CAS:1527-97-5		
2-Butanone-DNPH	CAS:958-60-1	Methacrolein-DNPH	CAS:5077-73-6		
Solvent : Acetonitrile		ampoule 1 ml	1 µg/ml	Ref : F128661	Price : FC

Aldehydes Mixtures (556 & 556.1) - 14 components**(EPA 556.1)**

Acetaldehyde	CAS:75-07-0	Ethanedial	CAS:107-22-2	Nonanal	CAS:124-19-6
Benzaldehyde	CAS:100-52-7	Formaldehyde	CAS:50-00-0	Octanal	CAS:124-13-0
Butanal	CAS:123-72-8	Heptanal	CAS:111-71-7	Pentanal	CAS:110-62-3
Cyclohexanone	CAS:108-94-1	Hexanal	CAS:66-25-1	Propanal	CAS:123-38-6
Decanal	CAS:112-31-2	Methylglyoxal	CAS:78-98-8		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F114771	Price : EF

Aldehydes Mixtures(556) - 15 components**(EPA 556)**

Acetaldehyde	CAS:75-07-0	Decanal	CAS:112-31-2	Methylglyoxal	CAS:78-98-8
Benzaldehyde	CAS:100-52-7	Ethanedial	CAS:107-22-2	Nonanal	CAS:124-19-6
Butanal	CAS:123-72-8	Formaldehyde	CAS:50-00-0	Octanal	CAS:124-13-0
Crotonaldehyde	CAS:123-73-9	Heptanal	CAS:111-71-7	Pentanal	CAS:110-62-3
Cyclohexanone	CAS:108-94-1	Hexanal	CAS:66-25-1	Propanal	CAS:123-38-6
Solvent : Acetonitrile		ampoule 1 ml	1 000 µg/ml	Ref : F114681	Price : EF

Carbonyl Compounds Standard Solution - 2 components**(EPA 603)**

Acrolein (2-Propenal)	CAS:107-02-8	Acrylonitrile	CAS:107-13-1		
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F114931	Price : CF
		ampoule 1 ml	2 000 µg/ml	Ref : F114941	Price : CF
Solvent : Water		ampoule 1 ml	1 000 µg/ml	Ref : F114951	Price : CF
		ampoule 1 ml	10 000 µg/ml	Ref : F114961	Price : CF

Carbonyl Compounds Standard Solution - 2 components**(EPA 8030A)**

Acrolein (2-Propenal)	CAS:107-02-8	Acrylonitrile	CAS:107-13-1		
Solvent : Water		ampoule 1 ml	100 µg/ml	Ref : F116441	Price : CF
		ampoule 1 ml	2 000 µg/ml	Ref : F116421	Price : CF

Carbonyl Compounds Standard Solution - 2 components**(EPA 8315A)**

Acetaldehyde	CAS:75-07-0	Formaldehyde	CAS:50-00-0		
Solvent : Water		ampoule 1 ml	1 000 µg/ml	Ref : F120008	Price : BG

Carbonyl Compounds Standard Solution - 4 components

Acetaldehyde-DNPH	CAS:1019-57-4	Acrolein-DNPH	CAS:888-54-0	Formaldehyde-DNPH	CAS:1081-15-8
Acetone-DNPH	CAS:1567-89-1				
Solvent : Acetonitrile		ampoule 1 ml	5 µg/ml	Ref: F128691	Price : DF

Carbonyl Compounds Standard Solution - 6 components**(EPA 8315A)**

Decanal-DNPH	id : 4110	Hexanal-DNPH	CAS:1527-97-5	Octanal-DNPH	CAS:1726-77-8
Heptanal-DNPH	CAS:2074-05-7	Nonanal-DNPH	id : 4100	Pentanal-DNPH	CAS:2057-84-3
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F120011	Price : DF

Carbonyl Compounds Standard Solution - 12 components**(EPA 554)**

Acetaldehyde	CAS:75-07-0	Decanal	CAS:112-31-2	Nonanal	CAS:124-19-6
Butanal	CAS:123-72-8	Formaldehyde	CAS:50-00-0	Octanal	CAS:124-13-0
Crotonaldehyde	CAS:123-73-9	Heptanal	CAS:111-71-7	Pentanal	CAS:110-62-3
Cyclohexanone	CAS:108-94-1	Hexanal	CAS:66-25-1	Propanal	CAS:123-38-6
Solvent : Acetonitrile		ampoule 1 ml	1 000 µg/ml	Ref : F114641	Price : EF

Carbonyl Compounds Standard Solution - 13 components

Acetaldehyde-DNPH	CAS:1019-57-4	n-Butyraldehyde-DNPH	CAS:1527-98-6	Propionaldehyde-DNPH	CAS:725-00-8
Acetone-DNPH	CAS:1567-89-1	Crotonaldehyde-DNPH	CAS:1527-96-4	m-Tolualdehyde-DNPH	CAS:2880-05-9
Acrolein-DNPH	CAS:888-54-0	Formaldehyde-DNPH	CAS:1081-15-8	Valeraldehyde-DNPH	CAS:2057-84-3
Benzaldehyde-DNPH	CAS:1157-84-2	Hexaldehyde-DNPH	CAS:1527-97-5		
2-Butanone-DNPH	CAS:958-60-1	Methacrolein-DNPH	CAS:5077-73-6		
Solvent : Acetonitrile		ampoule 1 ml	15 µg/ml	Ref : F128641	Price : FC

Carbonyl Compounds Standard Solution - 15 components

Acetaldehyde-DNPH	CAS:1019-57-4	Crotonaldehyde-DNPH	CAS:1527-96-4	Propionaldehyde-DNPH	CAS:725-00-8
Acetone-DNPH	CAS:1567-89-1	2,5-Dimethylbenzaldehyde-DNPH	CAS:152477-96-8	o-Tolualdehyde-DNPH	CAS:1773-44-0
Acrolein-DNPH	CAS:888-54-0	Formaldehyde-DNPH	CAS:1081-15-8	p-Tolualdehyde-DNPH	CAS:2571-00-8
Benzaldehyde-DNPH	CAS:1157-84-2	Hexaldehyde-DNPH	CAS:1527-97-5	m-Tolualdehyde-DNPH	CAS:2880-05-9
n-Butyraldehyde-DNPH	CAS:1527-98-6	Isovaleraldehyde-DNPH	CAS:2256-01-1	Valeraldehyde-DNPH	CAS:2057-84-3
Solvent : Acetonitrile		ampoule 1 ml	5 µg/ml	Ref : F128681	Price : FC
		ampoule 1 ml	15 µg/ml	Ref : F128631	Price : FC

Carbonyl Compounds Standard Solution - 16 components

Acetaldehyde-DNPH	CAS:1019-57-4	Cyclohexanone-DNPH	CAS:1589-62-4	Propanal-DNPH	CAS:725-00-8
Acetone-DNPH	CAS:1567-89-1	2,5-Dimethylbenzaldehyde-DNPH	CAS:152477-96-8	m-Tolualdehyde-DNPH	CAS:2880-05-9
Acrolein-DNPH	CAS:888-54-0	Formaldehyde-DNPH	CAS:1081-15-8	o-Tolualdehyde-DNPH	CAS:1773-44-0
Benzaldehyde-DNPH	CAS:1157-84-2	Hexaldehyde-DNPH	CAS:1527-97-5	p-Tolualdehyde-DNPH	CAS:2571-00-8
Butanal-DNPH	CAS:1527-98-6	Isovaleraldehyde-DNPH	CAS:2256-01-1		
Crotonaldehyde-DNPH	CAS:1527-96-4	Pentanal-DNPH	CAS:2057-84-3		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F120001	Price : FH

Carbonyl Compounds / NCC Standard Solution - 3 components

Acrolein (2-Propenal)	CAS:107-02-8	Acrylamide	CAS:79-06-1	Acrylonitrile	CAS:107-13-1
Solvent : Water		ampoule 1 ml	1 000 µg/ml	Ref : 119911	Price : DA

Derivatised Carbonyl Compounds Mixture - 12 components

Acetaldehyde-DNPH	CAS:1019-57-4	Decanal-DNPH	id : 4110	Nonanal-DNPH	id : 4100
Butanal-DNPH	CAS:1527-98-6	Formaldehyde-DNPH	CAS:1081-15-8	Octanal-DNPH	CAS:1726-77-8
Crotonaldehyde-DNPH	CAS:1527-96-4	Heptanal-DNPH	CAS:2074-05-7	Pentanal-DNPH	CAS:2057-84-3
Cyclohexanone-DNPH	CAS:1589-62-4	Hexaldehyde-DNPH	CAS:1527-97-5	Propanal-DNPH	CAS:725-00-8
Solvent : Methanol/Acetonitrile (4/1)		ampoule 1 ml	1 000 µg/ml	Ref : F114651	Price : BCA

Ethers & Alcohols

Ethers Standard Solution - 4 components

Bis-(2-chloroethoxy)-methane	CAS:111-91-1	bis(2-chloroisopropyl) ether	CAS:108-60-1	4-Chlorophenyl phenyl ether	CAS:7005-72-3
Bis-(2-chloroethyl)-ether	CAS:111-44-4				
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F117551	Price : DF

Ethers and Alcohols Standard Solution - 5 components

Bis-(2-chloroethyl)-ether	CAS:111-44-4	2-(2-Chloroethoxy)ethanol	CAS:628-89-7	Ethyleneglycol	CAS:107211
2-Chloroethanol	CAS:107-07-3	Diethylene glycol	CAS:111-46-6		
Solvent : Water		ampoule 1 ml	1 000 µg/ml	Ref : F119991	Price : DF

QC Haloethers Mixture - 5 components

Bis-(2-chloroethyl)-ether	CAS:111-44-4	Bis-(2-chloro-1-methylethyl)ether	CAS:108-60-1	4-Chlorophenyl phenyl ether	CAS:7005-72-3
Bis-(2-chloroethoxy)-methane	CAS:111-91-1	4-bromophenyl phenyl ether	CAS:101-55-3		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F115561	Price : DA

Haloorganic Acids & Esters

Chlorinated Acids Mixture A - 8 components

Acifluorfen	CAS:50594-66-6	2,4-D	CAS:94-75-7	Picloram	CAS:1918-02-1
Bentazone	CAS:25057-89-0	Dicamba	CAS:1918-00-9	Silvex	CAS:93-72-1
Chloramben	CAS:133-90-4	Dichlorprop	CAS:120-36-5		
Solvent : Acetonitrile		ampoule 1 ml	1 000 µg/ml	Ref : F114661	Price : EC

Chlorinated Acids Mixture B - 8 components

2,4-DB	CAS:94-82-6	MCPA	CAS:94-74-6	Pentachlorophenol	CAS:87-86-5
3,5-Dichlorobenzoic acid	CAS:51-36-5	Mecoprop (MCP)	CAS:7085-19-0	2,4,5-T	CAS:93-76-5
Dinoseb	CAS:88-85-7	4-Nitrophenol	CAS:100-02-7		
Solvent : Acetonitrile		ampoule 1 ml	1 000 µg/ml	Ref : F114671	Price : EC

Haloacetic Acids Mixture - 8 components**(EPA 552)**

Bromochloroacetic acid	CAS:5589-96-8	2,4-Dichlorophenol	CAS:120-83-2	Trichloroacetic acid	CAS:76-03-9
Dibromoacetic acid	CAS:631-64-1	Monobromoacetic acid	CAS:79-08-3	2,4,6-Trichlorophenol	CAS:88-06-2
Dichloroacetic acid	CAS:79-43-6	Monochloroacetic acid	CAS:79-11-8		

Solvent : Methyl-tert butylether ampoule 1 ml 1 000 µg/ml Ref : F114241 Price : FA

Haloacetic Acids Mixture without Surrogate - 10 components**(EPA 552)**

Bromochloroacetic acid	CAS:5589-96-8	40 µg/ml	Dichloroacetic acid	CAS:79-43-6	60 µg/ml
Bromodichloroacetic acid	CAS:71133-14-7	40 µg/ml	Monobromoacetic acid	CAS:79-08-3	40 µg/ml
Chlorodibromoacetic acid	CAS:5278-95-5	100 µg/ml	Monochloroacetic acid	CAS:79-11-8	60 µg/ml
Dalapon	CAS:75-99-0	40 µg/ml	Tribromoacetic acid	CAS:75-96-7	200 µg/ml
Dibromoacetic acid	CAS:631-64-1	20 µg/ml	Trichloroacetic acid	CAS:76-03-9	20 µg/ml

Solvent : Methyl-tert butylether ampoule 1 ml Ref : F114571 Price : FA

Haloacetic Acids Mixture with Surrogate - 11 components**(EPA 552)**

Bromochloroacetic acid	CAS:5589-96-8	40 µg/ml	Dibromoacetic acid	CAS:631-64-1	20 µg/ml
Bromodichloroacetic acid	CAS:71133-14-7	40 µg/ml	2,3-Dibromopropionic acid	CAS:600-05-5	100 µg/ml
Chlorodibromoacetic acid	CAS:5278-95-5	100 µg/ml	Monobromoacetic Acid	CAS:79-08-3	40 µg/ml
Dalapon	CAS:75-99-0	40 µg/ml	Monochloroacetic Acid	CAS:79-11-8	60 µg/ml
Dibromoacetic acid	CAS:631-64-1	20 µg/ml	Tribromoacetic acid	CAS:75-96-7	200 µg/ml
Dichloroacetic Acid	CAS:79-43-6	60 µg/ml			

Solvent : Methyl-tert butylether ampoule 1 ml Ref : F114431 Price : FA

Methylated Haloacetic Acids Mixture - 8 components**(EPA 552)**

Bromoacetic acid-methyl ester	CAS:96-32-2	Dibromoacetic acid-methyl ester	CAS:6482-26-4	Trichloroacetic acid-methyl ester	CAS:598-99-2
Bromochloroacetic acid-methyl ester	CAS:20428-74-4	Dichloroacetic acid-methyl ester	CAS:116-54-1	2,4,6-Trichloroanisole	CAS:87-40-1
Chloroacetic acid methyl ester	CAS:96-34-4	2,4-Dichloroanisole	CAS:553-82-2		

Solvent : Methyl-tert butylether ampoule 1 ml 1 000 µg/ml Ref : F114251 Price : GA

Methylated Haloacetic Acids Mixture without Surrogate - 10 components**(EPA 552)**

Bromodichloroacetic acid-methyl ester	CAS:20428-76-6	40 µg/ml	Methyl bromochloroacetate	CAS:20428-74-4	40 µg/ml
Dalapon methyl ester	CAS:17640-02-7	40 µg/ml	Methyl chloroacetate	CAS:96-34-4	60 µg/ml
Dibromochloroacetic acid-methyl ester	CAS:20428-75-5	100 µg/ml	Methyl dibromoacetate	CAS:6482-26-4	20 µg/ml
Dichloroacetic acid-methyl ester	CAS:116-54-1	60 µg/ml	Methyl tribromoacetate	CAS:3222-05-7	200 µg/ml
Methyl bromoacetate	CAS:96-32-2	40 µg/ml	Methyl trichloroacetate	CAS:598-99-2	20 µg/ml

Solvent : Methyl-tert butylether ampoule 1 ml Ref : F114591 Price : FF

Methylated Haloacetic Acids Mixture with Surrogate - 11 components**(EPA 552)**

Bromodichloroacetic acid-methyl ester	CAS:20428-76-6	40 µg/ml	Methyl chloroacetate	CAS:96-34-4	60 µg/ml
Dalapon methyl ester	CAS:17640-02-7	40 µg/ml	Methyl dibromoacetate	CAS:6482-26-4	20 µg/ml
Dibromochloroacetic acid-methyl ester	CAS:20428-75-5	100 µg/ml	Methyl tribromoacetate	CAS:3222-05-7	200 µg/ml
Dichloroacetic acid-methyl ester	CAS:116-54-1	60 µg/ml	Methyl trichloroacetate	CAS:598-99-2	20 µg/ml
Methyl bromoacetate	CAS:96-32-2	40 µg/ml	Methyl 2,3-dibromopropionate	CAS:1729-67-5	100 µg/ml
Methyl bromochloroacetate	CAS:20428-74-4	40 µg/ml			

Solvent : Methyl-tert butylether ampoule 1 ml Ref : F114581 Price : FF

Standard Solution - 7 components**(ISO 23631:2006)**

Bromochloroacetic Acid	CAS:5589-96-8	Dichloroacetic Acid	CAS:79-43-6	Trichloroacetic acid	CAS:76-03-9
Dalapon	CAS:75-99-0	Monobromoacetic Acid	CAS:79-08-3		
Dibromoacetic acid	CAS:631-64-1	Monochloroacetic Acid	CAS:79-11-8		

Solvent : Methyl-tert.butylether ampoule 10 ml 10 µg/ml Ref : F129500 Price : EF

Standard Solution - 7 components**(ISO 23631:2006)**

Bromochloroacetic acid-methyl ester	CAS:20428-74-4	Dichloroacetic acid-methyl ester	CAS:116-54-1	Trichloroacetic acid-methyl ester	CAS:598-99-2
Dalapon methyl ester	CAS:17640-02-7	Monobromoacetic Acid	CAS:96-32-2		
Dibromoacetic acid-methyl ester	CAS:6482-26-4	Monochloroacetic Acid	CAS:96-34-4		

Solvent : Methyl-tert.butylether ampoule 10 ml 10 µg/ml Ref : F129510 Price : EF

Hydrocarbon Standards**Alkane mixtures****n-Alkane Standard Solution - 6 components**

n-Decane (C10)	CAS:124-18-5	n-Eicosane (C20)	CAS:112-95-8	n-Triacontane (C30)	CAS:638-68-6
n-Hexadecane (C16)	CAS:544-76-3	n-Tetracosane (C24)	CAS:646-31-1	n-Tetracontane (C40)	CAS:4181-95-7

Solvent : n-Hexane ampoule 1 ml 100 µg/ml Ref : F062005 Price : DA

ampoule 5 ml 100 µg/ml Ref : F109791 Price : KJ

n-Alkane Standard Solution -16 components C10÷C40 (all even)**(ISO 9377-2-Mod)**

n-Decane (C10)	CAS:124-18-5	n-Docosane (C22)	CAS:629-97-0	n-Tetatriacontane (C34)	CAS:14167-59-0
n-Dodecane (C12)	CAS:112-40-3	n-Tetracosane (C24)	CAS:646-31-1	n-Hexatriacontane (C36)	CAS:630-06-8
n-Tetradecane (C14)	CAS:629-59-4	n-Hexacosane (C26)	CAS:630-01-3	n-Octatriacontane (C38)	CAS:7194-85-6
n-Hexadecane (C16)	CAS:544-76-3	n-Octacosane (C28)	CAS:630-02-4	n-Tetracontane (C40)	CAS:4181-95-7
n-Octadecane (C18)	CAS:593-45-3	n-Triacontane (C30)	CAS:638-68-6		
n-Eicosane (C20)	CAS:112-95-8	n-Dotriacontane (C32)	CAS:544-85-4		

Solvent : n-Hexane/Petroleum ether (1/1)	ampoule 1 ml	100 µg/ml	Ref : P858130	Price : GA
	ampoule 5 ml	100 µg/ml	Ref : F109771	Price : BGA
Solvent : n-Hexane	ampoule 1 ml	50 µg/ml	Ref : P858170	Price : GA
	ampoule 5 ml	50 µg/ml	Ref : F109781	Price : BGA

n-Alkane Standard Solution - 19 components C7+C8+C9+C10÷C40 (all even)					(ISO 9377-2-Mod)
n-Heptane (C7)	CAS:142-82-5	n-Octadecane (C18)	CAS:593-45-3	n-Dotriacontane (C32)	CAS:544-85-4
n-Octane (C8)	CAS:111-65-9	n-Eicosane (C20)	CAS:112-95-8	n-Tetracontane (C34)	CAS:14167-59-0
n-Nonane (C9)	CAS:111-84-2	n-Docosane (C22)	CAS:629-97-0	n-Hexatriacontane (C36)	CAS:630-06-8
n-Decane (C10)	CAS:124-18-5	n-Tetracosane (C24)	CAS:646-31-1	n-Octatriacontane (C38)	CAS:7194-85-6
n-Dodecane (C12)	CAS:112-40-3	n-Hexacosane (C26)	CAS:630-01-3	n-Tetracontane (C40)	CAS:4181-95-7
n-Tetradecane (C14)	CAS:629-59-4	n-Octacosane (C28)	CAS:630-02-4		
n-Hexadecane (C16)	CAS:544-76-3	n-Triacontane (C30)	CAS:638-68-6		
Solvent : Iso-Octane	ampoule 1 ml	50 µg/ml	Ref : F062015	Price : GF	
	ampoule 5 ml	50 µg/ml	Ref : F109821	Price : BJA	

Calibration Mixture Oils (A+B) Set					(ISO 9377-2-Mod)
Ampoule 1 Blank	n-Hexane		Ampoule 4 Mineral Oil	CAS : 8042-47-5 + Diesel Oil	0.6 mg/ml each
Ampoule 2 Mineral Oil	CAS : 8042-47-5 + Diesel Oil	0.2 mg/ml each	Ampoule 5 Mineral Oil	CAS : 8042-47-5 + Diesel Oil	0.8 mg/ml each
Ampoule 3 Mineral Oil	CAS : 8042-47-5 + Diesel Oil	0.4 mg/ml each	Ampoule 6 Mineral Oil	CAS : 8042-47-5 + Diesel Oil	1.0 mg/ml each
Solvent : n-Hexane	6 ampoules X 1 ml		Ref : F061955	Price : EF	
	6 ampoules X 5 ml		Ref : F109681	Price : BDA	

Calibration Mix Solution 1 - 17 components					(ISO 3924:1999)
n-Hexane (C6)	CAS:110-54-3	600 ug/ml	n-Octadecane (C18)	CAS:593-45-3	500 ug/ml
n-Heptane (C7)	CAS:142-82-5	600 ug/ml	n-Eicosane (C20)	CAS:112-95-8	200 ug/ml
n-Octane (C8)	CAS:111-65-9	800 ug/ml	n-Tetracosane (C24)	CAS:646-31-1	200 ug/ml
n-Nonane (C9)	CAS:111-84-2	800 ug/ml	n-Octacosane (C28)	CAS:630-02-4	50 ug/ml
n-Decane (C10)	CAS:124-18-5	1 200 ug/ml	n-Dotriacontane (C32)	CAS:544-85-4	100 ug/ml
n-Undecane (C11)	CAS:1120-21-4	1 200 ug/ml	n-Hexatriacontane (C36)	CAS:630-06-8	100 ug/ml
n-Dodecane (C12)	CAS:112-40-3	1 200 ug/ml	n-Tetracontane (C40)	CAS:4181-95-7	100 ug/ml
n-Tetradecane (C14)	CAS:629-59-4	1 200 ug/ml	n-Tetratetracontane (C44)	CAS:7098-22-8	100 ug/ml
n-Hexadecane (C16)	CAS:544-76-3	1 000 ug/ml			
Solvent : Chloroform	ampoule 1 ml		Ref : F110061	Price : GF	

Calibration Mix Solution 2 – 20 components					(ISO 3924:1999)
n-Pentane (C5)	CAS:109-66-0	n-Dodecane (C12)	CAS:112-40-3	n-Tetracosane (C24)	CAS:646-31-1
n-Hexane (C6)	CAS:110-54-3	n-Tetradecane (C14)	CAS:629-59-4	n-Octacosane (C28)	CAS:630-02-4
n-Heptane (C7)	CAS:142-82-5	n-Pentadecane (C15)	CAS:629-62-9	n-Dotriacontane (C32)	CAS:544-85-4
n-Octane (C8)	CAS:111-65-9	n-Hexadecane (C16)	CAS:544-76-3	n-Hexatriacontane (C36)	CAS:630-06-8
n-Nonane (C9)	CAS:111-84-2	n-Heptadecane (C17)	CAS:629-78-7	n-Tetracontane (C40)	CAS:4181-95-7
n-Decane (C10)	CAS:124-18-5	n-Octadecane (C18)	CAS:593-45-3	n-Tetratetracontane (C44)	CAS:7098-22-8
n-Undecane (C11)	CAS:1120-21-4	n-Eicosane (C20)	CAS:112-95-8		
Solvent : Chloroform	ampoule 1 ml	0.1% (w/w)	Ref : F110071	Price : HF	

Calibration Mixture C5-C44 - 18 components					(ASTM D2887)
n-Pentane (C5)	CAS:109-66-0	n-Undecane (C11)	CAS:1120-21-4	n-Tetracosane (C24)	CAS:646-31-1
n-Hexane (C6)	CAS:110-54-3	n-Dodecane (C12)	CAS:112-40-3	n-Octacosane (C28)	CAS:630-02-4
n-Heptane (C7)	CAS:142-82-5	n-Tetradecane (C14)	CAS:629-59-4	n-Dotriacontane (C32)	CAS:544-85-4
n-Octane (C8)	CAS:111-65-9	n-Hexadecane (C16)	CAS:544-76-3	n-Hexatriacontane (C36)	CAS:630-06-8
n-Nonane (C9)	CAS:111-84-2	n-Octadecane (C18)	CAS:593-45-3	n-Tetracontane (C40)	CAS:4181-95-7
n-Decane (C10)	CAS:124-18-5	n-Eicosane (C20)	CAS:112-95-8	n-Tetratetracontane (C44)	CAS:7098-22-8
Solvent : Carbon disulfide	ampoule 1 ml	100 µg/ml	Ref : F110231	Price : GF	

Column Resolution Test Mix C20-C24 - 2 components					(ASTM D5442)
n-Eicosane (C20)	CAS:112-95-8	n-Tetracosane (C24)	CAS:646-31-1		
Solvent : Cyclohexane	ampoule 1 ml	0.05 % wt	Ref : F112011	Price : EF	

Column Test Mixture - 17 components					(ASTM D2887)
n-Hexane (C6)	CAS:110-54-3	6 % wt	n-Octadecane (C18)	CAS:593-45-3	5 % wt
n-Heptane (C7)	CAS:142-82-5	6 % wt	n-Eicosane (C20)	CAS:112-95-8	2 % wt
n-Octane (C8)	CAS:111-65-9	8 % wt	n-Tetracosane (C24)	CAS:646-31-1	2 % wt
n-Nonane (C9)	CAS:111-84-2	8 % wt	n-Octacosane (C28)	CAS:630-02-4	1 % wt
n-Decane (C10)	CAS:124-18-5	12 % wt	n-Dotriacontane (C32)	CAS:544-85-4	1 % wt
n-Undecane (C11)	CAS:1120-21-4	12 % wt	n-Hexatriacontane (C36)	CAS:630-06-8	1 % wt
n-Dodecane (C12)	CAS:112-40-3	12 % wt	n-Tetracontane (C40)	CAS:4181-95-7	1 % wt
n-Tetradecane (C14)	CAS:629-59-4	12 % wt	n-Tetratetracontane (C44)	CAS:7098-22-8	1 % wt
n-Hexadecane (C16)	CAS:544-76-3	10 % wt			
	ampoule 1 ml		Ref : F112171	Price : EEA	

Column Test Mixture -2 components					(ASTM D2887- ISO 3924:1999)
n-Hexadecane (C16)	CAS:544-76-3	n-Octadecane (C18)	CAS:593-45-3		
Solvent : n-Octane	ampoule 1 ml	1% (w/v)	Ref : F110051	Price : FA	

Column Resolution Test Mix C16-C18 - 2 components**(ASTM D5307)**

n-Hexadecane (C16)	CAS:544-76-3	n-Octadecane (C18)	CAS:593-45-3		
Solvent : n-Octane		ampoule 1 ml	1 % wt	Ref : F111951	Price : EF

Diesel range Mixture - 10 components

n-Decane (C10)	CAS:124-18-5	n-Octadecane (C18)	CAS:593-45-3	n-Hexacosane (C26)	CAS:630-01-3
n-Dodecane (C12)	CAS:112-40-3	n-Eicosane (C20)	CAS:112-95-8	n-Octacosane (C28)	CAS:630-02-4
n-Tetradecane (C14)	CAS:629-59-4	n-Docosane (C22)	CAS:629-97-0		
n-Hexadecane (C16)	CAS:544-76-3	n-Tetracosane (C24)	CAS:646-31-1		
Solvent : n-Hexane		ampoule 1 ml	1 000 µg/ml	Ref : F112371	Price : FA

Extraction Solvent Stock Solution - 2 components**(ISO 9377-2-Mod)**

n-Decane (C10)	CAS:124-18-5	20 µl/l	n-Tetracontane (C40)	CAS:4181-95-7	20 mg/l
Solvent : n-Hexane		ampoule 1 ml		Ref : P858570	Price : DA
		ampoule 10 ml		Ref : F109741	Price : FE
Solvent : n-Hexane/Petroleum ether (1/1)		ampoule 1 ml		Ref : P858530	Price : DA
		ampoule 10 ml		Ref : F109751	Price : FE

Florisol Cartridge QC Standard Mixture - 2 components**(ISO 9377-2-Mod)**

Mineral Oil	CAS:8042-47-5	Diesel Oil	N/A (id : 381)		
Solvent : n-Hexane		ampoule 1 ml	1 000 µg/ml each	Ref : P858090	Price : DA
		ampoule 5 ml	1 000 µg/ml each	Ref : F109651	Price : FF

Fuel Oil Degradation/Retention Time Mix C17-C20 – 4 components**(ISO 3924:1999 - ASTM D2887)**

n-Heptadecane (C17)	CAS:629-78-7	Phytane (C20)	CAS:638-36-8	Pristane (C19)	CAS:1921-70-6
n-Octadecane (C18)	CAS:593-45-3				
Solvent : Dichloromethane/Carbon disulfide (1/1)		ampoule 1 ml	2 000 µg/ml	Ref : F110261	Price : DF

Hydrocarbons Standard Solution - 4 components

n-Decane (C10)	CAS:124-18-5	n-Triacontane (C30)	CAS:638-68-6	n-Tetracontane (C40)	CAS:4181-95-7
n-Eicosane (C20)	CAS:112-95-8				
Solvent : n-Hexane		ampoule 1 ml	100 µg/ml	Ref : F112221	Price : DA
Solvent : n-Hexane/Petroleum ether (1/1)		ampoule 1 ml	100 µg/ml	Ref : F062305	Price : DA

Hydrocarbons Standard Solution - 6 components

n-Decane (C10)	CAS:124-18-5	n-Hexadecane (C16)	CAS:544-76-3	n-Eicosane (C20)	CAS:112-95-8
n-Tetradecane (C14)	CAS:629-59-4	n-Octadecane (C18)	CAS:593-45-3	n-Tetracontane (C40)	CAS:4181-95-7
Solvent : n-Hexane		ampoule 1 ml	50 µg/ml	Ref : F112241	Price : DA

Hydrocarbons Standard Solution - 31 components

n-Decane (C10)	CAS:124-18-5	n-Docosane (C22)	CAS:629-97-0	n-Tritriacontane (C33)	CAS:630-05-7
n-Undecane (C11)	CAS:1120-21-4	n-Tricosane (C23)	CAS:638-67-5	n-Octadecane (C18)	CAS:593-45-3
n-Dodecane (C12)	CAS:112-40-3	n-Tetracosane (C24)	CAS:646-31-1	n-Tetatriacontane (C34)	CAS:14167-59-0
n-Tridecane (C13)	CAS:629-50-5	n-Pentacosane (C25)	CAS:629-99-2	n-Pentatriacontane (C35)	CAS:630-07-9
n-Tetradecane (C14)	CAS:629-59-4	n-Hexacosane (C26)	CAS:630-01-3	n-Heptatriacontane (C37)	CAS:7194-84-5
n-Pentadecane (C15)	CAS:629-62-9	n-Octacosane (C28)	CAS:630-02-4	n-Nonatriacontane (C39)	CAS:7194-86-7
n-Hexadecane (C16)	CAS:544-76-3	n-Nonacosane (C29)	CAS:630-03-5	n-Hexatriacontane (C36)	CAS:630-06-8
n-Heptadecane (C17)	CAS:629-78-7	n-Triacontane (C30)	CAS:638-68-6	n-Octatriacontane (C38)	CAS:7194-85-6
n-Nonadecane (C19)	CAS:629-92-5	n-Hentriacontane (C31)	CAS:630-04-6	n-Tetracontane (C40)	CAS:4181-95-7
n-Eicosane (C20)	CAS:112-95-8	n-Dotriacontane (C32)	CAS:544-85-4		
n-Heneicosane (C21)	CAS:629-94-7	n-Heptacosane (C27)	CAS:593-49-7		
Solvent : n-Hexane		ampoule 1 ml	100 µg/ml	Ref : F112231	Price : KF

Hydrocarbon Window Defining Standard C8-C40 – 35 components**(ISO 3924:1999 ; ASTM D2887)**

n-Octane (C8)	CAS:111-65-9	Pristane (C19)	CAS:1921-70-6	n-Triacontane (C30)	CAS:638-68-6
n-Nonane (C9)	CAS:111-84-2	n-Eicosane (C20)	CAS:112-95-8	n-Hentriacontane (C31)	CAS:630-04-6
n-Decane (C10)	CAS:124-18-5	Phytane (C20)	CAS:638-36-8	n-Dotriacontane (C32)	CAS:544-85-4
n-Undecane (C11)	CAS:1120-21-4	n-Heneicosane (C21)	CAS:629-94-7	n-Tritriacontane (C33)	CAS:630-05-7
n-Dodecane (C12)	CAS:112-40-3	n-Docosane (C22)	CAS:629-97-0	n-Tetatriacontane (C34)	CAS:14167-59-0
n-Tridecane (C13)	CAS:629-50-5	n-Tricosane (C23)	CAS:638-67-5	n-Pentatriacontane (C35)	CAS:630-07-9
n-Tetradecane (C14)	CAS:629-59-4	n-Tetracosane (C24)	CAS:646-31-1	n-Hexatriacontane (C36)	CAS:630-06-8
n-Pentadecane (C15)	CAS:629-62-9	n-Pentacosane (C25)	CAS:629-99-2	n-Heptatriacontane (C37)	CAS:7194-84-5
n-Hexadecane (C16)	CAS:544-76-3	n-Hexacosane (C26)	CAS:630-01-3	n-Octatriacontane (C38)	CAS:7194-85-6
n-Heptadecane (C17)	CAS:629-78-7	n-Heptacosane (C27)	CAS:593-49-7	n-Nonatriacontane (C39)	CAS:7194-86-7
n-Octadecane (C18)	CAS:593-45-3	n-Octacosane (C28)	CAS:630-02-4	n-Tetracontane (C40)	CAS:4181-95-7
n-Nonadecane (C19)	CAS:629-92-5	n-Nonacosane (C29)	CAS:630-03-5		
Solvent : Carbon disulfide/Dichloromethane (3/1)		ampoule 1 ml	500 µg/ml	Ref : F061935	Price : KF
		ampoule 5 ml	500 µg/ml	Ref : F110241	Price : DJA

Mixture of Mineral Oils (A+B) (Diesel Fuel / Lubricating Oil) : 1/1(neat) - 2 components**(ISO 9377-2-Mod)**

Mineral Oil 50%	CAS:8042-47-5	Diesel Oil 50%	N/A (id : 381)		
		ampoule 1 ml		Ref : F109631	Price : DF
		ampoule 5 ml		Ref : F061945	Price : KJ

Petrochemical Internal Standard C14-C17 - 4 components**(ASTM D5307)**

1-Tetradecene	CAS:1120-36-1	n-Hexadecane (C16)	CAS:544-76-3	n-Heptadecane (C17)	CAS:629-78-7
n-Pentadecane (C15)	CAS:629-62-9				
		bottle 40 ml	25 % wt	Ref : F111941	Price : EEA

Petrochemical Quantitative Standard C10-C44 - 14 components**(ASTM D5307)**

n-Decane (C10)	CAS:124-18-5	n-Hexadecane (C16)	CAS:544-76-3	n-Dotriacontane (C32)	CAS:544-85-4
n-Undecane (C11)	CAS:1120-21-4	n-Heptadecane (C17)	CAS:629-78-7	n-Hexatriacontane (C36)	CAS:630-06-8
n-Dodecane (C12)	CAS:112-40-3	n-Octadecane (C18)	CAS:593-45-3	n-Tetracontane (C40)	CAS:4181-95-7
n-Tridecane (C13)	CAS:629-50-5	n-Eicosane (C20)	CAS:112-95-8	n-Tetratetracontane (C44)	CAS:7098-22-8
n-Tetradecane (C14)	CAS:629-59-4	n-Tetracosane (C24)	CAS:646-31-1		
n-Pentadecane (C15)	CAS:629-62-9	n-Octacosane (C28)	CAS:630-02-4		
Solvent : Carbon disulfide		ampoule 1 ml	5 000 µg/ml	Ref : F111931	Price : BGA

Quantitative Linearity Wax Mixture C16-C44 - 12 components**(ASTM 5442)**

n-Hexadecane (C16)	CAS:544-76-3	n-Tetracosane (C24)	CAS:646-31-1	n-Dotriacontane (C32)	CAS:544-85-4
n-Octadecane (C18)	CAS:593-45-3	n-Hexacosane (C26)	CAS:630-01-3	n-Hexatriacontane (C36)	CAS:630-06-8
n-Eicosane (C20)	CAS:112-95-8	n-Octacosane (C28)	CAS:630-02-4	n-Tetracontane (C40)	CAS:4181-95-7
n-Docosane (C22)	CAS:629-97-0	n-Triacontane (C30)	CAS:638-68-6	n-Tetratetracontane (C44)	CAS:7098-22-8
Solvent : Cyclohexane		ampoule 1 ml	0.01 % wt	Ref : F111991	Price : FA

Quality Control Standard of Mineral Oils - 2 components**(ISO 9377-2-Mod)**

Mineral Oil	CAS:8042-47-5	Diesel Oil	N/A (id : 381)		
Solvent : Acetone		ampoule 1 ml	500 µg/ml each	Ref : P858050	Price : DF
		ampoule 5 ml	500 µg/ml each	Ref : F109641	Price : KJ
		ampoule 1 ml	1 000 µg/ml each	Ref : F061965	Price : DF
		ampoule 5 ml	1 000 µg/ml each	Ref : F145091	Price : KJ
Solvent : n-Hexane		ampoule 1 ml	5 000 µg/ml each	Ref : P858010	Price : DF
		ampoule 5 ml	5 000 µg/ml each	Ref : F109661	Price : KJ
		ampoule 1 ml	10 000 µg/ml each	Ref : F061975	Price : DF
		ampoule 5 ml	10 000 µg/ml each	Ref : F109671	Price : KJ

Retention Time Wax Mixture C16-C44 - 12 components**(ASTM D5442)**

n-Hexadecane (C16)	CAS:544-76-3	n-Tetracosane (C24)	CAS:646-31-1	n-Dotriacontane (C32)	CAS:544-85-4
n-Octadecane (C18)	CAS:593-45-3	n-Hexacosane (C26)	CAS:630-01-3	n-Hexatriacontane (C36)	CAS:630-06-8
n-Eicosane (C20)	CAS:112-95-8	n-Octacosane (C28)	CAS:630-02-4	n-Tetracontane (C40)	CAS:4181-95-7
n-Docosane (C22)	CAS:629-97-0	n-Triacontane (C30)	CAS:638-68-6	n-Tetratetracontane (C44)	CAS:7098-22-8
Solvent : Cyclohexane		ampoule 1 ml	8.3 % wt	Ref : F112001	Price : DEA

BTEX mixtures**ACs Standard Solution 6 components****(ISO11423-2:1997)**

Benzene	CAS:71-43-2	Toluene	CAS:108-88-3	m-Xylene	CAS:108-38-3
Ethylbenzene	CAS:100-41-4	o-Xylene	CAS:95-47-6	p-Xylene	CAS:106-42-3
Solvent : Acetone		ampoule 5 ml	50 ml/l	Ref : F129001	Price : EF
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F129231	Price : CD
		ampoule 10 ml	1 µg/ml	Ref : F122271	Price : DF

BTEX Standard Solution 6 components**(ISO 9377-2-Mod)**

Benzene	CAS:71-43-2	Toluene	CAS:108-88-3	m-Xylene	CAS:108-38-3
Ethylbenzene	CAS:100-41-4	o-Xylene	CAS:95-47-6	p-Xylene	CAS:106-42-3
Solvent : n-Hexane		ampoule 1 ml	200 µg/ml	Ref : F061985	Price : DF
		ampoule 5 ml	200 µg/ml	Ref : F109801	Price : KJ
		ampoule 5 ml	2000 µg/ml	Ref : F061995	Price : DF
		ampoule 5 ml	2000 µg/ml	Ref : F109811	Price : KJ

Valve Timing Calibration Mix (without Internal Standard) - 4 components**(ASTM D5580)**

Benzene	CAS:71-43-2	5%wt	Ethylbenzene	CAS:100-41-4	10%wt
Toluene	CAS:108-88-3	5%wt	o-Xylene	CAS:95-47-6	10%wt
Solvent : Iso-Octane 70% wt		ampoule 1 ml		Ref : F062235	Price : EF
		ampoule 10 ml		Ref : F112121	Price : BEA

Valve Timing Calibration Mix (with Internal Standard) – 5 components**(ASTM D5580)**

Benzene	CAS:71-43-2	4.5 %wt	Toluene	CAS:108-88-3	4.5 %wt
Ethylbenzene	CAS:100-41-4	9.0 %wt	o-Xylene	CAS:95-47-6	9.0 %wt
2-Hexanone	CAS:591-78-6	10.0 %wt			
Solvent : Iso-Octane 63 %wt		ampoule 1 ml		Ref : F062245	Price : EF
		ampoule 10 ml		Ref : F112131	Price : BEA

n-Alkanes + BTEX mixtures**ACs Standard Solution - 10 components**

n-Heptane (C7)	CAS:142-82-5	Benzene	CAS:71-43-2	o-Xylene	CAS:95-47-6
n-Decane (C10)	CAS:124-18-5	Ethylbenzene	CAS:100-41-4	p-Xylene	CAS:106-42-3
n-Eicosane (C20)	CAS:112-95-8	Toluene	CAS:108-88-3		
n-Tetracontane (C40)	CAS:111-84-2	m-Xylene	CAS:108-38-3		
Solvent : n-Pentane		ampoule 1 ml	50 µg/ml	Ref : F062345	Price : CF
		ampoule 5 ml	50 µg/ml	Ref : F112291	Price : GF

ACs Standard Solution - 12 components

n-Heptane (C7)	CAS:142-82-5	n-Eicosane (C20)	CAS:112-95-8	Toluene	CAS:108-88-3
n-Octane (C8)	CAS:111-65-9	n-Tetracontane (C40)	CAS:4181-95-7	m-Xylene	CAS:108-38-3
n-Nonane (C9)	CAS:111-84-2	Benzene	CAS:71-43-2	o-Xylene	CAS:95-47-6
n-Decane (C10)	CAS:124-18-5	Ethylbenzene	CAS:100-41-4	p-Xylene	CAS:106-42-3
Solvent : n-Pentane		ampoule 1 ml	50 µg/ml	Ref : F062335	Price : CF
		ampoule 5 ml	50 µg/ml	Ref : F112321	Price : HA

n-Alkane Standard Solution C7+C8+C9+C10+C40+ BTEX - 25 components

n-Heptane (C7)	CAS:142-82-5	n-Docosane (C22)	CAS:629-97-0	n-Tetracontane (C40)	CAS:4181-95-7
n-Octane (C8)	CAS:111-65-9	n-Tetracosane (C24)	CAS:646-31-1	Benzene	CAS:71-43-2
n-Nonane (C9)	CAS:111-84-2	n-Hexacosane (C26)	CAS:630-01-3	Ethylbenzene	CAS:100-41-4
n-Decane (C10)	CAS:124-18-5	n-Octacosane (C28)	CAS:630-02-4	Toluene	CAS:108-88-3
n-Dodecane (C12)	CAS:112-40-3	n-Triacontane (C30)	CAS:638-68-6	m-Xylene	CAS:108-38-3
n-Tetradecane (C14)	CAS:629-59-4	n-Dotriacontane (C32)	CAS:544-85-4	o-Xylene	CAS:95-47-6
n-Hexadecane (C16)	CAS:544-76-3	n-Tetatriacontane (C34)	CAS:14167-59-0	p-Xylene	CAS:106-42-3
n-Octadecane (C18)	CAS:593-45-3	n-Hexatriacontane (C36)	CAS:630-06-8		
n-Eicosane (C20)	CAS:112-95-8	n-Octatriacontane (C38)	CAS:7194-85-6		
Solvent : n-Hexane (ISO 9377-2-Mod)		ampoule 1 ml	50 µg/ml	Ref : F062025	Price : GA
		ampoule 5 ml	50 µg/ml	Ref : F109851	Price : BHA
Solvent : n-Pentane		ampoule 1 ml	50 µg/ml	Ref : F062355	Price : GA
		ampoule 5 ml	50 µg/ml	Ref : F112331	Price : BHA

Column Resolution Check Mix - 13 components**(ASTM E1387-01)**

n-Hexane (C6)	CAS:110-54-3	n-Hexadecane (C16)	CAS:544-76-3	2-Ethyltoluene	CAS:611-14-3
n-Octane (C8)	CAS:111-65-9	n-Octadecane (C18)	CAS:593-45-3	3-Ethyltoluene	CAS:620-14-4
n-Decane (C10)	CAS:124-18-5	n-Eicosane (C20)	CAS:112-95-8	p-Xylene	CAS:106-42-3
n-Dodecane (C12)	CAS:112-40-3	Toluene	CAS:108-88-3		
n-Tetradecane (C14)	CAS:629-59-4	1,2,4-Trimethylbenzene	CAS:95-63-6		
Solvent : Dichloromethane		ampoule 1 ml	2 000 µg/ml	Ref : F128941	Price : FA

Hydrocarbons Calibration Standard - 16 components**(ASTM D3710)**

n-Pentane (C5)	CAS:109-66-0	8.0 % wt	n-Pentadecane (C15)	CAS:629-62-9	2.5 % wt
n-Hexane (C6)	CAS:110-54-3	6.0 % wt	n-Butylbenzene	CAS:104-51-8	4.0 % wt
n-Heptane (C7)	CAS:142-82-5	10.0 % wt	2,4-Dimethylpentane	CAS:108-08-7	6.0 % wt
n-Octane (C8)	CAS:111-65-9	6.0 % wt	2-Methylbutane	CAS:78-78-4	10.0 % wt
n-Decane (C10)	CAS:124-18-5	4.0 % wt	2-Methylpentane	CAS:107-83-5	6.0 % wt
n-Dodecane (C12)	CAS:112-40-3	3.5 % wt	n-Propylbenzene	CAS:103-65-1	5.0 % wt
n-Tridecane (C13)	CAS:629-50-5	2.5 % wt	Toluene	CAS:108-88-3	11.0 % wt
n-Tetradecane (C14)	CAS:629-59-4	2.5 % wt	p-Xylene	CAS:106-42-3	13.0 % wt
		ampoule 1 ml		Ref : F110421	Price : BJA

p-Xylene Impurity Standard - 10 components**(ASTM D3798)**

Benzene	CAS:71-43-2	0.15 % wt	n-Propylbenzene	CAS:103-65-1	0.15 % wt
Cumene	CAS:98-82-8	0.15 % wt	Toluene	CAS:108-88-3	0.15 % wt
Ethylbenzene	CAS:100-41-4	0.15 % wt	m-Xylene	CAS:108-38-3	0.15 % wt
n-Octane (C8)	CAS:111-65-9	0.15 % wt	o-Xylene	CAS:95-47-6	0.15 % wt
n-Pentane (C5)	CAS:109-66-0	0.15 % wt	p-Xylene	CAS:106-42-3	98.65 % wt
		ampoule 1 ml		Ref : F110431	Price : HF

Nitrogen containing compounds (Amines, Aniline derivatives, Nitroaromatics)

NCC Standard Solution - 3 components

N-Nitroso-dimethylamine	CAS:62-75-9	N-Nitroso-diphenylamine	CAS:86-30-6	N-Nitroso-di-n-propylamine	CAS:621-64-7
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F116691	Price : DH

NCC Standard Solution - 9 components

N-Nitrosodibutylamine	CAS:924-16-3	N-Nitrosomethylethylamine	CAS:10595-95-6	N-Nitrosomorpholine	CAS:59-89-2
N-Nitrosodiethylamine	CAS:55-18-5	N-Nitroso-diphenylamine	CAS:86-30-6	N-Nitrosopiperidine	CAS:100-75-4
N-Nitroso-dimethylamine	CAS:62-75-9	N-Nitroso-di-n-propylamine	CAS:621-64-7	N-Nitrosopyrrolidine	CAS:930-55-2
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119341	Price : EF

NCC Standard Solution - 12 components

Benzidine	CAS:92-87-5	250 µg/ml	3,3'-Dimethylbenzidine	CAS:119-93-7	350 µg/ml
Benzoylprop ethyl	CAS:22212-55-1	350 µg/ml	Diuron	CAS:330-54-1	450 µg/ml
Carbaryl	CAS:63-25-2	1 000 µg/ml	Linuron	CAS:330-55-2	1 300 µg/ml
o-Chlorophenyl thiourea	CAS:5344-82-1	750 µg/ml	Monuron (CMU)	CAS:150-68-5	400 µg/ml
3,3'-Dichlorobenzidine	CAS:91-94-1	250 µg/ml	Rotenone	CAS:83-79-4	3 200 µg/ml
3,3'-Dimethoxybenzidine	CAS:119-90-4	750 µg/ml	Siduron	CAS:1982-49-6	450 µg/ml
Solvent : Acetonitrile/Methanol (1/1)		ampoule 1 ml		Ref : F119961	Price : FH

NCC Standard Solution - 19 components

Aniline	CAS:62-53-3	2-Chloro-4,6-dinitroaniline	CAS:3531-19-9	2-Nitroaniline	CAS:88-74-4
4-Bromoaniline	CAS:106-40-1	2-Chloro-4-nitroaniline	CAS:121-87-9	3-Nitroaniline	CAS:99-09-2
2-Bromo-6-chloro-4-nitroaniline	CAS:99-29-6	4-Chloro-2-nitroaniline	CAS:89-63-4	4-Nitroaniline	CAS:100-01-6
2-Bromo-4,6-dinitroaniline	CAS:1817-73-8	2,6-Dibromo-4-nitroaniline	CAS:827-94-1	2,4,6-Trichloroaniline	CAS:634-93-5
2-Chloroaniline	CAS:95-51-2	3,4-Dichloroaniline	CAS:95-76-1	2,4,5-Trichloroaniline	CAS:636-30-6
3-Chloroaniline	CAS:108-42-9	Dicloran (CNA)	CAS:99-30-9		
4-Chloroaniline	CAS:106-47-8	2,4-Dinitroaniline	CAS:97-02-9		
Solvent : Toluene		ampoule 1 ml	1 000 µg/ml	Ref : F117631	Price : FJ

NCC Standard Solution - 32 components

4-Aminoazobenzene	CAS:60-09-3	p-Dimethylaminoazobenzene	CAS:60-11-7	5-Nitro-o-anisidine	CAS:99-59-2
4-Aminobiphenyl	CAS:92-67-1	Diphenylamine	CAS:122-39-4	5-Nitro-o-toluidine	CAS:99-55-8
3-Amino-9-ethylcarbazole	CAS:132-32-1	1,2-Diphenylhydrazine	CAS:122-66-7	4,4'-Oxydianiline	CAS:101-80-4
Aniline	CAS:62-53-3	4,4'-Methylene-bis(2-chloroanil)	CAS:101-14-4	Phenacetin	CAS:62-44-2
4-Chloroaniline	CAS:106-47-8	4,4'-Methylene bis(N,Ndimethylaniline)	CAS:101-61-1	1,4-Phenylenediamine	CAS:106-50-3
5-Chloro-2-methylaniline	CAS:95-79-4	1-Naphthylamine	CAS:134-32-7	2-Picoline	CAS:109-06-8
3-(Chloromethyl) pyridine hydrochloride	CAS:6959-48-4	2-Naphthylamine	CAS:91-59-8	Propyzamide	CAS:23950-58-5
4-Chloro-1,2-phenylenediamine	CAS:95-83-0	Nicotine	CAS:54-11-5	Pyridine	CAS:110-86-1
4-Chloro-1,3-phenylenediamine	CAS:5131-60-2	2-Nitroaniline	CAS:88-74-4	o-Toluidine	CAS:95-53-4
p-Cresidine	CAS:120-71-8	3-Nitroaniline	CAS:99-09-2	2,4,5-Trimethylaniline	CAS:137-17-7
2,4-Diaminotoluene	CAS:95-80-7	4-Nitroaniline	CAS:100-01-6		
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119331	Price : GC

RCRA Analytes - 6 components

1,4-Dinitrobenzene	CAS:100-25-4	2,6-Dinitrotoluene	CAS:606-20-2	Nitrobenzene	CAS:98-95-3
2,4-Dinitrotoluene	CAS:121-14-2	1,4-Naphthoquinone	CAS:130-15-4	Pentachloronitrobenzene	CAS:82-68-8
Solvent : Iso-Octane/Toluene (1/1)		ampoule 1 ml	1 000 µg/ml	Ref : F117371	Price : DJ

Pesticides

Carbamates Laboratory Performances Check Solution - 4 components

Aldicarb-sulfoxide	CAS:1646-87-3	100 µg/ml	Carbofuran-3-hydroxy	CAS:16655-82-6	2 µg/ml
BDMC	CAS:672-99-1	10 µg/ml	Methiocarb	CAS:2032-65-7	20 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F114181	Price : DF

Carbamates / NCC / OPP Standard Solution - 11 components

Aldicarb	CAS:116-06-3	Metazachlor	CAS:67129-08-2	Parathion-methyl	CAS:298-00-0
Chlorfenvinphos	CAS:470-90-6	Metolachlor	CAS:51218-45-2	Terbutryn	CAS:886-50-0
Chlorpropham	CAS:101-21-3	Prometryn	CAS:7287-19-6	Vinclozolin	CAS:50471-44-8
Crimidine	CAS:535-89-7	Parathion-ethyl	CAS:56-38-2		
Solvent : Methanol		ampoule 1 ml	10 µg/ml	Ref : F128801	Price : DF

Carbamates Standard Solution - 2 components

Napropamide	CAS:15299-99-7	Propanil	CAS:709-98-8		
Solvent : Acetonitrile/Acetone (1/1)		ampoule 1 ml	100 µg/ml	Ref : F115941	Price : DF

Carbamates Standard Solution - 6 components**(EPA 634)**

Butylate	CAS:2008-41-5	EPTC	CAS:759-94-4	Pebulate	CAS:1114-71-2
Cycloate	CAS:1134-23-2	Molinat	CAS:2212-67-1	Vernolate	CAS:1929-77-7
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F115951	Price : EF

Carbamates Standard Solution - 10 components**(EPA 531.1)**

Aldicarb	CAS:116-06-3	Carbaryl	CAS:63-25-2	Methomyl	CAS:16752-77-5
Aldicarb-sulfone	CAS:1646-88-4	Carbofuran	CAS:1563-66-2	Oxamyl	CAS:23135-22-0
Aldicarb-sulfoxide	CAS:1646-87-3	Carbofuran-3-hydroxy	CAS:16655-82-6		
Baygon	CAS:114-26-1	Methiocarb	CAS:2032-65-7		
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F114161	Price : FA

Carbamates Standard Solution - 10 components**(EPA 8318)**

Aldicarb	CAS:116-06-3	Dioxacarb	CAS:6988-21-2	Promecarb	CAS:2631-37-0
Aldicarb-sulfone	CAS:1646-88-4	3-Hydroxycarbofuran	CAS:16655-82-6	Propoxur	CAS:114-26-1
Carbaryl	CAS:63-25-2	Methiocarb	CAS:2032-65-7		
Carbofuran	CAS:1563-66-2	Methomyl	CAS:16752-77-5		
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F119921	Price : GA

Carbamates Standard Solution - 21 components**(EPA 632)**

Aminocarb	CAS:2032-59-9	Fenuron TCA	CAS:4482-55-7	Monuron TCA	CAS:140-41-0
Barban	CAS:101-27-9	Fluometuron	CAS:2164-17-2	Neburon	CAS:555-37-3
Carbaryl	CAS:63-25-2	Linuron	CAS:330-55-2	Oxamyl	CAS:23135-22-0
Carbofuran	CAS:1563-66-2	Methiocarb	CAS:2032-65-7	Propham	CAS:122-42-9
Chlorpropham	CAS:101-21-3	Methomyl	CAS:16752-77-5	Propoxur	CAS:114-26-1
Diuron	CAS:330-54-1	Mexacarbate	CAS:315-18-4	Siduron	CAS:1982-49-6
Fenuron	CAS:101-42-8	Monuron	CAS:150-68-5	Swep	CAS:1918-18-9
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F115921	Price : EE

Carbamates Mixture - 29 components**(EPA 8321A)**

Aldicarb	CAS:116-06-3	Carbofuran	CAS:1563-66-2	Mexacarbate	CAS:315-18-4
Aldicarb-sulfone	CAS:1646-88-4	Chloroxuron	CAS:1982-47-4	Monuron (CMU)	CAS:150-68-5
Aldicarb-sulfoxide	CAS:1646-87-3	Chlorpropham	CAS:101-21-3	Neburon	CAS:555-37-3
Aminocarb	CAS:2032-59-9	Diuron	CAS:330-54-1	Oxamyl	CAS:23135-22-0
Barban	CAS:101-27-9	Fenuron	CAS:101-42-8	Propachlor	CAS:1918-16-7
Benomyl	CAS:17804-35-2	Fluometuron	CAS:2164-17-2	Propham	CAS:122-42-9
Bromacil	CAS:314-40-9	3-Hydroxycarbofuran	CAS:16655-82-6	Propoxur	CAS:114-26-1
Bendiocarb	CAS:22781-23-3	Linuron	CAS:330-55-2	Siduron	CAS:1982-49-6
Carbaryl	CAS:63-25-2	Methiocarb	CAS:2032-65-7	Tebuthiuron	CAS:34014-18-1
Carbendazim	CAS:10605-21-7	Methomyl	CAS:16752-77-5		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F119951	Price : GF

Chlorinated Herbicides Acids Mixture - 7 components**(EPA 515)**

Acifluorfen	CAS:50594-66-6	100 µg/ml	Dicamba	CAS:1918-00-9	300 µg/ml
Bentazone	CAS:25057-89-0	1 000 µg/ml	Picloram	CAS:1918-02-1	300 µg/ml
2,4-D	CAS:94-75-7	100 µg/ml	2,4,5-TP	CAS:93-72-1	100 µg/ml
2,4-DB	CAS:94-82-6	1 000 µg/ml			
Solvent : Methanol		ampoule 1 ml		Ref : F113851	Price : EA

Chlorinated Herbicides Mixture - 6 components**(EPA 515)**

3,5-Dichlorobenzoic acid	CAS:51-36-5	500 µg/ml	Pentachlorophenol	CAS:87-86-5	100 µg/ml
Dichlorprop	CAS:120-36-5	100 µg/ml	2,4,5-T	CAS:93-76-5	100 µg/ml
Dinoseb	CAS:88-85-7	200 µg/ml	Tetrachloroterephthalic acid	CAS:2136-79-0	100 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F113831	Price : EA

Chlorinated Herbicides Mixture - 6 components**(EPA 515)**

DCPA methyl ester	CAS:1861-32-1	100 µg/ml	Dinoseb methyl ether	CAS:6099-79-2	200 µg/ml
3,5-Dichlorobenzoic acid methyl ester	CAS:2905-67-1	500 µg/ml	Pentachloroanisole	CAS:1825-21-4	100 µg/ml
Dichlorprop methyl ester	CAS:57153-17-0	100 µg/ml	2,4,5-T methyl ester	CAS:1928-37-6	100 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F113841	Price : FA

Chlorinated Herbicides Mixture - 7 components**(EPA 515)**

Acifluorfen methyl ester	id : 15801	200 µg/ml	Dicamba-methyl ester	CAS:6597-78-0	300 µg/ml
Bentazon methyl derivative	CAS:2225-40-3	1 000 µg/ml	Picloram methyl ester	id : 15841	300 µg/ml
2,4-D methyl ester	CAS:1928-38-7	100 µg/ml	Silvex methyl ester	CAS:4841-20-7	100 µg/ml
2,4-DB methyl ester	CAS:18625-12-2	1 000 µg/ml			
Solvent : Methanol		ampoule 1 ml		Ref : F113861	Price : GA

Chlorinated Herbicides Mixture - 10 components**(EPA 615)**

2,4-D	CAS:94-75-7	100 µg/ml	Dinoseb	CAS:88-85-7	100 µg/ml
Dalapon	CAS:75-99-0	100 µg/ml	MCPA	CAS:94-74-6	10 000 µg/ml
2,4-DB	CAS:94-82-6	100 µg/ml	Mecoprop (MCP)	CAS:7085-19-0	10 000 µg/ml
Dicamba	CAS:1918-00-9	100 µg/ml	2,4,5-TP	CAS:93-72-1	10 000 µg/ml
Dichlorprop	CAS:120-36-5	100 µg/ml	2,4,5-T	CAS:93-76-5	10 000 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F115621	Price : HF

Chlorinated Herbicides Standard Solution - 16 components**(EPA 515)**

Acifluorfen	CAS:50594-66-6	100 µg/ml	3,5-Dichlorobenzoic acid	CAS:51-36-5	100 µg/ml
Bentazone	CAS:25057-89-0	200 µg/ml	Dichlorprop	CAS:120-36-5	300 µg/ml
Chloramben	CAS:133-90-4	100 µg/ml	Dinoseb	CAS:88-85-7	200 µg/ml
2,4-D	CAS:94-75-7	200 µg/ml	4-Nitrophenol	CAS:100-02-7	100 µg/ml
Dalapon	CAS:75-99-0	1 300 µg/ml	Pentachlorophenol	CAS:87-86-5	100 µg/ml
2,4-DB	CAS:94-82-6	800 µg/ml	Picloram	CAS:1918-02-1	100 µg/ml
Tetrachloroterephthalic acid	CAS:2136-79-0	100 µg/ml	2,4,5-TP	CAS:93-72-1	100 µg/ml
Dicamba	CAS:1918-00-9	100 µg/ml	2,4,5-T	CAS:93-76-5	100 µg/ml
Solvent : Methel-tert butylether		ampoule 1 ml		Ref : F113711	Price : GA

Chlorinated Phenoxyacid Compounds Mixture - 14 components**(EPA 8321A)**

2,4-D	CAS:94-75-7	Dicamba	CAS:1918-00-9	Silvex	CAS:93-72-1
2,4-D butoxyethyl ester	CAS:1929-73-3	Dichlorprop	CAS:120-36-5	2,4,5-T	CAS:93-76-5
2,4-D 2-ethylhexyl ester	CAS:1928-43-4	Dinoseb	CAS:88-85-7	2,4,5-T-butoxyethyl ester	CAS:2545-59-7
Dalapon	CAS:75-99-0	MCPA	CAS:94-74-6	2,4,5-T n-butyl ester	CAS:93-79-8
2,4-DB	CAS:94-82-6	Mecoprop (MCP)	CAS:7085-19-0		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F119941	Price : FF

Degradation Calibration Mixture - 2 components**(EPA 508.1)**

4,4'-DDT	CAS:50-29-3	200 µg/ml	Endrin	CAS:72-20-8	100 µg/ml
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113561	Price : CF

Degradation Products Mixture - 4 components**(EPA 508.1)**

4,4'-DDD (TDE)	CAS:72-54-8	200 µg/ml	Endrin aldehyde	CAS:7421-93-4	100 µg/ml
4,4'-DDE	CAS:72-55-9	200 µg/ml	Endrin ketone	CAS:53494-70-5	100 µg/ml
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113571	Price : EA

Industrial Chemicals & Triazine Herbicides Mixture - 2 components**(EPA 8141A)**

Atrazine	CAS:1912-24-9	Simazine (CAT)	CAS:122-34-9		
Solvent : Acetone		ampoule 1 ml	200 µg/ml	Ref : F119081	Price : DA

Laboratory Performance Check Solution - 5 components**(EPA 515)**

4,4'-Dibromooctafluorobiphenyl	CAS:10386-84-2	2 500 µg/ml	Dinoseb	CAS:88-85-7	40 µg/ml
3,5-Dichlorobenzoic acid	CAS:51-36-5	6 000 µg/ml	4-Nitrophenol	CAS:100-02-7	16 000 µg/ml
2,4-Dichlorophenylacetic acid	CAS:19719-28-9	5 000 µg/ml			
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113791	Price : FA

Methylated Chlorinated Herbicides Mixture - 8 components**(EPA 8151A)**

2,4-D methyl ester	CAS:1928-38-7	Dicamba-methyl ester	CAS:6597-78-0	Silvex methyl ester	CAS:4841-20-7
Dalapon methyl ester	CAS:17640-02-7	Dichlorprop methyl ester	CAS:57153-17-0	2,4,5-T methyl ester	CAS:1928-37-6
2,4-DB methyl ester	CAS:18625-12-2	Dinoseb methyl ether	CAS:6099-79-2		
Solvent : n-Hexane		ampoule 1 ml	20 µg/ml	Ref : F119161	Price : EA

Methylated Chlorinated Herbicides Mixture - 10 components**(EPA 615)**

2,4-D methyl ester	CAS:1928-38-7	100 µg/ml	Dinoseb methyl ether	CAS:6099-79-2	100 µg/ml
Dalapon methyl ester	CAS:17640-02-7	100 µg/ml	MCPA methyl ester	CAS:2436-73-9	10 000 µg/ml
2,4-DB methyl ester	CAS:18625-12-2	100 µg/ml	Mecoprop-methyl ester	CAS:23844-56-6	10 000 µg/ml
Dicamba-methyl ester	CAS:6597-78-0	100 µg/ml	2,4,5-TP	CAS:4841-20-7	10 000 µg/ml
Dichlorprop methyl ester	CAS:57153-17-0	100 µg/ml	2,4,5-T methyl ester	CAS:1928-37-6	10 000 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F115631	Price : JF

Methyl Derivatives Solution - 5 components**(EPA 515)**

4,4'-Dibromooctafluorobiphenyl	CAS:10386-84-2	250 µg/l	Methyl 3-amino-2,5-dichlorobenzoate	CAS:7286-84-2	600 µg/l
2,4-Dichlorophenylacetic acid methyl ester	id : 15851	500 µg/l	4-Nitroanisole	CAS:100-17-4	1 600 µg/l
Dinoseb methyl ether	CAS:6099-79-2	4 µg/l			
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113811	Price : EA

Methyl Derivatives Solution - 5 components**(EPA 515)**

4,4'-Dibromooctafluorobiphenyl	CAS:10386-84-2	250 µg/ml	Dinoseb	CAS:88-85-7	4 µg/ml
3,5-Dichlorobenzoic acid	CAS:51-36-5	600 µg/ml	4-Nitrophenol	CAS:100-02-7	16 000 µg/ml
2,4-Dichlorophenylacetic acid	CAS:19719-28-9	500 µg/ml			
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113821	Price : DF

Methyl Derivatives Standard Solution - 16 components**(EPA 515)**

Acifluorfen methyl ester	id : 15801	DCPA	CAS:1861-32-1	Picloram methyl ester	id : 15841
Bentazon methyl derivative	CAS:2225-40-3	Dicamba-methyl ester	CAS:6597-78-0	Silvex methyl ester	CAS:4841-20-7
Chloramben methyl ester	CAS:7286-84-2	Dichlorprop methyl ester	CAS:57153-17-0	2,4,5-T methyl ester	CAS:1928-37-6
2,4-D methyl ester	CAS:1928-38-7	Dinoseb methyl ether	CAS:6099-79-2	3,5-Dichlorobenzoic acid methyl ester	CAS:2905-67-1
Dalapon methyl ester	CAS:17640-02-7	4-Nitroanisole	CAS:100-17-4		
2,4-DB methyl ester	CAS:18625-12-2	Pentachloroanisole	CAS:1825-21-4		
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113751	Price : BCA

Methyl Derivatives Standard Solution - 16 components**(EPA 515)**

Acifluorfen methyl ester	id : 15801	100 µg/ml	3,5-Dichlorobenzoic acid methyl ester	CAS:2905-67-1	100 µg/ml
Bentazon methyl derivative	CAS:2225-40-3	200 µg/ml	Dichlorprop methyl ester	CAS:57153-17-0	300 µg/ml
Chloramben methyl ester	CAS:7286-84-2	100 µg/ml	Dinoseb methyl ether	CAS:6099-79-2	200 µg/ml
2,4-D methyl ester	CAS:1928-38-7	200 µg/ml	4-Nitroanisole	CAS:100-17-4	100 µg/ml
Dalapon methyl ester	CAS:17640-02-7	1 300 µg/ml	Pentachloroanisole	CAS:1825-21-4	100 µg/ml
2,4-DB methyl ester	CAS:18625-12-2	800 µg/ml	Picloram methyl ester	id : 15841	100 µg/ml
DCPA	CAS:1861-32-1	100 µg/ml	Silvex methyl ester	CAS:4841-20-7	100 µg/ml
Dicamba-methyl ester	CAS:6597-78-0	100 µg/ml	2,4,5-T methyl ester	CAS:1928-37-6	100 µg/ml
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113721	Price : HA

NCC Standard Solution - 7 components (EPA 633)

Bromacil	CAS:314-40-9	Metribuzin	CAS:21087-64-9	Tricyclazole	CAS:41814-78-2
Deet	CAS:134-62-3	Terbacil	CAS:5902-51-2		
Hexazinone	CAS:51235-04-2	Triadimefon	CAS:43121-43-3		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F115931	Price : DF

NCC Standard Solution - 8 components**(EPA 505)**

Atrazine	CAS:1912-24-9	Heptachlor	CAS:76-44-8	trans-Nonachlor	CAS:39765-80-5
cis-Chlordane	CAS:5103-71-9	Hexachlorocyclopentadiene	CAS:77-47-4	Simazine (CAT)	CAS:122-34-9
trans-Chlordane	CAS:5103-74-2	cis-Nonachlor	CAS:5103-73-1		
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F113171	Price : FF

NCC Standard Solution - 8 components

Atrazine	CAS:1912-24-9	Cyanazine	CAS:21725-46-2	Simazine (CAT)	CAS:122-34-9
Atrazine-desethyl	CAS:6190-65-4	Propazine	CAS:139-40-2	Terbuthylazine	CAS:5915-41-3
Atrazine-desisopropyl	CAS:1007-28-9	Sebuthylazine	CAS:7286-69-3		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : 128781	Price : DF

NCC Standard Solution - 11 components**(EPA 619)**

Ametryn	CAS:834-12-8	Prometryn	CAS:7287-19-6	Simazine (CAT)	CAS:122-34-9
Atraton	CAS:1610-17-9	Propazine	CAS:139-40-2	Terbuthylazine	CAS:5915-41-3
Atrazine	CAS:1912-24-9	Sebumeton	CAS:26259-45-0	Terbutryn	CAS:886-50-0
Prometon	CAS:1610-18-0	Simetryn	CAS:1014-70-6		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F115751	Price : EF

NCC Standard Solution - 16 components**(EPA 505)**

Alachlor	CAS:15972-60-8	10 µg/ml	Heptachlor	CAS:76-44-8	1 µg/ml
Aldrin	CAS:309-00-2	1 µg/ml	Hexachlorobenzene	CAS:118-74-1	1 µg/ml
Atrazine	CAS:1912-24-9	250 µg/ml	Hexachlorocyclopentadiene	CAS:77-47-4	1 µg/ml
alpha-Chlordane	CAS:5103-71-9	1 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	1 µg/ml
Gamma-Chlordane	CAS:5103-74-2	1 µg/ml	cis-Nonachlor	CAS:5103-73-1	1 µg/ml
Dieldrin	CAS:60-57-1	1 µg/ml	trans-Nonachlor	CAS:39765-80-5	1 µg/ml
Endrin	CAS:72-20-8	1 µg/ml	Simazine (CAT)	CAS:122-34-9	250 µg/ml
Gamma-HCH (Lindane)	CAS:58-89-9	1 µg/ml	cis-Heptachlorepoxide -, exo-isomer B)	CAS:1024-57-3	1 µg/ml
Solvent : Acetone		ampoule 1 ml		Ref : F113181	Price : GF

NCC Standard Solution - 17 components

Atrazine	CAS:1912-24-9	Hexazinone	CAS:51235-04-2	Monolinuron	CAS:1746-81-2
Atrazine-desethyl	CAS:6190-65-4	Isoproturon	CAS:34123-59-6	Propazine	CAS:139-40-2
Atrazine-desisopropyl	CAS:1007-28-9	Linuron	CAS:330-55-2	Sebuthylazine	CAS:7286-69-3
Chlortoluron	CAS:15545-48-9	Methabenzthiazuron	CAS:18691-97-9	Simazine (CAT)	CAS:122-34-9
Cyanazine	CAS:21725-46-2	Metobromuron	CAS:3060-89-7	Terbuthylazine	CAS:5915-41-3
Diuron	CAS:330-54-1	Metoxuron	CAS:19937-59-8		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F127651	Price : FC

NCC Standard Solution - 17 components**(ISO 11369:1997)**

Atrazine	CAS:1912-24-9	Isoproturon	CAS:34123-59-6	Metoxuron	CAS:19937-59-8
Atrazine-desethyl	CAS:6190-65-4	Linuron	CAS:330-55-2	Monolinuron	CAS:1746-81-2
Chlortoluron	CAS:15545-48-9	Metazachlor	CAS:67129-08-2	Sebuthylazine	CAS:7286-69-3
Cyanazine	CAS:21725-46-2	Methabenzthiazuron	CAS:18691-97-9	Simazine (CAT)	CAS:122-34-9
Diuron	CAS:330-54-1	Metobromuron	CAS:3060-89-7	Terbuthylazine	CAS:5915-41-3
Hexazinone	CAS:51235-04-2	Metolachlor	CAS:51218-45-2		
Solvent : Acetonitrile		ampoule 10 ml	10 µg/ml	Ref : F129021	Price : EF

NCC Standard Solution - 17 components**(EPA 8270C)**

Anilazine	CAS:101-05-3	3,3'-Dimethoxybenzidine	CAS:119-90-4	6-Propyl-2-thiouracil	CAS:51-52-5
o-Anisidine	CAS:90-04-0	3,3'-Dimethylbenzidine	CAS:119-93-7	Safrole	CAS:94-59-7
Barban	CAS:101-27-9	Fluchloralin	CAS:33245-39-5	Strychnine	CAS:57-24-9
Benazide	CAS:92-87-5	Isosafrole	CAS:120-58-1	Toluene-2,4-diisocyanate	CAS:584-84-9
Bromoxynil	CAS:1689-84-5	Mestranol	CAS:72-33-3	Trifluralin	CAS:1582-09-8
3,3'-Dichlorobenzidine	CAS:91-94-1	Mexacarbate	CAS:315-18-4		
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119721	Price : FF

NCC / OPP Laboratory performance check mixture - 6 components**(EPA 507)**

Atrazine	CAS:1912-24-9	15 µg/ml	Prometon	CAS:1610-18-0	30 µg/m
Bromacil	CAS:314-40-9	500 µg/m	Triphenyl phosphate	CAS:115-86-6	250 µg/m
1,3-Dimethyl-2-nitrobenzene	CAS:81-20-9	250 µg/m	Vernolate	CAS:1929-77-7	5 µg/m
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113431	Price : GF

NCC / OPP Standard Solution - 2 components**(EPA 507)**

Paraoxon-methyl	CAS:950-35-6	Simazine (CAT)	CAS:122-34-9		
Solvent : Acetone		ampoule 1 ml	1 000 µg/ml	Ref : F113401	Price : CF

NCC / OPP Standard Solution - 6 components**(EPA 507)**

Ametryn	CAS:834-12-8	Disulfoton	CAS:298-04-4	Merphos	CAS:150-50-5
Cycloate	CAS:1134-23-2	Fenamiphos	CAS:22224-92-6	Prometon	CAS:1610-18-0
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113341	Price : KC

NCC / OPP Standard Solution - 8 components**(EPA 507)**

Dichlorvos	CAS:62-73-7	Napropamide	CAS:15299-99-7	Tebuthiuron	CAS:34014-18-1
Fenarimol	CAS:60168-88-9	Pebulate	CAS:1114-71-2	Terbacil	CAS:5902-51-2
Fluridone	CAS:59756-60-4	Simetryn	CAS:1014-70-6		
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113381	Price : EF

NCC / OPP Standard Solution - 8 components**(EPA 507)**

Benfluralin	CAS:1861-40-1	Oxyfluorfen	CAS:42874-03-3	Propachlor	CAS:1918-16-7
Isopropalin	CAS:33820-53-0	Pendimethalin	CAS:40487-42-1	Trifluralin	CAS:1582-09-8
Oxadiazon	CAS:19666-30-9	Profluralin	CAS:26399-36-0		
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113391	Price : EF

NCC / OPP Standard Solution - 9 components**(EPA 507)**

Butachlor	CAS:23184-66-9	Metolachlor	CAS:51218-45-2	Norflurazon	CAS:27314-13-2
Carboxin	CAS:5234-68-4	Metribuzin	CAS:21087-64-9	Terbufos	CAS:13071-79-9
Diazinon	CAS:333-41-5	MGK 264	CAS:113-48-4	Vernolate	CAS:1929-77-7
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113361	Price : FA

NCC / OPP Standard Solution - 9 components**(EPA 507)**

Atrazine	CAS:1912-24-9	Ethoprophos	CAS:13194-48-4	Propazine	CAS:139-40-2
Diphenamid	CAS:957-51-7	Mevinphos	CAS:7786-34-7	Terbutryn	CAS:886-50-0
EPTC	CAS:759-94-4	Prometryn	CAS:7287-19-6	Triadimefon	CAS:43121-43-3
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113351	Price : EF

NCC / OPP Standard Solution - 10 components**(EPA 507)**

Alachlor	CAS:15972-60-8	Chlorpropham	CAS:101-21-3	trans-Tetrachlorvinphos	CAS:22248-79-9
Atraton	CAS:1610-17-9	Hexazinone	CAS:51235-04-2	Tricyclazole	CAS:41814-78-2
Bromacil	CAS:314-40-9	Molinate	CAS:2212-67-1		
Butylate	CAS:2008-41-5	Pronamide	CAS:23950-58-5		
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113371	Price : FA

NCC / OPP Standard Solution - 12 components**(ISO 10695:2000)**

Atrazine	CAS:1912-24-9	Parathion-methyl	CAS:298-00-0	Simazine (CAT)	CAS:122-34-9
Cyanazine	CAS:21725-46-2	Pendimethalin	CAS:40487-42-1	Terbutylazine	CAS:5915-41-3
Metazachlor	CAS:67129-08-2	Propazine	CAS:139-40-2	Trifluralin	CAS:1582-09-8
Parathion-ethyl	CAS:56-38-2	Sebutylazine	CAS:7286-69-3	Vinclozolin	CAS:50471-44-8
Solvent : Acetone		ampoule 10 ml	10 µg/ml	Ref : F128971	Price : EF

NCC / OPP Standard Solution - 21 components

Acetochlor	CAS:34256-82-1	Diazinon	CAS:333-41-5	Metolachlor	CAS:51218-45-2
Alachlor	CAS:15972-60-8	Dichlobenil	CAS:1194-65-6	Prometryn	CAS:7287-19-6
Atrazine	CAS:1912-24-9	Dimethoate	CAS:60-51-5	Propachlor	CAS:1918-16-7
Atrazine-desethyl	CAS:6190-65-4	Hexazinone	CAS:51235-04-2	Propazine	CAS:139-40-2
Chlorfenvinphos	CAS:470-90-6	Lenacil	CAS:2164-08-1	Simazine (CAT)	CAS:122-34-9
Cyanazine	CAS:21725-46-2	Metamitron	CAS:41394-05-2	Terbutylazine	CAS:5915-41-3
Desmetryn	CAS:1014-69-3	Metazachlor	CAS:67129-08-2	Terbutryn	CAS:886-50-0
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F127591	Price : GF

OCs Decomposition Solution - 2 components

p,p'-DDT	CAS:50-29-3	2 µg/ml	Endrin	CAS:72-20-8	1 µg/ml
Solvent : Acetone		ampoule 1 ml		Ref : F116821	Price : CF

OCs Laboratory Performance Check Solution - 4 components

Chlorothalonil	CAS:1897-45-6	50 µg/ml	Delta-HCH	CAS:319-86-8	40 µg/ml
Chlorpyrifos	CAS:2921-88-2	2 µg/ml	Chlorthal-dimethyl	CAS:1861-32-1	50 µg/ml
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113531	Price : DF

OCs Laboratory Performance Check Solution - 5 components

DBOB	CAS:10386-84-2	250 µg/ml	4-Nitrophenol	CAS:100-02-7	1 600 µg/ml
3,5-Dichlorobenzoic acid	CAS:51-36-5	600 µg/ml	Dinoseb	CAS:88-85-7	4 µg/ml
2,4-Dichlorophenylacetic acid	CAS:19719-28-9	500 µg/ml			
Solvent : Iso-octane		ampoule 1 ml		Ref : F119181	Price : DF

OCs / NCC Standard Solution - 12 components

Alachlor	CAS:15972-60-8	50 µg/ml	Gamma-HCH (Lindane)	CAS:58-89-9	20 µg/ml
Aldrin	CAS:309-00-2	20 µg/ml	Heptachlor	CAS:76-44-8	20 µg/ml
Atrazine	CAS:1912-24-9	500 µg/ml	Hexachlorobenzene	CAS:118-74-1	10 µg/ml
cis-Heptachlorepoxyde (cis-, exo-, isomer B)	CAS:1024-57-3	20 µg/ml	Hexachlorocyclopentadiene	CAS:77-47-4	20 µg/ml
Dieldrin	CAS:60-57-1	20 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	200 µg/ml
Endrin	CAS:72-20-8	20 µg/ml	Simazine (CAT)	CAS:122-34-9	100 µg/ml
Solvent : Acetone		ampoule 1 ml		Ref : F113161	Price : EF

OCs / NCS Standard Solution - 19 components

Alachlor	CAS:15972-60-8	Chlorothalonil	CAS:1897-45-6	Metribuzin	CAS:21087-64-9
Atrazine	CAS:1912-24-9	Cyanazine	CAS:21725-46-2	Permethrin	CAS:52645-53-1
Butachlor	CAS:23184-66-9	DCPA	CAS:1861-32-1	Propachlor	CAS:1918-16-7
alpha-Chlordane	CAS:5103-71-9	Etridiazole	CAS:2593-15-9	Simazine (CAT)	CAS:122-34-9
Gamma-Chlordane	CAS:5103-74-2	Hexachlorobenzene	CAS:118-74-1	Trifluralin	CAS:1582-09-8
Chlorobenzilate	CAS:510-15-6	Hexachlorocyclopentadiene	CAS:77-47-4		
Chloroneb	CAS:2675-77-6	Metolachlor	CAS:51218-45-2		
Solvent : Ethyl Acetate		ampoule 1 ml	100 µg/ml	Ref : F113551	Price : GF
			(Permethrin: 200 µg/ml)		

OCs / Phenoxyacetic Herbicides Standard Solution - 6 components

Acifluorfen	CAS:50594-66-6	Chloramben	CAS:133-90-4	3,5-Dichlorobenzoic acid	CAS:51-36-5
Bentazone	CAS:25057-89-0	DCPA diacid	CAS:2136-79-0	Picloram	CAS:1918-02-1
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F119120	Price : EF

OCs / Phenoxyacetic Herbicides Standard Solution - 12 components

2,4-D	CAS:94-75-7	Dichlorprop	CAS:120-36-5	4-Nitrophenol	CAS:100-02-7
Dalapon	CAS:75-99-0	Dinoseb	CAS:88-85-7	Pentachlorophenol	CAS:87-86-5
2,4-DB	CAS:94-82-6	MCPA	CAS:94-74-6	Silvex	CAS:93-72-1
Dicamba	CAS:1918-00-9	Mecoprop (MCPP)	CAS:7085-19-0	2,4,5-T	CAS:93-76-5
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F119110	Price : FA

OCs QC Check Sample Concentrate - 16 components

Aldrin	CAS:309-00-2	20 µg/ml	Dieldrin	CAS:60-57-1	20 µg/ml
Alpha-HCH	CAS:319-84-6	20 µg/ml	Endosulfan I	CAS:959-98-8	20 µg/ml
Beta-HCH	CAS:319-85-7	20 µg/ml	Endosulfan II	CAS:33213-65-9	100 µg/ml
Gamma-HCH (Lindane)	CAS:58-89-9	20 µg/ml	Endosulfan-total (sulfate)	CAS:1031-07-8	100 µg/ml
Delta-HCH	CAS:319-86-8	20 µg/ml	Endrin	CAS:72-20-8	100 µg/ml
4,4'-DDD (TDE)	CAS:72-54-8	100 µg/ml	Endrin aldehyde	CAS:7421-93-4	20 µg/ml
4,4'-DDE	CAS:72-55-9	20 µg/ml	Heptachlor	CAS:76-44-8	20 µg/ml
4,4'-DDT	CAS:50-29-3	100 µg/ml	Heptachlor epoxide	CAS:1024-57-3	20 µg/ml
Solvent : Iso-Octane		ampoule 1 ml		Ref : F115171	Price : FF

OCs QC Standard Solution - 17 components**(EPA 8080A)**

Aldrin	CAS:309-00-2	0.02 µg/ml	Endosulfan I	CAS:959-98-8	0.02 µg/ml
Alpha-HCH	CAS:319-84-6	0.02 µg/ml	Endosulfan II	CAS:33213-65-9	0.10 µg/ml
Beta-HCH	CAS:319-85-7	0.02 µg/ml	Endosulfan-total (sulfate)	CAS:1031-07-8	0.10 µg/ml
Gamma-HCH (Lindane)	CAS:58-89-9	0.02 µg/ml	Endrin	CAS:72-20-8	0.10 µg/ml
Delta-HCH	CAS:319-86-8	0.02 µg/ml	Endrin aldehyde	CAS:7421-93-4	0.02 µg/ml
4,4'-DDD (TDE)	CAS:72-54-8	0.10 µg/ml	Heptachlor	CAS:76-44-8	0.02 µg/ml
4,4'-DDE	CAS:72-55-9	0.02 µg/ml	Heptachlor epoxide	CAS:1024-57-3	0.02 µg/ml
4,4'-DDT	CAS:50-29-3	0.10 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	0.02 µg/ml
Dieldrin	CAS:60-57-1	0.02 µg/ml			

Solvent : Acetone

ampoule 1 ml

Ref : F116711

Price : GF

OCs Standard Solution - 3 components**(EPA 8080A)**

2,4'-DDD	CAS:53-19-0	2,4'-DDE	CAS:3424-82-6	2,4'-DDT	CAS:789-02-6
----------	-------------	----------	---------------	----------	--------------

Solvent : Iso-Octane

ampoule 1 ml

250 µg/ml

Ref : F116811

Price : DF

OCs Standard Solution - 5 components**(EPA 608)**

Chlorothalonil	CAS:1897-45-6	1 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	40 µg/ml
DCPA	CAS:1861-32-1	3 µg/ml	Permethrin	CAS:52645-53-1	400 µg/ml
Dicloran (CNA)	CAS:99-30-9	2 µg/ml			

Solvent : n-Hexane

ampoule 1 ml

Ref : F115211

Price : DF

OCs Standard Solution - 5 components**(EPA 608)**

Chlorothalonil	CAS:1897-45-6	100 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	100 µg/ml
DCPA	CAS:1861-32-1	100 µg/ml	Permethrin	CAS:52645-53-1	200 µg/ml
Dicloran (CNA)	CAS:99-30-9	100 µg/ml			

Solvent : Iso-Octane

ampoule 1 ml

Ref : F115191

Price : DF

OCs Standard Solution - 6 components**(EPA 8081A)**

Chlorobenzilate	CAS:510-15-6	1,2-Dibromo-3-chloropropane	CAS:96-12-8	Hexachlorocyclopentadiene	CAS:77-47-4
Diallate	CAS:2303-16-4	Hexachlorobenzene	CAS:118-74-1	Isodrin	CAS:465-73-6

Solvent : n-Hexane/Toluene (1/1)

ampoule 1 ml

1 000 µg/ml

Ref : F117231

Price : EF

OCs Standard Solution - 7 components**(EPA 508)**

Alachlor	CAS:15972-60-8	Hexachlorocyclopentadiene	CAS:77-47-4	Simazine (CAT)	CAS:122-34-9
Atrazine	CAS:1912-24-9	Metolachlor	CAS:51218-45-2		
Cyanazine	CAS:21725-46-2	Metribuzin	CAS:21087-64-9		

Solvent : Ethyl Acetate

ampoule 1 ml

500 µg/ml

Ref : F113461

Price : EA

OCs Standard Solution - 7 components**(EPA 608)**

Chlorobenzilate	CAS:510-15-6	1,2-Dibromo-3-chloropropane	CAS:96-12-8	Propachlor	CAS:1918-16-7
Chloroneb	CAS:2675-77-6	Etridiazole	CAS:2593-15-9		
Chloropropylate	CAS:5836-10-2	Pentachloronitrobenzene	CAS:82-68-8		

Solvent : Iso-Octane

ampoule 1 ml

1 000 µg/ml

Ref : F115181

Price : FA

OCs Standard Solution - 7 components**(EPA 608)**

Chlorobenzilate	CAS:510-15-6	200 µg/ml	Etridiazole	CAS:2593-15-9	40 µg/ml
Chloroneb	CAS:2675-77-6	40 µg/ml	Pentachloronitrobenzene	CAS:82-68-8	60 µg/ml
Chloropropylate	CAS:5836-10-2	200 µg/ml	Propachlor	CAS:1918-16-7	1 000 µg/ml
1,2-Dibromo-3-chloropropane	CAS:96-12-8	40 µg/ml			

Solvent : Iso-Octane

ampoule 1 ml

Ref : F115201

Price : EA

OCs Standard Solution - 9 components**(EPA 617)**

Captan	CAS:133-06-2	Dicofol	CAS:115-32-2	PCNB	CAS:82-68-8
Carbophenothion	CAS:786-19-6	Isodrin	CAS:465-73-6	Perthane	CAS:72-56-0
Dicloran (CNA)	CAS:99-30-9	Mirex	CAS:2385-85-5	Trifluralin	CAS:1582-09-8

Solvent : Iso-Octane

ampoule 1 ml

1 000 µg/ml

Ref : F115651

Price : EF

OCs Standard Solution - 9 components**(EPA 508)**

Chlorobenzilate	CAS:510-15-6	Chlorthal-dimethyl	CAS:1861-32-1	Permethrin	CAS:52645-53-1
Chloroneb	CAS:2675-77-6	Etridiazole	CAS:2593-15-9	Propachlor	CAS:1918-16-7
Chlorothalonil	CAS:1897-45-6	Hexachlorobenzene	CAS:118-74-1	Trifluralin	CAS:1582-09-8

Solvent : Methyl-tert butylether

ampoule 1 ml

1 000 µg/ml

Ref : F113481

Price : FA

OCs Standard Solution - 11 components

oxy-Chlordane	CAS:27304-13-8	Dieldrin	CAS:60-57-1	Hexachlorobenzene	CAS:118-74-1
4,4'-DDE	CAS:72-55-9	Beta-HCH	CAS:319-85-7	Mirex	CAS:2385-85-5
2,4'-DDT	CAS:789-02-6	Gamma-HCH (Lindane)	CAS:58-89-9	trans-Nonachlor	CAS:39765-80-5
4,4'-DDT	CAS:50-29-3	Heptachlor-exo-epoxide	CAS:1024-57-3		

Solvent : n-Nonane

ampoule 1 ml

10 µg/ml

Ref : F128701

Price : FC

OCs Standard Solution - 12 components**(EPA 508)**

alpha-Chlordane	CAS:5103-71-9	Chlorothalonil	CAS:1897-45-6	Hexachlorobenzene	CAS:118-74-1
Gamma-Chlordane	CAS:5103-74-2	Chlorpyrifos	CAS:2921-88-2	Propachlor	CAS:1918-16-7
Chlorobenzilate	CAS:510-15-6	DCPA	CAS:1861-32-1	Trifluralin	CAS:1582-09-8
Chloroneb	CAS:2675-77-6	Etridiazole	CAS:2593-15-9	Permethrin	CAS:52645-53-1
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113451	Price : GH
<i>(Etridiazole: 2 000 µg/ml)</i>					

OCs Standard Solution - 12 components**(EPA 525.1)**

Alachlor	CAS:15972-60-8	Gamma-Chlordane	CAS:5103-74-2	Heptachlor-exo-epoxide	CAS:1024-57-3
Aldrin	CAS:309-00-2	Gamma-HCH (Lindane)	CAS:58-89-9	Methoxychlor (DMTD)	CAS:72-43-5
Atrazine	CAS:1912-24-9	Endrin	CAS:72-20-8	trans-Nonachlor	CAS:39765-80-5
alpha-Chlordane	CAS:5103-71-9	Heptachlor	CAS:76-44-8	Simazine (CAT)	CAS:122-34-9
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F113991	Price : GA
		ampoule 1 ml	500 µg/ml	Ref : F114001	Price : BCA

OCs Standard Solution - 16 components**(EPA 608)**

Aldrin	CAS:309-00-2	Endosulfan II	CAS:33213-65-9	Gamma-HCH (Lindane)	CAS:58-89-9
4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	Delta-HCH	CAS:319-86-8
4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8
4,4'-DDT	CAS:50-29-3	Endrin aldehyde	CAS:7421-93-4	Heptachlor epoxide	CAS:1024-57-3
Dieldrin	CAS:60-57-1	Alpha-HCH	CAS:319-84-6		
Endosulfan I	CAS:959-98-8	Beta-HCH	CAS:319-85-7		
Solvent : Iso-Octane		ampoule 1 ml	20 µg/ml	Ref : F115051	Price : FC

OCs Standard Solution - 16 components**(EPA 508)**

4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	Delta-HCH	CAS:319-86-8
4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8
4,4'-DDT	CAS:50-29-3	Endrin aldehyde	CAS:7421-93-4	Heptachlor-exo-epoxide	CAS:1024-57-3
Dieldrin	CAS:60-57-1	Alpha-HCH	CAS:319-84-6	Methoxychlor (DMTD)	CAS:72-43-5
Endosulfan-alpha	CAS:959-98-8	Beta-HCH	CAS:319-85-7		
Endosulfan-beta	CAS:33213-65-9	Gamma-HCH (Lindane)	CAS:58-89-9		
Solvent : Methyl-tert butylether		ampoule 1 ml	1 000 µg/ml	Ref : F113441	Price : JD

OCs Standard Solution - 17 components**(EPA 8080A)**

Aldrin	CAS:309-00-2	Endosulfan II	CAS:33213-65-9	Gamma-HCH (Lindane)	CAS:58-89-9
4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	Delta-HCH	CAS:319-86-8
4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8
4,4'-DDT	CAS:50-29-3	Endrin aldehyde	CAS:7421-93-4	Heptachlor epoxide	CAS:1024-57-3
Dieldrin	CAS:60-57-1	Alpha-HCH	CAS:319-84-6	Methoxychlor (DMTD)	CAS:72-43-5
Endosulfan I	CAS:959-98-8	Beta-HCH	CAS:319-85-7		
Solvent : Acetone		ampoule 1 ml	2 000 µg/ml	Ref : F116701	Price : BEA

OCs Standard Solution - 17 components**(EPA 508.1)**

Aldrin	CAS:309-00-2	Endosulfan-beta	CAS:33213-65-9	Gamma-HCH (Lindane)	CAS:58-89-9
4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	Delta-HCH	CAS:319-86-8
4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8
4,4'-DDT	CAS:50-29-3	Endrin aldehyde	CAS:7421-93-4	Heptachlor-exo-epoxide	CAS:1024-57-3
Dieldrin	CAS:60-57-1	Alpha-HCH	CAS:319-84-6	Methoxychlor (DMTD)	CAS:72-43-5
Endosulfan-alpha	CAS:959-98-8	Beta-HCH	CAS:319-85-7		
Solvent : Ethyl Acetate		ampoule 1 ml	100 µg/ml	Ref : F113541	Price : GF

OCs Standard Solution - 18 components**(EPA 617)**

Aldrin	CAS:309-00-2	4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8
alpha-BHC	CAS:319-84-6	4,4'-DDT	CAS:50-29-3	Endrin ketone	CAS:53494-70-5
beta-BHC	CAS:319-85-7	Dieldrin	CAS:60-57-1	Endrin aldehyde	CAS:7421-93-4
gamma-BHC	CAS:58-89-9	Endosulfan I	CAS:959-98-8	Heptachlor	CAS:76-44-8
delta-BHC	CAS:319-86-8	Endosulfan II	CAS:33213-65-9	Methoxychlor (DMTD)	CAS:72-43-5
4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	cis-Heptachlorepoxyde (cis-, exo-, isomer B)	CAS:1024-57-3
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F115641	Price : BCF

OCs Standard Solution - 20 components**(EPA 8081A)**

Aldrin	CAS:309-00-2	4,4'-DDD (TDE)	CAS:72-54-8	Endrin	CAS:72-20-8
alpha-BHC	CAS:319-84-6	4,4'-DDE	CAS:72-55-9	Endrin aldehyde	CAS:7421-93-4
beta-BHC	CAS:319-85-7	4,4'-DDT	CAS:50-29-3	Endrin ketone	CAS:53494-70-5
gamma-BHC	CAS:58-89-9	Dieldrin	CAS:60-57-1	Heptachlor	CAS:76-44-8
delta-BHC	CAS:319-86-8	Endosulfan-alpha	CAS:959-98-8	Methoxychlor (DMTD)	CAS:72-43-5
alpha-Chlordane	CAS:5103-71-9	Endosulfan-beta	CAS:33213-65-9	cis-Heptachlorepoxyde (cis-, exo-, isomer B)	CAS:1024-57-3
Gamma-Chlordane	CAS:5103-74-2	Endosulfan-total (sulfate)	CAS:1031-07-8		
Solvent : n-Hexane/Toluene (1/1)		ampoule 1 ml	1 000 µg/ml	Ref : F117221	Price : BFA

OCs Pesticide Standard Solution - 20 components**(ISO 6468:1996)**

Aldrin	CAS:309-00-2	Dieldrin	CAS:60-57-1	Delta-HCH	CAS:319-86-8
2,4'-DDE	CAS:3424-82-6	Endosulfan-alpha	CAS:959-98-8	epsilon-HCH	CAS:6108-10-7
4,4'-DDE	CAS:72-55-9	Endosulfan-beta	CAS:33213-65-9	Heptachlor	CAS:76-44-8
2,4'-DDD (o,p'-TDE)	CAS:53-19-0	Endrin	CAS:72-20-8	Heptachlor-endo-epoxide	CAS:28044-83-9
4,4'-DDD (TDE)	CAS:72-54-8	Alpha-HCH	CAS:319-84-6	Heptachlor-exo-epoxide	CAS:1024-57-3
2,4'-DDT	CAS:789-02-6	Beta-HCH	CAS:319-85-7	Methoxychlor (DMTD)	CAS:72-43-5
4,4'-DDT	CAS:50-29-3	Gamma-HCH (Lindane)	CAS:58-89-9		
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F147410	Price : EG

OCs Standard Solution - 20 components

Aldrin	CAS:309-00-2	4,4'-DDT	CAS:50-29-3	Endrin ketone	CAS:53494-70-5
alpha-BHC	CAS:319-84-6	Dieldrin	CAS:60-57-1	Heptachlor	CAS:76-44-8
beta-BHC	CAS:319-85-7	Endosulfan I	CAS:959-98-8	Heptachlor epoxide	CAS:1024-57-3
delta-BHC	CAS:319-86-8	Endosulfan II	CAS:33213-65-9	Isodrin	CAS:465-73-6
gamma-BHC	CAS:58-89-9	Endosulfan-total (sulfate)	CAS:1031-07-8	Methoxychlor (DMTD)	CAS:72-43-5
4,4'-DDD (TDE)	CAS:72-54-8	Endrin	CAS:72-20-8	Mirex	CAS:2385-85-5
4,4'-DDE	CAS:72-55-9	Endrin aldehyde	CAS:7421-93-4		
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119351	Price :FF

OCs Standard Solution - 23 components

Aldrin	CAS:309-00-2	2,4'-DDT	CAS:789-02-6	Heptachlor	CAS:76-44-8
cis-Chlordane	CAS:5103-71-9	4,4'-DDT	CAS:50-29-3	Heptachlor-endo-epoxide	CAS:28044-83-9
trans-Chlordane	CAS:5103-74-2	Dieldrin	CAS:60-57-1	Heptachlor-exo-epoxide	CAS:1024-57-3
oxy-Chlordane	CAS:27304-13-8	Endrin	CAS:72-20-8	Hexachlorobenzene	CAS:118-74-1
2,4'-DDD	CAS:53-19-0	Alpha-HCH	CAS:319-84-6	Mirex	CAS:2385-85-5
4,4'-DDD (TDE)	CAS:72-54-8	Beta-HCH	CAS:319-85-7	cis-Nonachlor	CAS:5103-73-1
2,4'-DDE	CAS:3424-82-6	Gamma-HCH (Lindane)	CAS:58-89-9	trans-Nonachlor	CAS:39765-80-5
4,4'-DDE	CAS:72-55-9	Delta-HCH	CAS:319-86-8		
Solvent : n-Nonane		ampoule 1 ml	200 µg/ml	Ref : F128711	Price : KG
		ampoule 1 ml	2 000 µg/ml	Ref : F128731	Price : DAA

OCs Standard Solution - 24 components

Aldrin	CAS:309-00-2	4,4'-DDT	CAS:50-29-3	Beta-HCH	CAS:319-85-7
cis-Chlordane	CAS:5103-71-9	Endosulfan-alpha	CAS:959-98-8	Gamma-HCH (Lindane)	CAS:58-89-9
trans-Chlordane	CAS:5103-74-2	Dieldrin	CAS:60-57-1	Delta-HCH	CAS:319-86-8
2,4'-DDD	CAS:53-19-0	Endosulfan-beta	CAS:33213-65-9	Heptachlor	CAS:76-44-8
4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	Heptachlor-endo-epoxide	CAS:28044-83-9
2,4'-DDE	CAS:3424-82-6	Endrin	CAS:72-20-8	Heptachlor-exo-epoxide	CAS:1024-57-3
4,4'-DDE	CAS:72-55-9	Endrin aldehyde	CAS:7421-93-4	Hexachlorobenzene	CAS:118-74-1
2,4'-DDT	CAS:789-02-6	Alpha-HCH	CAS:319-84-6	Vinclozolin	CAS:50471-44-8
Solvent : Iso-Octane		ampoule 5 ml	20 µg/ml	Ref : F127681	Price GA

OCs Standard Solution - 29 components

Alachlor	CAS:15972-60-8	Chlorothalonil	CAS:1897-45-6	Endrin aldehyde	CAS:7421-93-4
Aldrin	CAS:309-00-2	DCPA	CAS:1861-32-1	Etridiazole	CAS:2593-15-9
Atrazine	CAS:1912-24-9	4,4'-DDD (TDE)	CAS:72-54-8	Gamma-HCH (Lindane)	CAS:58-89-9
alpha-BHC	CAS:319-84-6	4,4'-DDE	CAS:72-55-9	Heptachlor	CAS:76-44-8
beta-BHC	CAS:319-85-7	4,4'-DDT	CAS:50-29-3	Heptachlor-exo-epoxide	CAS:1024-57-3
delta-BHC	CAS:319-86-8	Dieldrin	CAS:60-57-1	Methoxychlor (DMTD)	CAS:72-43-5
cis-Chlordane	CAS:5103-71-9	Endosulfan-alpha	CAS:959-98-8	trans-Nonachlor	CAS:39765-80-5
Gamma-Chlordane	CAS:5103-74-2	Endosulfan-beta	CAS:33213-65-9	Permethrin	CAS:52645-53-1
Chlorobenzilate	CAS:510-15-6	Endosulfan-total (sulfate)	CAS:1031-07-8	Simazine (CAT)	CAS:122-34-9
Chloroneb	CAS:2675-77-6	Endrin	CAS:72-20-8		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F128751	Price : KC

OPP Standard Solution - 3 components**(EPA 622)**

Dichlorvos	CAS:62-73-7	10 µg/ml	cis-Tetrachlorvinphos	CAS:961-11-5	500 µg/ml
Phosdrin	CAS:7786-34-7	30 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F115781	Price: CH

OPP Standard Solution - 4 components**(EPA 614.1)**

Dioxathion	CAS:78-34-2	Ethion	CAS:563-12-2	Terbufos	CAS:13071-79-9
EPN	CAS:2104-64-5				
Solvent : n-Hexane		ampoule 1 ml	200 µg/ml	Ref : F115601	Price : DF

OPP Standard Solution - 4 components**(EPA 614.1)**

Dioxathion	CAS:78-34-2	10 µg/ml	Ethion	CAS:563-12-2	100 µg/ml
EPN	CAS:2104-64-5	200 µg/ml	Terbufos	CAS:13071-79-9	4 µg/ml
Solvent : n-Hexane		ampoule 1 ml		Ref : F115611	Price : DF

OPP Standard Solution - 7 components**(EPA 622)**

Ethoprop	CAS:13194-48-4	25 µg/ml	Chlorpyrifos-ethyl	CAS:2921-88-2	30 µg/ml
Parathion-methyl	CAS:298-00-0	30 µg/ml	Merphos	CAS:150-50-5	25 µg/ml
Ronnel	CAS:299-84-3	30 µg/ml	Diazinon	CAS:333-41-5	60 µg/ml
Chlorpyrifos methyl	CAS:5598-13-0	30 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml	Ref : F115791		Price : EF

OPP Standard Solution - 7 components**(EPA 8141A)**

Dimethoate	CAS:60-51-5	Monocrotophos	CAS:6923-22-4	TEPP	CAS:107-49-3
EPN	CAS:2104-64-5	Parathion-ethyl	CAS:56-38-2		
Malathion	CAS:121-75-5	Sulfotep	CAS:3689-24-5		
Solvent : n-Hexane/Acetone (9/1)		ampoule 1 ml	200 µg/ml	Ref : F119031	Price : EF

OPP Standard Solution - 8 components**(EPA 614)**

Azinphos-methyl	CAS:86-50-0	Disulfoton	CAS:298-04-4	Parathion-ethyl	CAS:56-38-2
Demeton (O+S)	CAS:8065-48-3	Ethion	CAS:563-12-2	Parathion-methyl	CAS:298-00-0
Diazinon	CAS:333-41-5	Malathion	CAS:121-75-5		
Solvent : Acetone		ampoule 1 ml	200 µg/ml	Ref : F115591	Price : EA

OPP Standard Solution - 9 components**(EPA 622)**

Azinphos-methyl	CAS:86-50-0	150 µg/ml	Fenthion	CAS:55-38-9	10 µg/ml
Coumaphos	CAS:56-72-4	150 µg/ml	Phorate	CAS:298-02-2	15 µg/ml
Demeton (O+S)	CAS:8065-48-3	25 µg/ml	Tokuthion	CAS:34643-46-4	50 µg/ml
Disulfoton	CAS:298-04-4	20 µg/ml	Trichloronat	CAS:327-98-0	15 µg/ml
Fensulfothion	CAS:115-90-2	150 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml	Ref : F115771		Price : EF

OPP Standard Solution - 9 components**(EPA 8141A)**

Aspon	CAS:3244-90-4	Dichlofenthion	CAS:97-17-6	Fonofos	CAS:944-22-9
Chlorpyrifos methyl	CAS:5598-13-0	Dicrotophos	CAS:141-66-2	Thionazin	CAS:297-97-2
Crotoxyphos	CAS:7700-17-6	Fenitrothion	CAS:122-14-5	Trichlorfon	CAS:52-68-6
Solvent : n-Hexane		ampoule 1 ml	200 µg/ml	Ref : F119051	Price : FA

OPP Standard Solution - 10 components**(EPA 8141A)**

Azinphos-ethyl	CAS:2642-71-9	Famphur	CAS:52-85-7	Terbufos	CAS:13071-79-9
Chlorfenvinphos	CAS:470-90-6	Leptophos	CAS:21609-90-5	Trithion	CAS:786-19-6
Dioxathion	CAS:78-34-2	Phosmet	CAS:732-11-6		
Ethion	CAS:563-12-2	Phosphamidon	CAS:13171-21-6		
Solvent : n-Hexane/Acetone (9/1)		ampoule 1 ml	200 µg/ml	Ref : F119041	Price : FA

OPP Standard Solution - 19 components**(EN 12918)**

Azinphos-ethyl	CAS:2642-71-9	Dichlorvos	CAS:62-73-7	Parathion-methyl	CAS:298-00-0
Azinphos-methyl	CAS:86-50-0	Dimethoate	CAS:60-51-5	Phosalone	CAS:2310-17-0
Bromophos	CAS:2104-96-3	Fenitrothion	CAS:122-14-5	Propetamphos	CAS:31218-83-4
Chlorfenvinphos	CAS:470-90-6	Fenthion	CAS:55-38-9	Triadimefon	CAS:43121-43-3
Chlorpyrifos-ethyl	CAS:2921-88-2	Malathion	CAS:121-75-5	Triazophos	CAS:24017-47-8
Chlorpyrifos methyl	CAS:5598-13-0	Mevinphos	CAS:7786-34-7		
Diazinon	CAS:333-41-5	Parathion-ethyl	CAS:56-38-2		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F128951	Price : GF
Solvent : Methanol		ampoule 10 ml	10 µg/ml	Ref : F131371	Price : GF

OPP Standard Solution - 20 components**(EPA 8141A)**

Azinphos-methyl	CAS:86-50-0	Disulfoton	CAS:298-04-4	Phorate	CAS:298-02-2
Bolstar	CAS:35400-43-2	Ethoprophos	CAS:13194-48-4	Phosdrin	CAS:7786-34-7
Chlorpyrifos	CAS:2921-88-2	Fensulfothion	CAS:115-90-2	Ronnel	CAS:299-84-3
Coumaphos	CAS:56-72-4	Fenthion	CAS:55-38-9	Stirofos	CAS:22248-79-9
Demeton (O+S)	CAS:8065-48-3	Merphos	CAS:150-50-5	Tokuthion	CAS:34643-46-4
Diazinon	CAS:333-41-5	Naled	CAS:300-76-5	Trichloronat	CAS:327-98-0
Dichlorvos	CAS:62-73-7	Parathion-methyl	CAS:298-00-0		
Solvent : n-Hexane/Acetone (9/1)		ampoule 1 ml	200 µg/ml	Ref : F119021	Price : KG

OPP Standard Solution - 21 components**(EPA 622)**

Azinphos-methyl	CAS:86-50-0	Dichlorvos	CAS:62-73-7	Naled	CAS:300-76-5
Bolstar	CAS:35400-43-2	Disulfoton	CAS:298-04-4	Parathion-methyl	CAS:298-00-0
Chlorpyrifos-ethyl	CAS:2921-88-2	Ethoprop	CAS:13194-48-4	Phorate	CAS:298-02-2
Chlorpyrifos methyl	CAS:5598-13-0	Fensulfothion	CAS:115-90-2	Ronnel	CAS:299-84-3
Coumaphos	CAS:56-72-4	Fenthion	CAS:55-38-9	Stirofos	CAS:22248-79-9
Demeton	CAS:8065-48-3	Merphos	CAS:150-50-5	Tokuthion	CAS:34643-46-4
Diazinon	CAS:333-41-5	Mevinphos	CAS:7786-34-7	Trichloronat	CAS:327-98-0
Solvent : Iso-Octane		ampoule 1 ml	100 µg/ml	Ref : F115761	Price : HF

OPP Standard Solution - 34 components**(EPA 8270C)**

Azinphos-methyl	CAS:86-50-0	EPN	CAS:2104-64-5	Phosmet	CAS:732-11-6
Benzenethiol	CAS:108-98-5	Ethion	CAS:563-12-2	Phosphamidon	CAS:13171-21-6
Captafol	CAS:2425-06-1	Fensulfothion	CAS:115-90-2	Sulfallate	CAS:95-06-7
Captan	CAS:133-06-2	Fenthion	CAS:55-38-9	Sulfoxide	CAS:120-62-7
Chlorfenvinphos	CAS:470-90-6	Hexamethylphosphoric triamide (HMPA)	CAS:680-31-9	Terbufos	CAS:13071-79-9
Coumaphos	CAS:56-72-4	Leptophos	CAS:21609-90-5	cis-Tetrachlorvinphos	CAS:961-11-5
Crotoxyphos	CAS:7700-17-6	Malathion	CAS:121-75-5	Tetraethyl dithiopyrophosphate	CAS:3689-24-5
Demeton-O	CAS:298-03-3	Mevinphos	CAS:7786-34-7	Tetraethylpyrophosphate	CAS:107-49-3
Demeton-S	CAS:126-75-0	Monocrotophos	CAS:6923-22-4	Tri-4-cresyl phosphate	CAS:78-32-0
Dichlorvos	CAS:62-73-7	Naled	CAS:300-76-5	Trimethyl phosphate	CAS:512-56-1
Dicrotophos	CAS:141-66-2	Phosalone	CAS:2310-17-0	Tris (2,3-dibromopropyl) phosphate	CAS:126-72-7
Dioxathion	CAS:78-34-2				

Solvent : Dichloromethane

ampoule 1 ml

100 µg/ml

Ref : F119371

Price : GE

Organochlorine Pesticides Mixture - 9 components**(EPA 505)**

Alachlor	CAS:15972-60-8	Endrin	CAS:72-20-8	Heptachlor-exo-epoxide	CAS:1024-57-3
Aldrin	CAS:309-00-2	Gamma-HCH (Lindane)	CAS:58-89-9	Hexachlorobenzene	CAS:118-74-1
Dieldrin	CAS:60-57-1	Heptachlor	CAS:76-44-8	Methoxychlor (DMTD)	CAS:72-43-5

Solvent : Toluene

ampoule 1 ml

100 µg/ml

Ref : F113151

Price : EF

Organochlorine Pesticide Mixture - 16 components**(EPA 625)**

Aldrin	CAS:309-00-2	Endosulfan II	CAS:33213-65-9	delta-HCH	CAS:319-86-8
4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	gamma-HCH	CAS:58-89-9
4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8
4,4'-DDT	CAS:50-29-3	Endrin aldehyde	CAS:7421-93-4	cis-Heptachlorepoxyde	CAS:1024-57-3
Dieldrin	CAS:60-57-1	alpha-HCH	CAS:319-84-6	(cis-,exo-, isomer B)	
Endosulfan I	CAS:959-98-8	beta-HCH	CAS:319-85-7		

Solvent : Hexane/Toluene (1/1)

ampoule 1 ml

2 000 µg/ml

Ref : F115891

Price : BCF

Organochlorine Pesticide Mixture - 17 components

Aldrin	CAS:309-00-2	Endosulfan II	CAS:33213-65-9	Gamma-HCH (Lindane)	CAS:58-89-9
4,4'-DDD (TDE)	CAS:72-54-8	Endosulfan-total (sulfate)	CAS:1031-07-8	Delta-HCH	CAS:319-86-8
4,4'-DDE	CAS:72-55-9	Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8
4,4'-DDT	CAS:50-29-3	Endrin aldehyde	CAS:7421-93-4	Heptachlor epoxide	CAS:1024-57-3
Dieldrin	CAS:60-57-1	Alpha-HCH	CAS:319-84-6	Methoxychlor (DMTD)	CAS:72-43-5
Endosulfan I	CAS:959-98-8	Beta-HCH	CAS:319-85-7		

Solvent : Acetone

ampoule 1 ml

2 000 µg/ml

Ref : F116701

Price : BEA

Organophosphorus Compounds Mixture - 15 components**(EPA 8321A)**

Asulam	CAS:3337-71-1	Fensulfothion	CAS:115-90-2	Parathion-methyl	CAS:298-00-0
Dichlorvos	CAS:62-73-7	Merphos	CAS:150-50-5	Phorate	CAS:298-02-2
Dimethoate	CAS:60-51-5	Monomyl	CAS:16752-77-5	Thiofanox	CAS:39196-18-4
Disulfoton	CAS:298-04-4	Monocrotophos	CAS:6923-22-4	Trichlorfon	CAS:52-68-6
Famphur	CAS:52-85-7	Naled	CAS:300-76-5	Tris(2,3-dibromopropyl)phosphate	CAS:126-72-7

Solvent : Acetonitrile

ampoule 1 ml

100 µg/ml

Ref : F119931

Price : GF

Pesticide Degradation Check Solution - 2 components**(EPA 508.1)**

4,4'-DDT	CAS:50-29-3	Endrin	CAS:72-20-8
----------	-------------	--------	-------------

Solvent : Ethyl Acetate

ampoule 1 ml

1 µg/ml

Ref : F113581

Price : CF

Pesticide Standard Solution - 2 components

Iprodione	CAS:36734-19-7	Methidathion	CAS:950-37-8
-----------	----------------	--------------	--------------

Solvent : Acetonitrile

ampoule 10 ml

10 µg/ml

Ref : F128851

Price : DF

Pesticide Standard Solution - 5 components**(EPA 627)**

Benfluralin	CAS:1861-40-1	Isopropalin	CAS:33820-53-0	Trifluralin	CAS:1582-09-8
Ethalfuralin	CAS:55283-68-6	Profluralin	CAS:26399-36-0		

Solvent : n-Hexane

ampoule 1 ml

1 000 µg/ml

Ref : F115911

Price : EF

Pesticide Standard Solution - 6 components

Chlorpyrifos	CAS:2921-88-2	Ethion	CAS:563-12-2	Parathion	CAS:56-38-2
Diazinon	CAS:333-41-5	Malathion	CAS:121-75-5	Parathion-methyl	CAS:298-00-0

Solvent : Methanol

ampoule 1 ml

100 µg/ml

Ref : F127671

Price : DF

Pesticide Standard Solution - 8 components

Atrazine	CAS:1912-24-9	Heptachlor	CAS:76-44-8	trans-Nonachlor	CAS:39765-80-5
cis-Chlordane	CAS:5103-71-9	Hexachlorocyclopentadiene	CAS:77-47-4	Simazine (CAT)	CAS:122-34-9
trans-Chlordane	CAS:5103-74-2	cis-Nonachlor	CAS:5103-73-1		

Solvent : Methanol

ampoule 1 ml

100 µg/ml

Ref : F113171

Price : FF

Pesticide Standard Solution (DIN 38407-14) – 8 components

2,4-D	CAS:94-75-7	Fenoprop	CAS:93-72-1	Mecoprop (MCP)	CAS:7085-19-0
2,4-DB	CAS:94-82-6	MCPA	CAS:94-74-6	2,4,5-T	CAS:93-76-5
Dichloroprop	CAS:120-36-5	MCPB	CAS:94-81-5		
Solvent : Methanol		ampoule 1 ml	10 µg/ml	Ref : F062355	Price : CF
		ampoule 10 ml	10 µg/ml	Ref : F128771	Price : BAA

Pesticide Standard Solution - 14 components**(EPA 8270C)**

Aramite	CAS:140-57-8	Dinoseb	CAS:88-85-7	Parathion-methyl	CAS:298-00-0
Chlorobenzilate	CAS:510-15-6	Disulfoton	CAS:298-04-4	Phorate	CAS:298-02-2
Diallate	CAS:2303-16-4	Famphur	CAS:52-85-7	Thionazin	CAS:297-97-2
Dimethoate	CAS:60-51-5	Kepone	CAS:143-50-0	O.O.O-Triethylphosphorothioate	CAS:126-68-1
Dinocap	CAS:39300-45-3	Parathion	CAS:56-38-2		
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119361	Price : GF

Pesticide Standard Solution - 16 components

Alachlor	CAS:15972-60-8	10 µg/ml	Heptachlor	CAS:76-44-8	1 µg/ml
Aldrin	CAS:309-00-2	1 µg/ml	Hexachlorobenzene	CAS:118-74-1	1 µg/ml
Atrazine	CAS:1912-24-9	250 µg/ml	Hexachlorocyclopentadiene	CAS:77-47-4	1 µg/ml
alpha-Chlordane	CAS:5103-71-9	1 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	5 µg/ml
Gamma-Chlordane	CAS:5103-74-2	1 µg/ml	cis-Nonachlor	CAS:5103-73-1	1 µg/ml
Dieldrin	CAS:60-57-1	1 µg/ml	trans-Nonachlor	CAS:39765-80-5	1 µg/ml
Endrin	CAS:72-20-8	1 µg/ml	Simazine (CAT)	CAS:122-34-9	250 µg/ml
Gamma-HCH (Lindane)	CAS:58-89-9	1 µg/ml	cis-Heptachlorepoxyde (cis-, exo-, isomer B)	CAS:1024-57-3	1 µg/ml
Solvent : Acetone		ampoule 1 ml		Ref : F113181	Price : FA

Pesticide Standard Solution - 16 components

Azinphos-ethyl	CAS:2642-71-9	Hexaconazole	CAS:79983-71-4	Propyzamide	CAS:23950-58-5
Bitertanol	CAS:55179-31-2	Kresoxim-methyl	CAS:143390-89-0	Resmethrin	CAS:10453-86-8
Bromopropylate	CAS:18181-80-1	Metalaxyl	CAS:57837-19-1	Triadimefon	CAS:43121-43-3
Cypermethrin	CAS:52315-07-8	Myclobutanil	CAS:88671-89-0	Triadimenol	CAS:55219-65-3
Dicofol	CAS:115-32-2	Prochloraz	CAS:67747-09-5		
Formothion	CAS:2540-82-1	Procymidone	CAS:32809-16-8		
Solvent : Acetonitrile		ampoule 1 ml	10 µg/ml	Ref : F128841	Price : EF

Pesticide Standard Solution - 18 components**(EPA 8270C)**

Acetophenone	CAS:98-86-2	Carbofuran	CAS:1563-66-2	Hydroquinone	CAS:123-31-9
2-Acetamidofluorene	CAS:53-96-3	Carbophenothion	CAS:786-19-6	Maleic anhydride	CAS:108-31-6
1-Acetyl-2-thiourea	CAS:591-08-2	Dichlone	CAS:117-80-6	1,4-Naphthoquinone	CAS:130-15-4
2-Aminoanthraquinone	CAS:117-79-3	Diethylstilbestrol	CAS:56-53-1	Nitrofen	CAS:1836-75-5
p-Benzoquinone	CAS:106-51-4	5,5-Diphenylhydantoin	CAS:57-41-0	4-Nitroquinoline-N-oxide	CAS:56-57-5
Carbaryl	CAS:63-25-2	Ethyl carbamate	CAS:51-79-6	Phthalic anhydride	CAS:85-44-9
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119381	Price : FF

Pesticide Standard Solution - 21 components

Azoxystrobin	CAS:131860-33-8	Demeton-S-methyl-sulfone	CAS:17040-19-6	Phosphamidon	CAS:13171-21-6
Captan	CAS:2425-06-1	Dichlorvos	CAS:62-73-7	Pirimiphos-methyl	CAS:29232-93-7
Captan	CAS:133-06-2	Famoxadone	CAS:131807-57-3	Prochloraz	CAS:67747-09-5
Carbaryl	CAS:63-25-2	Folpet	CAS:133-07-3	Quintozene	CAS:82-68-8
Chlorpyrifos	CAS:2921-88-2	Malathion	CAS:121-75-5	Triadimefon	CAS:43121-43-3
Chlorothalonil	CAS:1897-45-6	Malaoxon	CAS:1634-78-2	Triadimenol	CAS:55219-65-3
Deltamethrin	CAS:52918-63-5	Nitrofen	CAS:1836-75-5	Trichlorfon	CAS:52-68-6
Solvent : Acetonitrile		ampoule 1 ml	10 µg/ml	Ref : F128831	Price : FC

Pesticide Standard Solution - 25 components

Atrazine	CAS:1912-24-9	Ethofumesate	CAS:26225-79-6	Metolachlor	CAS:51218-45-2
Atrazine-desethyl	CAS:6190-65-4	Flurochloridone	CAS:61213-25-0	Metribuzin	CAS:21087-64-9
Atrazine-desisopropyl	CAS:1007-28-9	Flurtamone	CAS:96525-23-4	Prometryn	CAS:7287-19-6
Carbamazepin	CAS:298-46-4	Isoproturon	CAS:34123-59-6	Simazine (CAT)	CAS:122-34-9
Carbetamide	CAS:16118-49-3	Linuron	CAS:330-55-2	Terbutylazine	CAS:5915-41-3
Chloridazon	CAS:1698-60-8	Metamitron	CAS:41394-05-2	Terbutylazine-desethyl	CAS:30125-63-4
Chlortoluron	CAS:15545-48-9	Metazachlor	CAS:67129-08-2	Terbutryn	CAS:886-50-0
Dimefuron	CAS:34205-21-5	Methabenzthiazuron	CAS:18691-97-9		
Diuron	CAS:330-54-1	Metobromuron	CAS:3060-89-7		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F128861	Price : GF

Pesticide Standard Solution - 34 components

Alachlor	CAS:15972-60-8	Fenpropimorph	CAS:67564-91-4	Sebuthylazine	CAS:7286-69-3
Atrazine	CAS:1912-24-9	Flusilazole	CAS:85509-19-9	Simazine (CAT)	CAS:122-34-9
Atrazine-desethyl	CAS:6190-65-4	Folpet	CAS:133-07-3	Tebutam	CAS:35256-85-0
Atrazine-desisopropyl	CAS:1007-28-9	Metazachlor	CAS:67129-08-2	Terbumeton	CAS:33693-04-8
Carbofuran	CAS:1563-66-2	Metolachlor	CAS:51218-45-2	Terbumeton-desethyl	CAS:30125-64-5
Cyanazine	CAS:21725-46-2	Oxadiazon	CAS:19666-30-9	Terbuthylazine	CAS:5915-41-3
Deltamethrin	CAS:52918-63-5	Oxadixyl	CAS:77732-09-3	Terbuthylazine-desethyl	CAS:30125-63-4
Dimethenamid	CAS:87674-68-8	Parathion	CAS:56-38-2	Tri-allate	CAS:2303-17-5
Dimethoate	CAS:60-51-5	Parathion-methyl	CAS:298-00-0	Trifluralin	CAS:1582-09-8
Endosulfan-alpha	CAS:959-98-8	Pendimethalin	CAS:40487-42-1	Vinclozolin	CAS:50471-44-8
Endosulfan-beta	CAS:33213-65-9	Prometryn	CAS:7287-19-6		
Ethofumesate	CAS:26225-79-6	Propazine	CAS:139-40-2		
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F127621	Price : JA

Pesticide Standard Solution - 54 components

Azinphos-ethyl	CAS:2642-71-9	Diflufenican	CAS:83164-33-4	Hexaconazole	CAS:79983-71-4
Azinphos-methyl	CAS:86-50-0	Dimethoate	CAS:60-51-5	Malathion	CAS:121-75-5
Bromuconazole	CAS:116255-48-2	Disulfoton	CAS:298-04-4	Metazachlor	CAS:67129-08-2
Bromophos-ethyl	CAS:4824-78-6	Epoxiconazole	CAS:106325-08-0	Metconazole	CAS:125116-23-6
Bromophos-methyl	CAS:2104-96-3	Esfenvalerate	CAS:66230-04-4	Metolachlor	CAS:51218-45-2
Captan	CAS:133-06-2	Ethion	CAS:563-12-2	Oxadiazon	CAS:19666-30-9
Chlorfenvinphos	CAS:470-90-6	Ethoprophos	CAS:13194-48-4	Paclobutrazol	CAS:76738-62-0
Chlorothalonil	CAS:1897-45-6	Etrinfos	CAS:38260-54-7	Parathion	CAS:56-38-2
Chlorpyrifos	CAS:2921-88-2	Fenbuconazole	CAS:114369-43-6	Parathion-methyl	CAS:298-00-0
Chlorpyrifos methyl	CAS:5598-13-0	Fenitrothion	CAS:122-14-5	Permethrin	CAS:52645-53-1
lambda-Cyhalothrin	CAS:91465-08-6	Fenpropidin	CAS:67306-00-7	Procymidone	CAS:32809-16-8
Cypermethrin	CAS:52315-07-8	Fenpropimorph	CAS:67564-91-4	Propetamphos	CAS:31218-83-4
Cyproconazole	CAS:113096-99-4	Fenthion	CAS:55-38-9	Propiconazole	CAS:60207-90-1
Cyprodinil	CAS:121552-61-2	Fluquinconazole	CAS:136426-54-5	Tebuconazole	CAS:107534-96-3
Deltamethrin	CAS:52918-63-5	Flurochloridone	CAS:61213-25-0	Terbufos	CAS:13071-79-9
Diazinon	CAS:333-41-5	Flusilazole	CAS:85509-19-9	Triadimefon	CAS:43121-43-3
Dichlorvos	CAS:62-73-7	Folpet	CAS:133-07-3	Tri-allate	CAS:2303-17-5
Difenoconazole	CAS:119446-68-3	Formothion	CAS:2540-82-1	Trifluralin	CAS:1582-09-8
Solvent : n-Hexane		ampoule 10 ml	1 µg/ml	Ref : F127661	Price : BCF

Phenoxyacetic Herbicides Standard Solution - 3 components

2,4-D methyl ester	CAS:1928-38-7	Silvex methyl ester	CAS:4841-20-7	2,4,5-T methyl ester	CAS:1928-37-6
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F119171	Price : DA

Phenoxyacetic Herbicides Standard Solution - 10 components

Bentazone	CAS:25057-89-0	Dichlorprop	CAS:120-36-5	Mecoprop (MCP)	CAS:7085-19-0
Bromoxynil	CAS:1689-84-5	Fenoprop	CAS:93-72-1	2,4,5-Trichlorophenoxy acetic acid	CAS:93-76-5
2,4-D	CAS:94-75-7	MCPA	CAS:94-74-6		
2,4-DB	CAS:94-82-6	MCPB	CAS:94-81-5		
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F129031	Price : EF

Phenoxyalkanoic Carbonic Esters Standard Solution - 10 components

Bentazone methyl ester	id : 17751	Dichlorprop methyl ester	CAS:57153-17-0	Mecoprop methyl ester	CAS:23844-56-6
Bromoxynil methyl ester	id : 17761	Fenoprop methyl ester	CAS:4841-20-7	2,4,5-Trichlorophenoxy acetic acid	CAS:1928-37-6
2,4-D methyl ester	CAS:1928-38-7	MCPA methyl ester	CAS:2436-73-9		
2,4-DB methyl ester	CAS:18625-12-2	MCPB methyl ester	CAS:57153-18-1		
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F129041	Price : EE

SDWA Herbicides Mixture - 6 components

2,4-D	CAS:94-75-7	1 300 µg/ml	Pentachlorophenol	CAS:87-86-5	100 µg/ml
Dalapon	CAS:75-99-0	1 300 µg/ml	Picloram	CAS:1918-02-1	100 µg/ml
Dinoseb	CAS:88-85-7	200 µg/ml	Silvex	CAS:93-72-1	100 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F113871	Price : EA

SDWA Methylated Herbicides Mixture - 6 components

2,4-D methyl ester	CAS:1928-38-7	200 µg/ml	Pentachloroanisole	CAS:1825-21-4	100 µg/ml
Dalapon methyl ester	CAS:17640-02-7	1 300 µg/ml	Picloram methyl ester	id : 15841	100 µg/ml
Dinoseb methyl ether	CAS:6099-79-2	200 µg/ml	Silvex methyl ester	CAS:4841-20-7	100 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F113881	Price : FA

Semi Volatile Pesticide Stock Standard - 20 components

Aldrin	CAS:309-00-2	Endosulfan-alpha	CAS:959-98-8	Beta-HCH	CAS:319-85-7
alpha-Chlordane	CAS:5103-71-9	Endosulfan-beta	CAS:33213-65-9	Gamma-HCH (Lindane)	CAS:58-89-9
Gamma-Chlordane	CAS:5103-74-2	Endosulfan-total (sulfate)	CAS:1031-07-8	Delta-HCH	CAS:319-86-8
4,4'-DDD (TDE)	CAS:72-54-8	Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8
4,4'-DDE	CAS:72-55-9	Endrin aldehyde	CAS:7421-93-4	Heptachlor-exo-epoxide	CAS:1024-57-3
4,4'-DDT	CAS:50-29-3	Endrin ketone	CAS:53494-70-5	Methoxychlor (DMTD)	CAS:72-43-5
Dieldrin	CAS:60-57-1	Alpha-HCH	CAS:319-84-6		
Solvent : n-Hexane/Toluene (1/1)		ampoule 1 ml	2 000 µg/ml	Ref : F128761	Price : CKA

TCLP Pesticides Mixture - 5 components**(EPA 1311)**

Gamma-HCH (Lindane)	CAS:58-89-9	100 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	1 000 µg/ml
Endrin	CAS:72-20-8	100 µg/ml	cis-Heptachlorepoxyde	CAS:1024-57-3	100 µg/ml
Heptachlor	CAS:76-44-8	100 µg/ml	(cis, exo-, isomer B)		
Solvent : Methanol		ampoule 1 ml		Ref : F116111	Price : DF

TCLP Herbicides Spiking Mixture - 2 components**(EPA 1311)**

2,4-D	CAS:94-75-7	Silvex	CAS:93-72-1		
Solvent : Methanol		ampoule 1 ml	2 000 µg/ml	Ref : F116151	Price : DF

TCLP methylated Herbicides Spiking Mixture - 2 components**(EPA 1311)**

2,4-D methyl ester	CAS:1928-38-7	Silvex methyl ester	CAS:4841-20-7		
Solvent : Methanol		ampoule 1 ml	2 000 µg/ml	Ref : F116161	Price : DF

TCLP Pesticides Spiking Mixture - 2 components**(EPA 1311)**

Chlordane (technical)	CAS:57-74-9	2 000 µg/ml	Toxaphene	CAS:8001-35-2	4 000 µg/ml
Solvent : Methanol		ampoule 1 ml		Ref : F116141	Price : FA

TCLP Pesticides Spiking Mixture - 5 components**(EPA 1311)**

Endrin	CAS:72-20-8	Heptachlor	CAS:76-44-8	Methoxychlor (DMTD)	CAS:72-43-5
Gamma-HCH (Lindane)	CAS:58-89-9	cis-Heptachlorepoxyde	CAS:1024-57-3		
		(cis-, exo-, isomer B)			
Solvent : Methanol		ampoule 1 ml	2 000 µg/ml	Ref : F116131	Price : GA

TCLP Pesticides Spiking Mixture - 7 components**(EPA 1311)**

Chlordane (technical)	CAS:57-74-9	2 000 µg/ml	Methoxychlor (DMTD)	CAS:72-43-5	2 000 µg/ml
Endrin	CAS:72-20-8	2 000 µg/ml	Toxaphene	CAS:8001-35-2	4 000 µg/ml
Gamma-HCH (Lindane)	CAS:58-89-9	2 000 µg/ml	cis-Heptachlorepoxyde	CAS:1024-57-3	2 000 µg/ml
Heptachlor	CAS:76-44-8	2 000 µg/ml	(cis, exo-, isomer B)		
Solvent : Methanol		ampoule 1 ml		Ref : F116121	Price : GF

Triazine & Urea Pesticide Mixture - 29 components

Atrazine	CAS:1912-24-9	Dimefuron	CAS:34205-21-5	Metolachlor	CAS:51218-45-2
Atrazine-desethyl	CAS:6190-65-4	Diuron	CAS:330-54-1	Metoxuron	CAS:19937-59-8
Atrazine-desisopropyl	CAS:1007-28-9	Ethidimuron	CAS:30043-49-3	Monolinuron	CAS:1746-81-2
Bromacil	CAS:314-40-9	Ethofumesate	CAS:26225-79-6	Prometryn	CAS:7287-19-6
Carbetamide	CAS:16118-49-3	Isoproturon	CAS:34123-59-6	Propazine	CAS:139-40-2
Chloridazon	CAS:1698-60-8	Linuron	CAS:330-55-2	Simazine (CAT)	CAS:122-34-9
Chloroxuron	CAS:1982-47-4	Metamitron	CAS:41394-05-2	Terbutylazine	CAS:5915-41-3
Chlorpropham	CAS:101-21-3	Metazachlor	CAS:67129-08-2	Terbutylazine-desethyl	CAS:30125-63-4
Chlortoluron	CAS:15545-48-9	Methabenzthiazuron	CAS:18691-97-9	Terbutryn	CAS:886-50-0
Cyanazine	CAS:21725-46-2	Metobromuron	CAS:3060-89-7		
Solvent : Acetonitrile		ampoule 1 ml	100 µg/ml	Ref : F128821	Price : HA

Underivatized Analyte - 5 components**(EPA 515)**

2,4-Dichlorophenylacetic acid	CAS:19719-28-9	500 µg/ml	Dinoseb	CAS:88-85-7	4 µg/ml
4,4'-Dibromooctafluorobiphenyl	CAS:10386-84-2	250 µg/ml	4-Nitrophenol	CAS:100-02-7	1 600 µg/ml
3,5-Dichlorobenzoic acid	CAS:51-36-5	600 µg/ml			
Solvent : Methyl-tert butylether		ampoule 1 ml		Ref : F113801	Price : DF

Urea Pesticide Standard Solution - 10 components

Buturon	CAS:3766-60-7	Isoproturon	CAS:34123-59-6	Monolinuron	CAS:1746-81-2
Chlortoluron	CAS:15545-48-9	Methabenzthiazuron	CAS:18691-97-9	Monuron (CMU)	CAS:150-68-5
Diuron	CAS:330-54-1	Metobromuron	CAS:3060-89-7		
Fenuron	CAS:101-42-8	Metoxuron	CAS:19937-59-8		
Solvent : Methanol		ampoule 1 ml	10 µg/ml	Ref : F128791	Price : DF

Phenols**Chlorophenol Standard, Stock Solution - 15 components****(ISO14154:2005)**

2,3-Dichlorophenol	CAS:576-24-9	400 µg/ml	2,3,4,6-Tetrachlorophenol	CAS:58-90-2	600 µg/ml
2,4-Dichlorophenol	CAS:120-83-2	400 µg/ml	2,3,4-Trichlorophenol	CAS:15950-66-0	400 µg/ml
2,5-Dichlorophenol	CAS:583-78-8	400 µg/ml	2,3,5-Trichlorophenol	CAS:933-78-8	400 µg/ml
2,6-Dichlorophenol	CAS:87-65-0	400 µg/ml	2,3,6-Trichlorophenol	CAS:933-75-5	400 µg/ml
3,4-Dichlorophenol	CAS:95-77-2	400 µg/ml	2,4,5-Trichlorophenol	CAS:95-95-4	400 µg/ml
3,5-Dichlorophenol	CAS:591-35-5	400 µg/ml	2,4,6-Trichlorophenol	CAS:88-06-2	600 µg/ml
Pentachlorophenol	CAS:87-86-5	1 000 µg/ml	3,4,5-Trichlorophenol	CAS:609-19-8	200 µg/ml
2,3,4,5-Tetrachlorophenol	CAS:4901-51-3	200 µg/ml			
Solvent : Ethanol		ampoule 1 ml		Ref : F109881	Price : HF

Chlorophenyl Acetate Standard, Stock Solution - 15 components**(ISO14154:2005)**

2,3-Dichlorophenol acetate	CAS:61925-85-7	503 µg/ml	2,3,6-Trichlorophenol acetate	CAS:61925-87-9	485 µg/ml
2,4-Dichlorophenol acetate	CAS:6341-97-5	503 µg/ml	2,4,5-Trichlorophenol acetate	CAS:5393-75-9	485 µg/ml
2,5-Dichlorophenol acetate	CAS:30124-46-0	503 µg/ml	2,4,6-Trichlorophenol acetate	CAS:23399-90-8	728 µg/ml
2,6-Dichlorophenol acetate	CAS:28165-71-1	503 µg/ml	3,4,5-Trichlorophenol acetate	(id:15591)	243 µg/ml
3,4-Dichlorophenol acetate	CAS:17847-51-7	503 µg/ml	2,3,4,5-Tetrachlorophenol acetate	(id:14700)	250 µg/ml
3,5-Dichlorophenol acetate	CAS:61925-86-8	503 µg/ml	2,3,4,6-Tetrachlorophenol acetate	(id:14720)	750 µg/ml
2,3,4-Trichlorophenol acetate	(id:15561)	485 µg/ml	Pentachlorophenol acetate	CAS:1441-02-7	1 158 µg/ml
2,3,5-Trichlorophenol acetate	CAS:61925-88-0	485 µg/ml			
Solvent : Ethanol		ampoule 1 ml		Ref : F109901	Price : CBA

Phenols Internal Standard Solution - 2 components**(ISO 17495:2001)**

2,4-Dibromophenol	CAS:615-58-7	2,3,6-Trichlorophenol	CAS:933-75-5		
Solvent : Ethyl Acetate		ampoule 1 ml	10 µg/ml	Ref : F129061	Price : CF

Phenols Standard Solution - 3 components**(EPA 8041)**

2-Cyclohexyl-4,6-dinitrophenol	CAS:131-89-5	2,3,4,5-Tetrachlorophenol	CAS:4901-51-3	2,3,5,6-Tetrachlorophenol	CAS:935-95-5
Solvent : 2-Propanol		ampoule 1 ml	2 000 µg/ml	Ref : F116511	Price : HJ

Phenols Standard Solution - 9 components**(EPA 8041)**

4-Chloro-3-methylphenol	CAS:59-50-7	2-Methyl-4,6-dinitrophenol	CAS:534-52-1	Pentachlorophenol	CAS:87-86-5
o-Cresol	CAS:95-48-7	2-Nitrophenol	CAS:88-75-5	Phenol	CAS:108-95-2
2,4-Dichlorophenol	CAS:120-83-2	4-Nitrophenol	CAS:100-02-7	2,4,6-Trichlorophenol	CAS:88-06-2
Solvent : 2-Propanol		ampoule 1 ml	2 000 µg/ml	Ref : F116501	Price : EF

Phenols Standard Solution - 9 components**(EPA 8041)**

2-Chlorophenol	CAS:95-57-8	2,6-Dichlorophenol	CAS:87-65-0	Dinoseb	CAS:88-85-7
m-Cresol	CAS:108-39-4	2,4-Dimethylphenol	CAS:105-67-9	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
p-Cresol	CAS:106-44-5	2,4-Dinitrophenol	CAS:51-28-5	2,4,5-Trichlorophenol	CAS:95-95-4
Solvent : 2-Propanol		ampoule 1 ml	2 000 µg/ml	Ref : F116521	Price : FA

Phenols Standard Solution - 10 components

4-Chloro-3-methylphenol	CAS:59-50-7	2,4-Dimethylphenol	CAS:105-67-9	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
p-Cresol	CAS:106-44-5	2,6-Dimethylphenol	CAS:576-26-1	2,4,6-Trichlorophenol	CAS:88-06-2
2,4-Dichlorophenol	CAS:120-83-2	Pentachlorophenol	CAS:87-86-5		
2,6-Dichlorophenol	CAS:87-65-0	Phenol	CAS:108-95-2		
Solvent : Methanol		ampoule 1 ml	200 µg/ml	Ref : F128811	Price : EH

Phenols Standard Solution - 11 components**(EPA 604)**

4-Chloro-3-methylphenol	CAS:59-50-7	2,4-Dinitrophenol	CAS:51-28-5	Pentachlorophenol	CAS:87-86-5
2-Chlorophenol	CAS:95-57-8	2-Methyl-4,6-dinitrophenol	CAS:534-52-1	Phenol	CAS:108-95-2
2,4-Dichlorophenol	CAS:120-83-2	2-Nitrophenol	CAS:88-75-5	2,4,6-Trichlorophenol	CAS:88-06-2
2,4-Dimethylphenol	CAS:105-67-9	4-Nitrophenol	CAS:100-02-7		
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F114911	Price : EA

Phenols Standard Solution - 11 components**(EPA 625)**

4-Chloro-3-methylphenol	CAS:59-50-7	2,4-Dinitrophenol	CAS:51-28-5	Pentachlorophenol	CAS:87-86-5
2-Chlorophenol	CAS:95-57-8	2-Methyl-4,6-dinitrophenol	CAS:534-52-1	Phenol	CAS:108-95-2
2,4-Dichlorophenol	CAS:120-83-2	2-Nitrophenol	CAS:88-75-5	2,4,6-Trichlorophenol	CAS:88-06-2
2,4-Dimethylphenol	CAS:105-67-9	4-Nitrophenol	CAS:100-02-7		
Solvent : Dichloromethane		ampoule 1 ml	2 000 µg/ml	Ref : F115881	Price : FA

Phenols Standard Solution - 14 components**(ISO 17495:2001)**

2,4-Dichloro-6-nitrophenol	CAS:609-89-2	2,4-Dinitro-6-methylphenol	CAS:534-52-1	5-Methyl-2-nitrophenol	CAS:700-38-9
2,6-Dichloro-4-nitrophenol	CAS:618-80-4	2,4-Dinitrophenol	CAS:51-28-5	2-Nitrophenol	CAS:88-75-5
2,6-Dimethyl-4-nitrophenol	CAS:2423-71-4	4-Methyl-2-nitrophenol	CAS:119-33-5	3-Nitrophenol	CAS:554-84-7
2,5-Dinitrophenol	CAS:329-71-5	3-Methyl-2-nitrophenol	CAS:4920-77-8	4-Nitrophenol	CAS:100-02-7
2,6-Dinitrophenol	CAS:573-56-8	3-Methyl-4-nitrophenol	CAS:2581-34-2		
Solvent : Ethyl Acetate		ampoule 1 ml	500 µg/ml	Ref : F129051	Price : DJ

Phenols Standard Solution - 14 components

2,3-Dimethylphenol	CAS:526-75-0	3,5-Dimethylphenol	CAS:108-68-9	3-Methylphenol	CAS:108-39-4
2,4-Dimethylphenol	CAS:105-67-9	2-Ethylphenol	CAS:90-00-6	4-Methylphenol	CAS:106-44-5
2,5-Dimethylphenol	CAS:95-87-4	3-Ethylphenol	CAS:620-17-7	Phenol	CAS:108-95-2
2,6-Dimethylphenol	CAS:576-26-1	4-Ethylphenol	CAS:123-07-9	Thymol	CAS:89-83-8
3,4-Dimethylphenol	CAS:95-65-8	2-Methylphenol	CAS:95-48-7		
Solvent : Dichloromethane		ampoule 1 ml	1 000 µg/ml	Ref : F112411	Price : FD

Phenols Standard Solution - 17 components**(EPA 8041)**

4-Chloro-3-methylphenol	CAS:59-50-7	2,6-Dichlorophenol	CAS:87-65-0	Pentachlorophenol	CAS:87-86-5
2-Chlorophenol	CAS:95-57-8	2,4-Dimethylphenol	CAS:105-67-9	Phenol	CAS:108-95-2
o-Cresol	CAS:95-48-7	2,4-Dinitrophenol	CAS:51-28-5	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
m-Cresol	CAS:108-39-4	2-Methyl-4,6-dinitrophenol	CAS:534-52-1	2,4,5-Trichlorophenol	CAS:95-95-4
p-Cresol	CAS:106-44-5	2-Nitrophenol	CAS:88-75-5	2,4,6-Trichlorophenol	CAS:88-06-2
2,4-Dichlorophenol	CAS:120-83-2	4-Nitrophenol	CAS:100-02-7		
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F116571	Price : FJ

Phenols Standard Solution - 19 components

2-Amino-4-chlorophenol	CAS:95-85-2	2,4-Dichlorophenol	CAS:120-83-2	4-n-Octylphenol	CAS:1806-26-4
4-tert-Butylphenol	CAS:98-54-4	2,5-Dichlorophenol	CAS:583-78-8	4-tert-Octylphenol	CAS:140-66-9
4-Chloro-3-methylphenol	CAS:59-50-7	2,6-Dichlorophenol	CAS:87-65-0	Pentachlorophenol	CAS:87-86-5
2-Chlorophenol	CAS:95-57-8	3,4-Dichlorophenol	CAS:95-77-2	2,4,5-Trichlorophenol	CAS:95-95-4
3-Chlorophenol	CAS:108-43-0	3,5-Dichlorophenol	CAS:591-35-5	2,4,6-Trichlorophenol	CAS:88-06-2
4-Chlorophenol	CAS:106-48-9	Nonylphenol	CAS:25154-52-3		
2,3-Dichlorophenol	CAS:576-24-9	4-n-nonylphenol	CAS:104-40-5		
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F062365	Price : DF
		ampoule 5 ml	100 µg/ml	Ref : F112431	Price : DF

Phenols Standard Solution - 19 components**(EN 12673)**

2-Chlorophenol	CAS:95-57-8	3 000 µg/l	2,3,4,5-Tetrachlorophenol acetate	id : 14700	200 µg/l
3-Chlorophenol	CAS:108-43-0	3 000 µg/l	2,3,4,6-Tetrachlorophenol acetate	id : 14 720	200 µg/l
4-Chlorophenol	CAS:106-48-9	3 000 µg/l	2,3,5,6-Tetrachlorophenol Acetate	id : 14730	200 µg/l
2,3-Dichlorophenol	CAS:576-24-9	400 µg/l	2,3,4-Trichlorophenol acetate	id : 15561	300 µg/l
2,4-Dichlorophenol	CAS:120-83-2	400 µg/l	2,3,5-Trichlorophenol acetate	CAS:61925-88-0	300 µg/l
2,5-Dichlorophenol	CAS:583-78-8	400 µg/l	2,3,6-Trichlorophenol acetate	CAS:61925-87-9	300 µg/l
2,6-Dichlorophenol	CAS:87-65-0	400 µg/l	2,4,5-Trichlorophenol acetate	CAS:5393-75-9	300 µg/l
3,4-Dichlorophenol	CAS:95-77-2	400 µg/l	2,4,6-Trichlorophenol acetate	CAS:23399-90-8	300 µg/l
3,5-Dichlorophenol	CAS:591-35-5	400 µg/l	3,4,5-Trichlorophenol acetate	id : 15591	300 µg/l
Pentachlorophenol	CAS:87-86-5	100 µg/l			
Solvent : Ethanol		ampoule 1 ml		Ref : F131321	Price : BCA

Phenols Standard Solution - 20 components

2-Chlorophenol	CAS:95-57-8	2,6-Dichlorophenol	CAS:87-65-0	2,3,4-Trichlorophenol	CAS:15950-66-0
3-Chlorophenol	CAS:108-43-0	3,4-Dichlorophenol	CAS:95-77-2	2,3,5-Trichlorophenol	CAS:933-78-8
4-Chlorophenol	CAS:106-48-9	3,5-Dichlorophenol	CAS:591-35-5	2,3,6-Trichlorophenol	CAS:933-75-5
4-Chloro-3-methylphenol	CAS:59-50-7	Pentachlorophenol	CAS:87-86-5	2,4,5-Trichlorophenol	CAS:95-95-4
2,3-Dichlorophenol	CAS:576-24-9	2,3,4,5-Tetrachlorophenol	CAS:4901-51-3	2,4,6-Trichlorophenol	CAS:88-06-2
2,4-Dichlorophenol	CAS:120-83-2	2,3,4,6-Tetrachlorophenol	CAS:58-90-2	3,4,5-Trichlorophenol	CAS:609-19-8
2,5-Dichlorophenol	CAS:583-78-8	2,3,5,6-Tetrachlorophenol	CAS:935-95-5		
Solvent : Dichloromethane		ampoule 1 ml	1 000 µg/ml	Ref : F112421	Price : BAK

Phenols Standard Solution - 21 components**(EPA 8041)**

4-Chloro-3-methylphenol	CAS:59-50-7	2,6-Dichlorophenol	CAS:87-65-0	Pentachlorophenol	CAS:87-86-5
2-Chlorophenol	CAS:95-57-8	2,4-Dimethylphenol	CAS:105-67-9	Phenol	CAS:108-95-2
o-Cresol	CAS:95-48-7	Dinoseb	CAS:88-85-7	2,3,4,5-Tetrachlorophenol	CAS:4901-51-3
m-Cresol	CAS:108-39-4	2,4-Dinitrophenol	CAS:51-28-5	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
p-Cresol	CAS:106-44-5	2-Methyl-4,6-dinitrophenol	CAS:534-52-1	2,3,5,6-Tetrachlorophenol	CAS:935-95-5
2-Cyclohexyl-4,6-dinitrophenol	CAS:131-89-5	2-Nitrophenol	CAS:88-75-5	2,4,6-Trichlorophenol	CAS:88-06-2
2,4-Dichlorophenol	CAS:120-83-2	4-Nitrophenol	CAS:100-02-7	2,4,5-Trichlorophenol	CAS:95-95-4
Solvent : 2-Propanol		ampoule 1 ml	100 µg/ml	Ref : F116561	Price : HF

Phenols Standard Stock Solution - 26 components**(ISO 8165-2:1999)**

2-Benzylphenol	CAS:28994-41-4	2,5-Dichlorophenol	CAS:583-78-8	Pentachlorophenol	CAS:87-86-5
2-Chlorophenol	CAS:95-57-8	2,6-Dichlorophenol	CAS:87-65-0	2,4,5-Trichlorophenol	CAS:95-95-4
3-Chlorophenol	CAS:108-43-0	2,4-Dimethylphenol	CAS:105-67-9	2,4,6-Trichlorophenol	CAS:88-06-2
4-Chlorophenol	CAS:106-48-9	4-Ethylphenol	CAS:123-07-9	2,3,5-Trichlorophenol	CAS:933-78-8
4-Chloro-2-methylphenol	CAS:1570-64-5	2-Methylphenol	CAS:95-48-7	2,3,6-Trichlorophenol	CAS:933-75-5
4-Chloro-3-methylphenol	CAS:59-50-7	3-Methylphenol	CAS:108-39-4	2,3,4,5-Tetrachlorophenol	CAS:4901-51-3
2,4-Dichloro-3,5-dimethylphenol	CAS:133-53-9	4-Methylphenol	CAS:106-44-5	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
2,3-Dichlorophenol	CAS:576-24-9	Phenol	CAS:108-95-2	2,3,5,6-Tetrachlorophenol	CAS:935-95-5
2,4-Dichlorophenol	CAS:120-83-2	2-Phenylphenol	CAS:90-43-7		
Solvent : Methanol		ampoule 1 ml	30 µg/ml	Ref : F129181	Price : HF

Phenols Standard Stock Solution - 31 components**(ISO 8165-2:1999)**

2-Benzyl-4-chlorophenol	CAS:120-32-1	2,4-Dichloro-3,5-dimethylphenol	CAS:133-53-9	Pheno	CAS:108-95-2
2-Benzylphenol	CAS:28994-41-4	2,3-Dichlorophenol	CAS:576-24-9	2-Phenylphenol	CAS:90-43-7
4-tert-Butyl-2-chlorophenol	CAS:98-28-2	2,4-Dichlorophenol	CAS:120-83-2	2,3,4,5-Tetrachlorophenol	CAS:4901-51-3
2,6-Di-tert-butyl-4-methylphenol	CAS:128-37-0	2,5-Dichlorophenol	CAS:583-78-8	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
2-Chlorophenol	CAS:95-57-8	2,6-Dichlorophenol	CAS:87-65-0	2,3,5,6-Tetrachlorophenol	CAS:935-95-5
2-Chloro-5-methylphenol	CAS:615-74-7	2,4-Dimethylphenol	CAS:105-67-9	2,3,5-Trichlorophenol	CAS:933-78-8
3-Chlorophenol	CAS:108-43-0	4-Ethylphenol	CAS:123-07-9	2,3,6-Trichlorophenol	CAS:933-75-5
4-Chlorophenol	CAS:106-48-9	2-Methylphenol	CAS:95-48-7	2,4,5-Trichlorophenol	CAS:95-95-4
4-Chloro-2-isopropyl-5-methylphenol	CAS:89-68-9	3-Methylphenol	CAS:108-39-4	2,4,6-Trichlorophenol	CAS:88-06-2
4-Chloro-2-methylphenol	CAS:1570-64-5	4-Methylphenol	CAS:106-44-5		
4-Chloro-3-methylphenol	CAS:59-50-7	Pentachlorophenol	CAS:87-86-5		
Solvent : Methanol		ampoule 1 ml	300 µg/ml	Ref : F109300	Price : BBA

Phenol Standard Stock Solution - 33 components**(ISO 8165-1:1992)**

Phenol	CAS:108-95-2	2,6-Di-tert-butyl-4-methylphenol	CAS:128-37-0	1-Naphthol	CAS:90-15-3
2-Benzylphenol	CAS:28994-41-4	2,4-Dichloro-3,5-dimethylphenol	CAS:133-53-9	2-Naphthol	CAS:135-19-3
2-Benzyl-4-chlorophenol	CAS:120-32-1	2,3-Dichlorophenol	CAS:576-24-9	Pentachlorophenol	CAS:87-86-5
4-tert-Butyl-2-chlorophenol	CAS:98-28-2	2,4-Dichlorophenol	CAS:120-83-2	2-Phenylphenol	CAS:90-43-7
2-Chlorophenol	CAS:95-57-8	2,5-Dichlorophenol	CAS:583-78-8	2,3,4,5-Tetrachlorophenol	CAS:4901-51-3
3-Chlorophenol	CAS:108-43-0	2,6-Dichlorophenol	CAS:87-65-0	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
4-Chlorophenol	CAS:106-48-9	2,4-Dimethylphenol	CAS:105-67-9	2,3,5,6-Tetrachlorophenol	CAS:935-95-5
4-Chloro-2-isopropyl-5-methylphenol	CAS:89-68-9	4-Ethylphenol	CAS:123-07-9	2,3,5-Trichlorophenol	CAS:933-78-8
4-Chloro-2-methylphenol	CAS:1570-64-5	2-Methylphenol	CAS:95-48-7	2,3,6-Trichlorophenol	CAS:933-75-5
4-Chloro-3-methylphenol	CAS:59-50-7	3-Methylphenol	CAS:108-39-4	2,4,5-Trichlorophenol	CAS:95-95-4
2-Chloro-5-methylphenol	CAS:615-74-7	4-Methylphenol	CAS:106-44-5	2,4,6-Trichlorophenol	CAS:88-06-2
Solvent : Methanol		ampoule 1 ml	100 µg/ml	Ref : F109270	Price : KA

Phenols Surrogate Standard Mixture - 2 components**(EPA 8041)**

2-Fluorophenol	CAS:367-12-4	2,4,6-Tribromophenol	CAS:118-79-6		
Solvent : 2-Propanol		ampoule 1 ml	2 000 µg/ml	Ref : F116541	Price : DA

TCLP Acids Mixture - 6 components**(EPA 1311)**

2-Methylphenol	CAS:95-48-7	4-Methylphenol	CAS:106-44-5	2,4,5-Trichlorophenol	CAS:95-95-4
3-Methylphenol	CAS:108-39-4	Pentachlorophenol	CAS:87-86-5	2,4,6-Trichlorophenol	CAS:88-06-2
Solvent : Methanol		ampoule 1 ml	1 000 µg/ml	Ref : F116081	Price : EA

VOC / Phenols Standard Solution - 25 components

Benzoic acid	CAS:65-85-0	4,6-Dinitro-2-methylphenol	CAS:534-52-1	4-Nitrophenol	CAS:100-02-7
Benzyl alcohol	CAS:100-51-6	2,4-Dinitrophenol	CAS:51-28-5	Pentachlorophenol	CAS:87-86-5
4-Chloro-3-methylphenol	CAS:59-50-7	Ethyl methanesulfonate	CAS:62-50-0	Phenol	CAS:108-95-2
2-Chlorophenol	CAS:95-57-8	Hexachlorophen	CAS:70-30-4	Resorcinol	CAS:108-46-3
2-Cyclohexyl-4,6-dinitrophenol	CAS:131-89-5	Methyl methanesulfonate	CAS:66-27-3	2,3,4,6-Tetrachlorophenol	CAS:58-90-2
2,4-Dichlorophenol	CAS:120-83-2	2-Methylphenol	CAS:95-48-7	2,4,5-Trichlorophenol	CAS:95-95-4
2,6-Dichlorophenol	CAS:87-65-0	3-Methylphenol	CAS:108-39-4	2,4,6-Trichlorophenol	CAS:88-06-2
Diethyl sulfate	CAS:64-67-5	4-Methylphenol	CAS:106-44-5		
2,4-Dimethylphenol	CAS:105-67-9	2-Nitrophenol	CAS:88-75-5		
Solvent : Dichloromethane/Methanol P & T (1/1)		ampoule 1 ml	100 µg/ml	Ref : F119321	Price : GA

VOC / PAH / Phthalates Standard Solution - 58 components**(EPA 8270C)**

Acenaphthene	CAS:83-32-9	Dibenzo(a,e)pyrene	CAS:192-65-4	Hexachloroethane	CAS:67-72-1
Acenaphthylene	CAS:208-96-8	1,2-Dibromo-3-chloropropane	CAS:96-12-8	Hexachloropropene	CAS:1888-71-7
Anthracene	CAS:120-12-7	Di-n-butylphthalate	CAS:84-74-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Benzo(a)anthracene	CAS:56-55-3	1,2-Dichlorobenzene	CAS:95-50-1	Isophorone	CAS:78-59-1
Benzo(b)fluoranthene	CAS:205-99-2	1,3-Dichlorobenzene	CAS:541-73-1	3-Methylcholanthrene	CAS:56-49-5
Benzo(k)fluoranthene	CAS:207-08-9	1,4-Dichlorobenzene	CAS:106-46-7	2-Methylnaphthalene	CAS:91-57-6
Benzo(g,h,i)perylene	CAS:191-24-2	Diethyl phthalate	CAS:84-66-2	Naphthalene	CAS:91-20-3
Benzo(a)pyrene	CAS:50-32-8	7,12-Dimethylbenzo(a)anthracene	CAS:57-97-6	5-Nitroacenaphthene	CAS:602-87-9
Bis-(2-chloroethoxy)-methane	CAS:111-91-1	Dimethyl phthalate	CAS:131-11-3	Nitrobenzene	CAS:98-95-3
Bis-(2-chloroethyl)-ether	CAS:111-44-4	1,2-Dinitrobenzene	CAS:528-29-0	4-Nitrobiphenyl	CAS:92-93-3
Bis-(2-chloro-1-methylethyl)ether	CAS:108-60-1	1,3-Dinitrobenzene	CAS:99-65-0	Pentachlorobenzene	CAS:608-93-5
4-Bromophenyl-phenyl ether	CAS:101-55-3	1,4-Dinitrobenzene	CAS:100-25-4	Pentachloronitrobenzene	CAS:82-68-8
Butyl benzyl phthalate	CAS:85-68-7	2,4-Dinitrotoluene	CAS:121-14-2	Phenanthrene	CAS:85-01-8
1-Chloronaphthalene	CAS:90-13-1	2,6-Dinitrotoluene	CAS:606-20-2	Pyrene	CAS:129-00-0
2-Chloronaphthalene	CAS:91-58-7	Di-n-octylphthalate	CAS:117-84-0	1,2,4,5-Tetrachlorobenzene	CAS:95-94-3
4-Chlorophenyl phenyl ether	CAS:7005-72-3	Fluoranthene	CAS:206-44-0	1,2,4-Trichlorobenzene	CAS:120-82-1
Chrysene	CAS:218-01-9	Fluorene	CAS:86-73-7	1,3,5-Trinitrobenzene	CAS:99-35-4
Dibenz(a,j)acridine	CAS:224-42-0	Hexachlorobenzene	CAS:118-74-1	Phthalic acid,bis-2-ethylhexylester	CAS:117-81-7
Dibenzo(a,h)anthracene	CAS:53-70-3	Hexachlorobutadiene	CAS:87-68-3		
Dibenzofuran	CAS:132-64-9	Hexachlorocyclopentadiene	CAS:77-47-4		
Solvent : Dichloromethane/Benzene (3/1)		ampoule 1 ml	100 µg/ml	Ref : F119311	Price : BCA

Phthalates**Extractables Mixture - 9 components****(EPA 525.1)**

Bis(2-ethylhexyl)adipate	CAS:103-23-1	Di-n-butylphthalate	CAS:84-74-2	HCB	CAS:118-74-1
Bis(2-ethylhexyl)phthalate	CAS:117-81-7	Diethyl phthalate	CAS:84-66-2	Hexachlorocyclopentadiene	CAS:77-47-4
Butyl benzyl phthalate	CAS:85-68-7	Dimethyl phthalate	CAS:131-11-3	Pentachlorophenol	CAS:87-86-5
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F114011	Price : EF
			(Pentachlorophenol : 400 µg/ml)		
		ampoule 1 ml	500 µg/ml	Ref : F114021	Price : EF
			(Pentachlorophenol : 2 000 µg/ml)		

Phthalates Matrix Spike Solution - 2 components**(EPA 8061)**

Bis(2-ethylhexyl)phthalate	CAS:117-81-7	Butyl benzyl phthalate	CAS:85-68-7		
Solvent : Acetone		ampoule 1 ml	2 000 µg/ml	Ref : F116651	Price : CF

Phthalates QC Check Sample Concentrate - 6 components**(EPA 606)**

Bis(2-ethylhexyl)phthalate	CAS:117-81-7	50 µg/ml	Diethyl phthalate	CAS:84-66-2	25 µg/ml
Butyl benzyl phthalate	CAS:85-68-7	10 µg/ml	Dimethyl phthalate	CAS:131-11-3	25 µg/ml
di-n-butyl phthalate	CAS:84-74-2	25 µg/ml	di-n-octyl phthalate	CAS:117-84-0	50 µg/ml
Solvent : Acetone		ampoule 1 ml		Ref : F115031	Price : DF

Phthalates Standard Solution - 2 components

Phthalic acid, bis-2-ethylhexylester	CAS:117-81-7	Phthalic acid, bis-1-octyl ester	CAS:117-84-0		
Solvent : Toluene		ampoule 5 ml	100 µg/ml	Ref : F112391	Price : EA

Phthalates Standard Solution - 6 components**(EPA 606)**

Bis(2-ethylhexyl)phthalate	CAS:117-81-7	di-n-butyl phthalate	CAS:84-74-2	Dimethyl phthalate	CAS:131-11-3
Butyl benzyl phthalate	CAS:85-68-7	Diethyl phthalate	CAS:84-66-2	di-n-octyl phthalate	CAS:117-84-0
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F115021	Price : DF

Phthalates Standard Solution - 6 components**(EPA 8061)**

Phthalic acid, bis-2-ethylhexylester	CAS:117-81-7	Di-n-butylphthalate	CAS:84-74-2	Dimethyl phthalate	CAS:131-11-3
Phthalic acid, benzylbutyl ester	CAS:85-68-7	Diethyl phthalate	CAS:84-66-2	Di-n-octylphthalate	CAS:117-84-0
Solvent : n-Hexane		ampoule 1 ml	100 µg/ml	Ref : F116631	Price : EF
		ampoule 1 ml	1 000 µg/ml	Ref : F116621	Price : EF

Phthalates Standard Solution - 7 components**(EPA 506)**

Adipic acid-bis-2-ethylhexyl ester	CAS:103-23-1	Phthalic acid, bis-ethyl ester	CAS:84-66-2	Phthalic acid, bis-1-octyl ester	CAS:117-84-0
Phthalic acid, benzylbutyl ester	CAS:85-68-7	Phthalic acid, bis-2-ethylhexylester	CAS:117-81-7		
Phthalic acid, bis-butyl ester	CAS:84-74-2	Phthalic acid, bis-methyl ester	CAS:131-11-3		
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F113281	Price : EA

Phthalates Standard Solution - 7 components

Adipic acid-bis-2-ethylhexyl ester	CAS:103-23-1	Phthalic acid, bis-butyl ester	CAS:84-74-2	Phthalic acid, bis-1-octyl ester	CAS:117-84-0
Diisononyl phthalate	CAS:28553-12-0	Phthalic acid, bis-ethyl ester	CAS:84-66-2		
Phthalic acid, benzylbutyl ester	CAS:85-68-7	Phthalic acid, bis-2-ethylhexylester	CAS:117-81-7		
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F112381	Price : EC

Phthalates Standard Solution - 7 components**(EPA 506)**

Adipic acid-bis-2-ethylhexyl ester	CAS:103-23-1	1 200 µg/ml	Phthalic acid, bis-2-ethylhexylester	CAS:117-81-7	250 µg/ml
Phthalic acid, benzylbutyl ester	CAS:85-68-7	250 µg/ml	Phthalic acid, bis-methyl ester	CAS:131-11-3	100 µg/ml
Phthalic acid, bis-butyl ester	CAS:84-74-2	100 µg/ml	Phthalic acid, bis-1-octyl ester	CAS:117-84-0	650 µg/ml
Phthalic acid, bis-ethyl ester	CAS:84-66-2	100 µg/ml			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F113311	Price : EF

Phthalates Standard Solution - 11 components**(ISO 18856:2004)**

Butyl benzyl phthalate	CAS:85-68-7	Di-2-ethylhexylphthalate	CAS:117-81-7	Di-n-octylphthalate	CAS:117-84-0
Dibutylphthalate	CAS:84-74-2	Diethyl phthalate	CAS:84-66-2	Diundecyl phthalate	CAS:3648-20-2
Dicyclohexyl phthalate	CAS:84-61-7	Di-iso-butylphthalate	CAS:84-69-5	Phthalic acid, bis-propyl ester	CAS:131-16-8
Didecyl phthalate	CAS:84-77-5	Dimethyl phthalate	CAS:131-11-3		
Solvent : Ethyl Acetate		ampoule 1 ml	1 000 µg/ml	Ref : F129191	Price : DC

Phthalates Surrogate Standard Solution - 3 components**(EPA 8061)**

Dibenzyl phthalate	CAS:523-31-9	Diphenyl isophthalate	id : 16211	Diphenyl phthalate	CAS:84-62-8
Solvent : Acetone		ampoule 1 ml	50 µg/ml	Ref : F116681	Price : DA
		ampoule 1 ml	500 µg/ml	Ref : F116671	Price : DA

Polychlorinated Biphenyls (PCBs) Congeners Mixtures**PCBs Standard Solution - 6 components**

PCB 28	CAS:7012-37-5	PCB 101	CAS:37680-73-2	PCB 153	CAS:35065-27-1
PCB 52	CAS:35693-99-3	PCB 138	CAS:35065-28-2	PCB 180	CAS:35065-29-3
Solvent : Iso-Octane		ampoule 10 ml	10 µg/ml	Ref : F128561	Price : FG

PCBs Standard Solution - 7 components

PCB 28	CAS:7012-37-5	PCB 101	CAS:37680-73-2	PCB 153	CAS:35065-27-1
PCB 52	CAS:35693-99-3	PCB 138	CAS:35065-28-2	PCB 180	CAS:35065-29-3
PCB 209	CAS:2051-24-3				
Solvent : Iso-Octane		ampoule 10 ml	10 µg/ml	Ref : F128571	Price : FG

PCBs Standard Solution - 7 components

PCB 77	CAS:32598-13-3	PCB 114	CAS:74472-37-0	PCB 126	CAS:57465-28-8
PCB 81	CAS:70362-50-4	PCB 118	CAS:31508-00-6	PCB 156	CAS:38380-08-4
PCB 105	CAS:32598-14-4				
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F112631	Price : EF

PCBs Standard Solution - 8 components

PCB 28	CAS:7012-37-5	PCB 118	CAS:31508-00-6	PCB 180	CAS:35065-29-3
PCB 52	CAS:35693-99-3	PCB 138	CAS:35065-28-2	PCB 194	CAS:35694-08-7
PCB 101	CAS:37680-73-2	PCB 153	CAS:35065-27-1		
Solvent : n-Hexane		ampoule 1 ml	100 µg/ml	Ref : F112661	Price : FA

PCBs Standard Solution - 8 components

PCB 28	CAS:7012-37-5	PCB 118	CAS:31508-00-6	PCB 180	CAS:35065-29-3
PCB 52	CAS:35693-99-3	PCB 138	CAS:35065-28-2	PCB 194	CAS:35694-08-7
PCB 101	CAS:37680-73-2	PCB 153	CAS:35065-27-1		
Solvent : n-Hexane		ampoule 5 ml	1 µg/ml	Ref : F112651	Price : FC

PCBs Standard Solution - 15 components

PCB 18	CAS:37680-65-2	PCB 52	CAS:35693-99-3	PCB 149	CAS:38380-04-0
PCB 20	CAS:38444-84-7	PCB 101	CAS:37680-73-2	PCB 153	CAS:35065-27-1
PCB 28	CAS:7012-37-5	PCB 105	CAS:32598-14-4	PCB 170	CAS:35065-30-6
PCB 31	CAS:16606-02-3	PCB 118	CAS:31508-00-6	PCB 180	CAS:35065-29-3
PCB 44	CAS:41464-39-5	PCB 138	CAS:35065-28-2	PCB 194	CAS:35694-08-7
Solvent : n-Hexane		ampoule 5 ml	10 µg/ml	Ref : F128581	Price : BAC

PCBs Standard Solution - 19 components

PCB 1	CAS:2051-60-7	PCB 87	CAS:38380-02-8	PCB 170	CAS:35065-30-6
PCB 4	CAS:16605-91-7	PCB 101	CAS:37680-73-2	PCB 180	CAS:35065-29-3
PCB 18	CAS:37680-65-2	PCB 110	CAS:38380-03-9	PCB 181	CAS:74472-47-2
PCB 31	CAS:16606-02-3	PCB 138	CAS:35065-28-2	PCB 187	CAS:52663-68-0
PCB 44	CAS:41464-39-5	PCB 141	CAS:52712-04-6	PCB 206	CAS:40186-72-9
PCB 52	CAS:35693-99-3	PCB 151	CAS:52663-63-5		
PCB 66	CAS:32598-10-0	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 1 ml	100 µg/ml	Ref : F127821	Price : BCF

WHO Dioxine like PCB mixtures**ISS PCB Mixture – 18 components**

PCB 28	CAS:7012-37-5	PCB 110	CAS:38380-03-9	PCB 153	CAS:35065-27-1
PCB 52	CAS:35693-99-3	PCB 118	CAS:31508-00-6	PCB 170	CAS:35065-30-6
PCB 95	CAS:38379-99-6	PCB 138	CAS:35065-28-2	PCB 177	CAS:52663-70-4
PCB 99	CAS:38380-01-7	PCB 146	CAS:51908-16-8	PCB 180	CAS:35065-29-3
PCB 101	CAS:37680-73-2	PCB 149	CAS:38380-04-0	PCB 183	CAS:52663-69-1
PCB 105	CAS:32598-14-4	PCB 151	CAS:52663-63-5	PCB 187	CAS:52663-68-0
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F127881	Price : BDA

PCBs Standard Solution - 12 components

PCB 77	CAS:32598-13-3	PCB 118	CAS:31508-00-6	PCB 157	CAS:69782-90-7
PCB 81	CAS:70362-50-4	PCB 123	CAS:65510-44-3	PCB 167	CAS:52663-72-6
PCB 105	CAS:32598-14-4	PCB 126	CAS:57465-28-8	PCB 169	CAS:32774-16-6
PCB 114	CAS:74472-37-0	PCB 156	CAS:38380-08-4	PCB 189	CAS:39635-31-9
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F127871	Price : HF
		ampoule 10 ml	2 µg/ml	Ref : F112671	Price : JC

WHO/ISS PCB Mixtures**PCB Calibration mixture (acc. Guideline for Waste Fuels 2007/0697/A) - 7 components**

PCB 28	CAS:7012-37-5	PCB 118	CAS:31508-00-6	PCB 180	CAS:35065-29-3
PCB 52	CAS:35693-99-3	PCB 138	CAS:35065-28-2		
PCB 101	CAS:37680-73-2	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F128531	Price : EJ

PCBs Calibration Mixture without PCB 170 (acc.method 2) - 13 components

PCB 18	CAS:37680-65-2	PCB 101	CAS:37680-73-2	PCB180	CAS:35065-29-3
PCB 28	CAS:7012-37-5	PCB 118	CAS:31508-00-6	PCB194	CAS:35694-08-7
PCB 31	CAS:16606-02-3	PCB 138	CAS:35065-28-2	PCB 209	CAS:2051-24-3
PCB 44	CAS:41464-39-5	PCB 149	CAS:38380-04-0		
PCB 52	CAS:35693-99-3	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F131470	Price : JA

PCBs Calibration Mixture - 14 components**(EN12766 ; EN 61619)**

PCB 18	CAS:37680-65-2	PCB 101	CAS:37680-73-2	PCB 170	CAS:35065-30-6
PCB 28	CAS:7012-37-5	PCB 118	CAS:31508-00-6	PCB180	CAS:35065-29-3
PCB 31	CAS:16606-02-3	PCB 138	CAS:35065-28-2	PCB194	CAS:35694-08-7
PCB 44	CAS:41464-39-5	PCB 149	CAS:38380-04-0	PCB 209	CAS:2051-24-3
PCB 52	CAS:35693-99-3	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F127901	Price : JC

Dry Color Manufacturer's Association (DCMA) PCB Mixture – 10 components

PCB 1	CAS:2051-60-7	100 µg/ml	PCB 136	CAS:38411-22-2	10 µg/ml
PCB 11	CAS:2050-67-1	100 µg/ml	PCB 185	CAS:52712-05-7	5 µg/ml
PCB 29	CAS:15862-07-4	10 µg/ml	PCB 194	CAS:35694-08-7	5 µg/ml
PCB 47	CAS:2437-79-8	10 µg/ml	PCB 206	CAS:40186-72-9	5 µg/ml
PCB 117	CAS:68194-11-6	10 µg/ml	PCB 209	CAS:2051-24-3	5 µg/ml
Solvent : n-Hexane		ampoule 1 ml		Ref : F128551	Price : JA

Multicomponent Stock Solution of PCBs - 7 components**(ISO 10382:2002 - EN 15308)**

PCB 28	CAS:7012-37-5	PCB 118	CAS:31508-00-6	PCB 180	CAS:35065-29-3
PCB 52	CAS:35693-99-3	PCB 138	CAS:35065-28-2		
PCB 101	CAS:37680-73-2	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 1 ml	100 µg/ml	Ref : F208910	Price : ED
		ampoule 5 ml	100 µg/ml	Ref : F109831	Price : BAH
Solvent : n-Pentane		ampoule 1 ml	100 µg/ml	Ref : F131391	Price : ED

Multicomponent Stock Solution of PCBs & Organochlorine pesticides - 24 components**(ISO10382:2002)**

PCB 28	CAS:7012-37-5	Alpha-HCH	CAS:319-84-6	Heptachlor-exo-epoxide	CAS:1024-57-3
PCB 52	CAS:35693-99-3	Beta-HCH	CAS:319-85-7	Endosulfan-alpha	CAS:959-98-8
PCB 101	CAS:37680-73-2	Gamma-HCH (Lindane)	CAS:58-89-9	4,4'-DDE	CAS:72-55-9
PCB 118	CAS:31508-00-6	Aldrin	CAS:309-00-2	2,4'-DDD	CAS:53-19-0
PCB 138	CAS:35065-28-2	Dieldrin	CAS:60-57-1	2,4'-DDT	CAS:789-02-6
PCB 153	CAS:35065-27-1	Endrin	CAS:72-20-8	4,4'-DDD (TDE)	CAS:72-54-8
PCB 180	CAS:35065-29-3	Heptachlor	CAS:76-44-8	2,4'-DDE	CAS:3424-82-6
Hexachlorobenzene	CAS:118-74-1	Heptachlor-endo-epoxide	CAS:28044-83-9	4,4'-DDT	CAS:50-29-3
Solvent : Iso-Octane		ampoule 1 ml	100 µg/ml	Ref : F109841	Price : KF

Precision & recovery PCBs Standard Solution - 13 components**(ISO 17858:2007)**

PCB 77	CAS:32598-13-3	PCB 126	CAS:57465-28-8	PCB 170	CAS:35065-30-6
PCB 81	CAS:70362-50-4	PCB 156	CAS:38380-08-4	PCB 180	CAS:35065-29-3
PCB 105	CAS:32598-14-4	PCB 157	CAS:69782-90-7	PCB 189	CAS:39635-31-9
PCB 114	CAS:74472-37-0	PCB 167	CAS:52663-72-6		
PCB 123	CAS:65510-44-3	PCB 169	CAS:32774-16-6		
Solvent : n-Nonane		ampoule 1 ml	20 µg/ml	Ref : F062125	Price : GH
		ampoule 5 ml	20 µg/ml	Ref : F110031	Price : BKA

Quality Control Mix – 5 components**(EN12766)**

Arochlor 1242	CAS : 53469-21-9	0.50 µg/l	PCB 30	CAS : 35693-92-6	0.02 µg/l
Arochlor 1254	CAS : 11097-69-1	0.25 µg/l	PCB 209	CAS : 2051-24-3	0.02 µg/l
Arochlor 1260	CAS : 11096-82-5	0.25 µg/l			
Solvent : Iso-Octane		ampoule 1 ml		Ref : F131520	Price : EA

PCBs Standard Solution - 2 components**(EPA 8082)**

PCB 1016	CAS:12674-11-2	PCB 1260	CAS:11096-82-5		
Solvent : Iso-Octane		ampoule 1 ml	1 000 µg/ml	Ref : F117301	Price : DF

PCBs Standard Solution - 8 components**(EPA 525.1)**

PCB 1	CAS:2051-60-7	PCB 47	CAS:2437-79-8	PCB 171	CAS:52663-71-5
PCB 5	CAS:16605-91-7	PCB 98	CAS:60233-25-2	PCB 200	CAS:40186-71-8
PCB 29	CAS:15862-07-4	PCB 154	CAS:60145-22-4		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F114031	Price : KC
		ampoule 1 ml	500 µg/ml	Ref : F114041	Price : DAA

PCBs Standard Solution - 8 components**(EN 12766)**

PCB 1	CAS:2051-60-7	PCB 47	CAS:2437-79-8	PCB 171	CAS:52663-71-5
PCB 5	CAS:16605-91-7	PCB 89	CAS:73575-57-2	PCB 201	CAS:40186-71-8
PCB 29	CAS:15862-07-4	PCB 154	CAS:60145-22-4		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F127921	Price : JJ

PCBs Standard Solution - 19 components

PCB 1	CAS:2051-60-7	PCB 52	CAS:35693-99-3	PCB 180	CAS:35065-29-3
PCB 11	CAS:2050-67-1	PCB 66	CAS:32598-10-0	PCB 187	CAS:52663-68-0
PCB 18	CAS:37680-65-2	PCB 101	CAS:37680-73-2	PCB 194	CAS:35694-08-7
PCB 26	CAS:38444-81-4	PCB 118	CAS:31508-00-6	PCB 206	CAS:40186-72-9
PCB 31	CAS:16606-02-3	PCB 128	CAS:38380-07-3	PCB 209	CAS:2051-24-3
PCB 44	CAS:41464-39-5	PCB 138	CAS:35065-28-2		
PCB 49	CAS:41464-40-8	PCB 170	CAS:35065-30-6		
Solvent : Dichloromethane		ampoule 1 ml	100 µg/ml	Ref : F119851	Price : BAF

PCBs Standard Solution - 19 components**(EPA 8082)**

PCB 1	CAS:2051-60-7	PCB 87	CAS:38380-02-8	PCB 170	CAS:35065-30-6
PCB 5	CAS:16605-91-7	PCB 101	CAS:37680-73-2	PCB 180	CAS:35065-29-3
PCB 18	CAS:37680-65-2	PCB 110	CAS:38380-03-9	PCB 183	CAS:52663-69-1
PCB 31	CAS:16606-02-3	PCB 138	CAS:35065-28-2	PCB 187	CAS:52663-68-0
PCB 44	CAS:41464-39-5	PCB 141	CAS:52712-04-6	PCB 206	CAS:40186-72-9
PCB 52	CAS:35693-99-3	PCB 151	CAS:52663-63-5		
PCB 66	CAS:32598-10-0	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 1 ml	100 µg/ml	Ref : F117291	Price : BDA

PCBs Standard Solution - 32 components

PCB 18	CAS:37680-65-2	PCB 110	CAS:38380-03-9	PCB 156	CAS:38380-08-4
PCB 28	CAS:7012-37-5	PCB 114	CAS:74472-37-0	PCB 157	CAS:69782-90-7
PCB 31	CAS:16606-02-3	PCB 118	CAS:31508-00-6	PCB 167	CAS:52663-72-6
PCB 44	CAS:41464-39-5	PCB 123	CAS:65510-44-3	PCB 169	CAS:32774-16-6
PCB 52	CAS:35693-99-3	PCB 126	CAS:57465-28-8	PCB 170	CAS:35065-30-6
PCB 77	CAS:32598-13-3	PCB 128	CAS:38380-07-3	PCB 177	CAS:52663-70-4
PCB 81	CAS:70362-50-4	PCB 138	CAS:35065-28-2	PCB 180	CAS:35065-29-3
PCB 95	CAS:38379-99-6	PCB 146	CAS:51908-16-8	PCB 183	CAS:52663-69-1
PCB 99	CAS:38380-01-7	PCB 149	CAS:38380-04-0	PCB 187	CAS:52663-68-0
PCB 101	CAS:37680-73-2	PCB 151	CAS:52663-63-5	PCB 189	CAS:39635-31-9
PCB 105	CAS:32598-14-4	PCB 153	CAS:35065-27-1		
Solvent : Iso-Octane		ampoule 1 ml	10 µg/ml	Ref : F127891	Price : CFA

Polycyclic Aromatic Hydrocarbons (PAH)**PAH Calibration Mixture - 16 components**

Acenaphthene	CAS:83-32-9	20 µg/ml	Chrysene	CAS:218-01-9	3.5 µg/ml
Acenaphthylene	CAS:208-96-8	15 µg/ml	Dibenzo(a,h)anthracene	CAS:53-70-3	3.5 µg/ml
Anthracene	CAS:120-12-7	0.8 µg/ml	Fluoranthene	CAS:206-44-0	8 µg/ml
Benzo(a)anthracene	CAS:56-55-3	4 µg/ml	Fluorene	CAS:86-73-7	5 µg/ml
Benzo(a)pyrene	CAS:50-32-8	5 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	4.5 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	4 µg/ml	Naphthalene	CAS:91-20-3	20 µg/ml
Benzo(g,h,i)perylene	CAS:191-24-2	3.5 µg/ml	Phenanthrene	CAS:85-01-8	3.5 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	4.5 µg/ml	Pyrene	CAS:129-00-0	8.5 µg/ml
Solvent : Acetonitrile		ampoule 1 ml		Ref : F128591	Price : FJ

PAH Internal & Surrogate Standard fortification Solution - 3 components**(EPA 525.2)**

Acenaphthene D10	CAS:15067-26-2	Chrysene D12	CAS:1719-03-5	Phenanthrene D10	CAS:1517-22-2
Solvent : Acetone		ampoule 1 ml	500 µg/ml	Ref : F114111	Price : DF

PAH Internal Standard Mixture - 6 components

Acenaphthene D10	CAS:15067-26-2	1,4-Dichlorobenzene D4	CAS:3855-82-1	Perylene D12	CAS:1520-96-3
Chrysene D12	CAS:1719-03-5	Naphthalene D8	CAS:1146-65-2	Phenanthrene D10	CAS:1517-22-2
Solvent : Dichloromethane		ampoule 1 ml	4 000 µg/ml	Ref : F119871	Price : BBC

Multiple Compound Standard Solution - 15 components**(ISO 17993:2002 - ISO 15753:2006)**

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Fluorene	CAS:86-73-7
Anthracene	CAS:120-12-7	Benzo(k)fluoranthene	CAS:207-08-9	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Benzo(a)anthracene	CAS:56-55-3	Chrysene	CAS:218-01-9	Naphthalene	CAS:91-20-3
Benzo(a)pyrene	CAS:50-32-8	Dibenzo(a,h)anthracene	CAS:53-70-3	Phenanthrene	CAS:85-01-8
Benzo(b)fluoranthene	CAS:205-99-2	Fluoranthene	CAS:206-44-0	Pyrene	CAS:129-00-0
Solvent : Acetonitrile		ampoule 1 ml	200 µg/ml	Ref : F110011	Price : HF
		ampoule 5 ml	10 µg/ml	Ref : F110041	Price : GF
Solvent : Acetonitrile/THF/Methanol (1/1/1)		ampoule 10 ml	200 µg/ml	Ref : F146591	Price : CGA
Solvent : Toluene		ampoule 1 ml	100 µg/ml	Ref : F110001	Price : GF

(WHO/EU) PAH in drinking water - 6 components**(ISO 7981-1:2005)**

Fluoranthene	CAS:206-44-0	10 µg/ml	Benzo(g,h,i)perylene	CAS:191-24-2	2 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	2 µg/ml	Benzo(a)pyrene	CAS:50-32-8	2 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	2 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	2 µg/ml
Solvent : Acetonitrile		ampoule 5 ml		Ref : F109250	Price : EF

PAH /ACs Standard Solution - 22 components

Acenaphthene	CAS:83-32-9	4-Bromobiphenyl ether	CAS:101-55-3	Hexachlorobenzene	CAS:118-74-1
Acenaphthylene	CAS:208-96-8	1-Chloronaphthalene	CAS:90-13-1	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Naphthalene	CAS:91-20-3
Benzo(a)anthracene	CAS:56-55-3	Dibenzofuran	CAS:132-64-9	Phenanthrene	CAS:85-01-8
Benzo(a)pyrene	CAS:50-32-8	Dibenzo(a,h)anthracene	CAS:53-70-3	Pyrene	CAS:129-00-0
Benzo(b)fluoranthene	CAS:205-99-2	Dibenzothiophene	CAS:132-65-0	1,2,4-Trichlorobenzene	CAS:120-82-1
Benzo(g,h,i)perylene	CAS:191-24-2	Fluoranthene	CAS:206-44-0		
Benzo(k)fluoranthene	CAS:207-08-9	Fluorene	CAS:86-73-7		

Solvent : Dichloromethane/Toluene (3/1) ampoule 1 ml 100 µg/ml Ref : F119841 Price : GA

PAH QC Mixture - 16 components**(EPA 610, 8100, 8310)**

Acenaphthene	CAS:83-32-9	100 µg/ml	Chrysene	CAS:218-01-9	10 µg/ml
Acenaphthylene	CAS:208-96-8	100 µg/ml	Dibenzo(a,h)anthracene	CAS:53-70-3	10 µg/ml
Anthracene	CAS:120-12-7	100 µg/ml	Fluoranthene	CAS:206-44-0	10 µg/ml
Benzo(a)anthracene	CAS:56-55-3	10 µg/ml	Fluorene	CAS:86-73-7	100 µg/ml
Benzo(a)pyrene	CAS:50-32-8	10 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	10 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	10 µg/ml	Naphthalene	CAS:91-20-3	100 µg/ml
Benzo(g,h,i)perylene	CAS:191-24-2	10 µg/ml	Phenanthrene	CAS:85-01-8	100 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	5 µg/ml	Pyrene	CAS:129-00-0	10 µg/ml

Solvent : Acetonitrile ampoule 1 ml Ref : F117531 Price : GF

PAH Standard Solution - 4 components

Benzo(g,h,i)perylene D12	CAS:93951-66-7	Fluoranthene D10	CAS:93951-69-0	Naphthalene D8	CAS:1146-65-2
Benzo(a)pyrene D12	CAS:63466-71-7				

Solvent : Acetonitrile ampoule 1 ml 100 µg/ml Ref : F112481 Price : KE

PAH Standard Solution - 6 components

Benzo(b)fluoranthene	CAS:205-99-2	Benzo(g,h,i)perylene	CAS:191-24-2	Fluoranthene	CAS:206-44-0
Benzo(k)fluoranthene	CAS:207-08-9	Benzo(a)pyrene	CAS:50-32-8	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5

Solvent : Acetonitrile ampoule 5 ml 10 µg/ml Ref : F109260 Price : EF

ampoule 1 ml 100 µg/ml Ref : F112191 Price : EF

ampoule 1 ml 1 000 µg/ml Ref : F112181 Price : JC

Solvent : Methanol ampoule 1 ml 100 µg/ml Ref : F112641 Price : EF

PAH Standard Solution - 6 components

Benzo(b)fluoranthene	CAS:205-99-2	2 µg/ml	Benzo(a)pyrene	CAS:50-32-8	2 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	2 µg/ml	Fluoranthene	CAS:206-44-0	10 µg/ml
Benzo(g,h,i)perylene	CAS:191-24-2	2 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	2 µg/ml

Solvent : Acetonitrile ampoule 5 ml Ref : F109250 Price : EF

PAH Standard Solution - 6 components**(EPA 525.2)**

Acenaphthene D10	CAS:15067-26-2	1,3-Dimethyl-2-nitrobenzene	CAS:81-20-9	Phenanthrene D10	CAS:1517-22-2
Chrysene D12	CAS:1719-03-5	Perylene D12	CAS:1520-96-3	Triphenylphosphate	CAS:115-86-6

Solvent : Acetone ampoule 1 ml 500 µg/ml Ref : F114121 Price : EA

PAH Standard Solution - 7 components**(EPA 525.2)**

Acenaphthene D10	CAS:15067-26-2	Perylene D12	CAS:1520-96-3	Pyrene D10	CAS:1718-52-1
Chrysene D12	CAS:1719-03-5	Phenanthrene D10	CAS:1517-22-2	Triphenylphosphate	CAS:115-86-6
1,3-Dimethyl-2-nitrobenzene	CAS:81-20-9				

Solvent : Acetone ampoule 1 ml 500 µg/ml Ref : F114131 Price : FF

PAH Standard Solution - 7 components

Benzo(a)anthracene	CAS:56-55-3	Benzo(k)fluoranthene	CAS:207-08-9	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Benzo(b)fluoranthene	CAS:205-99-2	Benzo(a)pyrene	CAS:50-32-8		
Benzo(j)fluoranthene	CAS:205-82-3	Dibenzo(a,h)anthracene	CAS:53-70-3		

Solvent : Toluene ampoule 1 ml 10 µg/ml Ref : F112471 Price : EF

PAH Standard Solution - 8 components**(EPA 8100)**

Benzo(j)fluoranthene	CAS:205-82-3	7-H-Dibenzo(c,g)carbazole	CAS:194-59-2	Dibenzo(a,i)pyrene	CAS:189-55-9
Dibenz(a,h)acridine	CAS:226-36-8	Dibenzo(a,e)pyrene	CAS:192-65-4	3-Methylcholanthrene	CAS:56-49-5
Dibenz(a,j)acridine	CAS:224-42-0	Dibenzo(a,h)pyrene	CAS:189-64-0		

Solvent : Dichloromethane ampoule 1 ml 1 000 µg/ml Ref : F117521 Price : CKA

PAH Standard Solution - 13 components

(EPA 525.1)

Acenaphthylene	CAS:208-96-8	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Anthracene	CAS:120-12-7	Benzo(a)pyrene	CAS:50-32-8	Phenanthrene	CAS:85-01-8
Benzo(a)anthracene	CAS:56-55-3	Chrysene	CAS:218-01-9	Pyrene	CAS:129-00-0
Benzo(j)fluoranthene	CAS:205-82-3	Dibenzo(a,h)anthracene	CAS:53-70-3		
Benzo(k)fluoranthene	CAS:207-08-9	Fluorene	CAS:86-73-7		
Solvent : Acetone		ampoule 1 ml	100 µg/ml	Ref : F113971	Price : EF
Solvent : Toluene		ampoule 1 ml	500 µg/ml	Ref : F113981	Price : KF

PAH Standard Solution - 16 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Acenaphthylene	CAS:208-96-8	Benzo(a)pyrene	CAS:50-32-8	Naphthalene	CAS:91-20-3
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Phenanthrene	CAS:85-01-8
Benzo(a)anthracene	CAS:56-55-3	Dibenzo(a,h)anthracene	CAS:53-70-3	Pyrene	CAS:129-00-0
Benzo(b)fluoranthene	CAS:205-99-2	Fluoranthene	CAS:206-44-0		
Benzo(k)fluoranthene	CAS:207-08-9	Fluorene	CAS:86-73-7		
Solvent : Acetonitrile	(ISO13877:1998)	ampoule 5 ml	10 µg/ml	Ref : F109871	Price : GF
	(ISO 13877:1998, EPA 610, EPA8310)	ampoule 1 ml	100 µg/ml	Ref : F109861	Price : GF
	(EPA 550 & 550.1)	ampoule 1 ml	1 000 µg/ml	Ref : F114231	Price : BDF
Solvent : Cyclohexane		ampoule 1 ml	200 µg/ml	Ref : F131381	Price : JA
Solvent : Dichloromethane	(EPA 8100)	ampoule 1 ml	1 000 µg/ml	Ref : F117511	Price : BDF
Solvent : Dichloromethane/Benzene (1/1)		ampoule 1 ml	2 000 µg/ml	Ref : F128611	Price : BJF
Solvent : Dichloromethane /Methanol(1/1)	(EPA 610)	ampoule 1 ml	20 µg/ml	Ref : F115521	Price : GF
		ampoule 1 ml	100 µg/ml	Ref : F115511	Price : GF

PAH Standard Solution - 16 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Acenaphthylene	CAS:208-96-8	Benzo(a)pyrene	CAS:50-32-8	Naphthalene	CAS:91-20-3
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	Phenanthrene	CAS:85-01-8
Benzo(a)anthracene	CAS:56-55-3	Dibenzo(a,c)anthracene	CAS:215-58-7	Pyrene	CAS:129-00-0
Benzo(b)fluoranthene	CAS:205-99-2	Fluoranthene	CAS:206-44-0		
Benzo(k)fluoranthene	CAS:207-08-9	Fluorene	CAS:86-73-7		
Solvent : Acetonitrile		ampoule 5 ml	10 µg/ml	Ref : F112461	Price : GA
Solvent : Cyclohexane		ampoule 5 ml	10 µg/ml	Ref : F112451	Price : GA
Solvent : Toluene	(EPA 505)	ampoule 1 ml	100 µg/ml	Ref : F062285	Price : DF
		ampoule 5 ml	100 µg/ml	Ref : F112441	Price : BCA

PAH Standard Solution - 16 components

Acenaphthene D10	CAS:15067-26-2	Benzo(g,h,i)perylene D12	CAS:93951-66-7	Indeno(1,2,3-c,d)pyrene D12	CAS:203578-33-0
Acenaphthylene D8	CAS:93951-97-4	Benzo(a)pyrene D12	CAS:63466-71-7	Naphthalene D8	CAS:1146-65-2
Anthracene D10	CAS:1719-06-8	Chrysene D12	CAS:1719-03-5	Phenanthrene D10	CAS:1517-22-2
Benzo(a)anthracene D12	CAS:1718-53-2	Dibenz(a,h)anthracene D14	CAS:13250-98-1	Pyrene D10	CAS:1718-52-1
Benzo(b)Fluoranthene D12	CAS:93951-98-5	Fluoranthene D10	CAS:93951-69-0		
Benzo(k)Fluoranthene D12	CAS:93952-01-3	Fluorene D10	CAS:81103-79-9		
Solvent : Cyclohexane		ampoule 1 ml	100 µg/ml	Ref : F112491	Price : BKA

PAH Standard Solution - 16 components

(EPA 550 & 550.1)

Acenaphthene	CAS:83-32-9	1 000 µg/ml	Chrysene	CAS:218-01-9	625 µg/ml
Acenaphthylene	CAS:208-96-8	1 000 µg/ml	Dibenzo(a,h)anthracene	CAS:53-70-3	1.25 µg/ml
Anthracene	CAS:120-12-7	625 µg/ml	Fluoranthene	CAS:206-44-0	2.50 µg/ml
Benzo(a)anthracene	CAS:56-55-3	1 µg/ml	Fluorene	CAS:86-73-7	100 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	1 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	12.50 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	1.25 µg/ml	Naphthalene	CAS:91-20-3	1 000 µg/ml
Benzo(g,h,i)perylene	CAS:191-24-2	5 µg/ml	Phenanthrene	CAS:85-01-8	50 µg/ml
Benzo(a)pyrene	CAS:50-32-8	5 µg/ml	Pyrene	CAS:129-00-0	62.5 µg/ml
Solvent : Acetonitrile		ampoule 1 ml		Ref : F114221	Price : GF

PAH Standard Solution - 16 components

Acenaphthene	CAS:83-32-9	100 µg/ml	Chrysene	CAS:218-01-9	10 µg/ml
Acenaphthylene	CAS:208-96-8	100 µg/ml	Dibenzo(a,h)anthracene	CAS:53-70-3	10 µg/ml
Anthracene	CAS:120-12-7	100 µg/ml	Fluoranthene	CAS:206-44-0	10 µg/ml
Benzo(a)anthracene	CAS:56-55-3	10 µg/ml	Fluorene	CAS:86-73-7	100 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	10 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	10 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	5 µg/ml	Naphthalene	CAS:91-20-3	100 µg/ml
Benzo(g,h,i)perylene	CAS:191-24-2	10 µg/ml	Phenanthrene	CAS:85-01-8	100 µg/ml
Benzo(a)pyrene	CAS:50-32-8	10 µg/ml	Pyrene	CAS:129-00-0	10 µg/ml
Solvent : Acetonitrile		ampoule 1 ml		Ref : F128601	Price : FK

PAH Standard Solution - 16 components**(EPA 610)**

Acenaphthene	CAS:83-32-9	100 µg/ml	Chrysene	CAS:218-01-9	100 µg/ml
Acenaphthylene	CAS:208-96-8	200 µg/ml	Dibenzo(a,h)anthracene	CAS:53-70-3	200 µg/ml
Anthracene	CAS:120-12-7	100 µg/ml	Fluoranthene	CAS:206-44-0	200 µg/ml
Benzo(a)anthracene	CAS:56-55-3	100 µg/ml	Fluorene	CAS:86-73-7	200 µg/ml
Benzo(b)fluoranthene	CAS:205-99-2	200 µg/ml	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	100 µg/ml
Benzo(k)fluoranthene	CAS:207-08-9	100 µg/ml	Naphthalene	CAS:91-20-3	1 000 µg/ml
Benzo(g,h,i)perylene	CAS:191-24-2	200 µg/ml	Phenanthrene	CAS:85-01-8	100 µg/ml
Benzo(a)pyrene	CAS:50-32-8	100 µg/ml	Pyrene	CAS:129-00-0	100 µg/ml

Solvent : Methanol / Dichloromethane (1/1)

ampoule 1 ml

Ref: F115531

Price : HF

PAH Standard Solution - 18 components

Acenaphthene	CAS:83-32-9	Benzo(g,h,i)perylene	CAS:191-24-2	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5
Acenaphthylene	CAS:208-96-8	Benzo(a)pyrene	CAS:50-32-8	2-Methyl-Fluoranthene	CAS:33543-31-6
Anthracene	CAS:120-12-7	Chrysene	CAS:218-01-9	2-Methylnaphthalene	CAS:91-57-6
Benzo(a)anthracene	CAS:56-55-3	Dibenzo(a,h)anthracene	CAS:53-70-3	Naphthalene	CAS:91-20-3
Benzo(b)fluoranthene	CAS:205-99-2	Fluoranthene	CAS:206-44-0	Phenanthrene	CAS:85-01-8
Benzo(k)fluoranthene	CAS:207-08-9	Fluorene	CAS:86-73-7	Pyrene	CAS:129-00-0

Solvent : Methanol

ampoule 1 ml

100 µg/ml

Ref: F112201

Price : HA

PAH Standard Solution - 33 components

Acenaphthene	CAS:83-32-9	Benzo(e)pyrene	CAS:192-97-2	2-Methylnaphthalene	CAS:91-57-6
Acenaphthylene	CAS:208-96-8	Chrysene	CAS:218-01-9	1-Methylphenanthrene	CAS:832-69-9
Anthracene	CAS:120-12-7	Dibenzo(a,h)anthracene	CAS:53-70-3	2-Methylphenanthrene	CAS:2531-84-2
Benzo(a)anthracene	CAS:56-55-3	Dibenzofuran	CAS:132-64-9	1-Methylpyrene	CAS:2381-21-7
Benzo(b)fluoranthene	CAS:205-99-2	Dibenzothiophene	CAS:132-65-0	Naphthalene	CAS:91-20-3
Benzo(j)fluoranthene	CAS:205-82-3	1,6-Dimethylnaphthalene	CAS:575-43-9	Perylene	CAS:198-55-0
Benzo(k)fluoranthene	CAS:207-08-9	Fluoranthene	CAS:206-44-0	Phenanthrene	CAS:85-01-8
Benzo(b)naphtho(1,2-d)thiophene	id :479	3,6-Dimethylphenanthrene	CAS:1576-67-6	Picene	CAS:213-46-7
Benzo(b)naphtho(2,1-d)thiophene	CAS:239-35-0	Indeno(1,2,3-c,d)pyrene	CAS:193-39-5	Pyrene	CAS:129-00-0
Benzo(g,h,i)perylene	CAS:191-24-2	1-Methylchrysene	CAS:3351-28-8	2,3,5-Trimethylnaphthalene	CAS:2245-38-7
Benzo(a)pyrene	CAS:50-32-8	1-Methylnaphthalene	CAS:90-12-0	Triphenylene	CAS:217-59-4

Solvent : Iso-Octane

ampoule 1 ml

10 µg/ml

Ref: F112501

Price : BHF

Residual Solvents in Flexible packaging material**VOC Standard Solution - 10 components****(EN 13628)**

Acetic acid-isobutyl ester	CAS:110-19-0	Ethyl Acetate	CAS:141-78-6	Trichloroethene	CAS:79-01-6
Acetone	CAS:67-64-1	Methanol	CAS:67-56-1	Xylenes	CAS:1330-20-7
1-Butanol	CAS :71-36-3	Methylisobutylketone	CAS:108-10-1		
Ethanol	CAS:64-17-5	Toluene	CAS:108-88-3		

Solvent : Iso-Octane

ampoule 1 ml

200 µg/ml

Ref: F129221

Price : CJ

Residual Solvents in Medical Products

The analysis of residual solvents in pharmaceutical products has changed, particularly for products being sold into Europe. The International Conference on Harmonization (ICH) Guidelines for Residual Solvents is becoming the international standard and is being adopted by more pharmacopoeias, including the European and United States Pharmacopoeia, every year. The ICH method and compound list is more extensive than any method previously used and poses new challenges. Compounds in Class 1 are solvents considered to be of highest risk and to be avoided in pharmaceutical manufacturing. Use of Class 2 compounds is to be limited, as they pose a lower, but present, threat to health. Compounds in Class 3 pose the lowest toxic potential and may be used routinely in manufacturing.

Class 1 solvents**Mix recommended by the European Pharmacopoeia/ICH - 5 components**

Benzene	CAS:71-43-2	2 µg/ml	Tetrachloromethane	CAS:56-23-5	4 µg/ml
1,2-Dichloroethane	CAS:107-06-2	5 µg/ml	1,1,1-Trichloroethane	CAS:71-55-6	10 µg/ml
1,1-Dichloroethene	CAS:75-35-4	8 µg/ml			

Solvent : Dimethylsulphoxide/Water (9/1)

ampoule 1 ml

Ref: F131530

Price : DF

Class 2 solvents**European Pharmacopoeia/ICH Class 2 Mix A - 13 components**

Chlorobenzene	CAS:108-90-7	360 µg/ml	Methylcyclohexane	CAS:108-87-2	1 180 µg/ml
Cyclohexane	CAS:110-82-7	3 880 µg/ml	Toluene	CAS:108-88-3	890 µg/ml
cis-1,2-Dichloroethene	CAS:156-59-2	1 870 µg/ml	1,1,2-Trichloroethene	CAS:79-01-6	80 µg/ml
Dichloromethane	CAS:75-09-2	600 µg/ml	m-Xylene	CAS:108-38-3	1 302 µg/ml
N,N - Dimethylformamide	CAS:68-12-2	880 µg/ml	o-Xylene	CAS:95-47-6	195 µg/ml
Ethylbenzene	CAS:100-41-4	369 µg/ml	p-Xylene	CAS:106-42-3	304 µg/ml
n-Hexane (C6)	CAS:110-54-3	290 µg/ml			

Solvent : Dimethylsulphoxide/Water (9/1)

ampoule 1 ml

Ref: F131550

Price : EJ

European Pharmacopeia/ICH Class 2 Mix B - 10 components

Acetic acid-dimethylamide	CAS:127-19-5	1 090 µg/ml	2-Hexanone	CAS:591-78-6	50 µg/ml
Acetonitrile	CAS:75-05-8	410 µg/ml	Methanol	CAS:67-56-1	3 000 µg/ml
Chloroform	CAS:67-66-3	60 µg/ml	Nitromethane	CAS:75-52-5	50 µg/ml
1,2-Dimethoxyethane	CAS:110-71-4	100 µg/ml	Pyridine	CAS:110-86-1	200 µg/ml
Dioxan	CAS:123-91-1	380 µg/ml	1,2,3,4-Tetrahydronaphthalene	CAS:119-64-2	100 µg/ml
Solvent : Dimethylsulphoxide/Water (9/1)		ampoule 1 ml	Ref : F131670	Price : EF	

European Pharmacopeia/ICH Class 2 Mix C - 6 components

2-Ethoxyethanol	CAS:110-80-5	160 µg/ml	2-Methoxyethanol	CAS:109-86-4	50 µg/ml
Ethylene glycol	CAS:107211	620 µg/ml	N-Methyl Pyrrolidone	CAS:872-50-4	4840 µg/ml
Formamide	CAS:75-12-7	220 µg/ml	Sulfolan	CAS:126-33-0	160 µg/ml
Solvent : Water		ampoule 1 ml	Ref : F131540	Price : DF	

Revised USP 467 Residual Solvent Method

On July 1, 2008, the USP released a new major revision of its general chapter <467>: now entitled *Residual Solvent Analysis*. Closely resembling the International Conference on Harmonization's Q3C Guideline, the chapter applies to all compendial drug substances and products and provides concentration limits of Class 1, 2, and 3 residual solvents.

Residual Solvents – Class 1 – 5 components

Benzene	CAS:71-43-2	10 µg/ml	1,1-Dichloroethene	CAS:75-35-4	40 µg/ml
Carbon tetrachloride	CAS:56-23-5	20 µg/ml	Trichloroethene	CAS:79-01-6	50 µg/ml
1,2-Dichloroethane	CAS:107-06-2	25 µg/ml			
Solvent : Dimethylsulphoxide		ampoule 1 ml	Ref : F131680	Price : DF	

Residual Solvents – Class 2 – Mix A – 15 components

Acetonitrile	CAS:75-05-8	2.05 µg/ml	Methylcyclohexane	CAS:108-87-2	5.90 µg/ml
Chlorobenzene	CAS:108-90-7	1.8 µg/ml	Methylene chloride	CAS:75-09-2	3.00 µg/ml
Cyclohexane	CAS:110-82-7	19.4 µg/ml	Tetrahydrofuran	CAS:109-99-9	3.45 µg/ml
cis-1,2-Dichloroethene	CAS:156-59-2	4.7 µg/ml	Toluene	CAS:108-88-3	4.45 µg/ml
trans-1,2-Dichloroethene	CAS:156-60-5	4.7 µg/ml	m-Xylene	CAS:108-38-3	6.51 µg/ml
1,4-Dioxane	CAS:123-91-1	1.9 µg/ml	o-Xylene	CAS:95-47-6	0.98 µg/ml
Ethylbenzene	CAS:100-41-4	1.84 µg/ml	p-Xylene	CAS:106-42-3	1.52 µg/ml
Methanol	CAS:67-56-1	15 µg/ml			
Solvent : Dimethylsulphoxide		ampoule 1 ml	Ref : F131690	Price : EJ	

Residual Solvents – Class 2 – Mix B - 8 components

Chloroform	CAS:67-66-3	60 µg/ml	Nitromethane	CAS:75-52-5	50 µg/ml
1,2-Dimethoxyethane	CAS:110-71-4	100 µg/ml	Pyridine	CAS:110-86-1	200 µg/ml
n-Hexane (C6)	CAS:110-54-3	290 µg/ml	Tetralin	CAS:119-64-2	100 µg/ml
2-Hexanone	CAS:591-78-6	50 µg/ml	Trichloroethene	CAS:79-01-6	80 µg/ml
Solvent : Dimethylsulphoxide		ampoule 1 ml	Ref : F131700	Price : EF	

VOC Standard Solution - 1 component

Ethylene Oxide	CAS:75-21-8	1 000 µg/ml			
Solvent : Dimethylsulphoxide		ampoule 1 ml	1 000 µg/ml	Ref : F131770	Price : CF
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F131780	Price : CF

VOC Standard Solution - 4 components

Chloroform	CAS:67-66-3	60 µg/ml	Dioxan	CAS:123-91-1	380 µg/ml
Dichloromethane	CAS:75-09-2	600 µg/ml	Trichloroethene	CAS:79-01-6	80 µg/ml
Solvent : Dimethylsulphoxide		ampoule 1 ml	Ref : F131720	Price : DK	
Solvent : Methanol purge and trap		ampoule 1 ml	Ref : F131740	Price : DK	

VOC Standard Solution - 5 components

Benzene	CAS:71-43-2	2 µg/ml	Dioxan	CAS:123-91-1	380 µg/ml
Chloroform	CAS:67-66-3	60 µg/ml	Trichloroethene	CAS:79-01-6	80 µg/ml
Dichloromethane	CAS:75-09-2	600 µg/ml			
Solvent : Dimethylsulphoxide		ampoule 1 ml	Ref : F131710	Price : DF	
Solvent : Methanol purge & trap		ampoule 1 ml	Ref : F131730	Price : DF	

VOC Standard Solution - 5 components

Benzene	CAS:71-43-2	200 µg/ml	Dioxan	CAS:123-91-1	200 µg/ml
Chloroform	CAS:67-66-3	100 µg/ml	Trichloroethene	CAS:79-01-6	200 µg/ml
Dichloromethane	CAS:75-09-2	1 000 µg/ml			
Solvent : Dimethylsulphoxide		ampoule 1 ml	Ref : F131750	Price : DF	
Solvent : Methanol purge & trap		ampoule 1 ml	Ref : F131760	Price : DF	

Sulfur in Petroleum Products

A gravimetrically prepared calibration standard and NIST Traceable by weight, and have been verified by analytical methods.

Total Sulfur

Calibration Kit for Total sulfur – 0.125-2.55 w/wppm

Calibration range : Total Sulfur by weight from Di-N-Butylsulfide : 0.125-2.55 ppm

Ampoule A	Blank		Ampoule D	Total Sulfur from Di-N-Butylsulfide	0.50 w/w ppm
Ampoule B	Total Sulfur from Di-N-Butylsulfide	0.125 w/w ppm	Ampoule E	Total Sulfur from Di-N-Butylsulfide	1.00 w/w ppm
Ampoule C	Total Sulfur from Di-N-Butylsulfide	0.250 w/w ppm	Ampoule F	Total Sulfur from Di-N-Butylsulfide	2.50 w/w ppm
Solvent : Iso-octane		6 ampoules X 1 ml	Ref : F205024		Price : BCA

Check Standard for Total sulfur – 0.75 w/w ppm

Total Sulfur from Di-N-Butylsulfide	0.75 w/w ppm		
Solvent : Iso-octane	ampoule 1 ml	Ref : F205084	Price : CF

Calibration Kit for Total sulfur – 2.50 - 50 w/w ppm

Calibration range : Total Sulfur by weight from Di-N-Butylsulfide : 2.50 - 50 ppm

Ampoule A	Blank		Ampoule E	Total Sulfur from Di-N-Butylsulfide	15.00 w/w ppm
Ampoule B	Total Sulfur from Di-N-Butylsulfide	2.50 w/w ppm	Ampoule F	Total Sulfur from Di-N-Butylsulfide	20.00 w/w ppm
Ampoule C	Total Sulfur from Di-N-Butylsulfide	5.00 w/w ppm	Ampoule G	Total Sulfur from Di-N-Butylsulfide	25.00 w/w ppm
Ampoule D	Total Sulfur from Di-N-Butylsulfide	10.00 w/w ppm	Ampoule H	Total Sulfur from Di-N-Butylsulfide	50.00 w/w ppm
Solvent : Iso-octane		8 ampoules X 1 ml	Ref : F205094		Price : BGA

Check Standard for Total sulfur – 30 w/w ppm

Total Sulfur from Di-N-Butylsulfide	30 w/w ppm		
Solvent : Iso-octane	ampoule 1 ml	Ref : F205174	Price : CF

Calibration Kit for Total sulfur – 50 – 1 000 w/w ppm

Calibration range : Total Sulfur by weight from Di-N-Butylsulfide : 50 – 1 000 ppm

Ampoule A	Blank		Ampoule E	Total Sulfur from Di-N-Butylsulfide	250 w/w ppm
Ampoule B	Total Sulfur from Di-N-Butylsulfide	50 w/w ppm	Ampoule F	Total Sulfur from Di-N-Butylsulfide	500 w/w ppm
Ampoule C	Total Sulfur from Di-N-Butylsulfide	75 w/w ppm	Ampoule G	Total Sulfur from Di-N-Butylsulfide	1 000 w/w ppm
Ampoule D	Total Sulfur from Di-N-Butylsulfide	100 w/w ppm			
Solvent : Iso-octane		7 ampoules X 1 ml	Ref : F205184		Price : BEA

Check Standard for Total sulfur – 300 w/w ppm

Total Sulfur from Di-N-Butylsulfide	300 w/w ppm		
Solvent : Iso-octane	ampoule 1 ml	Ref : F205254	Price : CF

Calibration Kit for Total sulfur – 1 000 -6 000 w/wppm

Calibration range : Total Sulfur by weight from Di-N-Butylsulfide : 1 000- 6 000 ppm

Ampoule A	Blank		Ampoule D	Total Sulfur from Di-N-Butylsulfide	2 000 w/w ppm
Ampoule B	Total Sulfur from Di-N-Butylsulfide	1 000 w/w ppm	Ampoule E	Total Sulfur from Di-N-Butylsulfide	4 000 w/w ppm
Ampoule C	Total Sulfur from Di-N-Butylsulfide	1 500 w/w ppm	Ampoule F	Total Sulfur from Di-N-Butylsulfide	6 000 w/w ppm
Solvent : Iso-octane		6 ampoules X 1 ml	Ref : F205264		Price : BCA

Check Standard for Total sulfur – 3 000 w/w ppm

Total Sulfur from Di-N-Butylsulfide	3 000 w/w ppm		
Solvent : Iso-octane	ampoule 1 ml	Ref : F205324	Price : CF

Total Sulfur & Total Nitrogen

Calibration Kit for Total sulfur & Total Nitrogen– 0.125 -5.00 w/wppm

Calibration range : Total Sulfur by weight from Di-N-Butylsulfide : 1 000- 6 000 ppm

Ampoule A	Blank				
Ampoule B	Total Sulfur from Thiophene	0.125 w/w ppm	&	Total Nitrogen from Pyridine	0.125 w/w ppm
Ampoule C	Total Sulfur from Thiophene	0.25 w/w ppm	&	Total Nitrogen from Pyridine	0.25 w/w ppm
Ampoule D	Total Sulfur from Thiophene	0.50 w/w ppm	&	Total Nitrogen from Pyridine	0.50 w/w ppm
Ampoule E	Total Sulfur from Thiophene	1.00 w/w ppm	&	Total Nitrogen from Pyridine	1.00 w/w ppm
Ampoule F	Total Sulfur from Thiophene	2.50 w/w ppm	&	Total Nitrogen from Pyridine	2.50 w/w ppm
Ampoule G	Total Sulfur from Thiophene	5.00 w/w ppm	&	Total Nitrogen from Pyridine	5.00 w/w ppm
Solvent : Iso-octane		7 ampoules X 1 ml	Ref : F205334		Price : BHF

Calibration Kit for Total sulfur & Total Nitrogen – 5.00 – 50 w/wppm

Calibration range : Total Sulfur by weight from Di-N-Butylsulfide : 1 000- 6 000 ppm

Ampoule A Blank

Ampoule B	Total Sulfur from Thiophene	5 w/w ppm	&	Total Nitrogen from Pyridine	5 w/w ppm
Ampoule C	Total Sulfur from Thiophene	10 w/w ppm	&	Total Nitrogen from Pyridine	10 w/w ppm
Ampoule D	Total Sulfur from Thiophene	25 w/w ppm	&	Total Nitrogen from Pyridine	25 w/w ppm
Ampoule E	Total Sulfur from Thiophene	50 w/w ppm	&	Total Nitrogen from Pyridine	50 w/w ppm
Solvent : Iso-octane		5 ampoules X 1 ml		Ref : F205404	Price : BCF

Calibration Kit for Total sulfur & Total Nitrogen – 50 -1 000 w/wppm

Calibration range : Total Sulfur by weight from Di-N-Butylsulfide : 1 000- 6 000 ppm

Ampoule A Blank

Ampoule B	Total Sulfur from Thiophene	50 w/w ppm	&	Total Nitrogen from Pyridine	50 w/w ppm
Ampoule C	Total Sulfur from Thiophene	75 w/w ppm	&	Total Nitrogen from Pyridine	75 w/w ppm
Ampoule D	Total Sulfur from Thiophene	100 w/w ppm	&	Total Nitrogen from Pyridine	100 w/w ppm
Ampoule E	Total Sulfur from Thiophene	250 w/w ppm	&	Total Nitrogen from Pyridine	250 w/w ppm
Ampoule F	Total Sulfur from Thiophene	500 w/w ppm	&	Total Nitrogen from Pyridine	500 w/w ppm
Ampoule G	Total Sulfur from Thiophene	1 000 w/w ppm	&	Total Nitrogen from Pyridine	1 000 w/w ppm
Solvent : Iso-octane		7 ampoules X 1 ml		Ref : F205454	Price : BHF

Volatile & Semi-Volatile Organic Compounds**VOC / ACs Standard Solution - 8 components**

Benzene	CAS:71-43-2	Toluene	CAS:108-88-3	o-Xylene	CAS:95-47-6
Ethylbenzene	CAS:100-41-4	1,2,4-Trimethylbenzene	CAS:95-63-6	m-Xylene	CAS:108-38-3
Methyl-tert.butylether	CAS:1634-04-4	1,3,5-Trimethylbenzene	CAS:108-67-8		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112341	Price : EA

VOC / ACs Standard Solution - 9 components

Benzene	CAS:71-43-2	500 µg/ml	Toluene	CAS:108-88-3	1 500 µg/ml
Ethylbenzene	CAS:100-41-4	1 500 µg/ml	1,2,4-Trimethylbenzene	CAS:95-63-6	1 000 µg/ml
n-Heptane	CAS:142-82-5	500 µg/ml	o-Xylene	CAS:95-47-6	1 000 µg/ml
Iso-octane (2,2,4-trimethylpentane)	CAS:540-84-1	1 500 µg/ml	m-Xylene	CAS:108-38-3	1 000 µg/ml
2-Methylpentane	CAS:107-83-5	1 500 µg/ml			
Solvent : Methanol purge & trap		ampoule 1 ml		Ref : F112351	Price : EF

VOC / ACs Standard Solution - 10 components

Benzene	CAS:71-43-2	Toluene	CAS:108-88-3	m-Xylene	CAS:108-38-3
Ethylbenzene	CAS:100-41-4	1,2,4-Trimethylbenzene	CAS:95-63-6	p-Xylene	CAS:106-42-3
Methyl-tert.butylether	CAS:1634-04-4	1,3,5-Trimethylbenzene	CAS:108-67-8		
Naphthalene	CAS:91-20-3	o-Xylene	CAS:95-47-6		
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F112361	Price : EG

VOC / ACs Standard Solution - 16 components

Benzene	CAS:71-43-2	2.5 µg/ml	Tetrachloromethane	CAS:56-23-5	5 µg/ml
Bromodichloromethane	CAS:75-27-4	5 µg/ml	Toluene	CAS:108-88-3	2.5 µg/ml
Chloroform	CAS:67-66-3	5 µg/ml	Tribromomethane	CAS:75-25-2	5 µg/ml
Dibromochloromethane	CAS:124-48-1	5 µg/ml	1,1,1-Trichloroethane	CAS:71-55-6	5 µg/ml
1,2-Dichloroethane	CAS:107-06-2	15 µg/ml	Trichloroethene	CAS:79-01-6	5 µg/ml
Dichloromethane	CAS:75-09-2	100 µg/ml	m-Xylene	CAS:108-38-3	2.5 µg/ml
Ethylbenzene	CAS:100-41-4	2.5 µg/ml	o-Xylene	CAS:95-47-6	2.5 µg/ml
Tetrachloroethene	CAS:127-18-4	5 µg/ml	p-Xylene	CAS:106-42-3	2.5 µg/ml
Solvent : Methanol purge & trap		ampoule 1 ml		Ref : F127721	Price : FA

VOC / ACs Standard Solution - 19 components

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	Toluene	CAS:108-88-3
Bromodichloromethane	CAS:75-27-4	1,2-Dichloroethane	CAS:107-06-2	Vinylchloride	CAS:75-01-4
Chlorobenzene	CAS:108-90-7	Styrene	CAS:100-42-5	o-Xylene	CAS:95-47-6
Chloroform	CAS:67-66-3	Tetrachloromethane	CAS:56-23-5	m-Xylene	CAS:108-38-3
Dibromochloromethane	CAS:124-48-1	Tetrachloroethene	CAS:127-18-4	p-Xylene	CAS:106-42-3
1,2-Dichlorobenzene	CAS:95-50-1	Tribromomethane	CAS:75-25-2		
1,3-Dichlorobenzene	CAS:541-73-1	Trichloroethene	CAS:79-01-6		
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F127691	Price : GA

VOC / ACS Standard Solution - 22 components**(ISO 10301:1997)**

Bromochloromethane	CAS:74-97-5	1,1-Dichloroethene	CAS:75-35-4	Trichloroethene	CAS:79-01-6
Bromodichloromethane	CAS:75-27-4	cis-1,2-Dichloroethene	CAS:156-59-2	Tetrachloroethene	CAS:127-18-4
Chloroform	CAS:67-66-3	trans-1,2-Dichloroethene	CAS:156-60-5	Tetrachloromethane	CAS:56-23-5
Dibromochloromethane	CAS:124-48-1	Dichloromethane	CAS:75-09-2	Tribromomethane	CAS:75-25-2
1,2-Dibromoethane	CAS:106-93-4	1,2-Dichloropropane	CAS:78-87-5	1,1,1-Trichloroethane	CAS:71-55-6
1,1-Dichloroethane	CAS:75-34-3	1,3-Dichloropropane	CAS:142-28-9	1,1,2-Trichloroethane	CAS:79-00-5
1,2-Dichloroethane	CAS:107-06-2	cis-1,3-Dichloropropene	CAS:10061-01-5		
Dibromomethane	CAS:74-95-3	trans-1,3-Dichloropropene	CAS:10061-02-6		
Solvent : Iso-octane		ampoule 10 ml	10 µg/ml	Ref : F128961	Price : GF

VOC / ACS Standard Solution - 28 components

Benzene	CAS:71-43-2	1,2-Dichloroethane	CAS:107-06-2	Toluene	CAS:108-88-3
Bromodichloromethane	CAS:75-27-4	1,1-Dichloroethene	CAS:75-35-4	Tribromomethane	CAS:75-25-2
sec-Butylbenzene	CAS:135-98-8	cis-1,2-Dichloroethene	CAS:156-59-2	1,1,1-Trichloroethane	CAS:71-55-6
Chlorobenzene	CAS:108-90-7	trans-1,2-Dichloroethene	CAS:156-60-5	1,1,2-Trichloroethane	CAS:79-00-5
Chloroform	CAS:67-66-3	Dichloromethane	CAS:75-09-2	Trichloroethene	CAS:79-01-6
Dibromochloromethane	CAS:124-48-1	Ethylbenzene	CAS:100-41-4	o-Xylene	CAS:95-47-6
1,2-Dichlorobenzene	CAS:95-50-1	Fluorotrichloromethane	CAS:75-69-4	m-Xylene	CAS:108-38-3
1,3-Dichlorobenzene	CAS:541-73-1	Tetrachloromethane	CAS:56-23-5	p-Xylene	CAS:106-42-3
1,4-Dichlorobenzene	CAS:106-46-7	1,1,2,2-Tetrachloroethane	CAS:79-34-5		
1,1-Dichloroethane	CAS:75-34-3	Tetrachloroethene	CAS:127-18-4		
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F112281	Price : HC

VOC Purgeable Gas Solution - 5 components**(EPA 601)**

Bromomethane	CAS:74-83-9	Chloromethane	CAS:74-87-3	Vinylchloride	CAS:75-01-4
Chloroethane	CAS:75-00-3	Dichlorodifluoromethane	CAS:75-71-8		
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F114851	Price : EA

Haloalkanes Mixture - 35 components**(EPA 502, EPA 524.2)**

Bromochloromethane	CAS:74-97-5	Dichlorodifluoromethane	CAS:75-71-8	Hexachlorobutadiene	CAS:87-68-3
Bromodichloromethane	CAS:75-27-4	1,1-Dichloroethane	CAS:75-34-3	Methylene chloride	CAS:75-09-2
Bromoform	CAS:75-25-2	1,2-Dichloroethane	CAS:107-06-2	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromomethane	CAS:74-83-9	1,1-Dichloroethene	CAS:75-35-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Carbon tetrachloride	CAS:56-23-5	cis-1,2-Dichloroethene	CAS:156-59-2	Tetrachloroethene	CAS:127-18-4
Chloroethane	CAS:75-00-3	trans-1,2-Dichloroethene	CAS:156-60-5	1,1,1-Trichloroethane	CAS:71-55-6
Chloroform	CAS:67-66-3	1,2-Dichloropropane	CAS:78-87-5	1,1,2-Trichloroethane	CAS:79-00-5
Chloromethane	CAS:74-87-3	1,3-Dichloropropane	CAS:142-28-9	Trichloroethene	CAS:79-01-6
Dibromochloromethane	CAS:124-48-1	2,2-Dichloropropane	CAS:594-20-7	Trichlorofluoromethane	CAS:75-69-4
1,2-Dibromo-3-chloropropane	CAS:96-12-8	1,1-Dichloropropene	CAS:563-58-6	1,2,3-Trichloropropane	CAS:96-18-4
1,2-Dibromoethane	CAS:106-93-4	cis-1,3-Dichloropropene	CAS:10061-01-5	Vinyl chloride	CAS:75-01-4
Dibromomethane	CAS:74-95-3	trans-1,3-Dichloropropene	CAS:10061-02-6		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112991	Price : KF

Haloalkanes/ BTEX Calibration Solutions for headspace analysis - 21 components

Benzene	CAS:71-43-2	50 µg/ml	Tetrachloroethene	CAS:127-18-4	0.5 µg/ml
Chloroform	CAS:67-66-3	3 µg/ml	Tetrachloromethane	CAS:56-23-5	0.2 µg/ml
cis-1,2-Dichloroethene	CAS:156-59-2	120 µg/ml	Toluene	CAS:108-88-3	50 µg/ml
trans-1,2-Dichloroethene	CAS:156-60-5	40 µg/ml	Trichloroethene	CAS:79-01-6	2.5 µg/ml
Dichloromethane	CAS:75-09-2	40 µg/ml	1,2,3-Trimethylbenzene	CAS:526-73-8	50 µg/ml
Ethylbenzene	CAS:100-41-4	50 µg/ml	1,2,4-Trimethylbenzene	CAS:95-63-6	50 µg/ml
Isopropylbenzene	CAS:98-82-8	50 µg/ml	1,3,5-Trimethylbenzene	CAS:108-67-8	50 µg/ml
Naphthalene	CAS:91-20-3	100 µg/ml	m-Xylene	CAS:108-38-3	50 µg/ml
n-Propylbenzene	CAS:103-65-1	50 µg/ml	o-Xylene	CAS:95-47-6	50 µg/ml
Styrene	CAS:100-42-5	50 µg/ml	p-Xylene	CAS:106-42-3	50 µg/ml
1,1,1-Trichloroethane	CAS:71-55-6	1.5 µg/ml			
Solvent : Methanol purge & trap		ampoule 1 ml		Ref : F128871	Price : GA

Purgeable Halocarbon Mix - 28 components**(EPA 601)**

Bromodichloromethane	CAS:75-27-4	1,3-Dichlorobenzene	CAS:541-73-1	trans-1,3-Dichloropropene	CAS:10061-02-6
Bromoform	CAS:75-25-2	1,4-Dichlorobenzene	CAS:106-46-7	Methylene chloride	CAS:75-09-2
Bromomethane	CAS:74-83-9	Dichlorodifluoromethane	CAS:75-71-8	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Carbon tetrachloride	CAS:56-23-5	1,1-Dichloroethane	CAS:75-34-3	Tetrachloroethene	CAS:127-18-4
Chlorobenzene	CAS:108-90-7	1,2-Dichloroethane	CAS:107-06-2	1,1,1-Trichloroethane	CAS:71-55-6
Chloroethane	CAS:75-00-3	1,1-Dichloroethene	CAS:75-35-4	1,1,2-Trichloroethane	CAS:79-00-5
Chloroform	CAS:67-66-3	trans-1,2-Dichloroethene	CAS:156-60-5	Trichloroethene	CAS:79-01-6
Chloromethane	CAS:74-87-3	1,2-Dichloropropane	CAS:78-87-5	Trichlorofluoromethane	CAS:75-69-4
Dibromochloromethane	CAS:124-48-1	cis-1,3-Dichloropropene	CAS:10061-01-5	Vinylchloride	CAS:75-01-4
1,2-Dichlorobenzene	CAS:95-50-1				
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F114841	Price : HF

Volatile Halogenated Hydrocarbons - 18 components**(ISO 15009:2002)**

3-Chloropropene	CAS:107-05-1	Dichloromethane	CAS:75-09-2	Tetrachloromethane	CAS:56-23-5
1,2-Dichlorobenzene	CAS:95-50-1	1,2-Dichloropropane	CAS:78-87-5	1,1,1-Trichloroethane	CAS:71-55-6
1,1-Dichloroethane	CAS:75-34-3	cis-1,3-Dichloropropene	CAS:10061-01-5	1,1,2-Trichloroethane	CAS:79-00-5
1,2-Dichloroethane	CAS:107-06-2	trans-1,3-Dichloropropene	CAS:10061-02-6	Trichloroethene	CAS:79-01-6
cis-1,2-Dichloroethene	CAS:156-59-2	Monochlorobenzene	CAS:108-90-7	Trichloromethane	CAS:67-66-3
trans-1,2-Dichloroethene	CAS:156-60-5	Tetrachloroethene	CAS:127-18-4	1,2,3-Trichloropropene	CAS:96-18-4
Solvent : Methanol purge & trap		ampoule 1 ml	4 000 µg/ml	Ref : F062065	Price : BCA
		ampoule 5 ml	4 000 µg/ml	Ref : F109941	Price : EJA

Halomethanes Mixture - 12 components**(EPA 502, EPA 524.2)**

Bromochloromethane	CAS:74-97-5	Carbon tetrachloride	CAS:56-23-5	Dibromomethane	CAS:74-95-3
Bromodichloromethane	CAS:75-27-4	Chloroform	CAS:67-66-3	Dichlorodifluoromethane	CAS:75-71-8
Bromoform	CAS:75-25-2	Chloromethane	CAS:74-87-3	Methylene chloride	CAS:75-09-2
Bromomethane	CAS:74-83-9	Dibromochloromethane	CAS:124-48-1	Trichlorofluoromethane	CAS:75-69-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112961	Price : FA

Halomethanes Mixture - 14 components**(EPA 502, EPA 524.2)**

Chloroethane	CAS:75-00-3	cis-1,2-Dichloroethene	CAS:156-59-2	1,1,1-Trichloroethane	CAS:71-55-6
1,2-Dibromoethane	CAS:106-93-4	trans-1,2-Dichloroethene	CAS:156-60-5	1,1,2-Trichloroethane	CAS:79-00-5
1,1-Dichloroethane	CAS:75-34-3	1,1,1,2-Tetrachloroethane	CAS:630-20-6	Trichloroethene	CAS:79-01-6
1,2-Dichloroethane	CAS:107-06-2	1,1,2,2-Tetrachloroethane	CAS:79-34-5	Vinylchloride	CAS:75-01-4
1,1-Dichloroethene	CAS:75-35-4	Tetrachloroethene	CAS:127-18-4		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112971	Price : FA

Halopropanes Mixture - 9 components**(EPA 502, EPA 524.2)**

1,2-Dibromo-3-chloropropane	CAS:96-12-8	2,2-Dichloropropane	CAS:594-20-7	trans-1,3-Dichloropropene	CAS:10061-02-6
1,2-Dichloropropane	CAS:78-87-5	1,1-Dichloropropene	CAS:563-58-6	Hexachlorobutadiene	CAS:87-68-3
1,3-Dichloropropane	CAS:142-28-9	cis-1,3-Dichloropropene	CAS:10061-01-5	1,2,3-Trichloropropene	CAS:96-18-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112981	Price : FA

Semi-Volatiles Internal Standard Mixture - 6 components**(EPA 625)**

Acenaphthene D10	CAS:15067-26-2	1,4-Dichlorobenzene D4	CAS:3855-82-1	Perylene D12	CAS:1520-96-3
Chrysene D12	CAS:1719-03-5	Naphthalene D8	CAS:1146-65-2	Phenanthrene D10	CAS:1517-22-2
Solvent : Dichloromethane		ampoule 1 ml	2 000 µg/ml	Ref : F115901	Price : HF

VOC Internal Standard Solution - 2 components**(EPA 502 - EPA 8021B)**

2-Bromo-1-chloropropane	CAS:3017 95 6	Fluorobenzene	CAS:462-06-6		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112761	Price : DF

VOC Internal Standard Mixture - 2 components**(EPA 502)**

2-Bromo-1-chloropropane	CAS:3017-95-6	1-Chloro-2-fluorobenzene	CAS:348-51-6		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112771	Price : DF

VOC Internal Standard Mixture (Surrogate Standard Mixture) - 3 components**(EPA 8240)**

Bromochloromethane	CAS:74-97-5	Chlorobenzene D5	CAS:3114-55-4	1,4-Difluorobenzene	CAS:540-36-3
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F119221	Price : DF

VOC Purgeable Internal Standards - 3 components**(EPA 601)**

Bromochloromethane	CAS:74-97-5	2-Bromo-1-chloropropane	CAS:3017-95-6	1,4-Dichlorobutane	CAS:110-56-5
Solvent : Methanol purge & trap		ampoule 1 ml	20 µg/ml	Ref : F114861	Price : CF

VOC Internal & Surrogate Standard Solution - 3 components**(EPA 8010B)**

4-Bromochlorobenzene	CAS:106-39-8	Bromochloromethane	CAS:74-97-5	4-Bromofluorobenzene	CAS:460-00-4
Solvent : Methanol purge & trap		ampoule 1 ml	1 500 µg/ml	Ref : F116241	Price : CF

VOC Internal & Surrogate Standard Solution - 3 components**(EPA 624)**

4-Bromofluorobenzene	CAS:460-00-4	Fluorobenzene	CAS:462-06-6	Pentafluorobenzene	CAS:363-72-4
Solvent : Methanol purge & trap		ampoule 1 ml	20 000 µg/ml	Ref : F115861	Price : DF

VOC Internal & Surrogate Standard Solution - 3 components**(EPA 624)**

Bromochloromethane	CAS:74-97-5	2-Bromo-1-chloropropane	CAS:3017-95-6	1,4-Dichlorobutane	CAS:110-56-5
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F115841	Price : CF
		ampoule 1 ml	20 000 µg/ml	Ref : F115851	Price : FA

VOC Surrogate Standard Solution - 2 components**(EPA 8021B)**

4-Bromochlorobenzene	CAS:106-39-8	1,4-Dichlorobutane	CAS:110-56-5		
Solvent : Methanol purge & trap		ampoule 1 ml	1 500 µg/ml	Ref : F116361	Price : CF

VOC Surrogate Standard Mixture 3 components**(EPA 8240 ; EPA 8260)**

4-Bromofluorobenzene	CAS:460-00-4	1,2-Dichloroethane D4	CAS:17060-07-0	Toluene D8	CAS:2037-26-5
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F119231	Price : DF

VOC Stock Calibration Standard Solution 1- 14 components**(ISO 15680:2003)**

Bromochloromethane	CAS:74-97-5	1,1-Dichloroethane	CAS:75-34-3	2,2-Dichloropropane	CAS:594-20-7
Bromomethane	CAS:74-83-9	1,1-Dichloroethene	CAS:75-35-4	1,1,1-Trichloroethane	CAS:71-55-6
Chloroethane	CAS:75-00-3	cis-1,2-Dichloroethene	CAS:156-59-2	Trichlorofluoromethane	CAS:75-69-4
Chloromethane	CAS:74-87-3	trans-1,2-Dichloroethene	CAS:156-60-5	Trichloromethane	CAS:67-66-3
Dichlorodifluoromethane	CAS:75-71-8	Dichloromethane	CAS:75-09-2		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062085	Price : JF
		ampoule 5 ml	2 000 µg/ml	Ref : F109961	Price : DCA

VOC Stock Calibration Standard Solution 2 - 13 components**(ISO15680:2003)**

Benzene	CAS:71-43-2	1,2-Dichloropropane	CAS:78-87-5	Toluene	CAS:108-88-3
Bromodichloromethane	CAS:75-27-4	1,1-Dichloropropene	CAS:563-58-6	1,1,2-Trichloroethane	CAS:79-00-5
1,2-Dibromoethane	CAS:106-93-4	cis-1,3-Dichloropropene	CAS:10061-01-5	Trichloroethene	CAS:79-01-6
Dibromomethane	CAS:74-95-3	trans-1,3-Dichloropropene	CAS:10061-02-6		
1,2-Dichloroethane	CAS:107-06-2	Tetrachloromethane	CAS:56-23-5		
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F062095	Price : JA
		ampoule 5 ml	2 000 µg/ml	Ref : F109971	Price : CJA

TCLP Volatiles Mixture - 11 components**(EPA 1311)**

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	Tetrachloroethene	CAS:127-18-4
Carbon tetrachloride	CAS:56-23-5	1,2-Dichloroethane	CAS:107-06-2	Trichloroethene	CAS:79-01-6
Chlorobenzene	CAS:108-90-7	1,1-Dichloroethene	CAS:75-35-4	Vinylchloride	CAS:75-01-4
Chloroform	CAS:67-66-3	MEK	CAS:78-93-3		
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F116091	Price : EF

TCLP Semi-Volatile Spiking Mixture - 13 components**(EPA 1311)**

1,4-Dichlorobenzene	CAS:106-46-7	2-Methylphenol	CAS:95-48-7	Pyridine	CAS:110-86-1
2,4-Dinitrotoluene	CAS:121-14-2	3-Methylphenol	CAS:108-39-4	2,4,5-Trichlorophenol	CAS:95-95-4
Hexachlorobenzene	CAS:118-74-1	4-Methylphenol	CAS:106-44-5	2,4,6-Trichlorophenol	CAS:88-06-2
Hexachlorobutadiene	CAS:87-68-3	Nitrobenzene	CAS:98-95-3		
Hexachloroethane	CAS:67-72-1	Pentachlorophenol	CAS:87-86-5		
Solvent : Dichloromethane		ampoule 1 ml	1 000 µg/ml	Ref : F116101	Price : FF

Trihalomethanes Standard Solution - 4 components**(EPA 501)**

Bromodichloromethane	CAS:75-27-4	Dibromochloromethane	CAS:124-48-1	Chloroform	CAS:67-66-3
Bromoform	CAS:75-25-2				
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F112711	Price : CF
		ampoule 1 ml	200 µg/ml	Ref : F112721	Price : CF
		ampoule 1 ml	2 000 µg/ml	Ref : F112741	Price : CF
		ampoule 1 ml	5 000 µg/ml	Ref : F112731	Price : DA

VOC Performance Check Solution - 8 components**(EPA 601)**

Benzene	CAS:71-43-2	1,2-Dichloroethane	CAS:107-06-2	Trichloroethene	CAS:79-01-6
Carbontetrachloride	CAS:56-23-5	1,1-Dichloroethene	CAS:75-35-4	Vinylchloride	CAS:75-01-4
1,4-Dichlorobenzene	CAS:106-46-7	1,1,1-Trichloroethane	CAS:71-55-6		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F114871	Price : EA

VOC Standard Solution - 2 components**(EPA 8011)**

1,2-Dibromo-3-chloropropane	CAS:96-12-8	1,2-Dibromoethane	CAS:106-93-4		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F116251	Price : CF
		ampoule 1 ml	2 000 µg/ml	Ref : F116261	Price : CF

VOC Standard Solution - 4 components**(EPA 8015B)**

Diethylether	CAS:60-29-7	MEK	CAS:78-93-3	MIBK	CAS:108-10-1
Ethanol	CAS:64-17-5				
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F116271	Price : DF

VOC Standard Solution - 6 components**(ASTM D5580, EPA 502, EPA 524.2, EPA 8021B)**

Bromomethane	CAS:74-83-9	Chloromethane	CAS:74-87-3	Fluorotrichloromethane	CAS:75-69-4
Chloroethane	CAS:75-00-3	Dichlorodifluoromethane	CAS:75-71-8	Vinylchloride	CAS:75-01-4
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F112901	Price : FA
		ampoule 1 ml	2000 µg/ml	Ref : F112911	Price : HF

VOC Standard Solution - 8 components**(EPA 8260)**

Acrolein (2-Propenal)	CAS:107-02-8	Chloralhydrate	CAS:302-17-0	Paraldehyde	CAS:123-63-7
Acrylonitrile	CAS:107-13-1	2-Chloroethyl-vinylether	CAS:110-75-8	Vinylacetate	CAS:108-05-4
Bromoacetone	CAS:598-31-2	Chloroprene	CAS:126-99-8		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F119281	Price : EF

VOC Standard Solution - 12 components

Benzene	CAS:71-43-2	1,2-Dichloroethane	CAS:107-06-2	Tetrachloromethane	CAS:56-23-5
Bromodichloromethane	CAS:75-27-4	cis-1,2-Dichloroethene	CAS:156-59-2	Tribromomethane	CAS:75-25-2
Chloroform	CAS:67-66-3	Dichloromethane	CAS:75-09-2	1,1,1-Trichloroethane	CAS:71-55-6
Dibromochloromethane	CAS:123-91-1	Tetrachloroethene	CAS:127-18-4	Trichloroethene	CAS:79-01-6
Solvent : Methanol purge & trap		ampoule 1 ml	1 000 µg/ml	Ref : F127731	Price : EH

VOC Standard Solution - 12 components**(EPA 8015B)**

Acetonitrile	CAS:75-05-8	Ethyl alcohol	CAS:64-17-5	Methacrylonitrile	CAS:126-98-7
Acrylamide	CAS:79-06-1	Ethyl methacrylate	CAS:97-63-2	Methyl methacrylate	CAS:80-62-6
Diethylether	CAS:60-29-7	Isobutyl alcohol	CAS:78-83-1	MIBK	CAS:108-10-1
1,4-Dioxane	CAS:123-91-1	MEK	CAS:78-93-3	Propionitrile	CAS:107-12-0
Solvent : Methanol purge & trap		ampoule 1 ml	100 µg/ml	Ref : F116281	Price : EF

VOC Standard Solution - 14 components**(EPA 8260)**

Allyl chloride	CAS:107-05-1	Iodomethane	CAS:74-88-4	2-Picoline	CAS:109-06-8
Benzyl chloride	CAS:100-44-7	Nitrobenzene	CAS:98-95-3	n-Propyl amine	CAS:107-10-8
cis-1,4-Dichloro-2-butene	CAS:1476-11-5	2-Nitropropane	CAS:79-46-9	Pyridine	CAS:110-86-1
trans-1,4-Dichloro-2-butene	CAS:110-57-6	N-Nitrosodi-n-butylamine	CAS:924-16-3	o-Toluidine	CAS:95-53-4
Hexachloroethane	CAS:67-72-1	Pentachloroethane	CAS:76-01-7		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F119271	Price : FA

VOC Standard Solution - 15 components**(EPA 8260)**

Acetone	CAS:67-64-1	Dioxan	CAS:123-91-1	Methanol	CAS:67-56-1
1-Butanol	CAS:71-36-3	Ethanol	CAS:64-17-5	MIBK	CAS:108-10-1
2-Butanone	CAS:78-93-3	Ethyl Acetate	CAS:141-78-6	2-Pentanone	CAS:107-87-9
tert-butyl alcohol	CAS:75-65-0	2-Hexanone	CAS:591-78-6	1-Propanol	CAS:71-23-8
Diethylether	CAS:60-29-7	Isobutyl alcohol	CAS:78-83-1	2-Propanol	CAS:67-63-0
Solvent : Water		ampoule 1 ml	200 µg/ml	Ref : F119251	Price : FA

VOC Standard Solution - 16 components**(EPA 8240)**

Acrolein (2-Propenal)	CAS:107-02-8	2-Chloroethanol	CAS:107-07-3	Trichlorofluoromethane	CAS:75-69-4
Acrylonitrile	CAS:107-13-1	2-Chloroethyl-vinylether	CAS:110-75-8	1,2,3-Trichloropropane	CAS:96-18-4
Bromomethane	CAS:74-83-9	Chloromethane	CAS:74-87-3	Vinylacetate	CAS:108-05-4
Bromoacetone	CAS:598-31-2	Chloroprene	CAS:126-99-8	Vinylchloride	CAS:75-01-4
Chloralhydrate	CAS:302-17-0	Dichlorodifluoromethane	CAS:75-71-8		
Chloroethane	CAS:75-00-3	Ethylene Oxide	CAS:75-21-8		
Solvent : Methanol purge & trap/Xylene (1/traces)		ampoule 1 ml	200 µg/ml	Ref : F119211	Price : FF

VOC Standard Solution - 17 components**(EPA 8260)**

Acetonitrile	CAS:75-05-8	1,3-Dichloro-2-propanol	CAS:96-23-1	Malononitrile	CAS:109-77-3
Allyl alcohol	CAS:107-18-6	1,2,3,4-Diepoxybutane	CAS:1464-53-5	Methacrylonitrile	CAS:126-98-7
Carbon disulfide	CAS:75-15-0	Epichlorhydrin	CAS:106-89-8	Methyl methacrylate	CAS:80-62-6
2-Chloroethanol	CAS:107-07-3	Ethylcyanid	CAS:107-12-0	Propargyl alcohol	CAS:107-19-7
3-Chloropropionitrile	CAS:542-76-7	Ethyl methacrylate	CAS:97-63-2	b-Propiolactone	CAS:57-57-8
Crotonaldehyde	CAS:123-73-9	2-Hydroxypropionitrile	CAS:78-97-7		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F119261	Price : GA

VOC Standard Solution - 24 components**(EPA 524.2)**

Acetone	CAS:67-64-1	Diethylether	CAS:60-29-7	Methyl methacrylate	CAS:80-62-6
Acrylonitrile	CAS:107-13-1	Ethyl methacrylate	CAS:97-63-2	MIBK	CAS:108-10-1
Allyl chloride	CAS:107-05-1	Hexachloroethane	CAS:67-72-1	MTBE	CAS:1634-04-4
Carbon disulfide	CAS:75-15-0	2-Hexanone	CAS:591-78-6	Nitrobenzene	CAS:98-95-3
Chloroacetonitrile	CAS:107-14-2	MEK	CAS:78-93-3	2-Nitropropane	CAS:79-46-9
1-Chlorobutane	CAS:109-69-3	Methacrylonitrile	CAS:126-98-7	Pentachloroethane	CAS:76-01-7
trans-1,4-Dichloro-2-butene	CAS:110-57-6	Methyl acrylate	CAS:96-33-3	Propionitrile	CAS:107-12-0
1,1-Dichloro-2-propanone	CAS:513-88-2	Methyl iodide	CAS:74-88-4	THF	CAS:109-99-9
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F113951	Price : GF

VOC Standard Solution - 28 components**(EPA 8240)**

Acetonitrile	CAS:75-05-8	1,2,3,4-Diepoxybutane	CAS:1464-53-5	2-Picoline	CAS:109-06-8
Allyl alcohol	CAS:107-18-6	1,4-Dioxane	CAS:123-91-1	Propargyl alcohol	CAS:107-19-7
Allyl chloride	CAS:107-05-1	Epichlorohydrin	CAS:106-89-8	b-Propiolactone	CAS:57-57-8
Benzyl chloride	CAS:100-44-7	Ethyl methacrylate	CAS:97-63-2	Propionitrile	CAS:107-12-0
2-Chloroethanol	CAS:107-07-3	2-Hydroxypropionitrile	CAS:78-97-7	n- Propyl amine	CAS:107-10-8
3-Chloropropionitrile	CAS:542-76-7	Isobutyl alcohol	CAS:78-83-1	Pyridine	CAS:110-86-1
1,2-Dibromo-3-chloropropane	CAS:96-12-8	Malononitrile	CAS:109-77-3	1,1,1,2-Tetrachloroethane	CAS:630-20-6
1,2-Dibromoethane	CAS:106-93-4	Methacrylonitrile	CAS:126-98-7	1,2,3-Trichloropropane	CAS:96-18-4
Dibromomethane	CAS:74-95-3	Methyl methacrylate	CAS:80-62-6		
1,3-Dichloropropane-2-ol	CAS:96-23-1	Pentachloroethane	CAS:76-01-7		

Solvent : Methanol purge & trap

ampoule 1 ml

200 µg/ml

Ref : F119201

Price : GA

VOC Standard Solution - 28 components

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	1,2,3-Trichlorobenzene	CAS:87-61-6
Bromobenzene	CAS:108-86-1	Ethylbenzene	CAS:100-41-4	1,2,4-Trichlorobenzene	CAS:120-82-1
n-Butylbenzene	CAS:104-51-8	Hexachlorobutadiene	CAS:87-68-3	Trichloroethene	CAS:79-01-6
sec-Butylbenzene	CAS:135-98-8	Isopropylbenzene	CAS:98-82-8	1,2,4-Trimethylbenzene	CAS:95-63-6
tert-Butylbenzene	CAS:98-06-6	4-Isopropyltoluene	CAS:99-87-6	1,3,5-Trimethylbenzene	CAS:108-67-8
Chlorobenzene	CAS:108-90-7	Naphthalene	CAS:91-20-3	o-Xylene	CAS:95-47-6
2-Chlorotoluene	CAS:95-49-8	n-Propylbenzene	CAS:103-65-1	m-Xylene	CAS:108-38-3
4-Chlorotoluene	CAS:106-43-4	Styrene	CAS:100-42-5	p-Xylene	CAS:106-42-3
1,2-Dichlorobenzene	CAS:95-50-1	Tetrachloroethene	CAS:127-18-4		
1,3-Dichlorobenzene	CAS:541-73-1	Toluene	CAS:108-88-3		

Solvent : Methanol purge & trap

ampoule 1 ml

2 000 µg/ml

Ref : F113081

Price : HA

VOC Standard Solution - 31 components**(EPA 624)**

Benzene	CAS:71-43-2	1,2-Dichlorobenzene	CAS:95-50-1	Methylene chloride	CAS:75-09-2
Bromodichloromethane	CAS:75-27-4	1,3-Dichlorobenzene	CAS:541-73-1	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Bromoform	CAS:75-25-2	1,4-Dichlorobenzene	CAS:106-46-7	Tetrachloroethene	CAS:127-18-4
Bromomethane	CAS:74-83-9	1,1-Dichloroethane	CAS:75-34-3	Toluene	CAS:108-88-3
Carbon tetrachloride	CAS:56-23-5	1,2-Dichloroethane	CAS:107-06-2	1,1,1-Trichloroethane	CAS:71-55-6
Chlorobenzene	CAS:108-90-7	1,1-Dichloroethene	CAS:75-35-4	1,1,2-Trichloroethane	CAS:79-00-5
Chloroethane	CAS:75-00-3	trans-1,2-Dichloroethene	CAS:156-60-5	Trichloroethene	CAS:79-01-6
2-Chloroethyl-vinylether	CAS:110-75-8	1,2-Dichloropropane	CAS:78-87-5	Trichlorofluoromethane	CAS:75-69-4
Chloroform	CAS:67-66-3	cis-1,3-Dichloropropene	CAS:10061-01-5	Vinylchloride	CAS:75-01-4
Chloromethane	CAS:74-87-3	trans-1,3-Dichloropropene	CAS:10061-02-6		
Dibromochloromethane	CAS:124-48-1	Ethylbenzene	CAS:100-41-4		

Solvent : Methanol purge & trap

ampoule 1 ml

20 µg/ml

Ref : F115831

Price : FF

ampoule 1 ml

200 µg/ml

Ref : F115821

Price : GF

VOC Standard Solution - 33 components**(EPA 8010B)**

Benzylchloride	CAS:100-44-7	1,2-Dichlorobenzene	CAS:95-50-1	trans-1,3-Dichloropropene	CAS:10061-02-6
Bromobenzene	CAS:108-86-1	1,3-Dichlorobenzene	CAS:541-73-1	Methylene chloride	CAS:75-09-2
Bromoform	CAS:75-25-2	1,4-Dichlorobenzene	CAS:106-46-7	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromomethane	CAS:74-83-9	Dichlorobromomethane	CAS:75-27-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Carbon tetrachloride	CAS:56-23-5	Dichlorodifluoromethane	CAS:75-71-8	Tetrachloroethene	CAS:127-18-4
Chlorobenzene	CAS:108-90-7	1,1-Dichloroethane	CAS:75-34-3	1,1,1-Trichloroethane	CAS:71-55-6
Chloroethane	CAS:75-00-3	1,2-Dichloroethane	CAS:107-06-2	1,1,2-Trichloroethane	CAS:79-00-5
Chloroform	CAS:67-66-3	1,1-Dichloroethene	CAS:75-35-4	Trichloroethene	CAS:79-01-6
Chloromethane	CAS:74-87-3	trans-1,2-Dichloroethene	CAS:156-60-5	Trichlorofluoromethane	CAS:75-69-4
Dibromochloromethane	CAS:124-48-1	1,2-Dichloropropane	CAS:78-87-5	1,2,3-Trichloropropane	CAS:96-18-4
Dibromomethane	CAS:74-95-3	cis-1,3-Dichloropropene	CAS:10061-01-5	Vinylchloride	CAS:75-01-4

Solvent : Methanol purge & trap

ampoule 1 ml

200 µg/ml

Ref : F116181

Price : BAF

VOC Standard Solution - 34 components**(EPA 8240)**

Acetone	CAS:67-64-1	1,2-Dichloroethane	CAS:107-06-2	Styrene	CAS:100-42-5
Benzene	CAS:71-43-2	1,1-Dichloroethene	CAS:75-35-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Bromodichloromethane	CAS:75-27-4	trans-1,2-Dichloroethene	CAS:156-60-5	Tetrachloroethene	CAS:127-18-4
Bromoform	CAS:75-25-2	1,2-Dichloropropane	CAS:78-87-5	Toluene	CAS:108-88-3
2-Butanone	CAS:78-93-3	cis-1,3-Dichloropropene	CAS:10061-01-5	1,1,1-Trichloroethane	CAS:71-55-6
Carbon disulfide	CAS:75-15-0	trans-1,3-Dichloropropene	CAS:10061-02-6	1,1,2-Trichloroethane	CAS:79-00-5
Carbon tetrachloride	CAS:56-23-5	Ethanol	CAS:64-17-5	Trichloroethene	CAS:79-01-6
Chlorobenzene	CAS:108-90-7	Ethylbenzene	CAS:100-41-4	o-Xylene	CAS:95-47-6
Chlorodibromomethane	CAS:124-48-1	2-Hexanone	CAS:591-78-6	m-Xylene	CAS:108-38-3
Chloroform	CAS:67-66-3	Methyl iodide	CAS:74-88-4	p-Xylene	CAS:106-42-3
1,4-Dichloro-2-butene	CAS:764-41-0	Methylene chloride	CAS:75-09-2		
1,1-Dichloroethane	CAS:75-34-3	4-Methyl-2-pentanone	CAS:108-10-1		

Solvent : Methanol purge & trap

ampoule 1 ml

200 µg/ml

Ref : F119191

Price : HF

VOC Standard Solution - 40 components**(EPA 8010B)**

Allyl chloride	CAS:107-05-1	Dibromomethane	CAS:74-95-3	trans-1,3-Dichloropropene	CAS:10061-02-6
Bromobenzene	CAS:108-86-1	1,2-Dichlorobenzene	CAS:95-50-1	Ethylene dibromide	CAS:106-93-4
Bromoform	CAS:75-25-2	1,3-Dichlorobenzene	CAS:541-73-1	Methylene chloride	CAS:75-09-2
Bromomethane	CAS:74-83-9	1,4-Dichlorobenzene	CAS:106-46-7	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Carbon tetrachloride	CAS:56-23-5	Dichlorobromomethane	CAS:75-27-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Chlorobenzene	CAS:108-90-7	1,4-Dichloro-2-butene	CAS:764-41-0	Tetrachloroethene	CAS:127-18-4
Chloroethane	CAS:75-00-3	Dichlorodifluoromethane	CAS:75-71-8	1,1,1-Trichloroethane	CAS:71-55-6
2-Chloroethanol	CAS:107-07-3	1,1-Dichloroethane	CAS:75-34-3	1,1,2-Trichloroethane	CAS:79-00-5

Chloroform	CAS:67-66-3	1,2-Dichloroethane	CAS:107-06-2	Trichloroethene	CAS:79-01-6
1-Chlorohexane	CAS:544-10-5	1,1-Dichloroethene	CAS:75-35-4	Trichlorofluoromethane	CAS:75-69-4
Chloromethane	CAS:74-87-3	trans-1,2-Dichloroethene	CAS:156-60-5	1,2,3-Trichloropropane	CAS:96-18-4
4-Chlorotoluene	CAS:106-43-4	1,2-Dichloropropane	CAS:78-87-5	Vinylchloride	CAS:75-01-4
Dibromochloromethane	CAS:124-48-1	1,3-Dichloro-2-propanol	CAS:96-23-1		
1,2-Dibromo-3-chloropropane	CAS:96-12-8	cis-1,3-Dichloropropene	CAS:10061-01-5		
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F116201	Price : BBF

VOC Standard Solution - 48 components

(EPA 8260)

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	Naphthalene	CAS:91-20-3
Bromochloromethane	CAS:74-97-5	Dichlorodifluoromethane	CAS:75-71-8	Styrene	CAS:100-42-5
Bromodichloromethane	CAS:75-27-4	1,1-Dichloroethane	CAS:75-34-3	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromoform	CAS:75-25-2	1,2-Dichloroethane	CAS:107-06-2	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Bromomethane	CAS:74-83-9	1,1-Dichloroethene	CAS:75-35-4	Tetrachloroethene	CAS:127-18-4
Carbon tetrachloride	CAS:56-23-5	cis-1,2-Dichloroethene	CAS:156-59-2	Toluene	CAS:108-88-3
Chlorobenzene	CAS:108-90-7	trans-1,2-Dichloroethene	CAS:156-60-5	1,2,4-Trichlorobenzene	CAS:120-82-1
Chlorodibromomethane	CAS:124-48-1	1,2-Dichloropropane	CAS:78-87-5	1,1,1-Trichloroethane	CAS:71-55-6
Chloroethane	CAS:75-00-3	1,3-Dichloropropane	CAS:142-28-9	1,1,2-Trichloroethane	CAS:79-00-5
Chloroform	CAS:67-66-3	2,2-Dichloropropane	CAS:594-20-7	Trichloroethene	CAS:79-01-6
Chloromethane	CAS:74-87-3	cis-1,3-Dichloropropene	CAS:10061-01-5	Trichlorofluoromethane	CAS:75-69-4
1,2-Dibromo-3-chloropropane	CAS:96-12-8	trans-1,3-Dichloropropene	CAS:10061-02-6	1,2,3-Trichloropropane	CAS:96-18-4
1,2-Dibromoethane	CAS:106-93-4	Ethylbenzene	CAS:100-41-4	Vinylchloride	CAS:75-01-4
Dibromomethane	CAS:74-95-3	Hexachlorobutadiene	CAS:87-68-3	o-Xylene	CAS:95-47-6
1,2-Dichlorobenzene	CAS:95-50-1	Isopropylbenzene	CAS:98-82-8	m-Xylene	CAS:108-38-3
1,3-Dichlorobenzene	CAS:541-73-1	Methylene chloride	CAS:75-09-2	p-Xylene	CAS:106-42-3
Solvent : Methanol purge & trap		ampoule 1 ml	200 µg/ml	Ref : F119241	Price : KK

VOC Standard Solution - 52 components

Benzene	CAS:71-43-2	100 mg/l	2,5-Dichloronitrobenzene	CAS:89-61-2	100 mg/l
Bromodichloromethane	CAS:75-27-4	100 mg/l	3,4-Dichloronitrobenzene	CAS:99-54-7	100 mg/l
Chlorobenzene	CAS:108-90-7	100 mg/l	cis-1,3-Dichloropropene	CAS:10061-01-5	50 mg/l
1-Chloro-2,4-dinitrobenzene	CAS:97-00-7	100 mg/l	trans-1,3-Dichloropropene	CAS:10061-02-6	50 mg/l
2-Chloroethanol	CAS:107-07-3	100 mg/l	alpha-alpha-Dichlorotoluene	CAS:98-87-3	100 mg/l
Chloroform	CAS:67-66-3	100 mg/l	Ethylbenzene	CAS:100-41-4	100 mg/l
1-Chloro-2-nitrobenzene	CAS:88-73-3	100 mg/l	Hexachloro-1,3-butadiene	CAS:87-68-3	100 mg/l
1-Chloro-3-nitrobenzene	CAS:121-73-3	100 mg/l	Hexachloroethane	CAS:67-72-1	100 mg/l
1-Chloro-4-nitrobenzene	CAS:100-00-5	100 mg/l	Isopropylbenzene	CAS:98-82-8	100 mg/l
Chloroprene	CAS:126-99-8	100 mg/l	Nitrobenzene	CAS:98-95-3	100 mg/l
3-Chloropropene	CAS:107-05-1	100 mg/l	2-Nitrotoluene	CAS:88-72-2	100 mg/l
2-Chlorotoluene	CAS:95-49-8	100 mg/l	1,2,4,5-Tetrachlorobenzene	CAS:95-94-3	100 mg/l
3-Chlorotoluene	CAS:108-41-8	100 mg/l	1,1,2,2-Tetrachloroethane	CAS:79-34-5	100 mg/l
4-Chlorotoluene	CAS:106-43-4	100 mg/l	Tetrachloroethene	CAS:127-18-4	100 mg/l
alpha-Chlorotoluene	CAS:100-44-7	100 mg/l	Tetrachloromethane	CAS:56-23-5	100 mg/l
Dibromochloromethane	CAS:124-48-1	100 mg/l	Toluene	CAS:108-88-3	100 mg/l
1,2-Dichlorobenzene	CAS:95-50-1	100 mg/l	Tribromomethane	CAS:75-25-2	100 mg/l
1,3-Dichlorobenzene	CAS:541-73-1	100 mg/l	1,2,3-Trichlorobenzene	CAS:87-61-6	100 mg/l
1,4-Dichlorobenzene	CAS:106-46-7	100 mg/l	1,2,4-Trichlorobenzene	CAS:120-82-1	100 mg/l
1,1-Dichloroethane	CAS:75-34-3	100 mg/l	1,3,5-Trichlorobenzene	CAS:108-70-3	100 mg/l
1,2-Dichloroethane	CAS:107-06-2	100 mg/l	1,1,1-Trichloroethane	CAS:71-55-6	100 mg/l
1,1-Dichloroethene	CAS:75-35-4	100 mg/l	1,1,2-Trichloroethane	CAS:79-00-5	100 mg/l
cis-1,2-Dichloroethene	CAS:156-59-2	50 mg/l	Trichloroethene	CAS:79-01-6	100 mg/l
trans-1,2-Dichloroethene	CAS:156-60-5	50 mg/l	o-Xylene	CAS:95-47-6	100 mg/l
Dichloromethane	CAS:75-09-2	100 mg/l	m-Xylene	CAS:108-38-3	100 mg/l
2,3-Dichloronitrobenzene	CAS:3209-22-1	100 mg/l	p-Xylene	CAS:106-42-3	100 mg/l
Solvent : Methanol purge & trap		ampoule 10 ml		Ref : F112261	Price : KK

VOC Standard Solution - 54 components

(EPA 502, 8021B)

Benzene	CAS:71-43-2	1,3-Dichlorobenzene	CAS:541-73-1	Naphthalene	CAS:91-20-3
Bromobenzene	CAS:108-86-1	1,4-Dichlorobenzene	CAS:106-46-7	n-Propylbenzene	CAS:103-65-1
Bromochloromethane	CAS:74-97-5	1,1-Dichloroethane	CAS:75-34-3	Styrene	CAS:100-42-5
Bromodichloromethane	CAS:75-27-4	1,2-Dichloroethane	CAS:107-06-2	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromoform	CAS:75-25-2	1,1-Dichloroethene	CAS:75-35-4	1,1,2,2-Tetrachloroethane	CAS:79-34-5
n-Butylbenzene	CAS:104-51-8	cis-1,2-Dichloroethene	CAS:156-59-2	Tetrachloroethene	CAS:127-18-4
sec-Butylbenzene	CAS:135-98-8	trans-1,2-Dichloroethene	CAS:156-60-5	Toluene	CAS:108-88-3
tert-Butylbenzene	CAS:98-06-6	1,2-Dichloropropane	CAS:78-87-5	1,2,3-Trichlorobenzene	CAS:87-61-6
Carbon tetrachloride	CAS:56-23-5	1,3-Dichloropropane	CAS:142-28-9	1,2,4-Trichlorobenzene	CAS:120-82-1
Chlorobenzene	CAS:108-90-7	2,2-Dichloropropane	CAS:594-20-7	1,1,1-Trichloroethane	CAS:71-55-6
Chloroform	CAS:67-66-3	1,1-Dichloropropene	CAS:563-58-6	1,1,2-Trichloroethane	CAS:79-00-5
2-Chlorotoluene	CAS:95-49-8	cis-1,3-Dichloropropene	CAS:10061-01-5	Trichloroethene	CAS:79-01-6
4-Chlorotoluene	CAS:106-43-4	trans-1,3-Dichloropropene	CAS:10061-02-6	1,2,3-Trichloropropane	CAS:96-18-4
Dibromochloromethane	CAS:124-48-1	Ethylbenzene	CAS:100-41-4	1,2,4-Trimethylbenzene	CAS:95-63-6
1,2-Dibromo-3-chloropropane	CAS:96-12-8	Hexachlorobutadiene	CAS:87-68-3	1,3,5-Trimethylbenzene	CAS:108-67-8
1,2-Dibromoethane	CAS:106-93-4	Isopropylbenzene	CAS:98-82-8	o-Xylene	CAS:95-47-6
Dibromomethane	CAS:74-95-3	4-Isopropyltoluene	CAS:99-87-6	m-Xylene	CAS:108-38-3
1,2-Dichlorobenzene	CAS:95-50-1	Methylene chloride	CAS:75-09-2	p-Xylene	CAS:106-42-3
Solvent : Methanol purge & trap		ampoule 1 ml	2 000 µg/ml	Ref : F112781	Price : BDA

VOC Standard Solution - 60 components

Benzene	CAS:71-43-2	1,4-Dichlorobenzene	CAS:106-46-7	Styrene	CAS:100-42-5
Bromobenzene	CAS:108-86-1	1,1-Dichloroethane	CAS:75-34-3	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromochloromethane	CAS:74-97-5	1,2-Dichloroethane	CAS:107-06-2	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Bromodichloromethane	CAS:75-27-4	1,1-Dichloroethene	CAS:75-35-4	Tetrachloroethene	CAS:127-18-4
n-Butylbenzene	CAS:104-51-8	cis-1,2-Dichloroethene	CAS:156-59-2	Tetrachloromethane	CAS:56-23-5
sec-Butylbenzene	CAS:135-98-8	trans-1,2-Dichloroethene	CAS:156-60-5	Toluene	CAS:108-88-3
tert-Butylbenzene	CAS:98-06-6	Dichloromethane	CAS:75-09-2	Tribromomethane	CAS:75-25-2
Chlorobenzene	CAS:108-90-7	1,2-Dichloropropane	CAS:78-87-5	1,2,3-Trichlorobenzene	CAS:87-61-6
Chloroform	CAS:67-66-3	1,3-Dichloropropane	CAS:142-28-9	1,2,4-Trichlorobenzene	CAS:120-82-1
Chloroprene	CAS:126-99-8	1,1-Dichloropropene	CAS:563-58-6	1,3,5-Trichlorobenzene	CAS:108-70-3
3-Chloropropene	CAS:107-05-1	cis-1,3-Dichloropropene	CAS:10061-01-5	1,1,1-Trichloroethane	CAS:71-55-6
2-Chlorotoluene	CAS:95-49-8	trans-1,3-Dichloropropene	CAS:10061-02-6	1,1,2-Trichloroethane	CAS:79-00-5
3-Chlorotoluene	CAS:108-41-8	2,3-Dichloro-1-propene	CAS:78-88-6	Trichloroethene	CAS:79-01-6
4-Chlorotoluene	CAS:106-43-4	Ethylbenzene	CAS:100-41-4	1,2,3-Trimethylbenzene	CAS:526-73-8
Dibromochloromethane	CAS:124-48-1	Hexachloro-1,3-butadiene	CAS:87-68-3	1,2,4-Trimethylbenzene	CAS:95-63-6
1,2-Dibromo-3-chloropropane	CAS:96-12-8	Hexachloroethane	CAS:67-72-1	1,3,5-Trimethylbenzene	CAS:108-67-8
1,2-Dibromoethane	CAS:106-93-4	Isopropylbenzene	CAS:98-82-8	Vinylchloride	CAS:75-01-4
Dibromomethane	CAS:74-95-3	4-Isopropyltoluene	CAS:99-87-6	o-Xylene	CAS:95-47-6
1,2-Dichlorobenzene	CAS:95-50-1	Naphthalene	CAS:91-20-3	m-Xylene	CAS:108-38-3
1,3-Dichlorobenzene	CAS:541-73-1	n-Propylbenzene	CAS:103-65-1	p-Xylene	CAS:106-42-3

Solvent : Methanol purge & trap

ampoule 1 ml 200 µg/ml

Ref : F127701

Price : BGA

(Hexachloro-1,3-butadiene : 50 µg/ml)

VOC Standard Solution - 60 components

(EPA 502, 524.2, 8021.B)

Benzene	CAS:71-43-2	1,2-Dichlorobenzene	CAS:95-50-1	Naphthalene	CAS:91-20-3
Bromobenzene	CAS:108-86-1	1,3-Dichlorobenzene	CAS:541-73-1	n-Propylbenzene	CAS:103-65-1
Bromochloromethane	CAS:74-97-5	1,4-Dichlorobenzene	CAS:106-46-7	Styrene	CAS:100-42-5
Bromodichloromethane	CAS:75-27-4	Dichlorodifluoromethane	CAS:75-71-8	1,1,1,2-Tetrachloroethane	CAS:630-20-6
Bromoform	CAS:75-25-2	1,1-Dichloroethane	CAS:75-34-3	1,1,2,2-Tetrachloroethane	CAS:79-34-5
Bromomethane	CAS:74-83-9	1,2-Dichloroethane	CAS:107-06-2	Tetrachloroethene	CAS:127-18-4
n-Butylbenzene	CAS:104-51-8	1,1-Dichloroethene	CAS:75-35-4	Toluene	CAS:108-88-3
sec-Butylbenzene	CAS:135-98-8	cis-1,2-Dichloroethene	CAS:156-59-2	1,2,3-Trichlorobenzene	CAS:87-61-6
tert-Butylbenzene	CAS:98-06-6	trans-1,2-Dichloroethene	CAS:156-60-5	1,2,4-Trichlorobenzene	CAS:120-82-1
Carbon tetrachloride	CAS:56-23-5	1,2-Dichloropropane	CAS:78-87-5	1,1,1-Trichloroethane	CAS:71-55-6
Chlorobenzene	CAS:108-90-7	1,3-Dichloropropane	CAS:142-28-9	1,1,2-Trichloroethane	CAS:79-00-5
Chloroethane	CAS:75-00-3	2,2-Dichloropropane	CAS:594-20-7	Trichloroethene	CAS:79-01-6
Chloroform	CAS:67-66-3	1,1-Dichloropropene	CAS:563-58-6	Trichlorofluoromethane	CAS:75-69-4
Chloromethane	CAS:74-87-3	cis-1,3-Dichloropropene	CAS:10061-01-5	1,2,3-Trichloropropane	CAS:96-18-4
2-Chlorotoluene	CAS:95-49-8	trans-1,3-Dichloropropene	CAS:10061-02-6	1,3,5-Trimethylbenzene	CAS:108-67-8
4-Chlorotoluene	CAS:106-43-4	Ethylbenzene	CAS:100-41-4	1,2,4-Trimethylbenzene	CAS:95-63-6
Dibromochloromethane	CAS:124-48-1	Hexachlorobutadiene	CAS:87-68-3	Vinylchloride	CAS:75-01-4
1,2-Dibromo-3-chloropropane	CAS:96-12-8	Isopropylbenzene	CAS:98-82-8	o-Xylene	CAS:95-47-6
1,2-Dibromoethane	CAS:106-93-4	4-Isopropyltoluene	CAS:99-87-6	m-Xylene	CAS:108-38-3
Dibromomethane	CAS:74-95-3	Methylene chloride	CAS:75-09-2	p-Xylene	CAS:106-42-3

Solvent : Methanol purge & trap

ampoule 1 ml 2 000 µg/ml

Ref : F112751

Price : BGA

WHO/THM Mixture - 6 components

Bromodichloromethane	CAS:75-27-4	1 000 µg/ml	1,2-Dichloroethane	CAS:107-06-2	1 000 µg/ml
Chloroform	CAS:67-66-3	1 000 µg/ml	Trichloroethene	CAS:79-01-6	1 000 µg/ml
Dibromochloromethane	CAS:124-48-1	1 000 µg/ml	Tetrachloroethene	CAS:127-18-4	100 µg/ml

Solvent : Methanol purge & trap

ampoule 1 ml

Ref : F127751

Price : EC

USP test mixtures for GC analysis of organic volatile impurities

USP Test Mixture - 5 components

Benzene	CAS:71-43-2	1 000 µg/ml	Methylene chloride	CAS:75-09-2	1 000 µg/ml
Chloroform	CAS:67-66-3	500 µg/ml	Trichloroethene	CAS:79-01-6	1 000 µg/ml
1,4-Dioxane	CAS:123-91-1	1 000 µg/ml			

Solvent : Dimethylsulphoxide

ampoule 1 ml

Ref : F128881

Price : DJ

Solvent : Methanol purge & trap

ampoule 1 ml

Ref : F128891

Price : DJ

Revised USP 467 Calibration Mixture - 5 components

Benzene	CAS:71-43-2	2 µg/ml	Methylene chloride	CAS:75-09-2	600 µg/ml
Chloroform	CAS:67-66-3	60 µg/ml	Trichloroethene	CAS:79-01-6	80 µg/ml
1,4-Dioxane	CAS:123-91-1	380 µg/ml			

Solvent : Dimethylsulphoxide

ampoule 1 ml

Ref : F128901

Price : DF

Solvent : Methanol purge & trap

ampoule 1 ml

Ref : F128911

Price : DF

Revised USP 467 Calibration Mixture - 5 components

Benzene	CAS:71-43-2	1 000 µg/ml	Methylene chloride	CAS:75-09-2	5 000 µg/ml
Chloroform	CAS:67-66-3	500 µg/ml	Trichloroethene	CAS:79-01-6	1 000 µg/ml
1,4-Dioxane	CAS:123-91-1	1 000 µg/ml			

Solvent : Dimethylsulphoxide

ampoule 1 ml

Ref : F128921

Price : DJ

Solvent : Methanol purge & trap

ampoule 1 ml

Ref : F128931

Price : DJ

Single-component Solutions

ACENAPHTHENE			CAS 83-32-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P802380	BG
	100 µg/ml	1 ml	ref. P802390	BG
cyclohexane	10 µg/ml	10 ml	ref. P802400	BG
	100 µg/ml	1 ml	ref. P802410	BG
methanol	10 µg/ml	10 ml	ref. P802420	BG
	100 µg/ml	1 ml	ref. P802430	BG
ACENAPHTHENE-D10			CAS 15067-26-2	Price
cyclohexane	100 µg/ml	1 ml	ref. P802460	CB
toluene	100 µg/ml	1 ml	ref. P802470	CB
	2 000 µg/ml	1 ml	ref. F062395	CJ
ACENAPHTHYLENE			CAS 208-96-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P802500	BG
	100 µg/ml	1 ml	ref. P802510	BG
cyclohexane	10 µg/ml	10 ml	ref. P802520	BG
	100 µg/ml	1 ml	ref. P802530	BG
methanol	10 µg/ml	10 ml	ref. P802540	BG
	100 µg/ml	1 ml	ref. P802550	BG
ACENAPHTHYLENE-D8			CAS 93951-97-4	Price
cyclohexane	100 µg/ml	1 ml	ref. P802580	DC
n-nonane	100 µg/ml	1 ml	ref. P802590	DC
ACEPHATE			CAS 30560-19-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P802620	BG
	100 µg/ml	1 ml	ref. P802630	BG
cyclohexane	10 µg/ml	10 ml	ref. P802640	BG
	100 µg/ml	1 ml	ref. P802650	BG
ACETALDEHYDE			CAS 75-07-0	Price
ethanol (ethanol-CPA high purity)	5 µg/ml	10 ml	ref. P873820	BG
	10 µg/ml	10 ml	ref. P873830	BG
	20 µg/ml	10 ml	ref. P873840	BG
acetonitrile	1 000 µg/ml	1 ml	ref. P873850	BK
methanol	1 000 µg/ml	1 ml	ref. P873860	BK
water	1 000 µg/ml	5 ml	ref. P873870	CE
ACETALDEHYDE-DNPH			CAS 1019-57-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P835990	CC
	100 µg/ml	1 ml	ref. P836000	CC
ACETOCHLORE			CAS 34256-82-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P802760	BG
	100 µg/ml	1 ml	ref. P802770	BG
iso-octane	10 µg/ml	10 ml	ref. P802780	BG
	100 µg/ml	1 ml	ref. P802790	BG
methanol	10 µg/ml	10 ml	ref. P802800	BG
	100 µg/ml	1 ml	ref. P802810	BG
ACETONE-DNPH			CAS 1567-89-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P802920	BJ
	100 µg/ml	1 ml	ref. P802930	BJ
ACETONITRILE			CAS 75-05-8	Price
methanol	100 µg/ml	1 ml	ref. F116491	BG
ACETOPHENONE			CAS 98-86-2	Price
cyclohexane	100 µg/ml	1 ml	ref. P873890	BG
methanol	100 µg/ml	5 ml	ref. P873880	BG
ACETOPHENONE-DNPH			CAS 1677-87-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P802960	BJ
	100 µg/ml	1 ml	ref. P802970	BJ
ACIFLUORFEN			CAS 50594-66-6	Price
acetone	10 µg/ml	10 ml	ref. P803000	BG
	100 µg/ml	1 ml	ref. P803010	BG
ACLOINIFEN			CAS 74070-46-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P803160	BG
	100 µg/ml	1 ml	ref. P803170	BG
cyclohexane	10 µg/ml	10 ml	ref. P803180	BG
	100 µg/ml	1 ml	ref. P803190	BG
methanol	10 µg/ml	10 ml	ref. P803200	BG
	100 µg/ml	1 ml	ref. P803210	BG

ACRINATHRIN			CAS 101007-06-1	Price
cyclohexane	100 µg/ml	1 ml	ref. P873900	BH
	10 µg/ml	10 ml	ref. P873910	BH
ACROLEIN-DNPH			CAS 888-54-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P803320	CB
	100 µg/ml	1 ml	ref. P803330	CB
ACROLEIN (2-propenal)			CAS 107-02-8	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P803380	BG
	100 µg/ml	1 ml	ref. P803390	BG
water	10 µg/ml	10 ml	ref. P803360	BG
	100 µg/ml	1 ml	ref. P803370	BG
ACRYLAMIDE			CAS 79-06-1	Price
methanol	10 µg/ml	10 ml	ref. P803240	BG
	100 µg/ml	1 ml	ref. P803250	BG
water	1 000 µg/ml	1 ml	ref. F116471	BK
ACRYLIC ACID-METHYL ESTER			CAS 96-33-3	Price
methanol	100 µg/ml	5 ml	ref. P873920	BG
ACRYLONITRILE			CAS 107-13-1	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P803280	BG
	100 µg/ml	1 ml	ref. P803290	BG
	1 000 µg/ml	1 ml	ref. F116461	BK
ADIPIC ACID-BIS-2-ETHYLHEXYL ESTER			CAS 103-23-1	Price
acetone	10 µg/ml	10 ml	ref. P807230	BG
	100 µg/ml	1 ml	ref. P807240	BG
methanol	10 µg/ml	10 ml	ref. P807250	BG
	100 µg/ml	1 ml	ref. P807260	BG
ALACHLOR			CAS 15972-60-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P803450	BG
	100 µg/ml	1 ml	ref. P803460	BG
methanol	10 µg/ml	10 ml	ref. P803470	BG
	100 µg/ml	1 ml	ref. P803480	BG
ALDICARB			CAS 116-06-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P803550	BG
	100 µg/ml	1 ml	ref. P803560	BG
cyclohexane	10 µg/ml	10 ml	ref. P803570	BG
	100 µg/ml	1 ml	ref. P803580	BG
ALDICARB-SULFONE			CAS 1646-88-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P873940	BH
	100 µg/ml	1 ml	ref. P873930	BH
ALDICARB-SULFOXIDE			CAS 1646-87-3	Price
methyl-tert-butylether	10 µg/ml	10 ml	ref. P873960	BH
	100 µg/ml	1 ml	ref. P873950	BH
ALDRIN			CAS 309-00-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P803650	BG
	100 µg/ml	1 ml	ref. P803660	BG
methanol	10 µg/ml	10 ml	ref. P803670	BG
	100 µg/ml	1 ml	ref. P803680	BG
ALLETHRIN			CAS 584-79-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P873980	BG
iso-octane	10 µg/ml	10 ml	ref. P874350	BG
methanol	100 µg/ml	1 ml	ref. P873970	BG
ALLIDOCHLOR			CAS 93-71-0	Price
methanol	100 µg/ml	1 ml	ref. P879160	BG
ALLOXYDIN-SODIUM			CAS 66003-55-2	Price
acetonitrile	100 µg/ml	1 ml	ref. P882540	BG
AMETRYN			CAS 834-12-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P803750	BG
	100 µg/ml	1 ml	ref. P803760	BG
cyclohexane	10 µg/ml	10 ml	ref. P803770	BG
	100 µg/ml	1 ml	ref. P803780	BG
AMIDOSULFURON			CAS 120923-37-7	Price
acetone	10 µg/ml	10 ml	ref. P803810	BH
	100 µg/ml	1 ml	ref. P803820	BH

4-AMINOBIPHENYL			CAS 92-67-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P873470	BH
methanol	100 µg/ml	5 ml	ref. P873460	BH
AMINOCARB			CAS 2032-59-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P879180	BG
methanol	100 µg/ml	1 ml	ref. P879170	BG
2-AMINO-4-CHLOROPHENOL			CAS 95-85-2	Price
methanol	100 µg/ml	1 ml	ref. P872980	BG
2-AMINO-4,6-DINITROTOLUENE			CAS 35572-78-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P803850	CD
	100 µg/ml	1 ml	ref. P803860	CD
4-AMINO-2,6-DINITROTOLUENE			CAS 19406-51-0	Price
acetonitrile	1 000 µg/ml	1 ml	ref. F062405	HJ
	1 000 µg/ml	5 ml	ref. P873450	DBE
methanol	10 µg/ml	1 ml	ref. P873440	CC
AMINOMETHYL PHOSPHONIC ACID			CAS 1066-51-9	Price
water	10 µg/ml	10 ml	ref. P803890	BG
	100 µg/ml	1 ml	ref. P803900	BG
2-AMINO-4-NITROTOLUENE			CAS 99-55-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P832910	BG
	100 µg/ml	1 ml	ref. P832920	BG
methanol	10 µg/ml	10 ml	ref. P832930	BG
	100 µg/ml	1 ml	ref. P832940	BG
4-AMINOPHENYLETHER			CAS 101-80-4	Price
methanol	100 µg/ml	5 ml	ref. P873480	BG
2-AMINONAPHTHALENE			CAS 91-59-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P873010	CK
methanol	10 µg/ml	10 ml	ref. P873000	CK
	100 µg/ml	5 ml	ref. P872990	CK
AMITRAZ			CAS 33089-61-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P874040	BG
	100 µg/ml	1 ml	ref. P874030	BG
cyclohexane	100 µg/ml	1 ml	ref. P874020	BG
AMITROLE			CAS 61-82-5	Price
methanol	10 µg/ml	10 ml	ref. P804000	BG
	100 µg/ml	1 ml	ref. P804010	BG
5α-ANDROSTANE			CAS 438-22-2	Price
dichloromethane	1 000 µg/ml	5 ml	ref. P882210	FA
	2 000 µg/ml	5 ml	ref. P882220	JF
	2 000 µg/ml	1 ml	ref. P882230	DA
ANILINE			CAS 62-53-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P804160	BG
	100 µg/ml	1 ml	ref. P804170	BG
methanol	10 µg/ml	10 ml	ref. P804180	BG
	100 µg/ml	1 ml	ref. P804190	BG
ANILINE-D5			CAS 4165-61-1	Price
acetone	100 µg/ml	1 ml	ref. F146943	BJ
methyl-tert-butylether	100 µg/ml	1 ml	ref. P804230	BJ
	2 000 µg/ml	1 ml	ref. F062415	CH
ANTHRACENE			CAS 120-12-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P804260	BG
	100 µg/ml	1 ml	ref. P804270	BG
cyclohexane	10 µg/ml	10 ml	ref. P804280	BG
	100 µg/ml	1 ml	ref. P804290	BG
dichloromethane	10 µg/ml	10 ml	ref. P804300	BG
	100 µg/ml	1 ml	ref. P804310	BG
ANTHRACENE-D10			CAS 1719-06-8	Price
acetone	100 µg/ml	1 ml	ref. F146771	BJ
cyclohexane	100 µg/ml	1 ml	ref. P804340	BJ
methyl-tert-butylether	100 µg/ml	1 ml	ref. P804350	BJ
ANTHRAQUINONE			CAS 84-65-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P804380	BG
	100 µg/ml	1 ml	ref. P804390	BG
ARAMITE			CAS 140-57-8	Price
n-hexane	10 µg/ml	10 ml	ref. P804420	CK
	100 µg/ml	1 ml	ref. P804430	CK

AROCLOR 1016			CAS 12674-11-2	Price
dichloromethane	100 µg/ml	1 ml	ref. F119731	BH
n-hexane	10 µg/ml	10 ml	ref. P804460	BH
	100 µg/ml	1 ml	ref. F116741	BH
iso-octane	10 µg/ml	10 ml	ref. P804480	BH
	100 µg/ml	1 ml	ref. F115081	BH
	1 000 µg/ml	1 ml	ref. F115681	CA
methanol	100 µg/ml	1 ml	ref. F113211	BH
methyl-tert-butylether	1 000 µg/ml	1 ml	ref. F120061	CA
AROCLOR 1221			CAS 11104-28-2	Price
dichloromethane	100 µg/ml	1 ml	ref. F119741	BG
n-hexane	10 µg/ml	10 ml	ref. P804520	BG
	100 µg/ml	1 ml	ref. F116751	BG
iso-octane	10 µg/ml	10 ml	ref. P804540	BG
	100 µg/ml	1 ml	ref. F115091	BG
	1 000 µg/ml	1 ml	ref. F115691	BK
methanol	100 µg/ml	1 ml	ref. F113221	BG
methyl-tert-butylether	1 000 µg/ml	1 ml	ref. F120071	BK
AROCLOR 1232			CAS 11141-16-5	Price
dichloromethane	100 µg/ml	1 ml	ref. F119751	CC
n-hexane	10 µg/ml	10 ml	ref. P804580	CC
	100 µg/ml	1 ml	ref. F116761	CC
iso-octane	10 µg/ml	10 ml	ref. P804600	CC
	100 µg/ml	1 ml	ref. F115121	CC
	1 000 µg/ml	1 ml	ref. F115701	EF
methanol	100 µg/ml	1 ml	ref. F113231	CC
methyl-tert-butylether	1 000 µg/ml	1 ml	ref. F120081	EF
AROCLOR 1242			CAS 53469-21-9	Price
dichloromethane	100 µg/ml	1 ml	ref. F119761	BH
n-hexane	100 µg/ml	1 ml	ref. F116771	BH
iso-octane	100 µg/ml	1 ml	ref. F115131	BH
	1 000 µg/ml	1 ml	ref. F115711	CA
methanol	100 µg/ml	1 ml	ref. F113241	BH
methyl-tert-butylether	1 000 µg/ml	1 ml	ref. F120091	CA
AROCLOR 1248			CAS 12672-29-6	Price
dichloromethane	100 µg/ml	1 ml	ref. F119771	BJ
n-hexane	10 µg/ml	10 ml	ref. P804700	BJ
	100 µg/ml	1 ml	ref. F116781	BJ
iso-octane	10 µg/ml	10 ml	ref. P804720	BJ
	100 µg/ml	1 ml	ref. F115141	BJ
	1 000 µg/ml	1 ml	ref. F115721	CC
methanol	100 µg/ml	1 ml	ref. F113251	BJ
methyl-tert-butylether	1 000 µg/ml	1 ml	ref. F120101	CC
AROCLOR 1254			CAS 11097-69-1	Price
dichloromethane	100 µg/ml	1 ml	ref. F119781	BJ
n-hexane	10 µg/ml	10 ml	ref. P804760	BJ
	100 µg/ml	1 ml	ref. F116791	BJ
iso-octane	10 µg/ml	10 ml	ref. P804780	BJ
	100 µg/ml	1 ml	ref. F115151	BJ
	1 000 µg/ml	1 ml	ref. F115731	CC
methanol	100 µg/ml	1 ml	ref. F113261	BJ
methyl-tert-butylether	1 000 µg/ml	1 ml	ref. F120111	CC
AROCLOR 1260			CAS 11096-82-5	Price
dichloromethane	100 µg/ml	1 ml	ref. F119791	BJ
n-hexane	10 µg/ml	10 ml	ref. P804820	BJ
	100 µg/ml	1 ml	ref. F116801	BJ
iso-octane	10 µg/ml	10 ml	ref. P804840	BJ
	100 µg/ml	1 ml	ref. F115161	BJ
	1 000 µg/ml	1 ml	ref. F115741	CC
methanol	100 µg/ml	1 ml	ref. F113271	BJ
	1 000 µg/ml	1 ml	ref. F113631	CC
	5 000 µg/ml	1 ml	ref. F113641	CH
methyl-tert-butylether	1 000 µg/ml	1 ml	ref. F120121	CC
AROCLOR 1262			CAS 37324-23-5	Price
n-hexane	10 µg/ml	10 ml	ref. P804880	BG
	100 µg/ml	1 ml	ref. P804890	BG
AROCLOR 1268			CAS 11100-14-4	Price
hexane	10 µg/ml	10 ml	ref. P804920	BG
	100 µg/ml	1 ml	ref. P804930	BG
ASPON			CAS 3244-90-4	Price
methanol	100 µg/ml	1 ml	ref. P879190	BG

ASULAM			CAS 3337-71-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P874060	BG
methanol	100 µg/ml	5 ml	ref. P874050	BG
ATRATON			CAS 1610-17-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P874080	BH
methyl- <i>tert</i> -butylether	100 µg/ml	1 ml	ref. P874070	BH
ATRAZINE			CAS 1912-24-9	Price
acetone	10 µg/ml	10 ml	ref. P804960	BG
	100 µg/ml	1 ml	ref. P804970	BG
acetonitrile	10 µg/ml	10 ml	ref. P804980	BG
	100 µg/ml	1 ml	ref. P804990	BG
cyclohexane	10 µg/ml	10 ml	ref. P805000	BG
	100 µg/ml	1 ml	ref. P805010	BG
ATRAZINE-D5			CAS 163165-75-1	Price
acetone	100 µg/ml	1 ml	ref. P874100	BHG
acetonitrile	1 000 µg/ml	1 ml	ref. P874090	CBA
n-nonane	100 µg/ml	1 ml	ref. P805050	BHG
ATRAZINE-DESETHYL			CAS 6190-65-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P805200	BG
	100 µg/ml	1 ml	ref. P805210	BG
cyclohexane	10 µg/ml	10 ml	ref. P805220	BG
	100 µg/ml	1 ml	ref. P805230	BG
ATRAZINE-DESETHYL-DESIISOPROPYL			CAS 3397-62-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P805380	BG
	100 µg/ml	1 ml	ref. P805390	BG
cyclohexane	10 µg/ml	10 ml	ref. P805400	BG
	100 µg/ml	1 ml	ref. P805410	BG
ATRAZINE-DESIISOPROPYL			CAS 1007-28-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P805290	BG
	100 µg/ml	1 ml	ref. P805300	BG
cyclohexane	10 µg/ml	10 ml	ref. P805310	BG
	100 µg/ml	1 ml	ref. P805320	BG
ATRAZINE-2-HYDROXY			CAS 2163-68-0	Price
methanol	10 µg/ml	10 ml	ref. P805440	BG
	100 µg/ml	1 ml	ref. P805450	BG
AZACONAZOLE			CAS 60207-31-0	Price
acetone	10 µg/ml	10 ml	ref. P805520	BG
	100 µg/ml	1 ml	ref. P805530	BG
acetonitrile	10 µg/ml	10 ml	ref. P805540	BG
	100 µg/ml	1 ml	ref. P805550	BG
cyclohexane	10 µg/ml	10 ml	ref. P805560	BG
	100 µg/ml	1 ml	ref. P805570	BG
AZAMETHIPHOS			CAS 35575-96-3	Price
acetone	10 µg/ml	10 ml	ref. P874120	BH
dichloromethane	100 µg/ml	5 ml	ref. P874110	BH
AZINPHOS-ETHYL			CAS 2642-71-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P805600	BG
	100 µg/ml	1 ml	ref. P805610	BG
cyclohexane	10 µg/ml	10 ml	ref. P805620	BG
	100 µg/ml	1 ml	ref. P805630	BG
iso-octane	10 µg/ml	10 ml	ref. P805640	BG
	100 µg/ml	1 ml	ref. P805650	BG
AZINPHOS-METHYL			CAS 86-50-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P805680	BG
	100 µg/ml	1 ml	ref. P805690	BG
cyclohexane	10 µg/ml	10 ml	ref. P805700	BG
	100 µg/ml	1 ml	ref. P805710	BG
AZIPROTRYNE			CAS 4658-28-0	Price
acetonitrile	100 µg/ml	1 ml	ref. P879200	BG
AZOBENZENE			CAS 103-33-3	Price
iso-octane	10 µg/ml	10 ml	ref. P805780	BG
	100 µg/ml	1 ml	ref. P805790	BG
methanol	10 µg/ml	10 ml	ref. P805800	BG
	100 µg/ml	1 ml	ref. P872670	BG
AZOXYSTROBIN			CAS 131860-33-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P805840	BG
	100 µg/ml	1 ml	ref. P805850	BG
BARBAN			CAS 101-27-9	Price
acetonitrile	100 µg/ml	1 ml	ref. P879210	BG

BENALAXYL			CAS 71626-11-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P874140	BG
	100 µg/ml	1 ml	ref. P874130	BG
BENAZOLIN			CAS 3813-05-6	Price
methanol	100 µg/ml	1 ml	ref. P879220	BG
BENDIOCARB			CAS 22781-23-3	Price
acetonitrile	100 µg/ml	1 ml	ref. P874150	BG
iso-octane	10 µg/ml	10 ml	ref. P874160	BG
methanol	1000 µg/ml	1 ml	ref. F116001	BK
BENFLURALIN			CAS 1861-40-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P805920	BG
	100 µg/ml	1 ml	ref. P805930	BG
BENFURACARB			CAS 82560-54-1	Price
acetonitrile	100 µg/ml	1 ml	ref. P874170	BG
cyclohexane	10 µg/ml	10 ml	ref. P874180	BG
BENODANIL			CAS 15310-01-7	Price
methanol	10 µg/ml	10 ml	ref. P879240	BG
	100 µg/ml	1 ml	ref. P879230	BG
BENOMYL			CAS 17804-35-2	Price
methanol	10 µg/ml	10 ml	ref. P806090	BG
	100 µg/ml	1 ml	ref. P806100	BG
BENOXACOR			CAS 98730-04-2	Price
acetone	10 µg/ml	10 ml	ref. P806130	BG
	100 µg/ml	1 ml	ref. P806140	BG
BENSULFURON-METHYL			CAS 83055-99-6	Price
methanol	100 µg/ml	1 ml	ref. P874190	BG
BENSULIDE			CAS 741-58-2	Price
methanol	10 µg/ml	10 ml	ref. P879260	BG
	100 µg/ml	1 ml	ref. P879250	BG
	1 000 µg/ml	1 ml	ref. F115981	BK
BENSULTAP			CAS 17606-31-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P874210	BG
methanol	100 µg/ml	1 ml	ref. P874200	BG
BENTAZON			CAS 25057-89-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P806000	BG
	100 µg/ml	1 ml	ref. P806010	BG
	1 000 µg/ml	1 ml	ref. F116041	BK
methanol	10 µg/ml	10 ml	ref. P806020	BG
	100 µg/ml	1 ml	ref. P806030	BG
BENZALDEHYDE			CAS 100-52-7	Price
acetonitrile	1 000 µg/ml	5 ml	ref. P874220	CE
BENZALDEHYDE-DNPH			CAS 1157-84-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P806170	BG
	100 µg/ml	1 ml	ref. P806180	BG
BENZENE			CAS 71-43-2	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P806210	BG
	100 µg/ml	1 ml	ref. P874230	BG
BENZENE-D6			CAS 1076-43-3	Price
methanol purge & trap	100 µg/ml	1 ml	ref. P806250	BJ
	2 000 µg/ml	1 ml	ref. F062425	CH
4,4'-BENZIDINE			CAS 92-87-5	Price
methanol	10 µg/ml	10 ml	ref. P806280	BG
	100 µg/ml	1 ml	ref. P806290	BG
BENZO(a)ANTHRACENE			CAS 56-55-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P806320	BH
	100 µg/ml	1 ml	ref. P806330	BH
cyclohexane	10 µg/ml	10 ml	ref. P806340	BH
	100 µg/ml	1 ml	ref. P806350	BH
methanol	10 µg/ml	10 ml	ref. P806360	BH
	100 µg/ml	1 ml	ref. P806370	BH
BENZO(a)ANTHRACENE -D12			CAS 1718-53-2	Price
acetone	100 µg/ml	1 ml	ref. P806400	DD
cyclohexane	100 µg/ml	1 ml	ref. P806410	DD
methyl- <i>tert</i> -butylether	100 µg/ml	1 ml	ref. P806420	DD
BENZO(b)CHRYSENE			CAS 214-17-5	Price
acetonitrile	100 µg/ml	1 ml	ref. P874260	CG
cyclohexane	10 µg/ml	10 ml	ref. P874270	CG
toluene	100 µg/ml	1 ml	ref. P874250	CG

BENZO(b)FLUORANTHENE			CAS 205-99-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P806450	<i>BJ</i>
	100 µg/ml	1 ml	ref. P806460	<i>BJ</i>
cyclohexane	10 µg/ml	10 ml	ref. P806470	<i>BJ</i>
	100 µg/ml	1 ml	ref. P806480	<i>BJ</i>
BENZO(b)FLUORANTHENE-D12			CAS 93951-98-5	Price
acetonitrile	100 µg/ml	1 ml	ref. P806510	<i>DD</i>
cyclohexane	100 µg/ml	1 ml	ref. P806520	<i>DD</i>
BENZO(j)FLUORANTHENE			CAS 205-82-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P874290	<i>CC</i>
cyclohexane	10 µg/ml	10 ml	ref. P874300	<i>CC</i>
methanol	100 µg/ml	5 ml	ref. P874280	<i>EE</i>
BENZO(k)FLUORANTHENE			CAS 207-08-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P806550	<i>CA</i>
	100 µg/ml	1 ml	ref. P806560	<i>CA</i>
cyclohexane	10 µg/ml	10 ml	ref. P806570	<i>CA</i>
	100 µg/ml	1 ml	ref. P806580	<i>CA</i>
BENZOIC ACID			CAS 65-85-0	Price
methanol	10 µg/ml	10 ml	ref. P806640	<i>BG</i>
	100 µg/ml	1 ml	ref. P806650	<i>BG</i>
BENZO(g,h,i)PERYLENE			CAS 191-24-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P806680	<i>BJ</i>
	100 µg/ml	1 ml	ref. P806690	<i>BJ</i>
cyclohexane	10 µg/ml	10 ml	ref. P806700	<i>BJ</i>
	100 µg/ml	1 ml	ref. P806710	<i>BJ</i>
dichloromethane	10 µg/ml	10 ml	ref. P806720	<i>BJ</i>
	100 µg/ml	1 ml	ref. P806730	<i>BJ</i>
BENZO(g,h,i)PERYLENE-D12			CAS 93951-66-7	Price
cyclohexane	100 µg/ml	1 ml	ref. P806760	<i>GH</i>
BENZO(a)PYRENE			CAS 50-32-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P806790	<i>BJ</i>
	100 µg/ml	1 ml	ref. P806800	<i>BJ</i>
cyclohexane	10 µg/ml	10 ml	ref. P806810	<i>BJ</i>
	100 µg/ml	1 ml	ref. P806820	<i>BJ</i>
BENZO(a)PYRENE-D12			CAS 63466-71-7	Price
acetone	100 µg/ml	1 ml	ref. P806850	<i>BAB</i>
cyclohexane	100 µg/ml	1 ml	ref. P806860	<i>BAB</i>
BENZO(e)PYRENE			CAS 192-97-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P806890	<i>CA</i>
	100 µg/ml	1 ml	ref. P806900	<i>CA</i>
cyclohexane	10 µg/ml	10 ml	ref. P806910	<i>CA</i>
	100 µg/ml	1 ml	ref. P806920	<i>CA</i>
BENZOXIMATE			CAS 29104-30-1	Price
cyclohexane	100 µg/ml	1 ml	ref. P879270	<i>BG</i>
BENZTHIAZURON			CAS 1929-88-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P879290	<i>BG</i>
dichloromethane	100 µg/ml	1 ml	ref. P879280	<i>BG</i>
	100 µg/ml	5 ml	ref. P882890	<i>BG</i>
BENZYL BENZOATE			CAS 120-51-4	Price
n-hexane	10 µg/ml	10 ml	ref. P806950	<i>BG</i>
	100 µg/ml	1 ml	ref. P806960	<i>BG</i>
	5 000 µg/ml	1 ml	ref. F116661	<i>CE</i>
BENZYL SALICYLATE			CAS 118-58-1	Price
acetonitrile	1 000 µg/ml	1 ml	ref. P874310	<i>BH</i>
BIFENOX			CAS 42576-02-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P874330	<i>BG</i>
cyclohexane	10 µg/ml	10 ml	ref. P874340	<i>BG</i>
methanol	100 µg/ml	1 ml	ref. P874320	<i>BG</i>
BIFENTHRIN			CAS 82657-04-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P806990	<i>BG</i>
	100 µg/ml	1 ml	ref. P807000	<i>BG</i>
cyclohexane	10 µg/ml	10 ml	ref. P807010	<i>BG</i>
	100 µg/ml	1 ml	ref. P807020	<i>BG</i>
iso-octane	10 µg/ml	10 ml	ref. P807030	<i>BG</i>
	100 µg/ml	1 ml	ref. P807040	<i>BG</i>
methanol	10 µg/ml	10 ml	ref. P807050	<i>BG</i>
	100 µg/ml	1 ml	ref. P807060	<i>BG</i>
BIORESMETHRIN			CAS 28434-01-7	Price
iso-octane	10 µg/ml	10 ml	ref. P879310	<i>BG</i>
methanol	100 µg/ml	1 ml	ref. P879300	<i>BG</i>

BIPHENYL			CAS 92-52-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P807130	<i>BG</i>
	100 µg/ml	1 ml	ref. P807140	<i>BG</i>
n-hexane	1 000 µg/ml	1 ml	ref. F116021	<i>BK</i>
methanol	10 µg/ml	10 ml	ref. P807150	<i>BG</i>
	100 µg/ml	1 ml	ref. P807160	<i>BG</i>
BIS-(4-AMINOPHENYL)-METHANE			CAS 101-77-9	Price
acetonitrile	100 µg/ml	1 ml	ref. P874430	<i>BG</i>
toluene	100 µg/ml	1 ml	ref. P874420	<i>BG</i>
BIS-(2-CHLOROETHOXY)-METHANE			CAS 111-91-1	Price
dichloromethane	100 µg/ml	1 ml	ref. P874390	<i>BG</i>
BIS-(2-CHLOROETHYL)-ETHER			CAS 111-44-4	Price
acetone	100 µg/ml	5 ml	ref. P874400	<i>BG</i>
BIS-(2-CHLORO-1-METHYLETHYL)-ETHER			CAS 108-60-1	Price
n-hexane	100 µg/ml	5 ml	ref. P874360	<i>BH</i>
dichloromethane	100 µg/ml	1 ml	ref. P874370	<i>BH</i>
BIS(2-ETHYLHEXYL)PHTHALATE-D4			id : 21280	Price
acetone	100 µg/ml	1 ml	ref. P807330	<i>BCG</i>
nonane	100 µg/ml	1 ml	ref. P807340	<i>BCG</i>
BISPHENOL A			CAS 80-05-7	Price
methanol	1 000 µg/ml	1 ml	ref. P874440	<i>BH</i>
BITERTANOL			CAS 55179-31-2	Price
acetone	10 µg/ml	10 ml	ref. P807370	<i>BG</i>
	100 µg/ml	1 ml	ref. P807380	<i>BG</i>
acetonitrile	10 µg/ml	10 ml	ref. P807390	<i>BG</i>
	100 µg/ml	1 ml	ref. P807400	<i>BG</i>
cyclohexane	10 µg/ml	10 ml	ref. P807410	<i>BG</i>
	100 µg/ml	1 ml	ref. P807420	<i>BG</i>
BRODIFACOU			CAS 56073-10-0	Price
methanol	100 µg/ml	1 ml	ref. P874450	<i>BG</i>
BROMACIL			CAS 314-40-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P807450	<i>BG</i>
	100 µg/ml	1 ml	ref. P807460	<i>BG</i>
methanol	10 µg/ml	10 ml	ref. P807470	<i>BG</i>
	100 µg/ml	1 ml	ref. P807480	<i>BG</i>
BROMADIOLONE			CAS 28772-56-7	Price
methanol	100 µg/ml	1 ml	ref. P874470	<i>BG</i>
4-BROMOANILINE			CAS 106-40-1	Price
iso-octane	100 µg/ml	1 ml	ref. P873490	<i>BG</i>
toluene	100 µg/ml	5 ml	ref. P873500	<i>BG</i>
BROMOBENZENE			CAS 108-86-1	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P807510	<i>BG</i>
	100 µg/ml	1 ml	ref. P807520	<i>BG</i>
2-BROMOCHLOROBENZENE			CAS 694-80-4	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P807550	<i>BG</i>
	100 µg/ml	1 ml	ref. P807560	<i>BG</i>
4-BROMOCHLOROBENZENE			CAS 106-39-8	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P807590	<i>BG</i>
	100 µg/ml	1 ml	ref. P807600	<i>BG</i>
1-BROMO-2-CHLOROETHANE			CAS 107-04-0	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P807630	<i>BG</i>
	100 µg/ml	1 ml	ref. P807640	<i>BG</i>
BROMOCHLOROMETHANE			CAS 74-97-5	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P807670	<i>BG</i>
	100 µg/ml	1 ml	ref. P807680	<i>BG</i>
2-BROMO-1-CHLOROPROPANE			CAS 3017-95-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P807750	<i>BG</i>
	100 µg/ml	1 ml	ref. P807760	<i>BG</i>
	2 000 µg/ml	1 ml	ref. F112931	<i>DA</i>
BROMODICHLOROMETHANE			CAS 75-27-4	Price
methanol	10 µg/ml	10 ml	ref. P807710	<i>BG</i>
	100 µg/ml	1 ml	ref. P807720	<i>BG</i>
4-BROMO-3,5-DIMETHYLPHENYL-N-METHYLCARBAMATE (BDMC)			CAS 672-99-1	Price
methanol	10 µg/ml	10 ml	ref. P807800	<i>BG</i>
	100 µg/ml	1 ml	ref. P807810	<i>BG</i>
	100 µg/ml	10 ml	ref. F114171	<i>BG</i>

4-BROMO-FLUOROBENZENE			CAS 460-00-4	Price
acetone	10 µg/ml	10 ml	ref. P807840	BG
	100 µg/ml	1 ml	ref. P807850	BG
methanol purge & trap	10 µg/ml	10 ml	ref. P807860	BG
	100 µg/ml	1 ml	ref. P807870	BG
	2 000 µg/ml	1 ml	ref. F113921	CA
BROMOFORM (tribromomethane)			CAS 75-25-2	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P807900	BG
	100 µg/ml	1 ml	ref. P807910	BG
BROMOMETHANE			CAS 74-83-9	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P807940	CE
	100 µg/ml	1 ml	ref. P807950	CE
1-BROMO-2-NITROBENZENE			CAS 577-19-5	Price
acetone	10 µg/ml	10 ml	ref. P807980	BG
	100 µg/ml	1 ml	ref. P807990	BG
	1 000 µg/ml	1 ml	ref. F119061	BK
	5 000 µg/ml	1 ml	ref. F117261	CE
BROMOPHOS-ETHYL			CAS 4824-78-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P808040	BG
	100 µg/ml	1 ml	ref. P808050	BG
iso-octane	10 µg/ml	10 ml	ref. P808020	BG
	100 µg/ml	1 ml	ref. P808030	BG
methanol	10 µg/ml	10 ml	ref. P808060	BG
	100 µg/ml	1 ml	ref. P808070	BG
BROMOPHOS-METHYL			CAS 2104-96-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P808140	BG
	100 µg/ml	1 ml	ref. P808150	BG
cyclohexane	10 µg/ml	10 ml	ref. P808160	BG
	100 µg/ml	1 ml	ref. P808170	BG
2-BROMOPROPIONIC ACID			CAS 598-72-1	Price
methyl-tert-butylether	1 000 µg/ml	1 ml	ref. P881560	BH
	2 000 µg/ml	5 ml	ref. P881550	CA
BROMOPROPYLATE			CAS 18181-80-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P808200	BG
	100 µg/ml	1 ml	ref. P808210	BG
cyclohexane	10 µg/ml	10 ml	ref. P808220	BG
	100 µg/ml	1 ml	ref. P808230	BG
methanol	10 µg/ml	10 ml	ref. P808240	BG
	100 µg/ml	1 ml	ref. P808250	BG
BROMOXYNIL			CAS 1689-84-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P808360	BG
	100 µg/ml	1 ml	ref. P808370	BG
iso-octane	10 µg/ml	10 ml	ref. P808380	BG
	100 µg/ml	1 ml	ref. P808390	BG
BROMOXYNIL-OCTANOATE			CAS 1689-99-2	Price
methanol	10 µg/ml	10 ml	ref. P808420	BG
	100 µg/ml	1 ml	ref. P808430	BG
BROMUCONAZOLE			CAS 116255-48-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P808460	BH
	100 µg/ml	1 ml	ref. P808470	BH
cyclohexane	10 µg/ml	10 ml	ref. P808480	BH
	100 µg/ml	1 ml	ref. P808490	BH
methanol	10 µg/ml	10 ml	ref. P808500	BH
	100 µg/ml	1 ml	ref. P808510	BH
BUPIRIMATE			CAS 41483-43-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P808540	BG
	100 µg/ml	1 ml	ref. P808550	BG
cyclohexane	10 µg/ml	10 ml	ref. P808560	BG
	100 µg/ml	1 ml	ref. P808570	BG
methanol	10 µg/ml	10 ml	ref. P808580	BG
	100 µg/ml	1 ml	ref. P808590	BG
BUPROFEZIN			CAS 69327-76-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P878490	BG
	100 µg/ml	1 ml	ref. P874480	BG
BUTACHLOR			CAS 23184-66-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P808620	BG
	100 µg/ml	1 ml	ref. P808630	BG
1,3-BUTADIENE			CAS 106-99-0	Price
methanol	200 µg/ml	1 ml	ref. P872230	BG
	2 000 µg/ml	1 ml	ref. P872240	BK

BUTALBITAL			CAS 77-26-9	Price
methanol	10 µg/ml	10 ml	ref. P879330	BH
	1 000 µg/ml	1 ml	ref. P879320	CA
1-BUTANOL			CAS 71-36-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P808660	BG
	100 µg/ml	1 ml	ref. P808670	BG
2-BUTANOL			CAS 78-92-2	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P808700	BG
	100 µg/ml	1 ml	ref. P808710	BG
tert-BUTANOL			CAS 75-65-0	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P808740	BG
	100 µg/ml	1 ml	ref. P808750	BG
2-BUTANONE			CAS 78-93-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P808780	BG
	100 µg/ml	1 ml	ref. P808790	BG
2-BUTANONE-DNPH			CAS 958-60-1	Price
methanol	10 µg/ml	10 ml	ref. P808820	BG
	100 µg/ml	1 ml	ref. P808830	BG
BUTCARBOXIM			CAS 34681-10-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P879350	BG
	100 µg/ml	1 ml	ref. P879340	BG
BUTRALIN			CAS 33629-47-9	Price
iso-octane	10 µg/ml	10 ml	ref. P808860	BG
	100 µg/ml	1 ml	ref. P808870	BG
BUTURON			CAS 3766-60-7	Price
acetone	10 µg/ml	10 ml	ref. P808900	BG
	100 µg/ml	1 ml	ref. P808910	BG
cyclohexane	10 µg/ml	10 ml	ref. P808920	BG
	100 µg/ml	1 ml	ref. P808930	BG
BUTYLATE			CAS 2008-41-5	Price
methanol	100 µg/ml	1 ml	ref. P874510	BG
	100 µg/ml	5 ml	ref. P874500	BG
n-BUTYLBENZENE			CAS 104-51-8	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P808960	BG
	100 µg/ml	1 ml	ref. P808970	BG
sec-BUTYLBENZENE			CAS 135-98-8	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P809000	BG
	100 µg/ml	1 ml	ref. P809010	BG
tert-BUTYLBENZENE			CAS 98-06-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P809040	BG
	100 µg/ml	1 ml	ref. P809050	BG
tert-BUTYL-4-HYDROXYANISOLE (mixture of 2-and 3-i)			CAS 25013-16-5	Price
methanol	100 µg/ml	1 ml	ref. P878720	BG
	100 µg/ml	1 ml	ref. P878720	BG
n-BUTYL MERCAPTAN			CAS 109-79-5	Price
toluene	2 000 µg/ml	1 ml	ref. P877100	BH
	2 000 µg/ml	1 ml	ref. P877100	BH
4-tert-BUTYL PHENOL			CAS 98-54-4	Price
ethanol	1 000 µg/ml	1 ml	ref. P873660	BH
	1 000 µg/ml	1 ml	ref. P873660	BH
BUTYRALDEHYDE			CAS 123-72-8	Price
acetonitrile	1 000 µg/ml	5 ml	ref. P874520	CF
	100 µg/ml	1 ml	ref. P874530	BH
n-BUTYRALDEHYDE-DNPH			CAS 1527-98-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P809080	BG
	100 µg/ml	1 ml	ref. P809090	BG
ε-CAPROLACTAM			CAS 105-60-2	Price
dichloromethane	100 µg/ml	10 ml	ref. P809240	BG
	100 µg/ml	1 ml	ref. P809250	BG
CAPTAFOL			CAS 2425-06-1	Price
cyclohexane	100 µg/ml	1 ml	ref. P874540	BG
	10 µg/ml	10 ml	ref. P874550	BG
CAPTAN			CAS 133-06-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P809320	BG
	100 µg/ml	1 ml	ref. P809330	BG
cyclohexane	10 µg/ml	10 ml	ref. P809340	BG
	100 µg/ml	1 ml	ref. P809350	BG
CARBARYL			CAS 63-25-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P809410	BG
	100 µg/ml	1 ml	ref. P809420	BG

CARBAZOLE			CAS 86-74-8	Price
iso-octane	10 µg/ml	10 ml	ref. P809520	BG
	100 µg/ml	1 ml	ref. P809530	BG
methanol	10 µg/ml	10 ml	ref. P809540	BG
	100 µg/ml	1 ml	ref. P832400	BG
	1000 µg/ml	1 ml	ref. F115961	BK
dichloromethane	10 µg/ml	10 ml	ref. P809560	BG
	100 µg/ml	1 ml	ref. P809570	BG
CARBENDAZIM			CAS 10605-21-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P809600	BG
	100 µg/ml	1 ml	ref. P809610	BG
methanol	10 µg/ml	10 ml	ref. P809620	BG
	100 µg/ml	1 ml	ref. P809630	BG
CARBETAMIDE			CAS 16118-49-3	Price
acetone	100 µg/ml	1 ml	ref. P874560	BG
acetonitrile	10 µg/ml	10 ml	ref. P874570	BG
CARBOFURAN			CAS 1563-66-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P809790	BG
	100 µg/ml	1 ml	ref. P809800	BG
cyclohexane	10 µg/ml	10 ml	ref. P809810	BG
	100 µg/ml	1 ml	ref. P809820	BG
methanol	10 µg/ml	10 ml	ref. P809830	BG
	100 µg/ml	1 ml	ref. P809840	BG
CARBOFURAN D3			id : 7530	Price
acetone	100 µg/ml	1 ml	ref. P874580	GC
CARBOFURAN-3-HYDROXY			CAS 16655-82-6	Price
ethyl acetate	10 µg/ml	10 ml	ref. P874600	CC
	100 µg/ml	1 ml	ref. P874590	CC
CARBON DISULFIDE			CAS 75-15-0	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P809910	BG
	100 µg/ml	1 ml	ref. P809920	BG
CARBOPHENOTHION			CAS 786-19-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P810000	BG
	100 µg/ml	1 ml	ref. P810010	BG
iso-octane	10 µg/ml	10 ml	ref. P810020	BG
	100 µg/ml	1 ml	ref. P810030	BG
CARBOSULFAN			CAS 55285-14-8	Price
iso-octane	10 µg/ml	10 ml	ref. P874620	BG
	100 µg/ml	1 ml	ref. P874610	BG
CARBOXIN			CAS 5234-68-4	Price
acetone	10 µg/ml	10 ml	ref. P874650	BG
cyclohexane	10 µg/ml	10 ml	ref. P874640	BG
methanol	100 µg/ml	1 ml	ref. P874630	BG
CARFENTRAZONE-ETHYL			CAS 128639-02-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P874660	BH
cyclohexane	10 µg/ml	10 ml	ref. P874670	BH
CHIMETHIONAT			CAS 2439-01-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P874690	BG
	100 µg/ml	1 ml	ref. P874680	BG
CHLORALHYDRATE			CAS 302-17-0	Price
acetonitrile	1 000 µg/ml	1 ml	ref. P879360	BK
methanol	1 000 µg/ml	1 ml	ref. P879370	BK
CHLORAMBEN			CAS 133-90-4	Price
methanol	100 µg/ml	1 ml	ref. P874700	BG
CHLORBENSIDE			CAS 103-17-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P879440	BJ
methanol	100 µg/ml	1 ml	ref. P879430	BJ
	1 000 µg/ml	1 ml	ref. P879420	CC
CHLORBROMURON			CAS 13360-45-7	Price
acetone	10 µg/ml	10 ml	ref. P810130	BG
	100 µg/ml	1 ml	ref. P810140	BG
acetonitrile	10 µg/ml	10 ml	ref. P810150	BG
	100 µg/ml	1 ml	ref. P810160	BG
CHLORBUFAM			CAS 1967-16-4	Price
ethyl acetate	10 µg/ml	10 ml	ref. P879460	BG
methanol	100 µg/ml	1 ml	ref. P879450	BG

CHLORDANE (TECHNICAL)			CAS 57-74-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P810190	BG
	100 µg/ml	1 ml	ref. P810200	BG
dichloromethane	100 µg/ml	1 ml	ref. F119801	BG
n-hexane	10 µg/ml	10 ml	ref. P810210	BG
	100 µg/ml	1 ml	ref. F116721	BG
iso-octane	10 µg/ml	10 ml	ref. P810230	BG
	100 µg/ml	1 ml	ref. F115071	BG
	1 000 µg/ml	1 ml	ref. F115661	BK
methanol	10 µg/ml	10 ml	ref. P810250	BG
	100 µg/ml	1 ml	ref. F113191	BG
methyl-tert-butylether	2 500 µg/ml	1 ml	ref. F120051	CA
cis-CHLORDANE			CAS 5103-71-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P810290	BK
	100 µg/ml	1 ml	ref. P810300	BK
cyclohexane	10 µg/ml	10 ml	ref. P810310	BK
	100 µg/ml	1 ml	ref. P810320	BK
methanol	10 µg/ml	10 ml	ref. P810330	BK
	100 µg/ml	1 ml	ref. P810340	BK
oxy-CHLORDANE			CAS 27304-13-8	Price
acetone	10 µg/ml	1 ml	ref. P877510	CH
cyclohexane	10 µg/ml	1 ml	ref. P877520	CH
	10 µg/ml	10 ml	ref. P877540	CH
	100 µg/ml	1 ml	ref. P877530	EE
methanol	10 µg/ml	1 ml	ref. P877560	CH
	100 µg/ml	1 ml	ref. P877550	EE
trans-CHLORDANE			CAS 5103-74-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P810410	BH
	100 µg/ml	1 ml	ref. P810420	BH
cyclohexane	10 µg/ml	10 ml	ref. P810430	BH
	100 µg/ml	1 ml	ref. P810440	BH
methanol	10 µg/ml	10 ml	ref. P810450	BH
	100 µg/ml	1 ml	ref. P810460	BH
CHLORDECONE HYDRATE			CAS 143-50-0	Price
iso-octane	10 µg/ml	10 ml	ref. P874720	BG
methanol	100 µg/ml	1 ml	ref. P874710	BG
CHLORDENE			CAS 3734-48-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P879480	BK
methanol	100 µg/ml	1 ml	ref. P879470	BK
CHLORDIMEFORM FREE BASE			CAS 6164-98-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P879500	BG
methanol	100 µg/ml	1 ml	ref. P879490	BG
CHLORENDIC ACID DIBUTYL ESTER			CAS 1770-80-5	Price
acetone	10 µg/ml	10 ml	ref. P817540	BH
	100 µg/ml	1 ml	ref. P817550	BH
acetonitrile	10 µg/ml	10 ml	ref. P817560	BH
	100 µg/ml	1 ml	ref. P817570	BH
cyclohexane	10 µg/ml	10 ml	ref. P817580	BH
	100 µg/ml	1 ml	ref. P817590	BH
CHLORFENAPYR			CAS 122453-73-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P879520	BG
methanol	100 µg/ml	1 ml	ref. P879510	BG
CHLORFENPROP-METHYL			CAS 14437-17-3	Price
iso-octane	10 µg/ml	10 ml	ref. P879530	BG
CHLORFENSON			CAS 80-33-1	Price
iso-octane	10 µg/ml	10 ml	ref. P874740	BG
methanol	100 µg/ml	1 ml	ref. P874730	BG
CHLORFENVINPHOS			CAS 470-90-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P810490	BG
	100 µg/ml	1 ml	ref. P810500	BG
cyclohexane	10 µg/ml	10 ml	ref. P810510	BG
	100 µg/ml	1 ml	ref. P810520	BG
CHLORFLUAZURON			CAS 71422-67-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P874750	BG
CHLORIDAZON			CAS 1698-60-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P810580	BG
	100 µg/ml	1 ml	ref. P810590	BG
CHLORIMURON-ETHYL			CAS 90982-32-4	Price
acetonitrile	100 µg/ml	1 ml	ref. P879540	BH
cyclohexane	10 µg/ml	10 ml	ref. P879550	BH

CHLORMEPHOS			CAS 24934-91-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P810620	BG
	100 µg/ml	1 ml	ref. P810630	BG
CHLORMEQUAT CHLORIDE			CAS 999-81-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P810660	BG
	100 µg/ml	1 ml	ref. P810670	BG
methanol	10 µg/ml	10 ml	ref. P810680	BG
	100 µg/ml	1 ml	ref. P810690	BG
CHLOROACETONITRILE			CAS 107-14-2	Price
methyl- <i>tert</i> -butylether	10 µg/ml	10 ml	ref. P874460	BG
	100 µg/ml	1 ml	ref. P810770	BG
2-CHLOROANILINE			CAS 95-51-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P810800	BG
	100 µg/ml	1 ml	ref. P810810	BG
iso-octane	10 µg/ml	10 ml	ref. P810820	BG
	100 µg/ml	1 ml	ref. P810830	BG
3-CHLOROANILINE			CAS 108-42-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P810860	BG
	100 µg/ml	1 ml	ref. P810870	BG
iso-octane	10 µg/ml	10 ml	ref. P810880	BG
	100 µg/ml	1 ml	ref. P810890	BG
methanol	10 µg/ml	10 ml	ref. P810900	BG
	100 µg/ml	1 ml	ref. P810910	BG
4-CHLOROANILINE			CAS 106-47-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P810940	BG
	100 µg/ml	1 ml	ref. P810950	BG
cyclohexane	10 µg/ml	10 ml	ref. P810960	BG
	100 µg/ml	1 ml	ref. P810970	BG
dichloromethane	10 µg/ml	10 ml	ref. P810980	BG
	100 µg/ml	1 ml	ref. P810990	BG
CHLOROBENZENE			CAS 108-90-7	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P811020	BG
	100 µg/ml	1 ml	ref. P811030	BG
CHLOROBENZENE [D5]			CAS 3114-55-4	Price
methanol purge & trap	100 µg/ml	1 ml	ref. P811060	BK
CHLOROBENZILATE			CAS 510-15-6	Price
n-hexane	10 µg/ml	10 ml	ref. P811090	BG
	100 µg/ml	1 ml	ref. P811100	BG
methanol	10 µg/ml	10 ml	ref. P811110	BG
	100 µg/ml	1 ml	ref. P811120	BG
2-CHLOROBENZOIC ACID			CAS 118-91-2	Price
water	1 104 µg/ml	1 ml	ref. F128991	BH
1-CHLORO-2,4-DINITROBENZENE			CAS 97-00-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P811210	BG
	100 µg/ml	1 ml	ref. P811220	BG
methanol	10 µg/ml	10 ml	ref. P811230	BG
	100 µg/ml	1 ml	ref. P811240	BG
CHLOROETHANE			CAS 75-00-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P811270	BG
	100 µg/ml	1 ml	ref. P811280	BG
2-CHLOROETHANOL			CAS 107-07-3	Price
methanol purge & trap	1000 µg/ml	1 ml	ref. P811310	BK
	100 µg/ml	1 ml	ref. P811320	BG
2-CHLOROETHYL VINYLETHER			CAS 110-75-8	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P811350	BG
	100 µg/ml	1 ml	ref. F116211	BG
	5 000 µg/ml	1 ml	ref. F116221	CE
1-CHLORO-2-FLUOROBENZENE			CAS 348-51-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P811430	BG
	100 µg/ml	1 ml	ref. P811440	BG
	2 000 µg/ml	1 ml	ref. F112951	DA
1-CHLORO-4-FLUOROBENZENE			CAS 352-33-0	Price
methanol purge & trap	2 500 µg/ml	1 ml	ref. P879380	CB
CHLOROFORM			CAS 67-66-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P811470	BG
	100 µg/ml	1 ml	ref. P811480	BG
CHLOROMETHANE			CAS 74-87-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P811510	BG
	100 µg/ml	1 ml	ref. P811520	BG

3-CHLORO-2-METHYLANILINE			CAS 87-60-5	Price
acetonitrile	100 µg/ml	1 ml	ref. P873370	BG
	10 µg/ml	10 ml	ref. P873380	BG
3-CHLORO-4-METHYLANILINE			CAS 95-74-9	Price
methanol	10 µg/ml	10 ml	ref. P873390	BG
4-CHLORO-2-METHYLANILINE			CAS 95-69-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P873520	BG
methanol	100 µg/ml	1 ml	ref. P873510	BG
5-CHLORO-2-METHYLANILINE			CAS 95-79-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P873680	BG
methanol	100 µg/ml	5 ml	ref. P873670	BG
4-CHLORO-3-METHYLPHENOL			CAS 59-50-7	Price
methanol	10 µg/ml	10 ml	ref. P811550	BG
	100 µg/ml	1 ml	ref. P811560	BG
1-CHLORONAPHTHALENE			CAS 90-13-1	Price
methanol	100 µg/ml	5 ml	ref. P872430	BG
2-CHLORONAPHTHALENE			CAS 91-58-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P811590	CE
	100 µg/ml	1 ml	ref. P811600	CE
iso-octane	10 µg/ml	10 ml	ref. P811610	CE
	100 µg/ml	1 ml	ref. P811620	CE
methanol	10 µg/ml	10 ml	ref. P811630	CE
	100 µg/ml	1 ml	ref. P811640	CE
CHLORONEB			CAS 2675-77-6	Price
iso-octane	10 µg/ml	10 ml	ref. P874770	BG
methanol	100 µg/ml	1 ml	ref. P874760	BG
4-CHLORO-2-NITROANILINE			CAS 89-63-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P811670	BG
	100 µg/ml	1 ml	ref. P811680	BG
1-CHLORO-2-NITROBENZENE			CAS 88-73-3	Price
methanol	10 µg/ml	10 ml	ref. P811710	BG
	100 µg/ml	1 ml	ref. P811720	BG
1-CHLORO-3-NITROBENZENE			CAS 121-73-3	Price
n-hexane	10 µg/ml	10 ml	ref. P811750	BG
	100 µg/ml	1 ml	ref. P811760	BG
iso-octane	1 000 µg/ml	1 ml	ref. F117501	BK
methanol	10 µg/ml	10 ml	ref. P811770	BG
	100 µg/ml	1 ml	ref. P811780	BG
1-CHLORO-4-NITROBENZENE			CAS 100-00-5	Price
methanol	10 µg/ml	10 ml	ref. P811810	BG
	100 µg/ml	1 ml	ref. P811820	BG
4-CHLORO-3-NITROBENZOTRIFLUORIDE			CAS 121-17-5	Price
acetone	1 000 µg/ml	1 ml	ref. F119101	BH
iso-octane	1 000 µg/ml	5 ml	ref. P879410	BH
2-CHLORO-3-NITROTOLUENE			CAS 3970-40-9	Price
methanol	10 µg/ml	10 ml	ref. P873020	BG
2-CHLORO-4-NITROTOLUENE			CAS 121-86-8	Price
methanol	100 µg/ml	1 ml	ref. P873030	BG
2-CHLORO-6-NITROTOLUENE			CAS 83-42-1	Price
methanol	100 µg/ml	1 ml	ref. P873040	BG
4-CHLORO-2-NITROTOLUENE			CAS 89-59-8	Price
methanol	10 µg/ml	10 ml	ref. P873530	BG
1-CHLOROOCATADECANE			CAS 3386-33-2	Price
dichloromethane	10 µg/ml	10 ml	ref. P811850	BH
	100 µg/ml	1 ml	ref. P811860	BH
methanol	10 µg/ml	10 ml	ref. P811870	BH
	100 µg/ml	1 ml	ref. P811880	BH
1-CHLOROOCANE			CAS 111-85-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P811910	BH
	100 µg/ml	1 ml	ref. P811920	BH
CHLOROPHACINONE			CAS 3691-35-8	Price
cyclohexane	100 µg/ml	1 ml	ref. P874780	BG
methanol	100 µg/ml	1 ml	ref. P874790	BG
2-CHLOROPHENOL			CAS 95-57-8	Price
methanol	10 µg/ml	10 ml	ref. P811950	BG
	100 µg/ml	1 ml	ref. P811960	BG

3-CHLOROPHENOL			CAS 108-43-0	Price
methanol	10 µg/ml	10 ml	ref. P812020	BG
	100 µg/ml	1 ml	ref. P812030	BG
4-CHLOROPHENOL			CAS 106-48-9	Price
methanol	10 µg/ml	10 ml	ref. P812060	BG
	100 µg/ml	1 ml	ref. P812070	BG
water	725 µg/ml	1 ml	ref. F128981	CF
2-CHLOROPHENOL ACETATE			CAS 4525-75-1	Price
iso-octane	10 µg/ml	10 ml	ref. P873050	BJ
3-CHLOROPHENOL ACETATE			CAS 13031-39-5	Price
iso-octane	10 µg/ml	10 ml	ref. P873400	BJ
4-CHLOROPHENOL ACETATE			CAS 876-27-7	Price
iso-octane	10 µg/ml	10 ml	ref. P873540	BJ
4-CHLOROPHENYL PHENYL ETHER			CAS 7005-72-3	Price
methanol	10 µg/ml	10 ml	ref. P873650	BG
	100 µg/ml	5 ml	ref. P873550	BG
CHLOROPICRIN			CAS 76-06-2	Price
methanol	1 000 µg/ml	1 ml	ref. P874800	BH
CHLOROPRENE			CAS 126-99-8	Price
methanol	100 µg/ml	1 ml	ref. F116231	BG
3-CHLOROPROPENE			CAS 107-05-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P812160	BG
	100 µg/ml	1 ml	ref. P812170	BG
methanol purge & trap	10 µg/ml	10 ml	ref. P812180	BG
	100 µg/ml	1 ml	ref. P812190	BG
CHLOROPROPYLATE			CAS 5836-10-2	Price
iso-octane	10 µg/ml	10 ml	ref. P879570	BG
methanol	100 µg/ml	1 ml	ref. P879560	BG
CHLOROTHALONIL			CAS 1897-45-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P812220	BG
	100 µg/ml	1 ml	ref. P812230	BG
cyclohexane	10 µg/ml	10 ml	ref. P812240	BG
	100 µg/ml	1 ml	ref. P812250	BG
α-CHLOROTOLUENE			CAS 100-44-7	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P812320	BG
	100 µg/ml	1 ml	ref. P812330	BG
2-CHLOROTOLUENE			CAS 95-49-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P812360	BG
	100 µg/ml	1 ml	ref. P812370	BG
methanol	10 µg/ml	10 ml	ref. P812380	BG
	100 µg/ml	1 ml	ref. P812390	BG
3-CHLOROTOLUENE			CAS 108-41-8	Price
iso-octane	10 µg/ml	10 ml	ref. P812420	BG
	100 µg/ml	1 ml	ref. P812430	BG
methanol	10 µg/ml	10 ml	ref. P812440	BG
	100 µg/ml	1 ml	ref. P812450	BG
4-CHLOROTOLUENE			CAS 106-43-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P812480	BG
	100 µg/ml	1 ml	ref. P812490	BG
methanol purge & trap	10 µg/ml	10 ml	ref. P812500	BG
	100 µg/ml	1 ml	ref. P812510	BG
CHLOROXURON			CAS 1982-47-4	Price
acetone	10 µg/ml	10 ml	ref. P812540	BG
	100 µg/ml	1 ml	ref. P812550	BG
acetonitrile	10 µg/ml	10 ml	ref. P812560	BG
	100 µg/ml	1 ml	ref. P812570	BG
methanol	10 µg/ml	10 ml	ref. P812580	BG
	100 µg/ml	1 ml	ref. P812590	BG
CHLORPROPHARM			CAS 101-21-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P874830	BG
	100 µg/ml	1 ml	ref. P874820	BG
iso-octane	10 µg/ml	10 ml	ref. P874840	BG
CHLORPYRIPHOS			CAS 2921-88-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P812620	BG
	100 µg/ml	1 ml	ref. P812630	BG
cyclohexane	10 µg/ml	10 ml	ref. P812640	BG
	100 µg/ml	1 ml	ref. P812650	BG
methanol	10 µg/ml	10 ml	ref. P812660	BG
	100 µg/ml	1 ml	ref. P812670	BG

CHLORPYRIPHOS-D10 [Diethyl-D10]			CAS 285138-81-0	Price
acetone	100 µg/ml	1 ml	ref. P812700	GB
	100 µg/ml	1 ml	ref. P812710	GB
CHLORPYRIFOS-METHYL			CAS 5598-13-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P812740	BG
	100 µg/ml	1 ml	ref. P812750	BG
cyclohexane	10 µg/ml	10 ml	ref. P812760	BG
	100 µg/ml	1 ml	ref. P812770	BG
CHLORSULFURON			CAS 64902-72-3	Price
acetonitrile	100 µg/ml	1 ml	ref. P874850	BH
ethyl acetate	10 µg/ml	10 ml	ref. P874860	BH
CHLORTAL DIACID			CAS 2136-79-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P815350	CE
	100 µg/ml	1 ml	ref. P815360	CE
methanol	10 µg/ml	10 ml	ref. P815370	CE
	100 µg/ml	1 ml	ref. P815380	CE
CHLORTHAL-DIMETHYL			CAS 1861-32-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P815290	BG
	100 µg/ml	1 ml	ref. P815300	BG
methanol	10 µg/ml	10 ml	ref. P815310	BG
	100 µg/ml	1 ml	ref. P815320	BG
CHLORTHIAMID			CAS 1918-13-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P874880	BG
methanol	100 µg/ml	1 ml	ref. P874870	BG
CHLORTHION			CAS 500-28-7	Price
cyclohexane	100 µg/ml	1 ml	ref. P879590	BG
iso-octane	10 µg/ml	10 ml	ref. P879600	BG
methanol	100 µg/ml	1 ml	ref. P879580	BG
CHLORTOLURON			CAS 15545-48-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P812920	BG
	100 µg/ml	1 ml	ref. P812930	BG
methanol	10 µg/ml	10 ml	ref. P812940	BG
	100 µg/ml	1 ml	ref. P812950	BG
CHLOZOLINATE			CAS 84332-86-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P874900	DB
	100 µg/ml	1 ml	ref. P874890	DB
CHRYSENE			CAS 218-01-9	Price
acetone	10 µg/ml	10 ml	ref. P813010	BH
	100 µg/ml	1 ml	ref. P813020	BH
acetonitrile	10 µg/ml	10 ml	ref. P813030	BH
	100 µg/ml	1 ml	ref. P813040	BH
cyclohexane	10 µg/ml	10 ml	ref. P813050	BH
	100 µg/ml	1 ml	ref. P813060	BH
dichloromethane	10 µg/ml	10 ml	ref. P813070	BH
	100 µg/ml	1 ml	ref. P813080	BH
CHRYSENE-D12			CAS 1719-03-5	Price
acetonitrile	100 µg/ml	1 ml	ref. P813110	CE
cyclohexane	100 µg/ml	1 ml	ref. P813120	CE
CLODINAFOP-PROPARGYL ESTER			CAS 105512-06-9	Price
acetone	10 µg/ml	10 ml	ref. P813190	BG
	100 µg/ml	1 ml	ref. P813200	BG
acetonitrile	10 µg/ml	10 ml	ref. P813210	BG
	100 µg/ml	1 ml	ref. P813220	BG
CLOFENTEZINE			CAS 74115-24-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P874920	BG
methanol	100 µg/ml	1 ml	ref. P874910	BG
CLOMAZONE			CAS 81777-89-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P813330	BG
	100 µg/ml	1 ml	ref. P813340	BG
cyclohexane	10 µg/ml	10 ml	ref. P813350	BG
	100 µg/ml	1 ml	ref. P813360	BG
CLOPYRALID			CAS 1702-17-6	Price
acetone	10 µg/ml	10 ml	ref. P813390	BG
	100 µg/ml	1 ml	ref. P813400	BG
acetonitrile	10 µg/ml	10 ml	ref. P813410	BG
	100 µg/ml	1 ml	ref. P813420	BG
CLOQUINTOCET-1-METHYLHEXYL ESTER			CAS 99607-70-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P874930	BG
cyclohexane	10 µg/ml	10 ml	ref. P874940	BG

CORONENE			CAS 191-07-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P874950	CA
cyclohexane	10 µg/ml	10 ml	ref. P874960	CA
COUMAPHOS			CAS 56-72-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P813610	BG
	100 µg/ml	1 ml	ref. P813620	BG
iso-octane	10 µg/ml	10 ml	ref. P813630	BG
	100 µg/ml	1 ml	ref. P813640	BG
CRIMIDINE			CAS 535-89-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P874980	BG
methanol	100 µg/ml	1 ml	ref. P874970	BG
CROTONALDEHYDE			CAS 123-73-9	Price
methanol	1 000 µg/ml	1 ml	ref. P874990	BH
CROTONALDEHYDE-DNPH			CAS 1527-96-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P813750	BG
	100 µg/ml	1 ml	ref. P813760	BG
CYANAZINE			CAS 21725-46-2	Price
acetone	10 µg/ml	10 ml	ref. P813910	BG
	100 µg/ml	1 ml	ref. P813920	BG
acetonitrile	10 µg/ml	10 ml	ref. P813930	BG
	100 µg/ml	1 ml	ref. P813940	BG
cyclohexane	10 µg/ml	10 ml	ref. P813950	BG
	100 µg/ml	1 ml	ref. P813960	BG
CYANOFENPHOS			CAS 13067-93-1	Price
iso-octane	10 µg/ml	10 ml	ref. P879620	BG
methanol	100 µg/ml	1 ml	ref. P879610	BG
CYANOPHOS			CAS 2636-26-2	Price
iso-octane	10 µg/ml	10 ml	ref. P875010	BJ
methanol	100 µg/ml	1 ml	ref. P875000	BJ
CYCLOATE			CAS 1134-23-2	Price
methanol	100 µg/ml	1 ml	ref. P875020	BG
CYCLOHEXANE			CAS 110-82-7	Price
methanol	1 000 µg/ml	5 ml	ref. P875030	BH
CYCLOHEXANONE			CAS 108-94-1	Price
methanol	1 000 µg/ml	1 ml	ref. P875040	BG
CYCLOHEXANONE-DNPH			CAS 1589-62-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P814070	BG
	100 µg/ml	1 ml	ref. P814080	BG
CYCLOHEXIMIDE			CAS 66-81-9	Price
methanol	100 µg/ml	1 ml	ref. P879630	BG
2-CYCLOHEXYL-4,6-DINITROPHENOL			CAS 131-89-5	Price
methanol	100 µg/ml	1 ml	ref. P879640	BJ
CYCLOPENTA(c,d)PYRENE			id : 8110	Price
cyclohexane	10 µg/ml	10 ml	ref. P875070	DG
CYFLUTHRIN			CAS 68359-37-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P814110	BG
	100 µg/ml	1 ml	ref. P814120	BG
cyclohexane	10 µg/ml	10 ml	ref. P814130	BG
	100 µg/ml	1 ml	ref. P814140	BG
λ-CYHALOTHRIN			CAS 91465-08-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P814210	BH
	100 µg/ml	1 ml	ref. P814220	BH
cyclohexane	10 µg/ml	10 ml	ref. P814230	BH
	100 µg/ml	1 ml	ref. P814240	BH
CYMOXANIL			CAS 57966-95-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P814310	BG
	100 µg/ml	1 ml	ref. P814320	BG
CYPERMETHRIN			CAS 52315-07-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P875090	BG
cyclohexane	100 µg/ml	1 ml	ref. P875080	BG
iso-octane	10 µg/ml	10 ml	ref. P875100	BG
α-CYPERMETHRIN			CAS 67375-30-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P814390	BG
	100 µg/ml	1 ml	ref. P814400	BG
cyclohexane	10 µg/ml	10 ml	ref. P874010	BG
	100 µg/ml	1 ml	ref. P874000	BG
CYPERMETHRIN-D5 [Phenoxy-Ring-D5]			id : 21270	Price
acetone	100 µg/ml	1 ml	ref. P814350	BCG
nonane	100 µg/ml	1 ml	ref. P814360	BCG

CYPHENOTRIN			CAS 39515-40-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P875120	BH
methanol	100 µg/ml	1 ml	ref. P875110	BH
CYPROCONAZOLE			CAS 113096-99-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P814480	BH
	100 µg/ml	1 ml	ref. P814490	BH
cyclohexane	10 µg/ml	10 ml	ref. P814500	BH
	100 µg/ml	1 ml	ref. P814510	BH
methanol	10 µg/ml	10 ml	ref. P814520	BH
	100 µg/ml	1 ml	ref. P814530	BH
CYPRODINIL			CAS 121552-61-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P814560	BG
	100 µg/ml	1 ml	ref. P814570	BG
iso-octane	10 µg/ml	10 ml	ref. P814580	BG
methanol	100 µg/ml	1 ml	ref. P814590	BG
	10 µg/ml	10 ml	ref. P814660	BG
	100 µg/ml	1 ml	ref. P814670	BG
CYROMAZINE			CAS 66215-27-8	Price
ethyl acetate	10 µg/ml	10 ml	ref. P875140	BG
methanol	100 µg/ml	1 ml	ref. P875130	BG
2,4-D			CAS 94-75-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P872760	BG
	100 µg/ml	1 ml	ref. P872750	BG
2,4-D D3			CAS 202480-67-9	Price
acetone	100 µg/ml	1 ml	ref. P872770	HA
2,4-D-1-BUTYL ESTER			CAS 94-80-4	Price
cyclohexane	100 µg/ml	1 ml	ref. P879650	BG
2,4-D METHYL ESTER			CAS 1928-38-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P816030	BG
	100 µg/ml	1 ml	ref. P816040	BG
methanol	10 µg/ml	10 ml	ref. P816050	BG
	100 µg/ml	1 ml	ref. P816060	BG
DALAPON			CAS 75-99-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P814850	BG
	100 µg/ml	1 ml	ref. P814860	BG
methanol	10 µg/ml	10 ml	ref. P814870	BG
	100 µg/ml	1 ml	ref. P814880	BG
methyl-tert-butylether	10 µg/ml	10 ml	ref. P814890	BG
	100 µg/ml	1 ml	ref. P814900	BG
DALAPON METHYL ESTER			CAS 17640-02-7	Price
n-hexane	10 µg/ml	10 ml	ref. P814920	BH
	100 µg/ml	1 ml	ref. P814930	BH
methanol	10 µg/ml	10 ml	ref. P814940	BH
	100 µg/ml	1 ml	ref. P814950	BH
DAZOMET			CAS 533-74-4	Price
methanol	100 µg/ml	1 ml	ref. P879670	BG
2,4-DB			CAS 94-82-6	Price
acetone	10 µg/ml	10 ml	ref. P815050	BG
	100 µg/ml	1 ml	ref. P815060	BG
methanol	10 µg/ml	10 ml	ref. P815070	BG
	100 µg/ml	1 ml	ref. P815080	BG
2,4-DB METHYL ESTER			CAS 18625-12-2	Price
methanol	100 µg/ml	1 ml	ref. P879660	BG
DECABROMODIPHENYL ETHER			CAS 1163-19-5	Price
iso-octane/toluene (9/1)	10 µg/ml	10 ml	ref. P837320	BG
	100 µg/ml	1 ml	ref. P837330	BG
n-nonane	10 µg/ml	10 ml	ref. P837340	BG
	100 µg/ml	1 ml	ref. P837350	BG
DECAFLUOROBIPHENYL			CAS 434-90-2	Price
acetone	10 µg/ml	10 ml	ref. P816090	BG
	100 µg/ml	1 ml	ref. P816100	BG
acetonitrile	100 µg/ml	1 ml	ref. F119901	BG
dichloromethane	10 µg/ml	10 ml	ref. P816110	BG
	100 µg/ml	1 ml	ref. P816120	BG
DECAFLUOROTRIPHENYLPHOSPHINE			CAS 5074-71-5	Price
acetone	10 µg/ml	10 ml	ref. P816150	BG
	100 µg/ml	1 ml	ref. P816160	BG
dichloromethane	10 µg/ml	10 ml	ref. P816170	BG
	100 µg/ml	1 ml	ref. P816180	BG
	500 µg/ml	1 ml	ref. F114071	BG

DECAFLUOROTRIPHENYLPHOSPHINE OXYDE			id : 17201	Price
acetonitrile	100 µg/ml	1 ml	ref. F119971	CA
DECANAL			CAS 112-31-2	Price
methanol	1 000 µg/ml	1 ml	ref. P875150	BH
n-DECANE			CAS 124-18-5	Price
dichloromethane	100 µg/ml	5 ml	ref. P877110	BG
DECANOIC ACID METHYL ESTER			CAS 110-42-9	Price
n-hexane	10 000 µg/ml	1 ml	ref. P875160	BJ
n-DECYLBENZENE			CAS 104-72-3	Price
methanol	1 000 µg/ml	5 ml	ref. P877120	CG
2,4'-DDD			CAS 53-19-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P815530	BG
	100 µg/ml	1 ml	ref. P815540	BG
methanol	10 µg/ml	10 ml	ref. P815550	BG
	100 µg/ml	1 ml	ref. P815560	BG
4,4'-DDD (TDE)			CAS 72-54-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P815590	BG
	100 µg/ml	1 ml	ref. P815600	BG
methanol	10 µg/ml	10 ml	ref. P815610	BG
	100 µg/ml	1 ml	ref. P815620	BG
2,4'-DDE			CAS 3424-82-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P815690	BJ
	100 µg/ml	1 ml	ref. P815700	BJ
methanol	10 µg/ml	10 ml	ref. P815710	BJ
	100 µg/ml	1 ml	ref. P815720	BJ
4,4'-DDE			CAS 72-55-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P815750	BG
	100 µg/ml	1 ml	ref. P815760	BG
methanol	10 µg/ml	10 ml	ref. P815770	BG
	100 µg/ml	1 ml	ref. P815780	BG
4,4'-DDMU			CAS 1022-22-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P879690	BG
methanol	100 µg/ml	1 ml	ref. P879680	BG
2,4'-DDT			CAS 789-02-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P815850	BJ
	100 µg/ml	1 ml	ref. P815860	BJ
methanol	10 µg/ml	10 ml	ref. P815870	BJ
	100 µg/ml	1 ml	ref. P815880	BJ
4,4'-DDT			CAS 50-29-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P815940	BG
	100 µg/ml	1 ml	ref. P815950	BG
methanol	10 µg/ml	10 ml	ref. P815960	BG
	100 µg/ml	1 ml	ref. P815970	BG
DEET			CAS 134-62-3	Price
methanol	10 µg/ml	10 ml	ref. P879710	BG
	100 µg/ml	1 ml	ref. P879700	BG
DELTAMETHRIN			CAS 52918-63-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P816310	BG
	100 µg/ml	1 ml	ref. P816320	BG
cyclohexane	10 µg/ml	10 ml	ref. P816330	BG
	100 µg/ml	1 ml	ref. P816340	BG
DEMETON (O+S)			CAS 8065-48-3	Price
acetonitrile	100 µg/ml	1 ml	ref. P875190	BH
cyclohexane	10 µg/ml	10 ml	ref. P875200	BH
DEMETON (S)			CAS 126-75-0	Price
iso-octane	10 µg/ml	10 ml	ref. P875220	BH
methanol	100 µg/ml	1 ml	ref. P875210	BH
DEMETON-S-METHYL			CAS 919-86-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P816410	BG
	100 µg/ml	1 ml	ref. P816420	BG
iso-octane	10 µg/ml	10 ml	ref. P816430	BG
	100 µg/ml	1 ml	ref. P816440	BG
DEMETON-S-METHYL -SULFONE			CAS 17040-19-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P875240	BG
methanol	100 µg/ml	1 ml	ref. P875230	BG
DEMETON-S-METHYL-SULFOXIDE			CAS 301-12-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P836080	BG
	100 µg/ml	1 ml	ref. P836090	BG
cyclohexane/toluene (1/1)	10 µg/ml	10 ml	ref. P836100	BG
	100 µg/ml	1 ml	ref. P836110	BG

DESMETRYN			CAS 1014-69-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P816630	BG
	100 µg/ml	1 ml	ref. P816640	BG
iso-octane	10 µg/ml	10 ml	ref. P816650	BG
	100 µg/ml	1 ml	ref. P816660	BG
DIALIFOS			CAS 10311-84-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P879730	BG
	100 µg/ml	1 ml	ref. P879740	BG
DIALLATE			CAS 2303-16-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P875260	BH
methanol	100 µg/ml	1 ml	ref. P875250	BH
2,4-DIAMINOTOLUENE			CAS 95-80-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P872780	BG
toluene	100 µg/ml	5 ml	ref. P872790	BG
DIAZINON			CAS 333-41-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P816810	BG
	100 µg/ml	1 ml	ref. P816820	BG
cyclohexane	10 µg/ml	10 ml	ref. P816830	BG
	100 µg/ml	1 ml	ref. P816840	BG
DIAZINON-D10			CAS 100155-47-3	Price
acetone	100 µg/ml	1 ml	ref. P875270	GC
DIBENZ[a,h]ACRIDINE			CAS 226-36-8	Price
methanol	100 µg/ml	1 ml	ref. P875280	CA
DIBENZ[a,j]ACRIDINE			CAS 224-42-0	Price
toluene	100 µg/ml	1 ml	ref. P875360	CB
DIBENZO(a,h)ANTHRACENE			CAS 53-70-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P816910	BJ
	100 µg/ml	1 ml	ref. P816920	BJ
cyclohexane	10 µg/ml	10 ml	ref. P816930	BJ
	100 µg/ml	1 ml	ref. P816940	BJ
dichloromethane	10 µg/ml	10 ml	ref. P816950	BJ
	100 µg/ml	1 ml	ref. P816960	BJ
DIBENZO(a,h)ANTHRACENE-D14			CAS 13250-98-1	Price
cyclohexane	100 µg/ml	1 ml	ref. P816990	DF
DIBENZOFURAN			CAS 132-64-9	Price
iso-octane	10 µg/ml	10 ml	ref. P817020	BJ
	100 µg/ml	1 ml	ref. P817030	BJ
DIBENZO(a,e)PYRENE			CAS 192-65-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P875310	CC
toluene	100 µg/ml	1 ml	ref. P875290	CC
DIBENZO(a,h)PYRENE			CAS 189-64-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P875320	CJ
cyclohexane	10 µg/ml	10 ml	ref. P875330	CJ
DIBENZO(a,i)PYRENE			CAS 189-55-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P875340	EF
cyclohexane	10 µg/ml	10 ml	ref. P875350	EF
DIBENZOTHIOPHENE			CAS 132-65-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P817060	BG
	100 µg/ml	1 ml	ref. P817070	BG
iso-octane	10 µg/ml	10 ml	ref. P817080	BG
	100 µg/ml	1 ml	ref. P817090	BG
DIBENZYLPHTHALATE-D4 [Ring-D4]			id : 21310	Price
acetone	10 µg/ml	1 ml	ref. P817120	BBG
nonane	100 µg/ml	1 ml	ref. P817130	BBG
4,4'-DIBROMOBIPHENYL			CAS 92-86-4	Price
dichloromethane	10 µg/ml	10 ml	ref. P817180	BJ
	100 µg/ml	1 ml	ref. P817190	BJ
ethyl acetate	10 µg/ml	10 ml	ref. P817160	BJ
	100 µg/ml	1 ml	ref. P817170	BJ
	1 000 µg/ml	1 ml	ref. F113651	CC
n-hexane	100 µg/ml	1 ml	ref. F113691	BJ
	1 000 µg/ml	1 ml	ref. F147420	CC
iso-octane	1 000 µg/ml	1 ml	ref. F117561	CC
DIBROMOCHLOROMETHANE			CAS 124-48-1	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P817220	BG
	100 µg/ml	1 ml	ref. P817230	BG

1,2-DIBROMO-3-CHLOROPROPANE			CAS 96-12-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P817260	BG
	100 µg/ml	1 ml	ref. P817270	BG
methanol purge & trap	10 µg/ml	10 ml	ref. P817280	BG
	100 µg/ml	1 ml	ref. P817290	BG
1,2-DIBROMOETHANE			CAS 106-93-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P817320	BG
	100 µg/ml	1 ml	ref. P817330	BG
methanol purge & trap	10 µg/ml	10 ml	ref. P817340	BG
	100 µg/ml	1 ml	ref. P817350	BG
DIBROMOMETHANE			CAS 74-95-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P817380	BG
	100 µg/ml	1 ml	ref. P817390	BG
4,4'-DIBROMOOCTAFLUOROBIPHENYL			CAS 10386-84-2	Price
acetone	250 µg/ml	1 ml	ref. F119150	BG
n-hexane	10 µg/ml	10 ml	ref. P817420	BG
	100 µg/ml	1 ml	ref. P817430	BG
dichloromethane	10 µg/ml	10 ml	ref. P817440	BG
	100 µg/ml	1 ml	ref. P817450	BG
methanol	5 000 µg/ml	1 ml	ref. F113911	CE
methyl- <i>tert</i> -butylether	10 µg/ml	10 ml	ref. P817460	BG
	100 µg/ml	1 ml	ref. F113781	BG
2,4-DIBROMOPHENOL			CAS 615-58-7	Price
acetone	1 000 µg/ml	5 ml	ref. F109280	CF
ethanol	1 000 µg/ml	1 ml	ref. F131331	BJ
methanol	10 µg/ml	10 ml	ref. P872810	BJ
	10 µg/ml	5 ml	ref. F109310	BJ
	100 µg/ml	1 ml	ref. P872800	BJ
2-propanol	1 000 µg/ml	1 ml	ref. F116531	BJ
2,6-DIBROMOPHENOL			CAS 608-33-3	Price
ethanol	1 000 µg/ml	1 ml	ref. F131341	BH
1,2-DIBROMOPROPANE			CAS 78-75-1	Price
n-hexane	10 µg/ml	10 ml	ref. F114791	BG
	400 µg/ml	10 ml	ref. F114831	BG
	1 000 µg/ml	1 ml	ref. F114781	BH
2,3-DIBROMOPROPIONIC ACID			CAS 600-05-5	Price
methyl- <i>tert</i> -butylether	1 000 µg/ml	1 ml	ref. F114611	BH
	2 000 µg/ml	5 ml	ref. P884180	BJ
2,3-DIBROMOPROPIONIC ACID METHYL ESTER			CAS 1729-67-5	Price
methyl- <i>tert</i> -butylether	1 000 µg/ml	1 ml	ref. F114601	BH
2,5-DIBROMOTOLUENE			CAS 615-59-8	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P817500	BG
	100 µg/ml	1 ml	ref. P817510	BG
DI-<i>n</i>-BUTYL PHTHALATE-D4			id : 19521	Price
acetone	100 µg/ml	1 ml	ref. P817620	BCG
nonane	100 µg/ml	1 ml	ref. P817630	BCG
DIBUTYL TIN-DICHLORIDE			CAS 683-18-1	Price
methanol	10 µg/ml	10 ml	ref. P875390	BG
DICAMBA			CAS 1918-00-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P817660	BG
	100 µg/ml	1 ml	ref. P817670	BG
methanol	10 µg/ml	10 ml	ref. P817680	BG
	100 µg/ml	1 ml	ref. P817690	BG
DICAMBA METHYL ESTER			CAS 6597-78-0	Price
methanol	10 µg/ml	10 ml	ref. P817720	BH
	100 µg/ml	1 ml	ref. P817730	BH
DICAPTHON			CAS 2463-84-5	Price
methanol	100 µg/ml	1 ml	ref. P879760	BJ
	1 000 µg/ml	1 ml	ref. P879750	DA
toluene	10 µg/ml	10 ml	ref. P879770	BJ
DICHOLOBENIL			CAS 1194-65-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P817800	BG
	100 µg/ml	1 ml	ref. P817810	BG
DICHOLOFENTHION			CAS 97-17-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P817840	BG
	100 µg/ml	1 ml	ref. P817850	BG
DICHOFLUANID			CAS 1085-98-9	Price
iso-octane	10 µg/ml	10 ml	ref. P817880	BG
	100 µg/ml	1 ml	ref. P817890	BG

DICHLONE			CAS 117-80-6	Price
methanol	100 µg/ml	1 ml	ref. P879780	BG
DICHLORMID			CAS 37764-25-3	Price
acetonitrile	100 µg/ml	1 ml	ref. P875400	BG
cyclohexane	10 µg/ml	10 ml	ref. P875410	BG
2,3-DICHLOROANILINE			CAS 608-27-5	Price
methanol	100 µg/ml	10 ml	ref. P872720	BG
2,4-DICHLOROANILINE			CAS 554-00-7	Price
methanol	10 µg/ml	10 ml	ref. P817960	BG
	100 µg/ml	1 ml	ref. P817970	BG
toluene	10 µg/ml	10 ml	ref. P817980	BG
	100 µg/ml	1 ml	ref. P817990	BG
2,5-DICHLOROANILINE			CAS 95-82-9	Price
methanol	100 µg/ml	1 ml	ref. P872860	BG
2,6-DICHLOROANILINE			CAS 608-31-1	Price
methanol	100 µg/ml	1 ml	ref. P872890	BG
3,4-DICHLOROANILINE			CAS 95-76-1	Price
acetonitrile	100 µg/ml	1 ml	ref. P873290	BG
cyclohexane	10 µg/ml	10 ml	ref. P873300	BG
3,5-DICHLOROANILINE			CAS 626-43-7	Price
methanol	100 µg/ml	1 ml	ref. P873330	BG
2,6-DICHLOROBENZAMIDE			CAS 2008-58-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P872910	BG
	100 µg/ml	1 ml	ref. P872900	BG
1,2-DICHLOROBENZENE			CAS 95-50-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P818020	BG
	100 µg/ml	1 ml	ref. P818030	BG
cyclohexane	10 µg/ml	10 ml	ref. P818040	BG
	100 µg/ml	1 ml	ref. P818050	BG
methanol	10 µg/ml	10 ml	ref. P818060	BG
	100 µg/ml	1 ml	ref. P818070	BG
1,2-DICHLOROBENZENE-D4			CAS 2199-69-1	Price
methanol purge & trap	100 µg/ml	1 ml	ref. P818100	BK
	2 000 µg/ml	1 ml	ref. F113961	CJ
1,3-DICHLOROBENZENE			CAS 541-73-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P818130	BG
	100 µg/ml	1 ml	ref. P818140	BG
methanol	10 µg/ml	10 ml	ref. P818150	BG
	100 µg/ml	1 ml	ref. P818160	BG
1,4-DICHLOROBENZENE			CAS 106-46-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P818190	BG
	100 µg/ml	1 ml	ref. P818200	BG
iso-octane	10 µg/ml	10 ml	ref. P818210	BG
	100 µg/ml	1 ml	ref. P818220	BG
methanol	10 µg/ml	10 ml	ref. P818230	BG
	100 µg/ml	1 ml	ref. P818240	BG
1,4-DICHLOROBENZENE D4			CAS 3855-82-1	Price
dichloromethane	4 000 µg/ml	1 ml	ref. P880920	DD
methanol	2 000 µg/ml	1 ml	ref. P872380	CF
3,3'-DICHLOROBENZIDINE			CAS 91-94-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P818270	BG
	100 µg/ml	1 ml	ref. P818280	BG
dichloromethane	10 µg/ml	10 ml	ref. P818310	BG
	100 µg/ml	1 ml	ref. P818320	BG
methanol	10 µg/ml	10 ml	ref. P818290	BG
	100 µg/ml	1 ml	ref. P818300	BG
3,5-DICHLOROBENZOIC ACID			CAS 51-36-5	Price
methyl- <i>tert</i> -butylether	10 µg/ml	10 ml	ref. P818350	BG
	100 µg/ml	1 ml	ref. P818360	BG
	1 000 µg/ml	1 ml	ref. F114621	BK
3,5-DICHLOROBENZOIC ACID METHYL ESTER			CAS 2905-67-1	Price
methanol	1 000 µg/ml	1 ml	ref. P884190	EC
methyl- <i>tert</i> -butylether	1 000 µg/ml	1 ml	ref. F114631	EC
1,4-DICHLOROBUTANE			CAS 110-56-5	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P818440	BG
	100 µg/ml	1 ml	ref. P818450	BG
	2 000 µg/ml	1 ml	ref. F112941	DA
trans -1,4-DICHLORO-2-BUTENE			CAS 110-57-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P818480	BG
	100 µg/ml	1 ml	ref. P818490	BG

1,1-DICHLOROETHANE			CAS 75-34-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P818560	BG
	100 µg/ml	1 ml	ref. P818570	BG
1,2-DICHLOROETHANE			CAS 107-06-2	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P818600	BG
	100 µg/ml	1 ml	ref. P818610	BG
1,2-DICHLOROETHANE-D4			CAS 17060-07-0	Price
methanol purge & trap	100 µg/ml	1 ml	ref. P818640	CB
	2 000 µg/ml	1 ml	ref. F062435	CJ
1,1-DICHLOROETHENE			CAS 75-35-4	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P818670	BG
	100 µg/ml	1 ml	ref. P818680	BG
cis -1,2-DICHLOROETHENE			CAS 156-59-2	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P818710	BG
	100 µg/ml	1 ml	ref. P818720	BG
trans -1,2-DICHLOROETHENE			CAS 156-60-5	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P818750	BG
	100 µg/ml	1 ml	ref. P818760	BG
DICHLOROMETHANE			CAS 75-09-2	Price
methanol	10 µg/ml	10 ml	ref. P818790	BG
	100 µg/ml	1 ml	ref. P818800	BG
2,3-DICHLORONITROBENZENE			CAS 3209-22-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P818830	BG
	100 µg/ml	1 ml	ref. P818840	BG
methanol	10 µg/ml	10 ml	ref. P818850	BG
	100 µg/ml	1 ml	ref. P818860	BG
2,4-DICHLORONITROBENZENE			CAS 611-06-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P872820	BG
toluene	100 µg/ml	1 ml	ref. P872830	BG
2,5-DICHLORONITROBENZENE			CAS 89-61-2	Price
methanol	10 µg/ml	10 ml	ref. P818890	BG
	100 µg/ml	1 ml	ref. P818900	BG
3,4-DICHLORONITROBENZENE			CAS 99-54-7	Price
methanol	10 µg/ml	10 ml	ref. P818930	BG
	100 µg/ml	1 ml	ref. P818940	BG
3,5-DICHLORONITROBENZENE			CAS 618-62-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P873340	BG
DICHLOROPHEN			CAS 97-23-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P879800	BG
methanol	100 µg/ml	1 ml	ref. P879790	BG
2,3-DICHLOROPHENOL			CAS 576-24-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P818970	BG
	100 µg/ml	1 ml	ref. P818980	BG
methanol	10 µg/ml	10 ml	ref. P818990	BG
	100 µg/ml	1 ml	ref. P819000	BG
2,4-DICHLOROPHENOL			CAS 120-83-2	Price
methanol	10 µg/ml	10 ml	ref. P819030	BG
	100 µg/ml	1 ml	ref. P819040	BG
2,4-DICHLOROPHENOL-D3 [3,5,6-D3]			CAS 93951-74-7	Price
acetone	100 µg/ml	1 ml	ref. P819070	FJ
2,5-DICHLOROPHENOL			CAS 583-78-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P819100	BG
	100 µg/ml	1 ml	ref. P819110	BG
methanol	10 µg/ml	10 ml	ref. P819120	BG
	100 µg/ml	1 ml	ref. P819130	BG
2,6-DICHLOROPHENOL			CAS 87-65-0	Price
methanol	10 µg/ml	10 ml	ref. P872930	BG
	100 µg/ml	1 ml	ref. P872920	BG
3,4-DICHLOROPHENOL			CAS 95-77-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P819200	BG
	100 µg/ml	1 ml	ref. P819210	BG
methanol	10 µg/ml	10 ml	ref. P819220	BG
	100 µg/ml	1 ml	ref. P819230	BG
3,5-DICHLOROPHENOL			CAS 591-35-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P819260	BG
	100 µg/ml	1 ml	ref. P819270	BG
methanol	10 µg/ml	10 ml	ref. P819280	BG
	100 µg/ml	1 ml	ref. P819290	BG

2,4-DICHLOROPHENOL ACETATE			CAS 6341-97-5	Price
iso-octane	10 µg/ml	10 ml	ref. P872840	BJ
2,5-DICHLOROPHENOL ACETATE			CAS 30124-46-0	Price
iso-octane	10 µg/ml	10 ml	ref. P872870	BJ
2,6-DICHLOROPHENOL ACETATE			CAS 28165-71-1	Price
iso-octane	10 µg/ml	10 ml	ref. P872940	BJ
3,4-DICHLOROPHENOL ACETATE			CAS 17847-51-7	Price
iso-octane	10 µg/ml	10 ml	ref. P873310	BJ
3,5-DICHLOROPHENOL ACETATE			CAS 61925-86-8	Price
iso-octane	10 µg/ml	10 ml	ref. P873350	BJ
2,4-DICHLOROPHENYLACETIC ACID			CAS 19719-28-9	Price
acetone	100 µg/ml	1 ml	ref. F119130	BG
	2 000 µg/ml	1 ml	ref. F119140	CA
	2 000 µg/ml	5 ml	ref. P884220	DC
methanol	200 µg/ml	1 ml	ref. P884200	BG
	5 000 µg/ml	1 ml	ref. F113891	CE
methyl-tert-butylether	100 µg/ml	1 ml	ref. F113761	BG
2,4-DICHLOROPHENYLACETIC ACID METHYL ESTER			id : 15851	Price
methanol	5 000 µg/ml	1 ml	ref. F113901	CJ
methyl-tert-butylether	100 µg/ml	1 ml	ref. F113771	BG
1-(3,4-DICHLOROPHENYL)-3-METHYLUREA			CAS 3567-62-2	Price
acetone	10 µg/ml	10 ml	ref. P815410	BG
	100 µg/ml	1 ml	ref. P815420	BG
acetonitrile	10 µg/ml	10 ml	ref. P815430	BG
	100 µg/ml	1 ml	ref. P815440	BG
1,2-DICHLOROPROPANE			CAS 78-87-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P819320	BG
	100 µg/ml	1 ml	ref. P819330	BG
methanol purge & trap	10 µg/ml	10 ml	ref. P819340	BG
	100 µg/ml	1 ml	ref. P819350	BG
1,3-DICHLOROPROPANE			CAS 142-28-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P819380	BG
	100 µg/ml	1 ml	ref. P819390	BG
methanol	10 µg/ml	10 ml	ref. P819400	BG
	100 µg/ml	1 ml	ref. P819410	BG
2,2-DICHLOROPROPANE			CAS 594-20-7	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P819440	BG
	100 µg/ml	1 ml	ref. P819450	BG
1,3-DICHLOROPROPANE-2-OL			CAS 96-23-1	Price
methanol	5 000 µg/ml	1 ml	ref. P872270	CA
1,1-DICHLORO-2-PROPANONE			CAS 513-88-2	Price
acetone	100 µg/ml	5 ml	ref. P872040	BG
1,1-DICHLOROPROPENE			CAS 563-58-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P819480	BH
	100 µg/ml	1 ml	ref. P819490	BH
CIS -1,3-DICHLOROPROPENE			CAS 10061-01-5	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P819560	BJ
	100 µg/ml	1 ml	ref. P819570	BJ
TRANS -1,3-DICHLOROPROPENE			CAS 10061-02-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P819600	BJ
	100 µg/ml	1 ml	ref. P819610	BJ
1,3-DICHLOROPROPENE (cis + trans)			CAS 542-75-6	Price
methanol	200 µg/ml	1 ml	ref. P880830	BJ
	200 total µg/ml	1 ml	ref. P872300	BJ
	4 000 µg/ml	1 ml	ref. P872280	CG
2,3-DICHLOROPROPENE			CAS 78-88-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P819520	BG
	100 µg/ml	1 ml	ref. P819530	BG
2,3-DICHLOROPROPIONIC ACID METHYL ESTER			CAS 3674-09-7	Price
methyl-tert-butylether	1 000 µg/ml	1 ml	ref. P884230	BK
DICHLORPROP			CAS 120-36-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P819640	BG
	100 µg/ml	1 ml	ref. P819650	BG
methanol	10 µg/ml	10 ml	ref. P819660	BG
	100 µg/ml	1 ml	ref. P819670	BG

DICHLORPROP METHYL ESTER			CAS 57153-17-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P819740	BG
	100 µg/ml	1 ml	ref. P819750	BG
methanol	10 µg/ml	10 ml	ref. P819760	BG
	100 µg/ml	1 ml	ref. P819770	BG
DICHLORPROP-P			CAS 15165-67-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P819830	BG
	100 µg/ml	1 ml	ref. P819840	BG
DICHLORVOS			CAS 62-73-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P819870	BG
	100 µg/ml	1 ml	ref. P819880	BG
cyclohexane	10 µg/ml	10 ml	ref. P819890	BG
	100 µg/ml	1 ml	ref. P819900	BG
DICHLORVOS-D6 [Dimethyl-D6]			id : 8080	Price
cyclohexane	100 µg/ml	1 ml	ref. P819930	FK
n-nonane	100 µg/ml	1 ml	ref. P819940	FK
DICLOBUTRAZOL			CAS 75736-33-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P875440	BG
toluene	100 µg/ml	1 ml	ref. P875430	BG
DICLOFOP-METHYL			CAS 51338-27-3	Price
acetone	10 µg/ml	10 ml	ref. P817920	BG
	100 µg/ml	1 ml	ref. P817930	BG
DICLORAN			CAS 99-30-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P819970	BG
	100 µg/ml	1 ml	ref. P819980	BG
DICOFOL			CAS 115-32-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P820010	BG
	100 µg/ml	1 ml	ref. P820020	BG
cyclohexane	10 µg/ml	10 ml	ref. P820030	BG
	100 µg/ml	1 ml	ref. P820040	BG
iso-octane	10 µg/ml	10 ml	ref. P820050	BG
	100 µg/ml	1 ml	ref. P820060	BG
DICROTOPHOS			CAS 141-66-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P879820	BG
methanol	100 µg/ml	1 ml	ref. P879810	BG
DICYCLOHEXYLPHthalate-D4			id : 21250	Price
acetone	100 µg/ml	1 ml	ref. P820090	BKA
nonane	100 µg/ml	1 ml	ref. P820100	BKA
DIELDRIN			CAS 60-57-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P820130	BH
	100 µg/ml	1 ml	ref. P820140	BH
methanol	10 µg/ml	10 ml	ref. P820150	BH
	100 µg/ml	1 ml	ref. P820160	BH
DIENOCHLOR			CAS 2227-17-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P879830	BG
DIETHANOLAMINE			CAS 111-42-2	Price
methanol	100 µg/ml	1 ml	ref. P875450	BG
DIETHOFENCARB			CAS 87130-20-9	Price
acetone	10 µg/ml	10 ml	ref. P820390	BH
	100 µg/ml	1 ml	ref. P820400	BH
acetonitrile	10 µg/ml	10 ml	ref. P820410	BH
	100 µg/ml	1 ml	ref. P820420	BH
cyclohexane	10 µg/ml	10 ml	ref. P820430	BH
	100 µg/ml	1 ml	ref. P820440	BH
DIETHYLAMINE			CAS 109-89-7	Price
n-hexane	10 µg/ml	10 ml	ref. P820190	BG
	100 µg/ml	1 ml	ref. P820200	BG
1,2-DIETHYLBENZENE			CAS 135-01-3	Price
methanol	10 µg/ml	10 ml	ref. P820230	BG
	100 µg/ml	1 ml	ref. P820240	BG
1,3-DIETHYLBENZENE			CAS 141-93-5	Price
methanol	10 µg/ml	10 ml	ref. P820270	BG
	100 µg/ml	1 ml	ref. P820280	BG
1,4-DIETHYLBENZENE			CAS 105-05-5	Price
methanol	10 µg/ml	10 ml	ref. P820310	BG
	100 µg/ml	1 ml	ref. P820320	BG
DIETHYLETHER			CAS 60-29-7	Price
methanol	100 µg/ml	1 ml	ref. P875460	BG

DIFENOCONAZOLE			CAS 119446-68-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P820470	BG
	100 µg/ml	1 ml	ref. P820480	BG
cyclohexane	10 µg/ml	10 ml	ref. P820490	BG
	100 µg/ml	1 ml	ref. P820500	BG
DIFENOXURON			CAS 14214-32-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P879850	BG
methanol	100 µg/ml	1 ml	ref. P879840	BG
DIFLUBENZURON			CAS 35367-38-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P875480	BG
	100 µg/ml	1 ml	ref. P875470	BG
DIFLUFENICAN			CAS 83164-33-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P820530	BG
	100 µg/ml	1 ml	ref. P820540	BG
1,4-DIFLUOROBENZENE			CAS 540-36-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P820570	BG
	100 µg/ml	1 ml	ref. P820580	BG
DI-n-HEXYL PHTHALATE D4			id : 21240	Price
acetone	100 µg/ml	1 ml	ref. P820650	BCG
nonane	100 µg/ml	1 ml	ref. P820660	BCG
DIISOBUTYL PHTHALATE-D4			id : 21260	Price
acetone	100 µg/ml	1 ml	ref. P820730	BDF
nonane	100 µg/ml	1 ml	ref. P820740	BDF
DIMEFURON			CAS 34205-21-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P875500	BG
	100 µg/ml	1 ml	ref. P875490	BG
DIMETHACHLOR			CAS 50563-36-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P820800	BG
	100 µg/ml	1 ml	ref. P820810	BG
cyclohexane	10 µg/ml	10 ml	ref. P820820	BG
	100 µg/ml	1 ml	ref. P820830	BG
DIMETHENAMID			CAS 87674-68-8	Price
methanol	10 µg/ml	10 ml	ref. P820860	BH
	100 µg/ml	1 ml	ref. P820870	BH
DIMETHIPIN			CAS 55290-64-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P820900	BH
	100 µg/ml	1 ml	ref. P820910	BH
DIMETHIRIMOL			CAS 5221-53-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P820940	CJ
	100 µg/ml	1 ml	ref. P820950	CJ
DIMETHOATE			CAS 60-51-5	Price
acetone	10 µg/ml	10 ml	ref. P820980	BG
	100 µg/ml	1 ml	ref. P820990	BG
acetonitrile	10 µg/ml	10 ml	ref. P821000	BG
	100 µg/ml	1 ml	ref. P821010	BG
cyclohexane	10 µg/ml	10 ml	ref. P821020	BG
	100 µg/ml	1 ml	ref. P821030	BG
DIMETHOMORPH			CAS 110488-70-5	Price
acetone	10 µg/ml	10 ml	ref. P821090	BG
	100 µg/ml	1 ml	ref. P821100	BG
acetonitrile	10 µg/ml	10 ml	ref. P821110	BG
	100 µg/ml	1 ml	ref. P821120	BG
cyclohexane	10 µg/ml	10 ml	ref. P821130	BG
	100 µg/ml	1 ml	ref. P821140	BG
3,3'-DIMETHOXYBENZIDINE			CAS 119-90-4	Price
acetonitrile	100 µg/ml	1 ml	ref. P873230	BG
methanol	100 µg/ml	5 ml	ref. P873240	BG
DIMETHYLAMINE			CAS 124-40-3	Price
methanol	10 µg/ml	10 ml	ref. P821190	BG
	100 µg/ml	1 ml	ref. P821200	BG
water	10 µg/ml	10 ml	ref. P821170	BG
	100 µg/ml	1 ml	ref. P821180	BG
2,6-DIMETHYLANILINE			CAS 87-62-7	Price
iso-octane	100 µg/ml	1 ml	ref. P872950	BG
n,n-DIMETHYLANILINE			CAS 121-69-7	Price
dichloromethane	200 µg/ml	5 ml	ref. P877070	BH
2,5-DIMETHYLBENZALDEHYDE-DNPH			CAS 152477-96-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P821230	BG
	100 µg/ml	1 ml	ref. P821240	BG

3,3-DIMETHYLBENZIDINE			CAS 119-93-7	Price
methanol	100 µg/ml	5 ml	ref. P873250	FE
7,12-DIMETHYLBENZO(a)ANTHRACENE			CAS 57-97-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P873780	BK
toluene	100 µg/ml	5 ml	ref. P873770	DD
1,6-DIMETHYLNAPHTHALENE			CAS 575-43-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P821270	BJ
	100 µg/ml	1 ml	ref. P821280	BJ
cyclohexane	10 µg/ml	10 ml	ref. P821290	BJ
	100 µg/ml	1 ml	ref. P821300	BJ
1,3-DIMETHYL-2-NITROBENZENE			CAS 81-20-9	Price
methyl-tert-butylether	250 µg/ml	1 ml	ref. F113321	BH
	1 000 µg/ml	1 ml	ref. P872370	CA
3,6-DIMETHYLPHENANTRENE			CAS 1576-67-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P821330	CE
	100 µg/ml	1 ml	ref. P821340	CE
cyclohexane	10 µg/ml	10 ml	ref. P821350	CE
	100 µg/ml	1 ml	ref. P821360	CE
2,4-DIMETHYLPHENOL			CAS 105-67-9	Price
methanol	100 µg/ml	1 ml	ref. P872850	BG
2,4-DIMETHYLPHENOL-3,5,6-D3			CAS 93951-75-8	Price
n-nonane	100 µg/ml	1 ml	ref. P821390	EJ
2,5-DIMETHYLPHENOL			CAS 95-87-4	Price
methanol	100 µg/ml	1 ml	ref. P872880	BG
2,6-DIMETHYLPHENOL			CAS 576-26-1	Price
methanol	10 µg/ml	10 ml	ref. P821420	BG
	100 µg/ml	1 ml	ref. P821430	BG
3,4-DIMETHYLPHENOL			CAS 95-65-8	Price
methanol	100 µg/ml	1 ml	ref. P874410	BG
3,5-DIMETHYLPHENOL			CAS 108-68-9	Price
methanol	100 µg/ml	1 ml	ref. P873360	BG
DIMETHYL PHTHALATE-D4 [3,4,5,6-D4]			id : 21300	Price
acetone	100 µg/ml	1 ml	ref. P821460	BBG
nonane	100 µg/ml	1 ml	ref. P821470	BBG
DIMETHYL PHTHALATE-D6 [Methyl-D6]			id : 21320	Price
acetone	100 µg/ml	1 ml	ref. P821500	HC
nonane	100 µg/ml	1 ml	ref. P821510	HC
DIMETILAN			CAS 644-64-4	Price
cyclohexane	10 µg/ml	1 ml	ref. P875510	BG
DINICONAZOLE			CAS 83657-24-3	Price
acetone	10 µg/ml	10 ml	ref. P821540	BH
	100 µg/ml	1 ml	ref. P821550	BH
cyclohexane	10 µg/ml	10 ml	ref. P821560	BH
	100 µg/ml	1 ml	ref. P821570	BH
3,5-DINITROANILINE			CAS 618-87-1	Price
toluene	10 µg/ml	10 ml	ref. P821600	BG
	100 µg/ml	1 ml	ref. P821610	BG
1,2-DINITROBENZENE			CAS 528-29-0	Price
methanol	10 µg/ml	10 ml	ref. P821640	BG
	100 µg/ml	1 ml	ref. P821650	BG
1,3-DINITROBENZENE			CAS 99-65-0	Price
methanol	10 µg/ml	10 ml	ref. P821680	BG
	100 µg/ml	1 ml	ref. P821690	BG
2,4-DINITROPHENOL			CAS 51-28-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P821720	BG
	100 µg/ml	1 ml	ref. P821730	BG
methanol	10 µg/ml	10 ml	ref. P821740	BG
	100 µg/ml	1 ml	ref. P821750	BG
2,4-DINITROTOLUENE			CAS 121-14-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P821780	BG
	100 µg/ml	1 ml	ref. P821790	BG
methanol	10 µg/ml	10 ml	ref. P821800	BG
	100 µg/ml	1 ml	ref. P821810	BG
	100 µg/ml	5 ml	ref. P878220	CF
2,6-DINITROTOLUENE			CAS 606-20-2	Price
methanol	10 µg/ml	10 ml	ref. P821860	BG
	100 µg/ml	1 ml	ref. P821870	BG

3,4-DINITROTOLUENE			CAS 610-39-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P821900	BG
	100 µg/ml	1 ml	ref. P821910	BG
methanol	10 µg/ml	10 ml	ref. P821920	BG
	100 µg/ml	1 ml	ref. P821930	BG
DINOCAP			CAS 39300-45-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P821960	BG
	100 µg/ml	1 ml	ref. P821970	BG
DINOSEB			CAS 88-85-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P822000	BG
	100 µg/ml	1 ml	ref. P822010	BG
cyclohexane	10 µg/ml	10 ml	ref. P822020	BG
	100 µg/ml	1 ml	ref. P822030	BG
methanol	10 µg/ml	10 ml	ref. P822040	BG
	100 µg/ml	1 ml	ref. P822050	BG
DINOSEB METHYL ETHER			CAS 6099-79-2	Price
n-hexane	10 µg/ml	10 ml	ref. P822120	BH
	100 µg/ml	1 ml	ref. P822130	BH
methanol	10 µg/ml	10 ml	ref. P822140	BH
	100 µg/ml	1 ml	ref. P822150	BH
DINOTERB			CAS 1420-07-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P822180	BG
	100 µg/ml	1 ml	ref. P822190	BG
iso-octane	10 µg/ml	10 ml	ref. P822200	BG
	100 µg/ml	1 ml	ref. P822210	BG
DIOXACARB			CAS 6988-21-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P879870	BG
methanol	100 µg/ml	1 ml	ref. P879860	BG
DIOXAN			CAS 123-91-1	Price
dichloromethane	10 µg/ml	10 ml	ref. P822340	BG
	100 µg/ml	1 ml	ref. P822350	BG
methanol purge & trap	10 µg/ml	10 ml	ref. P822320	BG
	100 µg/ml	1 ml	ref. P822330	BG
1,4-DIOXANE-D8			CAS 17647-74-4	Price
methanol purge & trap	100 µg/ml	1 ml	ref. P822390	CA
DI-n-PENTHYL PHTHALATE-D4			id : 21290	Price
acetone	100 µg/ml	1 ml	ref. P822420	BCA
n-nonane	100 µg/ml	1 ml	ref. P822470	BCA
DIPHENAMID			CAS 957-51-7	Price
ethyl acetate	10 µg/ml	10 ml	ref. P875530	BG
methanol	100 µg/ml	1 ml	ref. P875520	BG
DIPHENYLAMINE			CAS 122-39-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P822500	BG
	100 µg/ml	1 ml	ref. P822510	BG
cyclohexane	10 µg/ml	10 ml	ref. P822520	BG
	100 µg/ml	1 ml	ref. P822530	BG
1,2-DIPHENYLHYDRAZINE			CAS 122-66-7	Price
methanol	10 µg/ml	10 ml	ref. P822560	BG
	100 µg/ml	1 ml	ref. P822570	BG
DI-n-PROPYL PHTHALATE-3,4,5,6-D4			CAS 358731-29-0	Price
acetone	100 µg/ml	1 ml	ref. P822600	DA
n-nonane	100 µg/ml	1 ml	ref. P822610	DA
DIQUAT DIBROMIDE HYDRATE			CAS 6385-62-2	Price
water	10 µg/ml	10 ml	ref. P822640	BG
	100 µg/ml	1 ml	ref. P822650	BG
DISULFOTON			CAS 298-04-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P822720	BG
	100 µg/ml	1 ml	ref. P822730	BG
cyclohexane	10 µg/ml	10 ml	ref. P822740	BG
	100 µg/ml	1 ml	ref. P822750	BG
DISULFOTON-SULFONE			CAS 2497-06-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P875550	BH
methanol	100 µg/ml	1 ml	ref. P875540	BH
DISULFOTON-SULFOXIDE			CAS 2497-07-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P875570	BH
methanol	100 µg/ml	1 ml	ref. P875560	BH
DITALIMFOS			CAS 5131-24-8	Price
n-hexane	100 µg/ml	5 ml	ref. P884250	BJ
cyclohexane	10 µg/ml	10 ml	ref. P884240	BJ

DITHIANON				CAS 3347-22-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P875580	BG	
DIURON				CAS 330-54-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P822820	BG	
	100 µg/ml	1 ml	ref. P822830	BG	
cyclohexane	10 µg/ml	10 ml	ref. P822840	BG	
	100 µg/ml	1 ml	ref. P822850	BG	
iso-octane	10 µg/ml	10 ml	ref. P822860	BG	
	100 µg/ml	1 ml	ref. P822870	BG	
methanol	10 µg/ml	10 ml	ref. P822880	BG	
	100 µg/ml	1 ml	ref. P822890	BG	
DIURON-D6				id :8150	Price
acetone	100 µg/ml	1 ml	ref. P822920	EJ	
nonane	100 µg/ml	1 ml	ref. P822930	EJ	
n-DOCOSANE				CAS 929-97-0	Price
dichloromethane	100 µg/ml	5 ml	ref. P877130	BG	
DOCOSANOIC ACID-METHYL ESTER				CAS 929-77-7	Price
n-hexane	10 000 µg/ml	1 ml	ref. P875600	CF	
n-DODECANE				CAS 112-40-3	Price
methanol	100 µg/ml	5 ml	ref. P877140	BG	
	5 000 µg/ml	1 ml	ref. P877150	CF	
DODECANOIC ACID-METHYL ESTER				CAS 111-82-0	Price
n-hexane	10 000 µg/ml	1 ml	ref. P875610	BJ	
EDIFENPHOS				CAS 17109-49-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P879910	BG	
methanol	100 µg/ml	1 ml	ref. P879900	BG	
n-EICOSANE				CAS 112-95-8	Price
dichloromethane	100 µg/ml	5 ml	ref. P877170	BG	
ELCOSANOIC ACID-METHYL ESTER				CAS 1120-28-1	Price
n-hexane	10 000 µg/ml	1 ml	ref. P875620	EC	
α-ENDOSULFAN				CAS 959-98-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P823060	BG	
	100 µg/ml	1 ml	ref. P823070	BG	
cyclohexane	10 µg/ml	10 ml	ref. P823080	BG	
	100 µg/ml	1 ml	ref. P823090	BG	
iso-octane	10 µg/ml	10 ml	ref. P823100	BG	
	100 µg/ml	1 ml	ref. P823110	BG	
methanol	10 µg/ml	10 ml	ref. P823120	BG	
	100 µg/ml	1 ml	ref. P823130	BG	
β-ENDOSULFAN				CAS 33213-65-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P823200	BG	
	100 µg/ml	1 ml	ref. P823210	BG	
cyclohexane	10 µg/ml	10 ml	ref. P823220	BG	
	100 µg/ml	1 ml	ref. P823230	BG	
iso-octane	10 µg/ml	10 ml	ref. P823240	BG	
	100 µg/ml	1 ml	ref. P823250	BG	
methanol	10 µg/ml	10 ml	ref. P823260	BG	
	100 µg/ml	1 ml	ref. P823270	BG	
(α+β)-ENDOSULFAN				CAS 115-29-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P823340	BG	
	100 µg/ml	1 ml	ref. P823350	BG	
cyclohexane	10 µg/ml	10 ml	ref. P823360	BG	
	100 µg/ml	1 ml	ref. P823370	BG	
methanol	10 µg/ml	10 ml	ref. P823380	BG	
	100 µg/ml	1 ml	ref. P823390	BG	
ENDOSULFAN-TOTAL (SULFATE)				CAS 1031-07-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P823420	BG	
	100 µg/ml	1 ml	ref. P823430	BG	
methanol	10 µg/ml	10 ml	ref. P823440	BG	
	100 µg/ml	1 ml	ref. P823450	BG	
ENDRIN				CAS 72-20-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P823480	BG	
	100 µg/ml	1 ml	ref. P823490	BG	
cyclohexane	10 µg/ml	10 ml	ref. P823500	BG	
	100 µg/ml	1 ml	ref. P823510	BG	
iso-octane	10 µg/ml	10 ml	ref. P823520	BG	
	100 µg/ml	1 ml	ref. P823530	BG	
methanol	10 µg/ml	10 ml	ref. P823540	BG	
	100 µg/ml	1 ml	ref. F113291	BG	

ENDRIN ALDEHYDE				CAS 7421-93-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P823580	CF	
	100 µg/ml	1 ml	ref. P823590	CF	
cyclohexane	10 µg/ml	10 ml	ref. P823600	CF	
	100 µg/ml	1 ml	ref. P823610	CF	
iso-octane	10 µg/ml	10 ml	ref. P823620	CF	
	100 µg/ml	1 ml	ref. P823630	CF	
methanol	10 µg/ml	10 ml	ref. P823640	CF	
	100 µg/ml	1 ml	ref. P823650	CF	
ENDRINE KETONE				CAS 53494-70-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P823720	CE	
	100 µg/ml	1 ml	ref. P823730	CE	
methanol	10 µg/ml	10 ml	ref. P823740	CE	
	100 µg/ml	1 ml	ref. F113301	CE	
EPICHLORHYDRIN				CAS 106-89-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P823780	BG	
	100 µg/ml	1 ml	ref. P823790	BG	
EPN				CAS 2104-64-5	Price
acetonitrile	100 µg/ml	1 ml	ref. P879920	BG	
cyclohexane	10 µg/ml	10 ml	ref. P879940	BG	
	100 µg/ml	1 ml	ref. P879930	BG	
EPOXICONAZOLE				CAS 106325-08-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P823820	BG	
	100 µg/ml	1 ml	ref. P823830	BG	
EPTC				CAS 759-94-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P875640	BG	
methanol	100 µg/ml	1 ml	ref. P875630	BG	
ESFENVALERATE				CAS 66230-04-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P823900	BG	
	100 µg/ml	1 ml	ref. P823910	BG	
cyclohexane	10 µg/ml	10 ml	ref. P823920	BG	
	100 µg/ml	1 ml	ref. P823930	BG	
ESTRIOL				CAS 50-27-1	Price
acetonitrile	100 µg/ml	1 ml	ref. P875650	BG	
ESTRONE				CAS 53-16-7	Price
acetonitrile	100 µg/ml	1 ml	ref. P875660	BG	
ETACONAZOLE				CAS 6027-93-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P879960	BG	
methanol	100 µg/ml	1 ml	ref. P879950	BG	
ETHALFLURALIN				CAS 55283-68-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P879980	BG	
methanol	100 µg/ml	1 ml	ref. P879970	BG	
ETHANOL				CAS 64-17-5	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P823980	BG	
	100 µg/ml	1 ml	ref. P823990	BG	
water	10 µg/ml	10 ml	ref. P823960	BG	
	100 µg/ml	1 ml	ref. P823970	BG	
ETHANOLAMINE				CAS 141-43-5	Price
methanol	100 µg/ml	1 ml	ref. P875670	BG	
ETHIDIMURON				CAS 30043-49-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P824060	BG	
	100 µg/ml	1 ml	ref. P824070	BG	
ETHIOFENCARB				CAS 29973-13-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P875680	BG	
	100 µg/ml	1 ml	ref. P875690	BG	
ETHION				CAS 563-12-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P824100	BG	
	100 µg/ml	1 ml	ref. P824110	BG	
cyclohexane	10 µg/ml	10 ml	ref. P824120	BG	
	100 µg/ml	1 ml	ref. P824130	BG	
methanol	10 µg/ml	10 ml	ref. P824140	BG	
	100 µg/ml	1 ml	ref. P824150	BG	
ETHOFUMESATE				CAS 26225-79-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P824220	BG	
	100 µg/ml	1 ml	ref. P824230	BG	
methanol	10 µg/ml	10 ml	ref. P824240	BG	
	100 µg/ml	1 ml	ref. P824250	BG	

ETHOPROPHOS			CAS 13194-48-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P824280	BG
	100 µg/ml	1 ml	ref. P824290	BG
cyclohexane	10 µg/ml	10 ml	ref. P824300	BG
	100 µg/ml	1 ml	ref. P824310	BG
ETHOXYQUIN			CAS 91-53-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P824340	BG
	100 µg/ml	1 ml	ref. P824350	BG
methanol	10 µg/ml	10 ml	ref. P824360	BG
	100 µg/ml	1 ml	ref. P824370	BG
ETHYLACETATE			CAS 141-78-6	Price
ethanol	10 µg/ml	10 ml	ref. P824400	BG
	100 µg/ml	1 ml	ref. P824410	BG
ETHYLAMINE			CAS 75-04-7	Price
methanol	10 000 µg/ml	5 ml	ref. P875720	BJ
ETHYLBENZENE			CAS 100-41-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P824440	BG
	100 µg/ml	1 ml	ref. P824450	BG
methanol purge & trap	10 µg/ml	10 ml	ref. P824460	BG
	100 µg/ml	1 ml	ref. P824470	BG
ETHYLBENZENE-D10			CAS 25837-05-2	Price
methanol purge & trap	100 µg/ml	1 ml	ref. P824530	BK
	2 000 µg/ml	1 ml	ref. F062445	CG
ETHYL-tert-BUTYLETHER (ETBE)			CAS 637-92-3	Price
methanol	100 µg/ml	5 ml	ref. P875740	BH
ETHYLENEGLYCOL			CAS 107-21-1	Price
methanol purge & trap	100 µg/ml	1 ml	ref. P824610	BH
	1 000 µg/ml	1 ml	ref. P824600	CA
ETHYLENEGLYCOL-MONOMETHYL ETHER			CAS 109-86-4	Price
methanol	100 µg/ml	1 ml	ref. P875730	BG
ETHYLENE OXIDE			CAS 75-21-8	Price
dimethylsulphoxide	1000 µg/ml	1 ml	ref. P824640	BH
methanol purge & trap	1000 µg/ml	1 ml	ref. P824650	BH
toluene	1000 µg/ml	1 ml	ref. P824660	BH
2-ETHYLHEXYL 2-HYDROXYBENZOATE			CAS 118-60-5	Price
cyclohexane	100 µg/ml	1 ml	ref. P873060	BG
ETHYL MERCAPTAN			CAS 75-08-1	Price
toluene	2 000 µg/ml	1 ml	ref. P875700	BH
ETHYL OCTANATE			CAS 106-32-1	Price
n-hexane	10 000 µg/ml	1 ml	ref. P875710	CD
2-ETHYLPHENOL			CAS 90-00-6	Price
methanol	10 µg/ml	10 ml	ref. P873080	BG
	100 µg/ml	1 ml	ref. P873070	BG
3-ETHYLPHENOL			CAS 620-17-7	Price
methanol	100 µg/ml	1 ml	ref. P873410	BG
4-ETHYLPHENOL			CAS 123-07-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P873580	BG
cyclohexane	100 µg/ml	1 ml	ref. P873570	BG
m-ETHYLTOLUENE			CAS 620-14-4	Price
methanol	100 µg/ml	1 ml	ref. P877020	BG
o-ETHYLTOLUENE			CAS 611-14-3	Price
methanol	100 µg/ml	1 ml	ref. P877420	BK
p-ETHYLTOLUENE			CAS 622-96-8	Price
methanol	100 µg/ml	1 ml	ref. P877790	BK
ETRIDIAZOLE			CAS 2593-15-9	Price
cyclohexane	100 µg/ml	1 ml	ref. P875750	BG
ETRIMFOS			CAS 38260-54-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P824750	BG
	100 µg/ml	1 ml	ref. P824760	BG
EUGENOL			CAS 97-53-0	Price
ethanol	1 000 µg/ml	1 ml	ref. P875770	BJ
FAMOXADONE			CAS 131807-57-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P875790	BH
	100 µg/ml	1 ml	ref. P875780	BH
FAMPHUR			CAS 52-85-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P880000	BG
methanol	100 µg/ml	1 ml	ref. P879990	BG

FENAMIDONE			CAS 161326-34-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P875820	BG
	100 µg/ml	1 ml	ref. P875830	BG
cyclohexane	10 µg/ml	10 ml	ref. P875800	BG
	100 µg/ml	1 ml	ref. P875810	BG
FENAMINOSULF			CAS 140-56-7	Price
methanol	100 µg/ml	1 ml	ref. P880020	BH
	1 000 µg/ml	1 ml	ref. P880010	CA
FENAMIPHOS			CAS 22224-92-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P875370	BG
	100 µg/ml	1 ml	ref. P824800	BG
cyclohexane	10 µg/ml	10 ml	ref. P824810	BG
	100 µg/ml	1 ml	ref. P824820	BG
FENAMIPHOS-SULFONE			CAS 31972-44-8	Price
ethyl acetate	10 µg/ml	10 ml	ref. P875840	BH
FENAMIPHOS-SULFOXIDE			CAS 31972-43-7	Price
methyl-tert-butylether	10 µg/ml	10 ml	ref. P875850	BH
FENARIMOL			CAS 60168-88-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P824850	BG
	100 µg/ml	1 ml	ref. P824860	BG
FENAZAFLOR			CAS 14255-88-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P824930	BK
	100 µg/ml	1 ml	ref. P824940	BK
FENAZAQUIN			CAS 120928-09-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P875880	BH
cyclohexane	10 µg/ml	10 ml	ref. P875860	BH
	100 µg/ml	1 ml	ref. P875870	BH
FENBUCONAZOLE			CAS 114369-43-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P824970	BH
	100 µg/ml	1 ml	ref. P824980	BH
methanol	10 µg/ml	10 ml	ref. P824990	BH
	100 µg/ml	1 ml	ref. P825000	BH
FENBUTATIN OXYDE			CAS 13356-08-6	Price
acetone	100 µg/ml	1 ml	ref. P875900	BG
cyclohexane	10 µg/ml	10 ml	ref. P875890	BG
FENCHLORPHOS			CAS 299-84-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P825070	BG
	100 µg/ml	1 ml	ref. P825080	BG
cyclohexane	10 µg/ml	10 ml	ref. P825090	BG
	100 µg/ml	1 ml	ref. P825100	BG
methanol	10 µg/ml	10 ml	ref. P825110	BG
	100 µg/ml	1 ml	ref. P825120	BG
FENHEXAMID			CAS 126833-17-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P875910	BG
methanol	100 µg/ml	1 ml	ref. P875920	BG
FENITROTHION			CAS 122-14-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P825190	BG
	100 µg/ml	1 ml	ref. P825200	BG
cyclohexane	10 µg/ml	10 ml	ref. P825210	BG
	100 µg/ml	1 ml	ref. P825220	BG
FENOBUARB			CAS 3766-81-2	Price
acetonitrile	100 µg/ml	1 ml	ref. P875970	BG
cyclohexane	10 µg/ml	10 ml	ref. P875980	BG
FENOPROP			CAS 93-72-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P854490	BG
	100 µg/ml	1 ml	ref. P854500	BG
ethyl acetate	10 µg/ml	10 ml	ref. P854510	BG
	100 µg/ml	1 ml	ref. P854520	BG
iso-octane	10 µg/ml	10 ml	ref. P854530	BG
	100 µg/ml	1 ml	ref. P854540	BG
methanol	10 µg/ml	10 ml	ref. P854550	BG
	100 µg/ml	1 ml	ref. P854560	BG
FENOPROP-METHYL ESTER			CAS 4841-20-7	Price
methanol	10 µg/ml	10 ml	ref. P854590	BG
	100 µg/ml	1 ml	ref. P854600	BG
FENOXAPROP			CAS 95617-09-7	Price
methanol	100 µg/ml	1 ml	ref. P875990	BG
FENOXAPROP-P-ETHYL			CAS 71283-80-2	Price
methanol	100 µg/ml	1 ml	ref. P876000	BG

FENOXYCARB			CAS 79127-80-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P825290	BH
	100 µg/ml	1 ml	ref. P825300	BH
FENPICLONIL			CAS 74738-17-3	Price
ethyl acetate	10 µg/ml	10 ml	ref. P880030	BG
FENPROPATHRIN			CAS 64257-84-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P876020	BG
iso-octane	10 µg/ml	10 ml	ref. P876030	BG
	100 µg/ml	1 ml	ref. P876010	BG
FENPROPIDIN			CAS 67306-00-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P825330	BG
	100 µg/ml	1 ml	ref. P825340	BG
FENPROIMORTH			CAS 67564-91-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P825370	BG
	100 µg/ml	1 ml	ref. P825380	BG
cyclohexane	10 µg/ml	10 ml	ref. P825390	BG
	100 µg/ml	1 ml	ref. P825400	BG
FENSON			CAS 80-38-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P876040	BG
methanol	100 µg/ml	1 ml	ref. P876050	BG
FENSULFOTHION			CAS 115-90-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P825430	BG
	100 µg/ml	1 ml	ref. P825440	BG
cyclohexane	10 µg/ml	10 ml	ref. P825450	BG
	100 µg/ml	1 ml	ref. P825460	BG
FENTHION			CAS 55-38-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P825490	BG
	100 µg/ml	1 ml	ref. P825500	BG
cyclohexane	10 µg/ml	10 ml	ref. P825510	BG
	100 µg/ml	1 ml	ref. P825520	BG
FENTHION-SULFONE			CAS 3761-42-0	Price
ethyl acetate	10 µg/ml	1 ml	ref. P876070	BH
	100 µg/ml	1 ml	ref. P876060	BH
FENTHION-SULFOXIDE			CAS 3761-41-9	Price
acetone	100 µg/ml	1 ml	ref. P876090	BH
acetonitrile	100 µg/ml	1 ml	ref. P876080	BH
FENTIN-ACETATE			CAS 900-95-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P825580	BG
	100 µg/ml	1 ml	ref. P825590	BG
FENTIN-CHLORIDE			CAS 639-58-7	Price
iso-octane	10 µg/ml	10 ml	ref. P825620	BG
	100 µg/ml	1 ml	ref. P825630	BG
FENURON			CAS 101-42-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P876120	BH
	100 µg/ml	1 ml	ref. P876110	BH
	1 000 µg/ml	5 ml	ref. P876100	CF
FENVALERATE			CAS 51630-58-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P825700	BG
	100 µg/ml	1 ml	ref. P825710	BG
iso-octane	10 µg/ml	10 ml	ref. P825720	BG
	100 µg/ml	1 ml	ref. P825730	BG
FERBAM			CAS 14484-64-1	Price
methanol	100 µg/ml	1 ml	ref. P876130	BG
FIPRONIL			CAS 120068-37-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P876150	BH
	100 µg/ml	1 ml	ref. P876140	BH
cyclohexane	10 µg/ml	10 ml	ref. P876160	BH
FLAMPROP-ISOPROPYL			CAS 52756-22-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P876170	BG
methanol	100 µg/ml	1 ml	ref. P876180	BG
FLAZASULFURON			CAS 104040-78-0	Price
acetone	10 µg/ml	10 ml	ref. P825840	BH
	100 µg/ml	1 ml	ref. P825850	BH
FLUAZIFOP			CAS 69335-91-7	Price
acetonitrile	100 µg/ml	1 ml	ref. P876190	BG
ethyl acetate	10 µg/ml	10 ml	ref. P876200	BG
FLUAZIFOP-BUTYL			CAS 69806-50-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P876210	BG
	100 µg/ml	1 ml	ref. P876220	BG

FLUAZIFOP-METYL			CAS 69335-90-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P876230	BJ
FLUAZINAM			CAS 79622-59-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P876240	BH
methanol	100 µg/ml	1 ml	ref. P876250	BH
FLUCHLORALIN			CAS 33245-39-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P880050	BG
methanol	100 µg/ml	1 ml	ref. P880040	BG
FLUCYTHRINATE			CAS 70124-77-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P876260	BH
iso-octane	100 µg/ml	1 ml	ref. P876270	BH
FLUFENACET			CAS 142459-58-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P876280	BG
cyclohexane	10 µg/ml	10 ml	ref. P876290	BG
FLUFENOXYURON			CAS 101463-69-8	Price
acetone	10 µg/ml	10 ml	ref. P825910	BG
	100 µg/ml	1 ml	ref. P825920	BG
cyclohexane	10 µg/ml	10 ml	ref. P825930	BG
	100 µg/ml	1 ml	ref. P825940	BG
FLUIDIOXONIL			CAS 131341-86-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P825970	BG
	100 µg/ml	1 ml	ref. P825980	BG
methanol	10 µg/ml	10 ml	ref. P825990	BG
	100 µg/ml	1 ml	ref. P826000	BG
FLUMETRALIN			CAS 62924-70-3	Price
acetone	100 µg/ml	5 ml	ref. P884970	BH
acetonitrile	100 µg/ml	1 ml	ref. P884980	BH
cyclohexane	10 µg/ml	10 ml	ref. P884990	BH
FLUMIOXAZIN			CAS 103361-09-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P826030	BJ
	100 µg/ml	1 ml	ref. P826040	BJ
cyclohexane	10 µg/ml	10 ml	ref. P826050	BJ
	100 µg/ml	1 ml	ref. P826060	BJ
FLUOMETURON			CAS 2164-17-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P876320	BG
	100 µg/ml	1 ml	ref. P876310	BG
	1 000 µg/ml	5 ml	ref. P876300	CE
FLUORANTHENE			CAS 206-44-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P826150	BH
	100 µg/ml	1 ml	ref. P826160	BH
cyclohexane	10 µg/ml	10 ml	ref. P826170	BH
	100 µg/ml	1 ml	ref. P826180	BH
methanol	10 µg/ml	10 ml	ref. P826190	BH
	100 µg/ml	1 ml	ref. P826200	BH
FLUORANTHENE-D10			CAS 93951-69-0	Price
acetone	100 µg/ml	1 ml	ref. P826230	DD
acetonitrile	100 µg/ml	1 ml	ref. P826240	DD
methanol purge & trap	100 µg/ml	1 ml	ref. P826250	DD
FLUORENE			CAS 86-73-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P826280	BH
	100 µg/ml	1 ml	ref. P826290	BH
cyclohexane	10 µg/ml	10 ml	ref. P826300	BH
	100 µg/ml	1 ml	ref. P826310	BH
methanol	10 µg/ml	10 ml	ref. P826320	BH
	100 µg/ml	1 ml	ref. P826330	BH
FLUORENE-D10			CAS 81103-79-9	Price
cyclohexane	100 µg/ml	1 ml	ref. P826360	CE
FLUOROBENZENE			CAS 462-06-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P826390	BG
	100 µg/ml	1 ml	ref. P826400	BG
	2 000 µg/ml	1 ml	ref. F112921	DA
2-FLUOROBIPHENYL			CAS 321-60-8	Price
dichloromethane	10 µg/ml	10 ml	ref. P826430	BG
	100 µg/ml	1 ml	ref. P826440	BG
FLUOROGLYCOFEN-ETHYL			CAS 77501-90-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P876330	BJ
1-FLUORONAPHTHALENE			CAS 321-38-0	Price
dichloromethane	10 µg/ml	10 ml	ref. P826470	BG
	100 µg/ml	1 ml	ref. P826480	BG

2-FLUOROPHENOL			CAS 367-12-4	Price
acetone	100 µg/ml	1 ml	ref. F146761	BH
dichloromethane	10 µg/ml	10 ml	ref. P826510	BH
	100 µg/ml	1 ml	ref. P826520	BH
methanol	2 000 µg/ml	1 ml	ref. F116591	CC
FLUOROTRICHLOROMETHANE			CAS 75-69-4	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P855980	CB
	100 µg/ml	1 ml	ref. P855990	CB
FLUQUINCONAZOLE			CAS 136426-54-5	Price
iso-octane	10 µg/ml	10 ml	ref. P826590	BH
	100 µg/ml	1 ml	ref. P826600	BH
FLURIDONE			CAS 59756-60-4	Price
methanol	100 µg/ml	1 ml	ref. P876340	BJ
FLUROCHLORIDONE			CAS 61213-25-0	Price
acetone	10 µg/ml	10 ml	ref. P826630	BG
	100 µg/ml	1 ml	ref. P826640	BG
acetonitrile	10 µg/ml	10 ml	ref. P826650	BG
	100 µg/ml	1 ml	ref. P826660	BG
cyclohexane	10 µg/ml	10 ml	ref. P826670	BG
	100 µg/ml	1 ml	ref. P826680	BG
FLUROXYPYR			CAS 69377-81-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P826710	BG
	100 µg/ml	1 ml	ref. P826720	BG
FLUROXYPYR-1- METHYLHEPTYL ESTER			CAS 81406-37-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P826750	BG
	100 µg/ml	1 ml	ref. P826760	BG
FLURPRIMIDOL			CAS 56425-91-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P876350	CE
FLURTAMONE			CAS 96525-23-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P876360	CC
FLUSILAZOLE			CAS 85509-19-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P875300	BG
	100 µg/ml	1 ml	ref. P826800	BG
FLUTOLANIL			CAS 66332-96-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P876370	CB
methanol	100 µg/ml	1 ml	ref. P876380	CB
FLUTRIAFOL			CAS 76674-21-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P826850	BG
	100 µg/ml	1 ml	ref. P826860	BG
toluene	10 µg/ml	10 ml	ref. P826870	BG
	100 µg/ml	1 ml	ref. P826880	BG
FOLPET			CAS 133-07-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P826910	BG
	100 µg/ml	1 ml	ref. P826920	BG
FOMESAFEN			CAS 72178-02-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P876390	BG
FONOFOS			CAS 994-22-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P826980	BG
	100 µg/ml	1 ml	ref. P826990	BG
cyclohexane	10 µg/ml	10 ml	ref. P827000	BG
	100 µg/ml	1 ml	ref. P827010	BG
FORAMSULFURON			CAS 173159-57-4	Price
acetonitrile	100 µg/ml	1 ml	ref. P876400	BG
FORMALDEHYDE			CAS 50-00-0	Price
water	10 µg/ml	10 ml	ref. P827080	BG
	100 µg/ml	1 ml	ref. P827090	BG
FORMALDEHYDE-DNPH			CAS 1081-15-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P827040	CC
	100 µg/ml	1 ml	ref. P827050	CC
FORMOTHION			CAS 2540-82-1	Price
acetone	10 µg/ml	10 ml	ref. P827120	BG
	100 µg/ml	1 ml	ref. P827130	BG
acetonitrile	10 µg/ml	10 ml	ref. P827140	BG
	100 µg/ml	1 ml	ref. P827150	BG

FOSTHIAZATE			CAS 98886-44-3	Price
acetone	10 µg/ml	10 ml	ref. P827180	CC
	100 µg/ml	1 ml	ref. P827190	CC
acetonitrile	10 µg/ml	10 ml	ref. P827200	CC
	100 µg/ml	1 ml	ref. P827210	CC
cyclohexane	10 µg/ml	10 ml	ref. P827220	CC
	100 µg/ml	1 ml	ref. P827230	CC
FUBERIDAZOLE			CAS 3878-19-1	Price
methanol	10 µg/ml	10 ml	ref. P880060	BG
FURALAXYL			CAS 57646-30-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P876410	BG
methanol	100 µg/ml	1 ml	ref. P876420	BG
FURALDEHYDE			CAS 98-01-1	Price
acetonitrile	1 000 µg/ml	1 ml	ref. P876430	BH
FURATHIOCARB			CAS 65907-30-4	Price
acetonitrile	100 µg/ml	1 ml	ref. P876440	BG
iso-octane	10 µg/ml	10 ml	ref. P876450	BG
GLUPHOSINATE AMMONIUM			CAS 77182-82-2	Price
water	10 µg/ml	10 ml	ref. P876470	BG
	100 µg/ml	1 ml	ref. P876460	BG
GLYPHOSATE			CAS 1071-83-6	Price
water	10 µg/ml	10 ml	ref. P827290	BG
	100 µg/ml	1 ml	ref. F114201	BG
HALOXYFOP			CAS 69806-34-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P876490	BG
	100 µg/ml	1 ml	ref. P876480	BG
HALOXYFOP-2-ETHOXYETHYL			CAS 87237-48-7	Price
acetone	10 µg/ml	10 ml	ref. P827390	BG
	100 µg/ml	1 ml	ref. P827400	BG
iso-octane	10 µg/ml	10 ml	ref. P827410	BG
	100 µg/ml	1 ml	ref. P827420	BG
α-HCH			CAS 319-84-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P816270	BG
	100 µg/ml	1 ml	ref. P816280	BG
methanol	10 µg/ml	10 ml	ref. P827620	BG
	100 µg/ml	1 ml	ref. P827480	BG
α-HCH-D6			CAS 86194-41-4	Price
cyclohexane	100 µg/ml	1 ml	ref. P827510	BFE
β-HCH			CAS 319-85-7	Price
acetone	10 µg/ml	10 ml	ref. P827540	BG
	100 µg/ml	1 ml	ref. P827550	BG
cyclohexane	10 µg/ml	10 ml	ref. P827560	BG
	100 µg/ml	1 ml	ref. P827570	BG
toluene	10 µg/ml	10 ml	ref. P827580	BG
	100 µg/ml	1 ml	ref. P827590	BG
δ-HCH			CAS 319-86-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P875180	BH
	100 µg/ml	1 ml	ref. P875170	BH
γ-HCH (lindane)			CAS 58-89-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P827690	BG
	100 µg/ml	1 ml	ref. P827700	BG
methanol	10 µg/ml	10 ml	ref. P827710	BG
	100 µg/ml	1 ml	ref. P827720	BG
HEPTACHLOR			CAS 76-44-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P827790	BH
	100 µg/ml	1 ml	ref. P827800	BH
cyclohexane	10 µg/ml	10 ml	ref. P827810	BH
	100 µg/ml	1 ml	ref. P827820	BH
methanol	10 µg/ml	10 ml	ref. P827830	BH
	100 µg/ml	1 ml	ref. P827840	BH
HEPTACHLOR -ENDO-EPOXIDE			CAS 28044-83-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P827910	CA
	100 µg/ml	1 ml	ref. P827920	CA
methanol	10 µg/ml	10 ml	ref. P827930	CA
	100 µg/ml	1 ml	ref. P827940	CA
HEPTACHLOR -EXO-EPOXIDE			CAS 1024-57-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P827970	CA
	100 µg/ml	1 ml	ref. P827980	CA
methanol	10 µg/ml	10 ml	ref. P827990	CA
	100 µg/ml	1 ml	ref. P828000	CA

n-HEPTADECANE			CAS 629-78-7	Price
ethylbenzene	100 000 µg/ml	2 ml	ref. P877190	GC
1-HEPTANAL			CAS 111-71-7	Price
acetonitrile	1 000 µg/ml	5 ml	ref. P872440	BG
n-HEPTANE			CAS 142-82-5	Price
methanol	100 µg/ml	5 ml	ref. P877200	BG
2-HEPTANONE			CAS 110-43-0	Price
methanol	1 000 µg/ml	5 ml	ref. P873090	BJ
3-HEPTANONE			CAS 106-35-4	Price
methanol/water (9/1)	10 µg/ml	10 ml	ref. P828100	BG
	100 µg/ml	1 ml	ref. P828110	BG
4-HEPTANONE			CAS 123-19-3	Price
methanol/water (9/1)	10 µg/ml	10 ml	ref. P828140	BG
	100 µg/ml	1 ml	ref. P828150	BG
HEPTENOPHOS			CAS 23560-59-0	Price
acetonitrile	100 µg/ml	1 ml	ref. P876520	BG
cyclohexane	10 µg/ml	10 ml	ref. P876530	BG
HEXABROMOBENZENE			CAS 87-82-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P876550	BG
	100 µg/ml	1 ml	ref. P876540	
HEXACHLOROBENZENE			CAS 118-74-1	Price
acetone	10 µg/ml	10 ml	ref. P828180	BG
	100 µg/ml	1 ml	ref. P828190	BG
acetonitrile	10 µg/ml	10 ml	ref. P828200	BG
	100 µg/ml	1 ml	ref. P828210	BG
cyclohexane	10 µg/ml	10 ml	ref. P828220	BG
	100 µg/ml	1 ml	ref. P828230	BG
iso-octane	10 µg/ml	10 ml	ref. P828240	BG
	100 µg/ml	1 ml	ref. P828250	BG
	1 000 µg/ml	10 ml	ref. F117491	BK
methanol	10 µg/ml	10 ml	ref. P828260	BG
	100 µg/ml	1 ml	ref. P828270	BG
HEXACHLOROBENZENE-13C6			CAS 93952-14-8	Price
acetone	100 µg/ml	1 ml	ref. P828300	BAF
n-nonane	100 µg/ml	1 ml	ref. P828310	BAF
HEXACHLORO-1,3-BUTADIENE			CAS 87-68-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P828340	BG
	100 µg/ml	1 ml	ref. P828350	BG
methanol	10 µg/ml	10 ml	ref. P828360	BG
	100 µg/ml	1 ml	ref. P828370	BG
HEXACHLOROCYCLOHEXANE (mix of isome)			CAS 608-73-1	Price
iso-octane	10 µg/ml	10 ml	ref. P876500	BG
methanol	100 µg/ml	1 ml	ref. P876510	BG
HEXACHLOROCYCLOPENTADIENE			CAS 77-47-4	Price
iso-octane	10 µg/ml	10 ml	ref. P828460	BG
	100 µg/ml	1 ml	ref. P828470	BG
methanol	10 µg/ml	10 ml	ref. P828480	BG
	100 µg/ml	1 ml	ref. P828490	BG
HEXACHLOROETHANE			CAS 67-72-1	Price
iso-octane	10 µg/ml	10 ml	ref. P828520	BG
	100 µg/ml	1 ml	ref. P828530	BG
methanol	10 µg/ml	10 ml	ref. P828540	BG
	100 µg/ml	1 ml	ref. P828550	BG
HEXACHLOROPHEN			CAS 70-30-4	Price
dichloromethane	10 µg/ml	10 ml	ref. P828600	BG
	100 µg/ml	1 ml	ref. P828610	BG
methanol	10 µg/ml	10 ml	ref. P828580	BG
	100 µg/ml	1 ml	ref. P828590	BG
HEXACONAZOLE			CAS 79983-71-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P828640	BG
	100 µg/ml	1 ml	ref. P828650	BG
cyclohexane	10 µg/ml	10 ml	ref. P828660	BG
	100 µg/ml	1 ml	ref. P828670	BG
iso-octane	10 µg/ml	10 ml	ref. P828680	BG
	100 µg/ml	1 ml	ref. P828690	BG
n-HEXACOSANE			CAS 630-01-3	Price
dichloromethane	100 µg/ml	5 ml	ref. P877210	CC

n-HEXADECANE			CAS 544-76-3	Price
acetone	10 µg/ml	10 ml	ref. P828720	BG
	100 µg/ml	1 ml	ref. P828730	BG
dichloromethane	100 µg/ml	5 ml	ref. P877220	CF
HEXAFLUMURON			CAS 86479-06-3	Price
acetone	10 µg/ml	10 ml	ref. P828760	BH
	100 µg/ml	1 ml	ref. P828770	BH
methyl-tert-butylether	10 µg/ml	10 ml	ref. P828780	BH
	100 µg/ml	1 ml	ref. P828790	BH
HEXALDEHYDE-DNPH			CAS 1527-95-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P828820	BH
	100 µg/ml	1 ml	ref. P828830	BH
2-HEXANONE			CAS 591-78-6	Price
methanol/water (9/1)	100 µg/ml	5 ml	ref. P873100	BG
HEXAZINONE			CAS 51235-04-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P828890	BG
	100 µg/ml	1 ml	ref. P828900	BG
cyclohexane	10 µg/ml	10 ml	ref. P828910	BG
	100 µg/ml	1 ml	ref. P828920	BG
HEXYTHIAZOX			CAS 78587-05-0	Price
acetonitrile	100 µg/ml	1 ml	ref. P876560	BG
toluene	10 µg/ml	10 ml	ref. P876570	BG
4-HYDROXYBENZOIC ACID-n-BUTYL ESTER			CAS 94-26-8	Price
methanol	100 µg/ml	1 ml	ref. P873620	BG
4-HYDROXYBENZOIC ACID ETHYL ESTER			CAS 120-47-8	Price
methanol	100 µg/ml	1 ml	ref. P873590	BG
4-HYDROXYBENZOIC ACID-ISOBUTYL ESTER			CAS 4274-02-3	Price
methanol	100 µg/ml	1 ml	ref. P873600	BG
4-HYDROXYBENZOIC ACID METHYL ESTER			CAS 99-76-3	Price
methanol	100 µg/ml	1 ml	ref. P873610	BG
4-HYDROXYBENZOIC ACID PROPYL ESTER			CAS 94-13-3	Price
methanol	100 µg/ml	1 ml	ref. P873630	BG
IMAZALIL			CAS 35554-44-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P828950	BG
	100 µg/ml	1 ml	ref. P828960	BG
cyclohexane	10 µg/ml	10 ml	ref. P828970	BG
	100 µg/ml	1 ml	ref. P828980	BG
IMAZAMETHABENZ-METHYL			CAS 81405-85-8	Price
acetone	10 µg/ml	10 ml	ref. P829090	BG
	100 µg/ml	1 ml	ref. P829100	BG
acetonitrile	10 µg/ml	10 ml	ref. P829110	BG
	100 µg/ml	1 ml	ref. P829120	BG
cyclohexane	10 µg/ml	10 ml	ref. P829130	BG
	100 µg/ml	1 ml	ref. P829140	BG
IMAZAPYR			CAS 81334-34-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P876580	BG
	100 µg/ml	1 ml	ref. P876590	BG
IMAZAQUIN			CAS 81335-37-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P880080	BG
methanol	100 µg/ml	1 ml	ref. P880070	BG
IMAZETHAPYR			CAS 81335-77-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P880100	BG
methanol	100 µg/ml	1 ml	ref. P880090	BG
IMIDACLOPRID			CAS 105827-78-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P829170	BG
	100 µg/ml	1 ml	ref. P829180	BG
INDAN			CAS 496-11-7	Price
methanol	10 µg/ml	10 ml	ref. P876610	BH
	100 µg/ml	5 ml	ref. P876600	BH
INDENE			CAS 95-13-6	Price
methanol	10 µg/ml	10 ml	ref. P876620	BG
toluene	50 µg/ml	1 ml	ref. P876630	BG
	100 µg/ml	5 ml	ref. P876640	BG
INDENO(1,2,3-c,d)PYRENE			CAS 193-39-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P829240	CE
	100 µg/ml	1 ml	ref. P829250	CE
cyclohexane	10 µg/ml	10 ml	ref. P829260	CE
	100 µg/ml	1 ml	ref. P829270	CE
dichloromethane	10 µg/ml	10 ml	ref. P829280	CE
	100 µg/ml	1 ml	ref. P829290	CE

IODOCARB			CAS 55406-53-6	Price
acetone	10 µg/ml	10 ml	ref. P829390	BH
	100 µg/ml	1 ml	ref. P829400	BH
IODOFENPHOS			CAS 18181-70-9	Price
cyclohexane	100 µg/ml	1 ml	ref. P876650	BG
iso-octane	10 µg/ml	10 ml	ref. P876660	BG
IODOMETHANE			CAS 74-88-4	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P829350	BG
	100 µg/ml	1 ml	ref. P829360	BG
IODOSULFURON-METHYL-SODIUM			CAS 144550-36-7	Price
acetone	10 µg/ml	10 ml	ref. P829430	BH
	100 µg/ml	1 ml	ref. P829440	BH
IOXYNIL			CAS 1689-83-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P829470	BG
	100 µg/ml	1 ml	ref. P829480	BG
cyclohexane	10 µg/ml	10 ml	ref. P829490	BG
	100 µg/ml	1 ml	ref. P829500	BG
IPIRODIONE			CAS 36734-19-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P829530	BG
	100 µg/ml	1 ml	ref. P829540	BG
cyclohexane	10 µg/ml	10 ml	ref. P829550	BG
	100 µg/ml	1 ml	ref. P829560	BG
IPROVALICARB			CAS 140923-17-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P876690	BG
	100 µg/ml	1 ml	ref. P876710	BG
cyclohexane	10 µg/ml	10 ml	ref. P876700	BG
IRGAROL 1051			CAS 28159-98-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P876720	BH
methanol	100 µg/ml	1 ml	ref. P876730	BH
ISAZOFOS			CAS 42509-80-8	Price
cyclohexane	100 µg/ml	1 ml	ref. P876740	BG
iso-octane	10 µg/ml	10 ml	ref. P876750	BG
ISOENZAN			CAS 297-78-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P829630	BK
	100 µg/ml	1 ml	ref. P829640	BK
methanol	10 µg/ml	10 ml	ref. P829650	BK
	100 µg/ml	1 ml	ref. P829660	BK
ISODRIN			CAS 465-73-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P829690	BG
	100 µg/ml	1 ml	ref. P829700	BG
iso-octane	10 µg/ml	10 ml	ref. P829710	BG
	100 µg/ml	1 ml	ref. P829720	BG
methanol	10 µg/ml	10 ml	ref. P829730	BG
	100 µg/ml	1 ml	ref. P829740	BG
ISOFENPHOS			CAS 25311-71-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P876770	BJ
	100 µg/ml	1 ml	ref. P876760	BJ
ISOFENPHOS-METHYL			CAS 99675-03-3	Price
acetonitrile	100 µg/ml	1 ml	ref. P879120	BJ
ISOPHORONE			CAS 78-59-1	Price
methanol	100 µg/ml	5 ml	ref. P876780	BG
	1 000 µg/ml	1 ml	ref. P876790	BK
ISOPROCARB			CAS 2631-40-5	Price
methanol	10 µg/ml	10 ml	ref. P880120	BG
	100 µg/ml	1 ml	ref. P880110	BG
ISOPROPALIN			CAS 33820-53-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P880140	BG
methanol	100 µg/ml	1 ml	ref. P880130	BG
ISOPROPYLBENZENE			CAS 98-82-8	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P829770	BG
	100 µg/ml	1 ml	ref. P829780	BG
4-ISOPROPYLTOLUENE			CAS 99-87-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P814270	BG
	100 µg/ml	1 ml	ref. P814280	BG
ISOPROTURON			CAS 34123-59-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P829810	BG
	100 µg/ml	1 ml	ref. P829820	BG

ISOPROTURON-D6 [Dimethylamino-D6]			CAS 217487-17-7	Price
acetone	100 µg/ml	1 ml	ref. P829850	EE
nonane	100 µg/ml	1 ml	ref. P829860	EE
ISOVALERALDEHYDE			CAS 590-86-3	Price
acetonitrile	1 000 µg/ml	1 ml	ref. P886800	BH
	1 000 µg/ml	5 ml	ref. P886810	CF
ISOVALERALDEHYDE-DNPH			CAS 2256-01-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P829930	BG
	100 µg/ml	1 ml	ref. P829940	BG
ISOXABEN			CAS 82558-50-7	Price
acetone	10 µg/ml	10 ml	ref. P829970	BG
	100 µg/ml	1 ml	ref. P829980	BG
iso-octane	10 µg/ml	10 ml	ref. P829990	BG
	100 µg/ml	1 ml	ref. P830000	BG
ISOXADIFEN-ETHYL			CAS 163520-33-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P830030	BH
	100 µg/ml	1 ml	ref. P830040	BH
cyclohexane	10 µg/ml	10 ml	ref. P830050	BH
	100 µg/ml	1 ml	ref. P830060	BH
ISOXAFLUTOLE			CAS 141112-29-0	Price
acetonitrile	100 µg/ml	1 ml	ref. P876820	BG
cyclohexane	10 µg/ml	10 ml	ref. P876830	BG
KRESOXIM-METHYL			CAS 143390-89-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P830130	BG
	100 µg/ml	1 ml	ref. P830140	BG
methanol	10 µg/ml	10 ml	ref. P830150	BG
	100 µg/ml	1 ml	ref. P830160	BG
LENACIL			CAS 2164-08-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P876850	BG
	100 µg/ml	1 ml	ref. P876840	BG
LEPTOPHOS			CAS 21609-90-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P880160	BG
methanol	100 µg/ml	1 ml	ref. P880150	BG
LETHANE 384			CAS 112-56-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P880180	BG
methanol	100 µg/ml	1 ml	ref. P880170	BG
LINALOOL			CAS 78-70-6	Price
ethanol	1 000 µg/ml	1 ml	ref. P876860	BH
LINURON			CAS 330-55-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P830190	BG
	100 µg/ml	1 ml	ref. P830200	BG
MALAOXON			CAS 1634-78-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P830340	BG
	100 µg/ml	1 ml	ref. P830350	BG
MALATHION			CAS 121-75-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P830380	BG
	100 µg/ml	1 ml	ref. P830390	BG
cyclohexane	10 µg/ml	10 ml	ref. P830400	BG
	100 µg/ml	1 ml	ref. P830410	BG
MALATHION-D6 [Dimethyl-D6]			id.: 8090	Price
acetone	100 µg/ml	1 ml	ref. P830440	FF
nonane	100 µg/ml	1 ml	ref. P830450	FF
MCPA			CAS 94-74-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P830560	BG
	100 µg/ml	1 ml	ref. P830570	BG
methanol	10 µg/ml	10 ml	ref. P830580	BG
	100 µg/ml	1 ml	ref. P830590	BG
MCPA-METHYL ESTER			CAS 2436-73-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P880200	BG
methanol	100 µg/ml	1 ml	ref. P880190	BG
MCPB			CAS 94-81-5	Price
acetone	10 µg/ml	10 ml	ref. P830750	BG
	100 µg/ml	1 ml	ref. P830760	BG
acetonitrile	10 µg/ml	10 ml	ref. P830770	BG
	100 µg/ml	1 ml	ref. P830780	BG
MECARBAM			CAS 2595-54-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P830840	BG
	100 µg/ml	1 ml	ref. P830850	BG
cyclohexane	10 µg/ml	10 ml	ref. P830860	BG
	100 µg/ml	1 ml	ref. P830870	BG

MECOPROP (MCPP)			CAS 7085-19-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P830940	BG
	100 µg/ml	1 ml	ref. P830950	BG
cyclohexane	10 µg/ml	10 ml	ref. P830960	BG
	100 µg/ml	1 ml	ref. P830970	BG
methanol	10 µg/ml	10 ml	ref. P830980	BG
	100 µg/ml	1 ml	ref. P830990	BG
MECOPROP-D3			CAS 352431-15-3	Price
acetone	100 µg/ml	1 ml	ref. P831020	FK
n-nonane	100 µg/ml	1 ml	ref. P831030	FK
MEFENACET			CAS 73250-68-7	Price
iso-octane	10 µg/ml	10 ml	ref. P876890	BG
methanol	100 µg/ml	1 ml	ref. P876900	BG
MEFENPYR-DIETHYL			CAS 135590-91-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P876910	BH
MEFLUIDIDE			CAS 53780-34-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P876920	BH
MEPRONIL			CAS 55814-41-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P879130	BG
methanol	100 µg/ml	1 ml	ref. P885630	BG
MERPHOS			CAS 150-50-5	Price
n-hexane	1 000 µg/ml	1 ml	ref. P876930	BAD
methanol	100 µg/ml	1 ml	ref. P876940	CF
METALAXYL			CAS 57837-19-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P831180	BG
	100 µg/ml	1 ml	ref. P831190	BG
cyclohexane	10 µg/ml	10 ml	ref. P831200	BG
	100 µg/ml	1 ml	ref. P831210	BG
methanol	10 µg/ml	10 ml	ref. P831220	BG
	100 µg/ml	1 ml	ref. P831230	BG
METALAXYL-M			CAS 70630-17-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P831260	BH
	100 µg/ml	1 ml	ref. P831270	BH
METALDEHYDE			CAS 9002-91-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P831340	BG
	100 µg/ml	1 ml	ref. P831350	BG
METAMITRON			CAS 41394-05-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P831380	BG
	100 µg/ml	1 ml	ref. P831390	BG
methanol	10 µg/ml	10 ml	ref. P831400	BG
	100 µg/ml	1 ml	ref. P831410	BG
METAZACHLOR			CAS 67129-08-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P831480	BG
	100 µg/ml	1 ml	ref. P831490	BG
cyclohexane	10 µg/ml	10 ml	ref. P831500	BG
	100 µg/ml	1 ml	ref. P831510	BG
METCONAZOLE			CAS 125116-23-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P831540	BG
	100 µg/ml	1 ml	ref. P831550	BG
cyclohexane	10 µg/ml	10 ml	ref. P831560	BG
	100 µg/ml	1 ml	ref. P831570	BG
METHABENZTHIAZURON			CAS 18691-97-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P831600	BG
	100 µg/ml	1 ml	ref. P831610	BG
methanol	10 µg/ml	10 ml	ref. P831620	BG
	100 µg/ml	1 ml	ref. P831630	BG
METHACRIFOS			CAS 62610-77-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P876960	BG
	100 µg/ml	1 ml	ref. P876950	BG
METHACRYLIC ACID -ETHYL-ESTER			CAS 97-63-2	Price
n-hexane	10 µg/ml	10 ml	ref. P824690	BH
	100 µg/ml	1 ml	ref. P824700	BH
methanol purge & trap	10 µg/ml	10 ml	ref. P824710	BH
	100 µg/ml	1 ml	ref. P824720	BH
	1 000 µg/ml	1 ml	ref. P876970	CA
METHACRYLIC ACID-METHYL ESTER			CAS 80-62-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P832710	BJ
	100 µg/ml	1 ml	ref. P832720	BJ

METHACRYLONITRILE			CAS 126-98-7	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P831700	BG
	100 µg/ml	1 ml	ref. P831710	BG
METHAMIDOPHOS			CAS 10265-92-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P831780	BG
	100 µg/ml	1 ml	ref. P831790	BG
methanol	10 µg/ml	10 ml	ref. P831800	BG
	100 µg/ml	1 ml	ref. P831810	BG
METHANOL			CAS 67-56-1	Price
dichloromethane	10 µg/ml	10 ml	ref. P831740	BG
	100 µg/ml	1 ml	ref. P831750	BG
METHIDATHION			CAS 950-37-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P831880	BG
	100 µg/ml	1 ml	ref. P831890	BG
cyclohexane	10 µg/ml	10 ml	ref. P831900	BG
	100 µg/ml	1 ml	ref. P831910	BG
METHIOCARB			CAS 2032-65-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P832020	BG
	100 µg/ml	1 ml	ref. P832030	BG
METHIOCARB SULFONE			CAS 2179-25-1	Price
methyl-tert-butylether	100 µg/ml	1 ml	ref. P876990	BK
METHIOCARB SULFOXIDE			CAS 2635-10-1	Price
methyl-tert-butylether	100 µg/ml	1 ml	ref. P877000	BK
METHOMYL			CAS 16752-77-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P832060	BG
	100 µg/ml	1 ml	ref. P832070	BG
METHOPRENE			CAS 40596-69-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P880220	BH
methanol	100 µg/ml	1 ml	ref. P880210	BH
METHOXYCHLOR (DMTD)			CAS 72-43-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P832140	BG
	100 µg/ml	1 ml	ref. P832150	BG
cyclohexane	10 µg/ml	10 ml	ref. P832160	BG
	100 µg/ml	1 ml	ref. P832170	BG
iso-octane	10 µg/ml	10 ml	ref. P832180	BG
	100 µg/ml	1 ml	ref. P832190	BG
methanol	10 µg/ml	10 ml	ref. P832200	BG
	100 µg/ml	1 ml	ref. P832210	BG
2-METHOXY-5-METHYLANILINE			CAS 120-71-8	Price
ethyl acetate	10 µg/ml	10 ml	ref. P873120	BG
methanol	100 µg/ml	5 ml	ref. P873110	BG
METHYL-tert-AMYLETHER			CAS 994-05-8	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P804120	CB
	100 µg/ml	1 ml	ref. P804130	CB
METHYL 2-BROMOBUTYRATE			CAS 3196-15-4	Price
methyl-tert-butylether	10 µg/ml	10 ml	ref. P832350	BG
	100 µg/ml	1 ml	ref. P832470	BG
METHYL 2-BROMOPROPIONATE			CAS 5445-17-0	Price
methyl-tert-butylether	10 µg/ml	10 ml	ref. P832310	BG
	100 µg/ml	1 ml	ref. P832320	BG
METHYL- tert-BUTYL ETHER			CAS 1634-04-4	Price
methanol	2 000 µg/ml	1 ml	ref. F116311	DC
1-METHYLCHRYSENE			CAS 3351-28-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P872470	DA
5-METHYLCHRYSENE			CAS 3697-24-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P873700	DF
cyclohexane	10 µg/ml	10 ml	ref. P873690	DF
6-METHYLCHRYSENE			CAS 1705-85-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P873720	CF
dichloromethane	100 µg/ml	5 ml	ref. P873710	FK
2-METHYL-4,6-DINITROPHENOL			CAS 534-52-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P832500	BG
	100 µg/ml	1 ml	ref. P832510	BG
methanol	10 µg/ml	10 ml	ref. P832520	BG
	100 µg/ml	1 ml	ref. P832530	BG
4,4'-METHYLENE-BIS(2-CHLOROANIL)			CAS 101-14-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P873430	BG
	100 µg/ml	1 ml	ref. P873420	BG

2--METHYL-FLUORANTHENE			CAS 33543-31-6	Price
dichloromethane	100 µg/ml	5 ml	ref. P873710	BJ
toluene	50 µg/ml	1 ml	ref. P873130	BJ
METHYL ISOBUTYL KETONE-DNPH			CAS 1655-42-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P832670	BJ
	100 µg/ml	1 ml	ref. P832680	BJ
METHYL MERCAPTAN			CAS 74-93-1	Price
methanol	1 000 µg/ml	10 ml	ref. P877010	BH
METHYL METHACRYLATE			CAS 80-62-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P832710	BJ
	100 µg/ml	1 ml	ref. P832720	BJ
1-METHYLNAPHTHALENE			CAS 90-12-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P832750	BG
	100 µg/ml	1 ml	ref. P832760	BG
cyclohexane	10 µg/ml	10 ml	ref. P832770	BG
	100 µg/ml	1 ml	ref. P832780	BG
methanol	10 µg/ml	10 ml	ref. P832790	BG
	100 µg/ml	1 ml	ref. P832800	BG
2-METHYLNAPHTHALENE			CAS 91-57-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P832830	BG
	100 µg/ml	1 ml	ref. P832840	BG
cyclohexane	10 µg/ml	10 ml	ref. P832850	BG
	100 µg/ml	1 ml	ref. P832860	BG
dichloromethane	10 µg/ml	10 ml	ref. P832870	BG
	100 µg/ml	1 ml	ref. P832880	BG
2-METHYL-NONANE			CAS 871-83-0	Price
dichloromethane	1 000 µg/ml	1 ml	ref. P873140	BH
2-METHYLPENTANE			CAS 107-83-5	Price
methanol	100 µg/ml	5 ml	ref. P873150	BG
4-METHYL-2-PENTANONE			CAS 108-10-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P832970	BG
	100 µg/ml	1 ml	ref. P832980	BG
methanol purge & trap	10 µg/ml	10 ml	ref. P832990	BG
	100 µg/ml	1 ml	ref. P833000	BG
1-METHYLPHENANTHRENE			CAS 832-69-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P833030	BH
	100 µg/ml	1 ml	ref. P833040	BH
cyclohexane	10 µg/ml	10 ml	ref. P833050	BH
	100 µg/ml	1 ml	ref. P833060	BH
2-METHYLPHENANTHRENE			CAS 2531-84-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P833090	CJ
	100 µg/ml	1 ml	ref. P833100	CJ
cyclohexane	10 µg/ml	10 ml	ref. P833110	CJ
	100 µg/ml	1 ml	ref. P833120	CJ
2-METHYLPHENOL			CAS 95-48-7	Price
methanol	100 µg/ml	1 ml	ref. P873160	BG
3-METHYLPHENOL			CAS 108-39-4	Price
methanol	10 µg/ml	10 ml	ref. P833150	BG
	100 µg/ml	1 ml	ref. P833160	BG
4-METHYLPHENOL			CAS 106-44-5	Price
cyclohexane	100 µg/ml	1 ml	ref. P873640	CC
2-METHYL-1-PROPANOL			CAS 78-83-1	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P833190	BG
	100 µg/ml	1 ml	ref. P833200	BG
1-METHYLPYRENE			CAS 2381-21-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P833230	CE
	100 µg/ml	1 ml	ref. P833240	CE
cyclohexane	10 µg/ml	10 ml	ref. P833250	CE
	100 µg/ml	1 ml	ref. P833260	CE
METOBROMURON			CAS 3060-89-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P833290	BG
	100 µg/ml	1 ml	ref. P833300	BG
cyclohexane	10 µg/ml	10 ml	ref. P833310	BG
	100 µg/ml	1 ml	ref. P833320	BG
METOLACHLOR			CAS 51218-45-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P833350	BG
	100 µg/ml	1 ml	ref. P833360	BG
cyclohexane	10 µg/ml	10 ml	ref. P833370	BG
	100 µg/ml	1 ml	ref. P833380	BG

METOSULAM			CAS 139528-85-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P877030	BH
METOXURON			CAS 19937-59-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P833480	BG
	100 µg/ml	1 ml	ref. P833490	BG
METRIBUZIN			CAS 21087-64-9	Price
acetone	10 µg/ml	10 ml	ref. P833520	BG
	100 µg/ml	1 ml	ref. P833530	BG
acetonitrile	10 µg/ml	10 ml	ref. P833540	BG
	100 µg/ml	1 ml	ref. P833550	BG
cyclohexane	10 µg/ml	10 ml	ref. P833560	BG
	100 µg/ml	1 ml	ref. P833570	BG
METSULFURON-METHYL			CAS 74223-64-6	Price
methyl-tert-butylether	10 µg/ml	10 ml	ref. P877050	BG
	100 µg/ml	1 ml	ref. P877040	BG
MEVINPHOS			CAS 7786-34-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P833670	BG
	100 µg/ml	1 ml	ref. P833680	BG
cyclohexane	10 µg/ml	10 ml	ref. P833690	BG
	100 µg/ml	1 ml	ref. P833700	BG
MGK 264			CAS 113-48-4	Price
methanol	100 µg/ml	1 ml	ref. P877060	BG
MIREX			CAS 2385-85-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P833730	BG
	100 µg/ml	1 ml	ref. P833740	BG
iso-octane	10 µg/ml	10 ml	ref. P833750	BG
	100 µg/ml	1 ml	ref. P833760	BG
MOLINATE			CAS 2212-67-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P833790	BG
	100 µg/ml	1 ml	ref. P833800	BG
MONOCROTAPHOS			CAS 6923-22-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P833860	BG
	100 µg/ml	1 ml	ref. P833870	BG
iso-octane	10 µg/ml	10 ml	ref. P833880	BG
	100 µg/ml	1 ml	ref. P833890	BG
MONOLINURON			CAS 1746-81-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P833920	BG
	100 µg/ml	1 ml	ref. P833930	BG
cyclohexane	10 µg/ml	10 ml	ref. P833940	BG
	100 µg/ml	1 ml	ref. P833950	BG
MONURON			CAS 150-68-5	Price
methanol	10 µg/ml	10 ml	ref. P833980	BG
	100 µg/ml	1 ml	ref. P833990	BG
MYCLOBUTANIL			CAS 88671-89-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P834050	BG
	100 µg/ml	1 ml	ref. P834060	BG
NALED			CAS 300-76-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P834090	BG
	100 µg/ml	1 ml	ref. P834100	BG
iso-octane	10 µg/ml	10 ml	ref. F115801	BG
NAPHTHALENE			CAS 91-20-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P834130	BG
	100 µg/ml	1 ml	ref. P834140	BG
cyclohexane	10 µg/ml	10 ml	ref. P834150	BG
	100 µg/ml	1 ml	ref. P834160	BG
methanol	10 µg/ml	10 ml	ref. P834170	BG
	100 µg/ml	1 ml	ref. P834180	BG
NAPHTHALENE -D8			CAS 1146-65-2	Price
cyclohexane	100 µg/ml	1 ml	ref. P834210	CC
	2 000 µg/ml	1 ml	ref. F062455	CK
dichloromethane	100 µg/ml	1 ml	ref. P834220	CC
	2 000 µg/ml	1 ml	ref. F062465	CK
methyl-tert-butylether	100 µg/ml	1 ml	ref. P834230	CC
	2 000 µg/ml	1 ml	ref. F062475	CK
NAPHTHALENE-1-ADELHYDE			CAS 66-77-3	Price
ethanol	1 000 µg/ml	1 ml	ref. P877080	CF
methanol	100 µg/ml	1 ml	ref. P877090	BH
NAPROPAMIDE			CAS 15299-99-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P834260	BG
	100 µg/ml	1 ml	ref. P834270	BG

NEBURON			CAS 555-37-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P834300	BG
	100 µg/ml	1 ml	ref. P834310	BG
cyclohexane	10 µg/ml	10 ml	ref. P834320	BG
	100 µg/ml	1 ml	ref. P834330	BG
NICLOSAMID			CAS 50-65-7	Price
methanol	100 µg/ml	1 ml	ref. P877240	BG
NICOSULFURON			CAS 111991-09-4	Price
acetone	10 µg/ml	10 ml	ref. P834360	BH
	100 µg/ml	1 ml	ref. P834370	BH
9-NITROANTHRACENE			CAS 602-60-8	Price
toluene	100 µg/ml	1 ml	ref. P873790	CA
NITROBENZENE			CAS 98-95-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P834400	BG
	100 µg/ml	1 ml	ref. P834410	BG
methanol	10 µg/ml	10 ml	ref. P834420	BG
	100 µg/ml	1 ml	ref. P834430	BG
NITROBENZENE-D5			CAS 4165-60-0	Price
acetone	1000 µg/ml	1 ml	ref. P834460	CA
dichloromethane	1000 µg/ml	1 ml	ref. P834470	CA
methyl- <i>tert</i> -butylether	1000 µg/ml	1 ml	ref. P834480	CA
6-NITROBENZENE(a)PYRENE			CAS 63041-90-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P873740	DB
cyclohexane	10 µg/ml	10 ml	ref. P873750	DB
toluene	100 µg/ml	1 ml	ref. P873730	DB
6-NITROCHRYSENE			CAS 7496-02-8	Price
toluene	100 µg/ml	1 ml	ref. P873760	CH
NITROFEN			CAS 1836-75-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P834510	BG
	100 µg/ml	1 ml	ref. P834520	BG
cyclohexane	10 µg/ml	10 ml	ref. P834530	BG
	100 µg/ml	1 ml	ref. P834540	BG
3-NITROFLUORANTHENE			CAS 892-21-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P834570	CA
	100 µg/ml	1 ml	ref. P834580	CA
2-NITROFLUORENE			CAS 607-57-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P834610	BG
	100 µg/ml	1 ml	ref. P834620	BG
NITROGUANIDINE			CAS 556-88-7	Price
methanol	10 µg/ml	10 ml	ref. P834650	BG
	100 µg/ml	1 ml	ref. P834660	BG
2-NITROMESITYLENE			CAS 603-71-4	Price
methanol	2 000 µg/ml	1 ml	ref. P881770	BH
1-NITRONAPHTHALENE			CAS 86-57-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P834690	BG
	100 µg/ml	1 ml	ref. P834700	BG
methanol	10 µg/ml	10 ml	ref. P834710	BG
	100 µg/ml	1 ml	ref. P834720	BG
2-NITRONAPHTHALENE			CAS 581-89-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P873170	CB
2-NITROPHENOL			CAS 88-75-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P873190	BG
methanol	100 µg/ml	1 ml	ref. P873180	BG
4-NITROPHENOL			CAS 100-02-7	Price
methanol	10 µg/ml	10 ml	ref. P834750	BG
	100 µg/ml	1 ml	ref. P834760	BG
4-NITROPHENOL-D4 [2,3,5,6-D4]			CAS 93951-79-2	Price
acetone	100 µg/ml	1 ml	ref. P834790	CA
2-NITROPROPANE			CAS 79-46-9	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P834820	BG
	100 µg/ml	1 ml	ref. P834830	BG
1-NITROPYRENE			CAS 5522-43-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P834860	BK
	100 µg/ml	1 ml	ref. P834870	BK
n-NITROSODIMETHYLAMINE			CAS 62-75-9	Price
methanol	10 µg/ml	10 ml	ref. P834900	BG
	100 µg/ml	1 ml	ref. P834910	BG

n-NITROSODIPHENYLAMINE			CAS 86-30-6	Price
methanol	10 µg/ml	10 ml	ref. P834940	BG
n-NITROSO-DI-n-PROPYLAMINE			CAS 621-64-7	Price
methanol	10 µg/ml	10 ml	ref. P834980	BG
	100 µg/ml	1 ml	ref. P834990	BG
n-NITROTHAL-ISOPROPYL			CAS 10552-74-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P880240	BG
methanol	100 µg/ml	1 ml	ref. P880230	BG
2-NITROTOLUENE			CAS 88-72-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P835020	BG
	100 µg/ml	1 ml	ref. P835030	BG
3-NITROTOLUENE			CAS 99-08-1	Price
acetone	10 µg/ml	10 ml	ref. P835060	BG
	100 µg/ml	1 ml	ref. P835070	BG
acetonitrile	10 µg/ml	10 ml	ref. P835080	BG
	100 µg/ml	1 ml	ref. P835090	BG
4-NITROTOLUENE			CAS 99-99-0	Price
acetone	10 µg/ml	10 ml	ref. P835120	BG
	100 µg/ml	1 ml	ref. P835130	BG
acetonitrile	10 µg/ml	10 ml	ref. P835140	BG
	100 µg/ml	1 ml	ref. P835150	BG
cis-NONACHLOR			CAS 5103-73-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P880280	CE
iso-octane	100 µg/ml	1 ml	ref. P880270	CE
methanol	100 µg/ml	1 ml	ref. P880260	CE
	1 000 µg/ml	1 ml	ref. P880250	KE
trans-NONACHLOR			CAS 39765-80-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P835180	CF
	100 µg/ml	1 ml	ref. P835190	CF
cyclohexane	10 µg/ml	10 ml	ref. P835200	CF
	100 µg/ml	1 ml	ref. P835210	CF
methanol	10 µg/ml	10 ml	ref. P835220	CF
	100 µg/ml	1 ml	ref. P835230	CF
n-NONADECANE			CAS 629-92-5	Price
dichloromethane	10 000 µg/ml	5 ml	ref. P877250	FE
NONANAL			CAS 124-19-6	Price
methanol	1 000 µg/ml	1 ml	ref. P877280	BH
n-NONANE			CAS 111-84-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P835260	BG
	100 µg/ml	1 ml	ref. P835270	BG
methanol	10 µg/ml	10 ml	ref. P835280	BG
	100 µg/ml	1 ml	ref. P835290	BG
n-NONATRIACONTANE			CAS 7194-86-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P835320	BG
	100 µg/ml	1 ml	ref. P835330	BG
NONYLPHENOL			CAS 25154-52-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P835360	BG
	100 µg/ml	1 ml	ref. P835370	BG
cyclohexane	10 µg/ml	10 ml	ref. P835380	BG
	100 µg/ml	1 ml	ref. P835390	BG
4-NONYLPHENOL (TECHNICAL)			CAS 84852-15-3	Price
toluene	100 µg/ml	1 ml	ref. F129201	BG
4-n-NONYLPHENOL			CAS 104-40-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P835420	BG
	100 µg/ml	1 ml	ref. P835430	BG
cyclohexane	10 µg/ml	10 ml	ref. P835440	BG
	100 µg/ml	1 ml	ref. P835450	BG
NORFLURAZON			CAS 27314-13-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P835510	BH
	100 µg/ml	1 ml	ref. P835520	BH
iso-octane	10 µg/ml	10 ml	ref. P835530	BH
	100 µg/ml	1 ml	ref. P835540	BH
NUARIMOL			CAS 63284-71-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P880330	BG
methanol	100 µg/ml	1 ml	ref. P880320	BG
OCTACHLORONAPHTHALENE			CAS 2234-13-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P877340	CA
cyclohexane	10 µg/ml	10 ml	ref. P877360	CA
iso-octane	10 µg/ml	10 ml	ref. P877350	CA
methanol	100 µg/ml	1 ml	ref. P877330	CA

n-OCTACOSANE			CAS 630-02-4	Price
dichloromethane	10 µg/ml	10 ml	ref. P835570	BG
	100 µg/ml	1 ml	ref. P835580	BG
OCTADECANOIC ACID			CAS 57-11-4	Price
n-hexane	1 000 µg/ml	1 ml	ref. P877370	BH
OCTADECANOIC ACID-METHYL ESTER			CAS 112-61-8	Price
n-hexane	10 000 µg/ml	1 ml	ref. P877380	BJ
OCTANAL			CAS 124-13-0	Price
acetonitrile	1 000 µg/ml	1 ml	ref. P877400	BH
	1 000 µg/ml	5 ml	ref. P877390	CF
n-OCTANE			CAS 111-65-9	Price
methanol	1 000 µg/ml	5 ml	ref. P877260	BG
OCTANOIC ACID METHYL ESTER			CAS 111-11-5	Price
n-hexane	10 µg/ml	10 ml	ref. P877410	BG
1-OCTANOL			CAS 111-87-5	Price
methanol	2 500 µg/ml	5 ml	ref. P872480	BJ
n-OCTATRIACONTANE			CAS 7194-85-6	Price
ethylbenzene	10 000 µg/ml	2 ml	ref. P877270	DJA
4-n-OCTYLPHENOL			CAS 1806-26-4	Price
iso-octane	10 µg/ml	10 ml	ref. P835640	BG
	100 µg/ml	1 ml	ref. P835650	BG
4-tert-OCTYLPHENOL			CAS 140-66-9	Price
iso-octane	10 µg/ml	10 ml	ref. P835680	BG
	100 µg/ml	1 ml	ref. P835690	BG
toluene	100 µg/ml	1 ml	ref. F129211	BG
OFURACE			CAS 58810-48-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P877440	BH
ethyl acetate	100 µg/ml	1 ml	ref. P877430	BH
OMETHOATE			CAS 1113-02-6	Price
acetone	10 µg/ml	10 ml	ref. P835720	BG
	100 µg/ml	1 ml	ref. P835730	BG
ORYZALIN			CAS 19044-88-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P835790	BG
	100 µg/ml	1 ml	ref. P835800	BG
methanol	1 000 µg/ml	1 ml	ref. F115991	BK
OXADIARGYL			CAS 39807-15-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P877500	BG
OXADIAZON			CAS 19666-30-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P835830	BG
	100 µg/ml	1 ml	ref. P835840	BG
iso-octane	10 µg/ml	10 ml	ref. P835850	BG
	100 µg/ml	1 ml	ref. P835860	BG
OXADIXYL			CAS 77732-09-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P835890	BG
	100 µg/ml	1 ml	ref. P835900	BG
cyclohexane	10 µg/ml	10 ml	ref. P835910	BG
	100 µg/ml	1 ml	ref. P835920	BG
iso-octane	10 µg/ml	10 ml	ref. P835930	BG
	100 µg/ml	1 ml	ref. P835940	BG
methanol	10 µg/ml	10 ml	ref. P835950	BG
	100 µg/ml	1 ml	ref. P835960	BG
OXAMYL			CAS 23135-22-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P836010	BG
	100 µg/ml	1 ml	ref. P836020	BG
OXYFLUORFEN			CAS 42874-03-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P877570	BG
methanol	100 µg/ml	1 ml	ref. P877580	BG
PACLOBUTRAZOL			CAS 76738-62-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P836140	BG
	100 µg/ml	1 ml	ref. P836150	BG
iso-octane	10 µg/ml	10 ml	ref. P836160	BG
	100 µg/ml	1 ml	ref. P836170	BG
methanol	10 µg/ml	10 ml	ref. P836180	BG
	100 µg/ml	1 ml	ref. P836190	BG
PARAFORMALDEHYDE			CAS 30525-89-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P836220	BG
	100 µg/ml	1 ml	ref. P836230	BG

PARAFUCHSIN HYDROCHLORIDE			CAS 569-61-9	Price
methanol	100 µg/ml	1 ml	ref. P877590	BG
PARAOXON-ETHYL			CAS 311-45-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P877600	BG
	100 µg/ml	1 ml	ref. P877610	BG
methanol	100 µg/ml	1 ml	ref. P877620	BG
	1 000 µg/ml	1 ml	ref. P877630	BK
PARAOXON-METHYL			CAS 950-35-6	Price
iso-octane	10 µg/ml	10 ml	ref. P877650	BG
methanol	100 µg/ml	1 ml	ref. P877640	BG
PARAQUAT DICHLORIDE			CAS 1910-42-5	Price
water	10 µg/ml	10 ml	ref. P836380	BG
	100 µg/ml	1 ml	ref. P836390	BG
PARATHION			CAS 56-38-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P836560	BG
	100 µg/ml	1 ml	ref. P836570	BG
cyclohexane	10 µg/ml	10 ml	ref. P836580	BG
	100 µg/ml	1 ml	ref. P836590	BG
PARATHION-METHYL			CAS 298-00-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P836660	BG
	100 µg/ml	1 ml	ref. P836670	BG
cyclohexane	10 µg/ml	10 ml	ref. P836680	BG
	100 µg/ml	1 ml	ref. P836690	BG
PARSOL 1789			CAS 70356-09-1	Price
methanol	100 µg/ml	1 ml	ref. P877660	BG
PBB 49 (2,2',4,5'-Tetrabromobiphenyl)			CAS 60044-24-8	Price
iso-octane	10 µg/ml	10 ml	ref. P836760	CD
	100 µg/ml	1 ml	ref. P836770	CD
PBB 103 (2,2',4,5',6-Pentabromobiphenyl)			CAS 59080-39-6	Price
iso-octane	10 µg/ml	10 ml	ref. P836840	DA
	100 µg/ml	1 ml	ref. P836850	DA
PBB 209 (Decabromobiphenyl)			CAS 13654-09-6	Price
iso-octane	10 µg/ml	10 ml	ref. P836920	BH
	100 µg/ml	1 ml	ref. P836930	BH
PBDE #3			CAS 101-55-3	Price
methanol	100 µg/ml	1 ml	ref. P877670	BG
PCB 1 (2-Chlorobiphenyl)			CAS 2051-60-7	Price
n-hexane	10 µg/ml	10 ml	ref. P837380	BG
	100 µg/ml	1 ml	ref. P837390	BG
iso-octane	10 µg/ml	10 ml	ref. P837400	BG
	100 µg/ml	1 ml	ref. P837410	BG
PCB 2 (3-Chlorobiphenyl)			CAS 2051-61-8	Price
n-hexane	10 µg/ml	10 ml	ref. P837440	BG
	100 µg/ml	1 ml	ref. P837450	BG
iso-octane	10 µg/ml	10 ml	ref. P837460	BG
	100 µg/ml	1 ml	ref. P837470	BG
PCB 3 (4-Chlorobiphenyl)			CAS 2051-62-9	Price
n-hexane	10 µg/ml	10 ml	ref. P837500	BG
	100 µg/ml	1 ml	ref. P837510	BG
iso-octane	10 µg/ml	10 ml	ref. P837520	BG
	100 µg/ml	1 ml	ref. P837530	BG
PCB 4 (2,2'-Dichlorobiphenyl)			CAS 13029-08-8	Price
n-hexane	10 µg/ml	10 ml	ref. P837560	BG
	100 µg/ml	1 ml	ref. P837570	BG
iso-octane	10 µg/ml	10 ml	ref. P837580	BG
	100 µg/ml	1 ml	ref. P837590	BG
PCB 5 (2,3-Dichlorobiphenyl)			CAS 16605-91-7	Price
n-hexane	10 µg/ml	10 ml	ref. P837620	BG
	100 µg/ml	1 ml	ref. P837630	BG
iso-octane	10 µg/ml	10 ml	ref. P837640	BG
	100 µg/ml	1 ml	ref. P837650	BG
PCB 6 (2,3'-Dichlorobiphenyl)			CAS 25569-80-6	Price
n-hexane	10 µg/ml	10 ml	ref. P837680	CE
	100 µg/ml	1 ml	ref. P837690	CE
PCB 7 (2,4-Dichlorobiphenyl)			CAS 33284-50-3	Price
n-hexane	10 µg/ml	10 ml	ref. P837720	BG
	100 µg/ml	1 ml	ref. P837730	BG
iso-octane	10 µg/ml	10 ml	ref. P837740	BG
	100 µg/ml	1 ml	ref. P837750	BG

PCB 8 (2,4'-Dichlorobiphenyl)				CAS 34883-43-7	Price
n-hexane	10 µg/ml	10 ml	ref. P837780	CE	
	100 µg/ml	1 ml	ref. P837790	CE	
iso-octane	10 µg/ml	10 ml	ref. P837800	CE	
	100 µg/ml	1 ml	ref. P837810	CE	
PCB 9 (2,5-Dichlorobiphenyl)				CAS 34883-39-1	Price
n-hexane	10 µg/ml	10 ml	ref. P837840	BG	
	100 µg/ml	1 ml	ref. P837850	BG	
iso-octane	10 µg/ml	10 ml	ref. P837860	BG	
	100 µg/ml	1 ml	ref. P837870	BG	
PCB 10 (2,6-Dichlorobiphenyl)				CAS 33146-45-1	Price
n-hexane	10 µg/ml	10 ml	ref. P837900	BG	
	100 µg/ml	1 ml	ref. P837910	BG	
PCB 11 (3,3'-dichlorobiphenyl)				CAS 2050-67-1	Price
n-hexane	10 µg/ml	10 ml	ref. P837940	BG	
	100 µg/ml	1 ml	ref. P837950	BG	
iso-octane	10 µg/ml	10 ml	ref. P837960	BG	
	100 µg/ml	1 ml	ref. P837970	BG	
PCB 12 (3,4-Dichlorobiphenyl)				CAS 2974-92-7	Price
n-hexane	10 µg/ml	10 ml	ref. P838000	CE	
	100 µg/ml	1 ml	ref. P838010	CE	
iso-octane	10 µg/ml	10 ml	ref. P838020	CE	
	100 µg/ml	1 ml	ref. P838030	CE	
PCB 13 (3,4'-Dichlorobiphenyl)				CAS 2974-90-5	Price
n-hexane	10 µg/ml	10 ml	ref. P838060	CE	
	100 µg/ml	1 ml	ref. P838070	CE	
PCB 14 (3,5-Dichlorobiphenyl)				CAS 34883-41-5	Price
n-hexane	10 µg/ml	10 ml	ref. P838100	BG	
	100 µg/ml	1 ml	ref. P838110	BG	
iso-octane	10 µg/ml	10 ml	ref. P838120	BG	
	100 µg/ml	1 ml	ref. P838130	BG	
PCB 15 (4,4'-Dichlorobiphenyl)				CAS 2050-68-2	Price
n-hexane	10 µg/ml	10 ml	ref. P838160	CE	
	100 µg/ml	1 ml	ref. P838170	CE	
iso-octane	10 µg/ml	10 ml	ref. P838180	CE	
	100 µg/ml	1 ml	ref. P838190	CE	
PCB 16 (2,2',3-Trichlorobiphenyl)				CAS 38444-78-9	Price
n-hexane	10 µg/ml	10 ml	ref. P838220	DC	
	100 µg/ml	1 ml	ref. P838230	DC	
PCB 17 (2,2',4-Trichlorobiphenyl)				CAS 37680-66-3	Price
n-hexane	10 µg/ml	10 ml	ref. P838260	DC	
	100 µg/ml	1 ml	ref. P838270	DC	
PCB 18 (2,2',5-Trichlorobiphenyl)				CAS 37680-65-2	Price
n-hexane	10 µg/ml	10 ml	ref. P838300	BG	
	100 µg/ml	1 ml	ref. P838310	BG	
iso-octane	10 µg/ml	10 ml	ref. P838320	BG	
	100 µg/ml	1 ml	ref. P838330	BG	
PCB 19 (2,2',6-Trichlorobiphenyl)				CAS 38444-73-4	Price
n-hexane	10 µg/ml	10 ml	ref. P838360	DC	
	100 µg/ml	1 ml	ref. P838370	DC	
PCB 20 (2,3,3'-Trichlorobiphenyl)				CAS 38444-84-7	Price
n-hexane	10 µg/ml	10 ml	ref. P838400	CE	
	100 µg/ml	1 ml	ref. P838410	CE	
PCB 21 (2,3,4-Trichlorobiphenyl)				CAS 55702-46-0	Price
n-hexane	10 µg/ml	10 ml	ref. P838440	CE	
	100 µg/ml	1 ml	ref. P838450	CE	
PCB 22 (2,3,4'-Trichlorobiphenyl)				CAS 38444-85-8	Price
n-hexane	10 µg/ml	10 ml	ref. P838480	DC	
	100 µg/ml	1 ml	ref. P838490	DC	
PCB 23 (2,3,5-Trichlorobiphenyl)				CAS 55720-44-0	Price
n-hexane	10 µg/ml	10 ml	ref. P838520	DC	
	100 µg/ml	1 ml	ref. P838530	DC	
PCB 24 (2,3,6-Trichlorobiphenyl)				CAS 55702-45-9	Price
n-hexane	10 µg/ml	10 ml	ref. P838560	CE	
	100 µg/ml	1 ml	ref. P838570	CE	
PCB 25 (2,3',4-Trichlorobiphenyl)				CAS 55712-37-3	Price
n-hexane	10 µg/ml	10 ml	ref. P838600	CE	
	100 µg/ml	1 ml	ref. P838610	CE	

PCB 26 (2,3',5-Trichlorobiphenyl)				CAS 38444-81-4	Price
n-hexane	10 µg/ml	10 ml	ref. P838640	CE	
	100 µg/ml	1 ml	ref. P838650	CE	
PCB 27 (2,3',6-Trichlorobiphenyl)				CAS 38444-76-7	Price
n-hexane	10 µg/ml	10 ml	ref. P838680	DC	
	100 µg/ml	1 ml	ref. P838690	DC	
PCB 28 (2,4,4'-Trichlorobiphenyl)				CAS 7012-37-5	Price
n-hexane	10 µg/ml	10 ml	ref. P838720	CE	
	100 µg/ml	1 ml	ref. P838730	CE	
iso-octane	10 µg/ml	10 ml	ref. P838740	CE	
	100 µg/ml	1 ml	ref. P838750	CE	
PCB 29 (2,4,5-Trichlorobiphenyl)				CAS 15862-07-4	Price
n-hexane	10 µg/ml	10 ml	ref. P838780	CE	
	100 µg/ml	1 ml	ref. P838790	CE	
	400 µg/ml	1 ml	ref. F131401	CE	
iso-octane	10 µg/ml	10 ml	ref. P838800	CE	
	100 µg/ml	1 ml	ref. P838810	CE	
PCB 30 (2,4,6 -Trichlorobiphenyl)				CAS 35693-92-6	Price
n-hexane	10 µg/ml	10 ml	ref. P838840	BG	
	100 µg/ml	1 ml	ref. P838850	BG	
	400 µg/ml	1 ml	ref. F131411	BG	
iso-octane	10 µg/ml	1 ml	ref. F131480	BG	
PCB 31 (2,4',5-Trichlorobiphenyl)				CAS 16606-02-3	Price
n-hexane	10 µg/ml	10 ml	ref. P838880	CE	
	100 µg/ml	1 ml	ref. P838890	CE	
iso-octane	10 µg/ml	10 ml	ref. P804660	CE	
	100 µg/ml	1 ml	ref. P838910	CE	
PCB 32 (2,4',6-Trichlorobiphenyl)				CAS 38444-77-4	Price
n-hexane	10 µg/ml	10 ml	ref. P838940	DC	
	100 µg/ml	1 ml	ref. P838950	DC	
PCB 33 (2',3,4-Trichlorobiphenyl)				CAS 38444-86-9	Price
n-hexane	10 µg/ml	10 ml	ref. P838980	CE	
	100 µg/ml	1 ml	ref. P838990	CE	
PCB 34 (2',3,5-Trichlorobiphenyl)				CAS 37680-68-5	Price
n-hexane	10 µg/ml	10 ml	ref. P839020	CE	
	100 µg/ml	1 ml	ref. P839030	CE	
PCB 35 (3,3',4-Trichlorobiphenyl)				CAS 37680-69-6	Price
n-hexane	10 µg/ml	10 ml	ref. P839060	CE	
	100 µg/ml	1 ml	ref. P839070	CE	
iso-octane	10 µg/ml	10 ml	ref. P839080	CE	
	100 µg/ml	1 ml	ref. P839090	CE	
PCB 36 (3,3',5-Trichlorobiphenyl)				CAS 38444-87-0	Price
n-hexane	10 µg/ml	10 ml	ref. P839120	DC	
	100 µg/ml	1 ml	ref. P839130	DC	
PCB 37 (3,4,4'-Trichlorobiphenyl)				CAS 38444-90-5	Price
n-hexane	10 µg/ml	10 ml	ref. P839160	DC	
	100 µg/ml	1 ml	ref. P839170	DC	
PCB 38 (3,4,5-Trichlorobiphenyl)				CAS 53555-66-1	Price
n-hexane	10 µg/ml	10 ml	ref. P839200	DC	
	100 µg/ml	1 ml	ref. P839210	DC	
PCB 39 (3,4',5-Trichlorobiphenyl)				CAS 38444-88-1	Price
n-hexane	10 µg/ml	10 ml	ref. P839240	DC	
	100 µg/ml	1 ml	ref. P839250	DC	
PCB 40 (2,2',3,3'-Tetrachlorobiphenyl)				CAS 38444-93-8	Price
n-hexane	10 µg/ml	10 ml	ref. P839280	BG	
	100 µg/ml	1 ml	ref. P839290	BG	
PCB 41 (2,2',3,4-Tetrachlorobiphenyl)				CAS 53663-59-9	Price
n-hexane	10 µg/ml	10 ml	ref. P839320	DC	
	100 µg/ml	1 ml	ref. P839330	DC	
PCB 42 (2,2',3,4')-Tetrachlorobiphenyl)				CAS 36559-22-5	Price
n-hexane	10 µg/ml	10 ml	ref. P839360	DC	
	100 µg/ml	1 ml	ref. P839370	DC	
PCB 43 (2,2',3,5-Tetrachlorobiphenyl)				CAS 70362-46-8	Price
n-hexane	10 µg/ml	10 ml	ref. P839400	DC	
	100 µg/ml	1 ml	ref. P839410	DC	
PCB 44 (2,2',3,5')-Tetrachlorobiphenyl)				CAS 41464-39-5	Price
n-hexane	10 µg/ml	10 ml	ref. P839440	CE	
	100 µg/ml	1 ml	ref. P839450	CE	
iso-octane	10 µg/ml	10 ml	ref. P839460	CE	
	100 µg/ml	1 ml	ref. P839470	CE	

PCB 45 (2,2',3,6-Tetrachlorobiphenyl)			CAS 70362-45-7	Price
n-hexane	10 µg/ml	10 ml	ref. P839500	DC
	100 µg/ml	1 ml	ref. P839510	DC
PCB 46 (2,2',3,6'-Tetrachlorobiphenyl)			CAS 41464-47-5	Price
n-hexane	10 µg/ml	10 ml	ref. P839540	DC
	100 µg/ml	1 ml	ref. P839550	DC
PCB 47 (2,2',4,4'-Tetrachlorobiphenyl)			CAS 2437-79-8	Price
n-hexane	10 µg/ml	10 ml	ref. P839580	BG
	100 µg/ml	1 ml	ref. P839590	BG
PCB 48 (2,2',4,5-Tetrachlorobiphenyl)			CAS 70362-47-9	Price
n-hexane	10 µg/ml	10 ml	ref. P839620	DC
	100 µg/ml	1 ml	ref. P839630	DC
PCB 49 (2,2',4,5'-Tetrachlorobiphenyl)			CAS 41464-40-8	Price
n-hexane	10 µg/ml	10 ml	ref. P839660	CE
	100 µg/ml	1 ml	ref. P839670	CE
PCB 50 (2,2',4,6-Tetrachlorobiphenyl)			CAS 62796-65-0	Price
n-hexane	10 µg/ml	10 ml	ref. P839700	CE
	100 µg/ml	1 ml	ref. P839710	CE
PCB 51 (2,2',4,6'-Tetrachlorobiphenyl)			CAS 68194-04-7	Price
n-hexane	10 µg/ml	10 ml	ref. P839740	EA
	100 µg/ml	1 ml	ref. P839750	EA
PCB 52 (2,2',5,5'-Tetrachlorobiphenyl)			CAS 35693-99-3	Price
n-hexane	10 µg/ml	10 ml	ref. P839780	CE
	100 µg/ml	1 ml	ref. P839790	CE
iso-octane	10 µg/ml	10 ml	ref. P839800	CE
	100 µg/ml	1 ml	ref. P839810	CE
PCB 53 (2,2',5,6-Tetrachlorobiphenyl)			CAS 41464-41-9	Price
n-hexane	10 µg/ml	10 ml	ref. P839840	BH
	100 µg/ml	1 ml	ref. P839850	BH
PCB 54 (2,2',6,6'-Tetrachlorobiphenyl)			CAS 15968-05-5	Price
n-hexane	10 µg/ml	10 ml	ref. P839880	BH
	100 µg/ml	1 ml	ref. P839890	BH
PCB 55 (2,3,3',4-Tetrachlorobiphenyl)			CAS 74338-24-2	Price
n-hexane	10 µg/ml	10 ml	ref. P839920	DC
	100 µg/ml	1 ml	ref. P839930	DC
PCB 56 (2,3,3',4'-Tetrachlorobiphenyl)			CAS 41464-43-1	Price
n-hexane	10 µg/ml	10 ml	ref. P839960	CE
	100 µg/ml	1 ml	ref. P839970	CE
PCB 57 (2,3,3',5-Tetrachlorobiphenyl)			CAS 70424-67-8	Price
n-hexane	10 µg/ml	10 ml	ref. P840000	EA
	100 µg/ml	1 ml	ref. P840010	EA
PCB 58 (2,3,3',5'-Tetrachlorobiphenyl)			CAS 41464-49-7	Price
n-hexane	10 µg/ml	10 ml	ref. P840040	DC
	100 µg/ml	1 ml	ref. P840050	DC
PCB 59 (2,3,3',6-Tetrachlorobiphenyl)			CAS 74472-33-6	Price
n-hexane	10 µg/ml	10 ml	ref. P840080	FA
	100 µg/ml	1 ml	ref. P840090	FA
PCB 60 (2,3,4,4'-Tetrachlorobiphenyl)			CAS 33025-41-1	Price
n-hexane	10 µg/ml	10 ml	ref. P840120	CE
	100 µg/ml	1 ml	ref. P840130	CE
PCB 61 (2,3,4,5-Tetrachlorobiphenyl)			CAS 33284-53-6	Price
n-hexane	10 µg/ml	10 ml	ref. P840160	CE
	100 µg/ml	1 ml	ref. P840170	CE
PCB 62 (2,3,4,6-Tetrachlorobiphenyl)			CAS 54230-22-7	Price
n-hexane	10 µg/ml	10 ml	ref. P840200	DC
	100 µg/ml	1 ml	ref. P840210	DC
PCB 63 (2,3,4',5-Tetrachlorobiphenyl)			CAS 74472-34-7	Price
n-hexane	10 µg/ml	10 ml	ref. P840240	EA
	100 µg/ml	1 ml	ref. P840250	EA
PCB 64 (2,3,4',6-Tetrachlorobiphenyl)			CAS 52663-58-8	Price
n-hexane	10 µg/ml	10 ml	ref. P840280	FA
	100 µg/ml	1 ml	ref. P840290	FA
PCB 65 (2,3,5,6-Tetrachlorobiphenyl)			CAS 33284-54-7	Price
n-hexane	10 µg/ml	10 ml	ref. P840320	CE
	100 µg/ml	1 ml	ref. P840330	CE

PCB 66 (2,3',4,4'-Tetrachlorobiphenyl)			CAS 32598-10-0	Price
n-hexane	10 µg/ml	10 ml	ref. P840360	CE
	100 µg/ml	1 ml	ref. P840370	CE
iso-octane	10 µg/ml	10 ml	ref. P840380	CE
	100 µg/ml	1 ml	ref. P840390	CE
PCB 67 (2,3',4,5-Tetrachlorobiphenyl)			CAS 73557-53-8	Price
n-hexane	10 µg/ml	10 ml	ref. P840420	FA
	100 µg/ml	1 ml	ref. P840430	FA
PCB 68 (2,3',4,5'-Tetrachlorobiphenyl)			CAS 73575-52-7	Price
n-hexane	10 µg/ml	10 ml	ref. P840460	DC
	100 µg/ml	1 ml	ref. P840470	DC
PCB 69 (2,3',4,6-Tetrachlorobiphenyl)			CAS 60233-24-1	Price
n-hexane	10 µg/ml	10 ml	ref. P840500	DC
	100 µg/ml	1 ml	ref. P840510	DC
PCB 70 (2,3',4',5-Tetrachlorobiphenyl)			CAS 32598-11-1	Price
n-hexane	10 µg/ml	10 ml	ref. P840540	CE
	100 µg/ml	1 ml	ref. P840550	CE
PCB 71 (2,3',4',6-Tetrachlorobiphenyl)			CAS 41464-46-4	Price
n-hexane	10 µg/ml	10 ml	ref. P840580	DC
	100 µg/ml	1 ml	ref. P840590	DC
PCB 72 (2,3',5,5'-Tetrachlorobiphenyl)			CAS 41464-42-0	Price
n-hexane	10 µg/ml	10 ml	ref. P840620	CE
	100 µg/ml	1 ml	ref. P840630	CE
PCB 73 (2,3',5',6-Tetrachlorobiphenyl)			CAS 74338-23-1	Price
n-hexane	10 µg/ml	10 ml	ref. P840660	EA
	100 µg/ml	1 ml	ref. P840670	EA
PCB 74 (2,4,4',5-Tetrachlorobiphenyl)			CAS 32690-93-0	Price
n-hexane	10 µg/ml	10 ml	ref. P840700	DC
	100 µg/ml	1 ml	ref. P840710	DC
PCB 75 (2,4,4',6-Tetrachlorobiphenyl)			CAS 32598-12-2	Price
n-hexane	10 µg/ml	10 ml	ref. P840740	DC
	100 µg/ml	1 ml	ref. P840750	DC
PCB 76 (2',3,4,5-Tetrachlorobiphenyl)			CAS 70362-48-0	Price
n-hexane	10 µg/ml	10 ml	ref. P840780	DC
	100 µg/ml	1 ml	ref. P840790	DC
PCB 77 (3,3',4,4'-Tetrachlorobiphenyl)			CAS 32598-13-3	Price
n-hexane	10 µg/ml	10 ml	ref. P840820	BG
	100 µg/ml	1 ml	ref. P840830	BG
iso-octane	10 µg/ml	10 ml	ref. P840840	BG
	100 µg/ml	1 ml	ref. P840850	BG
PCB 78 (3,3',4,5-Tetrachlorobiphenyl)			CAS 70362-49-1	Price
n-hexane	10 µg/ml	10 ml	ref. P840880	DC
	100 µg/ml	1 ml	ref. P840890	DC
PCB 79 (3,3',4,5'-Tetrachlorobiphenyl)			CAS 41464-48-6	Price
n-hexane	10 µg/ml	10 ml	ref. P840920	DC
	100 µg/ml	1 ml	ref. P840930	DC
PCB 80 (3,3',5,5'-Tetrachlorobiphenyl)			CAS 33284-52-5	Price
n-hexane	10 µg/ml	10 ml	ref. P840960	DC
	100 µg/ml	1 ml	ref. P840970	DC
PCB 81 (3,4,4',5-Tetrachlorobiphenyl)			CAS 70362-50-4	Price
n-hexane	10 µg/ml	10 ml	ref. P841000	CE
	100 µg/ml	1 ml	ref. P841010	CE
PCB 82 (2,2',3,3',4-Pentachlorobiphenyl)			CAS 52663-62-4	Price
n-hexane	10 µg/ml	10 ml	ref. P841040	DC
	100 µg/ml	1 ml	ref. P841050	DC
PCB 83 (2,2',3,3',5-Pentachlorobiphenyl)			CAS 60145-20-2	Price
n-hexane	10 µg/ml	10 ml	ref. P841080	DC
	100 µg/ml	1 ml	ref. P841090	DC
PCB 84 (2,2',3,3',6-Pentachlorobiphenyl)			CAS 52663-60-2	Price
n-hexane	10 µg/ml	10 ml	ref. P841120	DC
	100 µg/ml	1 ml	ref. P841130	DC
PCB 85 (2,2',3,4,4'-Pentachlorobiphenyl)			CAS 65510-45-4	Price
n-hexane	10 µg/ml	10 ml	ref. P841160	CE
	100 µg/ml	1 ml	ref. P841170	CE
PCB 86 (2,2',3,4,5-Pentachlorobiphenyl)			CAS 55312-69-1	Price
n-hexane	10 µg/ml	10 ml	ref. P841200	DC
	100 µg/ml	1 ml	ref. P841210	DC

PCB 87 (2,2',3,4,5'-Pentachlorobiphenyl)			CAS 38380-02-8	Price
n-hexane	10 µg/ml	10 ml	ref. P841240	CE
	100 µg/ml	1 ml	ref. P841250	CE
PCB 88 (2,2',3,4,6-Pentachlorobiphenyl)			CAS 55215-17-3	Price
n-hexane	10 µg/ml	10 ml	ref. P841280	DC
	100 µg/ml	1 ml	ref. P841290	DC
PCB 89 (2,2',3,4,6'-Pentachlorobiphenyl)			CAS 73575-57-2	Price
n-hexane	10 µg/ml	10 ml	ref. P841320	DC
	100 µg/ml	1 ml	ref. P841330	DC
PCB 90 (2,2',3,4',5-Pentachlorobiphenyl)			CAS 68194-07-0	Price
n-hexane	10 µg/ml	10 ml	ref. P841360	EA
	100 µg/ml	1 ml	ref. P841370	EA
PCB 91 (2,2',3,4',6-Pentachlorobiphenyl)			CAS 68194-05-8	Price
n-hexane	10 µg/ml	10 ml	ref. P841400	DC
	100 µg/ml	1 ml	ref. P841410	DC
PCB 92 (2,2',3,5,5'-Pentachlorobiphenyl)			CAS 52663-61-3	Price
n-hexane	10 µg/ml	10 ml	ref. P841440	DC
	100 µg/ml	1 ml	ref. P841450	DC
PCB 93 (2,2',3,5,6-Pentachlorobiphenyl)			CAS 73575-56-1	Price
n-hexane	10 µg/ml	10 ml	ref. P841480	DC
	100 µg/ml	1 ml	ref. P841490	DC
PCB 94 (2,2',3,5,6'-Pentachlorobiphenyl)			CAS 73575-55-0	Price
n-hexane	10 µg/ml	10 ml	ref. P841520	FA
	100 µg/ml	1 ml	ref. P841530	FA
PCB 95 (2,2',3,5',6-Pentachlorobiphenyl)			CAS 38379-99-6	Price
n-hexane	10 µg/ml	10 ml	ref. P841560	DC
	100 µg/ml	1 ml	ref. P841570	DC
PCB 96 (2,2',3,6,6'-Pentachlorobiphenyl)			CAS 73575-54-9	Price
n-hexane	10 µg/ml	10 ml	ref. P841590	DC
	100 µg/ml	1 ml	ref. P841600	DC
PCB 97 (2,2',3',4,5-Pentachlorobiphenyl)			CAS 41464-51-1	Price
n-hexane	10 µg/ml	10 ml	ref. P841630	CE
	100 µg/ml	1 ml	ref. P841640	CE
PCB 98 (2,2',3',4,6-Pentachlorobiphenyl)			CAS 60233-25-2	Price
n-hexane	10 µg/ml	10 ml	ref. P841670	DC
	100 µg/ml	1 ml	ref. P841680	DC
PCB 99 (2,2',4,4',5-Pentachlorobiphenyl)			CAS 38380-01-7	Price
n-hexane	10 µg/ml	10 ml	ref. P841710	DC
	100 µg/ml	1 ml	ref. P841720	DC
PCB 100 (2,2',4,4',6-Pentachlorobiphenyl)			CAS 39485-83-1	Price
n-hexane	10 µg/ml	10 ml	ref. P841750	DC
	100 µg/ml	1 ml	ref. P841760	DC
PCB 101 (2,2',4,5,5'-Pentachlorobiphenyl)			CAS 37680-73-2	Price
n-hexane	10 µg/ml	10 ml	ref. P841790	CE
	100 µg/ml	1 ml	ref. P841800	CE
iso-octane	10 µg/ml	10 ml	ref. P841810	CE
	100 µg/ml	1 ml	ref. P841820	CE
PCB 102 (2,2',4,5,6'-Pentachlorobiphenyl)			CAS 68194-06-9	Price
n-hexane	10 µg/ml	10 ml	ref. P841850	DC
	100 µg/ml	1 ml	ref. P841860	DC
PCB 103 (2,2',4,5',6-Pentachlorobiphenyl)			CAS 60145-21-3	Price
n-hexane	10 µg/ml	10 ml	ref. P841890	CE
	100 µg/ml	1 ml	ref. P841900	CE
PCB 104 (2,2',4,6,6'-Pentachlorobiphenyl)			CAS 56558-16-8	Price
n-hexane	10 µg/ml	10 ml	ref. P841930	DC
	100 µg/ml	1 ml	ref. P841940	DC
PCB 105 (2,3,3',4,4'-Pentachlorobiphenyl)			CAS 32598-14-4	Price
n-hexane	10 µg/ml	10 ml	ref. P841970	CE
	100 µg/ml	1 ml	ref. P841980	CE
iso-octane	10 µg/ml	10 ml	ref. P841990	CE
	100 µg/ml	1 ml	ref. P842000	CE
PCB 106 (2,3,3',4,5-Pentachlorobiphenyl)			CAS 70424-69-0	Price
n-hexane	10 µg/ml	10 ml	ref. P842030	DC
	100 µg/ml	1 ml	ref. P842040	DC
PCB 107 (2,3,3',4',5-Pentachlorobiphenyl)			CAS 70424-68-9	Price
n-hexane	10 µg/ml	10 ml	ref. P842070	DC
	100 µg/ml	1 ml	ref. P842080	DC

PCB 108 (2,3,3',4,5'-Pentachlorobiphenyl)			CAS 70362-41-3	Price
n-hexane	10 µg/ml	10 ml	ref. P842110	DC
	100 µg/ml	1 ml	ref. P842120	DC
PCB 109 (2,3,3',4,6-Pentachlorobiphenyl)			CAS 74472-35-8	Price
n-hexane	10 µg/ml	10 ml	ref. P842150	EA
	100 µg/ml	1 ml	ref. P842160	EA
PCB 110 (2,3,3',4',6-Pentachlorobiphenyl)			CAS 38380-03-9	Price
n-hexane	10 µg/ml	10 ml	ref. P842190	DC
	100 µg/ml	1 ml	ref. P842200	DC
PCB 111 (2,3,3',5,5'-Pentachlorobiphenyl)			CAS 39635-32-0	Price
n-hexane	10 µg/ml	10 ml	ref. P842230	FA
	100 µg/ml	1 ml	ref. P842240	FA
PCB 112 (2,3,3',5,6-Pentachlorobiphenyl)			CAS 74472-36-9	Price
n-hexane	10 µg/ml	10 ml	ref. P842270	DC
	100 µg/ml	1 ml	ref. P842280	DC
PCB 113 (2,3,3',5'6-Pentachlorobiphenyl)			CAS 68194-10-5	Price
n-hexane	10 µg/ml	10 ml	ref. P842310	FA
	100 µg/ml	1 ml	ref. P842320	FA
PCB 114 (2,3,4,4',5-Pentachlorobiphenyl)			CAS 74472-37-0	Price
n-hexane	10 µg/ml	10 ml	ref. P842350	CE
	100 µg/ml	1 ml	ref. P842360	CE
PCB 115 (2,3,4,4',6-Pentachlorobiphenyl)			CAS 74472-38-1	Price
n-hexane	10 µg/ml	10 ml	ref. P842390	DC
	100 µg/ml	1 ml	ref. P842400	DC
PCB 116 (2,3,4,5,6-Pentachlorobiphenyl)			CAS 18259-05-7	Price
n-hexane	10 µg/ml	10 ml	ref. P842430	CE
	100 µg/ml	1 ml	ref. P842440	CE
PCB 117 (2,3,4',5,6-Pentachlorobiphenyl)			CAS 68194-11-6	Price
n-hexane	10 µg/ml	10 ml	ref. P842470	DC
	100 µg/ml	1 ml	ref. P842480	DC
PCB 118 (2,3',4,4',5-Pentachlorobiphenyl)			CAS 31508-00-6	Price
n-hexane	10 µg/ml	10 ml	ref. P842510	CE
	100 µg/ml	1 ml	ref. P842520	CE
iso-octane	10 µg/ml	10 ml	ref. P842530	CE
	100 µg/ml	1 ml	ref. P842540	CE
PCB 119 (2,3',4,4',6-Pentachlorobiphenyl)			CAS 56558-17-9	Price
n-hexane	10 µg/ml	10 ml	ref. P842570	DC
	100 µg/ml	1 ml	ref. P842580	DC
PCB 120 (2,3',4,5,5'-Pentachlorobiphenyl)			CAS 68194-12-7	Price
n-hexane	10 µg/ml	10 ml	ref. P842610	DC
	100 µg/ml	1 ml	ref. P842620	DC
PCB 121 (2,3',4,5',6-Pentachlorobiphenyl)			CAS 56558-18-0	Price
n-hexane	10 µg/ml	10 ml	ref. P842650	DC
	100 µg/ml	1 ml	ref. P842660	DC
PCB 122 (2',3,3',4,5-Pentachlorobiphenyl)			CAS 76842-07-4	Price
n-hexane	10 µg/ml	10 ml	ref. P842690	DC
	100 µg/ml	1 ml	ref. P842700	DC
PCB 123 (2',3,4,4',5-Pentachlorobiphenyl)			CAS 65510-44-3	Price
n-hexane	10 µg/ml	10 ml	ref. P842730	DC
	100 µg/ml	1 ml	ref. P842740	DC
PCB 124 (2',3,4,5,5'-Pentachlorobiphenyl)			CAS 70424-70-3	Price
n-hexane	10 µg/ml	10 ml	ref. P842770	DC
	100 µg/ml	1 ml	ref. P842780	DC
PCB 125 (2',3,4,5,6'-Pentachlorobiphenyl)			CAS 74472-39-2	Price
n-hexane	10 µg/ml	10 ml	ref. P842810	DC
	100 µg/ml	1 ml	ref. P842820	DC
PCB 126 (3,3',4,4',5-Pentachlorobiphenyl)			CAS 57465-28-8	Price
n-hexane	10 µg/ml	10 ml	ref. P842850	CE
	100 µg/ml	1 ml	ref. P842860	CE
iso-octane	10 µg/ml	10 ml	ref. P842870	CE
	100 µg/ml	1 ml	ref. P842880	CE
PCB 127 (3,3',4,5,5'-Pentachlorobiphenyl)			CAS 39635-33-1	Price
n-hexane	10 µg/ml	10 ml	ref. P842910	DC
	100 µg/ml	1 ml	ref. P842920	DC
PCB 128 (2,2',3,3',4,4'-Hexachlorobiphenyl)			CAS 38380-07-3	Price
n-hexane	10 µg/ml	10 ml	ref. P842950	CE
	100 µg/ml	1 ml	ref. P842960	CE
iso-octane	10 µg/ml	10 ml	ref. P842970	CE
	100 µg/ml	1 ml	ref. P842980	CE

PCB 129 (2,2',3,3',4,5-Hexachlorobiphenyl)			CAS 55215-18-4	Price
n-hexane	10 µg/ml	10 ml	ref. P843010	DC
	100 µg/ml	1 ml	ref. P843020	DC
PCB 130 (2,2',3,3',4,5'-Hexachlorobiphenyl)			CAS 52663-68-8	Price
n-hexane	10 µg/ml	10 ml	ref. P843050	DC
	100 µg/ml	1 ml	ref. P843060	DC
PCB 131 (2,2',3,3',4,6-Hexachlorobiphenyl)			CAS 61798-70-7	Price
n-hexane	10 µg/ml	10 ml	ref. P843090	DC
	100 µg/ml	1 ml	ref. P843100	DC
PCB 132 (2,2',3,3',4,6'-Hexachlorobiphenyl)			CAS 38380-05-1	Price
n-hexane	10 µg/ml	10 ml	ref. P843130	DC
	100 µg/ml	1 ml	ref. P843140	DC
PCB 133 (2,2',3,3',5,5'-Hexachlorobiphenyl)			CAS 35694-04-3	Price
n-hexane	10 µg/ml	10 ml	ref. P843170	EA
	100 µg/ml	1 ml	ref. P843180	EA
PCB 134 (2,2',3,3',5,6-Hexachlorobiphenyl)			CAS 52704-70-8	Price
n-hexane	10 µg/ml	10 ml	ref. P843210	DC
	100 µg/ml	1 ml	ref. P843220	DC
PCB 135 (2,2',3,3',5,6'-Hexachlorobiphenyl)			CAS 52744-13-5	Price
n-hexane	10 µg/ml	10 ml	ref. P843250	FA
	100 µg/ml	1 ml	ref. P843260	FA
PCB 136 (2,2',3,3',6,6'-Hexachlorobiphenyl)			CAS 38411-22-2	Price
n-hexane	10 µg/ml	10 ml	ref. P843290	CE
	100 µg/ml	1 ml	ref. P843300	CE
PCB 137 (2,2',3,4,4',5-Hexachlorobiphenyl)			CAS 35694-06-5	Price
n-hexane	10 µg/ml	10 ml	ref. P843330	DC
	100 µg/ml	1 ml	ref. P843340	DC
PCB 138 (2,2',3,4,4',5'-Hexachlorobiphenyl)			CAS 35065-28-2	Price
n-hexane	10 µg/ml	10 ml	ref. P843370	CE
	100 µg/ml	1 ml	ref. P843380	CE
iso-octane	10 µg/ml	10 ml	ref. P843390	CE
	100 µg/ml	1 ml	ref. P843400	CE
PCB 139 (2,2',3,4,4',6-Hexachlorobiphenyl)			CAS 56030-56-9	Price
n-hexane	10 µg/ml	10 ml	ref. P843430	EA
	100 µg/ml	1 ml	ref. P843440	EA
PCB 140 (2,2',3,4,4',6'-Hexachlorobiphenyl)			CAS 59291-64-4	Price
n-hexane	10 µg/ml	10 ml	ref. P843470	DC
	100 µg/ml	1 ml	ref. P843480	DC
PCB 141 (2,2',3,4,5,5'-Hexachlorobiphenyl)			CAS 52712-04-6	Price
n-hexane	10 µg/ml	10 ml	ref. P843510	CE
	100 µg/ml	1 ml	ref. P843520	CE
PCB 142 (2,2',3,4,5,6-Hexachlorobiphenyl)			CAS 41411-61-4	Price
n-hexane	10 µg/ml	10 ml	ref. P843550	EA
	100 µg/ml	1 ml	ref. P843560	EA
PCB 143 (2,2',3,4,5,6'-Hexachlorobiphenyl)			CAS 68194-15-0	Price
n-hexane	10 µg/ml	10 ml	ref. P843590	DC
	100 µg/ml	1 ml	ref. P843600	DC
	400 µg/ml	1 ml	ref. F131421	FA
PCB 144 (2,2',3,4,5',6-Hexachlorobiphenyl)			CAS 68194-14-9	Price
n-hexane	10 µg/ml	10 ml	ref. P843630	EA
	100 µg/ml	1 ml	ref. P843640	EA
PCB 145 (2,2',3,4,6,6'-Hexachlorobiphenyl)			CAS 74472-40-5	Price
n-hexane	10 µg/ml	10 ml	ref. P843670	DC
	100 µg/ml	1 ml	ref. P843680	DC
PCB 146 (2,2',3,4',5,5'-Hexachlorobiphenyl)			CAS 51908-16-8	Price
n-hexane	10 µg/ml	10 ml	ref. P843710	EA
	100 µg/ml	1 ml	ref. P843720	EA
PCB 147 (2,2',3,4',5,6-Hexachlorobiphenyl)			CAS 68194-13-8	Price
n-hexane	10 µg/ml	10 ml	ref. P843750	DC
	100 µg/ml	1 ml	ref. P843760	DC
PCB 148 (2,2',3,4',5,6'-Hexachlorobiphenyl)			CAS 74472-41-6	Price
n-hexane	10 µg/ml	10 ml	ref. P843790	EA
	100 µg/ml	1 ml	ref. P843800	EA
PCB 149 (2,2',3,4',5',6-Hexachlorobiphenyl)			CAS 38380-04-0	Price
n-hexane	10 µg/ml	10 ml	ref. P843830	EA
	100 µg/ml	1 ml	ref. P843840	EA
iso-octane	10 µg/ml	10 ml	ref. P843850	EA
	100 µg/ml	1 ml	ref. P843860	EA

PCB 150 (2,2',3,4',6,6'-Hexachlorobiphenyl)			CAS 68194-08-1	Price
n-hexane	10 µg/ml	10 ml	ref. P843890	EA
	100 µg/ml	1 ml	ref. P843900	EA
PCB 151 (2,2',3,5,5',6-Hexachlorobiphenyl)			CAS 52663-63-5	Price
n-hexane	10 µg/ml	10 ml	ref. P843930	DC
	100 µg/ml	1 ml	ref. P843940	DC
PCB 152 (2,2',3,5,6,6'-Hexachlorobiphenyl)			CAS 68194-09-2	Price
n-hexane	10 µg/ml	10 ml	ref. P843970	DC
	100 µg/ml	1 ml	ref. P843980	DC
PCB 153 (2,2',4,4',5,5'-Hexachlorobiphenyl)			CAS 35065-27-1	Price
n-hexane	10 µg/ml	10 ml	ref. P844010	CE
	100 µg/ml	1 ml	ref. P844020	CE
iso-octane	10 µg/ml	10 ml	ref. P844030	CE
	100 µg/ml	1 ml	ref. P844040	CE
PCB 154 (2,2',4,4',5,6'-Hexachlorobiphenyl)			CAS 60145-22-4	Price
n-hexane	10 µg/ml	10 ml	ref. P844070	DC
	100 µg/ml	1 ml	ref. P844080	DC
PCB 155 (2,2',4,4',6,6'-Hexachlorobiphenyl)			CAS 33979-03-2	Price
n-hexane	10 µg/ml	10 ml	ref. P844110	CE
	100 µg/ml	1 ml	ref. P844120	CE
	400 µg/ml	1 ml	ref. F131431	CE
PCB 156 (2,3,3',4,4',5-Hexachlorobiphenyl)			CAS 38380-08-4	Price
n-hexane	10 µg/ml	10 ml	ref. P844150	CE
	100 µg/ml	1 ml	ref. P844160	CE
PCB 157 (2,3,3',4,4',5'-Hexachlorobiphenyl)			CAS 69782-90-7	Price
n-hexane	10 µg/ml	10 ml	ref. P844190	DC
	100 µg/ml	1 ml	ref. P844200	DC
PCB 158 (2,3,3',4,4',6-Hexachlorobiphenyl)			CAS 74472-42-7	Price
n-hexane	10 µg/ml	10 ml	ref. P844230	DC
	100 µg/ml	1 ml	ref. P844240	DC
PCB 159 (2,3,3',4,5,5'-Hexachlorobiphenyl)			CAS 39635-35-3	EA
n-hexane	10 µg/ml	10 ml	ref. P844270	EA
	100 µg/ml	1 ml	ref. P844280	EA
PCB 160 (2,3,3',4,5,6-Hexachlorobiphenyl)			CAS 41411-62-5	Price
n-hexane	10 µg/ml	10 ml	ref. P844310	DC
	100 µg/ml	1 ml	ref. P844320	DC
PCB 161 (2,3,3',4,5',6-Hexachlorobiphenyl)			CAS 74472-43-8	Price
n-hexane	10 µg/ml	10 ml	ref. P844350	DC
	100 µg/ml	1 ml	ref. P844360	DC
PCB 162 (2,3,3',4',5,5'-Hexachlorobiphenyl)			CAS 39635-34-2	Price
n-hexane	10 µg/ml	10 ml	ref. P844390	EA
	100 µg/ml	1 ml	ref. P844400	EA
PCB 163 (2,3,3',4',5,6-Hexachlorobiphenyl)			CAS 74472-44-9	Price
n-hexane	10 µg/ml	10 ml	ref. P844430	DC
	100 µg/ml	1 ml	ref. P844440	DC
PCB 164 (2,3,3',4',5',6-Hexachlorobiphenyl)			CAS 74472-45-0	Price
n-hexane	10 µg/ml	10 ml	ref. P844470	DC
	100 µg/ml	1 ml	ref. P844480	DC
PCB 165 (2,3,3',5,5',6-Hexachlorobiphenyl)			CAS 74472-46-1	Price
n-hexane	10 µg/ml	10 ml	ref. P844510	DC
	100 µg/ml	1 ml	ref. P844520	DC
PCB 166 (2,3,4,4',5,6-Hexachlorobiphenyl)			CAS 41411-63-6	Price
n-hexane	10 µg/ml	10 ml	ref. P844550	EA
	100 µg/ml	1 ml	ref. P844560	EA
PCB 167 (2,3',4,4',5,5'-Hexachlorobiphenyl)			CAS 52663-72-6	Price
n-hexane	10 µg/ml	10 ml	ref. P844590	CE
	100 µg/ml	1 ml	ref. P844600	CE
PCB 168 (2,3',4,4',5',6-Hexachlorobiphenyl)			CAS 59291-65-5	Price
n-hexane	10 µg/ml	10 ml	ref. P844630	DC
	100 µg/ml	1 ml	ref. P844640	DC
PCB 169 (3,3',4,4',5,5'-Hexachlorobiphenyl)			CAS 32774-16-6	Price
n-hexane	10 µg/ml	10 ml	ref. P844670	CE
	100 µg/ml	1 ml	ref. P844680	CE
PCB 170 (2,2',3,3',4,4',5-Heptachlorobiphenyl)			CAS 35065-30-6	Price
n-hexane	10 µg/ml	10 ml	ref. P844710	DC
	100 µg/ml	1 ml	ref. P844720	DC
iso-octane	10 µg/ml	10 ml	ref. P844730	DC
	100 µg/ml	1 ml	ref. P844740	DC

PCB 171 (2,2',3,3',4,4',6-Heptachlorobiphenyl)				CAS 52663-71-5	Price
n-hexane	10 µg/ml	10 ml	ref. P844770	DC	
	100 µg/ml	1 ml	ref. P844780	DC	
PCB 172 (2,2',3,3',4,5,5'-Heptachlorobiphenyl)				CAS 52663-74-8	Price
n-hexane	10 µg/ml	10 ml	ref. P844810	EA	
	100 µg/ml	1 ml	ref. P844820	EA	
PCB 173 (2,2',3,3',4,5,6-Heptachlorobiphenyl)				CAS 68194-16-1	Price
n-hexane	10 µg/ml	10 ml	ref. P844850	FA	
	100 µg/ml	1 ml	ref. P844860	FA	
PCB 174 (2,2',3,3',4,5,6'-Heptachlorobiphenyl)				CAS 38411-25-5	Price
n-hexane	10 µg/ml	10 ml	ref. P844890	DC	
	100 µg/ml	1 ml	ref. P844900	DC	
PCB 175 (2,2',3,3',4,5',6-Heptachlorobiphenyl)				CAS 40186-70-7	Price
n-hexane	10 µg/ml	10 ml	ref. P844930	EA	
	100 µg/ml	1 ml	ref. P844940	EA	
PCB 176 (2,2',3,3',4,6,6'-Heptachlorobiphenyl)				CAS 52663-65-7	Price
n-hexane	10 µg/ml	10 ml	ref. P844970	DC	
	100 µg/ml	1 ml	ref. P844980	DC	
PCB 177 (2,2',3,3',4',5,6-Heptachlorobiphenyl)				CAS 52663-70-4	Price
n-hexane	10 µg/ml	10 ml	ref. P845010	DC	
	100 µg/ml	1 ml	ref. P845020	DC	
PCB 178 (2,2',3,3',5,5',6-Heptachlorobiphenyl)				CAS 52663-67-9	Price
n-hexane	10 µg/ml	10 ml	ref. P845050	EA	
	100 µg/ml	1 ml	ref. P845060	EA	
PCB 179 (2,2',3,3',5,6,6'-Heptachlorobiphenyl)				CAS 52663-64-6	Price
n-hexane	10 µg/ml	10 ml	ref. P845090	DC	
	100 µg/ml	1 ml	ref. P845100	DC	
PCB 180 (2,2',3,4,4',5,5'-Heptachlorobiphenyl)				CAS 35065-29-3	Price
n-hexane	10 µg/ml	10 ml	ref. P845130	CE	
	100 µg/ml	1 ml	ref. P845140	CE	
iso-octane	10 µg/ml	10 ml	ref. P845150	CE	
	100 µg/ml	1 ml	ref. P845160	CE	
PCB 181 (2,2',3,4,4',5,6-Heptachlorobiphenyl)				CAS 74472-47-2	Price
n-hexane	10 µg/ml	10 ml	ref. P845190	EA	
	100 µg/ml	1 ml	ref. P845200	EA	
PCB 182 (2,2',3,4,4',5,6'-Heptachlorobiphenyl)				CAS 60145-23-5	Price
n-hexane	10 µg/ml	10 ml	ref. P845230	DC	
	100 µg/ml	1 ml	ref. P845240	DC	
PCB 183 (2,2',3,4,4',5',6-Heptachlorobiphenyl)				CAS 52663-69-1	Price
n-hexane	10 µg/ml	10 ml	ref. P845270	DC	
	100 µg/ml	1 ml	ref. P845280	DC	
PCB 184 (2,2',3,4,4',6,6'-Heptachlorobiphenyl)				CAS 74472-48-3	Price
n-hexane	10 µg/ml	10 ml	ref. P845310	DC	
	100 µg/ml	1 ml	ref. P845320	DC	
PCB 185 (2,2',3,4,5,5',6-Heptachlorobiphenyl)				CAS 52712-05-7	Price
n-hexane	10 µg/ml	10 ml	ref. P845350	DC	
	100 µg/ml	1 ml	ref. P845360	DC	
PCB 186 (2,2',3,4,5,6,6'-Heptachlorobiphenyl)				CAS 74472-49-4	Price
n-hexane	10 µg/ml	10 ml	ref. P845390	DC	
	100 µg/ml	1 ml	ref. P845400	DC	
PCB 187 (2,2',3,4',5,5',6-Heptachlorobiphenyl)				CAS 52663-68-0	Price
n-hexane	10 µg/ml	10 ml	ref. P845430	CE	
	100 µg/ml	1 ml	ref. P845440	CE	
iso-octane	10 µg/ml	10 ml	ref. P845450	CE	
	100 µg/ml	1 ml	ref. P845460	CE	
PCB 188 (2,2',3,4',5,6,6'-Heptachlorobiphenyl)				CAS 74487-85-7	Price
n-hexane	10 µg/ml	10 ml	ref. P845490	DC	
	100 µg/ml	1 ml	ref. P845500	DC	
PCB 189 (2,3,3',4,4',5,5'-Heptachlorobiphenyl)				CAS 39635-31-9	Price
n-hexane	10 µg/ml	10 ml	ref. P845530	DC	
	100 µg/ml	1 ml	ref. P845540	DC	
PCB 190 (2,3,3',4,4',5,6-Heptachlorobiphenyl)				CAS 41411-64-7	Price
n-hexane	10 µg/ml	10 ml	ref. P845570	DC	
	100 µg/ml	1 ml	ref. P845580	DC	

PCB 191 (2,3,3',4,4',5',6-Heptachlorobiphenyl)				CAS 74472-50-7	Price
n-hexane	10 µg/ml	10 ml	ref. P845610	DC	
	100 µg/ml	1 ml	ref. P845620	DC	
PCB 192 (2,3,3',4,5,5',6-Heptachlorobiphenyl)				CAS 74472-51-8	Price
n-hexane	10 µg/ml	10 ml	ref. P845650	DC	
	100 µg/ml	1 ml	ref. P845660	DC	
PCB 193 (2,3,3',4',5,5',6-Heptachlorobiphenyl)				CAS 69782-91-8	Price
n-hexane	10 µg/ml	10 ml	ref. P845690	DC	
	100 µg/ml	1 ml	ref. P845700	DC	
PCB 194 (2,2',3,3',4,4',5,5'-Octachlorobiphenyl)				CAS 35694-08-7	Price
n-hexane	10 µg/ml	10 ml	ref. P845730	CE	
	100 µg/ml	1 ml	ref. P845740	CE	
iso-octane	10 µg/ml	10 ml	ref. P845750	CE	
	100 µg/ml	1 ml	ref. P845760	CE	
PCB 195 (2,2',3,3',4,4',5,6-Octachlorobiphenyl)				CAS 52663-78-2	Price
n-hexane	10 µg/ml	10 ml	ref. P845790	DC	
	100 µg/ml	1 ml	ref. P845800	DC	
iso-octane	10 µg/ml	10 ml	ref. P845810	DC	
	100 µg/ml	1 ml	ref. P845820	DC	
PCB 196 (2,2',3,3',4,4',5,6'-Octachlorobiphenyl)				CAS 42740-50-1	Price
n-hexane	10 µg/ml	10 ml	ref. P845850	EA	
	100 µg/ml	1 ml	ref. P845860	EA	
PCB 197 (2,2',3,3',4,4',6,6'-Octachlorobiphenyl)				CAS 33091-17-7	Price
n-hexane	10 µg/ml	10 ml	ref. P845890	EA	
	100 µg/ml	1 ml	ref. P845900	EA	
PCB 198 (2,2',3,3',4,5,5',6-Octachlorobiphenyl)				CAS 68194-17-2	Price
n-hexane	10 µg/ml	10 ml	ref. P845930	EA	
	100 µg/ml	1 ml	ref. P845940	EA	
	400 µg/ml	1 ml	ref. F131441	EA	
PCB 199 (2,2',3,3',4,5,6,6'-Octachlorobiphenyl)				CAS 52663-73-7	Price
n-hexane	10 µg/ml	10 ml	ref. P845970	EA	
	100 µg/ml	1 ml	ref. P845980	EA	
PCB 200 (2,2',3,3',4,5',6,6'-Octachlorobiphenyl)				CAS 40186-71-8	Price
n-hexane	10 µg/ml	10 ml	ref. P846010	DC	
	100 µg/ml	1 ml	ref. P846020	DC	
PCB 201 (2,2',3,3',4,5,5',6-Octachlorobiphenyl)				CAS 52663-75-9	Price
n-hexane	10 µg/ml	10 ml	ref. P846050	DC	
	100 µg/ml	1 ml	ref. P846060	DC	
PCB 202 (2,2',3,3',5,5',6,6'-Octachlorobiphenyl)				CAS 2136-99-4	Price
n-hexane	10 µg/ml	10 ml	ref. P846090	DC	
	100 µg/ml	1 ml	ref. P846100	DC	
PCB 203 (2,2',3,4,4',5,5',6-Octachlorobiphenyl)				CAS 52663-76-0	Price
n-hexane	10 µg/ml	10 ml	ref. P846130	DC	
	100 µg/ml	1 ml	ref. P846140	DC	
PCB 204 (2,2',3,4,4',5,6,6'-Octachlorobiphenyl)				CAS 74472-52-9	Price
n-hexane	10 µg/ml	10 ml	ref. P846170	DC	
	100 µg/ml	1 ml	ref. P846180	DC	
PCB 205 (2,3,3',4,4',5,5',6-Octachlorobiphenyl)				CAS 74472-53-0	Price
n-hexane	10 µg/ml	10 ml	ref. P846210	DC	
	100 µg/ml	1 ml	ref. P846220	DC	
PCB 206 (2,2',3,3',4,4',5,5',6-Nonachlorobiphenyl)				CAS 40186-72-9	Price
n-hexane	10 µg/ml	10 ml	ref. P846250	DC	
	100 µg/ml	1 ml	ref. P846260	DC	
iso-octane	10 µg/ml	10 ml	ref. P846270	DC	
	100 µg/ml	1 ml	ref. P846280	DC	
PCB 207 (2,2',3,3',4,4',5,6,6'-Nonachlorobiphenyl)				CAS 52663-79-3	Price
n-hexane	10 µg/ml	10 ml	ref. P846310	DC	
	100 µg/ml	1 ml	ref. P846320	DC	
	400 µg/ml	1 ml	ref. F131451	FA	
PCB 208 (2,2',3,3',4,5,5',6,6'-Nonachlorobiphenyl)				CAS 52663-77-1	Price
n-hexane	10 µg/ml	10 ml	ref. P846350	DC	
	100 µg/ml	1 ml	ref. P846360	DC	

PCB 209 (2,2',3,3',4,4',5,5',6,6'-Decachlorobiphenyl)			CAS 2051-24-3	Price
acetone	10 µg/ml	10 ml	ref. P846390	BG
	100 µg/ml	1 ml	ref. P846400	BG
	500 µg/ml	1 ml	ref. F113681	BG
	1 000 µg/ml	1 ml	ref. F117311	BK
acetonitrile	10 µg/ml	10 ml	ref. P846410	BG
	100 µg/ml	1 ml	ref. P846420	BG
cyclohexane	10 µg/ml	10 ml	ref. P846430	BG
	100 µg/ml	1 ml	ref. P846440	BG
n-hexane	10 µg/ml	10 ml	ref. P846450	BG
	100 µg/ml	1 ml	ref. F113671	BG
	400 µg/ml	1 ml	ref. F131461	BG
iso-octane	10 µg/ml	10 ml	ref. P846470	BG
	100 µg/ml	1 ml	ref. P846480	BG
	10 µg/ml	1 ml	ref. F131500	BG
methyl- <i>tert</i> -butylether	5 000 µg/ml	1 ml	ref. F120141	FA
toluene	1 000 µg/ml	1 ml	ref. F113661	CE

PEBULATE			CAS 1114-71-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P877700	BG
methanol	100 µg/ml	1 ml	ref. P877690	BG

PENCONAZOLE			CAS 66246-88-6	Price
acetone	10 µg/ml	10 ml	ref. P846510	BG
	100 µg/ml	1 ml	ref. P846520	BG
acetonitrile	10 µg/ml	10 ml	ref. P846530	BG
	100 µg/ml	1 ml	ref. P846540	BG
cyclohexane	10 µg/ml	10 ml	ref. P846550	BG
	100 µg/ml	1 ml	ref. P846560	BG

PENCYCURON			CAS 66063-05-6	Price
iso-octane	10 µg/ml	10 ml	ref. P877720	BG
methanol	100 µg/ml	1 ml	ref. P877710	BG

PENDIMETHALIN			CAS 40487-42-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P846600	BG
	100 µg/ml	1 ml	ref. P846610	BG
cyclohexane	10 µg/ml	10 ml	ref. P846620	BG
	100 µg/ml	1 ml	ref. P846630	BG

PENTACHLOROANILINE			CAS 527-20-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P877740	BG
	100 µg/ml	1 ml	ref. P877750	BG

PENTACHLOROANISOLE			CAS 1825-21-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P846690	BG
	100 µg/ml	1 ml	ref. P846700	BG
methanol	10 µg/ml	10 ml	ref. P846710	BG
	100 µg/ml	1 ml	ref. P846720	BG

PENTACHLOROBENZENE			CAS 608-93-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P846790	BG
	100 µg/ml	1 ml	ref. P846800	BG
cyclohexane	10 µg/ml	10 ml	ref. P846810	BG
	100 µg/ml	1 ml	ref. P846820	BG

PENTACHLOROPHENOL			CAS 87-86-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P846890	BG
	100 µg/ml	1 ml	ref. P846900	BG
cyclohexane	10 µg/ml	10 ml	ref. P846910	BG
	100 µg/ml	1 ml	ref. P846920	BG
methanol	10 µg/ml	10 ml	ref. P846930	BG
	100 µg/ml	1 ml	ref. P846940	BG

PENTACHLOROPHENOL-13C6			CAS 85380-74-1	Price
acetone	100 µg/ml	1 ml	ref. P846970	GF
n-nonane	100 µg/ml	1 ml	ref. P846980	GF

PENTACHLOROPHENOL ACETATE			CAS 1441-02-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P877770	BH
cyclohexane	10 µg/ml	10 ml	ref. P877780	BH
	100 µg/ml	1 ml	ref. P877760	BH

PENTACONTANE			CAS 6596-40-3	Price
toluene	10 µg/ml	10 ml	ref. P847010	BH
	100 µg/ml	1 ml	ref. P847020	BH

n- PENTACOSANE			CAS 629-99-2	Price
dichloromethane	10 µg/ml	10 ml	ref. P847050	BH
	100 µg/ml	1 ml	ref. P847060	BH

PENTAFLUOROBENZENE			CAS 363-72-4	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P847160	BG
	100 µg/ml	1 ml	ref. P847170	BG

O-(2,3,4,5,6-PENTAFLUOROBENZYL) HYDROXIL-AMINE HYDROCHLORIDE			CAS 57981-02-9	Price
water	15 µg/ml	10 ml	ref. F114821	CD

PENTAFLUOROPHENOL			CAS 771-61-9	Price
acetone	100 µg/ml	1 ml	ref. F146781	BG
dichloromethane	10 µg/ml	10 ml	ref. P847200	BG
	100 µg/ml	1 ml	ref. P847210	BG
methanol	2 000 µg/ml	1 ml	ref. F116601	CA

n- PENTANE			CAS 109-66-0	Price
methanol	1 000 µg/ml	5 ml	ref. P877290	BH

2-PENTANONE			CAS 107-87-9	Price
methanol/water (9/1)	10 µg/ml	10 ml	ref. P847240	BG
	100 µg/ml	1 ml	ref. P847250	BG

3-PENTANONE			CAS 96-22-0	Price
methanol/water (9/1)	10 µg/ml	10 ml	ref. P847280	BG
	100 µg/ml	1 ml	ref. P847290	BG

PERMETHRIN			CAS 52645-53-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P847380	BG
	100 µg/ml	1 ml	ref. P847390	BG
cyclohexane	10 µg/ml	10 ml	ref. P847400	BG
	100 µg/ml	1 ml	ref. P847410	BG

cis- PERMETHRIN			CAS 61949-76-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P847480	CB
	100 µg/ml	1 ml	ref. P847490	CB
cyclohexane	10 µg/ml	10 ml	ref. P847500	CB
	100 µg/ml	1 ml	ref. P847510	CB

trans- PERMETHRIN			CAS 61949-77-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P847540	CA
	100 µg/ml	1 ml	ref. P847550	CA
cyclohexane	10 µg/ml	10 ml	ref. P847560	CA
	100 µg/ml	1 ml	ref. P847570	CA

PERTHANE			CAS 72-56-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P880350	BG
methanol	100 µg/ml	1 ml	ref. P880340	BG

PERYLENE			CAS 198-55-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P847630	BJ
	100 µg/ml	1 ml	ref. P847640	BJ
cyclohexane	10 µg/ml	10 ml	ref. P847650	BJ
	100 µg/ml	1 ml	ref. P847660	BJ

PERYLENE-D12			CAS 1520-96-3	Price
cyclohexane	100 µg/ml	1 ml	ref. P847690	CC
	2 000 µg/ml	1 ml	ref. F062485	CK
toluene	100 µg/ml	1 ml	ref. P847700	CC
	2 000 µg/ml	1 ml	ref. F062495	CK

PHENANTHRENE			CAS 85-01-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P847730	BG
	100 µg/ml	1 ml	ref. P847740	BG
cyclohexane	10 µg/ml	10 ml	ref. P847750	BG
	100 µg/ml	1 ml	ref. P847760	BG
methanol	10 µg/ml	10 ml	ref. P847770	BG
	100 µg/ml	1 ml	ref. P847780	BG

PHENANTHRENE-D10			CAS 1517-22-2	Price
acetone	100 µg/ml	1 ml	ref. P847810	CC
	2 000 µg/ml	1 ml	ref. F062505	CK
cyclohexane	100 µg/ml	1 ml	ref. P847820	CC
	2 000 µg/ml	1 ml	ref. F062515	CK
methyl- <i>tert</i> -butylether	100 µg/ml	1 ml	ref. P847830	CC
	2 000 µg/ml	1 ml	ref. F062525	CK

PHENMEDIPHAM			CAS 13684-63-4	Price
acetone	10 µg/ml	10 ml	ref. P847860	BG
	100 µg/ml	1 ml	ref. P847870	BG

PHENOL			CAS 108-95-2	Price
methanol	10 µg/ml	10 ml	ref. P847900	BG
	100 µg/ml	1 ml	ref. P847910	BG

PHENOL-D5			CAS 4165-62-2	Price
acetone	100 µg/ml	1 ml	ref. P847940	CA
methanol purge & trap	100 µg/ml	1 ml	ref. P847950	CA
	2 000 µg/ml	1 ml	ref. F116611	CG
n-nonane	100 µg/ml	1 ml	ref. P847960	CA

PHENOL-D6			CAS 13127-88-3	Price
methanol purge & trap	100 µg/ml	1 ml	ref. P847990	CC
	2 000 µg/ml	1 ml	ref. F062535	CK
methyl-tert-butylether	100 µg/ml	1 ml	ref. P848000	CC
	2 000 µg/ml	1 ml	ref. F062545	CK
PHENOTHIN			CAS 26002-80-2	Price
iso-octane	10 µg/ml	10 ml	ref. P877810	BG
methanol	100 µg/ml	1 ml	ref. P877800	BG
PHENTHOATE			CAS 2597-03-7	Price
iso-octane	10 µg/ml	10 ml	ref. P877830	BG
methanol	100 µg/ml	1 ml	ref. P877820	BG
2-PHENYLPHENOL			CAS 90-43-7	Price
acetone	100 µg/ml	1 ml	ref. P873200	BG
cyclohexane	10 µg/ml	10 ml	ref. P873210	BG
methanol	1 000 µg/ml	1 ml	ref. F116031	BK
PHORATE			CAS 298-02-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P848030	BG
	100 µg/ml	1 ml	ref. P848040	BG
PHORATE-SULFONE			CAS 2588-04-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P877840	BH
n-hexane	100 µg/ml	1 ml	ref. P877850	BH
PHORATE-SULFOXIDE			CAS 2588-03-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P877860	BG
PHOSALONE			CAS 2310-17-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P848110	BG
	100 µg/ml	1 ml	ref. P848120	BG
cyclohexane	10 µg/ml	10 ml	ref. P848130	BG
	100 µg/ml	1 ml	ref. P848140	BG
PHOSMET			CAS 732-11-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P848200	BG
	100 µg/ml	1 ml	ref. P848210	BG
cyclohexane	10 µg/ml	10 ml	ref. P848220	BG
	100 µg/ml	1 ml	ref. P848230	BG
PHOSPHAMIDON			CAS 13171-21-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P848290	BG
	100 µg/ml	1 ml	ref. P848300	BG
PHOXIM			CAS 14816-18-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P848330	BG
	100 µg/ml	1 ml	ref. P848340	BG
PHTHALIC ACID, BENZYL BUTYL ESTER			CAS 85-68-7	Price
methanol	100 µg/ml	1 ml	ref. P877880	BH
	1 000 µg/ml	1 ml	ref. P877870	CA
PHTHALIC ACID, BIS-BUTYL ESTER			CAS 84-74-2	Price
n-hexane	1 000 µg/ml	1 ml	ref. P877910	CA
methanol	100 µg/ml	1 ml	ref. P877920	BH
PHTHALIC ACID, BIS-ETHYL ESTER			CAS 84-66-2	Price
n-hexane	1 000 µg/ml	1 ml	ref. P877930	CA
methanol	100 µg/ml	1 ml	ref. P877940	BH
PHTHALIC ACID, BIS-ETHYL ESTER D4			CAS 93952-12-6	Price
acetone	100 µg/ml	1 ml	ref. P820350	BCA
n-nonane	100 µg/ml	1 ml	ref. P820360	BCA
PHTHALIC ACID, BIS-2-ETHYLHEXYLESTER			CAS 117-81-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P877990	BG
dichloromethane	10 µg/ml	10 ml	ref. P807290	BG
	100 µg/ml	1 ml	ref. P807300	BG
methanol	100 µg/ml	1 ml	ref. P878000	BG
	1 000 µg/ml	1 ml	ref. P878010	BK
PHTHALIC ACID, BIS-ISO-BUTYL ESTER			CAS 84-69-5	Price
acetone	2 500 µg/ml	5 ml	ref. P877950	CF
cyclohexane	10 µg/ml	10 ml	ref. P877960	BG
PHTHALIC ACID, BIS-METYL ESTER			CAS 131-11-3	Price
methanol	100 µg/ml	1 ml	ref. F116481	BH
	1 000 µg/ml	1 ml	ref. P877970	CA
PHTHALIC ACID, BIS-1-OCTYL ESTER			CAS 117-84-0	Price
methanol	10 µg/ml	10 ml	ref. P822240	BG
	100 µg/ml	1 ml	ref. P822250	BG
	5 000 µg/ml	1 ml	ref. P877890	CE

PHTHALIC ACID, BIS-1-OCTYL ESTER-D4			id : 21330	Price
acetone	100 µg/ml	1 ml	ref. P822280	BCG
nonane	100 µg/ml	1 ml	ref. P822290	BCG
PHYTANE			CAS 638-36-8	Price
iso-octane	10 µg/ml	10 ml	ref. P848370	BG
	100 µg/ml	1 ml	ref. P848380	BG
PICENE			CAS 213-46-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P848410	BG
	100 µg/ml	1 ml	ref. P848420	BG
cyclohexane	10 µg/ml	10 ml	ref. P848430	BG
	100 µg/ml	1 ml	ref. P848440	BG
PICLORAM			CAS 1918-02-1	Price
acetone	10 µg/ml	10 ml	ref. P848470	BG
	100 µg/ml	1 ml	ref. P848480	BG
acetonitrile	10 µg/ml	10 ml	ref. P848490	BG
	100 µg/ml	1 ml	ref. P848500	BG
	1 000 µg/ml	1 ml	ref. F116051	BK
methanol	10 µg/ml	10 ml	ref. P848510	BG
	100 µg/ml	1 ml	ref. P848520	BG
PICOXYSTROBIN			CAS 117428-22-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P878030	BG
PIPERONYL BUTOXIDE			CAS 51-03-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P848550	BG
	100 µg/ml	1 ml	ref. P848560	BG
cyclohexane	10 µg/ml	10 ml	ref. P848570	BG
	100 µg/ml	1 ml	ref. P848580	BG
methanol	10 µg/ml	10 ml	ref. P848590	BG
	100 µg/ml	1 ml	ref. P848600	BG
PIRIMICARB			CAS 23103-98-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P848630	BG
	100 µg/ml	1 ml	ref. P848640	BG
cyclohexane	10 µg/ml	10 ml	ref. P848650	BG
	100 µg/ml	1 ml	ref. P848660	BG
PIRIMIPHOS-ETHYL			CAS 23505-41-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P848720	BG
	100 µg/ml	1 ml	ref. P848730	BG
methanol	10 µg/ml	10 ml	ref. P848740	BG
	100 µg/ml	1 ml	ref. P848750	BG
PIRIMIPHOS-METHYL			CAS 29232-93-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P848810	BG
	100 µg/ml	1 ml	ref. P848820	BG
PRETILACHLOR			CAS 51218-49-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P878070	BG
methanol	100 µg/ml	1 ml	ref. P878080	BG
PROCHLORAZ			CAS 67747-09-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P848890	BG
	100 µg/ml	1 ml	ref. P848900	BG
cyclohexane	10 µg/ml	10 ml	ref. P848910	BG
	100 µg/ml	1 ml	ref. P848920	BG
PROCYMIDONE			CAS 32809-16-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P848950	BG
	100 µg/ml	1 ml	ref. P848960	BG
iso-octane	10 µg/ml	10 ml	ref. P848970	BG
	100 µg/ml	1 ml	ref. P848980	BG
methanol	10 µg/ml	10 ml	ref. P848990	BG
	100 µg/ml	1 ml	ref. P849000	BG
PROFENOFOS			CAS 41198-08-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P849070	BG
	100 µg/ml	1 ml	ref. P849080	BG
PROFLURALIN			CAS 26399-36-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P880380	BH
	100 µg/ml	1 ml	ref. P880370	BH
	1 000 µg/ml	1 ml	ref. P880360	CA
PROMECARB			CAS 2631-37-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P878100	BG
	100 µg/ml	1 ml	ref. P878090	BG
PROMETON			CAS 1610-18-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P849110	BG
	100 µg/ml	1 ml	ref. P849120	BG
cyclohexane	10 µg/ml	10 ml	ref. P849130	BG
	100 µg/ml	1 ml	ref. P849140	BG

PROMETRYN			CAS 7287-19-6	Price
acetone	10 µg/ml	10 ml	ref. P849170	BG
	100 µg/ml	1 ml	ref. P849180	BG
acetonitrile	10 µg/ml	10 ml	ref. P849190	BG
	100 µg/ml	1 ml	ref. P849200	BG
cyclohexane	10 µg/ml	10 ml	ref. P849210	BG
	100 µg/ml	1 ml	ref. P849220	BG
PROPACHLOR			CAS 1918-16-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P849280	BG
	100 µg/ml	1 ml	ref. P849290	BG
methanol	10 µg/ml	10 ml	ref. P849300	BG
	100 µg/ml	1 ml	ref. P849310	BG
PROPAMOCARB FREE BASE			CAS 24579-73-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P878110	BG
methanol	100 µg/ml	1 ml	ref. P878120	BG
PROPAMOCARB HYDROCHLORIDE			CAS 25606-41-1	Price
methanol	10 µg/ml	10 ml	ref. P849340	BG
	100 µg/ml	1 ml	ref. P849350	BG
PROPANIL			CAS 709-98-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P849410	BG
	100 µg/ml	1 ml	ref. P849420	BG
cyclohexane	10 µg/ml	10 ml	ref. P849430	BG
	100 µg/ml	1 ml	ref. P849440	BG
PROPAQUIZAFOP			CAS 111479-05-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P878130	BG
PROPARGITE			CAS 2312-35-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P849500	BG
	100 µg/ml	1 ml	ref. P849510	BG
cyclohexane	10 µg/ml	10 ml	ref. P849520	BG
	100 µg/ml	1 ml	ref. P849530	BG
iso-octane	10 µg/ml	10 ml	ref. P849540	BG
	100 µg/ml	1 ml	ref. P849550	BG
methanol	10 µg/ml	10 ml	ref. P849560	BG
	100 µg/ml	1 ml	ref. P849570	BG
PROPAZINE			CAS 139-40-2	Price
acetone	10 µg/ml	10 ml	ref. P849600	BG
	100 µg/ml	1 ml	ref. P849610	BG
acetonitrile	10 µg/ml	10 ml	ref. P849620	BG
	100 µg/ml	1 ml	ref. P849630	BG
cyclohexane	10 µg/ml	10 ml	ref. P849640	BG
	100 µg/ml	1 ml	ref. P849650	BG
PROPETAMPHOS			CAS 31218-83-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P849750	BH
	100 µg/ml	1 ml	ref. P849760	BH
cyclohexane	10 µg/ml	10 ml	ref. P849770	BH
	100 µg/ml	1 ml	ref. P849780	BH
PROPHAM			CAS 122-42-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P878160	BG
	100 µg/ml	1 ml	ref. P878150	BG
	1 000 µg/ml	5 ml	ref. P878140	CE
cyclohexane	10 µg/ml	10 ml	ref. P878170	BG
PROPICONAZOLE			CAS 60207-90-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P849810	BG
	100 µg/ml	1 ml	ref. P849820	BG
cyclohexane	10 µg/ml	10 ml	ref. P849830	BG
	100 µg/ml	1 ml	ref. P849840	BG
methanol	10 µg/ml	10 ml	ref. P849850	BG
	100 µg/ml	1 ml	ref. P849860	BG
PROPIONALDEHYDE			CAS 123-38-6	Price
acetonitrile	1 000 µg/ml	5 ml	ref. P878180	CF
methanol	1 000 µg/ml	1 ml	ref. P878190	BH
PROPIONALDEHYDE-DNPH			CAS 725-00-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P849890	BG
	100 µg/ml	1 ml	ref. P849900	BG
PROPIONITRILE			CAS 107-12-0	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P849930	BG
	100 µg/ml	1 ml	ref. P849940	BG
PROPOXUR			CAS 114-26-1	Price
acetonitrile	100 µg/ml	5 ml	ref. P878200	BH
cyclohexane	10 µg/ml	10 ml	ref. P878210	BH

PROPYLAMINE			CAS 107-10-8	Price
methanol	1 000 µg/ml	1 ml	ref. P878230	BH
n-PROPYLBENZENE			CAS 103-65-1	Price
methanol	100 µg/ml	1 ml	ref. P877310	BH
	1 000 µg/ml	1 ml	ref. P877300	CA
PROPYZAMIDE			CAS 23950-58-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P850010	BG
	100 µg/ml	1 ml	ref. P850020	BG
cyclohexane	10 µg/ml	10 ml	ref. P850030	BG
	100 µg/ml	1 ml	ref. P850040	BG
methanol	10 µg/ml	10 ml	ref. P850050	BG
	100 µg/ml	1 ml	ref. P850060	BG
PROSULFOCARB			CAS 52888-80-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P878260	BG
cyclohexane	10 µg/ml	10 ml	ref. P878250	BG
	100 µg/ml	1 ml	ref. P878240	BG
PROSULFURON			CAS 94125-34-5	Price
acetonitrile	100 µg/ml	1 ml	ref. P878270	BG
PROTHIOPHOS			CAS 34643-46-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P850130	BG
	100 µg/ml	1 ml	ref. P850140	BG
cyclohexane	10 µg/ml	10 ml	ref. P850150	BG
	100 µg/ml	1 ml	ref. P850160	BG
PYMETROZINE			CAS 123312-89-0	Price
acetonitrile	100 µg/ml	1 ml	ref. P878310	BG
methanol	10 µg/ml	10 ml	ref. P878320	BG
PYRACLOSTROBIN			CAS 175013-18-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P878330	BG
	100 µg/ml	1 ml	ref. P878340	BG
PYRAFLUFEN-ETHYL			CAS 129630-19-9	Price
cyclohexane	10 µg/ml	10 ml	ref. P878350	BJ
PYRAZOPHOS			CAS 13457-18-6	Price
acetone	10 µg/ml	10 ml	ref. P850190	BG
	100 µg/ml	1 ml	ref. P850200	BG
PYRENE			CAS 129-00-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P850230	BG
	100 µg/ml	1 ml	ref. P850240	BG
cyclohexane	10 µg/ml	10 ml	ref. P850250	BG
	100 µg/ml	1 ml	ref. P850260	BG
dichloromethane	10 µg/ml	10 ml	ref. P850290	BG
	100 µg/ml	1 ml	ref. P850300	BG
methanol	10 µg/ml	10 ml	ref. P850270	BG
	100 µg/ml	1 ml	ref. P850280	BG
PYRENE-D10			CAS 1718-52-1	Price
acetone	1 000 µg/ml	1 ml	ref. P850330	CH
	500 µg/ml	1 ml	ref. F114091	CE
acetonitrile	100 µg/ml	1 ml	ref. P850340	CC
cyclohexane	100 µg/ml	1 ml	ref. P850350	CC
PYRETHRIN			CAS 8003-34-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P850380	BH
	100 µg/ml	1 ml	ref. P850390	BH
cyclohexane	10 µg/ml	10 ml	ref. P850400	BH
	100 µg/ml	1 ml	ref. P850410	BH
methanol	100 µg/ml	1 ml	ref. P878360	BH
PYRIDABEN			CAS 96489-71-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P878400	BG
	100 µg/ml	1 ml	ref. P878390	BG
cyclohexane	10 µg/ml	10 ml	ref. P878410	BG
PYRIDAPHENTHION			CAS 119-12-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P850440	BH
	100 µg/ml	1 ml	ref. P850450	BH
cyclohexane	10 µg/ml	10 ml	ref. P850460	BH
	100 µg/ml	1 ml	ref. P850470	BH
methanol	10 µg/ml	10 ml	ref. P850480	BH
	100 µg/ml	1 ml	ref. P850490	BH
PYRIDATE			CAS 55512-33-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P850520	BG
	100 µg/ml	1 ml	ref. P850530	BG
cyclohexane	10 µg/ml	10 ml	ref. P850540	BG
	100 µg/ml	1 ml	ref. P850550	BG

PYRIDINE			CAS 110-86-1	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P850580	BG
	100 µg/ml	1 ml	ref. P850590	BG
PYRIDINE-D5			CAS 7291-22-7	Price
dichloromethane	100 µg/ml	1 ml	ref. P850620	CD
	2 000 µg/ml	1 ml	ref. F062555	CJ
methyl- <i>tert</i> -butylether	100 µg/ml	1 ml	ref. P850630	CD
	2 000 µg/ml	1 ml	ref. F062565	CJ
PYRIFENOX			CAS 88283-41-4	Price
iso-octane	10 µg/ml	10 ml	ref. P850660	BH
	100 µg/ml	1 ml	ref. P850670	BH
methanol	10 µg/ml	10 ml	ref. P850680	BH
	100 µg/ml	1 ml	ref. P850690	BH
PYRIMETHANIL			CAS 53112-28-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P850720	BG
	100 µg/ml	1 ml	ref. P850730	BG
methanol	10 µg/ml	10 ml	ref. P850740	BG
	100 µg/ml	1 ml	ref. P850750	BG
PYRIPROXYFEN			CAS 95737-68-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P878430	BG
	100 µg/ml	1 ml	ref. P878420	BG
cyclohexane	10 µg/ml	10 ml	ref. P878440	BG
QUINALPHOS			CAS 13593-03-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P850780	BG
	100 µg/ml	1 ml	ref. P850790	BG
methanol	10 µg/ml	10 ml	ref. P850800	BG
	100 µg/ml	1 ml	ref. P850810	BG
QUINCLORAC			CAS 84087-01-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P880410	BG
	100 µg/ml	1 ml	ref. P880400	BG
methanol	100 µg/ml	1 ml	ref. P880390	BG
QUINMERAC			CAS 90717-03-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P880430	BG
	100 µg/ml	1 ml	ref. P880420	BG
QUINOXYFEN			CAS 124495-18-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P850840	BH
	100 µg/ml	1 ml	ref. P850850	BH
cyclohexane	10 µg/ml	10 ml	ref. P850860	BH
	100 µg/ml	1 ml	ref. P850870	BH
methanol	10 µg/ml	10 ml	ref. P850880	BH
	100 µg/ml	1 ml	ref. P850890	BH
QUINTOZENE			CAS 82-68-8	Price
acetone	5 000 µg/ml	1 ml	ref. F117251	CF
acetonitrile	10 µg/ml	10 ml	ref. P850920	BH
	100 µg/ml	1 ml	ref. P850930	BH
cyclohexane	10 µg/ml	10 ml	ref. P850940	BH
	100 µg/ml	1 ml	ref. P850950	BH
ethyl acetate	1 000 µg/ml	1 ml	ref. F113611	CA
iso-octane	1 000 µg/ml	1 ml	ref. F113601	CA
methanol	100 µg/ml	1 ml	ref. F113591	BH
methyl- <i>tert</i> -butylether	1 000 µg/ml	1 ml	ref. F120131	CA
QUIZALOFOP-ETHYL			CAS 76578-14-8	Price
acetonitrile	100 µg/ml	1 ml	ref. P878460	BG
iso-octane	10 µg/ml	10 ml	ref. P878450	BG
RESMETHRIN			CAS 10453-86-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P851020	BG
	100 µg/ml	1 ml	ref. P851030	BG
cyclohexane	10 µg/ml	10 ml	ref. P851040	BG
	100 µg/ml	1 ml	ref. P851050	BG
iso-octane	10 µg/ml	10 ml	ref. P851080	BG
	100 µg/ml	1 ml	ref. P851090	BG
methanol	10 µg/ml	10 ml	ref. P851060	BG
	100 µg/ml	1 ml	ref. P851070	BG
RESORCINOL			CAS 108-46-3	Price
methanol	10 µg/ml	10 ml	ref. P878480	BG
	100 µg/ml	5 ml	ref. P878470	BG
RIMSULFURON			CAS 122931-48-0	Price
acetone	10 µg/ml	10 ml	ref. P851120	BH
	100 µg/ml	1 ml	ref. P851130	BH
ROTENONE			CAS 83-79-4	Price
methanol	1 000 µg/ml	1 ml	ref. F115971	BH

SALICYLIC ACID			CAS 69-72-7	Price
isopropanol (2-propanol)	2 000 µg/ml	1 ml	ref. F116061	CA
	2 000 µg/ml	5 ml	ref. P878490	DC
water	100 µg/ml	1 ml	ref. P878500	BG
SEBUTHYLAZINE			CAS 7286-69-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P851160	BG
	100 µg/ml	1 ml	ref. P851170	BG
cyclohexane	10 µg/ml	10 ml	ref. P851180	BG
	100 µg/ml	1 ml	ref. P851190	BG
SEBUTHYLAZINE-DESETHYL			CAS 37019-18-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P873810	BG
	100 µg/ml	1 ml	ref. P851260	BG
SECBUMETON			CAS 26259-45-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P851290	BG
	100 µg/ml	1 ml	ref. P851300	BG
cyclohexane	10 µg/ml	10 ml	ref. P851310	BG
	100 µg/ml	1 ml	ref. P851320	BG
iso-octane	10 µg/ml	10 ml	ref. P851330	BG
	100 µg/ml	1 ml	ref. P851340	BG
SIDURON			CAS 1982-49-6	Price
acetonitrile	1 000 µg/ml	1 ml	ref. P878510	CB
methanol	100 µg/ml	1 ml	ref. P878520	BJ
SIMAZINE			CAS 122-34-9	Price
acetone	10 µg/ml	10 ml	ref. P851400	BG
	100 µg/ml	1 ml	ref. P851410	BG
	1 000 µg/ml	1 ml	ref. F113411	BK
acetonitrile	10 µg/ml	10 ml	ref. P851420	BG
	100 µg/ml	1 ml	ref. P851430	BG
SIMAZINE-D5			CAS 220621-41-0	Price
acetone	100 µg/ml	1 ml	ref. P878530	EJ
SIMAZINE-2-HYDROXY			CAS 2599-11-3	Price
methanol	10 µg/ml	10 ml	ref. P878540	BJ
SIMETRYN			CAS 1014-70-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P878560	BG
cyclohexane	100 µg/ml	1 ml	ref. P878550	BG
SPIROXAMINE			CAS 118134-30-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P851560	BG
	100 µg/ml	1 ml	ref. P851570	BG
cyclohexane	10 µg/ml	10 ml	ref. P851580	BG
	100 µg/ml	1 ml	ref. P851590	BG
methanol	10 µg/ml	10 ml	ref. P851600	BG
	100 µg/ml	1 ml	ref. P851610	BG
STEARYL STEARATE			CAS 2778-96-3	Price
cyclohexane	2 000 µg/ml	10 ml	ref. P878570	JD
n-hexane	2 000 µg/ml	10 ml	ref. F109761	JD
STYRENE			CAS 100-42-5	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P851700	BG
	100 µg/ml	1 ml	ref. P851710	BG
SULCOTRIONE			CAS 99105-77-8	Price
acetone	10 µg/ml	10 ml	ref. P851740	BG
	100 µg/ml	1 ml	ref. P851750	BG
SULFOSULFURON			CAS 141776-32-1	Price
acetonitrile	10 µg/ml	1 ml	ref. P878580	BG
SULFOTEP			CAS 3689-24-5	Price
iso-octane	10 µg/ml	10 ml	ref. P878590	BG
SULPROFOS			CAS 35400-43-2	Price
cyclohexane	100 µg/ml	1 ml	ref. P878600	BG
SWEP			CAS 1918-18-9	Price
iso-octane	10 µg/ml	10 ml	ref. P880450	BG
methanol	100 µg/ml	1 ml	ref. P880440	BG
2,4,5-T			CAS 93-76-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P851900	BG
	100 µg/ml	1 ml	ref. P851910	BG
methanol	10 µg/ml	10 ml	ref. P851920	BG
	100 µg/ml	1 ml	ref. P851930	BG
2,4,5-T METHYL ESTER			CAS 1928-37-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P851960	BG
	100 µg/ml	1 ml	ref. P851970	BG
methanol	10 µg/ml	10 ml	ref. P851980	BG
	100 µg/ml	1 ml	ref. P851990	BG

TAU-FLUVALINATE			CAS 102851-06-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P852020	BG
	100 µg/ml	1 ml	ref. P852030	BG
cyclohexane	10 µg/ml	10 ml	ref. P852040	BG
	100 µg/ml	1 ml	ref. P852050	BG
methanol	10 µg/ml	10 ml	ref. P852060	BG
	100 µg/ml	1 ml	ref. P852070	BG

TBTO			CAS 56-35-9	Price
methanol	100 µg/ml	1 ml	ref. P878610	BH
	1 000 µg/ml	1 ml	ref. P878620	CA

TEBUCONAZOLE			CAS 107534-96-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P852100	BG
	100 µg/ml	1 ml	ref. P852110	BG
iso-octane	10 µg/ml	10 ml	ref. P852120	BG
	100 µg/ml	1 ml	ref. P852130	BG

TEBUFENOZIDE			CAS 112410-23-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P852160	BH
	100 µg/ml	1 ml	ref. P852170	BH

TEBUFENPYRAD			CAS 119168-77-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P878650	BG
cyclohexane	10 µg/ml	10 ml	ref. P878640	BG
	100 µg/ml	1 ml	ref. P878630	BG

TEBUTAM			CAS 35256-85-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P852200	BG
	100 µg/ml	1 ml	ref. P852210	BG

TEBUTHIURON			CAS 34014-18-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P878660	BG

TECNAZENE			CAS 117-18-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P878690	BG
iso-octane	100 µg/ml	1 ml	ref. P878680	BG

TEFLUBENZURON			CAS 83121-18-0	Price
acetone	10 µg/ml	10 ml	ref. P852240	BG
	100 µg/ml	1 ml	ref. P852250	BG
cyclohexane	10 µg/ml	10 ml	ref. P852260	BG
	100 µg/ml	1 ml	ref. P852270	BG

TEFLUTHRIN			CAS 79538-32-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P852300	BH
	100 µg/ml	1 ml	ref. P852310	BH
cyclohexane	10 µg/ml	10 ml	ref. P852320	BH
	100 µg/ml	1 ml	ref. P852330	BH

TEMEPHOS			CAS 3383-96-8	Price
acetone	10 µg/ml	10 ml	ref. P852360	BG
	100 µg/ml	1 ml	ref. P852370	BG
cyclohexane	10 µg/ml	10 ml	ref. P852380	BG
	100 µg/ml	1 ml	ref. P852390	BG

TERBACIL			CAS 5902-51-2	Price
iso-octane	10 µg/ml	10 ml	ref. P878700	BG

TERBUFOS			CAS 13071-79-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P852530	BG
	100 µg/ml	1 ml	ref. P852540	BG
cyclohexane	10 µg/ml	10 ml	ref. P852550	BG
	100 µg/ml	1 ml	ref. P852560	BG

TERBUFOS SULFONE			CAS 56070-16-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P852590	BH
	100 µg/ml	1 ml	ref. P852600	BH

TERBUFOS SULFOXIDE			CAS 10548-10-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P852630	BH
	100 µg/ml	1 ml	ref. P852640	BH

TERBUMETON			CAS 33693-04-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P852670	BG
	100 µg/ml	1 ml	ref. P852680	BG
methanol	10 µg/ml	10 ml	ref. P852690	BG
	100 µg/ml	1 ml	ref. P852700	BG

TERBUMETON-DESETHYL			CAS 30125-64-5	Price
methanol	10 µg/ml	10 ml	ref. P852730	BG
	100 µg/ml	1 ml	ref. P852740	BG

TERBUTHYLAZINE			CAS 5915-41-3	Price
acetone	10 µg/ml	10 ml	ref. P852770	BG
	100 µg/ml	1 ml	ref. P852780	BG
acetonitrile	10 µg/ml	10 ml	ref. P852790	BG
	100 µg/ml	1 ml	ref. P852800	BG
cyclohexane	10 µg/ml	10 ml	ref. P852810	BG
	100 µg/ml	1 ml	ref. P852820	BG

TERBUTHYLAZINE-DESETHYL			CAS 30125-63-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P852890	BG
	100 µg/ml	1 ml	ref. P852900	BG

TERBUTHYLAZINE-2-HYDROXY			CAS 66753-07-9	Price
acetone	10 µg/ml	10 ml	ref. P852930	BG
	100 µg/ml	1 ml	ref. P852940	BG

TERBUTRYN			CAS 886-50-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P852970	BG
	100 µg/ml	1 ml	ref. P852980	BG
cyclohexane	10 µg/ml	10 ml	ref. P852990	BG
	100 µg/ml	1 ml	ref. P853000	BG

o-TERPHENYL			CAS 84-15-1	Price
acetone	10 µg/ml	10 ml	ref. P853060	BG
	100 µg/ml	1 ml	ref. P853070	BG
dichloromethane	10 µg/ml	10 ml	ref. P853100	BG
	100 µg/ml	1 ml	ref. P853110	BG
n-hexane	10 µg/ml	10 ml	ref. P853080	BG
	100 µg/ml	1 ml	ref. P853090	BG

P-TERPHENYL			CAS 92-94-4	Price
dichloromethane	10 µg/ml	10 ml	ref. P853140	BG
	100 µg/ml	1 ml	ref. P853150	BG

P-TERPHENYL-D14			CAS 1718-51-0	Price
dichloromethane	500 µg/ml	1 ml	ref. F114101	BK
	2 000 µg/ml	1 ml	ref. P887010	FA
	2 000 µg/ml	5 ml	ref. P887000	CCJ

a-TERPINEOL			CAS 98-55-5	Price
dichloromethane	2 000 µg/ml	1 ml	ref. P882630	BH

1,2,3,4-TETRACHLOROBENZENE			CAS 634-66-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P872070	BG
methanol	100 µg/ml	1 ml	ref. P872080	BG

1,2,3,5-TETRACHLOROBENZENE			CAS 634-90-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P872150	BG
methanol	100 µg/ml	1 ml	ref. P872120	BG

1,2,4,5-TETRACHLOROBENZENE			CAS 95-94-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P853220	BG
	100 µg/ml	1 ml	ref. P853230	BG
cyclohexane	10 µg/ml	10 ml	ref. P853240	BG
	100 µg/ml	1 ml	ref. P853250	BG
methanol	10 µg/ml	10 ml	ref. P853260	BG
	100 µg/ml	1 ml	ref. P853270	BG

1,1,1,2-TETRACHLOROETHANE			CAS 630-20-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P853300	BG
	100 µg/ml	1 ml	ref. P853310	BG

1,1,2,2-TETRACHLOROETHANE			CAS 79-34-5	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P853340	BG
	100 µg/ml	1 ml	ref. P853350	BG

TETRACHLOROETHENE			CAS 127-18-4	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P853380	BG
	100 µg/ml	1 ml	ref. P853390	BG

TETRACHLOROMETHANE			CAS 56-23-5	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P809960	BG
	100 µg/ml	1 ml	ref. P809970	BG

2,3,4,5-TETRACHLOROPHENOL			CAS 4901-51-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P872550	CC
methanol	10 µg/ml	10 ml	ref. P872510	CC
	100 µg/ml	1 ml	ref. P872530	CC

2,3,4,6-TETRACHLOROPHENOL			CAS 58-90-2	Price
methanol	10 µg/ml	10 ml	ref. P853500	BH
	100 µg/ml	1 ml	ref. P853510	BH

2,3,5,6-TETRACHLOROPHENOL			CAS 935-95-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P872570	BK
	100 µg/ml	1 ml	ref. P872560	BK
isopropanol (2-propanol)	100 µg/ml	5 ml	ref. P872580	CK

2,4,5,6-TETRACHLORO <i>m</i>-XYLENE				CAS 877-09-8	Price
acetone	10 µg/ml	10 ml	ref. P853420	BG	
	100 µg/ml	1 ml	ref. P853430	BG	
	1 000 µg/ml	1 ml	ref. F117321	BK	
cyclohexane	10 µg/ml	10 ml	ref. P853440	BG	
	100 µg/ml	1 ml	ref. P853450	BG	
methanol	10 µg/ml	10 ml	ref. P853460	BG	
	100 µg/ml	1 ml	ref. P853470	BG	
TETRACHLORVINPHOS				CAS 22248-79-9	Price
iso-octane	10 µg/ml	10 ml	ref. P878730	BG	
TETRACONAZOLE				CAS 112281-77-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P853540	BH	
	100 µg/ml	1 ml	ref. P853550	BH	
iso-octane	10 µg/ml	10 ml	ref. P853580	BH	
	100 µg/ml	1 ml	ref. P853590	BH	
methanol	10 µg/ml	10 ml	ref. P853560	BH	
	100 µg/ml	1 ml	ref. P853570	BH	
<i>n</i>-TETRACOSANE				CAS 646-31-1	Price
dichloromethane	100 µg/ml	5 ml	ref. P877320	BG	
TETRACOSANOIC ACID				CAS 557-59-5	Price
cyclohexane	10 µg/ml	10 ml	ref. P878740	BG	
TETRACOSANOIC ACID-METHYL ESTER				CAS 2442-49-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P878750	BG	
<i>n</i>-TETRADECANE				CAS 629-59-4	Price
<i>n</i> -hexane	10 µg/ml	10 ml	ref. P852490	BJ	
	100 µg/ml	1 ml	ref. P852500	BJ	
TETRADECANOIC ACID-METHYL ESTER				CAS 124-10-7	Price
<i>n</i> -hexane	10 000 µg/ml	1 ml	ref. P878770	BJ	
methyl- <i>tert</i> -butylether	5 000 µg/ml	5 ml	ref. P878760	CF	
TETRADIFON				CAS 116-29-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P853620	BG	
	100 µg/ml	1 ml	ref. P853630	BG	
cyclohexane	10 µg/ml	10 ml	ref. P853640	BG	
	100 µg/ml	1 ml	ref. P853650	BG	
TETRAHYDROFURAN				CAS 109-99-9	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P853680	BG	
	100 µg/ml	1 ml	ref. P853690	BG	
1,2,3,4-TETRAHYDRONAPHTHALENE				CAS 119-64-2	Price
methanol	10 µg/ml	10 ml	ref. P872110	BG	
toluene	100 µg/ml	1 ml	ref. P872100	BG	
TETRAMETHRIN				CAS 7696-12-0	Price
acetonitrile	10 µg/ml	10 ml	ref. P853720	BG	
	100 µg/ml	1 ml	ref. P853730	BG	
cyclohexane	10 µg/ml	10 ml	ref. P853740	BG	
	100 µg/ml	1 ml	ref. P853750	BG	
1,2,3,5-TETRAMETHYLBENZENE				CAS 527-53-7	Price
methanol	10 µg/ml	10 ml	ref. P872170	BG	
1,2,4,5-TETRAMETHYLBENZENE				CAS 95-93-2	Price
cyclohexane	10 µg/ml	10 ml	ref. P872180	BG	
TETRAPROPYL TIN				CAS 2176-98-9	Price
dichloromethane	10 µg/ml	10 ml	ref. P853780	BG	
	100 µg/ml	1 ml	ref. P853790	BG	
TETRASUL				CAS 2227-13-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P878780	BH	
methanol	100 µg/ml	1 ml	ref. P878790	BH	
THIABENDAZOLE				CAS 148-79-8	Price
acetone	10 µg/ml	10 ml	ref. P854670	BG	
	100 µg/ml	1 ml	ref. P854680	BG	
acetonitrile	10 µg/ml	10 ml	ref. P853900	BG	
	100 µg/ml	1 ml	ref. P853910	BG	
methanol	10 µg/ml	10 ml	ref. P853920	BG	
	100 µg/ml	1 ml	ref. P853930	BG	
THIACHLOPRID				CAS 111988-49-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P880470	BG	
	100 µg/ml	1 ml	ref. P880460	BG	
THIAMETHOXAM				CAS 153719-23-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P880490	BG	
	100 µg/ml	1 ml	ref. P880480	BG	

THIAZAFLURON				CAS 25366-23-8	Price
acetone	10 µg/ml	10 ml	ref. P854000	BG	
	100 µg/ml	1 ml	ref. P854010	BG	
acetonitrile	10 µg/ml	10 ml	ref. P854020	BG	
	100 µg/ml	1 ml	ref. P854030	BG	
THIFENSULFURON-METHYL				CAS 79277-27-3	Price
methanol	100 µg/ml	1 ml	ref. P878800	BH	
THIOBENCARB				CAS 28249-77-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P878810	BG	
THIODICARB				CAS 59669-26-0	Price
acetone	10 µg/ml	10 ml	ref. P854060	BG	
	100 µg/ml	1 ml	ref. P854070	BG	
ethyl acetate	10 µg/ml	10 ml	ref. P854080	BG	
	100 µg/ml	1 ml	ref. P854090	BG	
THIOFANOX				CAS 39196-18-4	Price
cyclohexane	10 µg/ml	10 ml	ref. P878830	BG	
methanol	100 µg/ml	1 ml	ref. P878840	BG	
THIOFANOX-SULFOXIDE				CAS 39184-27-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P878850	BG	
methanol	100 µg/ml	1 ml	ref. P878860	BG	
THIOMETON				CAS 640-15-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P854120	BG	
	100 µg/ml	1 ml	ref. P854130	BG	
THIOPHANATE-ETHYL				CAS 23564-06-9	Price
acetonitrile	100 µg/ml	1 ml	ref. P878890	BH	
cyclohexane	10 µg/ml	10 ml	ref. P878870	BH	
methanol	100 µg/ml	1 ml	ref. P878880	BH	
THIOPHANATE-METHYL				CAS 23564-05-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P878900	BG	
methanol	100 µg/ml	1 ml	ref. P878910	BG	
THIRAM				CAS 137-26-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P854160	BG	
	100 µg/ml	1 ml	ref. P854170	BG	
cyclohexane	10 µg/ml	10 ml	ref. P854180	BG	
	100 µg/ml	1 ml	ref. P854190	BG	
TOLCLOFOS-METHYL				CAS 57018-04-9	Price
acetonitrile	100 µg/ml	1 ml	ref. P878920	BG	
cyclohexane	10 µg/ml	10 ml	ref. P878930	BG	
<i>p</i>-TOLUALDEHYDE				CAS 104-87-0	Price
acetonitrile	1 000 µg/ml	1 ml	ref. P878290	BH	
	1 000 µg/ml	5 ml	ref. P878280	CE	
<i>m</i>-TOLUALDEHYDE-2,4-DNPH				CAS 2880-05-9	Price
acetonitrile	10 µg/ml	10 ml	ref. P854250	BG	
	100 µg/ml	1 ml	ref. P854260	BG	
<i>o</i>-TOLUALDEHYDE-2,4-DNPH				CAS 1773-44-0	Price
acetonitrile	100 µg/ml	1 ml	ref. P877450	BG	
<i>p</i>-TOLUALDEHYDE-2,4-DNPH				CAS 2571-00-8	Price
acetonitrile	100 µg/ml	1 ml	ref. P878300	BG	
TOLUENE				CAS 108-88-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P854290	BG	
	100 µg/ml	1 ml	ref. P854300	BG	
TOLUENE-D8				CAS 2037-26-5	Price
acetone	50 µg/ml	1 ml	ref. F129011	BK	
methanol purge & trap	10 µg/ml	10 ml	ref. P854330	BK	
	100 µg/ml	1 ml	ref. P854340	BK	
<i>o</i>-TOLUIDINE				CAS 95-53-4	Price
acetonitrile	100 µg/ml	1 ml	ref. P877490	BG	
	1 000 µg/ml	1 ml	ref. P877480	BK	
cyclohexane	10 µg/ml	10 ml	ref. P877460	BG	
	100 µg/ml	1 ml	ref. P877470	BG	
TOLYLFLUANID				CAS 731-27-1	Price
cyclohexane	10 µg/ml	10 ml	ref. P854370	BG	
	100 µg/ml	1 ml	ref. P854380	BG	

TOXAPHENE			CAS 8001-35-2	Price
acetone	2 500 µg/ml	1 ml	ref. F113495	CA
dichloromethane	100 µg/ml	1 ml	ref. F119811	BG
n-hexane	10 µg/ml	10 ml	ref. P854410	BG
	100 µg/ml	1 ml	ref. F116731	BG
iso-octane	10 µg/ml	10 ml	ref. P854430	BG
	100 µg/ml	1 ml	ref. F115061	BG
	1 000 µg/ml	1 ml	ref. F115671	BK
methanol	10 µg/ml	10 ml	ref. P854450	BG
	100 µg/ml	1 ml	ref. F120031	BG
methyl- <i>tert</i> -butylether	2 500 µg/ml	1 ml	ref. F113491	CA

TRALOMETHRIN			CAS 66841-25-6	Price
iso-octane	10 µg/ml	10 ml	ref. P878940	BH
	100 µg/ml	1 ml	ref. P878950	BH

TRANSFLUTHRIN			CAS 118712-89-3	Price
cyclohexane	10 µg/ml	10 ml	ref. P854630	BG
	100 µg/ml	1 ml	ref. P854640	BG

<i>n</i> -TRIACONTANE			CAS 638-68-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P854750	BG
	100 µg/ml	1 ml	ref. P854760	BG
dichloromethane	10 µg/ml	10 ml	ref. P854770	BG
	100 µg/ml	1 ml	ref. P854780	BG

TRIADIMEFON			CAS 43121-43-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P854810	BG
	100 µg/ml	1 ml	ref. P854820	BG
cyclohexane	10 µg/ml	10 ml	ref. P854830	BG
	100 µg/ml	1 ml	ref. P854840	BG

TRIADIMENOL			CAS 55219-65-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P854870	BG
	100 µg/ml	1 ml	ref. P854880	BG
cyclohexane	10 µg/ml	10 ml	ref. P854890	BG
	100 µg/ml	1 ml	ref. P854900	BG
iso-octane	10 µg/ml	10 ml	ref. P854910	BG
	100 µg/ml	1 ml	ref. P854920	BG
methanol	10 µg/ml	10 ml	ref. P854930	BG
	100 µg/ml	1 ml	ref. P854940	BG

TRI-ALLATE			CAS 2303-17-5	Price
acetonitrile	10 µg/ml	10 ml	ref. P854970	BG
	100 µg/ml	1 ml	ref. P854980	BG
cyclohexane	10 µg/ml	10 ml	ref. P854990	BG
	100 µg/ml	1 ml	ref. P855000	BG

TRIASULFURON			CAS 82097-50-5	Price
acetone	10 µg/ml	10 ml	ref. P855030	BG
	100 µg/ml	1 ml	ref. P855040	BG

TRIAZOPHOS			CAS 24017-47-8	Price
cyclohexane	10 µg/ml	10 ml	ref. P855070	BG
	100 µg/ml	1 ml	ref. P855080	BG
methanol	10 µg/ml	10 ml	ref. P855090	BG
	100 µg/ml	1 ml	ref. P855100	BG

TRIBROMOACETIC ACID			CAS 75-96-7	Price
methyl- <i>tert</i> -butylether	100 µg/ml	1 ml	ref. P880510	BG
	1 000 µg/ml	1 ml	ref. P880500	BK

1,3,5-TRIBROMOBENZENE			CAS 626-39-1	Price
acetone	50 µg/ml	1 ml	ref. F117591	BG
n-hexane	1 000 µg/ml	1 ml	ref. F117581	BK

2,4,6-TRIBROMOPHENOL			CAS 118-79-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P855210	BG
	100 µg/ml	1 ml	ref. P855220	BG
dichloromethane	10 µg/ml	10 ml	ref. P855250	BG
	100 µg/ml	1 ml	ref. P855260	BG
ethanol	1.148 µg/ml	5 ml	ref. F109891	CC
	1 000 µg/ml	1 ml	ref. F131361	BK
methanol	10 µg/ml	10 ml	ref. P855230	BG
	100 µg/ml	1 ml	ref. P855240	BG
	2 000 µg/ml	1 ml	ref. F114921	CA

TRIBUFOS			CAS 78-48-8	Price
cyclohexane	100 µg/ml	1 ml	ref. P879720	BG
methyl- <i>tert</i> -butylether	1 000 µg/ml	1 ml	ref. F113421	BK

TRIBUTYL PHOSPHATE			CAS 126-73-8	Price
acetone	10 µg/ml	10 ml	ref. P855290	BG
	100 µg/ml	1 ml	ref. P855300	BG
cyclohexane	10 µg/ml	10 ml	ref. P855310	BG
	100 µg/ml	1 ml	ref. P855320	BG

TRIBUTYL TIN CHLORIDE TBTC			CAS 1461-22-9	Price
dichloromethane	10 µg/ml	10 ml	ref. P855370	BG
	100 µg/ml	1 ml	ref. P855380	BG
iso-octane	10 µg/ml	10 ml	ref. P855350	BG
	100 µg/ml	1 ml	ref. P855360	BG

TRICHLORFON			CAS 52-68-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P878960	BG
	100 µg/ml	1 ml	ref. P878970	BG

2,4,5-TRICHLOROANILINE			CAS 636-30-6	Price
cyclohexane	10 µg/ml	10 ml	ref. P872680	BG

2,4,6-TRICHLOROANISOLE			CAS 87-40-1	Price
iso-octane	10 µg/ml	10 ml	ref. P872730	BG

1,2,3-TRICHLOROBENZENE			CAS 87-61-6	Price
acetonitrile	10 µg/ml	10 ml	ref. P855490	BG
	100 µg/ml	1 ml	ref. P855500	BG
cyclohexane	10 µg/ml	10 ml	ref. P855510	BG
	100 µg/ml	1 ml	ref. P855520	BG
iso-octane	10 µg/ml	10 ml	ref. P855530	BG
	100 µg/ml	1 ml	ref. P855540	BG
methanol	10 µg/ml	10 ml	ref. P855550	BG
	100 µg/ml	1 ml	ref. P855560	BG

1,2,4-TRICHLOROBENZENE			CAS 120-82-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P855590	BG
	100 µg/ml	1 ml	ref. P855600	BG
n-hexane	10 µg/ml	10 ml	ref. P855610	BG
	100 µg/ml	1 ml	ref. P855620	BG
iso-octane	10 µg/ml	10 ml	ref. P855630	BG
	100 µg/ml	1 ml	ref. P855640	BG
methanol	10 µg/ml	10 ml	ref. P855650	BG
	100 µg/ml	1 ml	ref. P855660	BG
	100 µg/ml	5 ml	ref. P872210	BG

1,3,5-TRICHLOROBENZENE			CAS 108-70-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P855690	BG
	100 µg/ml	1 ml	ref. P855700	BG
cyclohexane	10 µg/ml	10 ml	ref. P855710	BG
	100 µg/ml	1 ml	ref. P855720	BG
methanol	10 µg/ml	10 ml	ref. P855730	BG
	100 µg/ml	1 ml	ref. P855740	BG

1,1,1-TRICHLOROETHANE			CAS 71-55-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P855800	BG
	100 µg/ml	1 ml	ref. P855810	BG

1,1,2-TRICHLOROETHANE			CAS 79-00-5	Price
methanol	10 µg/ml	10 ml	ref. P855860	BG
	100 µg/ml	1 ml	ref. P855870	BG

TRICHLOROETHENE			CAS 79-01-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P855940	BG
	100 µg/ml	1 ml	ref. P855950	BG

TRICHLORONAT			CAS 327-98-0	Price
cyclohexane	10 µg/ml	10 ml	ref. P856020	CA
	100 µg/ml	1 ml	ref. P856030	CA

2,3,5-TRICHLOROPHENOL			CAS 933-78-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P872600	BG
iso-octane	10 µg/ml	10 ml	ref. P872610	BG
isopropanol (2-propanol)	1 000 µg/ml	5 ml	ref. P872620	CK
methanol	100 µg/ml	1 ml	ref. P872590	BG

2,3,6-TRICHLOROPHENOL			CAS 933-75-5	Price
ethanol	1 000 µg/ml	1 ml	ref. F131351	BH
isopropanol (2-propanol)	1 000 µg/ml	5 ml	ref. P872650	CK
methanol	10 µg/ml	10 ml	ref. P872640	BH

2,4,5-TRICHLOROPHENOL			CAS 95-95-4	Price
acetone	10 µg/ml	10 ml	ref. P856060	BG
	100 µg/ml	1 ml	ref. P856070	BG
acetonitrile	10 µg/ml	10 ml	ref. P856080	BG
	100 µg/ml	1 ml	ref. P856090	BG
cyclohexane	10 µg/ml	10 ml	ref. P856100	BG
	100 µg/ml	1 ml	ref. P856110	BG
methanol	10 µg/ml	10 ml	ref. P856120	BG
	100 µg/ml	1 ml	ref. P856130	BG
2,4,6-TRICHLOROPHENOL			CAS 88-06-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P856160	BG
	100 µg/ml	1 ml	ref. P856170	BG
cyclohexane	10 µg/ml	10 ml	ref. P856180	BG
	100 µg/ml	1 ml	ref. P856190	BG
methanol	10 µg/ml	10 ml	ref. P856200	BG
	100 µg/ml	1 ml	ref. P856210	BG
3,4,5-TRICHLOROPHENOL			CAS 609-19-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P873270	CA
	100 µg/ml	1 ml	ref. P873260	CA
iso-octane	10 µg/ml	10 ml	ref. P873280	CA
2,4,5-TRICHLOROPHENOL ACETATE			CAS 5393-75-9	Price
iso-octane	10 µg/ml	10 ml	ref. P872690	CC
2,4,6-TRICHLOROPHENOL ACETATE			CAS 23399-90-8	Price
iso-octane	10 µg/ml	10 ml	ref. P872740	BJ
1,2,3-TRICHLOROPROPANE			CAS 96-18-4	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P856280	BG
	100 µg/ml	1 ml	ref. P856290	BG
methyl-tert-butylether	10 µg/ml	10 ml	ref. P856300	BG
	100 µg/ml	1 ml	ref. P856310	BG
α, 3,4-TRICHLOROTOLUENE			CAS 102-47-6	Price
dichloromethane	10 µg/ml	10 ml	ref. P856340	BG
	100 µg/ml	1 ml	ref. P856350	BG
1,1,1-TRICHLOROTRIFLUOROETHANE			CAS 354-58-5	Price
methanol	10 µg/ml	10 ml	ref. P856380	BG
	100 µg/ml	1 ml	ref. P856390	BG
1,1,2-TRICHLOROTRIFLUOROETHANE			CAS 76-13-1	Price
methanol	200 µg/ml	1 ml	ref. P872010	CC
	2 000 µg/ml	1 ml	ref. P872000	CK
TRICLOPYR			CAS 55335-06-3	Price
acetone	10 µg/ml	10 ml	ref. P856420	BG
	100 µg/ml	1 ml	ref. P856430	BG
TRICYCLAZOLE			CAS 41814-78-2	Price
methanol	100 µg/ml	1 ml	ref. P878990	BG
n-TRIDECANE			CAS 629-50-5	Price
methanol	10 µg/ml	10 ml	ref. P856460	BG
	100 µg/ml	1 ml	ref. P856470	BG
TRITIAZINE			CAS 1912-26-1	Price
acetonitrile	10 µg/ml	10 ml	ref. P887900	BG
methanol	100 µg/ml	1 ml	ref. P887890	BG
TRIFLOXYSTROBIN			CAS 141517-21-7	Price
acetonitrile	100 µg/ml	1 ml	ref. P879020	BG
cyclohexane	10 µg/ml	10 ml	ref. P879010	BG
TRIFLUMIZOLE			CAS 68694-11-1	Price
acetonitrile	100 µg/ml	1 ml	ref. P880520	BG
iso-octane	10 µg/ml	10 ml	ref. P880530	BG
TRIFLUMURON			CAS 64628-44-0	Price
acetone	10 µg/ml	10 ml	ref. P856530	BG
	100 µg/ml	1 ml	ref. P856540	BG
acetonitrile	10 µg/ml	10 ml	ref. P856550	BG
	100 µg/ml	1 ml	ref. P856560	BG
2',4',5'-TRIFLUOROACETOPHENONE			CAS 129322-83-4	Price
acetonitrile	2 000 µg/ml	1 ml	ref. F114811	CB
	10 000 µg/ml	1 ml	ref. F114801	CK

α,α,α-TRIFLUOROTOLUENE			CAS 98-08-8	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P856590	BG
	100 µg/ml	1 ml	ref. P856600	BG
	200 µg/ml	1 ml	ref. F113071	CF
	2 000 µg/ml	1 ml	ref. F113061	CF
TRIFLURALIN			CAS 1582-09-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P856630	BG
	100 µg/ml	1 ml	ref. P856640	BG
cyclohexane	10 µg/ml	10 ml	ref. P856650	BG
	100 µg/ml	1 ml	ref. P856660	BG
methanol	10 µg/ml	10 ml	ref. P856670	BG
	100 µg/ml	1 ml	ref. P856680	BG
TRIFLURALIN-D14			CAS 347841-79-6	Price
acetone	100 µg/ml	1 ml	ref. P856710	FJ
nonane	100 µg/ml	1 ml	ref. P856720	FJ
2,4,5-TRIMETHYLANILINE			CAS 137-17-7	Price
acetonitrile	10 µg/ml	10 ml	ref. P872710	CB
methanol	100 µg/ml	5 ml	ref. P872700	DC
1,2,3-TRIMETHYLBENZENE			CAS 526-73-8	Price
methanol	10 µg/ml	10 ml	ref. P856790	BG
	100 µg/ml	1 ml	ref. P856800	BG
1,2,4-TRIMETHYLBENZENE			CAS 95-63-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P856830	BG
	100 µg/ml	1 ml	ref. P856840	BG
1,3,5-TRIMETHYLBENZENE			CAS 108-67-8	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P856870	BG
	100 µg/ml	1 ml	ref. P856880	BG
3,3,5-TRIMETHYLCYCLOHEXYLSALICYLATE			CAS 118-56-9	Price
methanol	100 µg/ml	1 ml	ref. P873220	BJ
2,3,5-TRIMETHYLNAPHTHALENE			CAS 2245-38-7	Price
cyclohexane	10 µg/ml	10 ml	ref. P872630	CE
TRINEXAPAC-ETHYL			CAS 95266-40-3	Price
acetone	10 µg/ml	10 ml	ref. P856910	BG
	100 µg/ml	1 ml	ref. P856920	BG
cyclohexane	10 µg/ml	10 ml	ref. P856930	BG
	100 µg/ml	1 ml	ref. P856940	BG
TRIPHENYLENE			CAS 217-59-4	Price
acetonitrile	10 µg/ml	10 ml	ref. P857170	BH
	100 µg/ml	1 ml	ref. P857180	BH
cyclohexane	10 µg/ml	10 ml	ref. P857190	BH
	100 µg/ml	1 ml	ref. P857200	BH
TRIPHENYLPHOSPHATE			CAS 115-86-6	Price
ethyl acetate	10 µg/ml	10 ml	ref. P857250	BG
	100 µg/ml	1 ml	ref. P857260	BG
methyl-tert-butylether	10 µg/ml	10 ml	ref. P857230	BG
	100 µg/ml	1 ml	ref. P857240	BG
	500 µg/ml	1 ml	ref. F113331	BG
TRIPROPYLtin CHLORIDE			CAS 2279-76-7	Price
dichloromethane	2 000 µg/ml	1 ml	ref. P879030	CB
n-UNDECANE			CAS 1120-21-4	Price
methanol	10 µg/ml	10 ml	ref. P857330	BG
	100 µg/ml	1 ml	ref. P857340	BG
VALERALDEHYDE			CAS 110-62-3	Price
acetonitrile	1 000 µg/ml	5 ml	ref. P879040	CF
methanol	1 000 µg/ml	1 ml	ref. P879050	BH
VALERALDEHYDE-DNPH			CAS 2057-84-3	Price
acetonitrile	10 µg/ml	10 ml	ref. P857370	BG
	100 µg/ml	1 ml	ref. P857380	BG
VAMIDOTHION			CAS 2275-23-2	Price
acetonitrile	10 µg/ml	10 ml	ref. P857410	BG
	100 µg/ml	1 ml	ref. P857420	BG
VERNOLATE			CAS 1929-77-7	Price
methanol	100 µg/ml	1 ml	ref. P879060	BG

VINCLOZOLIN			CAS 50471-44-8	Price
acetonitrile	10 µg/ml	10 ml	ref. P857450	BG
	100 µg/ml	1 ml	ref. P857460	BG
cyclohexane	10 µg/ml	10 ml	ref. P857470	BG
	100 µg/ml	1 ml	ref. P857480	BG
iso-octane	10 µg/ml	10 ml	ref. P857490	BG
	100 µg/ml	1 ml	ref. P857500	BG
VINYLACETATE			CAS 108-05-4	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P857570	BG
	100 µg/ml	1 ml	ref. P857580	BG
VINYLCHLORIDE			CAS 75-01-4	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P857610	BG
	100 µg/ml	1 ml	ref. P857620	BG
WARFARIN			CAS 81-81-2	Price
acetonitrile	100 µg/ml	1 ml	ref. P880540	BG
m -XYLENE			CAS 108-38-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P857650	BG
	100 µg/ml	1 ml	ref. P857660	BG

o -XYLENE			CAS 95-47-6	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P857690	BG
	100 µg/ml	1 ml	ref. P857700	BG
P -XYLENE			CAS 106-42-3	Price
methanol purge & trap	10 µg/ml	10 ml	ref. P857730	BG
	100 µg/ml	1 ml	ref. P857740	BG
XYLENE (Mixed Isomers)			CAS 1330-20-7	Price
methanol purge & trap	100 µg/ml	1 ml	ref. P879080	BG
ZEARALENON			CAS 17924-92-4	Price
acetonitrile	100 µg/ml	5 ml	ref. P879090	BFK
ZIRAM			CAS 137-30-4	Price
acetonitrile	100 µg/ml	1 ml	ref. P879100	BG
methanol	100 µg/ml	1 ml	ref. P879110	BG



European Pharmacopoeia Reagents and Standards

Our Range of Standards *CPA-Pharm*



- Primary Opalescence Suspension
- Primary Colour Solutions
- Conductivity Standard Solutions
- Reference Buffer Solutions
- Reagents and Volumetric Solutions
- Ready to use Volumetric Solutions

The complete range is :

- manufactured in accordance with the European Pharmacopoeia,
- traceable to primary standards from NIST-CRM,
- provided with a Certificate of Analysis.

Quotation Request

Custom Made Mixture for GC, HPLC, GC-MS & LC-MS

Please send this form by fax or post. (photocopy for repeat use)

Name : Tel :

Company : Fax :

Address :

E-mail :

Mixture Description

Solvent :

Number of components : Quantity :

Concentration units : mg/ml ☐ µg/ml ☐ ng/ml ☐ other :

Ampoules

Volume per ampoule : 1 ml ☐ 2 ml ☐
5 ml ☐ 10 ml ☐

Headspace Ampoules

Volume per ampoule : 10 ml ☐ 20 ml ☐
Others :

Certain[®] Bottles

Volume per bottle : 1,5 ml ☐ 4,5 ml ☐
10 ml ☐

Amber bottles

Volume per bottle : 20 ml ☐ 40 ml ☐
100 ml ☐

Compound Name	CAS number	Conc.

Compound Name	CAS number	Conc.

Additionnals Comments :
.....
.....
.....

Certificate of Analysis

REFERENCE MATERIAL CALIBRATION

Reference Materials Calibration Laboratory

CPA Ltd

(*)

№ S24760

Date: 06.04.2010

**Description of the RM
subject to calibration:**

n-Alkane Standard Solution 6 components acc. ISO 9377-2-Mod
in n-Hexane

Catalog number:

9BAD.K1.ZT.1

Calibration method:

CRM's calibration procedure (MK-003)

**Result from calibration
(certified value)
/Uncertainty:**

Component	Chem. Formula	CAS No.	Certified Value and uncertainty (mg/l)*
n-Decane	C10H22	124-18-5	50.27 +/- 0.70
n-Hexadecane	C16H34	544-76-3	50.71 +/- 1.31
n-Eicosane	C20H42	112-95-8	49.97 +/- 1.68
n-Tetracosane	C24H50	646-31-1	50.06 +/- 1.18
n-Triacontane	C30H62	638-68-6	50.45 +/- 0.94
n-Tetracontane	C40H82	4181-95-7	50.09 +/- 2.01

* The expanded traceability of measurement is expressed as a mean quadratic uncertainty of measurement multiplied by the coverage factor of confidence interval $k = 2$ that corresponds to a confidence level of 95 % at a normal distribution. The mean quadratic uncertainty is determined in compliance with EA 4/02.

Metrological traceability:

n-Decane	ChemService, O-729
n-Hexadecane	ChemService, O-734
n-Eicosane	ChemService, F2183
n-Tetracosane	Dr Ehrenstorfer, C 17395800
n-Triacontane	ChemService, F2297
n-Tetracontane	Dr Ehrenstorfer, C 17395500

The metrological traceability is assured through calibration on GC/MS. The calibration curve is drawn using a series of standard solutions prepared from a certified reference material of Dr Ehrenstorfer, traceable to SI. All contributions in relation to the preparation of standard solutions are considered when evaluating the uncertainty.

This certified reference material is produced through weighing and dilution of the single compounds in Iso-octane.

Class A laboratory glassware is used.

The measurement results are traceable to SI. All analytical balances used for the preparation of the solution are periodically calibrated by an accredited calibration laboratory and are daily checked under an in-house procedure.

The results from temperature measurement are traceable to SI. The thermometers used for solution's calibration are calibrated from an ISO 17025 accredited laboratory. The ambient conditions are controlled with a hygrometer calibrated from an ISO 17025 accredited laboratory.

**Ambient conditions during
calibration:**

Temperature: 22.6 °C +/- 0.2 °C
Moisture: 48 % +/- 2 %



Certificate of Analysis

Starting material, purity, % (Lot N):	n-Decane	99.5	41054994
	n-Hexadecane	99	41048931
	n-Eicosane	99	41011256
	n-Tetracosane	99.5	41020807
	n-Triacontane	98.8	41011898
	n-Tetracontane	97.6	41011195

Certified density: 0.661 g/ml measured at 19.3°C

Calibration date (when applicable)

Expiration date: till 05.2012

Intended use:

This CRM is intended for:

- Calibration of TLC, GC/FID, GC/TCD, GC/ECD, GC/MS, HPLC/UV and HPLC/MS.
- Validation of analytical methods
- Preparation of "working reference samples"
- For detection limit and linearity studies

This certified reference material can be used directly or can be diluted in an appropriate high-purity matrix. Only a clean class A glassware should be used. Obtained concentration (in mg/l) after dilution is a result from the multiplication of certified value of CRM concentration and the CRM's volume used for dilution and divided into the flask's volume used for dilution.

Stability and storage:

This CRM is with a guaranteed stability for a period of 12 months. Stability is guaranteed provided that the solution is kept in its original packaging, at darkness and temperatures below 30°C.

Homogeneity:

This solution was mixed according to an in-house procedure and is guaranteed to be homogeneous.

This Certified reference Material was produced under ISO 9001:2008 Quality Control System.

The instructions of ISO Guide 34-2000 "General Requirements for the Competence of Reference Material Producers" were considered for the preparation of this solution.

Certificate of analysis was prepared according to ISO Guide 31:2000 „Reference materials - Contents of certificates and labels"

CALIBRATION PERFORMED BY:

Yordan Uzunov
(name, family name, signature)

MANAGER OF RMCL :

Krassimira Taralova
(name, family name, signature)

end

Additional information:

Gravimetric Data

Component	Purity %	Source Lot No	Weighed quantity, g	Final quantity, kg.10 ⁻³	Bulk/ Standard Solution lot No	Concentration mg/kg	Chemist ID
n-Decane	99.5	41054994	0.01246	2.6600	91095442	4660.8	AS
		91095442	1.0784	66.146	92076785	75.986	PA
n-Hexadecane	99	41048931	0.01188	2.4774	91095541	4747.4	AS
		91095541	1.0682	66.146	92076785	76.667	PA
n-Eicosane	99	41011256	0.01614	2.6545	91095459	6019.4	AS
		91095459	0.8303	66.146	92076785	75.559	PA
n-Tetracosane	99.5	41020807	0.0111	2.5824	91095558	4276.8	AS
		91095558	1.1706	66.146	92076785	75.689	PA
n-Triacontane	98.8	41011898	0.01436	2.7502	91095510	5158.8	AS
		91095510	0.978	66.146	92076785	76.275	PA
n-Tetracontane	97.6	41011195	0.01492	13.3611	91076670	1089.87	AS
		91076670	4.5966	66.146	92076785	75.737	PA

CAS Number Index

CAS	Page	CAS	Page	CAS	Page	CAS	Page	CAS	Page
50 - 00 - 0	128	86 - 30 - 6	133	103 - 33 - 3	113	121 - 73 - 3	117	527 - 20 - 8	140
50 - 27 - 1	125	86 - 50 - 0	113	103 - 65 - 1	142	121 - 75 - 5	130	527 - 53 - 7	145
50 - 29 - 3	120	86 - 57 - 7	133	104 - 40 - 5	133	121 - 86 - 8	117	528 - 29 - 0	124
50 - 32 - 8	114	86 - 73 - 7	127	104 - 51 - 8	115	122 - 14 - 5	126	533 - 74 - 4	119
50 - 65 - 7	133	86 - 74 - 8	116	104 - 72 - 3	120	122 - 34 - 9	143	534 - 52 - 1	131
51 - 03 - 6	141	87 - 40 - 1	146	104 - 87 - 0	145	122 - 39 - 4	124	535 - 89 - 7	119
51 - 28 - 5	124	87 - 60 - 5	117	105 - 05 - 5	123	122 - 42 - 9	142	540 - 36 - 3	123
51 - 36 - 5	121	87 - 61 - 6	146	105 - 60 - 2	115	122 - 66 - 7	124	541 - 73 - 1	121
52 - 68 - 6	146	87 - 62 - 7	123	105 - 67 - 9	124	123 - 07 - 9	126	542 - 75 - 6	122
52 - 85 - 7	126	87 - 65 - 0	122	106 - 32 - 1	126	123 - 19 - 3	129	544 - 76 - 3	129
53 - 16 - 7	125	87 - 68 - 3	129	106 - 35 - 4	129	123 - 38 - 6	142	554 - 00 - 7	121
53 - 19 - 0	120	87 - 82 - 1	129	106 - 39 - 8	114	123 - 72 - 8	115	555 - 37 - 3	133
53 - 70 - 3	120	87 - 86 - 5	140	106 - 40 - 1	114	123 - 73 - 9	119	556 - 88 - 7	133
55 - 38 - 9	127	88 - 06 - 2	147	106 - 42 - 3	148	123 - 91 - 1	124	557 - 59 - 5	145
56 - 23 - 5	144	88 - 72 - 2	133	106 - 43 - 4	118	124 - 10 - 7	145	563 - 12 - 2	125
56 - 35 - 9	144	88 - 73 - 3	117	106 - 44 - 5	132	124 - 13 - 0	134	563 - 58 - 6	122
56 - 38 - 2	134	88 - 75 - 5	133	106 - 46 - 7	121	124 - 18 - 5	120	569 - 61 - 9	134
56 - 55 - 3	113	88 - 85 - 7	124	106 - 47 - 8	117	124 - 19 - 6	133	575 - 43 - 9	124
56 - 72 - 4	119	89 - 59 - 8	117	106 - 48 - 9	118	124 - 40 - 3	123	576 - 24 - 9	122
57 - 11 - 4	134	89 - 61 - 2	122	106 - 89 - 8	125	124 - 48 - 1	120	576 - 26 - 1	124
57 - 74 - 9	116	89 - 63 - 4	117	106 - 93 - 4	121	126 - 73 - 8	146	577 - 19 - 5	115
57 - 97 - 6	124	90 - 00 - 6	126	106 - 99 - 0	115	126 - 75 - 0	120	581 - 89 - 5	133
58 - 89 - 9	128	90 - 12 - 0	132	107 - 02 - 8	111	126 - 98 - 7	131	583 - 78 - 8	122
58 - 90 - 2	144	90 - 13 - 1	117	107 - 04 - 0	114	126 - 99 - 8	118	584 - 79 - 2	111
59 - 50 - 7	117	90 - 43 - 7	141	107 - 05 - 1	118	127 - 18 - 4	144	590 - 86 - 3	130
60 - 29 - 7	123	91 - 20 - 3	132	107 - 06 - 2	122	129 - 00 - 0	142	591 - 35 - 5	122
60 - 51 - 5	123	91 - 53 - 2	126	107 - 07 - 3	117	131 - 11 - 3	141	591 - 78 - 6	129
60 - 57 - 1	123	91 - 57 - 6	132	107 - 10 - 8	142	131 - 89 - 5	119	594 - 20 - 7	122
61 - 82 - 5	112	91 - 58 - 7	117	107 - 12 - 0	142	132 - 64 - 9	120	598 - 72 - 1	115
62 - 53 - 3	112	91 - 59 - 8	112	107 - 13 - 1	111	132 - 65 - 0	120	600 - 05 - 5	121
62 - 73 - 7	123	91 - 94 - 1	121	107 - 14 - 2	117	133 - 06 - 2	115	602 - 60 - 8	133
62 - 75 - 9	133	92 - 52 - 4	114	107 - 21 - 1	126	133 - 07 - 3	128	603 - 71 - 4	133
63 - 25 - 2	115	92 - 67 - 1	112	107 - 83 - 5	132	133 - 90 - 4	116	606 - 20 - 2	124
64 - 17 - 5	125	92 - 86 - 4	120	107 - 87 - 9	140	134 - 62 - 3	120	607 - 57 - 8	133
65 - 85 - 0	114	92 - 87 - 5	113	108 - 05 - 4	148	135 - 01 - 3	123	608 - 27 - 5	121
66 - 77 - 3	132	92 - 94 - 4	144	108 - 10 - 1	132	135 - 98 - 8	115	608 - 31 - 1	121
66 - 81 - 9	119	93 - 71 - 0	111	108 - 38 - 3	148	137 - 17 - 7	147	608 - 33 - 3	121
67 - 56 - 1	131	93 - 72 - 1	126	108 - 39 - 4	132	137 - 26 - 8	145	608 - 73 - 1	129
67 - 66 - 3	117	93 - 76 - 5	143	108 - 41 - 8	118	137 - 30 - 4	148	608 - 93 - 5	140
67 - 72 - 1	129	94 - 13 - 3	129	108 - 42 - 9	117	139 - 40 - 2	142	609 - 19 - 8	147
69 - 72 - 7	143	94 - 26 - 8	129	108 - 43 - 0	118	140 - 56 - 7	126	610 - 39 - 9	124
70 - 30 - 4	129	94 - 74 - 6	130	108 - 46 - 3	143	140 - 57 - 8	112	611 - 06 - 3	122
71 - 36 - 3	115	94 - 75 - 7	119	108 - 60 - 1	114	140 - 66 - 9	134	611 - 14 - 3	126
71 - 43 - 2	113	94 - 80 - 4	119	108 - 67 - 8	147	141 - 43 - 5	125	615 - 58 - 7	121
71 - 55 - 6	146	94 - 81 - 5	130	108 - 70 - 3	146	141 - 66 - 2	123	615 - 59 - 8	121
72 - 20 - 8	125	94 - 82 - 6	119	108 - 86 - 1	114	141 - 78 - 6	126	618 - 62 - 2	122
72 - 43 - 5	131	95 - 13 - 6	129	108 - 88 - 3	145	141 - 93 - 5	123	618 - 87 - 1	124
72 - 54 - 8	120	95 - 47 - 6	148	108 - 90 - 7	117	142 - 28 - 9	122	620 - 14 - 4	126
72 - 55 - 9	120	95 - 48 - 7	132	108 - 94 - 1	119	142 - 82 - 5	129	620 - 17 - 7	126
72 - 56 - 0	140	95 - 49 - 8	118	108 - 95 - 2	140	143 - 50 - 0	116	621 - 64 - 7	133
74 - 83 - 9	115	95 - 50 - 1	121	109 - 66 - 0	140	148 - 79 - 8	145	622 - 96 - 8	126
74 - 87 - 3	117	95 - 51 - 2	117	109 - 79 - 5	115	150 - 50 - 5	131	626 - 39 - 1	146
74 - 88 - 4	130	95 - 53 - 4	145	109 - 86 - 4	126	150 - 68 - 5	132	626 - 43 - 7	121
74 - 93 - 1	132	95 - 57 - 8	117	109 - 89 - 7	123	156 - 59 - 2	122	629 - 50 - 5	147
74 - 95 - 3	121	95 - 63 - 6	147	109 - 99 - 9	145	156 - 60 - 5	122	629 - 59 - 4	145
74 - 97 - 5	114	95 - 65 - 8	124	110 - 42 - 9	120	189 - 55 - 9	120	629 - 78 - 7	129
75 - 00 - 3	117	95 - 69 - 2	117	110 - 43 - 0	129	189 - 64 - 0	120	629 - 92 - 5	133
75 - 01 - 4	148	95 - 74 - 9	117	110 - 56 - 5	121	191 - 07 - 1	119	629 - 99 - 2	140
75 - 04 - 7	126	95 - 76 - 1	121	110 - 57 - 6	121	191 - 24 - 2	114	630 - 01 - 3	129
75 - 05 - 8	111	95 - 77 - 2	122	110 - 62 - 3	147	192 - 65 - 4	120	630 - 02 - 4	134
75 - 07 - 0	111	95 - 79 - 4	117	110 - 75 - 8	117	192 - 97 - 2	114	630 - 20 - 6	144
75 - 08 - 1	126	95 - 80 - 7	120	110 - 82 - 7	119	193 - 39 - 5	129	634 - 66 - 2	144
75 - 09 - 2	122	95 - 82 - 9	121	110 - 86 - 1	143	198 - 55 - 0	140	634 - 90 - 2	144
75 - 15 - 0	116	95 - 85 - 2	112	111 - 11 - 5	134	205 - 82 - 3	114	636 - 30 - 6	146
75 - 21 - 8	126	95 - 87 - 4	124	111 - 42 - 2	123	205 - 99 - 2	114	637 - 92 - 3	126
75 - 25 - 2	115	95 - 93 - 2	145	111 - 44 - 4	114	206 - 44 - 0	127	638 - 36 - 8	141
75 - 27 - 4	114	95 - 94 - 3	144	111 - 65 - 9	134	207 - 08 - 9	114	638 - 68 - 6	146
75 - 34 - 3	122	95 - 95 - 4	147	111 - 71 - 7	129	208 - 96 - 8	111	639 - 58 - 7	127
75 - 35 - 4	122	96 - 12 - 8	121	111 - 82 - 0	125	213 - 46 - 7	141	640 - 15 - 3	145
75 - 65 - 0	115	96 - 18 - 4	147	111 - 84 - 2	133	214 - 17 - 5	113	644 - 64 - 4	124
75 - 69 - 4	128	96 - 22 - 0	140	111 - 85 - 3	117	217 - 59 - 4	147	646 - 31 - 1	145
75 - 96 - 7	146	96 - 23 - 1	122	111 - 87 - 5	134	218 - 01 - 9	118	672 - 99 - 1	114
75 - 99 - 0	119	96 - 33 - 3	111	111 - 91 - 1	114	224 - 42 - 0	120	683 - 18 - 1	121
76 - 06 - 2	118	97 - 00 - 7	117	112 - 31 - 2	120	226 - 36 - 8	120	694 - 80 - 4	114
76 - 13 - 1	147	97 - 17 - 6	121	112 - 40 - 3	125	297 - 78 - 9	130	709 - 98 - 8	142
76 - 44 - 8	128	97 - 23 - 4	122	112 - 56 - 1	130	298 - 00 - 0	134	725 - 00 - 8	142
77 - 26 - 9	115	97 - 53 - 0	126	112 - 61 - 8	134	298 - 02 - 2	141	731 - 27 - 1	145
77 - 47 - 4	129	97 - 63 - 2	131	112 - 95 - 8	125	298 - 04 - 4	124	732 - 11 - 6	141
78 - 48 - 8	146	98 - 01 - 1	128	113 - 48 - 4	132	299 - 84 - 3	126	741 - 58 - 2	113
78 - 59 - 1	130	98 - 06 - 6	115	114 - 26 - 1	142	300 - 76 - 5	132	759 - 94 - 4	125
78 - 70 - 6	130	98 - 08 - 8	147	115 - 29 - 7	125	301 - 12 - 2	120	771 - 61 - 9	140
78 - 75 - 1	121	98 - 54 - 4	115	115 - 32 - 2	123	302 - 17 - 0	116	786 - 19 - 6	116
78 - 83 - 1	132	98 - 55 - 5	144	115 - 86 - 6	147	309 - 00 - 2	111	789 - 02 - 6	120
78 - 87 - 5	122	98 - 82 - 8	130	115 - 90 - 2	127	311 - 45 - 5	134	832 - 69 - 9	132
78 - 88 - 6	122	98 - 86 - 2	111	116 - 06 - 3	111	314 - 40 - 9	114	834 - 12 - 8	111
78 - 92 - 2	115	98 - 95 - 3	133	116 - 29 - 0	145	319 - 84 - 6	128	871 - 83 - 0	132
78 - 93 - 3	115	99 - 08 - 1	133	117 - 18 - 0	144	319 - 85 - 7	128	876 - 27 - 7	118
79 - 00 - 5	146	99 - 30 - 9	123	117 - 80 - 6	121	319 - 86 - 8	128	877 - 09 - 8	145
79 - 01 - 6	146	99 - 54 - 7	122	117 - 81 - 7	141	321 - 38 - 0	127	886 - 50 - 0	144
79 - 06 - 1	111	99 - 55 - 8	112	117 - 84 - 0	141	321 - 60 - 8	127	888 - 54 - 0	111
79 - 34 - 5	144	99 - 65 - 0	124	118 - 56 - 9	147	327 - 98 - 0	146	892 - 21 - 7	133
79 - 46 - 9	133	99 - 76 - 3	129	118 - 58 - 1	114	330 - 54 - 1	125	900 - 95 - 8	127
80 - 05 - 7	114	99 - 87 - 6	130	118 - 60 - 5	126	330 - 55 - 2	130	919 - 86 - 8	120
80 - 33 - 1	116	99 - 99 - 0	133	118 - 74 - 1	129	333 - 41 - 5	120	929 - 77 - 7	125
80 - 38 - 6	127	100 - 00 - 5	117	118 - 79 - 6	146	348 - 51 - 6	117	929 - 97 - 0	125
80 - 62 - 6	131	100 - 02 - 7	133	118 - 91 - 2	117	352 - 33 - 0	117	933 - 75 - 5	146
80 - 62 - 6	132	100 - 41 - 4	126	119 - 12 - 0	142	354 - 58 - 5	147	933 - 78 - 8	146
81 - 20 - 9	124	100 - 42 - 5	143	119 - 64 - 2	145	363 - 72 - 4	140	935 - 95 - 5	144
81 - 81 - 2	148	100 - 44 - 7	118	119 - 90 - 4	123	367 - 12 - 4	128	950 - 35 - 6	134
82 - 68 - 8	143	100 - 52 - 7	113	119 - 93 - 7	124	434 - 90 - 2	119	950 - 37 - 8	131
83 - 32 - 9	111	101 - 14 - 4	131	120 - 12 - 7	112	438 - 22 - 2	112	957 - 51 - 7	124
83 - 42 - 1	117	101 - 21 - 3	118	120 - 36 - 5	122	460 - 00 - 4	115	958 - 60 - 1	115
83 - 79 - 4	143	101 - 27 - 9	113	120 - 47 - 8	129	462 - 06 - 6	127	959 - 98 - 8	125
84 - 15 - 1	144	101 - 42 - 8	127	120 - 51 - 4	114	465 - 73 - 6			

CAS Number Index

CAS	Page	CAS	Page	CAS	Page	CAS	Page	CAS	Page	CAS	Page
2425 - 06 - 1	115	10552 - 74 - 6	133	32598 - 11 - 1	136	41483 - 43 - 6	115	64902 - 72 - 3	118	81405 - 85 - 8	129
2436 - 73 - 9	130	10605 - 21 - 7	116	32598 - 12 - 2	136	41814 - 78 - 2	147	65510 - 44 - 3	137	81406 - 37 - 3	128
2437 - 79 - 8	136	11096 - 82 - 5	112	32598 - 13 - 3	136	42509 - 80 - 8	130	65510 - 45 - 4	136	81777 - 89 - 1	118
2439 - 01 - 2	116	11097 - 69 - 1	112	32598 - 14 - 4	137	42576 - 02 - 3	114	65907 - 30 - 4	128	82097 - 50 - 5	146
2442 - 49 - 1	145	11100 - 14 - 4	112	32690 - 93 - 0	136	42740 - 50 - 1	139	66003 - 55 - 2	111	82558 - 50 - 7	130
2463 - 84 - 5	121	11104 - 28 - 2	112	32774 - 16 - 6	138	42874 - 03 - 3	134	66063 - 05 - 6	140	82560 - 54 - 1	113
2497 - 06 - 5	124	11141 - 16 - 5	112	32809 - 16 - 8	141	43121 - 43 - 3	146	66215 - 27 - 8	119	82657 - 04 - 3	114
2497 - 07 - 6	124	12672 - 29 - 6	112	33025 - 41 - 1	136	50471 - 44 - 8	148	66230 - 04 - 4	125	83055 - 99 - 6	113
2531 - 84 - 2	132	12674 - 11 - 2	112	33089 - 61 - 1	112	50563 - 36 - 5	123	66246 - 88 - 6	140	83121 - 18 - 0	144
2540 - 82 - 1	128	13029 - 08 - 8	134	33091 - 17 - 7	139	50594 - 66 - 6	111	66332 - 96 - 5	128	83164 - 33 - 4	123
2571 - 00 - 8	145	13031 - 39 - 5	118	33146 - 45 - 1	135	51218 - 45 - 2	132	66753 - 07 - 9	144	83657 - 24 - 3	124
2588 - 03 - 6	141	13067 - 93 - 1	119	33213 - 65 - 9	125	51218 - 49 - 6	141	66841 - 25 - 6	146	84087 - 01 - 4	143
2588 - 04 - 7	141	13071 - 79 - 9	144	33245 - 39 - 5	127	51235 - 04 - 2	129	67129 - 08 - 2	131	84332 - 86 - 5	118
2593 - 15 - 9	126	13127 - 88 - 3	141	33284 - 50 - 3	134	51338 - 27 - 3	123	67306 - 00 - 7	127	84852 - 15 - 3	133
2595 - 54 - 2	130	13171 - 21 - 6	141	33284 - 52 - 5	136	51630 - 58 - 1	127	67375 - 30 - 8	119	85380 - 74 - 1	140
2597 - 03 - 7	141	13194 - 48 - 4	126	33284 - 53 - 6	136	51908 - 16 - 8	138	67564 - 91 - 4	127	85509 - 19 - 9	128
2599 - 11 - 3	143	13250 - 98 - 1	120	33284 - 54 - 7	136	52315 - 07 - 8	119	67747 - 09 - 5	141	86194 - 41 - 4	128
2631 - 37 - 0	141	13356 - 08 - 6	126	33543 - 31 - 6	132	52645 - 53 - 1	140	68194 - 04 - 7	132	86479 - 06 - 3	129
2631 - 40 - 5	130	13360 - 45 - 7	116	33629 - 47 - 9	115	52663 - 58 - 8	136	68194 - 05 - 8	137	87130 - 20 - 9	123
2635 - 10 - 1	131	13457 - 18 - 6	142	33693 - 04 - 8	144	52663 - 60 - 2	136	68194 - 06 - 9	137	87237 - 48 - 7	128
2636 - 26 - 2	119	13593 - 03 - 8	143	33820 - 53 - 0	130	52663 - 61 - 3	137	68194 - 07 - 0	137	87674 - 68 - 8	123
2642 - 71 - 9	113	13654 - 09 - 6	134	33979 - 03 - 2	138	52663 - 62 - 4	136	68194 - 08 - 1	138	88283 - 41 - 4	143
2675 - 77 - 6	117	13684 - 63 - 4	140	34014 - 18 - 1	144	52663 - 63 - 5	138	68194 - 09 - 2	138	88671 - 89 - 0	132
2778 - 96 - 3	143	14214 - 32 - 5	123	34123 - 59 - 6	130	52663 - 64 - 6	139	68194 - 10 - 5	137	90717 - 03 - 6	143
2880 - 05 - 9	145	14255 - 88 - 0	126	34205 - 21 - 5	123	52663 - 65 - 7	139	68194 - 11 - 6	137	90982 - 32 - 4	116
2905 - 67 - 1	121	14437 - 17 - 3	116	34256 - 82 - 1	111	52663 - 67 - 9	139	68194 - 12 - 7	137	91465 - 08 - 6	119
2921 - 88 - 2	118	14484 - 64 - 1	127	34643 - 46 - 4	142	52663 - 68 - 0	139	68194 - 13 - 8	138	93951 - 66 - 7	114
2974 - 90 - 5	135	14816 - 18 - 3	141	34681 - 10 - 2	115	52663 - 68 - 8	138	68194 - 14 - 9	138	93951 - 69 - 0	127
2974 - 92 - 7	135	15067 - 26 - 2	111	34883 - 39 - 1	135	52663 - 69 - 1	139	68194 - 15 - 0	138	93951 - 74 - 7	122
3017 - 95 - 6	114	15165 - 67 - 0	123	34883 - 41 - 5	135	52663 - 70 - 4	139	68194 - 16 - 1	139	93951 - 75 - 8	124
3060 - 89 - 7	132	15299 - 99 - 7	132	34883 - 43 - 7	135	52663 - 71 - 5	139	68194 - 17 - 2	139	93951 - 79 - 2	133
3114 - 55 - 4	117	15310 - 01 - 7	113	35065 - 27 - 1	138	52663 - 72 - 6	138	68359 - 37 - 5	119	93951 - 97 - 4	111
3196 - 15 - 4	131	15545 - 48 - 9	118	35065 - 28 - 2	138	52663 - 73 - 7	139	68694 - 11 - 1	147	93951 - 98 - 5	114
3209 - 22 - 1	122	15862 - 07 - 4	135	35065 - 29 - 3	139	52663 - 74 - 8	139	69327 - 76 - 0	115	93952 - 12 - 6	141
3244 - 90 - 4	112	15968 - 05 - 5	136	35065 - 30 - 6	138	52663 - 75 - 9	139	69335 - 90 - 6	127	93952 - 14 - 8	129
3337 - 71 - 1	113	15972 - 60 - 8	111	35256 - 85 - 0	144	52663 - 76 - 0	139	69335 - 91 - 7	127	94125 - 34 - 5	142
3347 - 22 - 6	125	16118 - 49 - 3	116	35367 - 38 - 5	123	52663 - 77 - 1	139	69377 - 81 - 7	128	95266 - 40 - 3	147
3351 - 28 - 8	131	16605 - 91 - 7	134	35400 - 43 - 2	143	52663 - 78 - 2	139	69782 - 90 - 7	138	95617 - 09 - 7	126
3383 - 96 - 8	144	16606 - 02 - 3	135	35554 - 44 - 0	129	52663 - 79 - 3	139	69782 - 91 - 8	139	95737 - 68 - 1	143
3386 - 33 - 2	117	16606 - 02 - 6	116	35572 - 78 - 2	112	52704 - 70 - 8	138	69806 - 34 - 4	128	96489 - 71 - 3	142
3397 - 62 - 4	113	16752 - 77 - 5	131	35575 - 96 - 3	113	52712 - 04 - 6	138	69806 - 50 - 4	127	96525 - 23 - 4	128
3424 - 82 - 6	120	17040 - 19 - 6	120	35693 - 92 - 6	135	52712 - 05 - 7	139	70124 - 77 - 5	127	98730 - 04 - 2	113
3567 - 62 - 2	122	17060 - 07 - 0	122	35693 - 99 - 3	136	52744 - 13 - 5	138	70356 - 09 - 1	134	98886 - 44 - 3	128
3674 - 09 - 7	122	17109 - 49 - 8	125	35694 - 04 - 3	138	52756 - 22 - 6	127	70362 - 41 - 3	137	99105 - 77 - 8	143
3689 - 24 - 5	143	17606 - 31 - 4	113	35694 - 06 - 5	138	52888 - 80 - 9	142	70362 - 45 - 7	136	99607 - 70 - 2	118
3691 - 35 - 8	117	17640 - 02 - 7	119	35694 - 08 - 7	139	52918 - 63 - 5	120	70362 - 46 - 8	135	99675 - 03 - 3	130
3697 - 24 - 3	131	17647 - 74 - 4	124	36559 - 22 - 5	135	53112 - 28 - 0	143	70362 - 47 - 9	136	100155 - 47 - 3	120
3734 - 48 - 3	116	17804 - 35 - 2	113	36734 - 19 - 7	130	53469 - 21 - 9	112	70362 - 48 - 0	136	101007 - 06 - 1	111
3761 - 41 - 9	127	17847 - 51 - 7	122	37019 - 18 - 4	143	53494 - 70 - 5	125	70362 - 49 - 1	136	101463 - 69 - 8	127
3761 - 42 - 0	127	17924 - 92 - 4	148	37324 - 23 - 5	112	53555 - 66 - 1	135	70362 - 50 - 4	136	102851 - 06 - 9	144
3766 - 60 - 7	115	18181 - 70 - 9	130	37680 - 65 - 2	135	53663 - 59 - 9	135	70424 - 67 - 8	136	103361 - 09 - 7	127
3766 - 81 - 2	126	18181 - 80 - 1	115	37680 - 66 - 3	135	53780 - 34 - 0	131	70424 - 68 - 9	137	104040 - 78 - 0	127
3813 - 05 - 6	113	18259 - 05 - 7	137	37680 - 68 - 5	135	54230 - 22 - 7	136	70424 - 69 - 0	137	105512 - 06 - 9	118
3855 - 82 - 1	121	18625 - 12 - 2	119	37680 - 69 - 6	135	55179 - 31 - 2	114	70424 - 70 - 3	137	105827 - 78 - 9	129
3878 - 19 - 1	128	18691 - 97 - 9	131	37680 - 73 - 2	137	55215 - 17 - 3	137	70630 - 17 - 0	131	106325 - 08 - 0	125
3970 - 40 - 9	117	19044 - 88 - 3	134	37764 - 25 - 3	121	55215 - 18 - 4	138	71283 - 80 - 2	126	107534 - 96 - 3	144
4165 - 60 - 0	133	19406 - 51 - 0	112	38260 - 54 - 7	126	55219 - 65 - 3	146	71422 - 67 - 8	116	110488 - 70 - 5	123
4165 - 61 - 1	112	19666 - 30 - 9	134	38379 - 99 - 6	137	55283 - 68 - 6	125	71626 - 11 - 4	113	111479 - 05 - 1	142
4165 - 62 - 2	140	19719 - 28 - 9	122	38380 - 01 - 7	137	55283 - 68 - 6	125	72178 - 02 - 0	128	111988 - 49 - 9	145
4274 - 02 - 3	129	19937 - 59 - 8	132	38380 - 02 - 8	137	55290 - 64 - 7	123	73250 - 68 - 7	131	111991 - 09 - 4	133
4525 - 75 - 1	118	21087 - 64 - 9	132	38380 - 03 - 9	137	55312 - 69 - 1	136	73557 - 53 - 8	136	112281 - 77 - 3	145
4658 - 28 - 0	113	21609 - 90 - 5	130	38380 - 04 - 0	138	55335 - 06 - 3	147	73575 - 52 - 7	136	112410 - 23 - 8	144
4824 - 78 - 6	115	21725 - 46 - 2	119	38380 - 05 - 1	138	55406 - 53 - 6	130	73575 - 54 - 9	137	113096 - 99 - 4	130
4841 - 20 - 7	126	22224 - 92 - 6	126	38380 - 07 - 3	137	55512 - 33 - 9	142	73575 - 55 - 0	137	114369 - 43 - 6	126
4901 - 51 - 3	144	22248 - 79 - 9	145	38380 - 08 - 4	138	55702 - 45 - 9	135	73575 - 56 - 1	137	116255 - 48 - 2	115
5074 - 71 - 5	119	22781 - 23 - 3	113	38411 - 22 - 2	138	55702 - 46 - 0	135	73575 - 57 - 2	137	117428 - 22 - 5	141
5103 - 71 - 9	116	23311 - 98 - 2	141	38411 - 25 - 5	139	55712 - 37 - 3	135	74070 - 46 - 5	111	118134 - 30 - 8	143
5103 - 73 - 1	133	23135 - 22 - 0	134	38444 - 73 - 4	135	55720 - 44 - 0	135	74115 - 24 - 5	118	118712 - 89 - 3	146
5103 - 74 - 2	116	23184 - 66 - 9	115	38444 - 76 - 7	135	55814 - 41 - 0	131	74223 - 64 - 6	132	119168 - 77 - 3	144
5131 - 24 - 8	124	23399 - 90 - 8	147	38444 - 77 - 4	135	56030 - 56 - 9	138	74338 - 23 - 1	136	119446 - 68 - 3	123
5221 - 53 - 4	123	23505 - 41 - 1	141	38444 - 78 - 9	135	56070 - 16 - 7	144	74338 - 24 - 2	136	120068 - 37 - 3	127
5234 - 68 - 4	116	23560 - 59 - 0	129	38444 - 81 - 4	135	56073 - 10 - 0	114	74472 - 33 - 6	136	120923 - 37 - 7	111
5393 - 75 - 9	147	23564 - 05 - 8	145	38444 - 84 - 7	135	56425 - 91 - 3	128	74472 - 34 - 7	136	120928 - 09 - 8	126
5445 - 17 - 0	131	23564 - 06 - 9	145	38444 - 85 - 8	135	56558 - 16 - 8	137	74472 - 35 - 8	137	121552 - 61 - 2	119
5522 - 43 - 0	133	23950 - 58 - 5	142	38444 - 86 - 9	135	56558 - 17 - 9	137	74472 - 36 - 9	137	122453 - 73 - 0	116
5598 - 13 - 0	118	24017 - 47 - 8	146	38444 - 87 - 0	135	56558 - 18 - 0	137	74472 - 37 - 0	137	122931 - 48 - 0	143
5836 - 10 - 2	118	24579 - 73 - 5	142	38444 - 88 - 1	135	57018 - 04 - 9	145	74472 - 38 - 1	137	123312 - 89 - 0	142
5902 - 51 - 2	144	24934 - 91 - 6	117	38444 - 90 - 5	135	57153 - 17 - 0	123	74472 - 39 - 2	137	124495 - 18 - 7	143
5915 - 41 - 3	144	25013 - 16 - 5	115	38444 - 93 - 8	135	57465 - 28 - 8	137	74472 - 40 - 5	138	125116 - 23 - 6	131
6027 - 93 - 4	125	25057 - 89 -									

Reference Number Index

Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page
F061935	5, 18, 74	F110231	18, 73	F113281	28, 94	F114821	39, 140	F116281	51, 107
F032345	76	F110241	5, 18, 74	F113291	27, 125	F114831	39, 121	F116291	52, 68
F061945	7, 75	F110261	18, 74	F113301	27, 125	F114841	40, 104	F116301	52, 68
F061955	7, 73	F110421	19, 76	F113311	28, 94	F114851	40, 104	F116311	52, 131
F061965	8, 75	F110431	19, 77	F113321	29, 124	F114861	40, 105	F116321	52
F061975	8, 75	F111801	20	F113331	29, 147	F114871	40, 106	F116331	52
F061985	8, 75	F111805	19	F113341	29, 81	F114881	40, 69	F116341	52
F061995	8, 75	F111806	19	F113351	29, 81	F114911	41, 91	F116361	53, 106
F062005	8	F111807	20	F113361	29, 81	F114921	41, 54, 146	F116421	53, 70
F062015	8, 73	F111851	20	F113371	29, 81	F114931	40, 70	F116441	53, 70
F062025	9, 76	F111861	20	F113381	29, 81	F114941	40, 70	F116461	53, 111
F062035	11, 69	F111871	20	F113391	29, 81	F114951	40, 53, 70	F116471	54, 111
F062045	11	F111881	20	F113401	29, 81	F114961	40, 53, 70	F116481	54, 141
F062055	11, 69	F111891	20	F113411	28, 143	F114971	41	F116491	54, 111
F062065	11, 105	F111901	20	F113421	29, 146	F115021	41, 94	F116501	54, 91
F062075	11, 68	F111911	19	F113431	29, 81	F115031	41, 94	F116511	54, 91
F062085	12, 106	F111921	21	F113441	30, 84	F115041	41	F116521	54, 91
F062095	12, 106	F111931	21, 75	F113451	30, 84	F115051	42, 84	F116531	54, 121
F062105	12	F111941	21, 75	F113461	30, 83	F115061	41, 146	F116541	55, 93
F062125	13, 96	F111951	21, 74	F113481	30, 83	F115071	41, 116	F116551	54
F062135	14	F111991	21, 75	F113491	30, 146	F115081	42, 57, 112	F116561	55, 92
F062145	16	F112001	21, 75	F113495	36, 146	F115091	42, 57, 112	F116571	54, 92
F062155	12, 17, 68	F112011	21, 73	F113531	30, 82	F115121	42, 57, 112	F116591	54, 128
F062175	20	F112021	22	F113541	31, 84	F115131	42, 57, 112	F116601	54, 140
F062185	20	F112071	22	F113551	31, 82	F115141	42, 57, 112	F116611	54, 140
F062195	20	F112121	22, 75	F113561	31, 79	F115151	42, 57, 112	F116621	55, 94
F062205	20	F112131	22, 76	F113571	31, 79	F115161	42, 57, 112	F116631	55, 94
F062215	20	F112141	22	F113581	31, 87	F115171	42, 82	F116651	55, 93
F062225	20	F112151	23	F113591	31, 143	F115181	42, 83	F116661	55, 114
F062235	22, 75	F112161	23	F113601	31, 143	F115191	42, 83	F116671	55, 94
F062245	22, 76	F112171	18, 73	F113611	31, 143	F115201	42, 83	F116681	55, 94
F062255	22	F112181	98	F113631	31, 112	F115211	42, 83	F116691	55, 77
F062265	23	F112191	98	F113641	31, 112	F115221	42	F116701	56, 84
F062275	23	F112201	100	F113651	31, 120	F115511	43, 99	F116711	56, 83
F062285	28, 99	F112221	74	F113661	32, 140	F115521	43, 99	F116721	55, 56, 116
F062305	74	F112231	74	F113671	32, 140	F115531	43, 100	F116731	55, 56, 146
F062325	76	F112241	74	F113681	32, 140	F115561	43, 71	F116741	55, 112
F062335	76	F112261	109	F113691	31, 120	F115571	43, 68	F116751	55, 112
F062355	88	F112281	104	F113711	32, 79	F115581	43, 68	F116761	55, 112
F062365	92	F112291	76	F113721	33, 80	F115591	44, 86	F116771	55, 112
F062395	111	F112321	76	F113751	33, 80	F115601	44, 85	F116781	55, 112
F062405	112	F112331	76	F113761	33, 122	F115611	44, 86	F116791	55, 112
F062415	112	F112341	103	F113771	33, 122	F115621	44, 79	F116801	55, 112
F062425	113	F112351	103	F113781	32, 121	F115631	44, 79	F116811	56, 83
F062435	122	F112361	103	F113791	32, 79	F115641	45, 84	F116821	55, 82
F062445	126	F112371	74	F113801	33, 90	F115651	44, 83	F116831	56
F062455	132	F112381	94	F113811	32, 79	F115661	44, 116	F117221	56, 84
F062465	132	F112391	94	F113821	33, 79	F115671	44, 146	F117231	56, 83
F062475	132	F112411	91	F113831	32, 78	F115681	45, 112	F117241	56
F062485	140	F112421	92	F113841	32, 78	F115691	45, 112	F117251	56, 143
F062495	140	F112431	92	F113851	32, 78	F115701	45, 112	F117261	56, 115
F062505	140	F112441	28, 99	F113861	32, 78	F115711	45, 112	F117291	57, 97
F062515	140	F112451	99	F113871	33, 89	F115721	45, 112	F117301	57, 96
F062525	140	F112461	99	F113881	33, 89	F115731	45, 112	F117311	57, 140
F062535	141	F112471	98	F113891	33, 122	F115741	45, 112	F117321	57, 145
F062545	141	F112481	98	F113901	33, 122	F115751	45, 80	F117331	57
F062555	143	F112491	99	F113911	32, 121	F115761	45, 87	F117371	57, 77
F062565	143	F112501	100	F113921	34, 115	F115771	45, 86	F117381	58
F109230	6	F112631	95	F113931	34	F115781	45, 85	F117491	57, 129
F109250	6, 97	F112641	98	F113941	34	F115791	45, 86	F117501	57, 117
F109260	6, 98	F112651	95	F113951	34, 107	F115801	45, 132	F117511	58, 99
F109270	6, 93	F112661	95	F113961	34, 121	F115811	46	F117521	58, 98
F109280	6, 121	F112671	95	F113971	36, 99	F115821	46, 108	F117531	43, 58, 65, 98
F109300	7, 92	F112711	24, 106	F113981	36, 99	F115831	46, 108	F117541	58
F109310	7, 121	F112721	24, 106	F113991	35, 84	F115841	46, 105	F117551	58, 71
F109631	7, 75	F112731	24, 106	F114001	35, 84	F115851	46, 105	F117561	58, 120
F109641	8, 75	F112741	24, 106	F114011	35, 93	F115861	46, 105	F117571	59
F109651	8, 74	F112751	26, 35, 53, 110	F114021	35, 93	F115871	46, 69	F117581	58, 146
F109661	8, 75	F112761	22, 25, 52, 105	F114031	36, 96	F115881	47, 91	F117591	59
F109671	8, 75	F112771	22, 25, 105	F114041	36, 96	F115891	46, 87	F117611	59
F109681	7, 73	F112781	26, 35, 53, 109	F114051	36	F115901	47, 105	F117621	59
F109741	8, 74	F112901	22, 25, 34, 53, 107	F114071	35, 119	F115911	47, 87	F117631	59, 77
F109751	8, 74	F112911	22, 25, 34, 53, 107	F114081	35	F115921	47, 78	F119021	60, 86
F109761	8, 143	F112921	25, 34, 127	F114091	35, 142	F115931	47, 80	F119031	59, 86
F109771	8, 73	F112931	25, 114	F114101	35, 144	F115941	47, 78	F119041	60, 86
F109781	8, 16, 73	F112941	25, 121	F114111	36, 97	F115951	48, 78	F119051	60, 86
F109791	8, 72	F112951	25, 117	F114121	36, 98	F115961	48, 116	F119061	59, 115
F109801	8, 75	F112961	25, 34, 105	F114131	36, 98	F115971	48, 143	F119071	59
F109811	8, 75	F112971	25, 34, 105	F114141	36	F115981	48, 113	F119081	59, 79
F109821	8, 73	F112981	25, 34, 105	F114151	36	F115991	48, 134	F119091	60
F109831	9, 96	F112991	25, 33, 104	F114161	37, 78	F116001	48, 113	F119101	60, 117
F109841	9, 96	F113001	24, 27, 69	F114171	37, 114	F116011	48	F119110	60, 82
F109851	9, 76	F113011	24, 68	F114181	37, 77	F116021	49, 114	F119120	60, 82
F109861	10, 43, 65, 99	F113021	24, 69	F114201	37, 128	F116031	49, 141	F119130	61, 122
F109871	10, 99	F113031	24, 69	F114211	37	F116041	49, 113	F119140	61, 122
F109881	11, 90	F113041	26, 68	F114221	37, 99	F116051	49, 141	F119150	60, 121
F109891	11, 146	F113051	24, 69	F114231	37, 99	F116061	49, 143	F119161	60, 79
F109901	11, 91	F113061	26, 40, 52, 147	F114241	38, 72	F116071	49	F119171	60, 89
F109911	11, 69	F113071	26, 40, 52, 147	F114251	38, 72	F116081	49, 93	F119181	60, 82
F109921	11	F113081	26, 68, 108	F114261	38, 72	F116091	50, 106	F119191	61, 108
F109931	11, 69	F113101	27, 69	F114271	38, 72	F116101	50, 106	F119201	61, 108
F109941	11, 105	F113111	27	F114281	38, 72	F116111	50, 90	F119211	61, 107
F109951	11, 68	F113121	27	F114291	38, 72	F116121	50, 90	F119221	61, 105
F109961	12, 106	F113131	27	F114301	38, 121	F116131	50, 90	F119231	61, 62, 106
F109971	12, 106	F113141	27	F114311	38, 121	F116141	50, 90	F119241	62, 109
F109981	12	F113151	28, 87	F114321	38, 121	F116151	49, 90	F119251	62, 107
F109991	12, 17, 68	F113161	28, 82	F114331	38, 121	F116161	49, 90	F119261	62, 107
F110001	12, 97	F113171	27, 80, 88	F114341	38, 70	F116171	50	F119271	62, 107
F110011	12, 97	F113181	27, 80, 88	F114351	39, 71	F116181	51, 108	F119281	62, 107
F110031	13, 96	F113191	28, 116	F114361	39, 71	F116201	51, 108	F119291	62
F110041	13, 97	F113211	28, 112	F114371	39, 71	F116211	50, 117	F119311	64, 93
F110051	5, 73, 18	F113221	28, 112	F114381	39, 70	F116221	50, 52, 117	F119321	64, 93
F110061	5, 73	F113231	28, 112	F114391	39, 121	F116231	51, 52, 118	F119331	63, 77
F110071	5, 73	F113241	28, 112	F114401	39, 121	F116241	51, 106	F119341	63, 77
F110101	16	F113251	28, 112	F114411	39, 147	F116251	51, 106	F119351	63, 85
F110211	18	F113261	28, 112	F114421	39, 147	F116261	51, 106	F119361	63, 88
F110221	18	F113271	28, 31, 112	F114431	39, 147	F116271	51, 106		

Reference Number Index

Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page
F119371	63. 87	F129041	13. 89	P803280	111	P805840	113	P807900	115
F119381	64. 88	F129051	13. 91	P803290	111	P805850	113	P807910	115
F119721	63. 81	F129061	13. 91	P803320	111	P805920	113	P807940	115
F119731	63. 112	F129181	7. 92	P803330	111	P805930	113	P807950	115
F119741	63. 112	F129191	14. 94	P803360	111	P806000	113	P807980	115
F119751	63. 112	F129201	14. 133	P803370	111	P806010	113	P807990	115
F119761	63. 112	F129211	14. 134	P803380	111	P806020	113	P808020	115
F119771	63. 112	F129221	16. 100	P803390	111	P806030	113	P808030	115
F119781	63. 112	F129231	10. 75	P803450	111	P806090	113	P808040	115
F119791	63. 112	F129500	14. 72	P803460	111	P806100	113	P808050	115
F119801	63. 116	F129510	14. 72	P803470	111	P806130	113	P808060	115
F119811	63. 146	F131321	15. 92	P803480	111	P806140	113	P808070	115
F119821	62	F131331	15. 121	P803550	111	P806170	113	P808140	115
F119831	62	F131341	15. 121	P803560	111	P806180	113	P808150	115
F119841	64. 98	F131351	15. 146	P803570	111	P806210	113	P808160	115
F119851	64. 97	F131361	15. 146	P803580	111	P806250	113	P808170	115
F119861	64	F131371	16. 86	P803650	111	P806280	113	P808200	115
F119871	65. 97	F131381	14. 99	P803660	111	P806290	113	P808210	115
F119901	65. 119	F131391	17. 96	P803670	111	P806320	113	P808220	115
F119911	66. 71	F131401	17. 135	P803680	111	P806330	113	P808230	115
F119921	66. 78	F131411	17. 135	P803750	111	P806340	113	P808240	115
F119931	66. 87	F131421	17. 138	P803760	111	P806350	113	P808250	115
F119941	66. 79	F131431	17. 138	P803770	111	P806360	113	P808360	115
F119951	66. 78	F131441	17. 139	P803780	111	P806370	113	P808370	115
F119961	67. 77	F131451	17. 139	P803810	111	P806400	113	P808380	115
F119971	67. 120	F131461	17. 140	P803820	111	P806410	113	P808390	115
F119981	67	F131470	15. 95	P803850	112	P806420	113	P808420	115
F119991	67. 71	F131480	15. 135	P803860	112	P806450	114	P808430	115
F120001	65. 71	F131500	15. 140	P803890	112	P806460	114	P808460	115
F120008	65. 70	F131520	15. 96	P803900	112	P806470	114	P808470	115
F120011	65. 70	F131530	100	P804000	112	P806480	114	P808480	115
F120031	27. 146	F131540	101	P804010	112	P806510	114	P808490	115
F120051	30. 116	F131550	100	P804120	131	P806520	114	P808500	115
F120061	30. 112	F131670	101	P804130	131	P806550	114	P808510	115
F120071	30. 112	F131680	101	P804160	112	P806560	114	P808540	115
F120081	30. 112	F131690	101	P804170	112	P806570	114	P808550	115
F120091	30. 112	F131700	101	P804180	112	P806580	114	P808560	115
F120101	30. 112	F131710	101	P804190	112	P806640	114	P808570	115
F120111	30. 112	F131720	101	P804230	112	P806650	114	P808580	115
F120121	30. 112	F131740	101	P804260	112	P806680	114	P808590	115
F120131	30. 143	F131750	101	P804270	112	P806690	114	P808620	115
F120141	31. 140	F131760	101	P804280	112	P806700	114	P808630	115
F122271	75	F131770	101	P804290	112	P806710	114	P808660	115
F127591	82	F131780	101	P804300	112	P806720	114	P808670	115
F127621	89	F145091	8. 75	P804310	112	P806730	114	P808700	115
F127651	80	F146591	12. 97	P804340	112	P806760	114	P808710	115
F127661	89	F146761	47. 128	P804350	112	P806790	114	P808740	115
F127671	87	F146771	47. 112	P804380	112	P806800	114	P808750	115
F127681	85	F146781	47. 140	P804390	112	P806810	114	P808780	115
F127691	103	F146943	47. 112	P804420	112	P806820	114	P808790	115
F127701	110	F147410	5. 85	P804430	112	P806850	114	P808820	115
F127721	103	F147420	31. 120	P804460	112	P806860	114	P808830	115
F127731	107	F147430	61	P804480	112	P806890	114	P808860	115
F127751	110	F147440	61	P804520	112	P806900	114	P808870	115
F127821	95	F205024	102	P804540	112	P806910	114	P808900	115
F127871	95	F205084	102	P804580	112	P806920	114	P808910	115
F127881	95	F205094	102	P804600	112	P806950	114	P808920	115
F127891	97	F205174	102	P804660	135	P806960	114	P808930	115
F127901	15. 96. 17	F205184	102	P804700	112	P806990	114	P808960	115
F127921	15. 96	F205254	102	P804720	112	P807000	114	P808970	115
F128531	95	F205264	102	P804760	112	P807010	114	P809000	115
F128551	96	F205324	102	P804780	112	P807020	114	P809010	115
F128561	94	F205334	102	P804820	112	P807030	114	P809040	115
F128571	94	F205404	103	P804840	112	P807040	114	P809050	115
F128581	95	F205454	103	P804880	112	P807050	114	P809080	115
F128591	97	F208910	96	P804890	112	P807060	114	P809090	115
F128601	99	LM12085	20	P804920	112	P807130	114	P809240	115
F128611	99	LM12095	20	P804930	112	P807140	114	P809250	115
F128631	71	P802380	111	P804960	113	P807150	114	P809320	115
F128641	71	P802390	111	P804970	113	P807160	114	P809330	115
F128651	70	P802400	111	P804980	113	P807230	111	P809340	115
F128661	70	P802410	111	P804990	113	P807240	111	P809350	115
F128671	69	P802420	111	P805000	113	P807250	111	P809410	115
F128681	71	P802430	111	P805010	113	P807260	111	P809420	115
F128691	70	P802460	111	P805050	113	P807290	141	P809520	116
F128701	83	P802470	111	P805200	113	P807300	141	P809530	116
F128711	85	P802500	111	P805210	113	P807330	114	P809540	116
F128731	85	P802510	111	P805220	113	P807340	114	P809560	116
F128741	87	P802520	111	P805230	113	P807370	114	P809570	116
F128751	85	P802530	111	P805290	113	P807380	114	P809600	116
F128761	89	P802540	111	P805300	113	P807390	114	P809610	116
F128771	88	P802550	111	P805310	113	P807400	114	P809620	116
F128781	80	P802580	111	P805320	113	P807410	114	P809630	116
F128791	90	P802590	111	P805380	113	P807420	114	P809790	116
F128801	77	P802620	111	P805390	113	P807450	114	P809800	116
F128811	91	P802630	111	P805400	113	P807460	114	P809810	116
F128821	90	P802640	111	P805410	113	P807470	114	P809820	116
F128831	88	P802650	111	P805440	113	P807480	114	P809830	116
F128841	88	P802760	111	P805450	113	P807510	114	P809840	116
F128851	87	P802770	111	P805520	113	P807520	114	P809910	116
F128861	88	P802780	111	P805530	113	P807550	114	P809920	116
F128871	104	P802790	111	P805540	113	P807560	114	P809960	144
F128881	110	P802800	111	P805550	113	P807590	114	P809970	144
F128891	110	P802810	111	P805560	113	P807600	114	P810000	116
F128901	110	P802920	111	P805570	113	P807630	114	P810010	116
F128911	110	P802930	111	P805600	113	P807640	114	P810020	116
F128921	110	P802960	111	P805610	113	P807670	114	P810030	116
F128931	110	P802970	111	P805620	113	P807680	114	P810130	116
F128941	23. 76	P803000	111	P805630	113	P807710	114	P810140	116
F128951	16. 86	P803010	111	P805640	113	P807720	114	P810150	116
F128961	9. 104	P803160	111	P805650	113	P807750	114	P810160	116
F128971	10. 81	P803170	111	P805680	113	P807760	114	P810190	116
F128981	9. 118	P803180	111	P805690	113	P807800	114	P810200	116
F128991	9. 117	P803190	111	P805700	113	P807810	114	P810210	116
F129001	10. 75	P803200	111	P805710	113	P807840	115	P810230	116
F129011	10. 145	P803210	111	P805780	113	P807850	115	P810250	116
F129021	10. 80	P803240	111	P805790	113	P807860	115	P810290	116
F129031	12. 89	P803250	111	P805800	113	P807870	115	P810300	116

Reference Number Index

Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page
P810310	116	P812450	118	P815320	118	P817630	121	P819660	122
P810320	116	P812480	118	P815350	118	P817660	121	P819670	122
P810330	116	P812490	118	P815360	118	P817670	121	P819740	123
P810340	116	P812500	118	P815370	118	P817680	121	P819750	123
P810410	116	P812510	118	P815380	118	P817690	121	P819760	123
P810420	116	P812540	118	P815410	122	P817720	121	P819770	123
P810430	116	P812550	118	P815420	122	P817730	121	P819830	123
P810440	116	P812560	118	P815430	122	P817800	121	P819840	123
P810450	116	P812570	118	P815440	122	P817810	121	P819870	123
P810460	116	P812580	118	P815530	120	P817840	121	P819880	123
P810490	116	P812590	118	P815540	120	P817850	121	P819890	123
P810500	116	P812620	118	P815550	120	P817880	121	P819900	123
P810510	116	P812630	118	P815560	120	P817890	121	P819930	123
P810520	116	P812640	118	P815590	120	P817920	123	P819940	123
P810580	116	P812650	118	P815600	120	P817930	123	P819970	123
P810590	116	P812660	118	P815610	120	P817960	121	P819980	123
P810620	117	P812670	118	P815620	120	P817970	121	P820010	123
P810630	117	P812700	118	P815690	120	P817980	121	P820020	123
P810660	117	P812710	118	P815700	120	P817990	121	P820030	123
P810670	117	P812740	118	P815710	120	P818020	121	P820040	123
P810680	117	P812750	118	P815720	120	P818030	121	P820050	123
P810690	117	P812760	118	P815750	120	P818040	121	P820060	123
P810770	117	P812770	118	P815760	120	P818050	121	P820090	123
P810800	117	P812920	118	P815770	120	P818060	121	P820100	123
P810810	117	P812930	118	P815780	120	P818070	121	P820130	123
P810820	117	P812940	118	P815850	120	P818100	121	P820140	123
P810830	117	P812950	118	P815860	120	P818130	121	P820150	123
P810860	117	P813010	118	P815870	120	P818140	121	P820160	123
P810870	117	P813020	118	P815880	120	P818150	121	P820190	123
P810880	117	P813030	118	P815940	120	P818160	121	P820200	123
P810890	117	P813040	118	P815950	120	P818190	121	P820230	123
P810900	117	P813050	118	P815960	120	P818200	121	P820240	123
P810910	117	P813060	118	P815970	120	P818210	121	P820270	123
P810940	117	P813070	118	P816030	119	P818220	121	P820280	123
P810950	117	P813080	118	P816040	119	P818230	121	P820310	123
P810960	117	P813110	118	P816050	119	P818240	121	P820320	123
P810970	117	P813120	118	P816060	119	P818270	121	P820350	141
P810980	117	P813190	118	P816090	119	P818280	121	P820360	141
P810990	117	P813200	118	P816100	119	P818290	121	P820390	123
P811020	117	P813210	118	P816110	119	P818300	121	P820400	123
P811030	117	P813220	118	P816120	119	P818310	121	P820410	123
P811060	117	P813330	118	P816150	119	P818320	121	P820420	123
P811090	117	P813340	118	P816160	119	P818350	121	P820430	123
P811100	117	P813350	118	P816170	119	P818360	121	P820440	123
P811110	117	P813360	118	P816180	119	P818440	121	P820470	123
P811120	117	P813390	118	P816270	128	P818450	121	P820480	123
P811210	117	P813400	118	P816280	128	P818480	121	P820490	123
P811220	117	P813410	118	P816310	120	P818490	121	P820500	123
P811230	117	P813420	118	P816320	120	P818560	122	P820530	123
P811240	117	P813610	119	P816330	120	P818570	122	P820540	123
P811270	117	P813620	119	P816340	120	P818600	122	P820570	123
P811280	117	P813630	119	P816410	120	P818610	122	P820580	123
P811310	117	P813640	119	P816420	120	P818640	122	P820650	123
P811320	117	P813750	119	P816430	120	P818670	122	P820660	123
P811350	117	P813760	119	P816440	120	P818680	122	P820730	123
P811430	117	P813910	119	P816630	120	P818710	122	P820740	123
P811440	117	P813920	119	P816640	120	P818720	122	P820800	123
P811470	117	P813930	119	P816650	120	P818750	122	P820810	123
P811480	117	P813940	119	P816660	120	P818760	122	P820820	123
P811510	117	P813950	119	P816810	120	P818790	122	P820830	123
P811520	117	P813960	119	P816820	120	P818800	122	P820860	123
P811550	117	P814070	119	P816830	120	P818830	122	P820870	123
P811560	117	P814080	119	P816840	120	P818840	122	P820900	123
P811590	117	P814110	119	P816910	120	P818850	122	P820910	123
P811600	117	P814120	119	P816920	120	P818860	122	P820940	123
P811610	117	P814130	119	P816930	120	P818890	122	P820950	123
P811620	117	P814140	119	P816940	120	P818900	122	P820980	123
P811630	117	P814210	119	P816950	120	P818930	122	P820990	123
P811640	117	P814220	119	P816960	120	P818940	122	P821000	123
P811670	117	P814230	119	P816990	120	P818970	122	P821010	123
P811680	117	P814240	119	P817020	120	P818980	122	P821020	123
P811710	117	P814270	130	P817030	120	P818990	122	P821030	123
P811720	117	P814280	130	P817060	120	P819000	122	P821090	123
P811750	117	P814310	119	P817070	120	P819030	122	P821100	123
P811760	117	P814320	119	P817080	120	P819040	122	P821110	123
P811770	117	P814350	119	P817090	120	P819070	122	P821120	123
P811780	117	P814360	119	P817120	120	P819100	122	P821130	123
P811810	117	P814390	119	P817130	120	P819110	122	P821140	123
P811820	117	P814400	119	P817160	120	P819120	122	P821170	123
P811850	117	P814480	119	P817170	120	P819130	122	P821180	123
P811860	117	P814490	119	P817180	120	P819200	122	P821190	123
P811870	117	P814500	119	P817190	120	P819210	122	P821200	123
P811880	117	P814510	119	P817220	120	P819220	122	P821230	123
P811910	117	P814520	119	P817230	120	P819230	122	P821240	123
P811920	117	P814530	119	P817260	121	P819260	122	P821270	124
P811950	117	P814560	119	P817270	121	P819270	122	P821280	124
P811960	117	P814570	119	P817280	121	P819280	122	P821290	124
P812020	118	P814580	119	P817290	121	P819290	122	P821300	124
P812030	118	P814590	119	P817320	121	P819320	122	P821330	124
P812060	118	P814660	119	P817330	121	P819330	122	P821340	124
P812070	118	P814670	119	P817340	121	P819340	122	P821350	124
P812160	118	P814850	119	P817350	121	P819350	122	P821360	124
P812170	118	P814860	119	P817380	121	P819380	122	P821390	124
P812180	118	P814870	119	P817390	121	P819390	122	P821420	124
P812190	118	P814880	119	P817420	121	P819400	122	P821430	124
P812220	118	P814890	119	P817430	121	P819410	122	P821460	124
P812230	118	P814900	119	P817440	121	P819440	122	P821470	124
P812240	118	P814920	119	P817450	121	P819450	122	P821500	124
P812250	118	P814930	119	P817460	121	P819480	122	P821510	124
P812320	118	P814940	119	P817500	121	P819490	122	P821540	124
P812330	118	P814950	119	P817510	121	P819520	122	P821550	124
P812360	118	P815050	119	P817540	116	P819530	122	P821560	124
P812370	118	P815060	119	P817550	116	P819560	122	P821570	124
P812380	118	P815070	119	P817560	116	P819570	122	P821600	124
P812390	118	P815080	119	P817570	116	P819600	122	P821610	124
P812420	118	P815290	118	P817580	116	P819610	122	P821640	124
P812430	118	P815300	118	P817590	116	P819640	122	P821650	124
P812440	118	P815310	118	P817620	121	P819650	122	P821680	124

Reference Number Index

Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page
P821690	124	P823730	125	P826040	127	P828260	129	P830440	130
P821720	124	P823740	125	P826050	127	P828270	129	P830450	130
P821730	124	P823780	125	P826060	127	P828300	129	P830560	130
P821740	124	P823790	125	P826150	127	P828310	129	P830570	130
P821750	124	P823820	125	P826160	127	P828340	129	P830580	130
P821780	124	P823830	125	P826170	127	P828350	129	P830590	130
P821790	124	P823900	125	P826180	127	P828360	129	P830750	130
P821800	124	P823910	125	P826190	127	P828370	129	P830760	130
P821810	124	P823920	125	P826200	127	P828460	129	P830770	130
P821860	124	P823930	125	P826230	127	P828470	129	P830780	130
P821870	124	P823960	125	P826240	127	P828480	129	P830840	130
P821900	124	P823970	125	P826250	127	P828490	129	P830850	130
P821910	124	P823980	125	P826280	127	P828520	129	P830860	130
P821920	124	P823990	125	P826290	127	P828530	129	P830870	130
P821930	124	P824060	125	P826300	127	P828540	129	P830940	131
P821960	124	P824070	125	P826310	127	P828550	129	P830950	131
P821970	124	P824100	125	P826320	127	P828580	129	P830960	131
P822000	124	P824110	125	P826330	127	P828590	129	P830970	131
P822010	124	P824120	125	P826360	127	P828600	129	P830980	131
P822020	124	P824130	125	P826390	127	P828610	129	P830990	131
P822030	124	P824140	125	P826400	127	P828640	129	P831020	131
P822040	124	P824150	125	P826430	127	P828650	129	P831030	131
P822050	124	P824220	125	P826440	127	P828660	129	P831180	131
P822120	124	P824230	125	P826470	127	P828670	129	P831190	131
P822130	124	P824240	125	P826480	127	P828680	129	P831200	131
P822140	124	P824250	125	P826510	128	P828690	129	P831210	131
P822150	124	P824280	126	P826520	128	P828720	129	P831220	131
P822180	124	P824290	126	P826590	128	P828730	129	P831230	131
P822190	124	P824300	126	P826600	128	P828760	129	P831260	131
P822200	124	P824310	126	P826630	128	P828770	129	P831270	131
P822210	124	P824340	126	P826640	128	P828780	129	P831340	131
P822240	141	P824350	126	P826650	128	P828790	129	P831350	131
P822250	141	P824360	126	P826660	128	P828820	129	P831380	131
P822280	141	P824370	126	P826670	128	P828830	129	P831390	131
P822290	141	P824400	126	P826680	128	P828890	129	P831400	131
P822320	124	P824410	126	P826710	128	P828900	129	P831410	131
P822330	124	P824440	126	P826720	128	P828910	129	P831480	131
P822340	124	P824450	126	P826750	128	P828920	129	P831490	131
P822350	124	P824460	126	P826760	128	P828950	129	P831500	131
P822390	124	P824470	126	P826800	128	P828960	129	P831510	131
P822420	124	P824530	126	P826850	128	P828970	129	P831540	131
P822470	124	P824600	126	P826860	128	P828980	129	P831550	131
P822500	124	P824610	126	P826870	128	P829090	129	P831560	131
P822510	124	P824640	126	P826880	128	P829100	129	P831570	131
P822520	124	P824650	126	P826910	128	P829110	129	P831600	131
P822530	124	P824660	126	P826920	128	P829120	129	P831610	131
P822560	124	P824690	131	P826980	128	P829130	129	P831620	131
P822570	124	P824700	131	P826990	128	P829140	129	P831630	131
P822600	124	P824710	131	P827000	128	P829170	129	P831700	131
P822610	124	P824720	131	P827010	128	P829180	129	P831710	131
P822640	124	P824750	126	P827040	128	P829240	129	P831740	131
P822650	124	P824760	126	P827050	128	P829250	129	P831750	131
P822720	124	P824800	126	P827080	128	P829260	129	P831780	131
P822730	124	P824810	126	P827090	128	P829270	129	P831790	131
P822740	124	P824820	126	P827120	128	P829280	129	P831800	131
P822750	124	P824850	126	P827130	128	P829290	129	P831810	131
P822820	125	P824860	126	P827140	128	P829350	130	P831880	131
P822830	125	P824930	126	P827150	128	P829360	130	P831890	131
P822840	125	P824940	126	P827180	128	P829390	130	P831900	131
P822850	125	P824970	126	P827190	128	P829400	130	P831910	131
P822860	125	P824980	126	P827200	128	P829430	130	P832020	131
P822870	125	P824990	126	P827210	128	P829440	130	P832030	131
P822880	125	P825000	126	P827220	128	P829470	130	P832060	131
P822890	125	P825070	126	P827230	128	P829480	130	P832070	131
P822920	125	P825080	126	P827290	128	P829490	130	P832140	131
P822930	125	P825090	126	P827390	128	P829500	130	P832150	131
P823060	125	P825100	126	P827400	128	P829530	130	P832160	131
P823070	125	P825110	126	P827410	128	P829540	130	P832170	131
P823080	125	P825120	126	P827420	128	P829550	130	P832180	131
P823090	125	P825190	126	P827480	128	P829560	130	P832190	131
P823100	125	P825200	126	P827510	128	P829630	130	P832200	131
P823110	125	P825210	126	P827540	128	P829640	130	P832210	131
P823120	125	P825220	126	P827550	128	P829650	130	P832310	131
P823130	125	P825290	127	P827560	128	P829660	130	P832320	131
P823200	125	P825300	127	P827570	128	P829690	130	P832350	131
P823210	125	P825330	127	P827580	128	P829700	130	P832400	116
P823220	125	P825340	127	P827590	128	P829710	130	P832470	131
P823230	125	P825370	127	P827620	128	P829720	130	P832500	131
P823240	125	P825380	127	P827690	128	P829730	130	P832510	131
P823250	125	P825390	127	P827700	128	P829740	130	P832520	131
P823260	125	P825400	127	P827710	128	P829770	130	P832530	131
P823270	125	P825430	127	P827720	128	P829780	130	P832670	132
P823340	125	P825440	127	P827790	128	P829810	130	P832680	132
P823350	125	P825450	127	P827800	128	P829820	130	P832710	131, 132
P823360	125	P825460	127	P827810	128	P829850	130	P832720	131, 132
P823370	125	P825490	127	P827820	128	P829860	130	P832750	132
P823380	125	P825500	127	P827830	128	P829930	130	P832760	132
P823390	125	P825510	127	P827840	128	P829940	130	P832770	132
P823420	125	P825520	127	P827910	128	P829970	130	P832780	132
P823430	125	P825580	127	P827920	128	P829980	130	P832790	132
P823440	125	P825590	127	P827930	128	P829990	130	P832800	132
P823450	125	P825620	127	P827940	128	P830000	130	P832830	132
P823480	125	P825630	127	P827970	128	P830030	130	P832840	132
P823490	125	P825700	127	P827980	128	P830040	130	P832850	132
P823500	125	P825710	127	P827990	128	P830050	130	P832860	132
P823510	125	P825720	127	P828000	128	P830060	130	P832870	132
P823520	125	P825730	127	P828100	129	P830130	130	P832880	132
P823530	125	P825840	127	P828110	129	P830140	130	P832910	112
P823540	125	P825850	127	P828140	129	P830150	130	P832920	112
P823580	125	P825910	127	P828150	129	P830160	130	P832930	112
P823590	125	P825920	127	P828180	129	P830190	130	P832940	112
P823600	125	P825930	127	P828190	129	P830200	130	P832970	132
P823610	125	P825940	127	P828200	129	P830340	130	P832980	132
P823620	125	P825970	127	P828210	129	P830350	130	P832990	132
P823630	125	P825980	127	P828220	129	P830380	130	P833000	132
P823640	125	P825990	127	P828230	129	P830390	130	P833030	132
P823650	125	P826000	127	P828240	129	P830400	130	P833040	132
P823720	125	P826030	127	P828250	129	P830410	130	P833050	132

Reference Number Index

Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page
P833060	132	P835090	133	P837690	134	P839620	136	P841710	137
P833090	132	P835120	133	P837720	134	P839630	136	P841720	137
P833100	132	P835130	133	P837730	134	P839660	136	P841750	137
P833110	132	P835140	133	P837740	134	P839670	136	P841760	137
P833120	132	P835150	133	P837750	134	P839700	136	P841790	137
P833150	132	P835180	133	P837780	135	P839710	136	P841800	137
P833160	132	P835190	133	P837790	135	P839740	136	P841810	137
P833190	132	P835200	133	P837800	135	P839750	136	P841820	137
P833200	132	P835210	133	P837810	135	P839780	136	P841850	137
P833230	132	P835220	133	P837840	135	P839790	136	P841860	137
P833240	132	P835230	133	P837850	135	P839800	136	P841890	137
P833250	132	P835260	133	P837860	135	P839810	136	P841900	137
P833260	132	P835270	133	P837870	135	P839840	136	P841930	137
P833290	132	P835280	133	P837900	135	P839850	136	P841940	137
P833300	132	P835290	133	P837910	135	P839880	136	P841970	137
P833310	132	P835320	133	P837940	135	P839890	136	P841980	137
P833320	132	P835330	133	P837950	135	P839920	136	P841990	137
P833350	132	P835360	133	P837960	135	P839930	136	P842000	137
P833360	132	P835370	133	P837970	135	P839960	136	P842030	137
P833370	132	P835380	133	P838000	135	P839970	136	P842040	137
P833380	132	P835390	133	P838010	135	P840000	136	P842070	137
P833480	132	P835420	133	P838020	135	P840010	136	P842080	137
P833490	132	P835430	133	P838030	135	P840040	136	P842110	137
P833520	132	P835440	133	P838060	135	P840050	136	P842120	137
P833530	132	P835450	133	P838070	135	P840080	136	P842150	137
P833540	132	P835510	133	P838100	135	P840090	136	P842160	137
P833550	132	P835520	133	P838110	135	P840120	136	P842190	137
P833560	132	P835530	133	P838120	135	P840130	136	P842200	137
P833570	132	P835540	133	P838130	135	P840160	136	P842230	137
P833670	132	P835570	134	P838160	135	P840170	136	P842240	137
P833680	132	P835580	134	P838170	135	P840200	136	P842270	137
P833690	132	P835640	134	P838180	135	P840210	136	P842280	137
P833700	132	P835650	134	P838190	135	P840240	136	P842310	137
P833730	132	P835680	134	P838220	135	P840250	136	P842320	137
P833740	132	P835690	134	P838230	135	P840280	136	P842350	137
P833750	132	P835720	134	P838260	135	P840290	136	P842360	137
P833760	132	P835730	134	P838270	135	P840320	136	P842390	137
P833790	132	P835790	134	P838300	135	P840330	136	P842400	137
P833800	132	P835800	134	P838310	135	P840360	136	P842430	137
P833860	132	P835830	134	P838320	135	P840370	136	P842440	137
P833870	132	P835840	134	P838330	135	P840380	136	P842470	137
P833880	132	P835850	134	P838360	135	P840390	136	P842480	137
P833890	132	P835860	134	P838370	135	P840420	136	P842510	137
P833920	132	P835890	134	P838400	135	P840430	136	P842520	137
P833930	132	P835900	134	P838410	135	P840460	136	P842530	137
P833940	132	P835910	134	P838440	135	P840470	136	P842540	137
P833950	132	P835920	134	P838450	135	P840500	136	P842570	137
P833980	132	P835930	134	P838480	135	P840510	136	P842580	137
P833990	132	P835940	134	P838490	135	P840540	136	P842610	137
P834050	132	P835950	134	P838520	135	P840550	136	P842620	137
P834060	132	P835960	134	P838530	135	P840580	136	P842650	137
P834090	132	P835990	111	P838560	135	P840590	136	P842660	137
P834100	132	P836000	111	P838570	135	P840620	136	P842690	137
P834130	132	P836010	134	P838600	135	P840630	136	P842700	137
P834140	132	P836020	134	P838610	135	P840660	136	P842730	137
P834150	132	P836080	120	P838640	135	P840670	136	P842740	137
P834160	132	P836090	120	P838650	135	P840700	136	P842770	137
P834170	132	P836100	120	P838680	135	P840710	136	P842780	137
P834180	132	P836110	120	P838690	135	P840740	136	P842810	137
P834210	132	P836140	134	P838720	135	P840750	136	P842820	137
P834220	132	P836150	134	P838730	135	P840780	136	P842850	137
P834230	132	P836160	134	P838740	135	P840790	136	P842860	137
P834260	132	P836170	134	P838750	135	P840820	136	P842870	137
P834270	132	P836180	134	P838780	135	P840830	136	P842880	137
P834300	133	P836190	134	P838790	135	P840840	136	P842910	137
P834310	133	P836220	134	P838800	135	P840850	136	P842920	137
P834320	133	P836230	134	P838810	135	P840880	136	P842950	137
P834330	133	P836380	134	P838840	135	P840890	136	P842960	137
P834360	133	P836390	134	P838850	135	P840920	136	P842970	137
P834370	133	P836560	134	P838880	135	P840930	136	P842980	137
P834400	133	P836570	134	P838890	135	P840960	136	P843010	138
P834410	133	P836580	134	P838910	135	P840970	136	P843020	138
P834420	133	P836590	134	P838940	135	P841000	136	P843050	138
P834430	133	P836660	134	P838950	135	P841010	136	P843060	138
P834460	133	P836670	134	P838980	135	P841040	136	P843090	138
P834470	133	P836680	134	P838990	135	P841050	136	P843100	138
P834480	133	P836690	134	P839020	135	P841080	136	P843130	138
P834510	133	P836760	134	P839030	135	P841090	136	P843140	138
P834520	133	P836770	134	P839060	135	P841120	136	P843170	138
P834530	133	P836840	134	P839070	135	P841130	136	P843180	138
P834540	133	P836850	134	P839080	135	P841160	136	P843210	138
P834570	133	P836920	134	P839090	135	P841170	136	P843220	138
P834580	133	P836930	134	P839120	135	P841200	136	P843250	138
P834610	133	P837320	119	P839130	135	P841210	136	P843260	138
P834620	133	P837330	119	P839160	135	P841240	137	P843290	138
P834650	133	P837340	119	P839170	135	P841250	137	P843300	138
P834660	133	P837350	119	P839200	135	P841280	137	P843330	138
P834690	133	P837380	134	P839210	135	P841290	137	P843340	138
P834700	133	P837390	134	P839240	135	P841320	137	P843370	138
P834710	133	P837400	134	P839250	135	P841330	137	P843380	138
P834720	133	P837410	134	P839280	135	P841360	137	P843390	138
P834750	133	P837440	134	P839290	135	P841370	137	P843400	138
P834760	133	P837450	134	P839320	135	P841400	137	P843430	138
P834790	133	P837460	134	P839330	135	P841410	137	P843440	138
P834820	133	P837470	134	P839360	135	P841440	137	P843470	138
P834830	133	P837500	134	P839370	135	P841450	137	P843480	138
P834860	133	P837510	134	P839400	135	P841480	137	P843510	138
P834870	133	P837520	134	P839410	135	P841490	137	P843520	138
P834900	133	P837530	134	P839440	135	P841520	137	P843550	138
P834910	133	P837560	134	P839450	135	P841530	137	P843560	138
P834940	133	P837570	134	P839460	135	P841560	137	P843590	138
P834980	133	P837580	134	P839470	135	P841570	137	P843600	138
P834990	133	P837590	134	P839500	136	P841590	137	P843630	138
P835020	133	P837620	134	P839510	136	P841600	137	P843640	138
P835030	133	P837630	134	P839540	136	P841630	137	P843670	138
P835060	133	P837640	134	P839550	136	P841640	137	P843680	138
P835070	133	P837650	134	P839580	136	P841670	137	P843710	138
P835080	133	P837680	134	P839590	136	P841680	137	P843720	138

Reference Number Index

Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page
P843750	138	P845790	139	P847900	140	P849940	142	P852030	144
P843760	138	P845800	139	P847910	140	P850010	142	P852040	144
P843790	138	P845810	139	P847940	140	P850020	142	P852050	144
P843800	138	P845820	139	P847950	140	P850030	142	P852060	144
P843830	138	P845850	139	P847960	140	P850040	142	P852070	144
P843840	138	P845860	139	P847990	141	P850050	142	P852100	144
P843850	138	P845890	139	P848000	141	P850060	142	P852110	144
P843860	138	P845900	139	P848030	141	P850130	142	P852120	144
P843890	138	P845930	139	P848040	141	P850140	142	P852130	144
P843900	138	P845940	139	P848110	141	P850150	142	P852160	144
P843930	138	P845970	139	P848120	141	P850160	142	P852170	144
P843940	138	P845980	139	P848130	141	P850190	142	P852200	144
P843970	138	P846010	139	P848140	141	P850200	142	P852210	144
P843980	138	P846020	139	P848200	141	P850230	142	P852240	144
P844010	138	P846050	139	P848210	141	P850240	142	P852250	144
P844020	138	P846060	139	P848220	141	P850250	142	P852260	144
P844030	138	P846090	139	P848230	141	P850260	142	P852270	144
P844040	138	P846100	139	P848290	141	P850270	142	P852300	144
P844070	138	P846130	139	P848300	141	P850280	142	P852310	144
P844080	138	P846140	139	P848330	141	P850290	142	P852320	144
P844110	138	P846170	139	P848340	141	P850300	142	P852330	144
P844120	138	P846180	139	P848370	141	P850330	142	P852360	144
P844150	138	P846210	139	P848380	141	P850340	142	P852370	144
P844160	138	P846220	139	P848410	141	P850350	142	P852380	144
P844190	138	P846250	139	P848420	141	P850380	142	P852390	144
P844200	138	P846260	139	P848430	141	P850390	142	P852490	145
P844230	138	P846270	139	P848440	141	P850400	142	P852500	145
P844240	138	P846280	139	P848470	141	P850410	142	P852530	144
P844270	138	P846310	139	P848480	141	P850440	142	P852540	144
P844280	138	P846320	139	P848490	141	P850450	142	P852550	144
P844310	138	P846350	139	P848500	141	P850460	142	P852560	144
P844320	138	P846360	139	P848510	141	P850470	142	P852590	144
P844350	138	P846390	140	P848520	141	P850480	142	P852600	144
P844360	138	P846400	140	P848550	141	P850490	142	P852630	144
P844390	138	P846410	140	P848560	141	P850520	142	P852640	144
P844400	138	P846420	140	P848570	141	P850530	142	P852670	144
P844430	138	P846430	140	P848580	141	P850540	142	P852680	144
P844440	138	P846440	140	P848590	141	P850550	142	P852690	144
P844470	138	P846450	140	P848600	141	P850580	143	P852700	144
P844480	138	P846470	140	P848630	141	P850590	143	P852730	144
P844510	138	P846480	140	P848640	141	P850620	143	P852740	144
P844520	138	P846510	140	P848650	141	P850630	143	P852770	144
P844550	138	P846520	140	P848660	141	P850660	143	P852780	144
P844560	138	P846530	140	P848720	141	P850670	143	P852790	144
P844590	138	P846540	140	P848730	141	P850680	143	P852800	144
P844600	138	P846550	140	P848740	141	P850690	143	P852810	144
P844630	138	P846560	140	P848750	141	P850720	143	P852820	144
P844640	138	P846600	140	P848810	141	P850730	143	P852890	144
P844670	138	P846610	140	P848820	141	P850740	143	P852900	144
P844680	138	P846620	140	P848890	141	P850750	143	P852930	144
P844710	138	P846630	140	P848900	141	P850780	143	P852940	144
P844720	138	P846690	140	P848910	141	P850790	143	P852970	144
P844730	138	P846700	140	P848920	141	P850800	143	P852980	144
P844740	138	P846710	140	P848950	141	P850810	143	P852990	144
P844770	139	P846720	140	P848960	141	P850840	143	P853000	144
P844780	139	P846790	140	P848970	141	P850850	143	P853060	144
P844810	139	P846800	140	P848980	141	P850860	143	P853070	144
P844820	139	P846810	140	P848990	141	P850870	143	P853080	144
P844850	139	P846820	140	P849000	141	P850880	143	P853090	144
P844860	139	P846890	140	P849070	141	P850890	143	P853100	144
P844890	139	P846900	140	P849080	141	P850920	143	P853110	144
P844900	139	P846910	140	P849110	141	P850930	143	P853140	144
P844930	139	P846920	140	P849120	141	P850940	143	P853150	144
P844940	139	P846930	140	P849130	141	P850950	143	P853220	144
P844970	139	P846940	140	P849140	141	P851020	143	P853230	144
P844980	139	P846970	140	P849170	142	P851030	143	P853240	144
P845010	139	P846980	140	P849180	142	P851040	143	P853250	144
P845020	139	P847010	140	P849190	142	P851050	143	P853260	144
P845050	139	P847020	140	P849200	142	P851060	143	P853270	144
P845060	139	P847050	140	P849210	142	P851070	143	P853300	144
P845090	139	P847060	140	P849220	142	P851080	143	P853310	144
P845100	139	P847160	140	P849280	142	P851090	143	P853340	144
P845130	139	P847170	140	P849290	142	P851120	143	P853350	144
P845140	139	P847200	140	P849300	142	P851130	143	P853380	144
P845150	139	P847210	140	P849310	142	P851160	143	P853390	144
P845160	139	P847240	140	P849340	142	P851170	143	P853420	145
P845190	139	P847250	140	P849350	142	P851180	143	P853430	145
P845200	139	P847280	140	P849410	142	P851190	143	P853440	145
P845230	139	P847290	140	P849420	142	P851260	143	P853450	145
P845240	139	P847380	140	P849430	142	P851290	143	P853460	145
P845270	139	P847390	140	P849440	142	P851300	143	P853470	145
P845280	139	P847400	140	P849450	142	P851310	143	P853500	144
P845310	139	P847410	140	P849510	142	P851320	143	P853510	144
P845320	139	P847480	140	P849520	142	P851330	143	P853540	145
P845350	139	P847490	140	P849530	142	P851340	143	P853550	145
P845360	139	P847500	140	P849540	142	P851400	143	P853560	145
P845390	139	P847510	140	P849550	142	P851410	143	P853570	145
P845400	139	P847540	140	P849560	142	P851420	143	P853580	145
P845430	139	P847550	140	P849570	142	P851430	143	P853590	145
P845440	139	P847560	140	P849600	142	P851560	143	P853620	145
P845450	139	P847570	140	P849610	142	P851570	143	P853630	145
P845460	139	P847630	140	P849620	142	P851580	143	P853640	145
P845490	139	P847640	140	P849630	142	P851590	143	P853650	145
P845500	139	P847650	140	P849640	142	P851600	143	P853680	145
P845530	139	P847660	140	P849650	142	P851610	143	P853690	145
P845540	139	P847690	140	P849750	142	P851700	143	P853720	145
P845570	139	P847700	140	P849760	142	P851710	143	P853730	145
P845580	139	P847730	140	P849770	142	P851740	143	P853740	145
P845610	139	P847740	140	P849780	142	P851750	143	P853750	145
P845620	139	P847750	140	P849810	142	P851900	143	P853780	145
P845650	139	P847760	140	P849820	142	P851910	143	P853790	145
P845660	139	P847770	140	P849830	142	P851920	143	P853900	145
P845690	139	P847780	140	P849840	142	P851930	143	P853910	145
P845700	139	P847810	140	P849850	142	P851960	143	P853920	145
P845730	139	P847820	140	P849860	142	P851970	143	P853930	145
P845740	139	P847830	140	P849890	142	P851980	143	P854000	145
P845750	139	P847860	140	P849900	142	P851990	143	P854010	145
P845760	139	P847870	140	P849930	142	P852020	144	P854020	145

Reference Number Index

Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page
P854030	145	P856060	147	P872380	121	P873610	129	P874730	116
P854060	145	P856070	147	P872430	117	P873620	129	P874740	116
P854070	145	P856080	147	P872440	129	P873630	129	P874750	116
P854080	145	P856090	147	P872470	131	P873640	132	P874760	117
P854090	145	P856100	147	P872480	134	P873650	118	P874770	117
P854120	145	P856110	147	P872510	144	P873660	115	P874780	117
P854130	145	P856120	147	P872530	144	P873670	117	P874790	117
P854160	145	P856130	147	P872550	144	P873680	117	P874800	118
P854170	145	P856160	147	P872560	144	P873690	131	P874820	118
P854180	145	P856170	147	P872570	144	P873700	131	P874830	118
P854190	145	P856180	147	P872580	144	P873710	131	P874840	118
P854250	145	P856190	147	P872590	146	P873710	132	P874850	118
P854260	145	P856200	147	P872600	146	P873720	131	P874860	118
P854290	145	P856210	147	P872610	146	P873730	133	P874870	118
P854300	145	P856280	147	P872620	146	P873740	133	P874880	118
P854330	145	P856290	147	P872630	147	P873750	133	P874890	118
P854340	145	P856300	147	P872640	146	P873760	133	P874900	118
P854370	145	P856310	147	P872650	146	P873770	124	P874910	118
P854380	145	P856340	147	P872670	113	P873780	124	P874920	118
P854410	146	P856350	147	P872680	146	P873790	133	P874930	118
P854430	146	P856380	147	P872690	147	P873810	143	P874940	118
P854450	146	P856390	147	P872700	147	P873820	111	P874950	119
P854490	126	P856420	147	P872710	147	P873830	111	P874960	119
P854500	126	P856430	147	P872720	121	P873840	111	P874970	119
P854510	126	P856460	147	P872730	146	P873850	111	P874980	119
P854520	126	P856470	147	P872740	147	P873860	111	P874990	119
P854530	126	P856530	147	P872750	119	P873870	111	P875000	119
P854540	126	P856540	147	P872760	119	P873880	111	P875010	119
P854550	126	P856550	147	P872770	119	P873890	111	P875020	119
P854560	126	P856560	147	P872780	120	P873900	111	P875030	119
P854590	126	P856590	147	P872790	120	P873910	111	P875040	119
P854600	126	P856600	147	P872800	121	P873920	111	P875070	119
P854630	146	P856630	147	P872810	121	P873930	111	P875080	119
P854640	146	P856640	147	P872820	122	P873940	111	P875090	119
P854670	145	P856650	147	P872830	122	P873950	111	P875100	119
P854680	145	P856660	147	P872840	122	P873960	111	P875110	119
P854750	146	P856670	147	P872850	124	P873970	111	P875120	119
P854760	146	P856680	147	P872860	121	P873980	111	P875130	119
P854770	146	P856710	147	P872870	122	P874000	119	P875140	119
P854780	146	P856720	147	P872880	124	P874010	119	P875150	120
P854810	146	P856790	147	P872890	121	P874020	112	P875160	120
P854820	146	P856800	147	P872900	121	P874030	112	P875170	128
P854830	146	P856830	147	P872910	121	P874040	112	P875180	128
P854840	146	P856840	147	P872920	122	P874050	113	P875190	120
P854870	146	P856870	147	P872930	122	P874060	113	P875200	120
P854880	146	P856880	147	P872940	122	P874070	113	P875210	120
P854890	146	P856910	147	P872950	123	P874080	113	P875220	120
P854900	146	P856920	147	P872980	112	P874090	113	P875230	120
P854910	146	P856930	147	P872990	112	P874100	113	P875240	120
P854920	146	P856940	147	P873000	112	P874110	113	P875250	120
P854930	146	P857170	147	P873010	112	P874120	113	P875260	120
P854940	146	P857180	147	P873020	117	P874130	113	P875270	120
P854970	146	P857190	147	P873030	117	P874140	113	P875280	120
P854980	146	P857200	147	P873040	117	P874150	113	P875290	120
P854990	146	P857230	147	P873050	118	P874160	113	P875300	128
P855000	146	P857240	147	P873060	126	P874170	113	P875310	120
P855030	146	P857250	147	P873070	126	P874180	113	P875320	120
P855040	146	P857260	147	P873080	126	P874190	113	P875330	120
P855070	146	P857330	147	P873090	129	P874200	113	P875340	120
P855080	146	P857340	147	P873100	129	P874210	113	P875350	120
P855090	146	P857370	147	P873110	131	P874220	113	P875360	120
P855100	146	P857380	147	P873120	131	P874230	113	P875370	126
P855210	146	P857410	147	P873130	132	P874250	113	P875390	121
P855220	146	P857420	147	P873140	132	P874260	113	P875400	121
P855230	146	P857450	148	P873150	132	P874270	113	P875410	121
P855240	146	P857460	148	P873160	132	P874280	114	P875430	123
P855250	146	P857470	148	P873170	133	P874290	114	P875440	123
P855260	146	P857480	148	P873180	133	P874300	114	P875450	123
P855290	146	P857490	148	P873190	133	P874310	114	P875460	123
P855300	146	P857500	148	P873200	141	P874320	114	P875470	123
P855310	146	P857570	148	P873210	141	P874330	114	P875480	123
P855320	146	P857580	148	P873220	147	P874340	114	P875490	123
P855350	146	P857610	148	P873230	123	P874350	111	P875500	123
P855360	146	P857620	148	P873240	123	P874360	114	P875510	124
P855370	146	P857650	148	P873250	124	P874370	114	P875520	124
P855380	146	P857660	148	P873260	147	P874390	114	P875530	124
P855490	146	P857690	148	P873270	147	P874400	114	P875540	124
P855500	146	P857700	148	P873280	147	P874410	124	P875550	124
P855510	146	P857730	148	P873290	121	P874420	114	P875560	124
P855520	146	P857740	148	P873300	121	P874430	114	P875570	124
P855530	146	P858010	8 . 75	P873310	122	P874440	114	P875580	125
P855540	146	P858050	8 . 75	P873330	121	P874450	114	P875600	125
P855550	146	P858090	8 . 74	P873340	122	P874460	117	P875610	125
P855560	146	P858130	8 . 73	P873350	122	P874470	114	P875620	125
P855590	146	P858170	8 . 73	P873360	124	P874480	115	P875630	125
P855600	146	P858450	8	P873370	117	P874500	115	P875640	125
P855610	146	P858490	8	P873380	117	P874510	115	P875650	125
P855620	146	P858500	8	P873390	117	P874520	115	P875660	125
P855630	146	P858530	8 . 74	P873400	118	P874530	115	P875670	125
P855640	146	P858570	8 . 74	P873410	126	P874540	115	P875680	125
P855650	146	P872000	147	P873420	131	P874550	115	P875690	125
P855660	146	P872010	147	P873430	131	P874560	116	P875700	126
P855690	146	P872040	122	P873440	112	P874570	116	P875710	126
P855700	146	P872070	144	P873450	112	P874580	116	P875720	126
P855710	146	P872080	144	P873460	112	P874590	116	P875730	126
P855720	146	P872100	145	P873470	112	P874600	116	P875740	126
P855730	146	P872110	145	P873480	112	P874610	116	P875750	126
P855740	146	P872120	144	P873490	114	P874620	116	P875770	126
P855800	146	P872150	144	P873500	114	P874630	116	P875780	126
P855810	146	P872170	145	P873510	117	P874640	116	P875790	126
P855860	146	P872180	145	P873520	117	P874650	116	P875800	126
P855870	146	P872210	146	P873530	117	P874660	116	P875810	126
P855940	146	P872230	115	P873540	118	P874670	116	P875820	126
P855950	146	P872240	115	P873550	118	P874680	116	P875830	126
P855980	128	P872270	122	P873570	126	P874690	116	P875850	130
P855990	128	P872280	122	P873580	126	P874700	116	P875840	126
P856020	146	P872300	122	P873590	129	P874710	116	P875850	126
P856030	146	P872370	124	P873600	129	P874720	116	P875860	126

Reference Number Index

Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page	Ref N°	Page
P875870	126	P876900	131	P877880	141	P878890	145	P879870	124
P875880	126	P876910	131	P877890	141	P878900	145	P879900	125
P875890	126	P876920	131	P877910	141	P878910	145	P879910	125
P875900	126	P876930	131	P877920	141	P878920	145	P879920	125
P875910	126	P876940	131	P877930	141	P878930	145	P879930	125
P875920	126	P876950	131	P877940	141	P878940	146	P879940	125
P875970	126	P876960	131	P877950	141	P878950	146	P879950	125
P875980	126	P876970	131	P877960	141	P878960	146	P879960	125
P875990	126	P876990	131	P877970	141	P878970	146	P879970	125
P876000	126	P877000	131	P877990	141	P878990	147	P879980	125
P876010	127	P877010	132	P878000	141	P879010	147	P879990	126
P876020	127	P877020	126	P878010	141	P879020	147	P880000	126
P876030	127	P877030	132	P878030	141	P879030	147	P880010	126
P876040	127	P877040	132	P878070	141	P879040	147	P880020	126
P876050	127	P877050	132	P878080	141	P879050	147	P880030	127
P876060	127	P877060	132	P878090	141	P879060	147	P880040	127
P876070	127	P877070	123	P878100	141	P879080	148	P880050	127
P876080	127	P877080	132	P878110	142	P879090	148	P880060	128
P876090	127	P877090	132	P878120	142	P879100	148	P880070	129
P876100	127	P877100	115	P878130	142	P879110	148	P880080	129
P876110	127	P877110	120	P878140	142	P879120	130	P880090	129
P876120	127	P877120	120	P878150	142	P879130	131	P880100	129
P876130	127	P877130	125	P878160	142	P879160	111	P880110	130
P876140	127	P877140	125	P878170	142	P879170	112	P880120	130
P876150	127	P877150	125	P878180	142	P879180	112	P880130	130
P876160	127	P877170	125	P878190	142	P879190	112	P880140	130
P876170	127	P877190	129	P878200	142	P879200	113	P880150	130
P876180	127	P877200	129	P878210	142	P879210	113	P880160	130
P876190	127	P877210	129	P878220	124	P879220	113	P880170	130
P876200	127	P877220	129	P878230	142	P879230	113	P880180	130
P876210	127	P877240	133	P878240	142	P879240	113	P880190	130
P876220	127	P877250	133	P878250	142	P879250	113	P880200	130
P876230	127	P877260	134	P878260	142	P879260	113	P880210	131
P876240	127	P877270	134	P878270	142	P879270	114	P880220	131
P876250	127	P877280	133	P878280	145	P879280	114	P880230	133
P876260	127	P877290	140	P878290	145	P879290	114	P880240	133
P876270	127	P877300	142	P878300	145	P879300	114	P880250	133
P876280	127	P877310	142	P878310	142	P879310	114	P880260	133
P876290	127	P877320	145	P878320	142	P879320	115	P880270	133
P876300	127	P877330	133	P878330	142	P879330	115	P880280	133
P876310	127	P877340	133	P878340	142	P879340	115	P880290	133
P876320	127	P877350	133	P878350	142	P879350	115	P880300	133
P876330	127	P877360	133	P878360	142	P879360	116	P880340	140
P876340	128	P877370	134	P878390	142	P879370	116	P880350	140
P876350	128	P877380	134	P878400	142	P879380	117	P880360	141
P876360	128	P877390	134	P878410	142	P879410	117	P880370	141
P876370	128	P877400	134	P878420	143	P879420	116	P880380	141
P876380	128	P877410	134	P878430	143	P879430	116	P880390	143
P876390	128	P877420	126	P878440	143	P879440	116	P880400	143
P876400	128	P877430	134	P878450	143	P879450	116	P880410	143
P876410	128	P877440	134	P878460	143	P879460	116	P880420	143
P876420	128	P877450	145	P878470	143	P879470	116	P880430	143
P876430	128	P877460	145	P878480	143	P879480	116	P880440	143
P876440	128	P877470	145	P878490	115	P879490	116	P880450	143
P876450	128	P877480	145	P878490	143	P879500	116	P880460	145
P876460	128	P877490	145	P878500	143	P879510	116	P880470	145
P876470	128	P877500	134	P878510	143	P879520	116	P880480	145
P876480	128	P877510	116	P878520	143	P879530	116	P880490	145
P876490	128	P877520	116	P878530	143	P879540	116	P880500	146
P876500	129	P877530	116	P878540	143	P879550	116	P880510	146
P876510	129	P877540	116	P878550	143	P879560	118	P880520	147
P876520	129	P877550	116	P878560	143	P879570	118	P880530	147
P876530	129	P877560	116	P878570	143	P879580	118	P880540	148
P876540	129	P877570	134	P878580	143	P879590	118	P880550	148
P876550	129	P877580	134	P878590	143	P879600	118	P880560	148
P876560	129	P877590	134	P878600	143	P879610	119	P880570	148
P876570	129	P877600	134	P878610	144	P879620	119	P880580	148
P876580	129	P877610	134	P878620	144	P879630	119	P880590	148
P876590	129	P877620	134	P878630	144	P879640	119	P880600	148
P876600	129	P877630	134	P878640	144	P879650	119	P880610	148
P876610	129	P877640	134	P878650	144	P879660	119	P880620	148
P876620	129	P877650	134	P878660	144	P879670	119	P880630	148
P876630	129	P877660	134	P878670	144	P879680	120	P880640	148
P876640	129	P877670	134	P878690	144	P879690	120	P880650	148
P876650	130	P877690	140	P878700	144	P879700	120	P880660	148
P876660	130	P877700	140	P878720	115	P879710	120	P880670	148
P876690	130	P877710	140	P878730	145	P879720	146	P880680	148
P876700	130	P877720	140	P878740	145	P879730	120	P880690	148
P876710	130	P877730	140	P878750	145	P879740	120	P880700	148
P876720	130	P877740	140	P878760	145	P879750	121	P880710	148
P876730	130	P877750	140	P878770	145	P879760	121	P880720	148
P876740	130	P877760	140	P878780	145	P879770	121	P880730	148
P876750	130	P877770	140	P878790	145	P879780	121	P880740	148
P876760	130	P877780	140	P878800	145	P879790	122	P880750	148
P876770	130	P877790	126	P878810	145	P879800	122	P880760	148
P876780	130	P877800	141	P878820	145	P879810	123	P880770	148
P876790	130	P877810	141	P878830	145	P879820	123	P880780	148
P876800	130	P877820	141	P878840	145	P879830	123	P880790	148
P876810	130	P877830	141	P878850	145	P879840	123	P880800	148
P876820	130	P877840	141	P878860	145	P879850	123	P880810	148
P876830	130	P877850	141	P878870	145	P879860	124	P880820	148
P876840	130	P877860	141	P878880	145			P880830	148
P876850	130	P877870	141					P880840	148
P876860	130							P880850	148
P876890	131							P880860	148



**The Experts in Custom-made Standards
Organic & Inorganic**

www.cpachem.com

